

# **uInProve**™ PLSV series GC injection valve

## THE MOST RELIABLE, DURABLE AND COMPACT GC VALVE TECHNOLOGY



### UNBEATABLE PERFORMANCE

- 1 million+ cycle in high purity trace gas application and still working without noticeable effect on chromatography
- Result never achieved before with a rotary valve
- Exceed diaphragm valve

### INNOVATIVE INSERT SLIDING VALVE TECHNOLOGY\*

- Reduced surface sealing area that improves leak integrity and reduces the required actuation torque
- Constant pressure to flow drop. Resolve diaphragm valve pressure/flow variation
- Slidable insert replaces the traditional rotor
- Valve inserts easily replaceable
- Inserts are offered in different materials
- High resistance to particulates

### MINIATURE DESIGN

- Smaller size than any other GC valve available on the market

### ACTUATION CONFIGURATION

- Compact electrical actuator
- Compact pneumatic actuator

## **uInProve**™ PLSV series GC valve

The uInProve PLSV series GC valve is a disruptive valve technology achieving both the lifetime of a diaphragm valve and the constant pressure drop and the simplicity of a rotary valve. This new valve technology is based on a reduced sealing surface area offered by the valve's inserts that replaces the traditional rotor. This revolutionary valve has been designed to respect our most elevated standards that we demand for.

### UNIQUE CHARACTERISTICS

- Available in 4, 6, 10, 12 ports (14 ports optional)
- Sample stream selection version
- Internal Sampling loop version
- Three temperature ranges version available
- Two pressure ranges version available
- Insert available in different materials





\*patent pending




## PURGED LIP SEALING VALVE (PLSV) TECHNOLOGY


### INNOVATIVE PURGE CONCEPT



### FLOW




### REDUCED SURFACE SEALING AREA




## ACTUATION PRINCIPLE

### TWO POSITIONS CONFIGURATION




**Position 1**




**Position 2**

\*patent pending

## INSERT CONFIGURATION



Purging groove to prevent cross port leaks




▼ This insert design, which is the standard configuration in this valve series, has a sealing surface which is 7 times less than a traditional rotary valve. This reduced surface, in combination with our innovative purge concept, result in an increase lifetime and better leak integrity. This unique sliding insert design is what provides a valve with unbeatable performance exceeding rotary and diaphragm valves.

▼ This insert is a variant of our standard insert design. It further reduces the sealing surface hence less friction, lower actuation torque and extended lifetime. It is used in the internal sampling loop and sample stream selection model.

## TYPICAL CHROMATOGRAPHIC CONFIGURATION

Here is a typical chromatographic method based on the ulnProve valve. Here, the valves are set in an heartcut method with a small restrictor for the purge. It is also to be noted that in this configuration, the sample loop is operated above ambient pressure. Our experience has demonstrated that this configuration, regardless of the valve technology, offers the best performance. It is however not a requirement, but a good design practice.





## ACTUATOR OPTION

### COMPACT PNEUMATIC ACTUATOR

For gas chromatographs that require a pneumatic actuation, this compact design is the perfect choice. Its compact size allows you to fit multiple valves in a single chassis for the most demanding applications.




- ▶ Allow easy replacement for diaphragm valve
- ▶ Compact dimension
- ▶ Operating temperature up to 180°C
- ▶ Stand-off extender available

### COMPACT ELECTRICAL ACTUATOR

This actuator has been designed for portable gas chromatographs or where no actuation air is available.

This is the most compact electrical actuator on the market, thanks to our valve design which requires much less torque compared to traditional rotary valves. Stand-off extenders are also offered for applications where the valve needs to be heated.



- ▶ The smallest electrical actuator on the market
- ▶ Perfect for portable gas chromatograph
- ▶ Used in ASDevices GCsSense
- ▶ Stand-off extender available




## **µInProve™** WITH INTERNAL SAMPLING LOOP

The same innovative sliding insert design can be used to offer valves with an internal sampling loop. In fact, one of the insert is used as the sampling loop. The insert internal volume is precisely machined to a specific volume. As the insert is made of rigid material, the volume remains constant regardless of the operating conditions, resolving the known issue of diaphragm valves with internal sampling loop. It is well known in the field that the internal sampling loop volume of diaphragm valves varies due to the elasticity of the flexible diaphragm membrane.

**Standard Sliding Insert**

**Sample Volume Insert**



0.1 µl

2.0 µl


Micro-Machined to a specific volume

Constant volume regardless of operating conditions


Available in different materials

Available between 0.1 to 2.0 µl  
Custom size available on request

**Sampling Position**







**Injection Position**



## PERFORMANCE TESTING – NEVER ACHIEVED BEFORE

The ulnProve chromatographic valve has been designed to offer a superior leak integrity and life time. In order to demonstrate its performance, the ulnProve valve has been subject to many tests and the ultimate one is a real trace measurement chromatographic application. Using the ultra sensitive ASDevices SePdd and iMov gas chromatographic platform, a GC method to measure trace ppb level N2 in H2 was setup and ran for over 1,000,000 actuations. This application, based on a heartcut method with two valves, was chosen as it is the best one to demonstrate leak integrity, life time, repeatability and constant flow characteristic over time.

With more than 1,000,000 actuations, the valves still offer the same outstanding performance.




Note: In order to achieve 1,000,000 cycles, the instrument was run continuously for one day on normal analysis mode with a cycle time of 5 minutes followed by an accelerated life testing mode where the valves were actuated every 5 seconds. This process, which was generating 104,256 actuations per week, was repeated 10 times for a total of 10 weeks.



## QUALITY TESTING

The ASDevices team has quite a long track record in developing high quality and high performance analytical and industrial valves. Consequently, we know what is required in term of design and process control to deliver the best quality valve.

Before entering the manufacturing process, every single components is thoroughly inspected for any defects.



### Cleanliness



During manufacturing, great care is taken to ensure that all parts are kept clean. Not only the wetted parts, but all parts. After being cleaned in an ultrasonic bath with water soluble detergents, all parts are cleaned with high pressure hot UHP water. Following this step, all wetted parts are subject to a proprietary inerting process. From that point, parts are continuously stored in clean packaging and handled with special gloves.

### Leak testing

Our unique patent pending purge system means our valves can't leak when operated within specification. However, it is still necessary during manufacturing to do a leak testing to make sure all parts are properly assembled and within specification.

Some manufacturers only do batch testing. We don't. We test every single valve with our ultra sensitive Enhanced Plasma Discharge technology\* and unique innovative leak testing method\*.




## HOW WE CATEGORISE VALVE

Chromatographic valves have always been categorised by their operating temperature with a fixed maximum operating pressure. The specified maximum operating pressure is also most of the time way too high for most chromatographic applications. Most chromatographic methods only involve pressures below 100 PSIG. Most capillary or packed column operate in the range of 10 to 50 PSIG. Pressure in the area of 300 PSIG are generally required for liquid sample only. Tuning a valve for high pressure operation means more stress on components hence an impact on its lifetime. Also, a valve tuned at 300 PSIG does not necessarily work better at 100 PSIG. In fact, its performance will be worst over time. This is why we offer not only different temperature ranges but also different pressure ranges. With this approach, we make sure that you have the best valve for your application.

\*patent pending






See AN-05 document for more configurations





## PHOTO GALLERY










# Specifications / Model


			PLSV-XX-XX-XX-LT-XX	PLSV-XX-XX-XX-MT-XX	PLSV-XX-XX-XX-HT-XX	
Standard maximum working pressure (PSI/kPa)						
SP: Standard pressure			1034/150	1034/150	1034/150	
HP: High pressure			2068/300	2068/300	2068/300	
Optional			Up to 700 Psig for liquid sample			
Maximum working temperature (°C/°F)			100/212	180/356	350/662	
Cross ports (Atm-cc/sec He)			Due to our patented design, it is impossible for our valve to develop a cross port leak when operated within specifications.			
In / Outboard (Atm-cc/sec He)						
Flow path to purge (Atm-cc/sec He)			3.0x10-12 6.0x10-13 (UHP version)	3.0x10-12	3.0x10-12	
Estimated actuations			1,000,000	750,000	500,000	
Valve cap material			300 Series stainless steel			
Cylinder body material			Anodized aluminium		Stainless steel	
Insert material			PEEK, Teflon, Vespel, other available			
Number of ports available			4, 6, 8, 10, 12, 14			
Sample stream configuration			Up to 10 ports			
Internal Sampling Loop Available			0.1 to 2 µL			
Port size (mm/in)			1/0.040			
Actuator type			Electrical / Pneumatic / Manual			
Actuation torque (lbs.in)			1 to 6 (# of Ports dependant)			
Pneumatic actuation pressure ( PSIG/kPa)			65 to 145 PSIG			
Surface cleaning procedure			Ultrasonic washing, O2 compatible			
Purged flow requirement			1 to 3 ml/min or no flow with miniature vaccum pump			
	4 Ports	6 Ports	8 Ports	10 Ports	12 Ports	14 Ports
Weight	205 g (0,45 lbs)	205 g (0,45 lbs)	296 g (0,65 lbs)	296 g (0,65 lbs)	296 g (0,65 lbs)	296 g (0,65 lbs)
Dimension (D x H)	39.2 x 66.7 mm (1,54 x 2,62 in)	39.2 x 66.7 mm (1,54 x 2,62 in)	48.2 x 66.7 mm (1,89 x 2,62 in)	48.2 x 66.7 mm (1,89 x 2,62 in)	48.2 x 66.7 mm (1,89 x 2,62 in)	48.2 x 66.7 mm (1,89 x 2,62 in)

## PLSV VALVE CONFIGURATION



### LEGAL DISCLAIMER

THE INFORMATION GIVEN IN THIS DOCUMENT (INCLUDING BUT NOT LIMITED TO CONTENTS OF REFERENCED WEBSITES) IS GIVEN AS A HINT FOR THE IMPLEMENTATION OUR PRODUCTS AND TECHNOLOGIES COMPONENT ONLY AND SHALL NOT BE REGARDED AS ANY DESCRIPTION OR WARRANTY OF A CERTAIN FUNCTIONALITY, CONDITION OR QUALITY OF THE ASDEVICES TECHNOLOGIES COMPONENT. THE RECIPIENT OF THIS DOCUMENT MUST VERIFY ANY FUNCTION DESCRIBED HEREIN IN THE REAL APPLICATION. ASDEVICES HEREBY DISCLAIMS ANY AND ALL WARRANTIES AND LIABILITIES OF ANY KIND (INCLUDING WITHOUT LIMITATION WARRANTIES OF NON-INFRINGEMENT OF INTELLECTUAL PROPERTY RIGHTS OF ANY THIRD PARTY) WITH RESPECT TO ANY AND ALL INFORMATION GIVEN IN THIS APPLICATION NOTE. THIS MATERIAL IS PROPRIETARY TO ANALYTICAL SENSING DEVICES AND CANNOT BE COPIED, REPRODUCED OR DISSEMINATED IN ANY WAY WITHOUT ITS PRIOR WRITTEN APPROVAL. THE TECHNOLOGY DESCRIBED HEREIN MAY BE SUBJECT TO PATENT PROTECTION OR OTHER FORMS OF INTELLECTUAL PROPERTY RIGHTS.



**Analytical Sensing Devices**

233, Jalbert St. W., Thetford Mines, QC, Canada G6G 7W1

Phone: 418.338.0299 Fax: 418.338.3485

Email: [info@asdevices.com](mailto:info@asdevices.com)

[www.asdevices.com](http://www.asdevices.com)