

Analysis of Permanent Gases:

More Challenging Than You Might Think

Mark Sinnott October 9, 2018


Analysis of Permanent Gases

- Sample Introduction Syringes and valves
- General discussion on PLOT columns
 - > The molesieve column is at the heart of permanent gas separations
- > Techniques when $CO_2 + C_2$'s, C_3 's, etc. is also needed
 - Column Isolation
 - Parallel columns
 - Cryogenic separations
 - Unique selectivity packed columns
- Fechniques for low level detection of hydrogen



- Aailent

Gas Sampling Valve



Gas sampling valve



Syringes for Gas injection

		Application	http://pconlab.net/Hamilton-Syringes.html
2	10–12° sharp, beveled, curved non-coring	Gas chromatography,	, septum piercing
3	Blunt, electro-polished	<u> </u>	uid chromatography (HPLC) injection, aphy (TLC), general liquid handling, ctions
3T	Blunt, electro-polished, coated with PTFE 19 mm from the tip	Thin layer chromatog	raphy (TLC) applications
4	Sharp 10-12° beveled needle	Life science/animal in	jections
5	Conical with side port for penetration without coring	Headspace, application minimal septum dama	ons prone to needle clogging, causes age
AS	Conical, non-coring designed to withstand multiple injections	Autosampler injection	n, pre-pierced septa



Gas Tight Syringes









PLOT columns

Column Type	Phase Type	Chromatographic Process	Stationary Phases
WCOT	Liquid	Gas - Liquid partitioning	Polysiloxanes PEG
PLOT	Solid	Gas - Solid adsorbtion	Porous Polymers Al_2O_3 , Zeolites, etc

WCOT = \underline{W} all \underline{C} oated \underline{O} pen \underline{T} ubular PLOT = \underline{P} orous \underline{L} ayer \underline{O} pen \underline{T} ubular

Why use a PLOT column?

Highly retentive


> (can resolve gases at non-cryo temperatures)

Unique selectivity

Permits higher initial oven temperatures than WCOT




PLOT Columns






WCOT Ethylene Analysis





PLOT Ethylene Analysis





JW Column PLOT column Portfolio- DB, HP, CP

HP-PLOT Al₂O₃ **GS-Alumina GS-Alumina KCI** HP-PLOT Al₂O₃ S HP-PLOT Al₂O₃ M CP-Al₂O₃/KCL $CP-AI_2O_3/Na2SO4$ **CP-SilcaPLOT CP-CarboBOND CP-CarboPLOT P7 GS-CarbonPLOT CP-PoraPLOT Amines CP-PoraPLOT S**


11

CP-PoraBOND Q **CP-PoraPLOT Q CP-PoraPLOT Q-HT** HP-PLOT Q **CP-PoraBOND U** HP-PLOT U GS-Q **CP-PoraPLOT U HP-PLOT** Molesieve **CP-Molsieve 5A GS-GasPro** ShinCarbon ST Select Permanent Gas Column



Selectivity Interactions in PLOT Phases

Shape / Size		Surface		
Zeolites	Bonded Carbon / Molecular Sieves	Porous Polymers	Bonded Silica	Al ₂ O ₃
Molesieve	Molesieve CarbonPLOT	PLOT-Q/U	GasPro SilicaPLOT	Alumina



Size Selectivity

> Vapor pressure always plays a leading role in solute interactions

Graphic of type KFI molecular sieve obtained from the International Zeolite Association web site (www.izo-online.org).



Molesieve Column

Separation based mostly on molecular size and shape

- Excellent at what it does but very limited
 - ≻H₂, O₂, N₂, CH₄, CO
 - ➢Noble Gases
 - ≻NO₂
 - ≻SF₆

>Limitations...cannot use for...

≻CO₂

13


October 9, 2018

- ≻ Water
- >C₂ HC's and larger



Molesieve Separations

Conditions



GC:

Sample:

Carrier:

Column:


6890

250 ml, split (75:1)

HP-PLOT/MoleSieve, 30 m x 0.32 mm x 12 µm


Helium, 2 µl/min





7.5

Noble Gases



Pub No.: 5964-2129E

2.5

5.0


0.0



Tips for using a Molesieve column

- Molesive material is commonly used as a moisture trap
- Samples must be reasonably dry
- Separation will change as column absorbs water over time
 - Resolution loss between CH₄ and CO indicates when column should be reconditioned (300 ° c 8 hours +)

Permanent Gases






But I also need to analyze for CO2....its getting complicated...

- \succ Recall that CO₂ cannot be done on a molesieve column (absorbed)
- > Alternative options:
 - > 2 unique injections on 2 columns
 - Molesieve + PLOT-Q
 - > 1 injection on 1 column at Cryo temps (-80 °c!)
 - GasPro
 - 1 Injection on 1 column ShinCarbon ST (packed column)
 - > 1 injection with 2 columns + valve to "direct traffic"
 - Column Isolation Molesieve + PLOT-Q
 - > 1 injection with 2 columns in parallel
 - Select Permanent Gas Column Molesieve + PLOT-Q ...did I mention this would get complicated?? ⁽ⁱ⁾









GasPro - Cryo


Permanent Gases in Hydrocarbon Blend















1 Injection 2 Columns + ValveColumn Isolation




Column isolation: Setting up valve timing




🔆 Agilent

Column isolation






Column Isolation – Flexibility of Elution Order....






1 Injection 2 Columns in Parallel...Select Permanent Gas Column System Set up






Select Permanent Gas Column





Select Permanent Gas Column






So you think that was complicated... what if I need to analyze for hydrogen too...

The trouble with hydrogen
 Must use TCD
 Sensitivity is carrier gas dependent
 Creates unique issues related sensitivity
 Recall the negative hydrogen peak from earlier



TCD

- When using He carrier sensitivity for hydrogen is on the order of only ~10%
- No problem right? I'll use Argon or Nitrogen carrier and get down to ppm levels for hydrogen...
- This is true, however now sensitivity for all the other compounds is now very poor...a real "catch 22"...
- > Lets see why this is the case...



"C'mon, c'mon - it's either one or the other."



Thermal Conductivity Detector - TCD


Name	Thermal Conductivity @ 400 K (mW m ⁻¹ K ⁻¹)			
Argon	22.4			
Hydrogen	230.9			
Helium	189.6			
Nitrogen	32.8			
Carbon Monoxide	32.3			
Carbon Dioxide	25.2			
Acetylene	33.3			
Ethylene	34.7			
Ethane	36.0			
Propane	31.0			
Butane	28.3			
Pentane	24.9			
Hexane	23.4			

		Thermal conductivity in mW m ⁻¹ K ⁻¹						
		100 K	200 K	300 K	400 K	500 K	600 K	Ref.
	Air	9.5	18.5	26.4	33.5	39.9	46.0	1
Ar	Argon $(P = 0)$	6.3	12.4	17.7	22.4	26.5	30.3	2, 3°
BF ₃	Boron trifluoride			19.0	24.6			4
CIH	Hydrogen chloride		9.2	14.5	19.5	24.0	28.1	4
F ₆ S	Sulfur hexafluoride $(P = 0)$			13.0	20.6	27.5	33.8	5
H ₂	Normal hydrogen (P = 0)	68.2	132.8	186.6	230.9	270.9	309.1	6
H ₂ O	Water $(P = 0)$			18.6	26.1	35.6	46.2	7
D ₂ O	Deuterium oxide $(P = 0)$			18.2	26.6	36.3	47.6	8
H ₂ S	Hydrogen sulfide			14.6	20.5	26.4	32.4	4
H ₃ N	Ammonia			25.1	37.2	53.1	68.6	9
He	Helium (P = 0)	74.7	118.3	155.7	189.6	221.4	251.6	10
Kr	Krypton $(P = 0)$		6.5	9.5	12.3	14.8	17.1	11
NO	Nitric oxide		17.8	25.9	33.1	39.6	46.2	4
N ₂	Nitrogen	9.4	18.3	26.0	32.8	39.0	44.8	1
N ₂ O	Nitrous oxide		9.8	17.4	26.0	34.1	41.8	4
Ne	Neon $(P = 0)$	22.3	37.4	49.4	59.9	69.5	78.5	12
O2	Oxygen	9.1	18.2	26.5	34.0	41.0	47.7	1
O ₂ S	Sulfur dioxide			9.6	14.3	20.0	25.6	4
Xe	Xenon $(P = 0)$		3.7	5.5	7.2	8.8	10.3	3°, 11
CCl ₂ F ₂	Dichlorodifluoromethane			9.9	15.0	20.1	25.2	13
CF ₄	Tetrafluoromethane $(P = 0)$			16.0	24.1	32.2	39.9	5
CO	Carbon monoxide ($P = 0$)			25.0	32.3	39.2	45.7	14
CO	Carbon dioxide		9.6	16.8	25.2	33.5	41.6	15
CHCl ₃	Trichloromethane			7.5	11.1	15.1		4
CH_4	Methane $(P = 0)$	10.4	21.8	34.4	50.0	68.4	88.6	16
CH₄O	Methanol				26.2	38.6	53.0	4
C ₂ Cl ₂ F ₄	1,2-Dichloro-1,1,2,2-tetrafluoroethane			10.3	15.7	21.1		13
C2Cl3F3	1,1,2-Trichloro-1,2,2-trifluoroethane			9.0	13.6	18.3		13
C ₂ H ₂	Acetylene			21.4	33.3	45.4	56.8	4
C2H4	Ethylene		11.3	20.6	34.7	49.9	68.6	17
C ₂ H ₆	Ethane		10.7	21.2	36.0	53.8	73.3	18
C ₂ H ₆ O	Ethanol			14.4	25.8	38.4	53.2	4
C ³ H ⁶ O	Acetone			11.5	20.2	30.6	42.7	4
C ₃ H ₈	Propane			18.5	31.0	46.4	64.6	19
C4F	Perfluorocyclobutane			12.5	19.5			13
C4H10	Butane			16.7	28.3	43.0	60.9	20
C4H10	Isobutane			17.1	28.9	43.2	60.2	21
C4H10O	Diethyl ether			15.1	25.0	37.1		4
C _s H ₁₂	Pentane				24.9	37.8	52.7	4
C ₆ H ₁₄	Hexane				23.4	35.4	48.7	4

' More accurate data covering a restricted temperature range.

https://ws680.nist.gov/publication/get_pdf.cfm?pub_id=907540







Workarounds for low level hydrogen detection

➤ 2 unique injections:

- 1 for hydrogen using Argon or N2 carrier 2nd injection for all other gases using helium carrier
- > 1 injection using column isolation and argon carrier + methanizer
 - TCD coupled to FID
 - Methanizer placed after TCD but before FID to convert CO/CO₂ to methane for enhanced detection of these gases by FID

> Sabatier Reaction (1897):
$$CO_2 + 4H_2 \xrightarrow[Ni]{400 \circ c} CH_4 + 2H_2O$$

Transformer Oil Gas Analyzer (TOGA)




Conversion of CO/CO2 to methane

Discovered by French chemist Paul Sabatier in 1897 (Sabatier Reaction)


$$CO_{2} + 4 H_{2} \xrightarrow{400 \circ c} CH_{4} + 2H_{2}O$$
Ni
$$Methanizer (H2 + Ni) \\ 350 \circ c FID$$








Transformer Oil Gas Analyzer (TOGA) Chromatograms






Transformer Oil Gas Analyzer (TOGA) Chromatograms





Transformer Oil Gas Analyzer (TOGA) Chromatograms





Summary

Analytes	Column	Technique
H ₂ /O ₂ /N ₂ /CH ₄ /CO	Molesieve	 1 Injection 1 Column
$H_2/O_2/N_2/CH_4/CO + Ar$	Molesieve (thick film)	 1 Injection 1 Column
$H_2/O_2/N_2/CH_4/CO + CO_2 + C2's$ etc	Molesieve + PLOT-Q	 2 Injections 2 Columns 1 Injection + Valve (Column Isolation) 1 Injection Parallel Columns (Select Perm Gas Column)
$H_2/O_2/N_2/CH_4/CO + CO_2 + C2$'s etc	GasPro	• Cryo (-80 ° c)
$H_2/O_2/N_2/CH_4/CO + CO_2 + C2$'s etc	ShinCarbon	Packed only
H ₂ /O ₂ /N ₂ /CH ₄ /CO + CO ₂ + C2's etc + Low level Hydrogen	Molesieve + PLOT-Q	 Single injection onto molesieve for H2 detection only Argon carrier + Methanizer (TOGA)



Contact Agilent Chemistries and Supplies Technical Support

We are always here to help!



1-800-227-9770 Option 3, Option 3:
Option 1 for GC/GCMS Columns and Supplies
Option 2 for LC/LCMS Columns and Supplies
Option 3 for Sample Preparation, Filtration and QuEChERS
Option 4 for Spectroscopy Supplies
Available in the USA 8-5 all time zones



gc-column-support@Agilent.com lc-column-support@agilent.com spp-support@agilent.com spectro-supplies-support@agilent.com



Outline

- Plot columns General discussion PLOT vs WCOT
- Molesieve
 - Very narrow functionality
 - Great for resolving O2/N2 at non-cryo temps
 - Traps CO2 and water
- The Problem with CO2 and HC larger than C1
 - Option 1: Make 2 unique injections
 - Option 2: Cryo + GasPro
 - Option 3: Column isolation
 - Option 4: ShinCarbon
- Detection issues when there is a need for H2
 - Explanation of TCD (catch 22)
 - Work arounds:
 - 2 unique injections
 - Argon carrier + Methanizer (TCD + FID in series)



