



## SELECTION GUIDE

---

# Capillary GC Column Selection Guide

Restek offers fused silica or stainless steel capillary columns in more than 900 stock combinations of stationary phase, stationary phase film thickness, column ID, and column length. If you are not sure about how to choose the best combination for your application, column selection can be a frustrating, hit-or-miss decision. The information in this guide can help you choose the proper column for your particular need. Further, it can help you to determine whether a column you already are using is the best choice, or if you might improve resolution, speed of analysis, and/or analyte quantification by using a different column.

## Parameters

- **Tubing Material**
- **Stationary Phase**
- **Stationary Phase Film Thickness**
- **Column Internal Diameter (ID)**
- **Column Length**



As always, your satisfaction with Restek chromatography columns is guaranteed. Please contact our **Technical Service Group** (or call 800-356-1688 or 814-353-1300, ext. 4), or your **Restek**





## SELECTION GUIDE

---

# Restek Capillary GC Column Selection Guide

## Tubing Material

Siltek Tubing

In either fused silica or stainless steel format, Restek columns offer excellent inertness, consistent column-to-column performance and, when installed and operated according to recommendations, long column lifetimes.

Fused silica Rxi®, Rtx®, or Rt™ columns are your best choice for most applications. They offer the highest coating efficiencies, ensuring the best resolution of closely eluting

compounds. Also, many Restek fused silica columns can be ordered with an [Integra-Guard™ integral guard column](#).

Rxi® columns, in particular, provide unmatched performance and exceptionally reliable column-to-column consistency. In developing these columns, we focused on achieving unsurpassed inertness, and the lowest bleed and the highest reproducibility possible. To achieve these goals, we hired the world's best polymer chemists and built a new state-of-the-art research facility. We established rigorous controls on tubing dimensions and surface activity, and we treat this highly uniform tubing with a unique deactivation chemistry, producing a consistent, inert surface on which to apply a polymer. We reformulated our polymers, ensuring neutrality and fine tuning

selectivity for retention time locking.

The resulting exceptionally low-bleed columns are ideal for trace-level GC/MS analysis: with an Rxi® column, you can chromatograph sub-nanogram levels of active acidic or basic compounds on the same column — often under the same conditions.

To ensure reproducibility, we redeveloped our manufacturing process to maximize column-to-column consistency. Every Rxi® column is individually tested and proven to meet our stringent specifications for coating efficiency, selectivity, film thickness, inertness, and bleed. We guarantee every Rxi® column you receive will be the best column you have ever used.

**Under harsh operating conditions, choose stainless steel MXT® columns:**

- Rough handling (e.g., field instruments or process GC).
- Portable instruments / other small ovens requiring tightly coiled columns.
- High temperature chromatography.

When the potential for breakage is high, MXT® columns are your best choice — they present little risk of spontaneous breakage when used at high temperatures or when coiled into small diameters. While fused silica columns generally cannot be used above 360°C, because the polyimide outer coating becomes brittle over time at high temperatures, MXT® columns function well at temperatures exceeding 450°C. With an MXT® column, the only limitation to oven temperature is the operating limit of the stationary phase.

To prepare MXT® columns, we treat the internal surface of the stainless steel tubing with our exclusive Siltek® surface treatment, making the surface as inert as deactivated fused silica. The Siltek® layer permeates the surface, rather than simply coating it, making the layer

exceptionally flexible, so the tubing can be coiled to very small diameters. Coating efficiencies for MXT® columns are slightly lower than for fused silica columns, because the metal tubing has a larger surface area. We offer MXT® columns with a wide variety of stationary phases, in 0.18mm, 0.25mm, 0.32mm, and 0.53mm ID. The minimum coil diameter for 0.25mm ID or 0.53mm ID MXT® columns is 1.5 inches or 2.5 inches, respectively.

# Restek Capillary GC Column Selection Guide

## Stationary Phase

The stationary phase is the single most important consideration when you are choosing a column. The interactions between the analytes and the functional groups of the stationary phase contribute more to the overall results of the analysis than any other factor.

**Table 1** summarizes the characteristics, chemical structures and, broadly, suggested uses for Restek general-purpose stationary phases. Change selectivity by choosing a stationary phase with a different percentage of substitution of a particular functional group (e.g., by switching from a 5% diphenyl/95% dimethyl polysiloxane stationary phase to a 20% diphenyl/80% dimethyl polysiloxane stationary phase) or by choosing a stationary phase with different functional groups (e.g., by switching from a diphenyl/dimethyl polysiloxane stationary phase to a polyethylene glycol stationary phase). Note that a stationary phase's selectivity for sample components follows the general chemical principle of "like prefers like": a nonpolar stationary phase, such as the Rxi®-1 methyl polysiloxane stationary phase, will preferentially retain and separate nonpolar compounds, such as straight-chain hydrocarbons, relative to polar compounds, such as alcohols. As methyl groups are replaced with more polar functionalities, such as phenyl or cyanopropyl groups, selectivity shifts toward more polar compounds. Polar phases, such as polyethylene glycol phases (e.g., Stabilwax®), are highly selective toward alcohols or other polar compounds.

**Table 2** lists Restek special-purpose stationary phases. Columns with these phases are our first recommendations for the applications noted.


**Table 3** lists stationary phases we recommend for environmental analyses by US EPA methods.

**Table 4** lists retention indices for test compounds on the stationary phases characterized in Table 1. A retention index is a mathematical derivation indicating the elution position of a compound with respect to normal (straight chain) hydrocarbons. For example, a retention index of 650 for benzene on a particular stationary phase indicates benzene will elute mid-way between n-hexane (RI=600) and n-heptane (RI=700). The longer a particular compound is held by a stationary phase, the greater the retention index will be for that compound. Similarly, the greater the separation between two compounds, the greater the difference between their retention indices. To review retention indices for a wide variety of compounds on a range of Restek stationary phases, see the [retention index tables](#) in our on-line Expert Center.

Note that if you will be using a selective detector, background levels in your chromatograms can be high if the stationary phase contains elements the detector is designed to respond to. For example, avoid using a cyanopropyl (CN)-containing stationary phase with a nitrogen-phosphorus detector, or a fluorine-containing stationary phase with an electron capture detector.

**Table 1** Stationary Phase Structures and Properties.

Structure Diagram	Properties
-------------------	------------





**Rxi@-1ms/Rtx@-1MS**

100% dimethyl polysiloxane

Stable to 360°C

Polarity: non-polar

Uses: solvents, petroleum products, pharmaceutical samples, waxes




**Rxi@-5ms/Rtx@-5/Rtx@-5MS**

5% diphenyl - 95% dimethyl polysiloxane

Stable to 360°C

Polarity: slightly polar

Uses: flavors, environmental samples, aromatic hydrocarbons





**Rxi@-5Sil MS/Rtx@-5Sil MS**

proprietary

Stable to 360°C

Polarity: slightly polar

Uses: flavors, environmental samples, pesticides, PCBs, aromatic hydrocarbons





**Rtx@-20**

20% diphenyl - 80% dimethyl polysiloxane

Stable to 310°C

Polarity: slightly polar

Uses: volatile compounds, alcohols





**Rtx@-1301, Rtx@-624, Rtx@-G43**

6% cyanopropylphenyl - 94% dimethyl polysiloxane

Stable to 280°C

Polarity: slightly polar

Uses: volatile compounds, insecticides, residue solvents in pharmaceutical products



**Rtx@-35**


35% diphenyl - 65% dimethyl polysiloxane

Stable to 300°C

Polarity: intermediately polar

Uses: pesticides, Aroclor PCBs, amines, nitrogen-containing herbicides






**Rtx@-200**

trifluoropropylmethyl polysiloxane

Stable to 360°C

Polarity: selective for lone pair electrons

Uses: environmental samples, solvents, Freon® gases, drugs, ketones, alcohols





**Rtx@-50**

100% methylphenyl polysiloxane

Stable to 340°C

Polarity: intermediately polar

Uses: FAMES, carbohydrates





**Rxi@-17**

50% diphenyl - 50% dimethyl polysiloxane

Stable to 300°C

Polarity: intermediately polar

Uses: triglycerides, phthalate esters, steroids, phenols





**Rtx@-1701**

14% cyanopropylphenyl - 86% dimethyl polysiloxane

Stable to 280°C

Polarity: intermediately polar

Uses: pesticides, Aroclor PCBs, alcohols, oxygenates





**Rtx@-65TG**

65% diphenyl - 35% dimethyl polysiloxane

Stable to 370°C

Polarity: intermediately polar

Uses: triglycerides, rosin acids, free fatty acids





**Rtx®-225**

50% cyanopropylmethyl - 50% phenylmethyl polysiloxane

Stable to 260°C

Polarity: polar

Uses: FAMES, carbohydrates




**Rf™-2330**

90% biscyanopropyl - 10% cyanopropylphenyl polysiloxane

Stable to 275°C

Polarity: very polar

Uses: *cis/trans* FAMES, dioxin isomers, rosin acids




**Stabilwax®**

Carbowax® PEG

Stable to 250°C

Polarity: polar

Uses: FAMES, flavors, acids, amines, solvents, xylene isomers



**Rtx®-Wax**

Carbowax® PEG

Stable to 250°C

Polarity: polar

Uses: FAMES, solvents, BTEX aromatics, flavors

[^ Back](#)

**Table 2** Restek special-purpose stationary phases.

Application	Stationary Phase
acids (underivatized)	Stabilwax®-DA
amines (underivatized); other basic compounds	Stabilwax®-DB
amines; other basic compounds	Rtx®-5 Amine Rtx®-35 Amine

blood alcohol; glycols	Rtx®-BAC1 Rtx®-BAC2
chiral compounds	Rt™-βDEX Rt™-γDEX
detailed hydrocarbon analysis (ASTM/CGSB)	Rtx®-IPONA
dioxin and furan congeners	Rtx®-Dioxin Rtx®-Dioxin2
fatty acid methyl esters (FAMEs)	FAMEWAX
flavor and fragrance components	Rtx®-1 F&F Rt™-CW20M F&F
nitroaromatic explosives e.g.: US EPA method 8095	Rtx®-TNT Rtx®-TNT2
organochlorine pesticides e.g.: US EPA methods 8081, 608, CLP Pesticides	Rtx®-CLPesticides Rtx®-CLPesticides2 Stx™-CLPesticides Stx™-CLPesticides2
organophosphorus pesticides e.g.: US EPA method 8141A	Rtx®-OPPesticides Rtx®-OPPesticides2
PCB congeners	Rtx®-PCB Stx™-500
residual solvents in pharmaceuticals	Rtx®-G27 Rtx®-G43
simulated distillation	MXT®-500 Sim Dist
simulated distillation: ASTM test method D2887	Rtx®-2887
volatile organic compounds e.g.: US EPA methods 502.2, 524.2, 601, 602, 624, 8010, 8020, 8260	Rtx®-VGC Rtx®-VMS Rtx®-VRX Rtx®-Volatiles Rtx®-502.2 Rtx®-624

[^ Back](#)

**Table 3** Stationary phases we recommend for environmental analyses.

For example chromatograms, click on the name of the stationary phase.

US EPA Method	Stationary Phase
502.2 (volatile organics)	<a href="#">Rtx®-502.2</a>
504.1 (dibromoethane/dibromochloropropane)	<a href="#">Rtx®-CLPesticides</a> <a href="#">Rtx®-CLPesticides2</a>
506 (phthalate & adipate esters)	<a href="#">Rxi®-5Sil MS/Rtx®-5Sil MS</a> <a href="#">Rxi®-1ms</a>
515/515.1 (chlorophenoxyacid herbicides)	<a href="#">Rtx®-440</a> <a href="#">Rtx®-50</a> <a href="#">Rtx®-CLPesticides</a> <a href="#">Rtx®-CLPesticides2</a>
524.2, Rev. IV (volatile organics)	<a href="#">Rtx®-VMS</a>
525.2 (semivolatile organics)	<a href="#">Rxi®-5ms</a>
526 (semivolatile organics)	<a href="#">Rxi®-5Sil MS/Rtx®-5Sil MS</a>
528 (phenols)	<a href="#">Rxi®-5Sil MS/Rtx®-5Sil MS</a>
551.1 (chlorinated disinfection byproducts)	<a href="#">Rtx®-5</a> <a href="#">Rtx®-200</a> <a href="#">Rtx®-1301</a>

552.2 (haloacetic acids)	<a href="#">Rtx@-5</a> <a href="#">Rtx@-200</a> <a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a>
601 (volatile organics)	<a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a> <a href="#">Stx™-CLPesticides</a> <a href="#">Stx™-CLPesticides2</a>
602 (volatile organics)	<a href="#">Rtx@-VMS</a> <a href="#">Rtx@-VGC</a>
604 (phenols)	<a href="#">Rtx@-5</a> <a href="#">Rtx@-50</a> <a href="#">Rtx@-200</a>
608 (organochlorine pesticides & PCBs)	<a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a> <a href="#">Rtx@-PCB</a> <a href="#">Stx™-500 (PCBs)</a>
610 (polynuclear aromatic hydrocarbons)	<a href="#">Rxi@-5ms</a> <a href="#">Rxi@-5Sil MS/Rtx@-5Sil MS</a> <a href="#">Rtx@-440</a>
615 (chlorophenoxyacid herbicides)	<a href="#">Rtx@-35</a> <a href="#">Rtx@-5</a>
619 (organonitrogen/organophosphorus pesticides)	<a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a> <a href="#">Rtx@-50</a> <a href="#">Rtx@-200</a>
624 (volatile organics)	<a href="#">Rtx@-VMS</a> , <a href="#">Rtx@-VGC</a>
1671 (volatile organics)	<a href="#">Stabilwax®</a> <a href="#">Stabilwax®-DB</a>
8010 (volatile organics)	<a href="#">Rtx@-VMS</a> <a href="#">Rtx@-VGC</a>
8020 (volatile organics)	<a href="#">Rtx@-VMS</a> <a href="#">Rtx@-VGC</a>
8021 (volatile organics)	<a href="#">Rtx@-VGC</a> <a href="#">Rtx@-VRX</a>
8081/8081A (organochlorine pesticides)	<a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a> <a href="#">Stx™-CLPesticides</a> <a href="#">Stx™-CLPesticides2</a> <a href="#">Rtx@-440</a> <a href="#">Rtx@-XLB</a>
8095 (nitroaromatic explosives)	<a href="#">Rtx@-TNT</a> <a href="#">Rtx@-TNT2</a>
8100 (polynuclear aromatic hydrocarbons)	<a href="#">Rtx@-5</a>
8140/8141/8141A (organophosphorus pesticides)	<a href="#">Rtx@-OPPesticides</a> <a href="#">Rtx@-OPPesticides2</a> <a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a>
8151/8151A (organochlorine pesticides)	<a href="#">Rtx@-CLPesticides</a> <a href="#">Rtx@-CLPesticides2</a> <a href="#">Rtx@-440</a>
8240 (volatile organics)	<a href="#">Rtx@-VMS</a> <a href="#">Rtx@-VGC</a>
8260/8260B (volatile organics)	<a href="#">Rtx@-VMS</a>
8270/8270D (semivolatile organics)	<a href="#">Rxi@-5Sil MS/Rtx@-5Sil MS</a> <a href="#">Rxi@-5ms</a> <a href="#">Rtx@-XLB</a>
Appendix IX (semivolatile organics)	<a href="#">Rxi@-5Sil MS/Rtx@-5Sil MS</a> <a href="#">MS</a>
CLP (semivolatile organics)	<a href="#">Rtx@-5ms</a>

TO-14/TO-15 (air toxins)	<a href="#">Rtx®-1</a> <a href="#">Rtx®-502.2</a>
--------------------------	--

[^ Back](#)

**Table 4** Retention indices for Restek general purpose stationary phases.

Phase	Benzene	Butanol	Pentanone	Nitropropane
Rtx®-1	651	651	667	705
Rtx®-5/Rtx®-5MS	667	667	689	743
Rtx®-20	711	704	740	820
Rtx®-1301/Rtx®-624	689	729	739	816
Rtx®-35	746	733	773	867
Rtx®-200	738	758	884	980
Rtx®-50	778	769	813	921
Rtx®-1701	721	778	784	881
Rtx®-65TG	794	779	825	938
Rtx®-225	847	937	958	958
Stabilwax®	963	1158	998	1230

[Back](#)

Restek Capillary GC Column Selection Guide

[< prev](#)[next >](#)

## CONTENTS

- Tubing Material
- Stationary Phase
- Stationary Phase Film Thickness
- Column Internal Diameter (ID)
- Column Length

## ALSO OF INTEREST

- GC Retention Time Indexes
- Capillary GC Column Installation Guide
- Capillary GC Column Preventative Maintenance Guide
- Capillary GC Column Troubleshooting Guide

# Restek Capillary GC Column Selection Guide

## Stationary Phase Film Thickness

Stationary phase film thickness affects the retention and elution temperature for each compound in the sample. A thicker film retains compounds longer, increasing the length of time each compound spends in the column ([Table 1](#)). A thinner film retains compounds less, reducing the length of time each compound spends in the column. Very volatile compounds should be analyzed on a thick film column, to increase the time they spend in the column and allow them to separate. High molecular weight compounds must be analyzed on a thin film column, to reduce analysis time to a practical interval, and help minimize bleed at the higher temperatures required to elute such compounds.


A comparative analysis of low boiling compounds on a 0.25 $\mu$ m, a 1.0 $\mu$ m, and a 5.0 $\mu$ m film of stationary phase, with all other variables held constant, shows the influence of film thickness ([Figure 1](#)). The 0.25 $\mu$ m film does not resolve butanol from benzene (peaks 1 & 2). The 1.0 $\mu$ m film provides about 80% resolution of this pair, but retention times are more than double those for the 0.25 $\mu$ m film. The 5.0 $\mu$ m film does not improve resolution between butanol and benzene, relative to the 1.0 $\mu$ m film, and retention times are increased by a factor of six relative to the 0.25 $\mu$ m film. So, for this particular sample, the 1.0 $\mu$ m phase film is best: resolution is suitable for quantifying the analytes, analysis time is acceptable, and a thicker film does not offer notable improvements. On the other hand, if we wanted to resolve very volatile C2 or C3

compounds, which would elute prior to peak 1, the 5.0 $\mu$ m film would be required.


Both sample capacity and bleed increase as stationary phase film thickness is increased.

Changes in the column ID/stationary phase film thickness ratio,  $\beta$ , are an important consideration when you want to make a change in column internal diameter or in stationary phase film thickness. Analyte retention increases as column internal diameter is decreased, and analyte retention decreases as stationary phase film thickness is decreased. When other column parameters and analysis conditions are held constant, a column with a smaller  $\beta$  value will be more retentive for a given analyte. To assure similar retention and resolution when you increase the column ID, you also must increase the stationary phase film thickness. Similarly, if you wish to reduce the column ID, but keep retention and resolution similar, you also must reduce the stationary phase film thickness. **Table 2** lists  $\beta$  values for common combinations of column ID and stationary phase film thickness; where  $\beta$  values for different combinations are similar, elution patterns will be similar. To find  $\beta$  values for other column ID/stationary phase film combinations, use the **Beta Calculator**.

**Figure 1** A sample containing low boiling components shows the differences in resolution among 0.25, 1.0, and 5.0 $\mu$ m columns. The 1.0 $\mu$ m column offers better resolution than the 0.25 $\mu$ m column, and the 5.0 $\mu$ m column does not offer any further improvements for compounds eluting after C6.




0.25µm, Rtx®-1



1.0µm, Rtx®-1

5.0µm, Rtx®-1





1. butanol
2. benzene
3. 2-pentanone
4. C7
5. 1-nitropropane
6. pyridine
7. C8
8. C9
9. C10

(Peak 9 elutes @ 117 minutes on the 5.0 $\mu$ m)

[^ Back](#)

**Table 1** Characteristics of thick film and thin film stationary phases.

Thick Films	Thin Films
higher sample capacity	lower sample capacity
less efficient	more efficient
longer retention times	shorter retention times
effectively retain lower molecular weight analytes	efficiently release higher molecular weight analytes
higher bleed	lower bleed

[^ Back](#)

**Table 2** Column ID: phase film ratio ( $\beta$ ) value calculation for film thickness vs. column IDs

phase film thickness (df) / $\beta$ value							
Column ID	0.10 $\mu$ m	0.25 $\mu$ m	0.50 $\mu$ m	1.00 $\mu$ m	1.50 $\mu$ m	3.00 $\mu$ m	5.00 $\mu$ m
0.18mm	450	180	90	45	30	15	9
0.25mm	625	250	125	63	42	21	13
0.32mm	800	320	160	80	53	27	16
0.53mm	1325	530	265	128	88	43	27

[^ Back](#)

Restek Capillary GC Column Selection Guide

< prevnext >

## CONTENTS

- Tubing Material
- Stationary Phase
- Stationary Phase Film Thickness
- Column Internal Diameter (ID)
- Column Length

## ALSO OF INTEREST

- GC Retention Time Indexes
- Capillary GC Column Installation Guide
- Capillary GC Column Preventative Maintenance Guide
- Capillary GC Column Troubleshooting Guide



## SELECTION GUIDE

---

# Restek Capillary GC Column Selection Guide

## Column Internal Diameter (ID)

Your need for **resolution**, the **concentrations of analytes** you anticipate in your samples, and the **sample introduction and analyte detection** instrumentation you are using will influence your choice of column internal diameter.

When all other column parameters and analysis conditions are held constant, analyte retention increases, and resolution improves, as column internal diameter is decreased. Sample capacity increases as internal diameter is increased. Bleed also increases as column diameter is increased and film thickness is held constant, because there is more stationary phase in the column, however for most applications the difference in bleed will be insignificant.

Columns with an ID of 0.25mm or less offer the highest column efficiencies, and therefore the greatest resolution, but have limited sample capacities, relative to wider ID columns. If concentrations of analytes in your samples exceed the capacity of the column, the analysis will be characterized by peak distortion, poor resolution, and poor reproducibility. For many applications a 0.32mm ID column offers the best balance of resolution and sample capacity.

0.53mm ID columns are best for high flow situations, such as accommodating incoming samples

from a purge and trap unit. Alternatively, columns with an ID of 0.25mm or less are the best choice for GC/MS applications — they provide optimal performance with small flows of carrier gas, and thus can be connected directly to the vacuum source of a mass spectrometer, eliminating the complications associated with a jet separator.

Use **Table 1** to compare typical column characteristics according to column ID, then select the ID that best matches your needs.

**Table 1** Characteristics of thick film and thin film stationary phases.

Select Column ID:

Column Parameters	Values
Typical Spitless Purge Time	
He split vent flow rate (50:1 split ratio)	
H2 split vent flow rate (50:1 split ratio)	
Column OD (Ferrule ID)	
Column Flow Rate (He @ 20cm/sec.)	
Column Flow Rate (H2 @ 40cm/sec.)	
Approximate Sample Capacity	
Typical Effective Plates/m (80% CE)	

Approximate column head pressure (He or H2 carrier gas)	
Column (m):	6101215202530406075105150
Head Pressure (psig):	

Always optimize the column flow rate by using linear velocity.

## Phase Ratio Calculations

Changes in the column ID/stationary phase film thickness ratio, beta, are an important consideration when you want to make a change in column internal diameter or in stationary phase film thickness. Analyte retention increases as column internal diameter is decreased, and analyte retention decreases as stationary phase film thickness is decreased. When other column parameters and analysis conditions are held constant, a column with a smaller beta value will be more retentive for a given analyte. To assure similar retention and resolution when you increase

the column ID, you also must increase the stationary phase film thickness. Similarly, if you wish to reduce the column ID, but keep retention and resolution similar, you also must reduce the stationary phase film thickness. **Table 2** lists beta values for common combinations of column ID and stationary phase film thickness; where beta values for different combinations are similar, elution patterns will be similar. To find beta values for other column ID/stationary phase film combinations, use the [Beta Calculator](#).

**Table 2** Column ID: phase film ratio (beta) values for commonly used column dimensions.

Column ID	phase film thickness (df) / beta value						
	0.10µm	0.25µm	0.50µm	1.00µm	1.50µm	3.00µm	5.00µm
0.18mm	450	180	90	45	30	15	9
0.25mm	625	250	125	63	42	21	13
0.32mm	800	320	160	80	53	27	16
0.53mm	1325	530	265	128	88	43	27

Restek Capillary GC Column Selection Guide

< prevnext >

## CONTENTS

- Tubing Material
- Stationary Phase
- Stationary Phase Film Thickness
- Column Internal Diameter (ID)
- Column Length

## ALSO OF INTEREST

- GC Retention Time Indexes
- Capillary GC Column Installation Guide
- Capillary GC Column Preventative Maintenance Guide

# Restek Capillary GC Column Selection Guide

## Column Length

Longer columns provide more resolving power, but increase analysis time and column purchase costs. As you consider whether the increase in resolution is worth the extra time and expense, remember that the benefit of using a longer column is much greater if you are performing a temperature programmed analysis than if you are performing an isothermal analysis. In an **isothermal analysis**, retention time varies directly with column length: if column length is doubled, analysis time will double as well. The increase in resolution will be only approximately 40%, however, because resolution is related to the square root of column length, as shown in the equation below. In a **temperature programmed analysis**, retention times are more dependent on temperature than on column length. As column length is increased, the increase in resolution will be the same as for an isothermal analysis, but there will be only a marginal increase in analysis time.

# Length Effects

Length affects resolution and speed of analysis.

$$\text{Resolution} = \frac{1}{4} \sqrt{\frac{L}{h}} \times \frac{k}{k+1} \times \frac{\alpha-1}{\alpha}$$

L = length    h = HETP    k = capacity factor     $\alpha$  = selectivity