





# VICI Jour

Products for Chromatography  
and Fluid/Gas Transfer

2005/06





## Tubing for Chromatography

PAGE 1–25




## High Pressure Polymeric & Stainless Steel Fittings Nanovolume™ Fittings, Tools

PAGE 26–47



## Low Pressure Polymeric Fittings

PAGE 48–63



## Filters & Mobile Phase Filters

PAGE 64–75



## Mobile Phase Handling Systems

PAGE 76–81



## Frits

PAGE 82–87



**PEEK Columns & Guard Columns**

PAGE 88–93

**Sample Loops**

PAGE 94–101

**Accessories/Miscellaneous**

PAGE 102–105

**Physical Data & Info**

PAGE 106–109

**Properties of Metals & Polymers**

PAGE 110–111

**Glossary**

PAGE 112–115

**Part Number Index**

PAGE 116–120

**Trademarks and Registered Trademarks**

PAGE 120

**International Distributors**

PAGE 122



C0392

## INDEX

PEEK Tubing	5
PEEKclad-FST™	8
Premium Grade Stainless Steel Tubing	10
Stainless Steel Tubing	11
Electroformed Nickel Tubing	12
Pre-Cut Premium Stainless Steel Tubing	13
PTFE & PTFE-Like Tubing	14
No-Ox Tubing	18
Thermo-Formed Polymeric Tubing	19
Knitted PTFE Reaction Coils	20
PEEK Tubing Elbows	21
Tubing Clip – The LC Tubing Organizer	21
Clean-Cut Polymer Tubing Cutter	22
Guillotine Polymer Tubing Cutter	22
Stainless Steel Tubing Cutter	23
Stainless Steel Tubing Pliers	23
Easy-Flange & Easy-Flange Combi-Kits	24

### Disclaimer:

The maximum holding pressure for any type of connection involving tubing and a ferrule, varies considerably with the tubing material, the ferrule material, the clearance between tubing OD and ferrule ID and the shape of the fitting detail.



## PEEK Tubing

- 1/32", 1/16", 1/8" and 1/4" OD Tubing Available
- Premium Grade Tubing with Tighter Tolerances
- Easy to Cut
- Biocompatible

PEEK tubing has the strength required to withstand continuous use at HPLC pressure without swelling or bursting. The dense polymer structure of VICI Jour PEEK tubing eliminates the permeability to organic solvents that causes other polymer tubing, such as ETFE to "sweat".

For easy identification of ID VICI Jour tubing is color coded. VICI recommends "Striped Color-Coded" tubing. It is manufactured of virgin, natural PEEK and therefore ensures maximum chemical resistance and biocompatibility. The color dye is applied only to the outer surface and thus not in contact with the fluid stream. "Solid Color-Coded" tubing however is made of pre-dyed raw material.

Remember PEEK tubing is not affected by halide salts, high strength buffers or other aggressive mobile phases that corrode stainless steel. The polymer

surface will not leach metal ions into the eluent or extract metal sensitive components from the sample. By far the most outstanding property is its excellent burst pressure.

PEEK is compatible with almost any of the solvents commonly used in HPLC. Dichloromethane, THF and DMSO may cause swelling in PEEK. Concentrated nitric and sulphuric acid will attack PEEK.

### Premium Grade PEEK Tubing – Dash-Striped Color-Coded

The new VICI Jour dash-striped premium grade PEEK tubing is the best choice for chromatography applications in which tube volume is of importance. The tighter tolerance of  $\pm 0.025$  mm ensures a more consistent internal volume of the transfer lines.



### Premium Grade PEEK Tubing – Dash-Striped Color-Coded

Part No.	OD	ID (mm)	Qty/pkg	Color	bar*	psi*
JR-TP-5998-M3	1/16"	0.064	3 m	Pink	435	6300
JR-TP-5998-M10	1/16"	0.064	10 m	Pink	435	6300
JR-TP-5999-M3	1/16"	0.13	3 m	Red	420	6100
JR-TP-5999-M10	1/16"	0.13	10 m	Red	420	6100
JR-TP-6000-M3	1/16"	0.18	3 m	Yellow	400	5800
JR-TP-6000-M10	1/16"	0.18	10 m	Yellow	400	5800
JR-TP-6001-M3	1/16"	0.25	3 m	Blue	386	5600
JR-TP-6001-M10	1/16"	0.25	10 m	Blue	386	5600

\* = Recommended maximum permanent working pressure (measured with ACN/Water 1:1 at room temperature)  
Other dimensions on request

## SPECS

**Material**  
PEEK

**Dimensions**  
OD and ID, see chart

**Tolerances**  
Premium Grade Tubing  
OD:  $\pm 0.025$  mm (.001")  
ID:  $\pm 0.025$  mm (.001")

**Pressure rating**  
See chart

**Chemical resistance**  
See chart on page 106

Max. recommended working temp.:  $< 100$  °C (continuous)  
for 1/16" OD tubing with ID up to 0.75 mm

**Spares & Tools**  
A clean burr-free perpendicular cut can be achieved with the VICI Jour Clean-Cut tubing cutter JR-797 see page 22

To bend PEEK tubing at the optimum radius, use our Tubing Elbows on page 21

We recommend PEEK finger-tight and flangeless fittings for your applications  
Fingertights see pages 27–30  
Flangeless see pages 49–53

VICI Jour Sample Loops for analytical applications are made of Premium Grade Tubing (see pages 96–100)

**Length Conversions**  
0.064 mm  $\approx$  .0025"  
0.13 mm  $\approx$  .005"  
0.18 mm  $\approx$  .007"  
0.25 mm  $\approx$  .010"

# Tubing for Chromatography

## PEEK Tubing

### SPECS

**Material**  
PEEK

**Dimensions**  
OD and ID, see chart

**Tolerances**  
For OD 1/16" PEEK tubing  
OD: +/- 0.05 mm (.002")  
ID: +/- 0.05 mm (.002")

For OD 1/8" PEEK tubing  
OD: +/- 0.10 mm (.004")  
ID: +/- 0.10 mm (.004")

**Pressure rating**  
See chart

**Chemical resistance**  
See chart on page 106

Max. recommended working  
temp.: < 100 °C (continuous)  
for 1/16" OD tubing with ID  
up to 0.75 mm

**Spares & Tools**  
A clean burr-free perpendicular  
cut can be achieved with  
the VICI Jour Clean-Cut tubing  
cutter JR 797 see page 22

To bend PEEK tubing at the  
optimum radius, use our  
Tubing Elbows on page 21

We recommend PEEK finger-  
tight and flangeless fittings  
for your applications  
Fingertights see pages 27-30  
Flangeless see pages 49-53

VICI Jour Sample Loops  
for analytical applications  
are made of Premium  
Grade Tubing  
(see pages 96-100)

### PEEK Tubing – Striped Color-Coded

Part No.	OD	ID (mm)	Qty/pkg	Color	bar*	psi*
JR-T-5999-M3	1/16"	0.13	3 m	Red	420	6100
JR-T-5999-M10	1/16"	0.13	10 m	Red	420	6100
JR-T-6000-M3	1/16"	0.18	3 m	Yellow	400	5800
JR-T-6000-M10	1/16"	0.18	10 m	Yellow	400	5800
JR-T-6001-M3	1/16"	0.25	3 m	Blue	386	5600
JR-T-6001-M10	1/16"	0.25	10 m	Blue	386	5600
JR-T-6002-M3	1/16"	0.50	3 m	Orange	350	4500
JR-T-6002-M10	1/16"	0.50	10 m	Orange	350	4500
JR-T-6003-M3	1/16"	0.75	3 m	Green	240	3500
JR-T-6003-M10	1/16"	0.75	10 m	Green	240	3500
JR-T-60031-M3	1/16"	1.00	3 m	Grey	165	2400
JR-T-60031-M10	1/16"	1.00	10 m	Grey	165	2400
JR-T-60032-M3	1/16"	1.40	3 m	Black	52	750
JR-T-60032-M10	1/16"	1.40	10 m	Black	52	750

\* = Recommended maximum permanent working pressure (measured with ACN/Water 1:1 at room temperature)  
Other pack sizes and tubing lengths are available on request.  
Please contact your local distributor or VICI directly.

### PEEK Tubing – Solid Color-Coded

Part No.	OD	ID (mm)	Qty/pkg	Color	bar*	psi*
JR-T-6007-M3	1/16"	0.13	3 m	Red	420	6100
JR-T-6007-M10	1/16"	0.13	10 m	Red	420	6100
JR-T-6008-M3	1/16"	0.18	3 m	Yellow	400	5800
JR-T-6008-M10	1/16"	0.18	10 m	Yellow	400	5800
JR-T-6009-M3	1/16"	0.25	3 m	Blue	386	5600
JR-T-6009-M10	1/16"	0.25	10 m	Blue	386	5600
JR-T-6010-M3	1/16"	0.50	3 m	Orange	350	4500
JR-T-6010-M10	1/16"	0.50	10 m	Orange	350	4500
JR-T-6011-M3	1/16"	0.75	3 m	Green	240	3500
JR-T-6011-M10	1/16"	0.75	10 m	Green	240	3500

\* = Recommended maximum permanent working pressure (measured with ACN/Water 1:1 at room temperature)  
Other pack sizes and tubing lengths are available on request.  
Please contact your local distributor or VICI directly.

### Length Conversions


0.13 mm ≈ .005"  
0.18 mm ≈ .007"  
0.25 mm ≈ .010"  
0.50 mm ≈ .020"  
0.75 mm ≈ .030"  
1.00 mm ≈ .039"  
1.40 mm ≈ .055"



### PEEK Tubing – Natural

Part No.	OD	ID (mm)	Qty/pkg	bar*	psi*
JR-T-6060-M3	0.36 mm	0.075	3 m	358	5200
JR-T-6060-M10	0.36 mm	0.075	10 m	358	5200
JR-T-5993-M3	1/32"	0.13	3 m	379	5500
JR-T-5993-M10	1/32"	0.13	10 m	379	5500
JR-T-5994-M3	1/32"	0.18	3 m	351	5100
JR-T-5994-M10	1/32"	0.18	10 m	351	5100
JR-T-5995-M3	1/32"	0.25	3 m	310	4500
JR-T-5995-M10	1/32"	0.25	10 m	310	4500
JR-T-60041-M3	1/8"	0.75	3 m	345	5000
JR-T-60041-M10	1/8"	0.75	10 m	345	5000
JR-T-6004-M3	1/8"	1.59	3 m	224	3250
JR-T-6004-M10	1/8"	1.59	10 m	224	3250
JR-T-60042-M3	1/8"	2.00	3 m	165	2400
JR-T-60042-M10	1/8"	2.00	10 m	165	2400
JR-T-6006	1/4"	3.17	per m	227	3300

\* = Recommended maximum permanent working pressure (measured with ACN/Water 1:1) at room temperature.  
Other pack sizes and tubing lengths are available on request.  
Please contact your local distributor or VICI directly.



### SPECS

**Material**  
PEEK

**Dimensions**  
OD and ID, see chart

**Tolerances**  
For OD 1/16" PEEK tubing  
OD: +/- 0.05 mm (.002")  
ID: +/- 0.05 mm (.002")

For OD 1/8" PEEK tubing  
OD: +/- 0.10 mm (.004")  
ID: +/- 0.10 mm (.004")

**Pressure rating**  
See chart

**Chemical resistance**  
See chart on page 106

Max. recommended working temp.: < 100 °C (continuous)  
for 1/16" OD tubing with ID up to 0.75 mm

**Spares & Tools**  
A clean burr-free perpendicular cut can be achieved with the VICI Jour Clean-Cut tubing cutter JR-797 see page 22

To bend PEEK tubing at the optimum radius, use our Tubing Elbows on page 21

We recommend PEEK finger-tight and flangeless fittings for your applications  
Fingertights see pages 27–30  
Flangeless see pages 49–53

VICI Jour Sample Loops for analytical applications are made of Premium Grade Tubing (see pages 96–100)

### Length Conversions

0.064 mm	≈ .0025"
0.075 mm	≈ .0030"
0.13 mm	≈ .005"
0.18 mm	≈ .007"
0.25 mm	≈ .010"
0.36 mm	≈ .014"
0.75 mm	≈ .030"
1.59 mm	≈ .062"
2.00 mm	≈ .079"
3.17 mm	≈ .125"

# Tubing for Chromatography

## PEEKclad-FST™ – PEEK-Clad Fused Silica Tubing

### SPECS

#### Material

PEEK, Fused Silica

#### Dimensions

See chart

#### Tolerances

For ID 50 µm +/- 3 µm,  
For IDs 100–150 µm +/- 5 µm

#### Pressure Rating

586 bar (8500 psi)

#### Spares & Tools

PEEK fingertights 1/16"  
see pages 28–30  
Nanovalome® fittings 1/32"  
see pages 38–39

We recommend PEEKclad-FST™ together with our Nanovalome® fittings on pages 38–39

#### Tech Tip

To ensure Zero Dead Volume connections, special cutting techniques and equipment are required. For this reason PEEKclad-FST™ is available only in a variety of pre-cut lengths.

## PEEKclad-FST™ – PEEK-Clad Fused Silica Tubing

- For Zero Dead Volume Connections
- Perfect for Capillary LC
- Inert Smooth Flow Path
- Available with 1/32" and 1/16" OD and 25 up to 300 µm ID

This product combines the inert and smooth inner surface of fused silica with very precise IDs, and the flexibility and easy handling of PEEK. Due to its mechanical strength it can be used both with polymeric or metal ferrules.

The ends are cut perfectly square with sophisticated equipment and polished for Zero Dead Volume connections in capillary LC. In contrast to a polymer sleeve/fused silica connection there is no possibility of liquid filling the gap between the sleeve and the tubing.

All this may result in a lower carry over or cross contamination between samples which can lead to improved reproducibility. The smooth inner wall surface

gives lower band broadening and therefore higher efficiency and resolution.


PEEKclad-FST™ is compatible with most organic solvents and strong acids (pH 0–10) but not compatible with hydrofluoric acid.

In fitting details the PEEK outside surface may be wetted and thus the chemical compatibility of PEEK has to be considered.

PEEK coated fused silica tubing can be manufactured in different proprietary processes. Trade names are PEEKclad-FST™ by Valco International Co. Inc. and VICI International and PEEKsil™ by SGE International Pty. Ltd.

## PEEKclad-FST™ – PEEK-Clad Fused Silica Tubing

Part No.	OD	ID	Length/mm	Color	Qty/pkg
JR-TS-3215	1/32"	25 µm	50	Orange	2
JR-TS-32110	1/32"	25 µm	100	Orange	2
JR-TS-32115	1/32"	25 µm	150	Orange	2
JR-TS-32125	1/32"	25 µm	250	Orange	2
JR-TS-32150	1/32"	25 µm	500	Orange	2
JR-TS-3225	1/32"	50 µm	50	Natural	2
JR-TS-32210	1/32"	50 µm	100	Natural	2
JR-TS-32215	1/32"	50 µm	150	Natural	2
JR-TS-32225	1/32"	50 µm	250	Natural	2
JR-TS-32250	1/32"	50 µm	500	Natural	2
JR-TS-3235	1/32"	75 µm	50	Black	2
JR-TS-32310	1/32"	75 µm	100	Black	2
JR-TS-32315	1/32"	75 µm	150	Black	2
JR-TS-32325	1/32"	75 µm	250	Black	2
JR-TS-32350	1/32"	75 µm	500	Black	2
JR-TS-3245	1/32"	100 µm	50	Red	2
JR-TS-32410	1/32"	100 µm	100	Red	2
JR-TS-32415	1/32"	100 µm	150	Red	2
JR-TS-32425	1/32"	100 µm	250	Red	2
JR-TS-32450	1/32"	100 µm	500	Red	2
JR-TS-3265	1/32"	150 µm	50	Purple	2
JR-TS-32610	1/32"	150 µm	100	Purple	2
JR-TS-32615	1/32"	150 µm	150	Purple	2
JR-TS-32625	1/32"	150 µm	250	Purple	2
JR-TS-32650	1/32"	150 µm	500	Purple	2





Part No.	OD	ID	Length/mm	Color	Qty/pkg
JR-TS-1615	1/16"	25 µm	50	Orange	5
JR-TS-16110	1/16"	25 µm	100	Orange	5
JR-TS-16115	1/16"	25 µm	150	Orange	5
JR-TS-16120	1/16"	25 µm	200	Orange	5
JR-TS-16150	1/16"	25 µm	500	Orange	2
JR-TS-1625	1/16"	50 µm	50	Natural	5
JR-TS-16210	1/16"	50 µm	100	Natural	5
JR-TS-16215	1/16"	50 µm	150	Natural	5
JR-TS-16220	1/16"	50 µm	200	Natural	5
JR-TS-16250	1/16"	50 µm	500	Natural	2
JR-TS-1635	1/16"	75 µm	50	Black	5
JR-TS-16310	1/16"	75 µm	100	Black	5
JR-TS-16315	1/16"	75 µm	150	Black	5
JR-TS-16320	1/16"	75 µm	200	Black	5
JR-TS-16350	1/16"	75 µm	500	Black	2
JR-TS-1645	1/16"	100 µm	50	Red	5
JR-TS-16410	1/16"	100 µm	100	Red	5
JR-TS-16415	1/16"	100 µm	150	Red	5
JR-TS-16420	1/16"	100 µm	200	Red	5
JR-TS-16450	1/16"	100 µm	500	Red	2
JR-TS-1665	1/16"	150 µm	50	Purple	5
JR-TS-16610	1/16"	150 µm	100	Purple	5
JR-TS-16615	1/16"	150 µm	150	Purple	5
JR-TS-16620	1/16"	150 µm	200	Purple	5
JR-TS-16650	1/16"	150 µm	500	Purple	2

Other inner diameters on request.

### SPECS

#### Material

PEEK, Fused Silica

#### Dimensions

See chart

#### Tolerances

For ID 50 µm +/- 3 µm,  
For IDs 100–150 µm +/- 5 µm

#### Pressure Rating

586 bar (8500 psi)

#### Spares & Tools

PEEK fingertights 1/16"  
see pages 28–30  
Nanovalome® fittings 1/32"  
see pages 38–39

We recommend PEEKclad-FST™ together with our Nanovalome® fittings on pages 38–39

#### Tech Tip

To ensure Zero Dead Volume connections, special cutting techniques and equipment are required. For this reason PEEKclad-FST™ is available only in a variety of pre-cut lengths.

# Tubing for Chromatography

## Premium Grade Stainless Steel Tubing

### SPECS

#### Material

SS316  
seamless

#### Dimensions

See chart

#### Tolerances

OD: +/- 0.05 mm (.002")

ID: +/- 0.025 mm (.001")

#### Pressure rating

For 1/16" OD tubing:  
varies with ID, min. 365 bar  
(5300 psi)

For 1/8" OD tubing:  
varies with ID, min. 315 bar  
(4600 psi)

#### Tech Tip

We recommend Pre-Cut  
Stainless Steel Tubing for Zero-  
Dead Volume Connections  
(see page 13)

#### Note

ID dimensions in mm are

## Premium Grade Stainless Steel Tubing

- Premium Grade - Seamless
- Variety of Sizes for HPLC and GC Applications
- Smooth Internal Surface
- Soft Annealed OD for Easy Ferrule Connections

### Premium Grade Stainless Steel Tubing

Part No.	OD	ID (mm)	Qty/pkg
JR-TSS.505-M3	1/32"	0.13	3 m
JR-TSS.505-M10	1/32"	0.13	10 m
JR-TSS.507-M3	1/32"	0.18	3 m
JR-TSS.507-M10	1/32"	0.18	10 m
JR-TSS.510-M3	1/32"	0.25	3 m
JR-TSS.510-M10	1/32"	0.25	10 m
JR-TSS.520-M3	1/32"	0.50	3 m
JR-TSS.520-M10	1/32"	0.50	10 m
JR-TSS105-M3	1/16"	0.13	3 m
JR-TSS105-M10	1/16"	0.13	10 m
JR-TSS110-M3	1/16"	0.25	3 m
JR-TSS110-M10	1/16"	0.25	10 m
JR-TSS115-M3	1/16"	0.38	3 m
JR-TSS115-M10	1/16"	0.38	10 m
JR-TSS120-M3	1/16"	0.50	3 m
JR-TSS120-M10	1/16"	0.50	10 m
JR-TSS130-M3	1/16"	0.75	3 m
JR-TSS130-M10	1/16"	0.75	10 m
JR-TSS140-M3	1/16"	1.00	3 m
JR-TSS140-M10	1/16"	1.00	10 m
JR-TSS230-M3	1/8"	0.75	3 m
JR-TSS230-M10	1/8"	0.75	10 m
JR-TSS240-M3	1/8"	1.00	3 m
JR-TSS240-M10	1/8"	1.00	10 m
JR-TSS260-M3	1/8"	1.52	3 m
JR-TSS260-M10	1/8"	1.52	10 m
JR-TSS267-M3	1/8"	1.78	3 m
JR-TSS267-M10	1/8"	1.78	10 m
JR-TSS285-M3	1/8"	2.16	3 m
JR-TSS285-M10	1/8"	2.16	10 m

### Length Conversions

0.13 mm ≈ .005"  
0.18 mm ≈ .007"  
0.25 mm ≈ .010"  
0.38 mm ≈ .015"  
0.50 mm ≈ .020"  
0.75 mm ≈ .030"  
1.00 mm ≈ .039"  
1.52mm ≈ .060"  
1.78 mm ≈ .070"  
2.16 mm ≈ .085"



### Stainless Steel Tubing

- Variety of Sizes for HPLC and GC Applications
- Smooth Internal Surface
- Soft Annealed OD for Easy Ferrule Connections

### Stainless Steel Tubing

Part No.	OD	ID (mm)	Qty/pkg
JR-T-625-04-M3	1/16"	0.13	3 m
JR-T-625-04-M10	1/16"	0.13	10 m
JR-T-625-05-M3	1/16"	0.18	3 m
JR-T-625-05-M10	1/16"	0.18	10 m
JR-T-625-10-M3	1/16"	0.25	3 m
JR-T-625-10-M10	1/16"	0.25	10 m
JR-T-625-20-M3	1/16"	0.50	3 m
JR-T-625-20-M10	1/16"	0.50	10 m
JR-T-625-30-M3	1/16"	0.75	3 m
JR-T-625-30-M10	1/16"	0.75	10 m
JR-T-625-40-M3	1/16"	1.00	3 m
JR-T-625-40-M10	1/16"	1.00	10 m
JR-T-626-00-M3	1/8"	2.10	3 m
JR-T-626-00-M10	1/8"	2.10	10 m
JR-T-628-00	1/4"	4.65	per m

### SPECS

#### Material

SS316  
welded under inert gas

#### Dimensions

See chart

#### Tolerances

OD:  $\pm 0.05$  mm (.002")  
ID:  $\pm 0.05$  mm (.002")

#### Pressure rating

For 1/16" and 1/8" OD tubing:  
varies with ID, min. 315 bar  
(4600 psi)

#### Tech Tip

We recommend Pre-Cut  
Stainless Steel Tubing for Zero-  
Dead Volume Connections  
(see page 13)



### Length Conversions

0.13 mm  $\approx$  .005"  
0.18 mm  $\approx$  .007"  
0.25 mm  $\approx$  .010"  
0.50 mm  $\approx$  .020"  
0.75 mm  $\approx$  .030"  
1.00 mm  $\approx$  .039"  
2.10 mm  $\approx$  .083"  
4.65 mm  $\approx$  .183"

# Tubing for Chromatography

## Electroformed Nickel Tubing

### SPECS

#### Material

Electroformed Nickel (EFNI)

#### Dimensions

See chart

#### Tolerances

OD:  $\pm 0.025$  mm

ID:  $\pm 0.025$  mm

#### Note

ID dimensions in mm are converted from nominal inch dimensions

## Electroformed Nickel Tubing

- Variety of Sizes for HPLC and GC Applications
- Extremely Smooth Internal Surface
- Ideal for Transfer Lines, Loops and Columns

This small bore tubing is made by electroplating nickel over a diamond drawn mandrel in a continuous process. When the mandrel is removed from the tubing, what's left is in an extremely inert and smooth interior surface – an incredible 1–2 microinch finish. The mirrorlike interior means that EFNI can be used instead of fused silica glass-lined tubing, or silica-lined tubing.

The highly inactive, non-adsorptive surface makes EFNI tubing ideal for transfer lines, sample loops, or columns, or any application which requires minimum carryover potential.

The ID of commercial drawn tubing is typically  $\pm 10\%$  of the nominal size. This variance added to the voids, pits, and striations in the wall leads to a large uncertainty in calculated volumes. Our EFNI tubing is held to  $\pm 0.025$  mm of the nominal size, and its microsmooth surface means that the absolute volume can be calculated very accurately.

Electroformed nickel tubing is electrolytically cut, electropolished, and steam cleaned, ready for use.

## Electroformed Nickel Tubing

Part No.	OD	ID (mm)	Qty/pkg
JR-TEFNI.502-M1	1/32"	0.05	1 m
JR-TEFNI.504-M1	1/32"	0.10	1 m
JR-TEFNI.505-M1	1/32"	0.13	1 m
JR-TEFNI.510-M3	1/32"	0.25	3 m
JR-TEFNI.510-M10	1/32"	0.25	10 m
JR-TEFNI.515-M3	1/32"	0.38	3 m
JR-TEFNI.515-M10	1/32"	0.38	10 m
JR-TEFNI.520-M3	1/32"	0.50	3 m
JR-TEFNI.520-M10	1/32"	0.50	10 m
JR-TEFNI130-M3	1/16"	0.75	3 m
JR-TEFNI130-M10	1/16"	0.75	10 m
JR-TEFNI140-M3	1/16"	1.00	3 m
JR-TEFNI140-M10	1/16"	1.00	10 m

### Length Conversions

0.05 mm  $\approx$  .0020"

0.10 mm  $\approx$  .0040"

0.13 mm  $\approx$  .005"

0.25 mm  $\approx$  .010"

0.38 mm  $\approx$  .015"

0.50 mm  $\approx$  .020"

0.75 mm  $\approx$  .030"

1.00 mm  $\approx$  .040"





### Pre-Cut Premium Stainless Steel Tubing

- Zero Dead Volume Connections
- Ready to Use
- Cleaned
- Soft Annealed OD for Easy Ferrule Connections

VICI offers pre-cut stainless steel tubing in lengths most commonly required in HPLC systems. These tubes are cut to square, burr-free ends for Zero Dead Volume connections.

VICI pre-cut tubing is electrolytically cut and specially steam cleaned to remove both organic and inorganic contaminants. Softener-free Polyethylene caps are used to avoid contamination of the tubing.



### Pre-cut Stainless Steel Tubing

Part No.	OD	ID (mm)	Color	Length (mm)
JR-T-97005	1/16"	0.13	Red	50
JR-T-97010	1/16"	0.13	Red	100
JR-T-97015	1/16"	0.13	Red	200
JR-T-97020	1/16"	0.13	Red	300
JR-T-97025	1/16"	0.18	Yellow	50
JR-T-97030	1/16"	0.18	Yellow	100
JR-T-97035	1/16"	0.18	Yellow	200
JR-T-97040	1/16"	0.18	Yellow	300
JR-T-97045	1/16"	0.25	Blue	50
JR-T-97055	1/16"	0.25	Blue	100
JR-T-97060	1/16"	0.25	Blue	200
JR-T-97065	1/16"	0.25	Blue	300
JR-T-97075	1/16"	0.50	Orange	50
JR-T-97080	1/16"	0.50	Orange	100
JR-T-97085	1/16"	0.50	Orange	200
JR-T-97090	1/16"	0.50	Orange	300
JR-T-97095	1/16"	0.75	Green	50
JR-T-97100	1/16"	0.75	Green	100
JR-T-97105	1/16"	0.75	Green	200
JR-T-97110	1/16"	0.75	Green	300

Australian Distributors: **TECH** *no* **logy**   
 HROM *alytic*   
 www.chromtech.net.au   
 sales@chromtech.net.au   
 ABN 14 643 445 058 **PTY LTD**   
 Tel: (03) 9762 2034 Fax: 61 3 9761 1169   
 ... for ALL Your chromatography supplies !

### SPECS

#### Material

Tubing: SS316L,  
seamless  
Caps: Polyethylene

#### Dimensions

OD: 1/16", ID see chart

#### Tolerances

OD:  $\pm 0.05$  mm (.002")  
ID:  $\pm 0.025$  mm (.001")

#### Pressure rating

For 1/16" OD tubing:  
varies with ID, min. 365 bar  
(5300 psi)

For 1/8" OD tubing:  
varies with ID, min. 315 bar  
(4600 psi)

#### Spares & Tools

We recommend our Valco  
Stainless Steel Fittings on  
page 40.

### Length Conversions

0.13 mm  $\approx$  .005"  
 0.18 mm  $\approx$  .007"  
 0.25 mm  $\approx$  .010"  
 0.50 mm  $\approx$  .020"  
 0.75 mm  $\approx$  .030"



### PTFE & PTFE-Like Tubing

- Low and Mid Pressure Applications
- Biocompatible

#### PTFE Tubing

PTFE tubing (PolyTetraFluoroEthylene) is chemically inert and suitable for lower pressure applications. Typically 1/16" OD tubing holds up to 62 bar (900 psi) and 1/8" up to 35 bar (500 psi).

For easy identification of the ID and to differentiate the transfer lines, PTFE tubing is available in a color-coded version.

#### Black PTFE Tubing

PTFE is too translucent and for certain light sensitive applications the use of black PTFE tubing is recommended. Typical applications:

- Light sensitive enzymatic reactions
- Luminescent reactions
- Connections to the flow-cell
- Radioactive HPLC

#### PFA Tubing

PFA tubing (PerFluoroAlkoxy) has excellent chemical stability and mechanical strength. Always use PFA instead of PTFE if gas permeability is an issue. 1/16" OD high purity tubing has a maximum pressure rating of up to 100 bar (1450 psi), while 1/8" OD maximum 72 bar (1050 psi).

#### FEP Tubing

FEP (FluorinatedEthylenePropylene) is chemically inert to most HPLC solvents and is recommended for low pressure applications. Typically 1/16" OD holds up to 117 bar (1700 psi) and 1/8" holds up to 70 bar (1000 psi). We recommend FEP for ion chromatography applications. Max. operating temperature for FEP is 80 °C.

#### ETFE Tubing

ETFE tubing has some outstanding properties which make it particularly suitable for HPLC since it has a higher burst pressure when compared to PTFE. The most common dimension is 1/16" x 0.25 mm which will withstand up to 186 bar (2700 psi) using water as mobile phase. It is important, however, to note that when using ETFE with organic solvents it may swell and will reduce the max. pressure to 70 bar (1000 psi). It is the ideal polymer for applications with chloride salts, which typically corrode stainless steel.



### PTFE Tubing – Natural

Part No.	OD	ID (mm)	Qty/pkg	bar*	psi*
JR-T-6805-M3	1/16"	0.18	3 m	62	900
JR-T-6805-M10	1/16"	0.18	10 m	62	900
JR-T-6805-M25	1/16"	0.18	25 m	62	900
JR-T-4011-M3	1/16"	0.25	3 m	55	800
JR-T-4011-M10	1/16"	0.25	10 m	55	800
JR-T-4011-M25	1/16"	0.25	25 m	55	800
JR-T-4183-M3	1/16"	0.50	3 m	50	700
JR-T-4183-M10	1/16"	0.50	10 m	50	700
JR-T-4183-M25	1/16"	0.50	25 m	50	700
JR-T-4036-M3	1/16"	0.75	3 m	37	550
JR-T-4036-M10	1/16"	0.75	10 m	37	550
JR-T-4036-M25	1/16"	0.75	25 m	37	550
JR-T-6807-M3	1/16"	1.00	3 m	25	350
JR-T-6807-M10	1/16"	1.00	10 m	25	350
JR-T-6807-M25	1/16"	1.00	25 m	25	350
JR-T-6800-M3	1/8"	1.59	3 m	35	500
JR-T-6800-M10	1/8"	1.59	10 m	35	500
JR-T-6800-M25	1/8"	1.59	25 m	35	500
JR-T-4037-M3	1/8"	2.40	3 m	18	250
JR-T-4037-M10	1/8"	2.40	10 m	18	250
JR-T-4037-M25	1/8"	2.40	25 m	18	250
JR-T-6801-M3	2.00 mm	1.70	3 m	10	150
JR-T-6801-M10	2.00 mm	1.70	10 m	10	150
JR-T-6801-M25	2.00 mm	1.70	25 m	10	150
JR-T-4039	4.00 mm	3.00	per m	17	250
JR-T-4041	4.76 mm	3.76	per m	14	200
JR-T-6810	1/4"	4.75	per m	17	250

### SPECS

#### Material PTFE

#### Dimensions See chart

#### Tolerances

For OD 1/16" tubing  
OD: +/- 0.05 mm (.002")  
ID: +/- 0.05 mm (.002")

For OD 1/8" tubing  
OD: +/- 0.10 mm (.004")  
ID: +/- 0.10 mm (.004")

#### Pressure rating

See chart on page 106

#### Chemical resistance

See chart on page 22

#### Special Info

Low pressure applications:  
Check carefully for gas permeability

#### Spares & Tools

The VICI Jour Clean-Cut tubing cutter gives you burr-free perpendicular cuts JR-797 page 21

We recommend PEEK finger-tight and flangeless fittings for your applications  
Fingertights pages 27–30  
Flangeless pages 49–53

#### Length Conversions

0.18 mm	≈ .007"
0.25 mm	≈ .010"
0.50 mm	≈ .020"
0.75 mm	≈ .030"
1.00 mm	≈ .039"
1.59 mm	≈ .062"
2.00 mm	≈ .079"
2.40 mm	≈ .094"
3.00 mm	≈ .118"
3.76 mm	≈ .148"
4.00 mm	≈ .157"
4.75 mm	≈ .187"

C0583





# Tubing for Chromatography

## PTFE & PTFE-Like Tubing

### SPECS

**Material**  
PTFE, PFA

**Dimensions**  
See chart

**Tolerances**  
For OD 1/16" tubing  
OD: +/- 0.05 mm (.002")  
ID: +/- 0.05 mm (.002")

For OD 1/8" tubing  
OD: +/- 0.10 mm (.004")  
ID: +/- 0.10 mm (.004")

**Pressure rating**  
See chart on page 106

**Chemical resistance**  
See chart on page 22

**Special Info**  
Low pressure applications:  
Check carefully for gas permeability

**Spares & Tools**  
The VICI Jour Clean-Cut tubing cutter gives you burr-free perpendicular cuts JR-797 page 21

We recommend PEEK fingertight and flangeless fittings for your applications  
Fingertights pages 27-30  
Flangeless pages 49-53

### PTFE Tubing – Solid Color Coded

Part No.	OD	ID (mm)	Qty/pkg	Color	bar*	psi*
JR-T-4011C-M3	1/16"	0.25	3 m	Blue	55	800
JR-T-4011C-M10	1/16"	0.25	10 m	Blue	55	800
JR-T-4011C-M25	1/16"	0.25	25 m	Blue	55	800
JR-T-4183C-M3	1/16"	0.50	3 m	Orange	50	700
JR-T-4183C-M10	1/16"	0.50	10 m	Orange	50	700
JR-T-4183C-M25	1/16"	0.50	25 m	Orange	50	700
JR-T-4036C-M3	1/16"	0.75	3 m	Green	38	550
JR-T-4036C-M10	1/16"	0.75	10 m	Green	38	550
JR-T-4036C-M25	1/16"	0.75	25 m	Green	38	550
JR-T-6808-M3	1/16"	0.25	3 m	Black	55	800
JR-T-6808-M10	1/16"	0.25	10 m	Black	55	800
JR-T-6808-M25	1/16"	0.25	25 m	Black	55	800
JR-T-6811-M3	1/16"	0.50	3 m	Black	50	700
JR-T-6811-M10	1/16"	0.50	10 m	Black	50	700
JR-T-6811-M25	1/16"	0.50	25 m	Black	50	700

\* = Recommended maximum permanent working pressure  
Other pack sizes and tubing lengths are available on request. Please contact your local distributor or VICI directly.

### PFA Tubing

Part No.	OD	ID (mm)	Qty/qty	bar*	psi*
JR-T-4001-M3	1/16"	0.50	3 m	100	1450
JR-T-4001-M10	1/16"	0.50	10 m	100	1450
JR-T-4001-M25	1/16"	0.50	25 m	100	1450
JR-T-4002-M3	1/16"	0.75	3 m	76	1100
JR-T-4002-M10	1/16"	0.75	10 m	76	1100
JR-T-4002-M25	1/16"	0.75	25 m	76	1100
JR-T-4007-M3	1/16"	1.00	3 m	55	800
JR-T-4007-M10	1/16"	1.00	10 m	55	800
JR-T-4007-M25	1/16"	1.00	25 m	55	800
JR-T-4003-M3	1/8"	1.59	3 m	72	1050
JR-T-4003-M10	1/8"	1.59	10 m	72	1050
JR-T-4003-M25	1/8"	1.59	25 m	72	1050

\* = Recommended maximum permanent working pressure  
Other pack sizes and tubing lengths are available on request. Please contact your local distributor or VICI directly.

### Length Conversions

0.25 mm = .010"  
0.36 mm = .014"  
0.50 mm = .020"  
0.75 mm = .030"  
1.00 mm = .039"  
1.59 mm = .062"



### FEP Tubing

Part No.	OD	ID (mm)	Qty/pkg	bar*	psi*
JR-T-6812-M3	1/16"	0.25	3 m	117	1700
JR-T-6812-M10	1/16"	0.25	10 m	117	1700
JR-T-6812-M25	1/16"	0.25	25 m	117	1700
JR-T-6802-M3	1/16"	0.50	3 m	96	1400
JR-T-6802-M10	1/16"	0.50	10 m	96	1400
JR-T-6802-M25	1/16"	0.50	25 m	96	1400
JR-T-6803-M3	1/16"	0.75	3 m	76	1100
JR-T-6803-M10	1/16"	0.75	10 m	76	1100
JR-T-6803-M25	1/16"	0.75	25 m	76	1100
JR-T-6806-M3	1/16"	1.00	3 m	52	750
JR-T-6806-M10	1/16"	1.00	10 m	52	750
JR-T-6806-M25	1/16"	1.00	25 m	52	750
JR-T-6804-M3	1/8"	1.59	3 m	69	1000
JR-T-6804-M10	1/8"	1.59	10 m	69	1000
JR-T-6804-M25	1/8"	1.59	25 m	69	1000
JR-T-6809	1/4"	4.35	per m	45	650

\* = Recommended maximum permanent working pressure

Other pack sizes and tubing lengths are available on request. Please contact your local distributor or VICI directly.

### ETFE Tubing

Part No.	OD	ID (mm)	Qty/pkg	bar*	psi*
JR-T-084-M3	1/32"	0.25	3 m	152	2200
JR-T-084-M10	1/32"	0.25	10 m	152	2200
JR-T-084-M25	1/32"	0.25	25 m	152	2200
JR-T-078-M3	1/16"	0.17	3 m	200	2900
JR-T-078-M10	1/16"	0.17	10 m	200	2900
JR-T-078-M25	1/16"	0.17	25 m	200	2900
JR-T-080-M3	1/16"	0.25	3 m	186	2700
JR-T-080-M10	1/16"	0.25	10 m	186	2700
JR-T-080-M25	1/16"	0.25	25 m	186	2700
JR-T-082-M3	1/16"	0.50	3 m	152	2200
JR-T-082-M10	1/16"	0.50	10 m	152	2200
JR-T-082-M25	1/16"	0.50	25 m	152	2200
JR-T-083-M3	1/16"	0.75	3 m	117	1700
JR-T-083-M10	1/16"	0.75	10 m	117	1700
JR-T-083-M25	1/16"	0.75	25 m	117	1700
JR-T-085-M3	1/16"	1.00	3 m	83	1200
JR-T-085-M25	1/16"	1.00	10 m	83	1200
JR-T-085-M10	1/16"	1.00	25 m	83	1200
JR-T-086-M3	1/8"	1.59	3 m	110	1600
JR-T-086-M10	1/8"	1.59	10 m	110	1600
JR-T-086-M25	1/8"	1.59	25 m	110	1600

\* = Recommended maximum permanent working pressure

Other pack sizes and tubing lengths are available on request. Please contact your local distributor or VICI directly.

### SPECS

**Material**  
FEP, ETFE

**Dimensions**  
See chart

**Tolerances**  
For OD 1/16" tubing  
OD: +/- 0.05 mm (.002")  
ID: +/- 0.05 mm (.002")

For OD 1/8" tubing  
OD: +/- 0.10 mm (.004")  
ID: +/- 0.10 mm (.004")

**Pressure rating**  
See chart on page 106

**Chemical resistance**  
See chart on page 22

**Special Info**  
Low pressure applications:  
Check carefully for gas permeability

**Spares & Tools**  
The VICI Jour Clean-Cut tubing cutter gives you burr-free perpendicular cuts JR-797 page 21

We recommend PEEK fingertight and flangeless fittings for your applications  
Fingertights pages 27-30  
Flangeless pages 49-53

### Length Conversions

0.17 mm ≈ .0067"  
0.25 mm ≈ .010"  
0.50 mm ≈ .020"  
0.75 mm ≈ .030"  
1.00 mm ≈ .039"  
1.59 mm ≈ .062"  
4.35 mm ≈ .171"



## No-Ox Tubing

### SPECS

#### Material

Tubing: FEP + PVDF

Insert: KEL-F

Ferrule: ETFE

Fitting: PEEK

#### Dimensions

OD: 1/8"

ID: 1.65 mm

#### Tolerances

OD: +/- 0.05 mm (.002")

ID: +/- 0.05 mm (.002")

#### Pressure rating

Low pressure only

#### Special info

Use the special designed KEL-F insert (P/N JR-6141-10) to prevent solvent contact of the PVDF layer

Due to the limited chemical resistance of the outer PVDF layer, the No-Ox tubing should not be used inside a solvent reservoir.

## No-Ox Tubing

- Prevents 'Re-gassing'
- Excellent Chemical Resistance

Most HPLC installations use PTFE tubing as internal low-pressure transfer lines. PTFE is easy to use and chemically extremely inert. But most PTFE and PTFE like polymers have micro pores which allow gases to diffuse freely.

Most solvent degassers are based on this effect. Gases dissolved in solvents are removed by diffusion in applying a vacuum on the outside of the tubing.

A contrary effect can be observed when degassed solvent will pass through PTFE tubing with standard air pressure on the outside. A 're-gassing' will be observed, gases will diffuse through the PTFE tubing back into the newly and costly degassed solvent.

We therefore created the No-Ox tubing. A FEP (PTFE like) tubing is coated on the outside with PVDF tubing. FEP has excellent chemical resistance, but micro

pores are present in the tubing wall. PVDF on the other hand is a very tight polymer, no micro pores can be observed, but it lacks chemical resistance.

Co-extruded No-Ox tubing will combine the two features of FEP and PVDF in the most positive manner, ensuring that degassed solvent will stay degassed until it reaches the HPLC pump. The PVDF-layer can easily be removed (using e.g. VICI Jour P/N 797 Clean-Cutter or a razor blade) and stripped from the FEP so only PTFE will be immersed by solvent i.e. in the solvent reservoir.

The No-Ox tubing is connected with a special KEL-F tubing insert (P/N 6141) to prevent solvent contact of the PVDF layer at the tubing end. The ferrule fits into the ID of the tubing and seals against the bottom of the port.



## No-Ox Tubing

Part No.	Description	Qty/pkg
JR-T-6130-M3	Tubing, No-Ox 1/8" x 1.65 mm ID	3 m
JR-T-6130-M10	Tubing, No-Ox 1/8" x 1.65 mm ID	10 m
JR-6140	No-Ox fitting Kit	1
	contains: 2 pcs KEL-F inserts, 2 pcs PEEK nuts and 2 pcs ETFE Ferrules	

## Spare Parts


Part No.	Description	Qty/pkg
JR-6141-10	Insert, KEL-F, for No-Ox fitting	10
JR-55071-10	Nut, PPS, flangeless 1/8", 1/4"-28	10
JR-55083-10	Adapter, PP, fingertight sleeve, black	10
JR-051-10	Ferrule, ETFE, 1/8"	10

### Thermo-Formed Polymeric Tubing


- Various Polymeric Tubing
- Telephone Cables and Customer Specific Items
- Axial Flexibility

VICI offers thermo-formed tubing in various of geometric shapes. Straight tubing, spirals or even more complex forms, just define the form and material on a drawing.

Contact you local distributor or VICI directly for a quotation.



PDF 810 000





# Tubing for Chromatography

## Knitted PTFE Reaction coils

### SPECS

**Material**  
PTFE

**Dimensions**  
OD: 1/16"  
ID: see chart

**Tolerances**  
OD:  $\pm 0.05$  mm (.002")  
ID:  $\pm 0.05$  mm (.002")

**Pressure rating**  
< 55 bar (< 800 psi)

**Recommended Volumetric Flow Rate:**

ID	Lin. Velocity	Vol. Flow
mm	cm/s	ml/min
0.25	10	0.29
0.50	10	1.2
0.75	10	2.7

### Special Info

Low pressure applications:  
Check carefully for gas permeability

## Knitted PTFE Reaction Coils

- Tortuous Flow Path of Optimal Flow Dynamic Design made from Inert PTFE
- Minimized Axial Dispersion combined with Efficient Radial Mixing
- Available in Different Lengths and IDs



### Preserved chromatographic peak shape

The Knitted PTFE Reaction coils (KRC) are made from tubing which has been knitted into a tortuous path to force the fast moving centre of the liquid stream to mix radially with the slower moving boundary liquid layer, thereby minimizing the axial dispersion. The shape of a chromatographic peak entering the reactor is thus preserved very efficiently. To ensure a stable radial mixing within the KRC, a linear flow rate of 10 cm/s or more is recommended. The KRC is usually the most optimal delay element for use in analytical flow systems, and is typically inserted in the flow path to create a delay line, so that a reaction or other event that requires a certain time can take place. Different delay or reaction times are accomplished by changing the inner diameter and length of the KRC, taking the flow rate through the KRC into consideration.

### Example application areas

KRC can be applied as a delay line in chromatographic separations. One example is the parallel coupling of a MS detector and a fraction collector in preparative separations (see left figure below). The delay time induced by the KRC allows time for detection and therefore an intelligent decision between fraction collection or waste without loss of chromatographic efficiency.

The KRC reactors can also be used in post-column reaction detection in HPLC (see right figure below). In this set-up the KRC also ensures thorough mixing between the column effluent and the added reagents.



## Knitted PTFE Reaction coils


Part No.	ID (mm)	Tubing length
JR-T-3000-123	0.25	1 m
JR-T-3000-223	0.25	2 m
JR-T-3000-423	0.25	4 m
JR-T-3000-133	0.25	10 m
JR-T-3000-233	0.25	20 m
JR-T-3000-125	0.50	1 m
JR-T-3000-225	0.50	2 m
JR-T-3000-425	0.50	4 m
JR-T-3000-135	0.50	10 m
JR-T-3000-235	0.50	20 m
JR-T-3000-128	0.75	1 m
JR-T-3000-228	0.75	2 m
JR-T-3000-428	0.75	4 m
JR-T-3000-138	0.75	10 m
JR-T-3000-238	0.75	20 m

### Length Conversions

0.25 mm  $\approx$  .010"  
0.50 mm  $\approx$  .020"  
0.75 mm  $\approx$  .030"

### PEEK Tubing Elbows

- Ideal for Routing 1/16" OD PEEK Tubing
- Optimum Radius, Prevents Kinking
- 90° and 180° Elbows



PEEK Tubing Elbows

Part No.	Description
JR-357090	Tubing Elbow 90°
JR-357180	Tubing Elbow 180°

### SPECS

**Material**  
Polypropylene

**Dimension**  
For OD 1/16" PEEK Tubing

### Tubing Clip – The LC Tubing Organizer

- Fixes 1/16" and 1/8" to Beakers, Flasks or Bottles

Here is a convenient small device for the busy chromatographer. The Tubing Clip, which holds both your 1/16" and 1/8" polymer tubing precisely where they should be in your beakers, flasks and bottles etc. up to 4 mm wall thickness.

The stainless steel spring helps to keep your Tubing Clip working for a long period of time and helps to keep your solvent line in the solvent.



### Tubing Clip

Part No.	Description
JR-9001	Tubing Clip

### SPECS

**Material**  
Polypropylene  
Stainless Steel spring

**Dimensions**  
For OD 1/16" and 1/8" tubing

**Special Info**  
Clip fitted with a SS spring suitable for glassware up to 4 mm thickness



### SPECS

#### Note

The VICI Jour Clean-Cut Tubing Cutter Part No. JR-797 is not for sale in Japan. For an alternative please see part No. JR-794 Guillotine Cutter.

### Clean-Cut Polymer Tubing Cutter

- Cuts PEEK, PTFE, ETFE and other Polymeric Tubing
- No Distortion of OD or Closing of ID
- Safety Lock Secures the Blade

It is difficult to obtain burr-free perpendicular cuts on polymeric tubing, but right angles and clean cuts are essential in LC to avoid leaks. The Clean-Cut is a tool for the most demanding LC applications which will do the job without distorting the outside diameter and/or closing the inside diameter of the tubing.

It is designed to cut PTFE, ETFE and polymer tubing in general and PEEK capillary tubing in particular. The compact design allows every chromatographer to have it in his pocket. A unique safety locking mechanism secures the blade when not in use.




### Clean-Cut Polymer Tubing Cutter

Part No.	Description
JR-797	Clean-Cut Tubing Cutter
JR-798	Replacement blade for JR-797

### Guillotine Polymer Tubing Cutter

- Cuts PEEK, PTFE, ETFE and other Polymeric Tubing



### Guillotine Polymer Tubing Cutter


Part No.	Description
JR-794	Guillotine Cutter for polymer tubing
JR-795	Replacement blade for JR-794

### Stainless Steel Tubing Cutter (for non-critical connections)

- For 1/16" and 1/8" SS Tubing
- Leaves Smooth, Uniform Cuts
- Minimum Deburring and/or Reaming Required
- Easy Replacement of the Cutting Wheel

This is the ideal tool for cutting 1/16" and 1/8" stainless steel tubing, with an ID > 0.5 mm, leaving smooth, uniform cuts with a minimum of deburring and/or

reaming for non-critical connections. The easily replaced cutting wheel scores tubing.



### SPECS

#### Tech Tip

Between sample injector, column and detector electrolytically cut and steam cleaned SS tubing should be used

SS Pre-Cut tubing (see page 13)

Part No.	Description
JR-792	Tubing Cutter for metal tubing
JR-793	Replacement cutting wheel for JR-792

### Stainless Steel Tubing Pliers (for non-critical connections)

- For Easy Cutting 1/16" OD SS Tubing off a Coil
- Reaches "Hard-To-Get" places in a HPLC System

The plier type of cutter can be used for quick cuts for non-critical applications or at 'hard to get' places. However for scientific equipment the use of Pre-Cut and cleaned tubing is strongly recommended.



Part No.	Description
JR-796	Tubing Pliers for metal tubing



### SPECS

#### Spares & Tools

Flanged Fittings

see pages 50 and 54

Polymeric Tubing

see pages 5-9 and 14-18

### Easy-Flange & Easy-Flange Combi-Kits

- Handy Tool for Perfect Flanges
- No Heat or Electricity Required
- No Contamination

The Easy-Flange flange-rolling tool has been developed for fittings requiring a flange and is suitable for PTFE and PTFE like tubes with OD between 1/16" and 1/8".

The tool forms the flange just by applying a mechanical force and does not require electrical power. The quality of the flange is considerably improved since

the flanges are formed without stressing the tubing by heat and because the specially designed negative conical profile of the flange forming part gives a shape ideal for maximum sealing properties. The Easy-Flange-Combi Kit is an upgraded version of our standard Easy-Flange Kit, which comes complete with PTFE tubing, PPS nuts and Polypropylene O-rings.



### Easy-Flange Kit

Part No.	Description	Qty
JR-201540	Easy-Flange, kit in plastic case, complete with the following parts:	
JR-202235	Flanging disc with 0.5 mm SS pin for PEEK tubing	1 pc
JR-201541	Flanging disc with 0.8 mm polymer pin	1 pc
JR-201554	Flanging disc with 0.8 mm Titanium pin	1 pc
JR-201536	Flanging disc with 1.3 mm polymer pin	1 pc
JR-201537	Flanging disc with 1.3 mm Titanium pin	1 pc
JR-797	Tool, Clean-Cut tubing Cutter	1 pc
JR-T-4036-M3	Tubing, PTFE, 1/16" x 0.75 mm ID, 3m/pkg	1 pkg

### Easy-Flange Combi-Kit

Part No.	Description	Qty
JR-201539	Easy-Flange, combi kit in plastic case, complete with the following parts:	
JR-202235	Flanging disc with 0.5 mm SS pin for PEEK tubing	1 pc
JR-201541	Flanging disc with 0.8 mm polymer pin	1 pc
JR-201554	Flanging disc with 0.8 mm Titanium pin	1 pc
JR-201536	Flanging disc with 1.3 mm polymer pin	1 pc
JR-201537	Flanging disc with 1.3 mm Titanium pin	1 pc
JR-797	Tool, Clean-Cut tubing Cutter	1 pc
JR-55060-10	Nut, PPS, flanged 1/16", 1/4"-28, 10/pkg	1 pkg
JR-201926-10	Ring, PP, flanged 1/16", 10/pkg	1 pkg
JR-55061-10	Nut, PPS, flanged 1/8", 1/4"-28, 10/pkg	1 pkg
JR-201928-10	Ring, PP, flanged 1/8", 10/pkg	1 pkg
JR-5508X-24	Adapter, PP, fingertight sleeve, assorted 2 pcs. per color, 24/pkg	1 pkg
JR-T-4036-M3	Tubing, PTFE, 1/16" x 0.75 mm ID, 3m/pkg	1 pkg
JR-T-6800-M3	Tubing, PTFE, 1/8" x 1/16" ID, 3m/pkg	1 pkg

### Easy-Flange Spare Parts & Options

Part No.	Description
JR-201541	Flanging disc with 0.8 mm polymer pin
JR-201554	Flanging disc with 0.8 mm Titanium pin
JR-201662	Flanging disc with 1.0 mm polymer pin
JR-202235	Flanging disc with 0.5 mm SS pin for PEEK tubing
JR-201536	Flanging disc with 1.3 mm polymer pin
JR-201537	Flanging disc with 1.3 mm Titanium pin
JR-201663	Flanging disc with 1.5 mm polymer pin
JR-201531	Flanging clamp pair
JR-797	Tool, Clean-Cut tubing Cutter







C0578

### INDEX

PPS One-Piece Hex-Head/ Fingertight Fitting	27
PEEK Color Coded & Natural One-Piece Fingertight Fittings	28
PEEK One-Piece Hex-Head Fittings	30
PEEK Nuts for use with Double Ferrule	30
PEEK Nuts with Single Ferrule	31
Color-It Fingertight Adapters	32
PEEK Unions	33
PEEK Tees, Crosses and Manifolds High Pressure	34
PEEK Biocompatible Mixing Tee	35
PEEK Tees, Crosses and Manifolds High Pressure	36
PEEK Plugs & Caps	37
PEEK Starter Kit	37
Nanovolume® Fittings	38
PEEK Liners for Fused Silica Tubing Connections	39
Valco Stainless Steel Nuts and Ferrules	40
Stainless Steel Nuts and Ferrules Rheodyne® Type	41
Stainless Steel Unions	42
Stainless Steel Tees & Crosses	43
Stainless Steel Plugs & Caps	44
Stainless Steel Starter Kit	44
ChromBox Kit	45
Tools	46-47

#### Disclaimer:

The maximum holding pressure for any type of connection involving tubing and a ferrule, varies considerably with the tubing material, the ferrule material, the clearance between tubing OD and ferrule ID and the shape of the fitting detail.

## PPS One-Piece Hex-Head/Fingertight Fitting

- Economically Priced Alternative to PEEK
- Excellent Chemical Resistance
- One-Piece – No Ferrule required
- Slim Design – Fingertight with Adapter

The new VICI Jour PPS One-Piece Hex-Head/Fingertight Fitting offers excellent chemical and pressure resistance.



## PPS One-Piece Hexhead/Fingertight Fitting

Part No.	Description	Qty/pkg
JR-58000-5	Fitting, PPS, hexhead one-piece, black, 10-32	5

## Color-Sleeve – fingertight adapter

This adapter is useful to color code different solvent lines and it turns the new VICI Jour low pressure nut into a fingertight fitting in an instant. The Color-Sleeve is available in 12 colors.

## Color-Sleeve – fingertight adapter

Part No.	Description	Qty/pkg
JR-55080-5	Adapter, PP, fingertight sleeve, natural	5
JR-55081-5	Adapter, PP, fingertight sleeve, white	5
JR-55082-5	Adapter, PP, fingertight sleeve, dark grey	5
JR-55083-5	Adapter, PP, fingertight sleeve, black	5
JR-55084-5	Adapter, PP, fingertight sleeve, lavender	5
JR-55085-5	Adapter, PP, fingertight sleeve, red	5
JR-55086-5	Adapter, PP, fingertight sleeve, yellow	5
JR-55087-5	Adapter, PP, fingertight sleeve, orange	5
JR-55088-5	Adapter, PP, fingertight sleeve, brown	5
JR-55089-5	Adapter, PP, fingertight sleeve, green	5
JR-55090-5	Adapter, PP, fingertight sleeve, blue	5
JR-55091-5	Adapter, PP, fingertight sleeve, purple	5
JR-5508X-24	Adapter, PP, fingertight sleeve, assorted 2 pcs. per color	24
JR-5508X-12	Adapter, PP, fingertight sleeve, assorted 2 pcs. each blue, red, green, yellow, black, white	12



## SPECS

## Material

PPS – PolyPhenylene Sulphide, black

Adapter: Polypropylene

## Chemical Resistance PPS

Chemical Class	Resistance
Acids, organic	A
Acids, inorganic	A/C*
Aldehydes	A
Alcohols	A
Bases	A/B
Esters	A
Halogenated Org.	A/B
Hydrocarbons	A
Ketones	A

\*conc. Halogen Acids

A = suitable

B = Marginal – dependant on application

C = Not recommended

Detailed chemical resistance chart see page 106

## Dimensions

For 1/16" OD Tubing

## Threads

10-32

## Pressure Rating

< 350 bar (< 5000 psi)

Varies with tubing material and ID

## Spares &amp; Tools

We offer 1/16" OD tubing for High Pressure applications in various materials (see pages 5–13)



### SPECS

**Material**  
PEEK

**Dimensions**  
For OD 1/16" Tubing

**Threads**  
10-32

**Pressure rating**  
< 350 bar (< 5000 psi)  
Varies with tubing material  
and ID

**Special Info**  
Fingertight Fittings – works  
with all 1/16" OD tubing and  
female 10-32 ports

Color Coded – six colors to  
identify tubing lines

Biocompatible

The coloration pigments  
are carefully chosen not  
to contain any heavy metals  
and hazardous chemicals.

**Spares & Tools**  
We offer 1/16" OD tubing  
for High Pressure applications  
in various materials  
(see pages 5–13)

### PEEK Color Coded & Natural One-Piece Fingertight Fittings

- Fingertight up to 350 bar (5000 psi)
- Excellent Chemical Resistance
- Biocompatible

To identify various tubing lines the PEEK fingertight fittings are available in six colors. PEEK fingertights are rated to 350 bar (5000 psi) and can be used in virtually any type of HPLC fitting detail on the market with 10-32 thread.



### PEEK Color Coded One-Piece Fingertight Fittings Molded Version

Part No.	Description	Qty/pkg
JR-55020-5	Fitting, PEEK, one-piece natural, 10-32	5
JR-55021-5	Fitting, PEEK, one-piece black, 10-32	5
JR-55022-5	Fitting, PEEK, one-piece red, 10-32	5
JR-55023-5	Fitting, PEEK, one-piece yellow, 10-32	5
JR-55024-5	Fitting, PEEK, one-piece blue, 10-32	5
JR-55025-5	Fitting, PEEK, one-piece green, 10-32	5



Part No. JR-55025-5

## PEEK Natural One-Piece Fingertight Fittings – Machined Versions

Part No.	Description	Qty/pkg
JR-5502-5	Fitting, PEEK, one-piece, 10-32	5
JR-5504-5	Fitting, PEEK, one-piece fingertight short head, 10-32	5
JR-5507-5	Fitting, PEEK, one-piece fingertight long, 10-32	5
JR-5508-5	Fitting, PEEK, one-piece fingertight narrow hex-head, 10-32	5



JR-5502-5



JR-5504-5



JR-5507-5



JR-5508-5



## SPECS

**Material**  
PEEK

**Dimensions**  
For OD 1/16" Tubing

**Threads**  
10-32

**Pressure rating**  
< 350 bar (< 5000 psi)  
Varies with tubing material and ID

**Special Info**  
Fingertight Fittings – works with all 1/16" OD tubing and female 10-32 ports

**Spares & Tools**  
We offer 1/16" OD tubing for High Pressure applications in various materials (see pages 5–13)



## Polymeric & Stainless Fittings – High Pressure

### PEEK One-Piece Hex-Head Fittings & PEEK Nuts for use with Double Ferrule

#### SPECS

**Material**  
PEEK

**Dimensions**  
For OD 1/16" Tubing

**Threads**  
10-32

**Pressure rating**  
< 350 bar (< 5000 psi)  
Varies with tubing material  
and ID

**Tech Tip**  
For tightening hex-head  
fittings, we recommend  
our ValvTool JR-800  
(see pages 46 and 95)

#### PEEK One-Piece Hex-Head Fittings

- Ideal for use with PEEK Tubing
- Withstand Pressures up to 350 bar (5000 psi)
- For any High Pressure Fittings

PEEK one-piece hex-head fittings are convenient with their integrated ferrule and easy to use with the Valv-Tool and the Color-It Color-Coded Fingertight Adapters. The 10-32 thread will fit virtually all standard high-pressure female connections. Connections with

the reusable PEEK one-piece hex-head will withstand pressures of up to 350 bar (5000 psi). PEEK one-piece hex-head fittings are available both in a short and long version.



1/4" = 6.35 mm

#### PEEK One-Piece Hex-Head Fittings

Part No.	Description	Qty/pkg
JR-55100-5	Fitting, PEEK, one-piece hex-head short, 10-32	5
JR-55110-5	Fitting, PEEK, one-piece hex-head long, 10-32	5

#### SPECS

**Material**  
PEEK

**Dimensions**  
For OD 1/16" Tubing

**Threads**  
10-32

**Tolerances**  
+/- 0.05 mm (.002")

**Pressure rating**  
< 420 bar (< 6000 psi)  
Varies with tubing material  
and ID

**Special Info**  
Double-Ferrule grips Tubing  
twice

Works with all 1/16" OD  
Tubing

Fits any 10-32 Fitting

Fits all Female compression  
Fittings Ideal for use with  
PEEK Tubing

#### PEEK Nuts for use with Double Ferrule

- High-Pressure Rating up to 420 bar (6000 psi)
- Ideal for use with 1/16" PEEK Tubing
- Fits any 10-32 Port

These hex-head nuts are for use with VICI Jour Double Ferrules. The fingertight version may be hand-tightened for operation at up to 420 bar (6000 psi). The hex-head version can be used for connections that are difficult to reach or closely spaced. These fittings fit virtually any female 1/16" fitting.



#### PEEK Nuts with Double Ferrule

Part No.	Description	Qty/pkg
JR-5511-5	Nut, PEEK, hex-head short, 10-32	5
JR-5510-5	Nut, PEEK, hex-head long, 10-32	5
JR-5580-5	Nut, PEEK, fingertight, 10-32	5
JR-5503-5	Nut, PEEK, universal fingertight, 10-32	5
JR-5004-5	Ferrule, PEEK, double 1/16"	5

## PEEK Nuts with Single Ferrule

- Reuseable Nut, Replaceable Ferrule
- For 1/16" and 1/8" OD Tubing
- Suitable for High Pressure 10-32 and 5/16"-24 Ports

Fingertight nuts have a knurled surface designed to allow sufficient sealing force to be applied to the ferrule without wrenches. Hex style nuts allow wrench

tightening; however, since they are polymeric, they can break and are recommended for use only when space is limited and fingers will not fit.



Part No. JR-5570-5



Part No. JR-5003-5



Part No. JR-ZN2PK-5



Part No. JR-ZF2PK-5

## Two-Piece Fingertight Nuts

Part No.	Description	Qty/pkg
JR-5570-5	Nut, PEEK, fingertight, 1/16", 10-32	5
JR-5003-5	Ferrule, Single, PEEK, 1/16"	5
JR-ZN2PK-5	Nut, PEEK, hex-head, 1/8", 5/16"-24	5
JR-ZF2PK-5	Ferrule, PEEK, 1/8"	5

## SPECS

**Material**  
PEEK

**Dimensions**  
For 1/16", 1/8" Tubing

**Threads**  
See Chart

**Pressure rating**  
< 280 bar (< 4000 psi)  
Varies with tubing material and ID



### SPECS

#### Material

Polypropylene

#### Dimensions

For 1/4" Hex-head nuts

### Color-It Fingertight Adapters

- Easy Solvent Line Identification
- Fits All Hex-Head PEEK and SS Fittings for 1/16" OD Tubing
- Six Different Colors Available

We recommend the Color-It Snap-On fingertight-extension for color coding our 1/4" hex-head nuts. Color-It is useful to color code different solvent lines and it turns your hex-head nuts into a fingertight fitting in an instant. Color-It knobs are available in six

different colors (blue, yellow, green, black, white and red). They can be used in connection with our PEEK and Stainless Steel hex-head nuts for 1/16" OD tubing (see pages 30 and 40–41)



### Color-It Fingertight Adapters

Part No.	Description	Qty/pkg
JR-55010-5	Color-It blue	5
JR-55011-5	Color-It yellow	5
JR-55012-5	Color-It green	5
JR-55013-5	Color-It black	5
JR-55014-5	Color-It white	5
JR-55015-5	Color-It red	5
JR-55016-12	Color-It assorted, 2 pcs. each color	12

## PEEK Unions

- Easy Connection of Transfer/Solvent Lines
- Low Dead Volume
- Completely Inert and Biocompatible

VICI Jour PEEK unions are usable with most common HPLC solvents (organic and inorganic) without risking corrosion, metal contamination, or protein degradation that might occur with stainless steel. Small diameter bores are used to insure low dead volume connections. Connections to these PEEK unions can

be made with any 10-32 fitting that uses polymeric ferrules. PEEK unions are available in "body only" versions as well as complete with our PEEK One-Piece Fingertight Fittings ready to connect to any 1/16" OD tubing (PEEK, stainless steel, PTFE, ETFE, etc.).



Part No. JR-1061



Part No. JR-ZBU1FPK



Part No. JR-ZRU21FPK

## PEEK Unions

Part No.	Description
JR-1061	Union, PEEK, 0.25 mm bore, complete
JR-1066	Union, PEEK, 0.50 mm bore, complete
JR-1067	Union, PEEK, 0.75 mm bore, complete
JR-061	Union, PEEK, 0.25 mm bore, body only
JR-066	Union, PEEK, 0.50 mm bore, body only
JR-067	Union, PEEK, 0.75 mm bore, body only

## PEEK Bulkhead Unions

Part No.	Description
JR-ZBU1CFPK	Bulkhead Union, PEEK, 0.25 mm bore, complete
JR-ZBU1MFPK	Bulkhead Union, PEEK, 0.50 mm bore, complete
JR-ZBU1FPK	Bulkhead Union, PEEK, 0.75 mm bore, complete

## PEEK Reducing Unions

Part No.	Description
JR-ZRU21FPK	Union, PEEK, reducing, 1/8" to 1/16", 0.75 mm bore, complete

## Spare Parts

Part No.	Description	Qty/pkg
JR-55021-S	Fitting, PEEK, one-piece black, 10-32 replacement for all JR-106X	5
JR-ZN2PK-S	Nut, PEEK, hex-head, 1/8", 5/16"-24, replacement for JR-ZRU21FPK	5
JR-ZF2PK-S	Ferrule, PEEK, 1/8", replacement for JR-ZRU21FPK	5

## SPECS

**Material**  
PEEK

**Dimensions**  
For 1/16" OD Tubing and 1/16" to 1/8"  
Others: see illustrations

**Threads**  
10-32 for 1/16" OD tubing  
5/16"-24 for 1/8" OD tubing

**Pressure rating**  
< 350 bar (< 5000 psi)  
Varies with tubing material and ID

**Note**  
Complete unions will be supplied with PEEK one-piece fittings.

PEEK unions have the standardized Valco fitting details with smooth transition of the taper into the threaded part.

**Tech Tip**  
Bulkhead Unions require a min. 10 mm panel hole.



# Polymeric & Stainless Fittings – High Pressure

## PEEK Tees, Crosses and Manifolds High Pressure

### SPECS

**Material**  
PEEK

**Dimensions**  
Tees and Crosses  
Diameter (body): 25.0 mm  
Thickness: 12.0 mm

**Threads**  
10-32  
For OD 1/16" Tubing

**Pressure rating**  
< 350 bar (< 5000 psi)  
Varies with tubing material  
and ID

**Special Info**  
100% PEEK flowpath  
for inertness and biocompatibility

**Tech Tip**  
Use our color coded one-piece fingertight fittings for easy solvent line identification (see pages 27–28)

**Note**  
Complete Tees and Crosses will be supplied with PEEK one-piece fittings.

PEEK unions have the standardized Valco fitting details with smooth transition of the taper into the threaded part.

## PEEK Tees, Crosses and Manifolds High Pressure

- Easy Connection of Solvent Lines
- Low Dead Volume
- Completely Inert and Biocompatible

VICI Jour PEEK tees, crosses and manifolds are usable with most common HPLC solvents (organic and inorganic) without risking corrosion, metal contamination, or protein degradation that might occur with stainless steel. Small diameter thru-holes are used to insure low dead volume connections. Connections to these PEEK tees, crosses and manifolds can be made with any 10-32 fitting that uses polymeric fer-

rules. PEEK tees, crosses and manifolds are available in "body only" versions as well as complete with our PEEK One-Piece Fingertight Fittings ready to connect to any 1/16" OD tubing (PEEK, stainless steel, PTFE, ETFE, etc.).



### PEEK Tees High Pressure

Part No.	Description
JR-1030	Tee, PEEK, HP, 0.25 mm bore, complete with one-piece fittings
JR-1032	Tee, PEEK, HP, 0.50 mm bore, complete with one-piece fittings
JR-1033	Tee, PEEK, HP, 0.75 mm bore, complete with one-piece fittings
JR-2030	Tee, PEEK, HP, 0.25 mm bore, body only
JR-2032	Tee, PEEK, HP, 0.50 mm bore, body only
JR-2033	Tee, PEEK, HP, 0.75 mm bore, body only

### PEEK Crosses High Pressure

Part No.	Description
JR-1040	Cross, PEEK, HP, 0.25 mm bore, complete with one-piece fittings
JR-1042	Cross, PEEK, HP, 0.50 mm bore, complete with one-piece fittings
JR-1043	Cross, PEEK, HP, 0.75 mm bore, complete with one-piece fittings
JR-2040	Cross, PEEK, HP, 0.25 mm bore body only
JR-2042	Cross, PEEK, HP, 0.50 mm bore body only
JR-2043	Cross, PEEK, HP, 0.75 mm bore body only



## PEEK Manifolds High Pressure

Part No.	Description
JR-Z6M1PK	Manifold, PEEK, 6 inlets to 1 outlet, 0.25 mm bore complete with one-piece fittings
JR-B-Z6M1PK	Manifold, PEEK, 6 inlets to 1 outlet, 0.25 mm bore, body only

## Spare Parts


Part No.	Description	Qty/pkg
JR-55021-5	Fitting, PEEK, one-piece black, 10-32	5

## PEEK Biocompatible Mixing Tee

- High-Pressure Gradient Mixing
- Post Column Reactions
- Completely Biocompatible

A unique mixing is the result of the turbulence generated by the introduction angle of the solvents as well as the flow through our PAT filter element – the two solvent streams mix more thoroughly than with a conventional tee.

The biocompatible tee can be used for high-pressure gradient mixing as well as for post column reactions.



The two solvent streams are perfectly mixed in the vortex mixing chamber

## PEEK Biocompatible Mixing Tee

Part No.	Description
JR-9000-0665	Tee, PEEK, mixing, complete with one-piece fittings

## SPECS

## Dimensions

Diameter (body): 25.4 mm  
Thickness (body): 22.9 mm

## SPECS

**Material**  
PEEK, PAT

## Dimensions

Width (body): 25.4 mm  
Thickness: 9.5 mm

## Tolerances

+/- 0.05 mm (.002")

## Threads

10-32  
For OD 1/16" Tubing

## Pressure rating

< 350 bar (< 5000 psi)

## Special info

Metal-free for analytical gradient HPLC

Built-in protein friendly PAT filter (PEEK Alloyed with Teflon)

Porosity 5 µm

Internal volume incl. frit = 4 µL



## PEEK Adapters

### SPECS

**Material**  
PEEK

**Dimensions**  
See chart

#### Pressure rating

M6 - 1/4"-28:  
< 70 bar (< 1000 psi)  
Male and female luer:  
< 3.5 bar (< 50 psi)  
Varies with tubing material  
and ID

#### Special Info

**PEEK Port Adapters**  
These adapters connect  
standard LC fittings to  
standard LC ports

They are available with  
1/4"-28, M6 or 10-32 threads  
(configured for VALCO)

#### PEEK Luer Adapters

Male and female luer adapt-  
ers connect luer fittings, 10-  
32, 1/4"-28 and M6 threaded  
fittings

Luer Adapters not available  
as luer lock

Thru-hole = 0.40 mm

## PEEK Adapters

- Easy Connection of Solvent Lines
- Low Dead Volume
- Completely Inert and Biocompatible



Part No. JR-062



Part No. JR-5595



Part No. JR-5505

## PEEK Adapters

Part No.	Description
JR-062	Adapter, PEEK, 10-32 to 1/4"-28
JR-5595	Pharmacia Adapter M6 to 10-32
JR-5505	Adapter, PEEK, M6 to 10-32

## PEEK Port and Luer Adapters


Part No.	Inlet		Outlet
JR-0611	Waters Female		10-32 Male
JR-0612	1/4"-28 Female		10-32 Male
JR-0622	1/4"-28 Female		M6 Male
JR-0625	1/4"-28 Female		Female Luer
JR-0626	1/4"-28 Female		Male Luer
JR-0613	M6 Female		10-32 Male
JR-0619	M6 Female		1/4"-28 Male
JR-0617	10-32 Female		1/4"-28 Male
JR-0623	10-32 Female		Female Luer

## PEEK Plugs & Caps


- Plugging Ports – Capping Solvent Lines
- All Common Thread Sizes Available



Part No. JR-408




Part No. JR-409



Part No. JR-410



Part No. JR-ZC1PK



Part No. JR-CCPK-1

## PEEK Plugs & Caps

Part No.	Description
JR-408	Plug, PEEK, 10-32
JR-409	Plug, PEEK, 1/4"-28
JR-410	Plug, PEEK, M6
JR-ZC1PK	Cap, PEEK, female 10-32
JR-CCPK-1	Cap, PEEK, female 1/4"-28

## PEEK Starter Kit

- All Parts Biocompatible
- Ideal Kit for All HPLC Users

For a number of interesting LC applications involving proteins, peptides, nucleic acids or other samples of biological origin metal systems may interact with samples or release transition metals that will deactivate columns.

By replacing stainless steel parts, such as tubing, fittings, ferrules, mobile phase filters, in-line filters etc. a biocompatible and metal-free environment is created for samples and mobile phase.

The PEEK Starter Kit is a collection of items, which is a good start for everybody working in the field of biochromatography.

The Kit consists of fittings, ferrules, tubing and solvent filters as well as a union, our Clean-Cut Tubing Cutter and other useful tools like our Tubing Elbows. All together for a very competitive price!



## PEEK Starter Kit

Part No.	Description
JR-35P	PEEK Starter Kit



## Nanovolume® Fittings

### SPECS

**Material**  
PEEK

**Dimensions**  
See illustrations

**Threads**  
Nuts: 6-40  
Ext. Nuts 1/16": 10-32  
Ext. Nuts 1/32": 8-36

**Pressure rating**  
< 350 bar (5000 psi)  
Varies with tubing material and ID

**Note**  
Nanovolume® Unions will be supplied with PEEK one-piece fingertight fittings, ext. nuts and ferrules for 1/32" OD tubing (see chart Nanovolume® Nuts & Ferrules).

The Nanovolume® Reducing Union will additionally be supplied with a PEEK ext. nut and ferrule for 1/16" OD tubing (see chart Nanovolume® Nuts & Ferrules).

\*patents pending


## Nanovolume® Fittings

- 100 and 150 µm Bores Available
- Nuts, Ferrules, Unions, Tees, Crosses and Ys
- For PEEK, PEEKclad-FST™ and Electro Formed Nickel Tubing


VICI Nanovolume® Fittings are designed for Capillary LC used in HTS techniques, Proteomics, Genomics and other newer developments in industrial routine research. Small sample size, reduction of analysis time and higher expectations on instrumentation as well as on accessories are required. VICI Nanovolume®

connectors include our new one-piece fingertight fittings for 1/32 OD tubing with a patented\* collapsible ferrule.

The VICI Nanovolume® Fittings are as easy to handle as our standard analytical fittings and unions. The only difference is their size.



Part No. JR-C-NNFFPK



Part No. JR-C-EN.5FPKB



Part No. JR-ZGF.5PK



Part No. JR-EN1FPKB



Part No. JR-ZGF1PK

## Nanovolume® Nuts & Ferrules

Part No.	Description
JR-C-NNFFPK	Fitting, PEEK, one-piece, grooved ferrule, 1/32"
JR-C-EN.5FPKB	Nut, PEEK, External, 1/32"
JR-ZGF.5PK	Ferrule, PEEK, grooved, 1/32"
JR-EN1FPKB	Nut, PEEK, External, 1/16"
JR-ZGF1PK	Ferrule, PEEK, grooved, 1/16"
JR-C-NPFPK	Plug, PEEK, 1/32"



Part No.  
JR-C-NEU.5XFPK, JR-C-NEU.5FPK

## Nanovolume® Unions

Part No.	Description
JR-C-NEU.5XFPK	Union, PEEK, 100 µm bore, Int. 1/32" – Ext. 1/32" ZDVF
JR-C-NEU.5FPK	Union, PEEK, 150 µm bore, Int. 1/32" – Ext. 1/32" ZDVF



Part No. JR-C-NEU1FPK

## Nanovolume® Reducing Union

Part No.	Description
JR-C-NEU1FPK	Union, PEEK, 150 µm bore, Reducing, Ext. 1/16" – Int. 1/32"



## Nanovolume® Tees

Part No.	Description
JR-C-NTXFPK	Tee, PEEK, 100 µm bore, 1/32", complete
JR-C-NTFPK	Tee, PEEK, 150 µm bore, 1/32", complete



Part No.  
JR-C-NTXFPK, JR-C-NTFPK

## Nanovolume® Crosses

Part No.	Description
JR-C-NXXFPK	Cross, PEEK, 100 µm bore ID, 1/32", complete
JR-C-NXFPK	Cross, PEEK, 150 µm bore, 1/32", complete



Part No.  
JR-C-NXXFPK, JR-C-NXFPK

## Nanovolume® Ys

Part No.	Description
JR-C-NYXFPK	Y, PEEK, 100 µm bore, 1/32", complete
JR-C-NYFPK	Y, PEEK, 150 µm bore, 1/32", complete



Part No.  
JR-C-NYXFPK, JR-C-NYFPK

## SPECS

**Material**  
PEEK

**Dimensions**  
Diameter (body): 21.6 mm  
Thickness: 9.5 mm

**Threads**  
6-40

**Pressure rating**  
< 350 bar (5000 psi)  
Varies with tubing material and ID

**Note**  
Nanovolume® Tees, Crosses and Ys will be supplied with PEEK one-piece fingertight fittings for 1/32" OD tubing (see chart Nanovolume® Nuts & Ferrules).

## PEEK Liners for Fused Silica Tubing Connections

- Available for the Most Common Fused Silica Tubing Sizes
- Two Different Lengths
- Liners with 1µm Screen

## PEEK Liners

Part No.	OD	for Fused Silica OD	Length (mm)	Qty/pkg
JR-C-NL.15L-5	1/32"	125–175 µm	27	5
JR-C-NL.15S-5	1/32"	125–175 µm	10	5
JR-C-NL.20L-5	1/32"	175–225 µm	27	5
JR-C-NL.20S-5	1/32"	175–225 µm	10	5
JR-C-NL.25L-5	1/32"	225–275 µm	27	5
JR-C-NL.25S-5	1/32"	225–275 µm	10	5
JR-C-NL.30L-5	1/32"	275–325 µm	27	5
JR-C-NL.30S-5	1/32"	275–325 µm	10	5
JR-C-NL.35L-5	1/32"	325–375 µm	27	5
JR-C-NL.35S-5	1/32"	325–375 µm	10	5

## PEEK Liners with 1µm Screen

Part No.	OD	for Fused Silica OD	Length (mm)
JR-C-NLS1.15	1/32"	125–175 µm	28
JR-C-NLS1.20	1/32"	175–225 µm	28
JR-C-NLS1.25	1/32"	225–275 µm	28
JR-C-NLS1.30	1/32"	275–325 µm	28
JR-C-NLS1.35	1/32"	325–375 µm	28



## SPECS

**Material**  
Liners: PEEK  
Screens: SS316L

**Dimensions**  
See chart

**Tech Tip**  
Use 10 mm Liners with External nut JR-C-EN.SFPKB  
Use 27/28 mm Liners with One-Piece PEEK Fitting JR-C-NNFPK.



### SPECS

#### Material

Nuts: SS303

Ferrules: SS316

#### Dimensions

OD 1/32", 1/16", 1/8"

#### Threads

See chart

#### Pressure rating

< 500 bar (< 7250 psi)

Varies with tubing material and ID

#### Spares & Tools

We recommend our ValvTool and wrenches for tightening SS Nuts (see page 46).

#### Stainless Steel Tubing

see pages 10-13

## Valco Stainless Steel Nuts and Ferrules

- Excellent Fit to Valco Fitting Details
- Different Lengths Available
- For 1/32", 1/16" and 1/8" OD Tubing

These nuts are for use with all standard Valco internal fittings and most valves. They may be used with fittings from other manufacturers as well. The medium and long types are ideal for metal Valco and Chemin-

ert valves with more than 8 ports. In using nuts with different length more space is created for easier access with a wrench.



Part No. JR-ZN.5-5



Part No. JR-ZF.5S6-5



Part No. JR-ZN1-5



Part No. JR-MZN1-5



Part No. JR-LZN1-5



Part No. ZF1S6



Part No. JR-ZN2-5



Part No. JR-ZF2S6-5

## Valco Stainless Steel Nuts and Ferrules

Part No.	Description	Qty/pkg
JR-ZN.5-5	Nut, SS, 1/32", 6-40	5
JR-ZF.5S6-5	Ferrule, SS, 1/32"	5
JR-ZN1-5	Nut, SS, 1/16", 10-32	5
JR-MZN1-5	Nut, SS, 1/16" Medium, 10-32	5
JR-LZN1-5	Nut, SS, 1/16" Long, 10-32	5
JR-ZF1S6-5	Ferrule, SS, 1/16"	5
JR-ZN2-5	Nut, SS, 1/8", 5/16"-24	5
JR-ZF2S6-5	Ferrule, SS, 1/8"	5

## Stainless Steel Nuts and Ferrules Rheodyne® Type

- Excellent Fit to Rheodyne® Fitting Details
- Different Lengths Available
- For 1/16" OD Tubing

These nuts and ferrules are for use with metal Rheodyne® valves. Three different lengths are available in order to accommodate any possible connecting situation. The standard nut is suitable for easy accessible

connections (couplers, detectors, columns etc.). For very tight and difficult to reach areas we recommend the long and extra long nut.



Part No. JR-S6-5



Part No. JR-S7-5



Part No. JR-S8-5



Part No. JR-S9-5



Part No. JR-60-5

## Stainless Steel Nuts and Ferrules Rheodyne® Type

Part No.	Description	Qty/pkg
JR-S6-5	Nut, SS, short, Rheodyne® type, 10-32	5
JR-S7-5	Nut, SS, long head, Rheodyne® type, 10-32	5
JR-S8-5	Nut, SS, x-long, Rheodyne® type, 10-32	5
JR-S9-5	Ferrule, SS, 1/16"	5
JR-60-5	Ferrule, SS, 1/16", 2-step Rheodyne® type	5

### SPECS

**Material**  
SS316

**Dimensions**  
For OD 1/16" Tubing

**Threads**  
See chart

**Pressure rating**  
< 500 bar (< 7250 psi)  
Varies with tubing material and ID

**Spares & Tools**  
We recommend our ValvTool and wrenches for tightening SS Nuts (see page 46).

Stainless Steel Tubing  
see pages 10–13



### SPECS

#### Material

Body: SS 316  
Nuts: SS303  
Ferrules: SS316

#### Dimensions

See illustrations

#### Pressure rating

< 500 bar (< 7250 psi)  
Varies with tubing material  
and ID


### Stainless Steel Unions

- Easy Connection of Solvent Lines
- Low Dead Volume
- Different Bores Available
- For 1/16" and 1/8" OD SS Tubing

VICI Stainless Steel Unions are made for joining two pieces of 1/16" OD tubing or 1/8" with 1/16" OD tubing. They are precision machined with a focus on the concentricity of the bore as well as on the finish of the surfaces to give an absolute leak-free connection.

VICI SS Unions have internal threads, which we rec-

ommend for most applications. Internal threads make stronger connections and internal fitting details provide the lowest dead volume for high performance instrumentation.



### Stainless Steel Unions

Part No.	Description	Drawing
JR-ZU1XCS6	Union, SS, bore 0.15 mm, 1/16", 10-32, complete	A
JR-ZU1CS6	Union, SS, bore 0.25 mm, 1/16", 10-32, complete	A
JR-ZU1MS6	Union, SS, bore 0.50 mm, 1/16", 10-32, complete	A
JR-ZU1S6	Union, SS, bore 0.75 mm, 1/16", 10-32, complete	A
JR-ZRU21	Union, SS, reducing 0.75 mm bore, 1/8" (5/16"-24) to 1/16" (10-32) complete	B

(All above SS Unions are complete with nuts and ferrules.)

### Spare Parts

Part No.	Description	Qty/pkg
JR-ZN1-5	Nut, SS, 1/16", 10-32, repl. for JR-ZU1x series + JR-ZRU21, 1/16" nut	5
JR-ZF1S6-5	Ferrule, SS, 1/16", repl. for JR-ZU1 series + JR-ZRU21, 1/16" ferrule	5
JR-ZN2-5	Nut, SS, 1/8", 5/16"-24, repl. for JR-ZRU21, 1/8" nut	5
JR-ZF2S6-5	Ferrule, SS, 1/8", repl. for JR-ZRU21, 1/8" ferrule	5



## Stainless Steel Tees & Crosses

- Easy Connection of Solvent Lines
- Low Dead Volume
- Different Bores Available
- For 1/16" OD Tubing

VICI Stainless Steel Tees & Crosses are suitable for connecting 3 or 4 solvent lines. They are precision machined with a focus on the concentricity of the bore as well as on the finish of the surfaces to give an absolute leak free connection.

VICI SS Tees & Crosses have internal threads, which we recommend for most applications. Internal threads make stronger connections and internal fitting details provide the lowest dead volume for high performance instrumentation.



## SPECS

### Material

Bodies: SS316  
Nuts: SS303  
Ferrules: SS316

### Dimensions

Diameter (body): 23.4 mm  
Thickness: 12.7 mm

### Threads

See chart

### Pressure rating

< 350 bar (< 5000 psi)  
Varies with tubing material and ID

## Stainless Steel Tees

Part No.	Description
JR-ZT1C	Tee, SS, 0.25 mm bore, 1/16", 10-32, complete
JR-ZT1M	Tee, SS, 0.50 mm bore, 1/16", 10-32, complete
JR-ZT1	Tee, SS, 0.75 mm bore, 1/16", 10-32, complete

(All above SS Tees are complete with nuts and ferrules.)

## Stainless Steel Crosses

Part No.	Description
JR-ZX1C	Cross, SS, 0.25 mm bore, 1/16", 10-32, complete
JR-ZX1M	Cross, SS, 0.50 mm bore, 1/16", 10-32, complete
JR-ZX1	Cross, SS, 0.75 mm bore, 1/16", 10-32, complete

(All above SS Crosses are complete with nuts and ferrules.)

## Spare Parts

Part No.	Description	Qty/pkg
JR-ZN1-5	Nut, SS, 1/16", 10-32	5
JR-ZF1S6-5	Ferrule, SS, 1/16"	5



### SPECS

#### Material

Plug 1/16": SS316

Plug 1/8": SS303

Caps; SS303

#### Dimensions

See illustration

#### Thread size

Plug & Cap 1/16": 10-32


Plug & Cap 1/8": 5/16"-24

#### Spares & Tools

We recommend our ValvTool and wrenches for tightening SS Plugs & Caps (see page 46).

### Stainless Steel Plugs & Caps

- Plugging Ports – Capping Solvent Lines
- For 1/16" and 1/8" OD Tubing



Part No. JR-ZP1



Part No. JR-ZP2



Part No. JR-ZC1



Part No. JR-ZC2

### Stainless Steel Plugs & Caps

Part No.	Description
JR-ZP1	Plug, SS, 1/16"
JR-ZP2	Plug, SS, 1/8"
JR-ZC1	Cap, SS, female 1/16"
JR-ZC2	Cap, SS, female 1/8"

### Stainless Steel Starter Kit

- Ideal Kit for HPLC Users
- For High Pressure Applications

The Kit consists of fittings, ferrules, tubing and solvent filters as well as a union, our Stainless Steel Tubing Cutter and other useful tools like our Valvtool. All together for a very competitive price!



### Stainless Steel Starter Kit

Part No.	Description
JR-355	Stainless Steel Starter Kit

## ChromBox Kit

– For Start-Up and Standard Use in all HPLC Systems

The ChromBox is an all-purpose analytical spare part kit for all our customers using HPLC. It contains the most often used spare parts made both of high quality SS 316L as well as PEEK and Teflon. All together for a very competitive price!



## ChromBox Kit

Part No.	Description
JR-34N	ChromBox Kit



### ValvTool

- For making connections in 'Hard-to-Reach' areas
- Ideal for Nuts with 1/4" Hex and 5/16" Heads

The new ValvTool is a time-saving device, which will provide easy access to many hard-to-reach areas. The unique design with its slotted wrench allows access to nuts where a loop or a capillary may otherwise

make it difficult. The ValvTool is the ideal for Rheodyne® and Valco valves and fittings as well as for most of the HPLC and PEEK fittings on the market with 1/4" and 5/16" heads.

#### ValvTool

Part No.	Description
JR-800	Tool, ValvTool 1/4" and 5/16" wrench




### File

Standard file

#### File

Part No.	Description
JR-802	Tool, File



### Wrench

- Open wrenches for nuts with 1/4", 5/16", 3/8" or 7/16" heads

#### Wrench

Part No.	Description
JR-804	Tool, Wrench 1/4" and 5/16" open wrench
JR-805	Tool, Wrench 3/8" and 7/16" open wrench



## Tweezer

- Ideal for picking small Fittings etc.

### Tweezers

Part No.	Description
JR-806	Tool, Tweezers



## Ferrule Removal Kit

- Removes Sticking Polymeric Ferrules out of Fitting Details
- Two Tools For Different Ferrule Sizes

### Ferrule Removal Kit

Part No.	Description
JR-FRK1	Ferrule Removal Kit




## Hex Key Set

- Fits Any Socket Head Screw on Any VICI Valve or/Actuator
- Key Sizes: 0.050", 1/16", 5/64", 3/32", 7/64", 1/8", 9/64", 5/32"

### Hex Key Set

Part No.	Description
JR-HKS	Hex Key Set




## Drill Set


- To Enlarge the ID of Fused Silica Adapters
- Useful Tool when Fused Silica Breaks in an Adapter
- 20 Drills Sized from 0.34 to 1 mm (0.0135" to 0.039")

### Drill Set

Part No.	Description
JR-PV	Drill set, incl. Pin Vise







0302

## INDEX

Low Pressure PPS nuts for 1/4"-28 Flat Bottom Fitting Details	49
Color-Sleeves – Fingertight Adapters	51
PEEK Color Coded Flangeless Nuts	52
Polyacetal Color Coded Flangeless Nuts	53
Polypropylene Flanged Fingertight Nuts	54
Flanged Tubing Connector	55
Unions	56
Bulkhead Unions	57
PEEK Tees, Crosses and Manifolds Low Pressure	58
PEEK Biocompatible Mixing Tees – Low Pressure	59
Manual Flow Selection Valves	60
Peristaltic Tubing Adapter Set	61
Low Pressure Polypropylene Adapters	62
Flanged Fittings Kit	63


### Disclaimer:

The maximum holding pressure for any type of connection involving tubing and a ferrule, varies considerably with the tubing material, the ferrule material, the clearance between tubing OD and ferrule ID and the shape of the fitting detail.

## Low Pressure PPS Nuts for 1/4"-28 Flat Bottom Fitting Detail

- Economically Priced Alternative to PEEK
- Excellent Chemical Resistance
- Fits 1/16" or 1/8" OD Tubing
- Slim Design – Fingertight with Adapter

The new VICI Jour PPS 1/4"-28 Fitting series offers excellent chemical and pressure resistance. There is a choice of three different sealing options: collapsible ferrule, inverted ferrule or flanged with washer.



### Flangeless type: Collapsible ferrule\* connection

Very effective sealing on narrow areas with no distortion of the tubing  
No flanging tools required.

Part No.	Description	Qty/pkg
JR-55050-10	Nut, PPS, for collapsible ferrule 1/16", 1/4"-28	10
JR-55051-10	Nut, PPS, for collapsible ferrule 1/8", 1/4"-28	10
JR-CFL-CB1KF	Ferrule, CTFE, collapsible, 1/16"	10
JR-CFL-CB2KF	Ferrule, CTFE, collapsible, 1/8"	10

## SPECS

### Material

Nuts: PPS – PolyPhenylene Sulphide, black

Collapsible Ferrules: CTFE

Inverted Ferrules: ETFE

Washers: Polypropylene

Adapters: Polypropylene

### Chemical Resistance PPS

Chemical Class	Resistance
Acids, organic	A
Acids, inorganic	A/C*
Aldehydes	A
Alcohols	A
Bases	A/B
Esters	A
Halogenated Org.	A/B
Hydrocarbons	A
Ketones	A

\*conc. Halogen Acids

A = suitable

B = Marginal –  
dependant on application

C = Not recommended

Detailed chemical resistance chart see page 106

### Dimensions

For 1/16" and 1/8" OD Tubing

Others: see illustration

### Threads

1/4"-28

### Pressure Rating

1/16" up to 105 bar (1500 psi)

1/8" up to 35 bar (500 psi)

Varies with tubing material and ID

### Tech Tip

For flanging tubing we recommend the VICI Jour Easy Flange Kit (see page 24)

\* Collapsible Ferrule

(see Glossary):

Patent No. 6,575,501



# Polymeric Fittings – Low Pressure

## Low Pressure PPS nuts for 1/4"-28 Flat Bottom Fitting Detail

### SPECS

#### Material

Nuts: PPS – PolyPhenylene Sulphide, black

Collapsible Ferrules: CTFE

Inverted Ferrules: ETFE

Washers: Polypropylene

Adapters: Polypropylene

#### Chemical Resistance PPS

Chemical Class	Resistance
Acids, organic	A
Acids, inorganic	A/C*
Aldehydes	A
Alcohols	A
Bases	A/B
Esters	A
Halogenated Org.	A/B
Hydrocarbons	A
Ketones	A

\*conc. Halogen Acids

A = suitable

B = Marginal –

dependant on application

C = Not recommended

Detailed chemical resistance chart see page 106

#### Dimensions

For 1/16" and 1/8" OD Tubing

Others: see illustration

#### Threads

1/4"-28

#### Pressure Rating

1/16" up to 105 bar (1500 psi)


1/8" up to 35 bar (500 psi)

Flanged with Washer: 30bar (435psi)

Varies with tubing material and ID

#### Tech Tip


For flanging tubing we recommend the VICI Jour Easy Flange Kit (see page 24)



#### Flangeless type: Inverted ferrule connection

No flanging tools required.



Part No.	Description	Qty/pkg
JR-55070-10	Nut, PPS, flangeless 1/16", 1/4"-28	10
JR-55071-10	Nut, PPS, flangeless 1/8", 1/4"-28	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10



#### Flanged type: Flange/washer connection

Solvent comes only into contact with the tubing.

Part No.	Description	Qty/pkg
JR-55060-10	Nut, PPS, flanged 1/16", 1/4"-28	10
JR-55061-10	Nut, PPS, flanged 1/8", 1/4"-28	10
JR-201926-10	Washer, PP, flanged, 1/16"	10
JR-201928-10	Washer, PP, flanged, 1/8"	10



### SPECS

#### Material

Adapters: Polypropylene

#### Color-Sleeve – fingertight adapter

This sleeve is useful to color code different solvent lines and it turns the new VICI Jour low pressure nut into a fingertight fitting in an instant. The Color-Sleeve is available in 12 colors.

Part No.	Description	Qty/pkg
JR-55080-10	Adapter, PP, fingertight sleeve, natural	10
JR-55081-10	Adapter, PP, fingertight sleeve, white	10
JR-55082-10	Adapter, PP, fingertight sleeve, dark grey	10
JR-55083-10	Adapter, PP, fingertight sleeve, black	10
JR-55084-10	Adapter, PP, fingertight sleeve, lavender	10
JR-55085-10	Adapter, PP, fingertight sleeve, red	10
JR-55086-10	Adapter, PP, fingertight sleeve, yellow	10
JR-55087-10	Adapter, PP, fingertight sleeve, orange	10
JR-55088-10	Adapter, PP, fingertight sleeve, brown	10
JR-55089-10	Adapter, PP, fingertight sleeve, green	10
JR-55090-10	Adapter, PP, fingertight sleeve, blue	10
JR-55091-10	Adapter, PP, fingertight sleeve, purple	10
JR-5508X-24	Adapter, PP, fingertight sleeve, assorted 2 pcs. per color	24
JR-5508X-12	Adapter, PP, fingertight sleeve, assorted 2 pcs. each blue, red, green, yellow, black, white	12



### SPECS

#### Material

Nuts: PEEK  
Ferrules: ETFE

#### Dimensions

For 1/16" and 1/8" OD tubing  
Others: see illustration

#### Threads

1/4"-28

#### Pressure rating

1/16" up to 105 bar (1500 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material  
and ID

## PEEK Color Coded Flangeless Nuts

- Flange-Free – Compatible with All 1/4"-28 Flat Bottom Ports
- Excellent Chemical Resistance
- Fits 1/16" or 1/8" OD Tubing
- True Fingertight – No Tools Required

PEEK Flangeless Nuts in combination with the inverted ferrule do not require flanging tools. A hand tightening procedure makes leak free seals. The inert ETFE ferrule is the only part, besides the tubing, which is in

contact with the fluid – the nut is not in direct contact with the fluid stream (see illustration) Six colors are available.




### PEEK Color Coded Flangeless Nuts 1/16"

Part No.	Description	Color	Qty/pkg
JR-20124-10	Nut, PEEK, flangeless	Natural	10
JR-20125-10	Nut, PEEK, flangeless	Yellow	10
JR-20126-10	Nut, PEEK, flangeless	Red	10
JR-20127-10	Nut, PEEK, flangeless	Blue	10
JR-20128-10	Nut, PEEK, flangeless	Green	10
JR-20129-10	Nut, PEEK, flangeless	Black	10
JR-041-10	Ferrule, ETFE, 1/16"		10



Part No. JR-20112-10

Flange-free 1-16



### PEEK Color Coded Flangeless Nuts 1/8"

Part No.	Description	Color	Qty/pkg
JR-20112-10	Nut, PEEK, flangeless	Green	10
JR-20113-10	Nut, PEEK, flangeless	Blue	10
JR-20114-10	Nut, PEEK, flangeless	Yellow	10
JR-20115-10	Nut, PEEK, flangeless	Red	10
JR-20116-10	Nut, PEEK, flangeless	Natural	10
JR-20117-10	Nut, PEEK, flangeless	Black	10
JR-051-10	Ferrule, ETFE, 1/8"		10

## Polyacetal Color Coded Flangeless Nuts

- Value for Money Alternative to PEEK
- Flange-Free – Compatible with All 1/4"-28 Flat Bottom Ports
- Fits 1/16" or 1/8" OD Tubing
- True Fingertight – No Tools Required

Polyacetal Flangeless Nuts in combination with the inverted ferrule do not require flanging tools. A hand tightening procedure makes leak free seals. The inert ETFE ferrule is the only part, besides the tubing,

which is in contact with the fluid – the nut is not in direct contact with the fluid stream (see illustration) Six colors are available.




### Polyacetal Color Coded Flangeless fingertight Nuts 1/16"

Part No.	Description	Color	Qty/pkg
JR-040GR-10	Nut, Polyacetal, flangeless 1/16"	green	10
JR-040BL-10	Nut, Polyacetal, flangeless 1/16"	blue	10
JR-040YL-10	Nut, Polyacetal, flangeless 1/16"	yellow	10
JR-040RE-10	Nut, Polyacetal, flangeless 1/16"	red	10
JR-040WH-10	Nut, Polyacetal, flangeless 1/16"	white	10
JR-040BK-10	Nut, Polyacetal, flangeless 1/16"	black	10
JR-040WL-10	Nut, Polyacetal, flangeless 1/16", 1" long	white	10
JR-041-10	Ferrule, ETFE, 1/16"		10



Flange-free 1-16



### Polyacetal Color Coded Flangeless fingertight Nuts 1/8"

Part No.	Description	Color	Qty/pkg
JR-050GR-10	Nut, Polyacetal, flangeless 1/8"	green	10
JR-050BL-10	Nut, Polyacetal, flangeless 1/8"	blue	10
JR-050YL-10	Nut, Polyacetal, flangeless 1/8"	yellow	10
JR-050RE-10	Nut, Polyacetal, flangeless 1/8"	red	10
JR-050WH-10	Nut, Polyacetal, flangeless 1/8"	white	10
JR-050BK-10	Nut, Polyacetal, flangeless 1/8"	black	10
JR-050WL-10	Nut, Polyacetal, flangeless 1/8", 1" long white		10
JR-051-10	Ferrule, ETFE, 1/8"		10



## SPECS

### Material

Nuts: Polyacetal  
Ferrules: ETFE

### Dimensions

For 1/16" and 1/8" OD Tubing  
Others: see illustration

### Threads

1/4"-28

### Pressure rating

1/16" up to 70 bar (1000 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material and ID

### Note

Flangeless Fittings for 1/4" OD tubing available: see spare parts for Unions on page 56



## Polymeric Fittings – Low Pressure

### Polypropylene Flanged Fingertight Nuts

#### SPECS

##### Material

Nuts: Fiber Glass reinforced Polypropylene  
Washers: Polypropylene

##### Dimensions

For OD 1/16", 1/8" Tubing  
Others: see illustration

##### Threads

1/4"-28

##### Pressure rating

< 30 bar (< 435 psi)  
Varies with tubing material and ID

##### Tech Tip

We recommend our Easy Flange Kit for flanging tubing (see page 24).

### Polypropylene Flanged Fingertight Nuts


- For Flanged Connections
- Compatible to All 1/4"-28 Flat Bottom Ports
- For 1/16" or 1/8" OD Tubing
- Fiber Glass Reinforced PP with PP Washers

A number of chromatography systems still require, or prefer flanged connections. Consequently, we offer a series of such nuts. By flanging the tubing, the washer never comes in contact with the mobile phase.



#### Flanged Fingertight Nuts 1/16"

Part No.	Description	Color	Qty/pkg
JR-201094BK-10	Nut, PP, flanged	black	10
JR-201094WH-10	Nut, PP, flanged	white	10
JR-202118-10	Nut, PP, flanged, x-long 1"	black	10
JR-202120-10	Nut, PP, flanged, x-long 1"	white	10
JR-202122-10	Nut, PP, flanged	blue	10
JR-202123-10	Nut, PP, flanged	green	10
JR-202124-10	Nut, PP, flanged	yellow	10
JR-202125-10	Nut, PP, flanged	red	10
JR-201926-10	Washer, PP, flanged 1/16"		10



#### Flanged Fingertight Nuts 1/8"


Part No.	Description	Color	Qty/pkg
JR-201580BK-10	Nut, PP, flanged	black	10
JR-201580WH-10	Nut, PP, flanged	white	10
JR-202119-10	Nut, PP, flanged, x-long 1"	white	10
JR-202121-10	Nut, PP, flanged, x-long 1"	black	10
JR-202142-10	Nut, PP, flanged	blue	10
JR-202143-10	Nut, PP, flanged	green	10
JR-202144-10	Nut, PP, flanged	yellow	10
JR-202145-10	Nut, PP, flanged	red	10
JR-201928-10	Washer, PP, flanged 1/8"		10



## Flanged Tubing Connector

- For Easy Connection of Flanged Tubing
- For 1/16" and 1/8" Tubing
- Zero Dead Volume

The VICI Flanged Tubing Connector fits all of our flanged fittings. With its female 1/4"-28 thread it is designed for use on tubing with a flanged end. We recommend them instead of a union to make a Zero Dead Volume connection.



## SPECS

**Material**  
PEEK

**Threads**  
1/4"-28

**Pressure rating**  
< 30 bar (< 435 psi)  
Varies with tubing material and ID

**Tech Tip**  
We recommend our Easy Flange Kit for flanging tubing (see page 24).

## Flanged Tubing Union

Part No.	Description
JR-CEN1PK	Flanged Tubing Connector, 1/16"
JR-CEN2PK	Flanged Tubing Connector, 1/8"



## Unions

### SPECS

#### Materials

Unions: PEEK, Polyamid  
Nuts: PEEK, Polyacetal  
Ferrules: ETFE, CTFE

#### Dimensions

See illustrations

#### Threads

Unions: 1/4"-28  
Reducing Union: 1/2"-20 to 1/4"-28

#### Pressure Rating of PEEK Unions

1/16" up to 105 bar (1500 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material and ID

#### Note

Complete versions with 0.50 and 0.75 bore are supplied with PPS flangeless nuts and ETFE ferrules for 1/16" OD tubing (see page 50).

Complete versions with 1.30 mm bore are supplied with PEEK flangeless nuts and ETFE ferrules for 1/8" OD tubing (see page 50).

The Reducing Union is supplied with Polyacetal flangeless nuts and ferrules (ETFE/CTFE) for 1/8" and 1/4" OD tubing

#### Tech Tip


VICI Unions have 1/4"-28 flat bottom ports. All of our three different sealing options (collapsible ferrule, inverted ferrule or flanged with washer) can be used (see pages 49-54).

To avoid diffusion of gases in your solvent line use our No-Ox tubing on page 18.

## Unions

- For Easy Connection of Low Pressure Solvent Lines
- Low Dead Volume
- Completely Inert and Biocompatible

These unions allow the connection of two pieces of 1/16", 1/8" or 1/8" to 1/4" OD tubing quickly and reliably. The unions are inert and biocompatible.



## Unions

Part No.	Description	Qty/pkg
JR-060-5	Union, Polyamid, through bore LP body only, 1/4"-28, for flanged tubing only	5
JR-1064	Union, PEEK, LP 0.5 mm bore, complete	1
JR-1068	Union, PEEK, LP 0.75 mm bore, complete	1
JR-1065	Union, PEEK, LP 1.3 mm bore, complete	1
JR-064	Union, PEEK, LP 0.5 mm bore, body only	1
JR-068	Union, PEEK, LP 0.75 mm bore, body only	1
JR-065	Union, PEEK, LP 1.3 mm bore, body only	1



## Reducing Union

Part No.	Description
JR-CU4LPK	Reducing Union, PEEK, LP 1.5 mm bore 1/4" to 1/8", complete

## Spare Parts

Part No.	Description	Qty/pkg
JR-20129-10	Nut, Peek, flangeless 1/16", black	10
JR-20117-10	Nut, Peek, flangeless 1/8", black	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10
R-050WH-10	Nut, Polyacetal, flangeless 1/8", white	10
JR-CFL-4D	Nut, Polyacetal, flangeless 1/4"	1
JR-CFL-CB4KF-S	Ferrule, CTFE, 1/4"	1



## Bulkhead Unions

- Panel mounting
- For Easy Connection of Low Pressure Solvent Lines
- Low Dead Volume
- Completely Inert and Biocompatible



### Bulkhead Unions

Part No.	Description
JR-CBUMPK	Bulkhead Union, PEEK, LP 0.75 mm bore, complete
JR-CBULPK	Bulkhead Union, PEEK, LP 1.50 mm bore, complete

### Spare Parts

Part No.	Description	Qty/pkg
JR-55070-10	Nut, PPS, flangeless 1/16"	10
JR-55071-10	Nut, PPS, flangeless 1/8"	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10
JR-55086-10	Adapter, PP, fingertight sleeve, yellow	10
JR-55090-10	Adapter, PP, fingertight sleeve, blue	10

Adapters available in different colors (see page 51)

## SPECS

### Material

Unions: PEEK  
Nuts: PPS  
Ferrules: ETFE  
Sleeves: Polypropylene

### Dimensions

See illustration

### Threads

Unions: 1/4"-28

### Pressure Rating

1/16" up to 105 bar (1500 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material and ID

### Note

0.75 mm bore version is supplied with PPS flangeless fittings and ETFE ferrules for 1/16" OD tubing (see page 50).

1.50 mm bore version is supplied with PPS flangeless fittings and ETFE ferrules for 1/8" OD tubing (see page 50).

### Tech Tip

Bulkhead Unions require a min. 10 mm panel hole.



# Polymeric Fittings – Low Pressure

## PEEK Tees, Crosses and Manifolds Low Pressure

### SPECS

#### Material

Body: PEEK  
Fittings: PPS  
Ferrules: ETFE  
Sleeves: Polypropylene

#### Dimensions

Tees and Crosses  
Diameter (body): 25.0 mm  
Thickness: 12.0 mm

Manifolds 5 ports  
Diameter (body): 22.2 mm  
Thickness (body): 17.1 mm

Manifolds 9 ports  
Diameter (body): 22.2 mm  
Thickness (body): 28.6 mm

#### Threads

1/4"-28

#### Pressure Rating

1/16" up to 105 bar (1500 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material  
and ID

#### Note

Complete versions with 0.50, 0.75 and 0.80 mm bore are supplied with PPS flangeless nuts and ETFE ferrules for 1/16" OD tubing (see page 50).

Complete versions with 1.50 mm bore are supplied with PPS flangeless nuts and ETFE ferrules for 1/8" OD tubing (see page 50).

#### Tech Tip

VICI PEEK Tees, Crosses & Manifolds have 1/4"-28 flat bottom ports. All of our three different sealing options (collapsible ferrule, inverted ferrule or flanged with washer) can be used (see pages 49–54).

To avoid diffusion of gases in your solvent line use our No-Ox tubing on page 18

## PEEK Tees, Crosses and Manifolds Low Pressure


- Biocompatible
- Different Bore Sizes Available
- For 1/16" and 1/8" OD Tubing
- 1/4"-28 Ports

VICI Tees, Crosses and Manifolds for low pressure applications are machined from PEEK. The bodies have 1/4"-28 female ports which accept any flanged or flangeless fitting for 1/16" and 1/8" OD tubing. Choose

0.80 mm and smaller thru-holes for applications where low flow rate or low dead volume is required. Tees, Crosses and Manifolds with 1.50 mm thru-holes are preferred for high-flow rate applications.


### PEEK Tees

Part No.	Description
JR-6029	Tee, PEEK, LP bore 0.50 mm, complete
JR-6030	Tee, PEEK, LP bore 0.8 mm, complete
JR-6031	Tee, PEEK, LP bore 1.5 mm, complete
JR-6034	Tee, PEEK, LP bore 0.50 mm, body only
JR-6032	Tee, PEEK, LP bore 0.8 mm, body only
JR-6033	Tee, PEEK, LP bore 1.5 mm, body only



### PEEK Crosses

Part No.	Description
JR-6039	Cross, PEEK, LP bore 0.50 mm, complete
JR-6040	Cross, PEEK, LP bore 0.8 mm, complete
JR-6041	Cross, PEEK, LP bore 1.5 mm, complete
JR-6044	Cross, PEEK, LP bore 0.50 mm, body only
JR-6042	Cross, PEEK, LP bore 0.8 mm, body only
JR-6043	Cross, PEEK, LP bore 1.5 mm, body only



### PEEK Manifolds

Part No.	Description
JR-C5M1PK	Manifold, PEEK, 5 ports LP 0.75 mm, complete
JR-C5M2PK	Manifold, PEEK, 5 ports LP 1.50 mm, complete
JR-C9M1PK	Manifold, PEEK, 9 ports LP 0.75 mm, complete
JR-C9M2PK	Manifold, PEEK, 9 ports LP 1.50 mm, complete



### Spare Parts

Part No.	Description	Qty/pkg
JR-55070-10	Nut, PPS, flangeless 1/16", 1/4"-28	10
JR-55071-10	Nut, PPS, flangeless 1/8", 1/4"-28	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10

Color-Sleeves – fingertight adapters: see page 51

## PEEK Biocompatible Mixing Tees – Low Pressure

- For Low Pressure Mixing
- Completely Biocompatible
- Mountable

A unique mixing is the result of the turbulence generated by the introduction angle of the solvents. The streams mix more thoroughly than with a conventional tee.



## PEEK Biocompatible Mixing Tees – Low Pressure

Part No.	Description
JR-CM1XPK	Mixing Tee, PEEK, 0.75 mm bore, complete
JR-CM2XPK	Mixing Tee, PEEK, 1.50 mm bore, complete

## Spare Parts

Part No.	Description	Qty/pkg
JR-55070-10	Nut, PPS, flangeless 1/16", 1/4"-28	10
JR-55071-10	Nut, PPS, flangeless 1/8", 1/4"-28	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10

Color-Sleeves – fingertight adapters: see page 51

## SPECS

## Material

Body: PEEK  
Fittings: PPS  
Ferrules: ETFE

## Dimensions

Length (body): 44.5 mm  
Thickness: 22.2 mm

## Threads

1/4"-28

## Pressure rating

1/16" up to 105 bar (1500 psi)  
1/8" up to 35 bar (500 psi)  
Varies with tubing material and ID

## Note

0.75 bore version is supplied with PPS flangeless fittings and ferrules for 1/16" OD tubing (see page 50).

1.50 mm bore version is supplied with PPS flangeless fittings and ferrules for 1/8" OD tubing (see page 50).



### SPECS

#### Material

Wetted parts: PTFE, CTFE  
Body: Polyacetal  
Fittings: PPS  
Ferrule: ETFE  
Rotor bar: SS

#### Dimensions

Bore: 0.8 mm  
Others: see illustration

#### Threads

1/4"-28

#### Pressure Rating

< 35 bar (< 500 psi)  
Varies with tubing material and ID

#### Note

Complete versions are supplied with PPS flangeless nuts and ETFE ferrules for 1/16" OD tubing.

Valves can also be used with 1/8" OD tubing. 1/8" nuts and ferrules: see spare parts.

#### Tech Tip


VICI Manual Flow Selection Valves have 1/4"-28 flat bottom ports. All of our three different sealing options (collapsible ferrule, inverted ferrule or flanged with washer) can be used (see pages 49-54).

To avoid diffusion of gases in your solvent line use our No-Ox tubing on page 18.

## Manual Flow Selection Valves

- Several Flow Patterns Available
- Biocompatible and Inert
- Pressure Resistant up to 35 bar (500 psi)
- Low Dead Volume

VICI Jour Manual Valves are the ideal tools for liquid flow selection. The valves are mountable and are pre-fitted with the new VICI Jour PPS 1/4"-28 Fittings and ETFE Ferrules for 1/16" OD tubing.



### Manual Flow Selection Valves

Part No.	Description	
JR-MSV2S	Valve, 2-port, manual, on/off, complete	A
JR-MSV4L	Valve, 4-port, manual, L-flow, complete	B
JR-MSV4T	Valve, 4-port, manual, T-flow, complete	B
JR-MSV4P	Valve, 4-port, manual, 2 pairs, complete	B
JR-MSV5D	Valve, 5-port, manual, distribution 1 inlet to 4 outlets, complete	C

### Spare Parts

Part No.	Description	Qty/pkg
JR-55070-10	Nut, PPS, flangeless 1/16", 1/4"-28	10
JR-55071-10	Nut, PPS, flangeless 1/8", 1/4"-28	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10

Color-Sleeves – fingertight adapters: see page 51



## Peristaltic Tubing Adapter Set

- For Connecting Soft Tubing to PTFE & PTFE Like Tubing
- Works With Peristaltic Tubing up to 4.50 mm OD
- For 1/16" and 1/8" PTFE & PTFE Like Tubing
- Biocompatible
- Low Dead Volume

The "One-Fits-All" set consists of everything to connect 1/16" and 1/8" OD PTFE & PTFE like tubing to soft tubing like silicone up to 4.50 mm OD. The adapter

comes with three different bushings (for the different soft tubing ODs) and two Fittings/Ferrules (one each for 1/16" and 1/8" OD tubing).



## Peristaltic Tubing Adapter Set

Part No.	Description
JR-PTAS	Adapter, PEEK, for Peristaltic to PTFE tubing connections, complete

## SPECS

### Material

Body and Bushings: PEEK  
Fittings: PPS  
Ferrules: ETFE

### Dimensions

Bore: 0.50 mm  
Length: 28 mm  
(without Fitting)  
Diameter: 12 mm

### Threads

Barbed to 1/4"-28

### Tech Tip

The set is equipped with fittings and ferrules for inverted ferrule connections. The adapter also works with our other sealing options (see pages 49-54)







### SPECS

**Material**  
Polypropylene

### Low Pressure Polypropylene Adapters

#### Low Pressures Polypropylene Adapters

Part No.	Description		Qty/pkg
JR-070211-10	Adapter, PP, male Luer to 1/16" barbed		10
JR-070212-10	Adapter, PP, male Luer to 3/32" barbed		10
JR-070213-10	Adapter, PP, male Luer to 1/8" barbed		10
JR-070111-10	Adapter, PP, female Luer to 1/16" barbed		10
JR-070112-10	Adapter, PP, female Luer to 3/32" barbed		10
JR-070113-10	Adapter, PP, female Luer to 1/8" barbed		10
JR-070101-10	Adapter, PP, female to female Luer		10
JR-070124-10	Adapter, PP, female Luer to 10-32 male		10
JR-070123-10	Adapter, PP, female Luer to 1/4"-28 male		10
JR-071123-10	Adapter, PP, 1/4"-28 male to 1/16" barbed		10
JR-071323-10	Adapter, PP, 1/4"-28 male to 1/8" barbed		10
JR-072350-10	Adapter, PP, plug 1/4"-28 male		10
JR-071124-10	Adapter, PP, 10-32 male to 1/16" barbed		10
JR-071324-10	Adapter, PP, 10-32 male to 1/8" barbed		10
JR-072450-10	Adapter, PP, plug 10-32 male		10
JR-070224-10	Adapter, PP, male Luer to 10-32 male		10
JR-070225-10	Adapter, PP, male Luer to 1/4"-28 female		10

## Flanged Fittings Kit

- A complete Kit of Flanged Fittings and accessories
- All Parts Biocompatible
- Excellent chemical resistance

The Flanged Fittings Kit contains all the materials that are required to make flanged connections for 1/16" and 1/8" OD tubing: flanged fittings, washers, unions,

4 different tubing sizes, our Clean-Cut Tubing Cutter, and a pair of tweezers.



## SPECS


### Tech Tip

Use our Easy Flange Kit for flanging your Tubing (see Page 24)

## Flanged Fittings Kit

Part No.	Description
JR-201529	Flanged Fitting Kit





C0577

## INDEX

Last Drop Mobile Phase Filter	65
Last Drop Filter/Spargers	66
Economy Last Drop Mobile Phase Filters	67
No-Met Biocompatible Mobile Phase Filter	68
SS Mobile Phase Filters & Helium Spargers	69
In-Line Low Pressure Filter Cartridges	70
PEEK In-Line Filter Kit – High Pressure	71
PEEK Pre-Column Filter Kit – High Pressure	72
Sure-Guard	73
SS High-Pressure In-Line Filter	74




## Last Drop Mobile Phase Filter

- No Loss of Mobile Phase
- Biocompatible PTFE Frits or SS Frits
- Three Different Porosities
- Two Connector Types

The Last Drop Mobile Phase Filter utilizes a flat filter element which sits parallel to the bottom of the reservoir. This design allows the filter to draw all but the last 2% of the mobile phase from the reservoir without drawing air into the system. Conventional cylindrical mobile phase filters begin to draw air into the

system when less than 10% of the solvent remains in the reservoir. We recommend the metal free PTFE version for sensitive biochromatography applications where metal surfaces may corrode and contaminate the solvent with ions.



### SPECS

#### Materials

Body: PTFE  
Frits: PTFE, SS316 (see chart)  
Tripod: PEEK  
Fitting: PEEK  
Ferrule: ETFE

#### Dimensions

Fitting: 1/4"-28 for 1/8" OD tubing  
Tripod: for 1.5, 2.2 and 3.5 mm ID tubing.  
Others: see illustration

#### Spares & Tools

For spare parts see chart on page 66

### Last Drop Mobile Phase Filters

Part No.	Description
JR-9000-0520	Filter, PTFE, Last Drop 2.5 $\mu$ m, stepped tubing connector
JR-9000-0520F	Filter, PTFE, Last Drop 2.5 $\mu$ m, fitting connector
JR-9000-0521	Filter, PTFE, Last Drop 5 $\mu$ m, stepped tubing connector
JR-9000-0521F	Filter, PTFE, Last Drop 5 $\mu$ m, fitting connector
JR-9000-0522	Filter, PTFE, Last Drop 10 $\mu$ m, stepped tubing connector
JR-9000-0522F	Filter, PTFE, Last Drop 10 $\mu$ m, fitting connector
JR-9000-0530	Filter, SS, Last Drop 2 $\mu$ m, stepped tubing connector
JR-9000-0530F	Filter, SS, Last Drop 2 $\mu$ m, fitting connector
JR-9000-0531	Filter, SS, Last Drop 5 $\mu$ m, stepped tubing connector
JR-9000-0531F	Filter, SS, Last Drop 5 $\mu$ m, fitting connector
JR-9000-0532	Filter, SS, Last Drop 10 $\mu$ m, stepped tubing connector
JR-9000-0532F	Filter, SS, Last Drop 10 $\mu$ m, fitting connector

## Last Drop Filter/Spargers

### SPECS

#### Materials

Body: PTFE  
Frits: PTFE, SS316 (see chart)  
Tripod: PEEK  
Fitting: PEEK  
Ferrule: ETFE

#### Dimensions

Fitting: 1/4"-28 for 1/8" OD tubing  
Tripod: for 1.5, 2.2 and 3.5 mm ID tubing.  
Others: see illustration

#### Tech Tip


We recommend our No-Ox tubing to prevent "regassing" of Helium degassed solvents (see page 18)

## Last Drop Filter/Spargers

- Parallel Filtering and Sparging
- Biocompatible PTFE Frits or SS Frits
- Three Different Porosities

This Filter/Sparger filter combines filtration and sparging in one single unit. The PTFE housing contains a mobile phase filter, either a stainless steel or a

PTFE filter element. The stepped PEEK connector is for the solvent line, the fitting connection for the Helium line.



## Last Drop Filter/Spargers

Part No.	Description
JR-9000-0602	Filter/Sparger, PTFE, Last Drop 2.5 µm filter, 10 µm sparger
JR-9000-0603	Filter/Sparger, PTFE, Last Drop 5 µm filter, 10 µm sparger
JR-9000-0604	Filter/Sparger, PTFE, Last Drop 10 µm filter, 10 µm sparger
JR-9000-0640	Filter/Sparger, SS, Last Drop 2 µm filter, 10 µm sparger
JR-9000-0641	Filter/Sparger, SS, Last Drop 5 µm filter, 10 µm sparger
JR-9000-0642	Filter/Sparger, SS, Last Drop 10 µm filter, 10 µm sparger

## Spare Parts

Part No.	Description	Qty/pkg
JR-20116-10	Nut, PEEK, flangeless, 1/8", natural	10
JR-051-10	Ferrule, ETFE, 1/8"	10
JR-8000-0485	Tripod Adapter, PEEK, universal	1





## Economy Last Drop Mobile Phase Filters

- Very Competitively Priced
- Biocompatible PTFE Frits or SS Frits
- Three Different Porosities
- Two Diameters

We can offer Economy Mobile Phase Filters in 2 sizes, depending upon the dimension of the bottleneck. The filter body is made of PTFE which has been specially selected to be resistant against virtually all com-

mon mobile phases. We recommend the 2  $\mu\text{m}$  porosity for flow rates less than 10 mL/min. For higher flow rates use the 10  $\mu\text{m}$  filters.



### Economy Last Drop Mobile Phase Filter

Part No.	Description
JR-4676-2.5TF	Filter, PTFE, Economy mobile phase 2,5 $\mu\text{m}$ , OD 19 mm
JR-4676-5TF	Filter, PTFE, Economy mobile phase 5 $\mu\text{m}$ , OD 19 mm
JR-4676-10TF	Filter, PTFE, Economy mobile phase 10 $\mu\text{m}$ , OD 19 mm
JR-4677-2.5TF	Filter, PTFE, Economy mobile phase 2,5 $\mu\text{m}$ , OD 22 mm
JR-4677-5TF	Filter, PTFE, Economy mobile phase 5 $\mu\text{m}$ , OD 22 mm
JR-4677-10TF	Filter, PTFE, Economy mobile phase 10 $\mu\text{m}$ , OD 22 mm
JR-4676-2	Filter, SS, Economy mobile phase 2 $\mu\text{m}$ , OD 19 mm
JR-4676-10	Filter, SS, Economy mobile phase 10 $\mu\text{m}$ , OD 19 mm
JR-4677-2	Filter, SS, Economy mobile phase 2 $\mu\text{m}$ , OD 22 mm
JR-4677-10	Filter, SS, Economy mobile phase 10 $\mu\text{m}$ , OD 22 mm

### SPECS

#### Material

Body: PTFE

Frits: PTFE/SS316 (see chart)

#### Dimension

For 1/8" OD tubing

Others: see illustration

#### Tech Tip

Easy to replace – easy to clean  
– use ultrasound or replace with new

We recommend our tubing JR-T-6800-M3 (1/8" OD x 1.60 mm ID) for connecting the filter (see page 15)

### SPECS

#### Material

Body: Polyethylene/PTFE

(see chart)

Adapter: PEEK

Fitting: PEEK

Ferrule: ETFE

#### Dimension

See illustration

#### Tech Tip

The Polyethylene filter is hydrophobic, therefore it may initially require some priming with methanol or acetonitrile.

#### Spares & Tools

We recommend our tubing

JR-T 6800-M3 (1/8" OD x

1.60 mm ID) for connecting

the filter (see page 15)

## No-Met Biocompatible Mobile Phase Filter


- Very Competitively Priced
- Inert and Biocompatible
- Replacement Filters
- Also Recommended for IC

Stainless steel is not an acceptable material for a growing number of applications involving the separation of labile macro biomolecules. It is of great importance that all of the biological activity of the samples is retained and that the biological specificity remains unchanged.

High salt buffer concentrations can corrode stainless steel and the metal ions released from these filters

may react with the biomolecules of interest or even contaminate the eluents in Ion chromatography.

The No-met filters are designed from inert polymeric components, which effectively eliminate metal contamination of the fluid path. The economy version can be easily slipped over 1/8" OD tubing, no fitting required.



### No-Met Mobile Phase Filter

Part No.	Description
JR-32171	Filter, PTFE, mobile phase No-Met 5 µm, 1/8"
JR-32172	Filter, PTFE, mobile phase replacement No-Met 5 µm
JR-32178	Filter, PE, mobile phase No-Met < 20 µm, 1/8"
JR-32179	Filter, PE, mobile phase replacement No-Met < 20 µm
JR-32174	Filter, PE, Economy No-Met 5 µm, 1/8"




## SS Mobile Phase Filters & Helium Spargers

- Ideal for Helium Sparging
- Three Different Porosities
- Fitting and Pipe Connectors

VICI Jour Mobile Phase Filters protect HPLC systems from small particles in the mobile phase. These filters are made from SS316 with PEEK or PTFE connectors and are suitable for most solvents.

The complete line has versions for both analytical and preparative applications.

VICI Jour Helium Spargers (2  $\mu$ m versions) offer an inexpensive way to prepare and maintain mobile phases free of dissolved gases. Connected to a regulated supply of Helium gas (0–400 mL/min) they effectively remove dissolved oxygen, nitrogen and other atmospheric gases from the mobile phase.



### Mobile Phase Filters & Helium Spargers

Part No.	Drawing	Porosity	Fitting/Pipe OD	Suggested Max. Flow Rate
JR-3675-2	A	2 $\mu$ m	1/8" Pipe OD	10 mL/min
JR-367008-2	B	2 $\mu$ m	1/8" Fitting ID	8 mL/min
JR-367008-10	B	10 $\mu$ m	1/8" Fitting ID	20 mL/min
JR-367008-20	B	20 $\mu$ m	1/8" Fitting ID	20 mL/min
JR-367016-2	C	2 $\mu$ m	1/16" Fitting ID	8 mL/min
JR-367016-10	C	10 $\mu$ m	1/16" Fitting ID	20 mL/min
JR-367016-20	C	20 $\mu$ m	1/16" Fitting ID	20 mL/min
JR-3678-2	D	2 $\mu$ m	1/8" Pipe OD	50 mL/min
JR-3678-25	D	25 $\mu$ m	1/8" Pipe OD	100 mL/min

### SPECS

#### Material

Body: SS316  
Pipe: SS316  
Pipe Adapter: PTFE  
Fitting Adapter: PEEK  
Fittings: PEEK  
Ferrules: ETFE

#### Dimension

See illustration

#### Tech Tip

We recommend our No-Ox tubing to prevent "regassing" of Helium degassed solvents (see page 18)



### SPECS

#### Material

Cartridge: PTFE/CTFE  
Filter Screen: SS316

#### Dimensions

For All 1/4"-28 fitting details  
OD: 5.20 mm  
Thickness: 2.03 mm  
Bore: 0.8 mm  
Filter Surface Diameter 2 mm

#### Tolerances

+/- 0.05 mm (.002")

#### Tech Tip

We recommend using these Filter Cartridges in combination with the VICI Unions on page 56. They can be used with other LP solvent line connectors like our Tees, Crosses & Manifolds on page 58 or our VICI Manual Flow Selection Valves on page 60.

### In-Line Low Pressure Filter Cartridges

- Easy to Replace
- Compact Design
- Fits All 1/4"-28 Fitting Details
- Three Different Porosities

These filters are convenient since they can be simply dropped into any 1/4"-28 fitting detail. Therefore, all LP solvent line connectors like unions, tees, crosses etc. can be used as cartridge holder. The filter is made

of a stainless steel screen pressed in an inert cartridge. The inner design of the cartridge ensures the equal distribution of the solvent to the screen. (Nut and union shown not included)



### In-Line Low/Medium Pressure Filter for LC

Part No.	Description	Qty/pkg
JR-CFE-S2-5	In-Line Filter Cartridge, 2 µm	5
JR-CFE-S10-5	In-Line Filter Cartridge, 10 µm	5
JR-CFE-S75-5	In-Line Filter Cartridge, 75 µm	5



## PEEK In-Line Filter Kit – High Pressure

- 100% Biocompatible
- Minimal Hold-up Volume
- Biocompatible PAT Frits or Ti Frits
- Different Porosities

The VICI Jour In-Line Filter traps fines and other particles from samples and mobile phases before they damage valuable instruments and columns.

This design is made entirely of PEEK for biocompatibility and chemical resistance. PAT (PEEK Alloyed with

Teflon) filter elements are used for complete biocompatibility. The design has virtually no hold-up volume and can be used in analytical applications with virtually no band broadening or loss of efficiency.



### In-Line Filter Kit – High Pressure

Part No.	Description	Internal Volume $\mu\text{L}$ *
JR-68250	Filter, PEEK, In-Line, PAT frit PEEK-encased 5 $\mu\text{m}$	12.4
JR-68251	Filter, PEEK, In-Line, PAT frit PEEK-encased 10 $\mu\text{m}$	12.4
JR-68247	Filter, PEEK, In-Line, Ti frit PEEK-encased 0.5 $\mu\text{m}$	4.4
JR-68253	Filter, PEEK, In-Line, Ti frit PEEK-encased 2 $\mu\text{m}$	5.8

### Spare Parts

Part No.	Description	Frit volume $\mu\text{L}$	Qty/pkg
JR-68152-5	Frit, PAT, PEEK-encased 5 $\mu\text{m}$	9.77	5
JR-9000-0460-5	Frit, PAT, PEEK-encased 10 $\mu\text{m}$	9.77	5
JR-1125-05P-5	Frit, Ti, PEEK-encased 0.5 $\mu\text{m}$	1.77	5
JR-1125-2P-5	Frit, Ti, PEEK-encased 2 $\mu\text{m}$	3.20	5

### SPECS

#### Material

Body: PEEK

Frit: see chart

#### Dimensions

Bore 0.25 mm/0.4 mm

See illustration

#### Threads

10-32 female

#### Pressure rating

< 350 bar (< 5000 psi)

#### Tech Tip

We recommend 2  $\mu\text{m}$  frits for columns with 5  $\mu\text{m}$  or larger particles and 0.5  $\mu\text{m}$  frits for smaller particles.

\* The internal volume is the total flushed volume between the ends of the connecting tubes, including the frit volume.

Calculation of frit volume: see Tech Info on page 109

### SPECS

#### Material

Body: PEEK  
Frit: see chart

#### Dimensions

Bore 0.25 mm/0.4 mm  
See illustration

#### Threads

10-32 female to 10-32 male

#### Pressure rating

< 350 bar (< 5000 psi)

#### Tech Tip

We recommend 2 µm frits for columns with 5 µm or larger particles and 0.5 µm frits for smaller particles.

\* The internal volume is the hold-up volume between the end of the connecting tube and the filter 10-32 connector outlet, including the frit volume.

Calculation of frit volume:  
see Tech Info on page 109

## PEEK Pre-Column Filter Kit – High Pressure

- 100% Biocompatible
- Minimal Hold-up Volume
- Biocompatible PAT Frits or Ti Frits
- Different Porosities

The VICI Jour Pre-Column Filter Kit has a standard 10-32 fitting and fits direct to most column types without introducing an additional dead volume.

The Filter Kit protects expensive columns against fines and particles, which may otherwise accumulate on the column frit, leading to split peaks and high backpressure.



### PEEK Pre-Column Filter Kit – High Pressure

Part No.	Description	Internal Volume µL*
JR-68260	Filter, PEEK, pre-column, PAT frit PEEK-encased 5 µm	14.4
JR-68261	Filter, PEEK, pre-column, PAT frit PEEK-encased 10 µm	14.4
JR-68258	Filter, PEEK, pre-column, Ti frit PEEK-encased 0.5 µm	6.4
JR-68262	Filter, PEEK, pre-column, Ti frit PEEK-encased 2 µm	7.8

### Spare Parts

Part No.	Description	Frit volume µL	Qty/pkg
JR-68152-5	Frit, PAT, PEEK-encased 5 µm	9.77	5
JR-9000-0460-5	Frit, PAT, PEEK-encased 10 µm	9.77	5
JR-1125-05P-5	Frit, Ti, PEEK-encased 0.5 µm	1.77	5
JR-1125-2P-5	Frit, Ti, PEEK-encased 2 µm	3.20	5





## Sure-Guard

- Minimized Hold-up Volume
- Disposable Filter Guard
- Biocompatible Ti and SS Frits
- Two Different Porosities

The VICI Jour Sure-Guard disposable In-Line filter offers an easy and inexpensive way of protecting valuable columns against fines and particles. It is easily connected directly to any column with an

inlet for 1/16" OD tubes and 10-32 threads. The VICI Jour Sure-Guard can be changed in seconds without tools.



## Sure-Guard

Part No.	Description	Frit volume $\mu\text{L}$	Qty/pkg
JR-0611-SS05-3	Sure-Guard, disposable In-Line Filter, SS frit 0.5 $\mu\text{m}$	0.45	3
JR-0611-SS2-3	Sure-Guard, disposable In-Line filter, SS frit 2 $\mu\text{m}$	0.61	3
JR-0611-Ti05-3	Sure-Guard, disposable In-Line filter, Ti frit 0.5 $\mu\text{m}$	0.45	3
JR-0611-Ti2-3	Sure-Guard, disposable In-Line filter, Ti frit 2 $\mu\text{m}$	0.61	3

## SPECS

### Material

Body: PEEK  
Filter: SS316, Titanium  
(see chart)

### Dimensions

Bore 0.4 mm  
See illustration

### Threads

10-32 female to 10-32 male

### Pressure rating

< 350 bar (< 5000 psi)

### Tech Tip

We recommend 2  $\mu\text{m}$  frits for columns with 5  $\mu\text{m}$  or larger particles and 0.5  $\mu\text{m}$  frits for smaller particles.

### Spares & Tools

We recommend our striped PEEK tubing for easy identification. (see pages 5-6)

Calculation of frit volume:  
see Tech Info on page 109

### SPECS

#### Material

Body: SS316,  
Frit: PEEK encased SS316

#### Dimensions

Bores: see chart  
Others: see illustrations

#### Filter Dimensions

See page 85

#### Threads

10-32 female to 10-32 female

#### Tolerances

+/- 0.05 mm (.002")

#### Pressure rating

< 350 bar (< 5000 psi)

#### Tech Tip

We recommend 2 µm frits for columns with 5 µm or larger particles and 0.5 µm frits for smaller particles.

\* The internal volume is the flushed volume between the ends of the connecting tubes, including the frit volume.

Calculation of frit volume:  
see Tech Info on page 109

## SS High-Pressure In-Line Filter

- Minimized Hold-up Volume
- High-Pressure Applications up to 350 bar (5000 psi)
- Two Different Bores for Analytical up to Semi Prep Flow Rates
- Three Different Porosities

The ViCl Jour SS In-Line Filter has virtually no hold-up volume and can be used in various applications without significant band broadening or loss of efficiency. The 0.25 mm bore versions are intended for analytical use between autosampler/injection valve and column. The 0.75 mm bore versions are suitable

for higher flow rates in semi prep systems or for solvent supply from pump to autosampler with minimal dead volume. These in-line filters are the ideal column protectors and available with frit porosities of 0.5, 2 or 5 µm.



### SS High-Pressure In-Line Filter

Part No.	Description	Internal Volume µL*
JR-68230-05	Filter, SS, In-Line, 0.25 mm bore, SS 0.5 µm	0.2
JR-68230-2	Filter, SS, In-Line, 0.25 mm bore, SS 2 µm	0.3
JR-68230-5	Filter, SS, In-Line, 0.25 mm bore, SS 5 µm	0.3
JR-68231-05	Filter, SS, In-Line, 0.75 mm bore, SS 0.5 µm	1.7
JR-68231-2	Filter, SS, In-Line, 0.75 mm bore, SS 2 µm	1.9
JR-68231-5	Filter, SS, In-Line, 0.75 mm bore, SS 5 µm	1.9

### Spare Parts

Part No.	Description	Frit volume µL	Qty/pkg
JR-1110-05P-5	Frit, SS, PEEK-encased 0.5 µm, for P/N 68230 series	0.14	5
JR-1110-2P-5	Frit, SS, PEEK-encased 2 µm, for P/N 68230 series	0.18	5
JR-1110-5P-5	Frit, SS, PEEK-encased 5 µm, for P/N 68230 series	0.20	5
JR-1111-05P-5	Frit, SS, PEEK-encased 0.5 µm, for P/N 68231 series	0.56	5
JR-1111-2P-5	Frit, SS, PEEK-encased 2 µm, for P/N 68231 series	0.75	5
JR-1111-5P-5	Frit, SS, PEEK-encased 5 µm, for P/N 68231 series	0.81	5







CC0172

### INDEX

Safety-Coated Glass Bottles	77
Opti-Cap	78
Eco-Cap	78
SolviFlex	79
Check Valve for Helium Sparging Line	80

## Safety-Coated Glass Bottles

- Special Design for Helium Sparging of LC Mobile Phases
- Safety Coating to Minimize Risk of Injury

VICI Jour Safety-Coated Glass Bottles are available in 1 and 2 L volume sizes.

Schott Duran® glass bottles with GL45 screw threading are coated on the outside with a layer of sturdy polymer. This reduces the likelihood of the bottle breaking if it is dropped. Even if the bottle breaks,

glass fragments will be kept within the polymer coating. Solvents will be kept within the coating for a short while.

The coating also protects solvents against UV light up to 385 nm.



C0487

## Safety-Coated Glass Bottles

Part No.	Description
JR-9000-0005	Bottle, Glass, safety-coated 1 L with solid cap
JR-9000-0010	Bottle, Glass, safety-coated 2 L with solid cap

## SPECS

### Material

Bottle: Duran® glass  
Cap: Polypropylene  
Pouring Ring: Polypropylene

### Pressure Rating

0.8 bar (11.6 psi)

### Special Info

The bottles are not suitable for use with pressures above 0.8 bar (11.6 psi), use only for continuous Helium sparging.

The safety coating is not intended to replace safe working conditions and precautions in laboratories.

### SPECS

#### Material

Screw Collar: Polyethylene  
Pouring Ring: Polypropylene  
Opti-Cap insert: PTFE  
Eco-Cap insert: ETFE  
O-Rings: EPDM  
Plugs: PEEK  
Luer Plugs: Polypropylene  
Others: see charts

#### Dimensions

For GL45 or S40 threads  
(see chart)

#### Bores

Opti-Cap: 2 x 1/8", 1 x 1/16"  
Eco-Cap: 3 x 1/8"

#### Threads

Opti-Cap: 1/4"-28

#### Special Info

These caps are not suitable for building up a Helium atmosphere within the solvent bottle, use only for continuous Helium sparging.

#### Tech Tip

Once the mobile phase is degassed, we recommend our No-Ox tubing between bottle cap and pump to prevent "regassing" of Helium degassed solvents (see page 18)

### Opti-Cap

- Ideal to Helium Sparge and Deliver Mobile Phase
- Also Available as Bottle Kit

The VICI Jour Opti-Cap is the most economical way to Helium sparge and deliver mobile phases. The Opti-Cap (PTFE insert with polyethylene collar) is available for standard GL45 threads or S40 threads used by Merck®. It has 3 ports for tubing insertion, 2 x 1/8" and 1 x 1/16", for different options (e.g. Helium sparging). The tubing ports are made so that you push the tubing through the hole or you can use any 1/4"-28 fittings. If required, just plug the port that you do not need.

The Kit consists of all needed parts to sparge and deliver LC mobile phases like fittings, plugs and tubing. Bottle Kits incl. additionally 1L or 2L Safety Coated Glass Bottles



### Opti-Cap

Part No.	Description
JR-9000-0001	Opti-Cap, PTFE, GL45, incl. 1x Plug
JR-9000-0006	Opti-Cap, PTFE, S40, incl. 1x Plug

### Opti-Cap Kits

JR-9000-0002	Opti-Cap Kit
JR-9000-0003	Opti-Cap Kit, incl. 1L Safety Coated Bottle
JR-9000-0004	Opti-Cap Kit, incl. 2L Safety Coated Bottle

Spare Parts see next page

### Eco-Cap

- Competitively Priced Alternative for Sparging and delivering Mobile Phases
- Also Available as Bottle Kit

The Eco-Cap is made of ETFE insert with polyethylene collar and has three holes for tubing insertion. The Eco-Cap has no threading in the inside of the tubing ports – tubing can easily put through.



### Eco-Cap

Part No.	Description
JR-9000-0007	Eco-Cap, ETFE, GL45, incl. 1 x Male Luer Plug





## Spare Parts


Part No.	Description
JR-9000-0005	Bottle, Glass, safety-coated 1L with solid cap
JR-9000-0010	Bottle, Glass, safety-coated 2L with solid cap
JR-367008-2	Filter, SS, 2 µm, 1/8" Fitting
JR-367008-10	Filter, SS, 10 µm, 1/8" Fitting
JR-T-6800-M3	Tubing, PTFE, 1/8" x 1/16" ID, 3m/pkg
JR-T-4037-M3	Tubing, PTFE, 1/8" x 2.40 mm ID, 3m/pkg
JR-20113-10	Nut, PEEK, flange-free blue 1/8", 1/4"-28, 10/pkg
JR-051-10	Ferrule, ETFE, 1/8", 10/pkg
JR-409	Plug, PEEK, 1/4"-28
JR-070350-10	Plug, PP, Male Luer, 10/pkg

## SolviFlex

- Be 'Mobile' with Mobile Phase Bottles
- Available for GL45 and S40 Thread Sizes
- Customized Versions on Request

The SolviFlex has a flexible telephone cable like PFA tubing which expands and contracts and consequently gives the flexibility to move mobile phase bottles on the bench. All SolviFlex assemblies consist of a biocompatible mobile phase filter and all the re-

quired parts to connect it to any 1/4"-28 port. Standard cap assemblies available for bottles with thread sizes GL45 or Merck S40. Customized versions are available on request – the cap has to be supplied by the customer.



## SPECS

## Material (wetted parts)

GL45 Cap: Polypropylene  
Tubing: PFA  
Filter: PE  
Ferrules: ETFE

## Dimension

1/8" Tubing OD

## Threads

See chart

## Special Info

Complete with fittings 1/4"-28 and PE filter 5 µm.

Color and design of S40 caps may differ from illustration

## SolviFlex

Part No.	Description
JR-6120-50030	Solviflex, assembly without cap
JR-6120-50031	Solviflex, assembly with GL45 cap
JR-6120-50032	Solviflex, assembly with S40 cap
JR-6120-50033	Solviflex, customized version, cap by customer

#### SPECS

##### Material

Body: Polypropylene  
Spring: Stainless Steel  
Unions: Nylon  
Adapters: Polypropylene  
Nuts: POM  
Ferrules: ETFE

##### Dimensions

For OD 1/16" or 1/8" Tubing

##### Threads

1/4"-28

##### Cracking pressure


0.1 bar (1.5 psi)

### Check Valve for Helium Sparging Line

- Prevents Backflow of Liquids
- Fits 1/16" or 1/8" OD Tubing

Low Pressure Check Valves are a simple and inexpensive way to prevent backflow from solvent reservoir into the sparging line. When a Helium supply is switched off, the Helium gas will diffuse through the tubing wall, allowing the solvent to be sucked into

the sparging line. The Check Valve which is located just outside the bottle protects gas line components such as regulators and valves from any liquid and also eliminates cross-contamination of the mobile phase reservoir.



#### Low pressure Check Valve

Part No.	Description
JR-9001-0620	Check Valve, Anti-Backflow, fittings for 1/16" tubing
JR-9001-0621	Check Valve, Anti-Backflow, fittings for 1/8" tubing







CM079

## INDEX

PAT Frits	83
Stainless Steel Frits – Unidense Type	84
Stainless Steel Frits – PEEK Encased	85
Titanium Frits – Unidense Type	86
Titanium Frits – PEEK Encased	87

## PAT Frits

- PAT - PEEK Alloyed with Teflon®
- Superior Chemical Resistance
- Excellent Mechanical Strength

PAT stands for PEEK Alloyed with Teflon®, a porous material for frits and filters in HPLC. This material is made by compressing PEEK and Teflon® particles under rigidly controlled conditions to form PAT porous sheets. PAT frits are suitable for use with all HPLC and IC mobile phases, including strong buffers and chloride salts which corrode common stainless steel frits. Compared with other porous polymer frits (Polyethylene, Polypropylene, PTFE), PAT frits exhibit superior chemical resistance and mechanical strength. PAT frits may be used with any HPLC column or filter but they are particularly well suited for separations involving

proteins, peptides, nucleic acids or other samples of biological origin. PAT frits are also ideal for IC systems where metal frits may interact with ionic samples or release transition metals that will deactivate ion exchange columns.

PAT Frits have a superior mechanical strength. The strength is derived from the thermally bonded structure of the PEEK and Teflon® alloy. Therefore, very thin PAT Filters are strong enough to be used without support rings and can easily be pressed into column ends.



## PAT Frits

Part No.	Diameter in Inches	Diameter in mm	Thickness in Inches	Thickness in mm	Porosity in µm	Frit volume in µL	Qty/pkg
JR-6145-3062-5	0.062	1.59	0.031	0.78	5	0.54	5
JR-6145-3125-5	0.125	3.17	0.031	0.78	5	2.18	5
JR-6145-6188-5	0.188	4.77	0.062	1.59	5	9.87	5
JR-6145-3500-5	0.500	12.70	0.031	0.78	5	34.91	5
JR-6146-3250-5	0.250	6.35	0.031	0.78	10	8.73	5
JR-6146-3500-5	0.500	12.70	0.031	0.78	10	34.91	5

## PAT Frits PEEK encased

Part No.	Ring OD in Inches	Ring OD in mm	Frit OD in Inches	Frit OD in mm	Thickness in Inches	Thickness in mm	Porosity in µm	Frit volume in µL	Qty/pkg
JR-9000-0477-5	0.2	5.08	0.102	2.59	0.062	1.59	5	2.91	5
JR-68152-5	0.25	6.35	0.187	4.77	0.062	1.59	5	9.77	5
JR-9000-0460-5	0.25	6.35	0.187	4.77	0.062	1.59	10	9.77	5
JR-67152-5	0.375	9.52	0.312	7.92	0.062	1.59	5	27.19	5

Other dimensions and porosities are available on request. Minimum order quantity: 100 pcs. Please contact your local distributor or VICI directly.

## SPECS

## Material

PAT (PEEK Alloyed with Teflon®)  
PAT + PEEK ring

## Dimensions

See chart

## Tolerances

+/- 0.05 mm (.002")

## Pressure rating

< 350 bar (< 5000 psi)

## Special Info

Compatible with most HPLC and IC Mobile Phases - Even those that corrode Stainless Steel. For further reference see chemical compatibility of PEEK.

## Note

Frit volumes are theoretical and are calculated by multiplying overall frit volume times the porosity proportion. PAT frits have a porosity proportion of 35%.

Stated frit porosities are only nominal and do not reflect the maximum pore size of a frit. The chart on page 109 will help to choose the right frit for your application.




## Stainless Steel Frits – Unidense Type

## Stainless Steel Frits – Unidense Type

For both standard frits and those produced to custom specifications, the frits offer unexcelled reproducibility and precision. Proprietary manufacturing processes ensure uniformity and precise, repeatable dimensional control. VICI Stainless Steel Frits are man-

ufactured of SS316L material which has extra low carbon content. The sintering process used results in a strong bonding of the sinter metal particles. Depending upon the grade of powder used, the porosity of the frit can be precisely controlled.



## Stainless Steel Frits – Unidense Type

Part No.	Diameter in Inches	Diameter in mm	Thickness in Inches	Thickness in mm	Porosity in $\mu\text{m}$	Frit volume in $\mu\text{L}$	Qty/pkg
JR-.5FR1-5	1/16"	1.59	0.030	0.76	0.5	0.39	5
JR-1FR1-5	1/16"	1.59	0.030	0.76	1	0.45	5
JR-2FR1-5	1/16"	1.59	0.030	0.76	2	0.53	5
JR-10FR1-5	1/16"	1.59	0.030	0.76	10	0.60	5
JR-.5FR2-5	1/8"	3.18	0.040	1.02	0.5	2.09	5
JR-1FR2-5	1/8"	3.18	0.040	1.02	1	2.41	5
JR-2FR2-5	1/8"	3.18	0.040	1.02	2	2.81	5
JR-10FR2-5	1/8"	3.18	0.040	1.02	10	3.22	5
JR-.5FR4-5	1/4"	6.35	0.040	1.02	0.5	8.37	5
JR-2FR4-5	1/4"	6.35	0.040	1.02	2	11.26	5
JR-10FR4-5	1/4"	6.35	0.040	1.02	10	12.87	5
JR-.5FR6-5	3/8"	9.53	0.040	1.02	0.5	18.82	5
JR-2FR6-5	3/8"	9.53	0.040	1.02	2	25.34	5
JR-10FR6-5	3/8"	9.53	0.040	1.02	10	28.96	5
JR-.5FR8-5	1/2"	12.70	0.040	1.02	0.5	33.46	5
JR-2FR8-5	1/2"	12.70	0.040	1.02	2	45.05	5
JR-10FR8-5	1/2"	12.70	0.040	1.02	10	51.48	5
JR-2FR1K-5	1"	25.40	0.060	1.52	2	270.28	5
JR-10FR1K-5	1"	25.40	0.060	1.52	10	308.89	5

Other dimensions and porosities are available on request. Minimum order quantity: 100 pcs.  
Please contact your local distributor or VICI directly.





## Stainless Steel Frits – PEEK Encased

- For HPLC Columns or In-Line Filters
- Low Dead Volume
- Various Porosities and Dimensions Available

VICI Jour PEEK encased frits with 1/4" ring OD are intended for HPLC columns with ID's of 3.2, 3.9 and 4.6 mm. By surrounding the frit with a PEEK ring the actual porous surface is reduced to the column ID.

The compressed PEEK ring forms a seal between the bottom of the fitting and the column end prevents the mobile phase and the sample from entering the previously accessible poorly swept areas.



## Stainless Steel Frits – PEEK Encased

Part No.	Ring OD in Inches	Ring OD in mm	Frit OD in Inches	Frit OD in mm	Thickness in Inches	Thickness in mm	Porosity in $\mu\text{m}$	Frit volume in $\mu\text{L}$	Qty/pkg
JR-1110-05P-5	0.118	3.00	0.038	0.97	0.028	0.70	0.5	0.14	5
JR-1110-2P-5	0.118	3.00	0.038	0.97	0.028	0.70	2	0.18	5
JR-1110-5P-5	0.118	3.00	0.038	0.97	0.028	0.70	5	0.2	5
JR-1111-05P-5	0.118	3.00	0.077	1.96	0.028	0.70	0.5	0.56	5
JR-1111-2P-5	0.118	3.00	0.077	1.96	0.028	0.70	2	0.75	5
JR-1111-5P-5	0.118	3.00	0.077	1.96	0.028	0.70	5	0.81	5
JR-1104-05P-5	1/4	6.35	0.078	2.00	0.062	1.59	0.5	1.26	5
JR-1104-2P-5	1/4	6.35	0.078	2.00	0.062	1.59	2	1.7	5
JR-1103-05P-5	1/4	6.35	0.181	4.60	0.032	0.80	0.5	3.51	5
JR-1103-2P-5	1/4	6.35	0.181	4.60	0.032	0.80	2	4.72	5
JR-1102-05P-5	1/4	6.35	0.125	3.20	0.062	1.59	0.5	3.24	5
JR-1102-2P-5	1/4	6.35	0.125	3.20	0.062	1.59	2	4.36	5
JR-1101-05P-5	1/4	6.35	0.153	3.90	0.062	1.59	0.5	4.86	5
JR-1101-2P-5	1/4	6.35	0.153	3.90	0.062	1.59	2	6.54	5
JR-1100-05P-5	1/4	6.35	0.181	4.60	0.062	1.59	0.5	6.8	5
JR-1100-2P-5	1/4	6.35	0.181	4.60	0.062	1.59	2	9.15	5

Other dimensions and porosities are available on request. Minimum order quantity: 100 pcs. Please contact your local distributor or VICI directly.

## SPECS

### Material

Frits: SS316L  
Ring: PEEK

### Dimensions

See chart

### Tolerances

$\pm 0.05$  mm (.002")

### Special Info

Eliminates Poorly Swept Volume – Improves Peak Symmetry.

### Note

Stated frit porosities are only nominal and do not reflect the maximum pore size of a frit. The chart on page 109 will help to choose the right frit for your application.

## SPECS

**Material**  
Titanium (sintered)

**Dimensions**  
See chart

**Tolerances**  
OD < 5 mm,  
+/- 0.05 mm (.002")  
OD 5–12.5 mm,  
+/- 0.08 mm (.003")  
OD 12.5–25.4 mm,  
+/- 0.20 mm (.008")

**Special Info**  
Standard frits are without  
chamfers – chamfered frits are  
available on request


**Note**  
Stated frit porosities are only  
nominal and do not reflect the  
maximum pore size of a frit.  
The chart on page 109 will help  
to choose the right frit for your  
application.

## Titanium Frits – Unidense Type

- Biocompatible – Recommended for Protein Analysis
- Higher Corrosion Resistance
- Various Pore Sizes and Dimensions

Titanium is preferred over stainless steel for applications related to the analysis of particularly sensitive substances in general and large biomolecules in particular. Evidently proteins have a tendency to adsorb

on the stainless steel (Iron) frit and consequently reduce recovery – or – may even decompose. Titanium also offers higher corrosion resistance and can be used with solvents containing halides.



## Titanium Frits – Unidense Type

Part No.	Diameter in Inches	Diameter in mm	Thickness in Inches	Thickness in mm	Porosity in µm	Frit volume in µL	Qty/pkg
JR-2FR1TI-5	1/16"	1.59	0.040	1.02	2	0.70	5
JR-2FR2TI-5	1/8"	3.18	0.040	1.02	2	2.81	5
JR-2FR4TI-5	1/4"	6.35	0.040	1.02	2	11.26	5
JR-2FR6TI-5	3/8"	9.53	0.040	1.02	2	25.34	5
JR-2FR8TI-5	1/2"	12.70	0.040	1.02	2	45.05	5
JR-10FR8TI-5	1/2"	12.70	0.040	1.02	10	51.48	5
JR-2FR1KTI-5	1"	25.40	0.040	1.02	2	270.28	5

Other dimensions and porosities are available on request. Minimum order quantity: 100 pcs.  
Please contact your local distributor or VICI directly.



## Titanium Frits – PEEK Encased

- For HPLC Columns or In-Line Filters
- Low Dead Volume
- Biocompatible – Recommended for Protein Analysis
- Higher Corrosion Resistance
- Various Pore Sizes and Dimensions

In terms of "biocompatibility" the frit is the critical spot for contributing metal ions to the mobile phase and interaction with protein samples. VICI Jour PEEK encased frits with 1/4" ring OD are intended for HPLC columns with ID's of 3.2 and 4.6 mm. By surrounding the frit with a PEEK ring the actual porous surface is

reduced to the column ID. The compressed PEEK ring forms a seal between the bottom of the fitting and the column end prevents the mobile phase and the sample from entering the previously accessible poorly swept areas.



## Titanium Frits – PEEK Encased

Part No.	Ring OD in Inches	Ring OD in mm	Frit OD in Inches	Frit OD in mm	Thickness in Inches	Thickness in mm	Porosity in $\mu\text{m}$	Frit volume in $\mu\text{L}$	Qty/ pkg
JR-1127-05P-5	1/4"	6.35	0.125	3.20	0.062	1.59	0.5	0.84	5
JR-1127-2P-5	1/4"	6.35	0.125	3.20	0.062	1.59	2	1.53	5
JR-1125-05P-5	1/4"	6.35	0.181	4.60	0.062	1.59	0.5	1.77	5
JR-1125-2P-5	1/4"	6.35	0.181	4.60	0.062	1.59	2	3.20	5

Other dimensions and porosities are available on request. Minimum order quantity: 100 pcs.  
Please contact your local distributor or VICI directly.

## SPECS

### Material

Frit: Titanium  
Ring: PEEK

### Dimensions

See chart

### Tolerances

+/- 0.05 mm (.002")

### Special Info

Eliminates Poorly Swept Volume – Improves Peak Symmetry

### Note

Stated frit porosities are only nominal and do not reflect the maximum pore size of a frit. The chart on page 109 will help to choose the right frit for your application.





C0584

### INDEX

PEEK Columns	89
In-Line Biocompatible Guard Column	92
Column Coupler	
One-Piece Fingertight	93
Column End Plugs	93

## PEEK Columns

- Biocompatible
- High Pressure Applications
- Packing Devices Available

### PEEK Columns 2.1 mm ID

We have designed this column for narrow-bore applications and where only limited amounts of sample are available for analysis. Because 2.1 mm ID columns operate at 1/5 of the flow rate of their 4.6 mm counterparts, the same sample mass injected on a 2.1 mm column will produce five times the detector signal. This greater mass sensitivity provides the same analysis as a 4.6 mm column with only 20% of the sample. Only PEEK and PAT contact the fluid stream.

Wetted Surfaces: PEEK and PAT frits  
Pressure rating: Max. 210 bar (3000 psi)

\* using the VICI Jour Packing Kit JR-67-250-4:  
max. pressure during packing 350 bar (5000 psi)

### PEEK Columns 7.5 mm ID

Ideal for Semi-prep Biochromatography. The basic column design is shown in the figure and consists of a length of threaded PEEK tubing and two PEEK end fittings with the packed bed held in place by two PAT frits. Connections are made through PEEK inlet and outlet fittings which accept 1/16" OD tubing and standard fingertight fittings. Max operating pressure 350 bar (5000 psi). Low Protein Binding ensures high recovery and accurate quantification.

Wetted Surfaces: PEEK and PAT frits  
Pressure rating: Up to 350 bar (5000 psi)

### PEEK Columns 4.6 mm & 4.0 mm ID

Today's HPLC and IC applications are more demanding than ever. The limitations of traditional stainless steel pose significant problems for a growing number of important biotechnology and ion chromatography applications.

Wetted Surfaces: PEEK and PAT (PEEK Alloyed with Teflon) frits  
Pressure Rating: Max. 350 bar (5000 psi) continuous. During packing Max. 560 bar (8000 psi)\*




### SPECS

**Material**  
PEEK, black and natural  
PAT

**Dimensions**  
See chart


**Tolerances**  
+/- 0.05 mm (.002")

**Pressure rating**  
See page 89




### PEEK Columns 2.1 mm ID

Part No.	Description
JR-65001	Column, PEEK, 30 x 2.1 mm ID, complete
JR-65002	Column, PEEK, 50 x 2.1 mm ID, complete
JR-65003	Column, PEEK, 100 x 2.1 mm ID, complete
JR-65005	Column, PEEK, 150 x 2.1 mm ID, complete
JR-65006	Column, PEEK, 250 x 2.1 mm ID, complete




### PEEK Columns 4.0 mm ID

Part No.	Description
JR-68181	Column, PEEK, 50 x 4.0 mm ID, complete
JR-68182	Column, PEEK, 100 x 4.0 mm ID, complete
JR-68183	Column, PEEK, 125 x 4.0 mm ID, complete
JR-68184	Column, PEEK, 150 x 4.0 mm ID, complete
JR-68185	Column, PEEK, 250 x 4.0 mm ID, complete
JR-68186	Column, PEEK, 300 x 4.0 mm ID, complete



### PEEK Columns 4.6 mm ID

Part No.	Description
JR-66175	Column, PEEK, 50 x 4.6 mm ID, complete
JR-66176	Column, PEEK, 100 x 4.6 mm ID, complete
JR-66180	Column, PEEK, 125 x 4.6 mm ID, complete
JR-66177	Column, PEEK, 150 x 4.6 mm ID, complete
JR-66178	Column, PEEK, 250 x 4.6 mm ID, complete
JR-66179	Column, PEEK, 300 x 4.6 mm ID, complete



### PEEK Columns 7.5 mm ID

Part No.	Description
JR-67202	Column, PEEK, 50 x 7.5 mm ID, complete
JR-67222	Column, PEEK, 75 x 7.5 mm ID, complete
JR-67204	Column, PEEK, 100 x 7.5 mm ID, complete
JR-67206	Column, PEEK, 150 x 7.5 mm ID, complete
JR-67208	Column, PEEK, 250 x 7.5 mm ID, complete
JR-67220	Column, PEEK, 300 x 7.5 mm ID, complete



### Packing Devices

Part No.	Description
JR-67-250-2	Packing Kit, PEEK, for 2.1 mm ID columns, complete
JR-67-250-4	Packing Kit, PEEK, for 4.0 and 4.6 mm ID columns, complete
JR-67-600	Packing Kit, PEEK, for 7.5 mm ID columns, complete

### Spare Parts

Part No.	Description
JR-062	Endfitting, PEEK, for 2.1 mm ID columns
JR-66140	Endfitting, PEEK, for 4.0 and 4.6 mm ID columns
JR-67140	Endfitting, PEEK, for 7.5 mm ID columns
JR-9000-0477	Frit, PAT, PEEK-encased, 5 µm for 2.1 mm ID columns
JR-68152	Frit, PAT, PEEK-encased, 5 µm for 4.0 and 4.6 mm ID columns
JR-67152	Frit, PAT, PEEK-encased, 5 µm for 7.5 mm ID columns
JR-67190	Frit spacer, PEEK, for 7.5 mm ID columns
JR-67140-1A	Packing adapter, PEEK, for 2.1 mm ID columns
JR-67140-2B	Packing adapter, PEEK, for 4.0 and 4.6 mm ID columns
JR-67140-3A	Packing adapter, PEEK, for 7.5 mm ID columns
JR-65011	Replacement packing seal, PEEK, for 2.1 mm ID columns
JR-66158	Replacement packing seal, PEEK, for 4.0 and 4.6 mm ID columns

### SPECS

**Material**  
PEEK, PAT

#### Pressure rating

For 2.1 mm ID columns: 210 bar (3000 psi)

For 4.0 to 7.5 mm ID columns: 350 bar (5000 psi)

#### Special note

Our Packing Kits are 'ready-to-use' tools to pack our PEEK columns.

Every Kit consists of the needed slurry reservoir, packing adapter, packing seals etc.

Further technical support material is available on our webpage

## In-Line Biocompatible Guard Column

### SPECS

**Material**  
PEEK, PAT

**Dimensions**  
See illustration

**Pressure rating**  
< 350 bar (< 5000 psi)

**Special note**  
Our Packing Kit is a 'ready-to-use' tool to pack our PEEK guard cartridges. Every Kit consists of the needed slurry reservoir, packing adapter, packing seals etc.

Further technical support material is available on our webpage

## In-Line Biocompatible Guard Column

- Biocompatible
- Easy to Fill
- Packing Device for Self Packing Available

VICI Jour In-Line PEEK Guard Column System extends the life of valuable analytical columns by trapping contaminants from the sample and mobile phase, thereby preventing loss of retention. Only PEEK and PAT contact the mobile phase and sample. This system is ideal for use with proteins or other biopoly-

mers that may be contaminated or absorbed by metal surfaces. VICI Jour Guard Columns can be easily dry packed or slurry packed. Guard columns may be changed in minutes without the use of tools - all connections may be made by hand.



## In-Line PEEK Guard Column

Part No.	Description
JR-68252	In-Line guard cartridge holder, PEEK, complete with 1 guard cartridge JR-68242
JR-68242	Column, PEEK, guard cartridge 8.0 x 3.0 mm ID, complete 2 pcs PAT-frit 5 µm

## Packing Device

Part No.	Description
JR-27305	Packing Kit, PEEK, for In-Line guard cartridge system, complete


## Spare Parts

Part No.	Description
JR-6145-6188	Replacement Frit, PAT, 5 µm for guard cartridge
JR-9	Packing seal, CTFE

## Column Coupler One-Piece Fingertight

- Zero Dead Volume Connections
- Fits all Column Designs
- Biocompatible
- Self Adjusting
- Four Bore Sizes

The unique feature of the VICI Jour Column Coupler is that its design fits all pilot lengths, such as Valco®, Waters®, Upchurch®, Rheodyne®, etc. It adapts automatically and bottoms out in any fitting detail and consequently creates only small void volumes.



### One-Piece Fingertight Column Coupler

Part No.	Description
JR-26501	Coupler, PEEK, universal 0.13 mm ID, red
JR-26502	Coupler, PEEK, universal 0.17 mm ID, yellow
JR-26503	Coupler, PEEK, universal 0.25 mm ID, blue
JR-26504	Coupler, PEEK, universal 0.50 mm ID, orange

## Column End Plugs

- Seals Analytical Columns
- One Size – Fits All
- With Wrench Flats

The column end plugs have a fingertight design. One size fits all columns with 10-32 threads. Easy handling is guaranteed by the additional 5/16" wrench flat, which allows loosening the plug with a wrench

in case it should be stuck. The long design and the threading down to the ferrule ensure a safe closing of valuable columns.



### Column End Plugs

Part No.	Color	Qty/pkg
JR-4071-10	Plug, Nylon, for HPLC columns, white, 10-32	10
JR-4072-10	Plug, Nylon, for HPLC columns, black, 10-32	10
JR-4073-10	Plug, Nylon, for HPLC columns, red, 10-32	10
JR-4074-10	Plug, Nylon, for HPLC columns, green, 10-32	10
JR-4075-10	Plug, Nylon, for HPLC columns, blue, 10-32	10

### SPECS

#### Material

Tubing: PEEK  
Fittings: PEEK  
Sleeve: Polypropylene

#### Dimensions

See illustration

#### Threads

10-32 male to 10-32 male

#### Pressure rating

< 350 bar (< 5000 psi)

### SPECS

#### Material

Nylon

#### Dimensions

For 10-32 fitting details





C0579

### INDEX

ValvTool	95
PEEK Sample Loops for Valco Valves	96
Stainless Steel Sample Loops for Valco Valves	97
PEEK Sample Loops for Cheminert Valves	98
Stainless Steel Sample Loops for Cheminert Valves	99
PEEK Sample Loops for Rheodyne® Valves	100
Stainless Steel Sample Loops for Rheodyne® Valves	101



## ValvTool

- Handy Tool for the Lab
- Tighten Nuts with a Tube or Loop installed
- Fits Nuts with 1/4" or 5/16" Heads

The new ValvTool is a time-saving device, which will provide easy access to many hard-to-reach areas. The unique design with its slotted wrench allows tightening nuts where a loop or a capillary may other-

wise make it difficult. The ValvTool is ideal for Valco, Cheminert and other valves and fittings as well as for most of the SS and PEEK fittings on the market with 1/4" and 5/16" heads.



## SPECS

### Dimensions

For nuts with 1/4" and 5/16" head

## ValvTool

Part No.	Description
JR-800	ValvTool, 1/4" and 5/16" open wrench

### SPECS

#### Material

Loops: PEEK  
Nuts: PEEK  
Ferrules: PEEK  
Caps: Polyethylene

#### Dimensions

See chart

Tolerances of tube ID  
Premium Grade PEEK Tubing  
(see page 5)  
+/- 0.025 mm (.001")

#### Note

VICI Jour PEEK Sample Loops are made of Premium Grade PEEK Tubing. The tighter ID tolerance of the tubing minimizes the tolerance of the loop volume.

All PEEK sample loops will be supplied with PEEK hex-head nuts and PEEK ferrules for reliable connections.

This design avoids "twisting" of the tube.



2x JR-5005 + 1 x JR-5510

## PEEK Sample Loops for Valco Valves

- Fits All Valco W and UW Type Valves
- Wide Range of Volumes Available
- Biocompatible
- Made from Premium Grade PEEK Tubing

VICI Jour Sample Loops are made from PEEK tubing and fittings. PEEK is a mechanically strong, chemically inert polymer ideal for HPLC applications. Since VICI Jour PEEK Sample Loops contain no metals or other extractable compounds, no foreign materials will be introduced into the system. Further softener free PE caps are used to avoid contamination of the tubing.

VICI Jour PEEK Sample Loops are particularly well suited for separations involving proteins, peptides nucleic acids or other samples of biological origin as well as in IC systems where metal traces may interact with ionic samples or release transition metals that may deactivate ion exchange columns.



### PEEK Sample Loops for Valco Valves 1/16" OD

Part No.	Description	bar*	psi*
JR-SLV-75-5	5 µL	400	5800
JR-SLV-75-10	10 µL	400	5800
JR-SLV-75-20	20 µL	400	5800
JR-SLV-75-50	50 µL	400	5800
JR-SLV-75-100	100 µL	400	5800
JR-SLV-75-250	250 µL	240	3500
JR-SLV-75-500	500 µL	240	3500
JR-SLV-75-1000	1 mL	240	3500

Other dimensions are available on request. Please contact your local distributor or VICI directly.  
\* = Recommended maximum working pressure at room temperature

### Spare Parts

Part No.	Description	Qty/pkg
JR-5510-5	Fitting, PEEK, Hex-head long, 10-32, 1/16"	5
JR-5005-5	Single Ferrule, PEEK, 1/16"	5





## Stainless Steel Sample Loops for Valco Valves

- Fits Valco W Type Valves
- Wide Range of Volumes Available

VICI Jour Stainless Steel Sample Loops are made of high quality SS316 and are cut using an electrolytic cutting process followed by ultrasonic and high pres-

sure steam cleaning to remove contaminants, dried and capped with softener free caps to avoid the re-contamination of the tubing.



### Stainless Steel Loop for Valco Valves 1/16" OD

Part No.	Description
JR-SLV-76-5	5 µL
JR-SLV-76-10	10 µL
JR-SLV-76-20	20 µL
JR-SLV-76-50	50 µL
JR-SLV-76-100	100 µL
JR-SLV-76-250	250 µL
JR-SLV-76-500	500 µL
JR-SLV-76-1000	1 mL

Other dimensions are available on request. Please contact your local distributor or VICI directly.

### Spare Parts

Part No.	Description	Qty/pkg
JR-ZN1-5	Nut, SS, 10-32, 1/16"	5
JR-ZF1S6-5	Ferrule, SS, 1/16"	5

## SPECS

### Material

Loops: SS316  
Nuts: SS303  
Ferrules: SS316  
Caps: Polyethylene

### Dimensions

See chart

### Tolerances of tube ID

Premium Grade SS Tubing  
+/- 0.025 mm (.001")

### Pressure rating

< 350 bar (5000 psi)

### Note

All Loops will be supplied with SS Hex-Head nuts and SS ferrules for reliable connections.

Not recommended for use with halides. High pressure resistance, low permeability

### Tech Tip

We offer nuts in three lengths to facilitate wrench access on valves. For loops frequently replaced we recommend the extra long nut – particularly when space is restricted (see page 40).

## Sample Loops

### PEEK Sample Loops for Cheminert Valves

#### SPECS

##### Material

Loops: PEEK  
Nuts: PEEK  
Ferrules: PEEK  
Caps: Polyethylene

##### Dimensions

See chart

##### Tolerances of tube ID

Premium Grade PEEK Tubing  
(see page 5)  
+/- 0.025 mm (.001")

##### Note

VICI Jour PEEK Sample Loops are made of Premium Grade PEEK Tubing. The tighter ID tolerance of the tubing minimizes the tolerance of the loop volume.

All PEEK sample loops will be supplied with PEEK hex-head nuts and PEEK ferrules for reliable connections.

This design avoids "twisting" of the tube.



2x JR-5005 + 1 x JR-5510

### PEEK Sample Loops for Cheminert Valves

- Fits All Cheminert C1, C2, C2V and C3 Type Valves
- Wide Range of Volumes Available
- Biocompatible
- Made from Premium Grade PEEK Tubing

VICI Jour Sample Loops are made from PEEK tubing and fittings. PEEK is a mechanically strong, chemically inert polymer ideal for HPLC applications. Since VICI Jour PEEK Sample Loops contain no metals or other extractable compounds, no foreign materials will be introduced into the system. Further softener free PE caps are used to avoid contamination of the tubing.

VICI Jour PEEK Sample Loops are particularly well suited for separations involving proteins, peptides nucleic acids or other samples of biological origin as well as in IC systems where metal traces may interact with ionic samples or release transition metals that may deactivate ion exchange columns.



### PEEK Sample Loops for Cheminert Valves 1/16" OD

Part No.	Description	bar*	psi*
JR-SLC-85-5	5 µL	400	5800
JR-SLC-85-10	10 µL	400	5800
JR-SLC-85-20	20 µL	400	5800
JR-SLC-85-50	50 µL	400	5800
JR-SLC-85-100	100 µL	400	5800
JR-SLC-85-250	250 µL	240	3500
JR-SLC-85-500	500 µL	240	3500
JR-SLC-85-1000	1 mL	240	3500

Other dimensions are available on request. Please contact your local distributor or VICI directly.

\* = Recommended maximum working pressure at room temperature

### Spare Parts

Part No.	Description	Qty/pkg
JR-5510-5	Fitting, PEEK, Hex-head long, 10-32, 1/16"	5
JR-5005-5	Single Ferrule, PEEK, 1/16"	5

## Stainless Steel Sample Loops for Cheminert Valves

- Fits All Cheminert C1, C2, C2V and C3 Type Valves
- Wide Range of Volumes Available

VICI Jour Stainless Steel Sample Loops are made of high quality SS316 and are cut using an electrolytic cutting process followed by ultrasonic and high pres-

sure steam cleaning to remove contaminants, dried and capped with softener free caps to avoid the re-contamination of the tubing.



### Stainless Steel Loop for Cheminert Valves 1/16" OD

Part No.	Description
JR-SLC-86-5	5 µL
JR-SLC-86-10	10 µL
JR-SLC-86-20	20 µL
JR-SLC-86-50	50 µL
JR-SLC-86-100	100 µL
JR-SLC-86-250	250 µL
JR-SLC-86-500	500 µL
JR-SLC-86-1000	1 mL

Other dimensions are available on request. Please contact your local distributor or VICI directly.

### Spare Parts

Part No.	Description	Qty/pkg
JR-ZN1-5	Nut, SS, 10-32, 1/16"	5
JR-ZF1S6-5	Ferrule, SS, 1/16"	5

### SPECS

#### Material

Loops: SS316  
Nuts: SS303  
Ferrules: SS316  
Caps: Polyethylene

#### Dimensions

See chart

#### Tolerances of tube ID

Premium Grade SS Tubing  
+/- 0.025 mm (.001")

#### Pressure rating

< 350 bar (5000 psi)

#### Note

All Loops will be supplied with SS Hex-Head nuts and SS ferrules for reliable connections.

Not recommended for use with halides. High pressure resistance, low permeability

#### Tech Tip

We offer nuts in three lengths to facilitate wrench access on valves. For loops frequently replaced we recommend the extra long nut – particularly when space is restricted (see page 40).



### SPECS

#### Material

Loops: PEEK  
Nuts: PEEK  
Ferrules: PEEK  
Caps: Polyethylene

#### Dimensions

See chart

#### Tolerances of tube ID

Premium Grade PEEK Tubing  
(see page 5)  
+/- 0.025 mm (.001")

#### Note

VICI Jour PEEK Sample Loops are made of Premium Grade PEEK Tubing. The tighter ID tolerance of the tubing minimizes the tolerance of the loop volume.

All PEEK sample loops will be supplied with PEEK hex-head nuts and PEEK ferrules for reliable connections.

This design avoids "twisting" of the tube.



2x JR-5005 + 1 x JR-5510

## PEEK Sample Loops for Rheodyne® Valves

- Fits All Rheodyne® Metal and Metal Free Valves
- Wide Range of Volumes Available
- Biocompatible
- Made from Premium Grade PEEK Tubing

VICI Jour Sample Loops are made from PEEK tubing and fittings. PEEK is a mechanically strong, chemically inert polymer ideal for HPLC applications. Since Our PEEK Sample Loops contain no metals or other extractable compounds, no foreign materials will be introduced into the system. Further softener free PE caps are used to avoid contamination of the tubing. Our PEEK Sample Loops are particularly well suited for separations involving proteins, peptides, nucleic acids or other samples of biological origin as well as in IC systems where metal traces may interact with ionic samples or release transition metals that may deactivate ion exchange columns.



### PEEK Sample Loops for Rheodyne® Valves 1/16" OD

Part No.	Description	bar*	psi*
JR-SLR-2190	1 µL	435	6300
JR-SLR-2191	2 µL	420	6100
JR-SLR-2192	3 µL	420	6100
JR-SLR-2193	4 µL	420	6100
JR-SLR-2194	5 µL	400	5800
JR-SLR-2196	10 µL	400	5800
JR-SLR-2195	20 µL	386	5600
JR-SLR-2200	30 µL	386	5600
JR-SLR-2201	40 µL	386	5600
JR-SLR-2202	50 µL	350	4500
JR-SLR-2204	100 µL	350	4500
JR-SLR-2206	200 µL	240	3500
JR-SLR-2207	250 µL	240	3500
JR-SLR-2210	500 µL	240	3500
JR-SLR-2222	1 mL	240	3500

Other dimensions are available on request. Please contact your local distributor or VICI directly.

\* = Recommended maximum working pressure at room temperature

### Spare Parts

Part No.	Description	Qty/pkg
JR-5510-5	Fitting, PEEK, Hex-head long, 10-32, 1/16"	5
JR-5005-5	Single Ferrule, PEEK, 1/16"	5



## Stainless Steel Sample Loops for Rheodyne® Valves

- Fits All Rheodyne Valves
- Special Loops for Valve 7725
- Wide Range of Volumes Available

VICI Jour Stainless Steel Sample Loops are made of high quality SS316 and are cut using an electrolytic cutting process followed by ultrasonic and high pres-

sure steam cleaning to remove contaminants, dried and capped with softener free caps to avoid the re-contamination of the tubing.



### Stainless Steel Loops for Rheodyne® Valves 1/16" OD

Part No.	Description
JR-SLR-111	5 µL
JR-SLR-112	10 µL
JR-SLR-113	20 µL
JR-SLR-114	50 µL
JR-SLR-115	100 µL
JR-SLR-116	200 µL
JR-SLR-117	500 µL
JR-SLR-118	1 mL

Other dimensions are available on request. Please contact your local distributor or VICI directly.

### Stainless Steel Loops for Rheodyne® Valves No 7725 1/16" OD

Part No.	Description
JR-SLR-111-77	5 µL
JR-SLR-112-77	10 µL
JR-SLR-113-77	20 µL
JR-SLR-114-77	50 µL
JR-SLR-115-77	100 µL
JR-SLR-116-77	200 µL
JR-SLR-117-77	500 µL
JR-SLR-118-77	1 mL

Other dimensions are available on request. Please contact your local distributor or VICI directly.

### Spare Parts

Part No.	Description	Qty/pkg
JR-57-5	Nut, SS, long head, Rheodyne® type, for Valve No. 7725	5
JR-58-5	Nut, SS, x-long, Rheodyne® type,	5
JR-59-5	Ferrule, SS, 1/16"	5
JR-60-5	Ferrule, SS, 1/16", 2-step Rheodyne® type	5

### SPECS

#### Material

Loops: SS316  
Nuts: SS316  
Ferrules: SS316  
Caps: Polyethylene

#### Dimensions

See chart

#### Tolerances of tube ID

Premium Grade SS Tubing  
+/- 0.025 mm (.001")

#### Pressure rating

< 350 bar (5000 psi)

#### Note

All Loops will be supplied with Rheodyne® type SS Hex-Head nuts and SS ferrules for reliable connections.

Not recommended for use with halides. High pressure resistance, low permeability

#### Tech Tip

We offer nuts in three lengths to facilitate wrench access on valves. For loops frequently replaced we recommend the extra long nut – particularly when space is restricted (see page 41).

All SS Loops for Rheodyne® are available with the 2-step SS Ferrule Rheodyne® type on request.



C0431

## INDEX

Backpressure Regulators	103
VICI Mininert™ Valves for Vials	104



## Backpressure Regulators

- Improves Baseline Stability by Preventing Bubble Formation in Flow Cell
- Adjustable Backpressure between 0–105 bar (1500 psi) without Disconnecting
- Low Dead Volume
- Biocompatible

Even properly degassed mobile phases may contain some dissolved gases that can release bubbles in the detector flow cell resulting in baseline noise and drift. The Backpressure Regulator provides a quick and convenient way to improve the stability of detector baseline by providing a constant backpressure on the flow cell which stops the release of dissolved gases from the mobile phase. A spring-loaded FFKM (Perfluorinated rubber) diaphragm provides the force that generates the backpressure. Spring tension may be adjusted using the backpressure adjustment screw

on the top of the regulator to vary the backpressure without disconnecting.

Because of the unique design, the Backpressure Regulator has only 0.6  $\mu$ L internal volume. Installed between two detectors in a multidetector system, or between a detector and fraction collector system, or between a detector and fraction collector it will not affect the peak shape. Typically, 90% or more of a given component's efficiency will be retained upon passing through the regulator. And the fluid path is completely biocompatible!



### Backpressure Regulators

Part No.	Description
JR-39021	Variable between 0–7 bar (0–100 psi)
JR-39022	Variable between 7–20 bar (100–300 psi)
JR-39023	Variable between 20–40 bar (300–600 psi)
JR-39024	Variable between 40–105 bar (600–1500 psi)

### SPECS

#### Material

PEEK  
FFKM (Perfluorinated rubber)

#### Dimensions

For OD 1/16" Tubing

#### Threads

10-32

#### Tolerances

+/- 0.05 mm (.002")

#### Pressure rating

See Chart

#### Special Info

6  $\mu$ L Internal Volume  
May be used between detector and other devices without affecting the performance

## Backpressure Regulators

- Improves Baseline Stability by Preventing Bubble Formation in Flow Cell
- Adjustable Backpressure between 0–105 bar (1500 psi) without Disconnecting
- Low Dead Volume
- Biocompatible

Even properly degassed mobile phases may contain some dissolved gases that can release bubbles in the detector flow cell resulting in baseline noise and drift. The Backpressure Regulator provides a quick and convenient way to improve the stability of detector baseline by providing a constant backpressure on the flow cell which stops the release of dissolved gases from the mobile phase. A spring-loaded FFKM (Perfluorinated rubber) diaphragm provides the force that generates the backpressure. Spring tension may be adjusted using the backpressure adjustment screw

on the top of the regulator to vary the backpressure without disconnecting.

Because of the unique design, the Backpressure Regulator has only 0.6  $\mu$ L internal volume. Installed between two detectors in a multidetector system, or between a detector and fraction collector system, or between a detector and fraction collector it will not affect the peak shape. Typically, 90% or more of a given component's efficiency will be retained upon passing through the regulator. And the fluid path is completely biocompatible!



### Backpressure Regulators

Part No.	Description
JR-39021	Variable between 0–7 bar (0–100 psi)
JR-39022	Variable between 7–20 bar (100–300 psi)
JR-39023	Variable between 20–40 bar (300–600 psi)
JR-39024	Variable between 40–105 bar (600–1500 psi)

### SPECS

#### Material

PEEK  
FFKM (Perfluorinated rubber)

#### Dimensions

For OD 1/16" Tubing

#### Threads

10-32

#### Tolerances

+/- 0.05 mm (.002")

#### Pressure rating

See Chart

#### Special Info

6  $\mu$ L Internal Volume  
May be used between detector and other devices without affecting the performance

## PEEK & Polymer Chemical Resistance

PEEK exhibits excellent resistance to a wide range of organic and inorganic chemicals. The compatibility of PEEK with many chemicals at 20 °C (68 °F) has been investigated and the results for unreinforced grades are given in the table below.

PEEK is compatible with almost any of the solvents used in HPLC. The only solvent which will attack PEEK are concentrated nitric acid and sulphuric acids.

However, PEEK tubing can safely withstand 20–30% nitric acid when passivating a system.

Methylene Chloride, DMSO and THF may cause swelling in PEEK. The highest temperature we recommend for PEEK is 100 °C. Up to this temperature the tubing will maintain the pressure rating stated.

A = Suitable

B = Marginal-dependent on application

C = Not recommended

### PEEK and Polymer chemical resistance

Resistance at 20 °C	PEEK	Poly-ethylene	Poly-propylene	PPS	PVDF	PTFE	ETFE
Acetaldehyde	A						
Acetic acid (20%)	A	A	A	A	A	A	A
Acetic acid (80%)	A	A	A	A	B	A	
Acetic acid (glacial)	A	A	A	A	A	A	A
Acetone	A	B	A	A	C	A	A
Acetonitrile	A			A		A	A
Acrylic acid	A			A			
Ammonia, anhydrous	A			A			
Ammonia (10%)	A	B	A	A	A	A	A
Ammonia (Liquid)	B			A			
Ammonium hydroxide	A	A	A	A	A	A	A
Aqua regia	C			C			
Aromatic hydrocarbons	A	B	C	A			
Benzene	A	B	C	A	A	A	A
Benzoic acid	A			A			
Benzaldehyde	A			A			
Bromine/dibromoethane	C			C			
Bromine (dry)	C			C			
Bromine (wet)	C			C			
Boric acid	A			A			
Butanol	A	A	A	A	A	A	A
Calcium hydroxide	A			A			
Carbon tetrachloride	A			A			
Chlorine (gas)	A			C			
Chlorine (liquid)	C			C			
Chloroacetic acid	A	B	B	A	A	A	A
Chlorobenzene	A			A			
Chloroform	A	B	B	A	A	A	A
Cyclohexane	A	B	C	A	A	A	A
Cyclohexanone	A	C	C	A	C	A	A
Diethylamine	A	C	A		C	A	A
Diethylether	A			A			
Diethylformamide	A	A	A	A	C	A	A
Dioxane	A			A			A
Ethanol	A	B	A	A		A	A
Ether	A	B	C	A	B	A	A
Ethyl acetate	A	B	A	A	C	A	A
Ethylene chloride		B	B	A	A	A	A
Ethylene glycol	A	A	A	A	A	A	A
Heptane	A	B	B	A	A	A	A
Hexane	A	B	B	A	A	A	A
Hydrobromic acid (100%)	C	B	B	A	A	A	

### INDEX

PEEK & Polymer Chemical Resistance 106–107

Conversions 108–109

Length Conversions  
& Some Physical Data  
about the Polymer  
& PEEK material we use 109



## General Conversion Factors

Inches to mm:	Inches x 25.4 = mm
Feet to Meter:	Feet x 0.3048 = Meter
Meter to Feet:	Meter x 3.28 = Feet
bar to psi:	bar x 14.5 = psi
psi to bar:	psi x 0.069 = bar
psi to KPa:	psi x 6.9 = KPa

## Pressure Conversions

bar	psi	KPa	MPa
0.1	1.5	10	0.01
0.2	2.9	20	0.02
0.3	4.4	30	0.03
0.4	5.8	40	0.04
0.5	7.3	50	0.05
0.6	8.7	60	0.06
0.7	10.2	70	0.07
0.8	11.6	80	0.08
0.9	13.1	90	0.09
1	14.5	100	0.1
2	29.0	200	0.2
3	43.5	300	0.3
4	58.0	400	0.4
5	72.5	500	0.5
6	87.0	600	0.6
7	101.5	700	0.7
8	116.0	800	0.8
9	130.5	900	0.9
10	145.0	1000	1
20	290.1	2000	2
30	435.1	3000	3
40	580.1	4000	4
50	725.2	5000	5
60	870.2	6000	6
70	1015.2	7000	7
80	1160.3	8000	8
90	1305.3	9000	9
100	1450.3	10000	10
200	2900.7	20000	20
300	4351.0	30000	30
400	5801.3	40000	40
500	7251.6	50000	50
600	8702.0	60000	60
700	10152.3	70000	70
800	11602.6	80000	80
900	13052.9	90000	90
1000	14503.3	100000	100

## Dimension Conversions

mm	Inches	mm	Inches	Inches	mm
0.064	0.0025"	3.00	0.118"	1/32"	0.79
0.13	0.005"	4.00	0.157"	1/16"	1.59
0.17	0.007"	4.60	0.181"	3/32"	2.38
0.25	0.010"	6.00	0.236"	1/8"	3.17
0.50	0.020"	6.35	0.250"	1/4"	6.35
0.75	0.030"	6.40	0.252"	5/16"	7.94
1.00	0.039"	7.00	0.276"	1/2"	12.70
1.50	0.059"	10.00	0.394"	1"	25.40
2.00	0.079"	12.70	0.500"	1ft	304.8

## Temperature Conversions

°C	°F	°C	°F	°C	°F
-40	-40	135	275	310	590
-35	-31	140	284	315	599
-30	-22	145	293	320	608
-25	-13	150	302	325	617
-20	-4	155	311	330	626
-15	5	160	320	335	635
-10	14	165	329	340	644
-5	23	170	338	345	653
0	32	175	347	350	662
5	41	180	356	375	707
10	50	185	365	400	752
15	59	190	374	425	797
20	68	195	383	450	842
25	77	200	392	475	887
30	86	205	401	500	932
35	95	210	410	525	977
40	104	215	419	550	1022
45	113	220	428	575	1067
50	122	225	437	600	1112
55	131	230	446	625	1157
60	140	235	455	650	1202
65	149	240	464	675	1247
70	158	245	473	700	1292
75	167	250	482	725	1337
80	176	255	491	750	1382
85	185	260	500	775	1427
90	194	265	509	800	1472
95	203	270	518	825	1517
100	212	275	527	850	1562
105	221	280	536	875	1607
110	230	285	545	900	1652
115	239	290	554	925	1697
120	248	295	563	950	1742
125	257	300	572	975	1787
130	266	305	581	1000	1832

## Length Conversions

mm	Inches	mm	Inches	Inches	mm
0.064	0.0025"	3.00	0.118"	1/32"	0.79
0.13	0.005"	4.00	0.157"	1/16"	1.59
0.17	0.007"	4.60	0.181"	3/32"	2.38
0.25	0.010"	6.00	0.236"	1/8"	3.17
0.50	0.020"	6.35	0.250"	1/4"	6.35
0.75	0.030"	6.40	0.252"	5/16"	7.94
1.00	0.039"	7.00	0.276"	1/2"	12.70
1.50	0.059"	10.00	0.394"	1"	25.40
2.00	0.079"	12.70	0.500"		

## Volume Chart

ID/mm	ID/inch	μl/cm	μl/inch	ID/mm	ID	μl/cm	μl/inch
0.050	.002"	0.02	0.05	1.00	.040"	7.85	20.59
0.064	.0025"	0.03	0.08	1.40	.055"	15.39	38.93
0.075	.003"	0.04	0.12	1.52	.060"	18.15	46.33
0.10	.004"	0.08	0.21	1.59	.062"	19.86	49.47
0.13	.005"	0.13	0.32	1.65	.065"	21.38	54.38
0.17	.0067"	0.23	0.58	1.70	.067"	22.70	57.78
0.18	.007"	0.25	0.63	1.78	.070"	24.88	63.06
0.25	.010"	0.49	1.29	2.00	.079"	31.42	80.32
0.38	.015"	1.13	2.90	2.10	.083"	34.64	88.66
0.50	.020"	1.96	5.15	2.16	.085"	36.64	92.99
0.75	.030"	4.42	11.58	2.40	.094"	45.24	113.72

This chart is only a guideline. Because of ID tolerances for our tubing (+/-0.025mm (.001") for premium grade tubing and +/- 0.050 mm (.002") for standard tubing) volumes have to be measured to get an exact result.

## Frit Volume

Internal frit volumes are theoretical and calculated by multiplying the total frit volume by the proportion of the pores in the frit (porosity proportion). Stated

frit porosities are only nominal and do not reflect the maximum pore size of a frit. The chart below will help you in selecting the right frit for your application.

## SS and Titanium Frits

Porosity in μm	Max. pore size in μm	Particles retained in μm	Porosity in %
0.2	1	1.0	22
0.5	2	1.5	26
1	3	3.0	30
2	5	3.5-5.0	35
5	15	15.0	38


## Standard Threads

Place fitting on the fitting detail below to identify the thread size.


6-40  
(for 1/32" OD tubing)



10-32  
(for 1/16" OD tubing)




5/16"-24  
(for 1/8" OD tubing)



1/4"-28  
(for 1/16" and 1/8" OD tubing)



1/2"-20  
(for 1/4" OD tubing)





### Properties of Metals

#### Electroformed nickel (EFNI)

We electroplate pure nickel over a diamond drawn mandrel in a continuous process. Then we carefully separate and remove the mandrel from the tubing. The result is an extremely inert and smooth interior surface - an incredible 1–2 microinch finish. It is widely used for transfer lines, since it minimizes the potential for carryover or cross contamination often found with mill drawn Nickel 200, due to its rough interior surface. Unlike glass or silica-lined stainless, EFNI can easily accept tight bends and cutting without heating, and does not release damaging glass fragments or silica particles. Electroformed nickel has more in common with fused silica than drawn nickel tubing in terms of surface inertness and smoothness.

#### Stainless steel, type 316

This is the standard tubing material for chromatography, suitable for a wide variety of applications. It is cold drawn seamless, not welded, with close tolerances held on both ID and OD. We neither recommend nor offer Type 304 stainless steel for analytical applications.

Austenitic stainless steels may be used for most chromatographic applications. Type 316 is most commonly used for HPLC because of its superior chloride ion resistance.

#### Stainless steel, type 303

Recommended for GC use and general purpose connections. It combines excellent machining characteristics with good resistance to corrosion and high temperature oxidation. However, it is susceptible to attack by chlorides, iodides, and bromides.

#### Titanium

Although it is more difficult to machine than common alloys containing aluminum and vanadium, Valco uses Grade 2 pure titanium in order to avoid possible contamination of the sample stream with these metals. Good for organic and inorganic salts except aluminum and calcium chlorides, and all alkalis except boiling concentrated potassium hydroxide.

Good with dilute, low temperature formic, lactic, sulfuric, hydrochloric, and phosphoric acids, but rapidly attacked by hydrofluoric acid. Good with dilute nitric acid at low temperatures; corrodes at high concentrations and temperatures. Can ignite with fuming nitric acid. Attacked by oxalic acid, concentrated phosphoric acid, hot trichloroacetic acid, and zinc chloride.



### Properties of Metals

#### Electroformed nickel (EFNI)

We electroplate pure nickel over a diamond drawn mandrel in a continuous process. Then we carefully separate and remove the mandrel from the tubing. The result is an extremely inert and smooth interior surface - an incredible 1–2 microinch finish. It is widely used for transfer lines, since it minimizes the potential for carryover or cross contamination often found with mill drawn Nickel 200, due to its rough interior surface. Unlike glass or silica-lined stainless, EFNI can easily accept tight bends and cutting without heating, and does not release damaging glass fragments or silica particles. Electroformed nickel has more in common with fused silica than drawn nickel tubing in terms of surface inertness and smoothness.

#### Stainless steel, type 316

This is the standard tubing material for chromatography, suitable for a wide variety of applications. It is cold drawn seamless, not welded, with close tolerances held on both ID and OD. We neither recommend nor offer Type 304 stainless steel for analytical applications.

Austenitic stainless steels may be used for most chromatographic applications. Type 316 is most commonly used for HPLC because of its superior chloride ion resistance.

#### Stainless steel, type 303

Recommended for GC use and general purpose connections. It combines excellent machining characteristics with good resistance to corrosion and high temperature oxidation. However, it is susceptible to attack by chlorides, iodides, and bromides.

#### Titanium

Although it is more difficult to machine than common alloys containing aluminum and vanadium, Valco uses Grade 2 pure titanium in order to avoid possible contamination of the sample stream with these metals. Good for organic and inorganic salts except aluminum and calcium chlorides, and all alkalis except boiling concentrated potassium hydroxide.

Good with dilute, low temperature formic, lactic, sulfuric, hydrochloric, and phosphoric acids, but rapidly attacked by hydrofluoric acid. Good with dilute nitric acid at low temperatures; corrodes at high concentrations and temperatures. Can ignite with fuming nitric acid. Attacked by oxalic acid, concentrated phosphoric acid, hot trichloroacetic acid, and zinc chloride.

## Properties of Polymers

### CTFE

Chlorotrifluoroethylene, is the generic name for the material produced as Kel-F®. It is very resistant to all chemicals except THF and some halogenated solvents, and is resistant to all inorganic corrosive liquids, including oxidizing acids. CTFE can be used at temperatures up to 100 °C. Swells in ketones.

### EPDM

Ethylene Propylene Diene Monomer. Ethylene-propylene synthetic rubbers offer excellent heat resistance, oxidation, ozone and weather aging due to their stable, saturated polymer backbone structure. As non-polar elastomers, they have good electrical resistivity, as well as resistance to polar solvents, such as water, acids, alkalies, phosphate esters and many ketones and alcohols.

### ETFE

Ethyltrifluoroethylene is the generic name for the material such as Tefzel®. A fluoropolymer used for sealing surfaces, it is resistant to most chemical attack; however, some chlorinated chemicals will cause a physical swelling of ETFE tubing.

### FEP

Fluorinated ethylene propylene is another member of the fluorocarbon family with similar chemical properties. It is generally more rigid than PTFE, with somewhat increased tensile strength. It is typically more transparent than PTFE, slightly less porous, and less permeable to oxygen. FEP is not as subject to compressive creep at room temperature as PTFE, and because of its slightly higher coefficient of friction is easier to retain in a compression fitting.

### FFKM (Perfluorinated rubber)

FFKM (Perfluorinated rubber) is a Perfluoroelastomer. This is a class of materials such as Kalrez®. It is chemically resistant to polar solvents, organic solvents, inorganic and organic acids and bases, fuels, oils, lubricants, inorganic salts, aldehydes, metal halogen compounds, chlorine, sodium hydroxide, aromatics, alcohols, steam, and strong oxidizing agents. The normal temperature service range is -40°C to 315°C and up to 343°C in intermittent service.

### PEEK

Considered relatively inert and biocompatible, polyetheretherketone tubing can withstand temperatures up to 100 °C. Under the right circumstances, .005" – .020" ID tubing can be used up to 5000 psi for a limited time, and .030" to 3000 psi. Larger IDs are typically good to 500 psi.

These limits will be substantially reduced at elevated temperatures and in contact with some solvents or acids.

Its mechanical properties allow PEEK to be used instead of stainless in many situations and in some environments where stainless would be too reactive. However, PEEK can be somewhat absorptive of solvents and analytes, notably methylene chloride,

DMSO, THF, and high concentrations of sulfuric and nitric acid. This tubing is highly prone to "kinking", or sealing off, if held in a sharp bend over time.

### PFA

Perfluoroalkoxy is a fluorocarbon with chemical and mechanical properties similar to FEP. More rigid than either PTFE or FEP. Commonly used for injection molded parts.

### PPS

Polyphenylene sulphide is the generic name for the material produced as Fortron®, Ryton®, and Techtron®. It is very resistant to all solvents, acids, and bases.

### PTFE

Polytetrafluoroethylene is the generic name for the class of materials such as Teflon®. It offers superior chemical resistance but is limited in pressure and temperature capabilities. Because it's so easy to handle, it is often used in low pressure situations where stainless steel might cause adsorption. PTFE tubing is relatively porous, and compounds of low molecular weight can diffuse through the tubing wall.

### Polyacetal

Polyacetal or Polyoxymethylene (POM) is the generic name for the material produced as Delrin®. It is an engineering polymer with high mechanical strength and rigidity, excellent dimensional stability and excellent resistance to moisture and solvents. It has a wide useful temperature range in air of -50° C to +90°C. It is not suitable for use with acids, alkalis and oxidizing agents and has limited resistance to dichloromethane and methyl ethyl ketone.

### Polyamide

Polyamide is a thermoplastic polymer also known as Nylon. Nylon has very good physical properties but moisture can have a significant effect. It exhibits very good heat resistance and excellent chemical and wear resistance.

### Polyethylene

Polyethylene is a semi-crystalline material with excellent chemical resistance, and good fatigue and wear resistance. Polyethylenes provide good resistance to organic solvents, staining, and have low moisture absorption rates.

### Polypropylene

Widely used polymer for non-wetted parts. Attacked by strong oxidizers, aromatic and chlorinated hydrocarbons.

### PVDF

PVDF, polyvinylidene fluoride, has excellent resistance to most mineral and organic acids, aliphatic and aromatic hydrocarbons, and halogenated solvents. Poor resistance to acetone, MEK, THF, and potassium and sodium hydroxide. Often supplied as Kynar®.



## Glossary

## A

**Adapter**

A type of fitting which provides a method of joining two components of differing thread types or systems.

**Analytical column**

A long narrow tube packed or coated with one of many available chemically diverse compounds that can separate the components in a sample according to their boiling point, polarity, molecular size, or combination thereof. A column of some kind is used with most chromatographic techniques.

## B

**Backflush**

The use of valving to reverse the flow through a column in order to "backflush" or purge heavier components from the column. See application demonstrating backflush.

**Biocompatibility**


Defines the materials used in a system (i.e. fittings, tubing, and valves) that do not change the bioactivity of the biological substances that come into contact with the surface of these materials. Note that in chromatographic systems, the tubing and column (frits, packing material and wall) contribute over 99% of the surface area and the valves and fittings are insignificant.

**Bore**

The diameter of the orifice through the fitting. See capillary bore, through-type bore, and large bore.

**Bulkhead fitting**

A type of fitting in which the fitting body is inserted through an instrument panel or mounting bracket, to which it is affixed with a mounting nut. The Valco fitting body is uniquely undercut so that it "bites" into the panel when the mounting nut is tightened, eliminating the need for a lock washer. Illustration

**Butt connection**

A type of connection in which the two tube ends are directly and squarely in contact, usually effected with a through-type union. Typically used with fused silica connections, or small bore metal tubing.

## C

**Cap**

A cap is used to dead-end a piece of tubing with a nut and ferrule attached.

**Capillary bore**

The smallest available standard orifice in a given fitting design (usually 0.25 mm). Typically denoted by suffix "C" in the product number.

**Collapsible ferrule**

A one-piece ferrule engineered to collapse as it is tightened. The collapse takes place in two very narrow areas, and results in a very effective seal with virtually no distortion of the tubing ID and no dead volume. \*Patent No. 6,575,501

**Compression fitting**


A style of fitting in which a threaded nut compresses a tapered ferrule onto tubing as the nut is tightened. Valco metal ferrules cut a ring into the tubing wall while polymer types rely on surface compression to form a seal and hold the tube by friction.

**Connecting volume**

The volume between two or more connections. This may be cleanly swept, thus not contributing to peak distortion, or may be "dead volume" such as that found in fittings with larger bores than the connecting tubing.

**Cross**

A type of distribution fitting which connects four pieces of tubing, arranging them in the pattern of a cross.

**Detail**

See fitting detail.



## Glossary, cont.

### Distribution fitting

A generic term for tees, crosses, and manifolds, used to provide multiple access points to "distribute" a gas or liquid through a system. CAUTION! Using a distribution fitting in reverse to merge multiple streams may create dead volume. Special manifolds are available for this application.

## F

### Female Nut

A type of compression fitting in which the nut has female threads.

### FIA

Flow Injection Analysis. A simple and versatile analytical technique for automating wet chemical analyses based on the manipulation of a sample zone formed from the injection of the sample into a continuous stream of fluid used as a carrier.

### Ferrule

One of the components of a compression fitting; the conical piece of metal or plastic that compresses onto the tube as it is forced into a tapered seat. Valco metal ferrules are unique in that they attach to and seal at the tube by cutting a shallow ring into it, instead of by actually swaging it. This is preferable since it introduces no flow restriction.

### Filter

A type of union or reducing union which traps the particulates in a stream. The filtering element is typically a mesh screen or sintered frit.

### Fingertight Fitting

A fingertight fitting is a knurled or winged the nut which is hand tightened to compress the ferrule. Note: Polymer type ferrules rely on surface compression to form a seal and hold the tube by friction.


### Fitting

A One-Piece combination of ferrule and nut with male threads.


### Fitting detail

A fitting detail has the thread to accept the nut, a tapered ferrule seat and a pilot for the tube, thus the counterpart where the nut, ferrule and tube make the seal and are held in place. The Valco fitting details have the perfect shape, with a smooth transition from the ferrule seat into the cylindrical and threaded part. Many other fitting details are poor in design where the ferrule is forced in a too shallow ferrule seat. The result is a 'knicked' ferrule which cannot apply its sealing force anymore to the front of the ferrule.

### Valco Design



### Other Style



### Flanged fitting

A type of fitting used with fluoropolymer tubing (PTFE, FEP) in which a flange is made at the tube end. Connections are made at the flange either by compressing the flange into a flat detail (typically 1/4"-28 threaded) or by butting two flanges together. A special flanging tool forms the flanges.

### Flangeless fitting

Similar in application to the flanged fitting, but the flange is not required. A ferrule system (collapsible ferrule or inverted ferrule) is used which grips/compresses the tube. This fitting type can be used with virtually any polymeric tubing since the tube end does not have to be formed, but simply square cut. Typically used in 1/4"-28 threaded fittings, it is usually interchangeable with flanged fittings.

### Frit

A filter element typically made of stainless, Hastelloy, Titanium, or polymers, usually 0.75 mm or 1 mm thick. Frits may provide better filtration than screens, but because they are thicker there is greater mixing potential, and they typically result in increased pressure drop.

## G

### GC

Gas Chromatography. An analytical method incorporating an injection system, analytical column, controlled temperature zone, and detector. An inert carrier gas moves the sample through the column, which separates the sample components into discrete bands which are measured as they pass through the detector.

### Guard column

A column used in series between the injector and analytical column to prevent certain types of components from entering the analytical column.

## Glossary, cont.

## H

**HPLC**

High Performance Liquid Chromatography. An analytical system consisting of an injector, pump, analytical column, and detector. Using a liquid mobile phase, the sample is pumped through the column, where it is separated into discrete sample component bands which are detected and measured as the bands elute from the column.

## I

**ID**

Internal diameter.

**Inert**

Technically, unreactive with other substances; however, in the instrumentation field, "inert" is a relative term. Often polymers are termed inert but are soluble in some fluids and can react with some compounds.

## L

**LC**

Liquid Chromatography. Any of a variety of low to medium pressure techniques which use a liquid mobile phase as the carrier to move sample. Similar to HPLC.

**Large bore**

A bore that is larger than the standard for a given fitting; a fitting ordered with a large bore will have a larger flow orifice than the standard or capillary bore fitting of the same design.

**Luer adapter**

An adapter that connects a tapered luer fitting (square nib) of a syringe to a tube or tube fitting.

## M

**Male Nut**

A type of compression fitting in which the nut has male threads.

**Make up**

The point at which a ferrule, nut, and tube are assembled in the fashion which will effect a leak-free seal. In most compression fittings, that is accomplished by compressing the tube with the small end of the ferrule. With Valco metal ferrules, the ferrule usually makes up on the tube by cutting a shallow ring in it.

**Manifold**

A type of distribution fitting in which a single source is directed to multiple outlets, or vice versa. Caution! Using a common distribution fitting in reverse to merge multiple streams may create dead volume. Special manifolds are available for this application.

**Microbore column**

A liquid chromatography column of narrow bore (typically 2 mm or less) for improved resolution.

## N

**Nanovolume**

Nanovolume fittings and valves are designed for minimum system volume, have typically bores of 100µm and are used in applications with extreme low flows e.g. 200 nL/min.

**NPT**

National Pipe Thread; a standardized tapered pipe fitting. See pipe thread.

**Nut**

The tensioning component of a compression fitting. As the threaded nut is tightened into the fitting detail, it pushes the ferrule forward into the tapered ferrule seat, causing it to make up on the tube.

## O

**OD**

Outside diameter.


**One Piece Fingertight**

A one-piece fingertight fitting is a knurled or winged nut with integrated ferrule. It is convenient to use since the ferrule is retrieved together with the nut from a port. Note: Polymer type ferrules rely on surface compression to form a seal and hold the tube by friction.

## P

**Pilot**

The tubing which extends beyond a ferrule in a made-up fitting. See also "Pilot depth", and "Zero dead volume" below.

**Pilot depth**

The length of the tubing diameter cavity beyond the tapered ferrule seat within a fitting detail. Valco fitting pilot depths are tightly controlled which allows interchanging of components without the risk of leaks or dead volume.

**Pipe thread**

The external or internal threads of a fitting designed conical thread



## Glossary, cont.

faces. This type of fitting does not "bottom out" in the detail. Typically used with Teflon tape or other compound to lubricate the threads; however, since the diffusion rate of air components through the Teflon tape is considerable, pipe fittings should not be used in systems where leakage rates are critical.

### Plug

A plug is used to block a fitting detail in a union, Tee, cross or valve.

## S

### Screen

A replaceable filter element generally made of Type 316 stainless steel, usually 0.003" thick. Screens clog less frequently than frits but they are less effective filters. Because they are thinner there is less hold-up volume and thus less mixing.

### Port

The connection, orifice, seal, or septum, etc. through which sample, mobile phase or eluents may be added or withdrawn.

## R

### Reducing ferrule

A ferrule which allows a smaller tube to be used in a fitting detail designed for a larger tube. Caution should be taken if standard reducing ferrules without integral pilots are used, since dead volume may be created in the fitting pilot depth.

### Reducing union

A fitting which joins two tubes of different ODs. The bore of the fitting should typically match the ID of the smaller tube.

## T

### Tee

A type of distribution fitting which connects three pieces of tubing, arranging them in the pattern of a "T".

### Through-type bore

A bore which is slightly larger than the OD of the tubing which is used with the given fitting. A union with a through-type bore allows the tube ends to butt directly together, or for one tube to run completely through the fitting. Denoted by suffix "T" in the product number. In order to assure correct pilot lengths, we recommend that ferrules be made up on the tubing in a standard union.

## U

### Union

A fitting for connecting two pieces of tubing of the same OD, or different OD's (see Reducing Union).

**Female Union** - A type of compression fitting in which the fitting body has female threads. Also referred to as internal union.

**Male Union** - A type of compression fitting in which the fitting body has male threads. Also referred to as external union.

### Unswept volume

The volume of any portion of a fitting which is in the flowpath but which is a different diameter than the primary flow orifice through the tubing/fitting assembly, or any area not directly swept by the fluid flow. This can also be known as "dead volume" if it is very poorly swept.

## W

### Wetted surfaces

The surfaces which are contacted by the sample stream.

## Y

### Y

A type of distribution fitting which connects three pieces of tubing, arranging them in the pattern of a "Y". Occasionally referred to as a "wye".


## Z

### ZDV

Zero dead volume. See definition below.

### Zero dead volume (ZDV)

Describes a connection which does not add volume to the system beyond what an extension of tubing would in its place.





## Part Number Index

Part. No.	Page	Part. No.	Page	Part. No.	Page
JR-.5FR1-5	84	JR-071124-10	62	JR-20125-10	52
JR-.5FR2-5	84	JR-071323-10	62	JR-20126-10	52
JR-.5FR4-5	84	JR-071324-10	62	JR-20127-10	52
JR-.5FR6-5	84	JR-072350-10	62	JR-20128-10	52
JR-.5FR8-5	84	JR-072450-10	62	JR-20129-10	52
JR-040BK-10	53	JR-1030	34	JR-201529	63
JR-040BL-10	53	JR-1032	34	JR-201531	25
JR-040GR-10	53	JR-1033	34	JR-201536	25
JR-040RE-10	53	JR-1040	34	JR-201537	25
JR-040WH-10	53	JR-1042	34	JR-201539	25
JR-040WL-10	53	JR-1043	34	JR-201540	24
JR-040YL-10	53	JR-1061	33	JR-201541	25
JR-041-10	50	JR-1064	56	JR-201554	25
JR-050BK-10	53	JR-1065	56	JR-201580BK-10	54
JR-050BL-10	53	JR-1066	33	JR-201580WH-10	54
JR-050GR-10	53	JR-1067	33	JR-201662	25
JR-050RE-10	53	JR-1068	56	JR-201663	25
JR-050WH-10	53	JR-10FR1-5	84	JR-201926-10	50
JR-050WL-10	53	JR-10FR1K-5	84	JR-201928-10	50
JR-050YL-10	53	JR-10FR2-5	84	JR-202118-10	54
JR-051-10	50	JR-10FR4-5	84	JR-202119-10	54
JR-060-5	56	JR-10FR6-5	84	JR-202120-10	54
JR-061	33	JR-10FR8-5	84	JR-202121-10	54
JR-0611	36	JR-10FR8TI-5	86	JR-202122-10	54
JR-0611-SS05-3	73	JR-1100-05P-5	85	JR-202123-10	54
JR-0611-SS2-3	73	JR-1100-2P-5	85	JR-202124-10	54
JR-0611-TI05-3	73	JR-1101-05P-5	85	JR-202125-10	54
JR-0611-TI2-3	73	JR-1101-2P-5	85	JR-202142-10	54
JR-0612	36	JR-1102-05P-5	85	JR-202143-10	54
JR-0613	36	JR-1102-2P-5	85	JR-202144-10	54
JR-0617	36	JR-1103-05P-5	85	JR-202145-10	54
JR-0619	36	JR-1103-2P-5	85	JR-202235	25
JR-062	36	JR-1104-05P-5	85	JR-2030	34
JR-0622	36	JR-1104-2P-5	85	JR-2032	34
JR-0623	36	JR-1110-05P-5	85	JR-2033	34
JR-0625	36	JR-1110-2P-5	85	JR-2040	34
JR-0626	36	JR-1110-5P-5	85	JR-2042	34
JR-064	56	JR-1111-05P-5	85	JR-2043	34
JR-065	56	JR-1111-2P-5	85	JR-26501	93
JR-066	33	JR-1111-5P-5	85	JR-26502	93
JR-067	33	JR-1125-05P-5	87	JR-26503	93
JR-068	56	JR-1125-2P-5	87	JR-26504	93
JR-070101-10	62	JR-1127-05P-5	87	JR-27305	92
JR-070111-10	62	JR-1127-2P-5	87	JR-2FR1-5	84
JR-070112-10	62	JR-1FR1-5	84	JR-2FR1K-5	84
JR-070113-10	62	JR-1FR2-5	84	JR-2FR1KTI-5	86
JR-070123-10	62	JR-201094BK-10	54	JR-2FR1TI-5	86
JR-070124-10	62	JR-201094WH-10	54	JR-2FR2-5	84
JR-070211-10	62	JR-20112-10	52	JR-2FR2TI-5	86
JR-070212-10	62	JR-20113-10	52	JR-2FR4-5	84
JR-070213-10	62	JR-20114-10	52	JR-2FR4TI-5	86
JR-070224-10	62	JR-20115-10	52	JR-2FR6-5	84
JR-070225-10	62	JR-20116-10	52	JR-2FR6TI-5	86
JR-070350-10	79	JR-20117-10	52	JR-2FR8-5	84
JR-071123-10	62	JR-20124-10	52	JR-2FR8TI-5	86



## Part Number Index, cont.

Part. No.	Page	Part. No.	Page	Part. No.	Page
JR-32171	68	JR-55024-5	28	JR-6032	58
JR-32172	68	JR-5502-5	29	JR-6033	58
JR-32174	68	JR-55025-5	28	JR-6034	58
JR-32178	68	JR-5503-5	30	JR-6039	58
JR-32179	68	JR-5504-5	29	JR-6040	58
JR-34N	45	JR-5505	36	JR-6041	58
JR-357090	21	JR-55050-10	49	JR-6042	58
JR-357180	21	JR-55051-10	49	JR-6043	58
JR-35P	37	JR-55060-10	50	JR-6044	58
JR-35S	44	JR-55061-10	50	JR-60-5	41
JR-367008-10	69	JR-55070-10	50	JR-6120-50030	79
JR-367008-2	69	JR-55071-10	50	JR-6120-50031	79
JR-367008-20	69	JR-5507-5	29	JR-6120-50032	79
JR-367016-10	69	JR-55080-10	51	JR-6120-50033	79
JR-367016-2	69	JR-55080-5	27	JR-6140	18
JR-367016-20	69	JR-55081-10	51	JR-6141-10	18
JR-3675-2	69	JR-55081-5	27	JR-6145-3062-5	83
JR-3678-2	69	JR-55082-10	51	JR-6145-3125-5	83
JR-3678-25	69	JR-55082-5	27	JR-6145-3500-5	83
JR-39021	103	JR-55083-10	51	JR-6145-6188-5	83
JR-39022	103	JR-55083-5	27	JR-6146-3250-5	83
JR-39023	103	JR-55084-10	51	JR-6146-3500-5	83
JR-39024	103	JR-55084-5	27	JR-65001	90
JR-4071-10	93	JR-5508-5	29	JR-65002	90
JR-4072-10	93	JR-55085-10	51	JR-65003	90
JR-4073-10	93	JR-55085-5	27	JR-65005	90
JR-4074-10	93	JR-55086-10	51	JR-65006	90
JR-4075-10	93	JR-55086-5	27	JR-65011	91
JR-408	37	JR-55087-10	51	JR-66140	91
JR-409	37	JR-55087-5	27	JR-66158	91
JR-410	37	JR-55088-10	51	JR-66175	90
JR-4676-10	67	JR-55088-5	27	JR-66176	90
JR-4676-10TF	67	JR-55089-10	51	JR-66177	90
JR-4676-2	67	JR-55089-5	27	JR-66178	90
JR-4676-2.5TF	67	JR-5508X-12	27	JR-66179	90
JR-4676-5TF	67	JR-5508X-24	27	JR-66180	90
JR-4677-10	67	JR-55090-10	51	JR-67140	91
JR-4677-10TF	67	JR-55090-5	27	JR-67140-1A	91
JR-4677-2	67	JR-55091-10	51	JR-67140-2B	91
JR-4677-2.5TF	67	JR-55091-5	27	JR-67140-3A	91
JR-4677-5TF	67	JR-55100-5	30	JR-67152-5	83
JR-5003-5	31	JR-5510-5	30	JR-67190	91
JR-5004-5	30	JR-55110-5	30	JR-67202	90
JR-5005-5	96	JR-5511-5	30	JR-67204	90
JR-55010-5	32	JR-5570-5	31	JR-67206	90
JR-55011-5	32	JR-5580-5	30	JR-67208	90
JR-55012-5	32	JR-5595	36	JR-67220	90
JR-55013-5	32	JR-56-5	41	JR-67222	90
JR-55014-5	32	JR-57-5	41	JR-67-250-2	91
JR-55015-5	32	JR-58000-5	27	JR-67-250-4	91
JR-55016-12	32	JR-58-5	41	JR-67-600	91
JR-55020-5	28	JR-59-5	41	JR-68152-5	83
JR-55021-5	28	JR-6029	58	JR-68181	90
JR-55022-5	28	JR-6030	58	JR-68182	90
JR-55023-5	28	JR-6031	58	JR-68183	90



## Part Number Index, cont.

Part. No.	Page	Part. No.	Page	Part. No.	Page
JR-68184	90	JR-9000-0602	66	JR-C-NYXFPK	39
JR-68185	90	JR-9000-0603	66	JR-CU4LPK	56
JR-68186	90	JR-9000-0604	66	JR-EN1FPKB	38
JR-68230-05	74	JR-9000-0640	66	JR-FRK1	47
JR-68230-2	74	JR-9000-0641	66	JR-HKS	47
JR-68230-5	74	JR-9000-0642	66	JR-LZN1-5	40
JR-68231-05	74	JR-9000-0665	35	JR-MSV2S	60
JR-68231-2	74	JR-9001	21	JR-MSV4L	60
JR-68231-5	74	JR-9001-0620	80	JR-MSV4P	60
JR-68242	92	JR-9001-0621	80	JR-MSV4T	60
JR-68247	71	JR-B-Z6M1PK	35	JR-MSV5D	60
JR-68250	71	JR-C5M1PK	58	JR-MZN1-5	40
JR-68251	71	JR-C5M2PK	58	JR-PTAS	61
JR-68252	92	JR-C9M1PK	58	JR-PV	47
JR-68253	71	JR-C9M2PK	58	JR-SLC-85-10	98
JR-68258	72	JR-CBULPK	57	JR-SLC-85-100	98
JR-68260	72	JR-CBUMPK	57	JR-SLC-85-1000	98
JR-68261	72	JR-CCPK-1	37	JR-SLC-85-20	98
JR-68262	72	JR-C-EN.5FPKB	38	JR-SLC-85-250	98
JR-792	23	JR-CEN1PK	55	JR-SLC-85-5	98
JR-793	23	JR-CEN2PK	55	JR-SLC-85-50	98
JR-794	22	JR-CFE-S10-5	70	JR-SLC-85-500	98
JR-795	22	JR-CFE-S2-5	70	JR-SLC-86-10	99
JR-796	23	JR-CFE-S75-5	70	JR-SLC-86-100	99
JR-797	22	JR-CFL-4D	56	JR-SLC-86-1000	99
JR-798	22	JR-CFL-CB1KF	49	JR-SLC-86-20	99
JR-800	46	JR-CFL-CB2KF	49	JR-SLC-86-250	99
JR-8000-0485	66	JR-CFL-CB4KF-5	56	JR-SLC-86-5	99
JR-802	46	JR-CM1XPK	59	JR-SLC-86-50	99
JR-804	46	JR-CM2XPK	59	JR-SLC-86-500	99
JR-805	46	JR-C-NERU1FPK	38	JR-SLR-111	101
JR-806	47	JR-C-NEU.5FPK	38	JR-SLR-111-77	101
JR-9	92	JR-C-NEU.5XFPK	38	JR-SLR-112	101
JR-9000-0001	78	JR-C-NL.15L-5	39	JR-SLR-112-77	101
JR-9000-0002	78	JR-C-NL.15S-5	39	JR-SLR-113	101
JR-9000-0003	78	JR-C-NL.20L-5	39	JR-SLR-113-77	101
JR-9000-0004	78	JR-C-NL.20S-5	39	JR-SLR-114	101
JR-9000-0005	77	JR-C-NL.25L-5	39	JR-SLR-114-77	101
JR-9000-0006	78	JR-C-NL.25S-5	39	JR-SLR-115	101
JR-9000-0007	78	JR-C-NL.30L-5	39	JR-SLR-115-77	101
JR-9000-0010	77	JR-C-NL.30S-5	39	JR-SLR-116	101
JR-9000-0460-5	83	JR-C-NL.35L-5	39	JR-SLR-116-77	101
JR-9000-0477-5	83	JR-C-NL.35S-5	39	JR-SLR-117	101
JR-9000-0520	65	JR-C-NLS1.15	39	JR-SLR-117-77	101
JR-9000-0520F	65	JR-C-NLS1.20	39	JR-SLR-118	101
JR-9000-0521	65	JR-C-NLS1.25	39	JR-SLR-118-77	101
JR-9000-0521F	65	JR-C-NLS1.30	39	JR-SLR-2190	100
JR-9000-0522	65	JR-C-NLS1.35	39	JR-SLR-2191	100
JR-9000-0522F	65	JR-C-NNFFPK	38	JR-SLR-2192	100
JR-9000-0530	65	JR-C-NPFPK	38	JR-SLR-2193	100
JR-9000-0530F	65	JR-C-NTFPK	39	JR-SLR-2194	100
JR-9000-0531	65	JR-C-NTXFPK	39	JR-SLR-2195	100
JR-9000-0531F	65	JR-C-NXFPK	39	JR-SLR-2196	100
JR-9000-0532	65	JR-C-NXXFPK	39	JR-SLR-2200	100
JR-9000-0532F	65	JR-C-NYFPK	39	JR-SLR-2201	100



## Part Number Index, cont.

Part. No.	Page	Part. No.	Page	Part. No.	Page
JR-SLR-2202	100	JR-T-4183-**	15	JR-T-97065	13
JR-SLR-2204	100	JR-T-4183C-**	16	JR-T-97075	13
JR-SLR-2206	100	JR-T-5993-**	7	JR-T-97080	13
JR-SLR-2207	100	JR-T-5994-**	7	JR-T-97085	13
JR-SLR-2210	100	JR-T-5995-**	7	JR-T-97090	13
JR-SLR-2222	100	JR-T-5999-**	6	JR-T-97095	13
JR-SLV-75-10	96	JR-T-6000-**	6	JR-T-97100	13
JR-SLV-75-100	96	JR-T-6001-**	6	JR-T-97105	13
JR-SLV-75-1000	96	JR-T-6002-**	6	JR-T-97110	13
JR-SLV-75-20	96	JR-T-6003-**	6	JR-TEFNI.502-M1	12
JR-SLV-75-250	96	JR-T-60031-**	6	JR-TEFNI.504-M1	12
JR-SLV-75-5	96	JR-T-60032-**	6	JR-TEFNI.505-M1	12
JR-SLV-75-50	96	JR-T-6004-**	7	JR-TEFNI.510-**	12
JR-SLV-75-500	96	JR-T-60041-**	7	JR-TEFNI.515-**	12
JR-SLV-76-10	97	JR-T-60042-**	7	JR-TEFNI.520-**	12
JR-SLV-76-100	97	JR-T-6006	7	JR-TEFNI130-**	12
JR-SLV-76-1000	97	JR-T-6007-**	6	JR-TEFNI140-**	12
JR-SLV-76-20	97	JR-T-6008-**	6	JR-TP-5998-**	5
JR-SLV-76-250	97	JR-T-6009-**	6	JR-TP-5999-**	5
JR-SLV-76-5	97	JR-T-6010-**	6	JR-TP-6000-**	5
JR-SLV-76-50	97	JR-T-6011-**	6	JR-TP-6001-**	5
JR-SLV-76-500	97	JR-T-6060-**	7	JR-TS-16110	9
JR-T-078-**	17	JR-T-6130-**	18	JR-TS-16115	9
JR-T-080-**	17	JR-T-625-04-**	11	JR-TS-16120	9
JR-T-082-**	17	JR-T-625-05-**	11	JR-TS-161210	9
JR-T-083-**	17	JR-T-625-10-**	11	JR-TS-161215	9
JR-T-084-**	17	JR-T-625-20-**	11	JR-TS-161225	9
JR-T-085-**	17	JR-T-625-30-**	11	JR-TS-16125	9
JR-T-086-**	17	JR-T-625-40-**	11	JR-TS-161250	9
JR-T-3000-123	20	JR-T-626-00-**	11	JR-TS-1615	9
JR-T-3000-125	20	JR-T-628-00	11	JR-TS-16150	9
JR-T-3000-128	20	JR-T-6800-**	15	JR-TS-16210	9
JR-T-3000-133	20	JR-T-6801-**	15	JR-TS-16215	9
JR-T-3000-135	20	JR-T-6802-**	17	JR-TS-16220	9
JR-T-3000-138	20	JR-T-6803-**	17	JR-TS-1625	9
JR-T-3000-223	20	JR-T-6804-**	17	JR-TS-16250	9
JR-T-3000-225	20	JR-T-6805-**	15	JR-TS-16310	9
JR-T-3000-228	20	JR-T-6806-**	17	JR-TS-16315	9
JR-T-3000-233	20	JR-T-6807-**	15	JR-TS-16320	9
JR-T-3000-235	20	JR-T-6808-**	16	JR-TS-1635	9
JR-T-3000-238	20	JR-T-6809	17	JR-TS-16350	9
JR-T-3000-423	20	JR-T-6810	15	JR-TS-16410	9
JR-T-3000-425	20	JR-T-6811-**	16	JR-TS-16415	9
JR-T-3000-428	20	JR-T-6812-**	17	JR-TS-16420	9
JR-T-4001-**	16	JR-T-97005	13	JR-TS-1645	9
JR-T-4002-**	16	JR-T-97010	13	JR-TS-16450	9
JR-T-4003-**	16	JR-T-97015	13	JR-TS-16610	9
JR-T-4007-**	16	JR-T-97020	13	JR-TS-16615	9
JR-T-4011-**	15	JR-T-97025	13	JR-TS-16620	9
JR-T-4011C-**	16	JR-T-97030	13	JR-TS-1665	9
JR-T-4036-**	15	JR-T-97035	13	JR-TS-16650	9
JR-T-4036C-**	16	JR-T-97040	13	JR-TS-16710	9
JR-T-4037-**	15	JR-T-97045	13	JR-TS-16715	9
JR-T-4039	15	JR-T-97055	13	JR-TS-16720	9
JR-T-4041	15	JR-T-97060	13	JR-TS-1675	9



### Part Number Index, cont.

Part. No.	Page	Part. No.	Page	Part. No.	Page
JR-TS-16750	9	JR-TS-3265	8	JR-ZGF1PK	38
JR-TS-16810	9	JR-TS-32650	8	JR-ZN.5-5	40
JR-TS-16815	9	JR-TSS.505-**	10	JR-ZN1-5	40
JR-TS-16820	9	JR-TSS.507-**	10	JR-ZN2-5	40
JR-TS-1685	9	JR-TSS.510-**	10	JR-ZN2PK-5	31
JR-TS-16850	9	JR-TSS.520-**	10	JR-ZP1	44
JR-TS-32110	8	JR-TSS105-**	10	JR-ZP2	44
JR-TS-32115	8	JR-TSS110-**	10	JR-ZRU21	42
JR-TS-32125	8	JR-TSS115-**	10	JR-ZRU21FPK	33
JR-TS-3215	8	JR-TSS120-**	10	JR-ZT1	43
JR-TS-32150	8	JR-TSS130-**	10	JR-ZT1C	43
JR-TS-32210	8	JR-TSS140-**	10	JR-ZT1M	43
JR-TS-32215	8	JR-TSS230-**	10	JR-ZU1CS6	42
JR-TS-32225	8	JR-TSS240-**	10	JR-ZU1MS6	42
JR-TS-3225	8	JR-TSS260-**	10	JR-ZU1S6	42
JR-TS-32250	8	JR-TSS267-**	10	JR-ZU1XCS6	42
JR-TS-32310	8	JR-TSS285-**	10	JR-ZX1	43
JR-TS-32315	8	JR-Z6M1PK	35	JR-ZX1C	43
JR-TS-32325	8	JR-ZBU1CFPK	33	JR-ZX1M	43
JR-TS-3235	8	JR-ZBU1FPK	33	PS-614158	104
JR-TS-32350	8	JR-ZBU1MFPK	33	PS-614160	104
JR-TS-32410	8	JR-ZC1	44	PS-614161	104
JR-TS-32415	8	JR-ZC1PK	37	PS-614163	104
JR-TS-32425	8	JR-ZC2	44	PS-614170	104
JR-TS-3245	8	JR-ZF.5S6-5	40	PS-614250	104
JR-TS-32450	8	JR-ZF1S6-5	40	PS-644350	104
JR-TS-32610	8	JR-ZF2PK-5	31	PS-644850	104
JR-TS-32615	8	JR-ZF2S6-5	40		
JR-TS-32625	8	JR-ZGF.5PK	38		

### Trademarks and Registered Trademarks

Cheminert® is a Registered Trademark of Valco Instruments Co. Inc. and VICI International

Delrin® is a Registered Trademark of E.I. du Pont de Nemours and Company

Duran® is a Registered Trademark of Schott AG

Fortron® is a Registered Trademark of Ticona

Kalrez® is a Registered Trademark of DuPont Dow Elastomers

Kel-F® is a Registered Trademark of 3M Company

Kynar® is a Registered Trademark of Elf Atochem North America, Inc.

Nanovolume® is a Registered Trademark of Valco Instruments Co. Inc.

PEEK™ is a Trademark of Victrex plc

PEEKclad-FST™ is a Trademark of Valco Instruments Co. Inc. and VICI International

PEEKsil™ is a Trademark of SGE International Pty. Ltd.

Rheodyne® is a Registered Trademark of Rheodyne LLC

Ryton® is a Registered Trademark of Chevron Phillips Chemical Company

Techtron® is a Registered Trademark of Quadrant EPP

Teflon® is a Registered Trademark of E.I. du Pont de Nemours and Company

Tefzel® is a Registered Trademark of E.I. du Pont de Nemours and Company

VICI® is a Registered Trademark of Valco Instruments Co. Inc. and VICI International