

Filters



Filters with a removable screen

These filters come with a removable 2μ screen. The standard screen can be replaced with any screen of the proper diameter, *but not by a frit*. These filters are suitable for streams with frequent contamination, since the filtering element is easily changed. Standard material is Type 316 series stainless.

Patent Numbers 4,281,679 and 4,173,363

P

Description	Bore	Standard		Bulkhead	
		Prod No	Price	Prod No	Price
1/32" to 1/32"	0.25 mm	ZUFR.5		ZBUFR.5	
1/16" to 1/32"	0.25 mm	ZRUFR1.5		ZBRUFR1.5	
1/16" to 1/16"	0.25 mm	ZUFR1C		ZBUFR1C	
	0.50 mm	ZUFR1		ZBUFR1	
1/8" to 1/16"	0.75 mm	ZRUFR21		ZBRUFR21	
1/8" to 1/8"	2.00 mm	ZUFR2		ZBUFR2	
1/4" to 1/16"	1.00 mm	ZRUFR41		ZBRUFR41	
1/4" to 1/8"	2.00 mm	ZRUFR42		ZBRUFR42	



TECH TIP

Should you use a filter with a frit or one with a screen?

Screens have much higher flow capacity (Cv), but frits are the best choice for maximum filtration or if your application requires Hastelloy C or titanium. However, since they are thicker than screens, frits allow more mixing, and the downside of their superior filtration is that they clog more often than screens.

Note! The difference in thickness also means that frits and screens cannot be used interchangeably in the same fitting body:

A frit must always be replaced with a frit.

A screen must always be replaced with a screen.

Replacement screens..... page 53

0.25 mm	=	.010"
0.50 mm	=	.020"
0.75 mm	=	.030"
1.0 mm	=	.040"
1.5 mm	=	.060"
2.0 mm	=	.080"
4.6 mm	=	.180"
6.0 mm	=	.236"
6.4 mm	=	.253"
7.0 mm	=	.275"
10.0 mm	=	.400"
27.0 mm	=	1.08"
1/32"	=	0.8 mm
1/16"	=	1.6 mm
1/8"	=	3.2 mm
1/4"	=	6.4 mm
3/8"	=	9.5 mm
1/2"	=	12.7 mm

5/16"	=	.312"	=	7.9 mm
3/8"	=	.375"	=	9.5 mm
7/16"	=	.437"	=	11.1 mm