

.....

Solenoid Operated Micro-Pumps



130SP Series Micro-Pump



INDEX

Page 3 **Micro-Pump Selection Guide**

Page 4 **030SP Series Micro-Pump**

Ported Micro-Pumps (1/4"-28 UNF) for a very small dispense volume (4µl)

Page 6 **120SP Series Micro-Pump**

Ported Micro-Pumps (1/4"-28 UNF) for precise dispense volumes from 10 to 60µl

Page 8 **130SP Series Micro-Pump**

Ported Micro-Pumps (1/4"-28 UNF) for precise dispense volumes from 10 to 60µl (inert body)

Page 10 **150SP Series Micro-Pump**

Ported Micro-Pumps (5/16"-24 UNF) for precise dispense volumes from 100 to 250µl

Page 12 **039SP Series Micro-Pump**

Manifold mounted Micro-Pumps for a very small dispense volume (4µl)


Page 14 **139SP Series Micro-Pump**

Manifold mounted Micro-Pumps for precise dispense volumes from 10 to 60µl

Page 16 **Manifolds, Controller and Tubing**

Page 17 **Fittings**

Page 18 **Micro-Pumps Tech Tips - operation and installation**




MICRO-PUMPS GENERAL INFORMATION

What is a Micro-Pump?

A Micro-Pump is a solenoid operated device designed to provide a precise, repeatable and discrete dispensed volume of fluid. The

flow path is isolated from the operating mechanism by a flexible diaphragm. When the solenoid is energized, the diaphragm is retracted creating a partial vacuum within the pump body. This pulls liquid through the inlet check valve (A) and simultaneously closes the outlet check valve (B). When the

solenoid is de-energized a spring pushes the diaphragm down, expelling a discrete volume of liquid through check valve B while simultaneously closing check valve A. Micro-Pumps require a complete on-off cycle for each discrete dispense. Repeatedly cycling the solenoid creates a pulsed flow (refer to "Accurate discrete dispense volumes" in next column).



Features of the Bio-Chem Valve™ Micro-Pump

Inert materials

Our pumps provide a non-metallic inert fluid path for the dispensing of high purity or aggressive fluids. There is a range of different materials available for all the wetted parts of the pumps - body, diaphragm and check valve. Material combinations can be chosen to suit the application (refer to individual product selection pages for standard combinations - custom combinations are available, refer to page 18).

Body materials: PPS, PTFE, PEEK™, POM

Diaphragm materials: EPDM, PTFE

Check valve materials: EPDM, FKM, FFKM

Self-priming

At start-up, the pump is able to draw air. The suction created by the larger pumps is sufficient to pull liquids from an unpressurized container located up to 4' 3" (1.3m) beneath the pump. Once the pump is primed, it is able to generate around 5psi (0.3bar) pressure, equating to 11' 6" (3.5m) of water.

Continuous duty

The pumps are capable of continuous duty. They are suitable for up to 20 million actuations, corresponding to nearly 3,000 hours of continuous use at a 2 Hz cycle rate.

Accurate discrete dispense volumes

Dispense volumes range from 4µl to 250µl per cycle. The pumps can be cycled at up to 4 Hz for the smallest version and 1.6 Hz for the largest. Pumps can be operated at less than the maximum cycle rate by increasing the length of the "off" time. The "on" time should remain unchanged to retain dispense accuracy.

Micro-Pump Selection Guide

1. Select pump style; either Ported or Manifold mount and work from the appropriate table:
 - Ported for direct connection with 1/4"-28 fittings (5/16"-24 for 150SP)
 - Manifold mount for use with manifolds (see page 16)

Then:

2. Locate the volumetric characteristics that best suit your needs
3. Choose your preferred body material depending on the level of chemical inertness you require
4. Turn to the pages indicated to see full details and ordering information for each pump.

Ported	Volumetric output		Body Material			
	Discrete Dispense Vol (μl)	Max flow rate (ml/min)	PTFE	PPS	PEEK™	POM
	4	0.96		030SP (pg. 4)		
	10	1.2				
	20	2.4				
	30	3.6	130SP (pg. 8)	120SP (pg. 6)	120SP (pg. 6)	130SP (pg. 8)
	40	4.8				
	50	6.0				
	60	7.2				
	100	9.6				
	125	12.0				
	150	14.4				
	175	16.8		150SP (pg. 10)	150SP (pg. 10)	
	200	19.2				
	225	21.6				
	250	24.0				

Manifold mounted	Volumetric output		Body Material			
	Discrete Dispense Vol (μl)	Max flow rate (ml/min)	PTFE	PPS	PEEK™	POM
	4	0.96		039SP (pg. 12)		
	10	1.2				
	20	2.4				
	30	3.6	139SP (pg. 14)		139SP (pg. 14)	139SP (pg. 14)
	40	4.8				
	50	6.0				
	60	7.2				

Polymers referenced in this brochure:

EPDM = ethylene-propylene-diene
 ETFE = ethylene tetrafluoroethylene
 FEP = fluorinated ethylene propylene
 FKM = fluorinated elastomer
 FFKM = perfluoro elastomer
 PEEK™ = polyetheretherketone
 POM = polyoxymethylene (Acetal resin)
 PPS = polyphenylene sulfide
 PTFE = polytetrafluoroethylene.

030SP SERIES MICRO-PUMP

For precise dispensing of 4µl and flow rates
up to 0.96 ml/min

- Self-priming
- 4µl discrete dispense volume
- 960µl/min maximum flow rate
- 1/4"-28 UNF threaded ports

The 030SP series Micro-Pumps are solenoid operated, with the operating mechanism isolated from the flow path by a diaphragm. Check valves situated at the inlet and outlet of the pump control the direction of flow.

Materials available for the wetted parts are:

- Body materials: PPS
- Diaphragm materials: PTFE
- Check valve materials: FKM

030SP series options

PART NO.	VDC	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 4µl dispense					
030SP124-4TV	12	4	PPS	PTFE	FKM
24 VDC; 4µl dispense					
030SP244-4TV	24	4	PPS	PTFE	FKM

ARRANGEMENT




SPECIFICATIONS

030SP Fluid Data	
Dispense Volume (µl)	4
Set-point accuracy	+/- 25%
Repeatability	+/- 5%
Max flow rate (µl/min)	960
Internal vol (µl)	130

030SP Electrical Data			
Voltage	Power @70°F (21°C)	Current @70°F (21°C)	Effective continuous power @ max cycle rate
12 VDC	1.9 Watts	0.22 amps	0.9 Watts
24 VDC	1.9 Watts	0.11 amps	0.9 Watts

030SP Cycle Rates		
Min "on" time	Min "off" time	Max cycle rate
125 msec	125 msec	4.0 Hz


030SP Micro-Pumps can be cycled at up to 4 Hz. To maintain pumping precision the voltage "on" time should remain fixed - the pumping rate can be changed by increasing the "off" time.



Recommended tubing for 030SP

Inlet & outlet, 1/32" (0.80mm) ID, hardwall tubing,
PART NO. 008T16-080

INSTALLATION DRAWING



120SP SERIES MICRO-PUMP

For precise dispensing between 10 and 60µl and flow rates up to 7.2 ml/min

- Self-priming
- 10-60µl discrete dispense volumes
- Up to 7.2 ml/min maximum flow rate
- 1/4"-28 UNF threaded ports

The 120SP series Micro-Pumps are solenoid operated, with the operating mechanism isolated from the flow path by a diaphragm. Check valves situated at the inlet and outlet of the pump control the direction of flow. The combination of materials for each component can be selected to best suit your specific application.

Materials available for the wetted parts are:


- Body materials: PPS, PEEK™
- Diaphragm materials: PTFE, EPDM
- Check valve materials: EPDM, FKM, FFKM

120SP series options

NOTE: For 24 VDC, replace 120SP12 with 120SP24 in any of the part numbers listed.


PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 10µl dispense (Note: PTFE diaphragm for all 10 µl options)				
120SP1210-4TE	10	PPS	PTFE	EPDM
120SP1210-4TV	10	PPS	PTFE	FKM
120SP1210-4TP	10	PPS	PTFE	FFKM
120SP1210-5TE	10	PEEK™	PTFE	EPDM
120SP1210-5TV	10	PEEK™	PTFE	FKM
120SP1210-5TP	10	PEEK™	PTFE	FFKM
12 VDC; 20µl dispense				
120SP1220-4EE	20	PPS	EPDM	EPDM
120SP1220-4TV	20	PPS	PTFE	FKM
120SP1220-4TP	20	PPS	PTFE	FFKM
120SP1220-5EE	20	PEEK™	EPDM	EPDM
120SP1220-5TV	20	PEEK™	PTFE	FKM
120SP1220-5TP	20	PEEK™	PTFE	FFKM
12 VDC; 30µl dispense				
120SP1230-4EE	30	PPS	EPDM	EPDM
120SP1230-4TV	30	PPS	PTFE	FKM
120SP1230-4TP	30	PPS	PTFE	FFKM
120SP1230-5EE	30	PEEK™	EPDM	EPDM
120SP1230-5TV	30	PEEK™	PTFE	FKM
120SP1230-5TP	30	PEEK™	PTFE	FFKM

ARRANGEMENT



PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 40µl dispense				
120SP1240-4EE	40	PPS	EPDM	EPDM
120SP1240-4TV	40	PPS	PTFE	FKM
120SP1240-4TP	40	PPS	PTFE	FFKM
120SP1240-5EE	40	PEEK™	EPDM	EPDM
120SP1240-5TV	40	PEEK™	PTFE	FKM
120SP1240-5TP	40	PEEK™	PTFE	FFKM
12 VDC; 50µl dispense				
120SP1250-4EE	50	PPS	EPDM	EPDM
120SP1250-4TV	50	PPS	PTFE	FKM
120SP1250-4TP	50	PPS	PTFE	FFKM
120SP1250-5EE	50	PEEK™	EPDM	EPDM
120SP1250-5TV	50	PEEK™	PTFE	FKM
120SP1250-5TP	50	PEEK™	PTFE	FFKM
12 VDC; 60µl dispense (Note: EPDM diaphragm for all 60 µl options)				
120SP1260-4EE	60	PPS	EPDM	EPDM
120SP1260-5EE	60	PEEK™	EPDM	EPDM

INSTALLATION DRAWING




SPECIFICATIONS

120SP Fluid Data						
Dispense Volume (µl)	10	20	30	40	50	60
Set-point accuracy	+/- 4%	+/- 4%	+/- 3%	+/- 3%	+/- 2%	+/- 2%
Repeatability	+/- 3%	+/- 3%	+/- 3%	+/- 2%	+/- 2%	+/- 2%
Max flow rate (µl/min)	1200	2400	3600	4800	6000	7200
Internal vol (µl)	105	105	105	105	105	105

120SP Electrical Data				120SP Cycle Rates		
Voltage	Power @70°F (21°C)	Current @70°F (21°C)	Effective continuous power @ max cycle rate	Min "on" time	Min "off" time	Max cycle rate
12 VDC	4.0 Watts	0.32 amps	1.2 Watts	150 msec	350 msec	2.0 Hz
24 VDC	4.0 Watts	0.16 amps	1.2 Watts			

Recommended tubing for 120SP
Inlet & outlet, 1/32" (0.80mm) ID, hardwall tubing, PART NO. 008T16-080

120SP Micro-Pumps can be cycled at up to 2 Hz. To maintain pumping precision the voltage "on" time should remain fixed - the pumping rate can be changed by increasing the "off" time.



130SP SERIES MICRO-PUMP

For precise dispensing between 10 and 60µl and flow rates up to 7.2 ml/min

- Self-priming
- 10-60µl discrete dispense volumes
- Up to 7.2 ml/min maximum flow rate
- 1/4"-28 UNF threaded ports
- Most inert body material for harshest applications

The 130SP series Micro-Pumps are solenoid operated, with the operating mechanism isolated from the flow path by a diaphragm. Check valves situated at the inlet and outlet of the pump control the direction of flow. The combination of materials for each component can be selected to best suit your specific application.

Materials available for the wetted parts are:


- Body materials: PTFE, POM
- Diaphragm materials: PTFE, EPDM
- Check valve materials: EPDM, FKM, FFKM

130SP series options


NOTE: For 24 VDC, replace 130SP12 with 130SP24 in any of the part numbers listed.

PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 10µl dispense (Note: PTFE diaphragm for all 10 µl options)				
130SP1210-1TP	10	PTFE	PTFE	FFKM
130SP1210-6TV	10	POM	PTFE	FKM
130SP1210-6TE	10	POM	PTFE	EPDM
12 VDC; 20µl dispense				
130SP1220-1TP	20	PTFE	PTFE	FFKM
130SP1220-6TV	20	POM	PTFE	FKM
130SP1220-6EE	20	POM	EPDM	EPDM
12 VDC; 30µl dispense				
130SP1230-1TP	30	PTFE	PTFE	FFKM
130SP1230-6TV	30	POM	PTFE	FKM
130SP1230-6EE	30	POM	EPDM	EPDM
12 VDC; 40µl dispense				
130SP1240-1TP	40	PTFE	PTFE	FFKM
130SP1240-6TV	40	POM	PTFE	FKM
130SP1240-6EE	40	POM	EPDM	EPDM
12 VDC; 50µl dispense				
130SP1250-1TP	50	PTFE	PTFE	FFKM
130SP1250-6TV	50	POM	PTFE	FKM
130SP1250-6EE	50	POM	EPDM	EPDM
12 VDC; 60µl dispense				
130SP1260-6EE	60	POM	EPDM	EPDM

ARRANGEMENT



INSTALLATION DRAWING




SPECIFICATIONS

130SP Volumetric Data						
Dispense Volume (µl)	10	20	30	40	50	60
Set-point accuracy	+/- 4%	+/- 4%	+/- 3%	+/- 3%	+/- 2%	+/- 2%
Repeatability	+/- 3%	+/- 3%	+/- 3%	+/- 2%	+/- 2%	+/- 2%
Max flow rate (µl/min)	1200	2400	3600	4800	6000	7200
Internal vol (µl)	105	105	105	105	105	105

130SP Electrical Data				130SP Cycle Rates		
Voltage	Power @70°F (21°C)	Current @70°F (21°C)	Effective continuous power @ max cycle rate	Min "on" time	Min "off" time	Max cycle rate
12 VDC	4.0 Watts	0.32 amps	1.2 Watts	150 msec	350 msec	2.0 Hz
24 VDC	4.0 Watts	0.16 amps	1.2 Watts			

Recommended tubing for 130SP
Inlet & outlet, 1/32" (0.80mm) ID, hardwall tubing, PART NO. 008T16-080

130SP Micro-Pumps can be cycled at up to 2 Hz. To maintain pumping precision the voltage "on" time should remain fixed - the pumping rate can be changed by increasing the "off" time.



150SP SERIES MICRO-PUMP

For precise dispensing between 100 and 250µl and flow rates up to 24 ml/min

- Self-priming
- 100-250µl discrete dispense volumes
- Up to 24 ml/min maximum flow rate
- 5/16"-24 UNF threaded ports

The 150SP series Micro-Pumps are solenoid operated, with the operating mechanism isolated from the flow path by a diaphragm. Check valves situated at the inlet and outlet of the pump control the direction of flow. The combination of materials for each component can be selected to best suit your specific application.

Materials available for the wetted parts are:

- Body materials: PPS, PEEK™
- Diaphragm materials: EPDM
- Check valve materials: EPDM

150SP series options


NOTE: For 24 VDC, replace 150SP12 with 150SP24 in any of the part numbers listed.

PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 100µl dispense				
150SP12100-4EE	100	PPS	EPDM	EPDM
150SP12100-5EE	100	PEEK™	EPDM	EPDM
12 VDC; 125µl dispense				
150SP12125-4EE	125	PPS	EPDM	EPDM
150SP12125-5EE	125	PEEK™	EPDM	EPDM
12 VDC; 150µl dispense				
150SP12150-4EE	150	PPS	EPDM	EPDM
150SP12150-5EE	150	PEEK™	EPDM	EPDM
12 VDC; 175µl dispense				
150SP12175-4EE	175	PPS	EPDM	EPDM
150SP12175-5EE	175	PEEK™	EPDM	EPDM
12 VDC; 200µl dispense				
150SP12200-4EE	200	PPS	EPDM	EPDM
150SP12200-5EE	200	PEEK™	EPDM	EPDM
12 VDC; 225µl dispense				
150SP12225-4EE	225	PPS	EPDM	EPDM
150SP12225-5EE	225	PEEK™	EPDM	EPDM
12 VDC; 250µl dispense				
150SP12250-4EE	250	PPS	EPDM	EPDM
150SP12250-5EE	250	PEEK™	EPDM	EPDM

ARRANGEMENT



INSTALLATION DRAWING




SPECIFICATIONS

150SP Fluid Data							
Dispense Volume (µl)	100	125	150	175	200	225	250
Set-point accuracy	+/- 5%	+/- 5%	+/- 4%	+/- 4%	+/- 4%	+/- 3%	+/- 3%
Repeatability	+/- 1%	+/- 1%	+/- 1%	+/- 1%	+/- 0.5%	+/- 0.5%	+/- 0.5%
Max flow rate (µl/min)	9600	12000	14400	16800	19200	21600	24000
Internal vol (µl)	710	710	710	710	710	710	710

150SP Electrical Data				150SP Cycle Rates		
Voltage	Power @70°F (21°C)	Current @70°F (21°C)	Effective continuous power @ max cycle rate	Min "on" time	Min "off" time	Max cycle rate
12 VDC	8.0 Watts	0.66 amps	3.2 Watts	200 msec	400 msec	1.6 Hz
24 VDC	8.0 Watts	0.33 amps	3.2 Watts			

Recommended tubing for 150SP	
Inlet & outlet, 1/8" (3.2mm) ID, hardwall tubing, PART NUMBER 008T47-032	

150SP Micro-Pumps can be cycled at up to 1.6 Hz. To maintain pumping precision the voltage "on" time should remain fixed - the pumping rate can be changed by increasing the "off" time.



039SP SERIES MICRO-PUMP

For precise dispensing of 4µl and flow rates up to 0.96 ml/min in a manifold mountable design

- Self-priming
- 4µl discrete dispense volume
- 960µl/min maximum flow rate
- Manifold mountable

This sibling to the 030SP Micro-Pump duplicates the performance characteristics but is supplied ready for mounting in your manifold. Please contact us if you would like us to supply the manifold (see page 16).


Materials available for the wetted parts of the pump are:

- Body materials: PPS
- Diaphragm materials: PTFE
- Check valve materials: FKM

039SP series options

PART NO.	VDC	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 4µl dispense					
039SP124-4TV	12	4	PPS	PTFE	FKM
24 VDC; 4µl dispense					
039SP244-4TV	24	4	PPS	PTFE	FKM

ARRANGEMENT




SPECIFICATIONS


039SP Volumetric Data	
Dispense Volume (µl)	4
Set-point accuracy	+/- 25%
Repeatability	+/- 5%
Max flow rate (µl/min)	960
Internal vol (µl)	130

039SP Electrical Data			
Voltage	Power @70°F (21°C)	Current @70°F (21°C)	Effective continuous power @ max cycle rate
12 VDC	1.9 Watts	0.22 amps	0.9 Watts
24 VDC	1.9 Watts	0.11 amps	0.9 Watts


039SP Cycle Rates		
Min "on" time	Min "off" time	Max cycle rate
125 msec	125 msec	4.0 Hz



INSTALLATION DRAWING



MANIFOLD INTERFACE DRAWING



139SP SERIES MICRO-PUMP

For precise dispensing between 10 and 60µl and flow rates up to 7.2 ml/min in a manifold mountable design

- Self-priming
- 10-60µl discrete dispense volumes
- Up to 7.2 ml/min maximum flow rate
- Manifold mountable

This sibling to the 130SP Micro-Pump duplicates the performance characteristics but is supplied ready for mounting in your manifold. Please contact us if you would like us to supply the manifold (see page 16). Materials available for the wetted parts are:


- Body materials: PTFE, POM, PEEK™
- Diaphragm materials: PTFE, EPDM
- Check valve materials: EPDM, FKM, FFKM

139SP series options

NOTE: For 24 VDC, replace 139SP12 with 139SP24 in any of the part numbers listed.

PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 10µl dispense (Note: PTFE diaphragm for all 10 µl options)				
139SP1210-1TP	10	PTFE	PTFE	FFKM
139SP1210-5TP	10	PEEK™	PTFE	FFKM
139SP1210-5TV	10	PEEK™	PTFE	FKM
139SP1210-5TE	10	PEEK™	PTFE	EPDM
139SP1210-6TV	10	POM	PTFE	FKM
139SP1210-6TE	10	POM	PTFE	EPDM
12 VDC; 20µl dispense				
139SP1220-1TP	20	PTFE	PTFE	FFKM
139SP1220-5TP	20	PEEK™	PTFE	FFKM
139SP1220-5TV	20	PEEK™	PTFE	FKM
139SP1220-5TE	20	PEEK™	PTFE	EPDM
139SP1220-6TV	20	POM	PTFE	FKM
139SP1220-6EE	20	POM	EPDM	EPDM
12 VDC; 30µl dispense				
139SP1230-1TP	30	PTFE	PTFE	FFKM
139SP1230-5TP	30	PEEK™	PTFE	FFKM
139SP1230-5TV	30	PEEK™	PTFE	FKM
139SP1230-5TE	30	PEEK™	PTFE	EPDM
139SP1230-6TV	30	POM	PTFE	FKM
139SP1230-6EE	30	POM	EPDM	EPDM

ARRANGEMENT




PART NO.	DISPENSE VOL (µL)	BODY MATERIAL	DIAPHRAGM MATERIAL	CHECK VALVE MATERIAL
12 VDC; 40µl dispense				
139SP1240-1TP	40	PTFE	PTFE	FFKM
139SP1240-5TP	40	PEEK™	PTFE	FFKM
139SP1240-5TV	40	PEEK™	PTFE	FKM
139SP1240-5TE	40	PEEK™	PTFE	EPDM
139SP1240-6TV	40	POM	PTFE	FKM
139SP1240-6EE	40	POM	EPDM	EPDM
12 VDC; 50µl dispense				
139SP1250-1TP	50	PTFE	PTFE	FFKM
139SP1250-5TP	50	PEEK™	PTFE	FFKM
139SP1250-5TV	50	PEEK™	PTFE	FKM
139SP1250-5TE	50	PEEK™	PTFE	EPDM
139SP1250-6TV	50	POM	PTFE	FKM
139SP1250-6EE	50	POM	EPDM	EPDM
12 VDC; 60µl dispense				
139SP1260-6EE	60	POM	EPDM	EPDM


SPECIFICATIONS

The 139SP has the same specifications as the 130SP (see page 7)

INSTALLATION DRAWING



MANIFOLD INTERFACE DRAWING



MANIFOLDS



Custom manifold for (1) 139SP Micro-Pump (shown) and (3) isolation valves (not shown). Blue lines indicate the fluid path; the red dots are ruby balls used as plugs.

Custom-built manifolds are used to organize multiple Micro-Pumps and other Fluid Control Devices such as Isolation Valves into an efficient, pre-assembled, space-saving module that is designed to meet your specific flow needs. Manifolds can range from simple blocks for two devices to complex shapes with intricate flow paths for many devices. Bio-Chem Fluidics has produced complex manifolds for as many as 84 Micro-Pumps on a single block.

Features:

- Reduction of internal equipment space requirements.
- Allows for the combining of valves, tubing, pumps and connectors into a single, pre-assembled component.
- Elimination of unsightly and unmanageable wiring and tubing.
- Helps to reduce inventory.
- Reduces production time and costs associated with testing, handling and assembling multiple components.
- Materials of construction to suit fluid characteristics including, but not limited to; PTFE, POM, PEEK™, acrylic and PPS.

Please contact your local Bio-Chem Fluidics facility to discuss your manifold requirements with one of our engineers.



Custom manifold for (2) 139SP Micro-Pumps (not shown).

FCD CONTROLLER

The Bio-Chem Fluidics' Fluid Control Device (FCD) Controller is designed to provide end-user programmed control signals to any combination of eight such devices, including the full Bio-Chem Valve™ range of Micro-Pumps and solenoid operated valves.

The controller uses intelligent part number recognition technology to take the guesswork out of programming – simply by entering the part number the controller will recognize the function of the FCD and will generate applicable control signals. For applications requiring more control, the FCD Controller has a user friendly symbol driven interface allowing for rapid programming of up to 999 steps. Built-in PC interface also allows for remote programming from any computer.

Features:

- Eight FCD's can be controlled, each can be programmed independently of the others.
- Supplied pre-loaded with control data for the complete range of Bio-Chem Valve™ valves and pumps. The standard USB interface ensures product data will be updated as needed.
- Stores up to 10 programs, each with up to 999 steps.
- LED's indicate the presence of a device and whether or not it is actuated.
- Adaptive power supply accepts either 115VAC or 230VAC
- RoHS compliant, CE marked.

FCD Controller to be released Summer 2009



CONTROLLER

PART NUMBER	DESCRIPTION	QTY
BBX8	8-position FCD controller	1 pk

TUBING

Inert PTFE tubing can be used with virtually all chemicals, solvents and corrosive materials, even at elevated temperatures. It can be sterilized in-line by steam, chemical methods, or autoclaving. This semi-rigid tubing is ideal for use with the Bio-Chem Valve™ range of Micro-Pumps.

PTFE TUBING				
PART NUMBER	OD	ID	LENGTH	QTY
008T16-080-20	1.6mm (1/16")	0.8mm (1/32")	20m	ea
008T16-080-200	1.6mm (1/16")	0.8mm (1/32")	200m	ea
008T47-032-10	4.7mm (3/16")	3.2mm (1/8")	10m	ea
008T47-032-100	4.7mm (3/16")	3.2mm (1/8")	100m	ea

Other sizes are available and we also stock flexible, Silicone tubing for other applications (for example, Pinch Valves). This PTFE tubing is the recommended tubing for use with our Omni-Lok™ fittings (see following page).

OMNI-LOK™ INVERTED CONE FITTINGS

Removable and reusable system for quick and convenient low-pressure connections

- Pressure rated up to 250psi (17 bar)
- For 1/16", 1/8" or 3/16" OD semi-rigid tubing e.g. PTFE, ETFE, FEP
- For flat-bottom 1/4"-28 UNF or 5/16"-24 UNF ports

Omni-Lok™ inverted cone fittings provide a simple, easy to use low-pressure connection. Only the ETFE cone and the tubing itself are in the fluid path.

No tools are required to assemble the flangeless fitting quickly and economically - just slip the fitting nut and the ETFE cone over the tubing and screw into the port. None of the parts are permanently attached to the tubing, so that the fitting nuts and inverted cones can easily be removed and re-used. A recess in the fitting nut houses the inverted cone. This allows maximum thread engagement with the port. The system seals up to 250psi (17 bar) pressure even in shallow PTFE ports. Note: The Omni-Lok™ inverted cone and fitting nut for 3/16" OD tubing and 5/16"-24 UNF flat-bottom ports is pressure rated up to 30 psi (2 bar).

Fitting nuts in robust, glass-filled polypropylene are available in a range of different colors for easy line identification. Nuts are also available in PEEK™ with standard and compact head designs (see the Omnifit® Fittings Systems Brochure for our full range).

For 1/16" OD Tubing

INVERTED CONES FOR 1/16" OD TUBING

PART NUMBER	DESCRIPTION	QTY
008CZ16	ETFE inverted cone	10pk

NUTS FOR 1/16" OD TUBING

PART NUMBER	MATERIAL	COLOR	THREAD	QTY
008NC16-YC5U	PP	Blue	1/4•28	10pk
008NC16-YC5G	PP	Green	1/4•28	10pk

For 1/8" OD Tubing

INVERTED CONES FOR 1/8" OD TUBING

PART NUMBER	DESCRIPTION	QTY
008CZ32	ETFE inverted cone	10pk

NUTS FOR 1/8" OD TUBING


PART NUMBER	MATERIAL	COLOR	THREAD	QTY
008NC32-YC5U	PP	Blue	1/4•28	10pk
008NC32-YC5G	PP	Green	1/4•28	10pk
008NC32-YC5N	PP	Orange	1/4•28	10pk
008NC32-YC5R	PP	Red	1/4•28	10pk
008NC32-YC5Y	PP	Yellow	1/4•28	10pk

TECH TIP:

Need more connection options?

See the Omnifit® Fitting Systems Brochure for our full range of threaded fittings, connectors and adaptors.

SPECIFICATIONS



ETFE inverted cone


For 3/16" OD Tubing

INVERTED CONE FOR 3/16" OD TUBING


PART NUMBER	DESCRIPTION	QTY
008CZ47	ETFE inverted cone	10pk

NUTS FOR 3/16" OD TUBING

PART NUMBER	MATERIAL	COLOR	THREAD	QTY
008NC47-YC7U	PP	Blue	5/16•24	10pk




008NC16-YC5U
Nut for 1/16" OD tubing
1/4•28, blue



008CZ16
Omni-Lok™ inverted
cone for 1/16" OD tubing

For use
with 030SP,
120SP &
130SP series
pumps



008NC47-YC7U
Nut for 3/16" OD tubing
5/16•24, blue



008CZ47
Omni-Lok™ inverted cone
for 3/16" OD tubing

For use with
150SP series
pumps

MICRO-PUMP TECH TIPS

OPERATING PARAMETERS

Design Specifications: Bio-Chem Fluidics’ Micro-Pumps are chemically compatible with a wide range of liquids, by virtue of the materials of construction. Specifications detailed in this brochure were determined via testing with distilled water under precise conditions. This means that the dispense rate for your pump may vary depending on your specific liquid. Other factors that can have an impact on operation include:

- Orientation
- Vertical distance between fluid reservoir and pump and then between pump and collection vessel
- Bore and length of inlet and outlet tubing
- Operating temperature

Please contact us to discuss your application and to get our recommendations for installation.

Pressure limits: Although Micro-Pumps are capable of producing outlet pressures of up to 5 psi (0.35 bar) while a dispense is taking place, for optimal dispense accuracy, the pressure on both the inlet and the outlet side of the pump should be kept between ± 0.5 psi (0.035 bar), equivalent to a head of ± 12” (300mm) water.


During the pump’s up-stroke, suction is created on the inlet. Positive pressure is generated at the outlet during the down-stroke. When the pump is not actuated, it will shut-off flow as long as the pressure on the inlet does not exceed the maximum holding pressure. To ensure correct operation, pressure on the inlet side should never exceed 2 psi (0.14 bar) even when the pump is in the closed position. The check valves in the pump prevent fluid from flowing against the intended flow direction.

Priming: Micro-Pumps must be fully primed prior to operation to ensure that all air is removed from the pump cavity. Priming is achieved by cycling the pump until no air bubbles are seen in the dispense. This normally takes 30-60 seconds. Excessive air bubbles in the dispense are generally caused by air leaks due to loose fittings - check all the fittings in the system and tighten accordingly.

Lead Wires: As a standard all lead wires are PTFE coated. Lead wires are provided with stripped ends for easy wiring into your control system - refer to drawings on product pages for more details. Different lengths and terminal connectors can be provided - refer to customization notes below.

INSTALLATION TIPS

Orientation: Pumps should be installed with the solenoid portion of the pump pointing upwards, downwards or in a horizontal position with the outlet on top. This ensures that any air in the system will be evacuated quickly and also minimizes the effects of a pressure head acting to keep the check elements open when they should be closed.



Preferred mounting positions

Tubing: Unlike centrifugal pumps where the outlet is normally larger than the inlet (to reduce the discharge head on the pump), our Micro-Pumps actually prefer to have the same sized tubing on the inlet and outlet. We recommend hardwall tubing for the connections and offer 1/16” OD x 1/32” ID (our part number 008T16-080) and 3/16” OD x 1/8” ID (008T47-032) PTFE tubing that can be installed using our Omni-Lok™ 1/4”-28 and 5/16”-24 fittings. For more details refer to pages 16 and 17.

Mounting options: The Micro-Pumps can be installed into your equipment with a variety of mounting options including mounting clips, rings and flanges. Some of the pumps can be mounted directly via mounting holes that are drilled into the pump body. For more details refer to the “Mounting Accessories & Options” spec sheet.



MC-100
Mounting clip



MR-100
Mounting ring

CUSTOMIZED SOLUTIONS

We understand that many applications require customized solutions. Our design and prototyping expertise enables us to offer simple modifications of standard products as well as completely customized designs. Over 90% of the Micro-Pumps we sell are customized to one extent or another. Customizable options include (but are not limited to):

- Materials of construction
- Operating voltage
- Dispense volume
- Mounting options
- Tagging / labeling
- Length and/or style of connecting leads
- Custom manifolds

We look forward to working with you to meet your design engineering objectives!


THE BIO-CHEM FLUIDICS BRAND FAMILY

Bio-Chem Fluidics is dedicated to providing instrument manufacturers and laboratories with the industry's best choice of inert, miniature fluid handling components.


Under our Bio-Chem Valve™ and Omnifit® brands we offer a complete fluid system solution for a wide range of industries including analytical chemistry, clinical diagnostics and medical device manufacturers as well as a world-class labware portfolio for the scientific community.



INERT SOLENOID VALVES AND PUMPS, ELECTRIC ROTARY VALVES



INERT FLUID HANDLING COMPONENTS AND LABWARE



Trademarks

PEEK™ is a registered trademark of Victrex plc.

Omnifit®, Omni-Lok™ are trademarks of Bio-Chem Fluidics Ltd.

Bio-Chem Valve™ is a trademark of Bio-Chem Fluidics Inc.



Bio-Chem Fluidics Inc
85 Fulton Street, Boonton NJ 07005 USA
t: (973) 263 3001 f: 973 263 2880 e: sales.us@biochemfluidics.com

Bio-Chem Fluidics Ltd
2 College Park, Coldhams Lane, Cambridge CB1 3HD UK
t: +44 (0) 1223 416642 f: +44 (0) 1223 416787 e: sales.eu@biochemfluidics.com