

# **High Speed PONA Analysis For Detailed Hydrocarbon Analysis Extended (DHAX) Using Hydrogen Carrier Gas For The Determination of Individual Components In Spark Ignition Fuels**

*Barry L. Burger, Neil Johansen, Valerie Gamble, Donald Rhoades*  
Restek Corporation 110 Benner Circle, Bellefonte, Pa. 16823 814.353.1300



**Turning Visions into Reality**  
[www.restek.com](http://www.restek.com)

## Abstract:

ASTM method D-6730-01 is specifically designed for the determination of the individual hydrocarbons present in spark ignition fuels, as well as fuel blends containing oxygenates such as MTBE, ETBE, t butanol, ethanol etc.

To maximize the resolution of these complex mixtures of spark ignition fuels, ASTM recommends a 100 meter x .025 mm id x 0.5um capillary column containing a deposition of 100% dimethyl polysiloxane liquid phase for the main analytical column. To control the selectivity of the aromatic compounds a short (typically 2-3 meters) tuning column containing a deposition of 5% diphenyl 95% dimethyl polysiloxane is connected to the inlet of the primary analytical column. To meet the demanding resolution and retention criteria in ASTM method D-6730-01, and in Canadian General Standards Board (CGSB) methodology for detailed hydrocarbon analysis (DHA), Restek Innovation chemists have reformulated our Rtx-1PONA column.



## Project Goal

ASTM method D-6730-01 suggests using helium as the carrier gas, at a linear velocity of 24cm/sec. (approximately 2.3 ml/min). The tridecane (C13) retention time this combination yields, approximately 146 minutes, greatly limits sample throughput per day. Restek's enhanced Rtx-1PONA column meets or exceeds all criteria in the method, but does so in 30% less time: retention time for C13 is 97 minutes, using helium as the carrier gas.

In most applications hydrogen is a better alternative to helium as the carrier gas, because it can be used at much higher linear velocities without compromising critical resolutions. A revision to ASTM D 6730-01 proposed by Neil Johansen Inc. (Aztec, New Mexico), in association with Envantage Analytical Software Inc. (Cleveland Ohio), has established optimal DHAX (detailed hydrocarbon analysis-extended) parameters, including specifying hydrogen as the carrier gas. Using these conditions analysis time is reduced to within 71 minutes (C13) – a 53% reduction versus using helium as the carrier gas. The method is also extended to include middle distillates having final boiling points up to 509 C / 948 F. (n-38).



# Introduction

Restek provided Neil Johansen Inc. with an enhanced Rtx 1 PONA column (100m x 0.25mm id x .5um df) and an Rtx 5 PONA tuning column for DHAX method development. Three (3) meters of the Rtx 5 tuning column was connected to the main analytical column, (Rtx 1 PONA) through a Universal Angled Press-tight Connector and was installed in a PerkinElmer Autosystem XL GC equipped with a flame ionization detector and programmable pneumatic control. The data system used was ChromPerfect Spirit (Justice Laboratory Software, Denville, NJ). The individual hydrocarbons were identified by using Dragon-DHA software, developed by Envantage Analytical Software Inc. in association with Neil Johansen Inc., which uses algorithms to process high-resolution chromatographic data.

The proposed DHAX method was optimized with hydrogen carrier gas at a rate of 3.62 mls/minute, constant flow, producing a linear velocity of 55 cm/sec. The column was briefly conditioned for 30 minutes at 300 C. The temperature was then lowered to 35 C. The dead time was adjusted to elute methane in 3.05 minutes  $\pm$  0.05 minutes, then a DHA Oxygenate setup blend was injected to determine the column's suitability for the method.



Figure # 1 illustrates the front end slice of the 30 component oxy set-up blend analysis. Key measurements such as: column efficiency, peak skewness, resolution and K' are calculated from this data. The calculated values shown here indicate the analytical column easily conforms to the ASTM D-6730-01 criteria.

Once established that the Rtx 1 PONA met method criteria, GC oven program parameters were entered and trial injections of the DHA/Oxy setup blend were begun. The analysis is permitted to run until the elution of C13 (tridecane). Based on the resolution of five critical pairs: t butanol / 2 methylbutene 2, 1 methylcyclopentene / benzene, 2,3,3 trimethylpentane / toluene, p xylene / 2,3 dimethylheptane, and 1 methylnaphthalene / C13, the tuning column length is incrementally reduced until the method criteria are met. In this example, the final tuning column length is 2.36 meters. Figure # 2 illustrates the complete 30 component Oxy set-up blend analysis. All critical pairs are sufficiently resolved, and meet method criteria. Tridecane (C13) retention time is less than 71 minutes. A 53% reduction in analysis time using hydrogen carrier gas, as compared to using helium. Figures # 3-6 illustrate the analysis of the PONA VI reference standard processed by Dragon DHA software. Note: The analysis of the 400-plus component PONA VI reference standard, plus a complete analytical report generated by Dragon-DHA software, can also be viewed on the Restek website: [www.restek.com/PONA](http://www.restek.com/PONA)




Fig. 1

## Rtx 1 PONA 100 meter x .25mm x .5um

Oven temperature 35 C Isothermal

Constant Flow H<sub>2</sub> @ 3.62 mls / min. ( 55 cm/sec.)

- C5 efficiency = 586,825 total theoretical plates
- C5 K' = 0.476
- t butanol skew = 2.10
- Resolution t butanol / 2methylbutene 2 = 5.39




Data courtesy of Neil Johansen in association with  
Envantage Analytical Software, Inc.

Fig. 2

**Rtx 1 PONA 100 meters x 0.25mm x .5um**  
**Tuning column; Rtx 5 2.36 meters x 0.25mm x 1um**

Data courtesy of Neil Johansen in association with  
 Envantage Analytical Software, Inc.


Peak #		
1.	ethanol	16. Toluene
2.	C5	17. C8
3.	t butanol	18. Ethylbenzene
4.	2 methylbutene 2	19. P xylene
5.	2,3 dimethylbutane	20. 2,3 DMH
6.	MTBE	21. C9
7.	C6	22. 5 methylnonane
8.	1 methylcyclopentene	23. 1,2 methylethylbenzene
9.	Benzene	24. C10
10.	cyclohexane	25. C11
11.	3 ethylpentane	26. 1,2,3,5,tetramethylbenzene
12.	1 t 2 dimethylcyclopentane	27. Naphthalene
13.	C7	28. C12
14.	2,2,3 trimethylpentane	29. 1 methyl naphthalene
15.	2,3,3 trimethylpentane	30. C13




**Column temperatrue program**


Initial temp..... 5C  
 Initial time..... 8.32 min (elution time of C5)  
 Program rate # 1...22 C / min.  
 1<sup>st</sup> hold temp .....48 C  
 1<sup>st</sup> hold time.....26.32 min. ( elution of ethylbenzene)  
 2<sup>nd</sup> program rate....3.20 C / min.  
 2<sup>nd</sup> temperature.....141 C (elution of dodecane)  
 2<sup>nd</sup> hold time..... 0 min.  
 3<sup>rd</sup> program rate... .1 C / min.  
 Final temp.....300 C.


Carrier gas.....H2  
 Constant Flow.....3.62 cc/min. ( 55 cm / sec )






















## Conclusion:

The benefits of hydrogen carrier gas are obvious: all critical components are resolved, per method D-6730-01, in the greatly reduced time of 70.5 minutes (C13), versus 146 minutes or 97 minutes using helium carrier gas. Relative to the results anticipated in the method, we virtually doubled sample throughput.

To assure that all Rtx 1 PONA columns will perform as demonstrated, each column is evaluated for film thickness, column efficiency, peak skewness, selectivity, resolution, and bleed to guarantee performance and reproducibility.

