

Selecting a GC Column

Strategic column choices can improve lab productivity by assuring that speed and resolution are optimized. While the number of choices available can be daunting, consideration of the resolution equation variables—separation factor, retention (capacity) factor, and efficiency—simplifies the decision. Separation factor determines which stationary phase is most appropriate. Once the phase has been chosen, physical dimensions (inner diameter, film thickness, length) can be selected based on retention factor and efficiency. Understanding how separation factor, retention factor, and efficiency influence separations allows analysts to make effective, informed choices and quickly select the best column for specific separations.

$$R = \frac{1}{4} \sqrt{N} \times \left(\frac{k}{k+1} \right) \times (\alpha - 1)$$

A measure of **Efficiency**.
This term is affected by:

- Length
- Inner diameter
- Carrier gas type and linear velocity

A measure of **Retention**.
This term is affected by:

- Inner diameter
- Film thickness
- Temperature

A measure of **Peak Separation**.
This term is affected by:

- Stationary phase composition
- Temperature

$N = L/H$ = Effective theoretical plate number
 L = Column length
 H = HETP = Height equivalent to a theoretical plate

k = Retention factor
 α = Separation factor
Baseline resolution ($R = 1.5$) is the goal.

Speed Up and Simplify GC Method Development with Restek's EZGC® Online Suite



- Model chromatograms
- Get column recommendations

- Translate methods **NEW!**
- Calculate flows **NEW!**

Download these FREE web apps today!

www.restek.com/ezgc

Separation Factor (α)

Choosing the right stationary phase is the first step toward optimizing your GC separation. It is the most important decision you will make because separation factor (α) has the greatest impact on resolution, and it is strongly affected by stationary phase polarity and selectivity.

Stationary phase polarity is determined by the type and amount of functional groups in the stationary phase. Structures for Restek stationary phases are presented in order of polarity on page 15. When choosing a column, consider the polarity of both the stationary phase and your target analytes. If the stationary phase and analyte polarities are similar, then the attractive forces are strong and more retention will result. Greater retention often results in increased resolution. Stationary phase polarity strongly influences column selectivity and separation factor, making it a useful consideration when selecting a column.

Stationary phase selectivity is defined by IUPAC as the extent to which other substances interfere with the determination of a given substance. Selectivity is directly related to stationary phase composition and how it interacts with target compounds through intermolecular forces (e.g., hydrogen bonding, dispersion, dipole-dipole interactions, and shape selectivity). As methyl groups in the stationary phase are replaced by different functionalities, such as phenyl or cyanopropyl pendant groups, compounds that are more soluble with those functional groups (e.g., aromatics or polar compounds, respectively) will interact more and be retained longer, often leading to better resolution and increased selectivity. In another example of the effect of stationary phase-analyte interactions, an Rtx®-200 stationary phase is highly selective for analytes containing lone pair electrons, such as halogen, nitrogen, or carbonyl groups, due to interactions with the fluorine pendant group in this phase. Selectivity can be approximated using existing applications or retention indices (Table I), making these useful tools for comparing phases and deciding which is most appropriate for a specific analysis.

Table I: Retention Indices for Restek Phases

Phase	Benzene	Butanol	Pentanone	Nitropropane
Rtx-1	651	651	667	705
Rtx-5/Rtx-5MS	667	667	689	743
Rtx-20	711	704	740	820
Rtx-1301/Rtx-624	689	729	739	816
Rtx-35	746	733	773	867
Rtx-200	738	758	884	980
Rtx-50	778	769	813	921
Rtx-1701	721	778	784	881
Rtx-65TG	794	779	825	938
Rtx-225	847	937	958	958
Stabilwax	963	1,158	998	1,230

Stationary phase polarity and selectivity also affect how much sample loading capacity the column will have for a particular analyte; the more soluble an analyte is in the stationary phase, the greater the sample loading capacity will be for that analyte. For example, a nonpolar stationary phase will have higher sample loading capacity for a nonpolar compound (e.g., pentane) than for a polar compound (e.g., ethanol).

Restek's Learning Network

Solutions to fit every level of expertise, budget, schedule, and analytical topic.

Learn more at www.restek.com/seminars



RESTEK
Distributor

CHROMALYTICS

in AUSTRALIA : Contact +81 3 9762 2034

SHOPPE

www.chromalytic.net.au

e-mail : sales @ chromtech.net.au

The relationship between polarity, selectivity, and sample loading capacity can be illustrated using an example. Consider the analysis of benzene and butanol, which have nearly the same boiling point, on an Rtx®-20 column (diphenyl dimethyl polysiloxane stationary phase). Since the benzene molecule is structurally more similar to the diphenyl phase than butanol is, benzene will solvate into the stationary phase more readily than butanol based on the concept that “like dissolves like.” Since benzene solvates more readily with the stationary phase, it has more interaction with the stationary phase as it elutes through the column and will be retained longer. Since butanol solvates less with the stationary phase, it has fewer interactions with the stationary phase and less will be retained. Therefore, the elution order of these two compounds on an Rtx®-20 column will be butanol first and benzene second. In addition, since benzene is more soluble in the diphenyl phase, the column has more capacity for benzene. This results in a more symmetrical peak shape for benzene than for butanol. A more polar column, such as a polyethylene glycol (PEG) column, will provide retention and better peak shape for butanol compared to benzene.

Due to their influence on separation factor, polarity and selectivity are primary considerations when selecting a column. However, temperature limits must also be considered. In general, highly polar stationary phases have lower maximum operating temperatures, so choosing a column with the appropriate maximum operating temperature, as well as optimal polarity and selectivity for the type of compounds being analyzed, is crucial.


Retention Factor (k)

The retention factor (k) of a column is based on the time an analyte spends in the stationary phase relative to the time it spends in the carrier gas. It is influenced primarily by column inner diameter (ID), phase film thickness, and temperature. Retention factor is sometimes referred to as capacity factor, which should not be confused with sample loading capacity. As a general rule, the thicker the film and the smaller the inner diameter, the more an analyte will be retained. Note that as temperature increases, k decreases; therefore, at higher temperatures analytes stay in the carrier gas longer and are less retained.

Chromatogram Search Tool

Search by **compound name**,
synonym, **CAS #**, or **keyword**

www.restek.com/chromatograms



When selecting column ID, consider the type of injection, the detector being used, and the concentration of sample (amount on-column). The injection technique is important because the column ID may need to be selected based on whether a split, splitless, direct, cool on-column injection, or other sample transfer method is being used. For example, 0.53 mm ID columns are ideal for cool on-column injections since the syringe needle (26 gauge) will fit into the large column ID. In addition to column ID, the detector and its flow requirements must be considered. For example, some MS detectors can only operate under column flow rates of up to 1.5 mL/min; therefore, a 0.53 mm ID column, which requires higher flows for proper chromatography, is not an option for MS work. Table II shows typical column characteristics for columns of various inner diameters.

Table II: General Column Characteristics Based on ID

Characteristic	Column Inner Diameter (mm)					
	0.10	0.15	0.18	0.25	0.32	0.53
Nitrogen flow (mL/min)	0.2	0.3	0.3	0.4	0.6	0.9
Helium flow (mL/min)	0.6	0.8	1.0	1.4	1.8	3.0
Hydrogen flow (mL/min)	0.7	1.1	1.3	1.8	2.3	3.7
Sample loading capacity (ng)	2.5	10	20	50	125	500
Theoretical plates/meter	11,000	7,000	6,000	4,000	3,000	2,000

Note: Flows listed are for maximum efficiency. Sample loading capacities are estimates only. Actual sample loading capacity varies with film thickness and analyte.

Fused Silica Capillary & PLOT Column Ferrule Guide

GC Column ID	Ferrule ID
0.15 mm	0.4
0.18 mm	0.4
0.25 mm	0.4
0.32 mm	0.5
0.53 mm	0.8

Film thickness has a direct effect on the retention and elution temperature for each sample component. Extremely volatile compounds should be analyzed on thick film columns to increase the time the compounds spend in the stationary phase, which allows them to better separate. High molecular weight compounds must be analyzed on thinner film columns. This reduces the length of time that the analytes stay in the column and minimizes phase bleed at higher elution temperatures. Film thickness also affects the amount of material that can be injected onto the column without overloading it. A thicker film column can be used for higher concentration samples.

Technical Service

For quick answers to commonly asked questions any time of the day, visit www.restek.com/answers or contact us directly:

In the U.S.

Phone: 1-800-356-1688, ext. 4

e-mail: support@restek.com

Hours of operation (Eastern Time):

Monday - Thursday, 8:00 a.m. to 6:00 p.m.

Friday, 8:00 a.m. to 5:00 p.m.

Outside the U.S.

Contact International Technical Service at intltechsupport@restek.com or find a local distributor at www.restek.com/distributor



Checking for leaks, using a Restek® electronic leak detector, is an easy way to protect your instrument and column from damage.

Film thickness also directly affects phase ratio (β), which must be accounted for when changing to a column with a different inner diameter. When inner diameter increases, film thickness (d_f) must also increase in order to provide comparable resolution and retention. Table III shows β values for common column dimensions; similar values indicate similar separations on different ID columns.

Table III: Phase Ratio (β) Values for Common Column Dimensions*

Column ID	Film Thickness (d_f) / β Value						
	0.10 μm	0.25 μm	0.50 μm	1.0 μm	1.5 μm	3.0 μm	5.0 μm
0.18 mm	450	180	90	45	30	15	9
0.25 mm	625	250	125	63	42	21	13
0.32 mm	800	320	160	80	53	27	16
0.53 mm	1325	530	265	128	88	43	27

* $\beta = r/2d_f$ (r =internal radius of tubing; d_f = phase film thickness)

Efficiency (N)

Column efficiency (N) is the column length divided by the height equivalent of a theoretical plate (HETP). The effective number of theoretical plates is affected by how well the phase has been coated onto the column walls, and it is measured by how narrow the peaks are when they elute out of the column. Higher column efficiency (N) results in greater resolution between peaks. Inner diameter also influences efficiency; a simple rule of thumb is the smaller the column ID, the more efficient the column.

Capillary columns are made in various lengths, typically 10, 15, 30, 60, and 105 meters. Longer columns provide more resolving power, but will also increase analysis time and cost more. When column length is doubled, analysis time will increase by as much as a factor of two. However, doubling the column length increases resolution by only approximately 40% since the column length term is under the square root function in the resolution equation. When selecting column length, the increase in resolution obtained in a longer column must be weighed against the increase in cost and analysis time.

Conclusion

A basic understanding of the resolution equation allows analysts to make more effective column choices. Phase choice should be influenced primarily by separation factor, which can be approximated by considering the structures of both the phase and the analyte, as well as by referencing retention indices or existing applications. Retention factor and efficiency also affect peak separations and should be considered when choosing column inner diameter, film thickness, and length. By better understanding these factors, analysts can simplify the column selection process, optimize separations, and increase lab productivity.

What Do the Temperature Limits Mean?

All Restek columns have published minimum and maximum operating temperatures that establish the working range for the stationary phase. Note that these ranges vary with the thickness of the coating.

Rxi®-5Sil MS Columns (fused silica)

ID	d_f (μm)	temp. limits
0.25 mm	0.25	-60 to 320/350 °C
0.32 mm	0.50	-60 to 320/350 °C
0.53 mm	1.50	-60 to 320/330 °C

The second temperature is the **maximum temperature-programmed operating temperature**, the temperature to which the column can be heated for short periods of time (i.e., during a temperature-programmed analysis). If only one temperature is listed, it is both the isothermal and the maximum temperature.


The **minimum operating temperature** defines the lowest usable temperature before the stationary phase solidifies. Operating the column below the minimum temperature will not harm the phase, but poor peak shape and other chromatography problems may occur.

Many phases list two maximum operating temperatures. The first temperature is the **maximum isothermal operating temperature**. This is the temperature to which the columns are guaranteed to meet the minimum bleed specification (i.e., lowest bleed level).

Structures, Polarities, Properties, and Uses for Restek® Capillary Column Phases, in Order of Increasing Polarity

Rxi®-1ms, Rxi®-1HT, Rtx®-1


Dimethyl polysiloxane



Similar to: (100%-methyl)-polysiloxane
Polarity: nonpolar
Uses: solvents, petroleum products, pharmaceutical samples, waxes
[G1] [G2] [G38]

Rxi®-5ms, Rxi®-5HT, Rtx®-5, Rtx®-5MS


Diphenyl dimethyl polysiloxane



Similar to: (5%-phenyl)-methylpolysiloxane
Polarity: slightly polar
Uses: flavors, environmental, aromatic hydrocarbons
[G27] [G36]

Rxi®-5Sil MS


1,4-bis(dimethylsiloxy)phenylene dimethyl polysiloxane



Similar to: (5%-phenyl)-methylpolysiloxane
Polarity: slightly polar
Uses: flavors, environmental, pesticides, PCBs, aromatic hydrocarbons

Rtx®-20


Diphenyl dimethyl polysiloxane



Similar to: (20%-phenyl)-methylpolysiloxane
Polarity: slightly polar
Uses: volatile compounds, alcohols
[G28] [G32]


Rtx®-1301, Rtx®-624, Rtx®-G43

Cyanopropylmethyl phenylmethyl polysiloxane



Similar to: (6%-cyanopropylphenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: volatile compounds, insecticides
[G43]


Rxi®-624Sil MS



Similar to: (6%-cyanopropylphenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: volatile compounds, insecticides, residue solvents in pharmaceutical products


Rtx®-35

Diphenyl dimethyl polysiloxane



Similar to: (35%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: pesticides, PCBs, amines, nitrogen-containing herbicides
[G42]


Rxi®-35Sil MS



Similar to: (35%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: pesticides, PCBs, amines, nitrogen-containing herbicides

Rtx®-200


Trifluoropropylmethyl polysiloxane



Similar to: (trifluoropropyl)-methylpolysiloxane
Polarity: selective for lone pair electrons
Uses: environmental, solvents, Freon® gases, drugs, ketones, alcohols
[G6]

Rtx®-50


Phenyl methyl polysiloxane



Similar to: (50%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: FAMES, carbohydrates
[G3]


Rxi®-17

Diphenyl dimethyl polysiloxane



Similar to: (50%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: triglycerides, phthalate esters, steroids, phenols
[G3]


Rxi®-17Sil MS



Similar to: (50%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: triglycerides, phthalate esters, steroids, phenols
[G3]

Rtx®-1701


Cyanopropylmethyl phenylmethyl polysiloxane



Similar to: (14%-cyanopropylphenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: pesticides, PCBs, alcohols, oxygenates
[G46]

Rtx®-65, Rtx®-65TG


Diphenyl dimethyl polysiloxane



Similar to: (65%-phenyl)-methylpolysiloxane
Polarity: intermediately polar
Uses: triglycerides, rosin acids, free fatty acids
[G17]

Rtx®-225


Cyanopropylmethyl phenylmethyl polysiloxane



Similar to: (50%-cyanopropylmethyl)-methylphenylpolysiloxane
Polarity: polar
Uses: FAMES, carbohydrates
[G7] [G19]

Rtx®-2330


Biscyanopropyl cyanopropylphenyl polysiloxane



Similar to: (95%-cyanopropyl)-phenyl polysiloxane
Polarity: polar
Uses: cis/trans FAMES, dioxin isomers, rosin acids
[G8] [G48]

Stabilwax®, Rtx®-Wax, Stabilwax®-MS

Polyethylene glycol



Polarity: polar
Uses: FAMES, flavors, acids, amines, solvents, xylene isomers
[G14] [G15]

note

Structures, polarities, and properties also apply to metal MXT® stationary phases.

Columns by Phase

Restek	Phase Description	USP Nomenclature*	Agilent	SGE	Phenomenex	Macherey- Nagel	Supelco	Alltech	Quadrex
Rtx-1 (p. 41) MXT-1 (p. 107)	dimethyl polysiloxane	G1, G2, G38	HP-1, DB-1, CP-Sil 5 CB	BP1	ZB-1	OPTIMA 1	SPB-1	AT-1, EC-1	007-1
Rxi-1HT (p. 39)	dimethyl polysiloxane		DB-1ht		ZB-1HTinferno			AT-1ht	
Rxi-1ms (p. 29)	dimethyl polysiloxane (low bleed)	G1, G2, G38	HP-1ms, HP-1msUI, DB-1ms, DB-1msUI, VF-1ms, Ultra 1	BP1	ZB-1, ZB-1ms	OPTIMA 1 MS, OPTIMA 1 MS Accent	SPB-1, Equity-1	AT-1ms	007-1
Rtx-5MS (p. 44) Rtx-5 (p. 42-43, 96) MXT-5 (p. 108)	diphenyl dimethyl polysiloxane	G27, G36	HP-5, DB-5, CP-Sil 8 CB	BP5	ZB-5	OPTIMA 5	SPB-5	EC-5, AT-5	007-5
Rxi-5HT (p. 39)	diphenyl dimethyl polysiloxane		DB-5ht, VF-5ht	HT5	ZB-5HTinferno	OPTIMA 5HT			
Rxi-5ms (p. 30)	diphenyl dimethyl polysiloxane (low bleed)	G27, G36	HP-5msSV, HP-5ms, HP-5msUI, DB-5, Ultra-2, CP-Sil 8 CB	BP5ms	ZB-5, ZB-5msi	OPTIMA 5 MS, OPTIMA 5 MS	SPB-5, Equity-5	AT-5ms	007-5
Rxi-5Sil MS (p. 32, 67, 75, 77, 85)	1,4-bis(dimethylsiloxy)phenylene dimethyl polysiloxane		DB-5ms, DB-5msUI, VF-5ms	BPX5	ZB-5ms, ZB- SemiVolatiles	OPTIMA 5MS Accent	SLB-5ms		007-5MS
Rxi-XLB (p. 31, 71)	proprietary phase		DB-XLB, VF-Xms		MR1, ZB-XLB	OPTIMA XLB			
Rtx-20 (p. 45)	diphenyl dimethyl polysiloxane	G28, G32					SPB-20	EC-20, AT-20	007-20
Rtx-35 (p. 46)	diphenyl dimethyl polysiloxane	G42	HP-35, DB-35	BPX35, BPX608	ZB-35		SPB-35, SPB-608	AT-35, AT-35ms	007-35
Rxi-35Sil MS (p. 34, 82)	proprietary phase		DB-35ms, DB35msUI, VF-35ms	BPX35	MR2	OPTIMA 35 MS			
Rtx-50 (p. 47) MXT-50 (p. 109)	phenyl methyl polysiloxane	G3	HP-50+, CP-Sil 24 CB				SPB-50	AT-50	007-17
Rxi-17 (p. 35)	diphenyl dimethyl polysiloxane	G3	HP-17, DB-17, DB-17ht, DB-608		ZB-50	OPTIMA 17	SPB-17		
Rxi-17Sil MS (p. 36, 76)	proprietary phase	G3	DB-17ms, VF-17ms	BPX50		OPTIMA 17 MS			
Rtx-65 (p. 47)	diphenyl dimethyl polysiloxane								007-65HT
Rxi-624Sil MS (p. 38, 81, 87, 94)	proprietary phase	G43	DB-624, VF-624ms, CP-Select 624 CB	BP624		OPTIMA 624 LB			
Rtx-1301 (p. 51) Rtx-624 (p. 52) MXT-1301 (p. 109)	cyanopropylmethyl phenylmethyl polysiloxane	G43	DB-1301, DB-624, DB-624UI, VF-1301ms, VF-624ms, CP-1301	BP624	ZB-624	OPTIMA 1301, OPTIMA 624	SPB-624	AT-624, AT-1301	007-1301, 007-624
Rtx-1701 (p. 53) MXT-1701 (p. 109)	cyanopropylmethyl phenylmethyl polysiloxane	G46	DB-1701P, DB-1701, CP-Sil 19 CB, VF-1701ms, VF-1701 Pesticides	BP10	ZB-1701, ZB-1701P	OPTIMA 1701	Equity-1701	AT-1701	007-1701
Rtx-200 (p. 48) MXT-200 (p. 110)	trifluoropropylmethyl polysiloxane	G6	DB-210, DB-200, VF-200ms			OPTIMA 210		AT-210	
Rtx-200MS (p. 48-49)	trifluoropropylmethyl polysiloxane (low bleed)	G6	VF-200ms						
Rtx-225 (p. 54)	cyanopropylmethyl phenylmethyl polysiloxane	G7, G19	DB-225ms, CP-Sil 43 CB	BP225		OPTIMA 225	SPB-225	AT-225	007-225
Rtx-440 (p. 50)	proprietary phase		Restek innovation						
Rtx-2330 (p. 55)	biscyanopropyl cyanopropylphenyl polysiloxane	G8, G48	VF-23ms	BPX70			SP-2330, SP-2331, SP-2380	AT-Silar90	007-23
Rt-2560 (p. 56, 83)	bicyanopropyl polysiloxane		HP-88, CP-Sil 88				SP-2560		
Rtx-Wax (p. 57)	polyethylene glycol	G14, G15, G16, G20, G39	DB-Wax, CP-Wax 52 CB	BP20	ZB-Wax	OPTIMA WAX		AT-WAXms, EC-WAX	007-CW
Stabilwax (p. 58, 95) Stabilwax-MS (p. 59) MXT-WAX (p. 110)	polyethylene glycol	G14, G15, G16, G20, G39	HP-INNOWax, CP-Wax 52 CB, VF-WAX MS		ZB-WAXplus	OPTIMA WAXplus	Supelcowax-10	AT-WAX	

See page 117 for Restek PLOT Column Phase Cross-Reference chart.

*See page 147 for our USP Liquid Phase and Solid Support Cross-Reference.



Application-Specific Columns by Industry

Restek	Applications	Agilent	Supelco	Macherey-Nagel	SGE	Alltech	Phenomenex
Chiral Columns							
Rt-βDEXm, Rt-βDEXsm, Rt-βDEXse, Rt-βDEXsp, Rt-βDEXsa, Rt-βDEXcst, Rt-γDEXsa (p. 103)	Chiral compounds						
Clinical, Forensic, & Toxicology							
Rtx-BAC Plus 1 (p. 65)	Blood alcohol testing	DB-ALC1					ZB-BAC1
Rtx-BAC Plus 2 (p. 65)		DB-ALC2					ZB-BAC2
Environmental							
Rxi-5Sil MS (p. 67, 75, 77)	Semivolatiles - EPA Methods 8270, 625, 525	DB-5ms, DB-5msUI, VF-5ms	SLB-5ms	OPTIMA 5MS Accent	BPX5		ZB-5ms, ZB-SemiVolatiles
Rtx-VMS (p. 78)	Volatiles - EPA Methods 8260, 624, 524	Restek innovation					
Rxi-624Sil MS (p. 81)	Volatiles - EPA Method 624	DB-624, VF-624ms, CP-Select 624 CB		OPTIMA 624 LB	BP624		
Rtx-502.2 (p. 80)	Volatiles - EPA Methods 8010, 8020, 502.2, 601, 602	DB-502.2	VOCOL			AT-502.2	
Rtx-Volatiles (p. 80)			VOCOL				
Rtx-VRX (p. 79)		DB-VRX					
Rtx-CLPesticides (p. 72)	Organochlorine pesticides -	DB-CLP1					ZB-CLP1
Rtx-CLPesticides2 (p. 72)	EPA Methods 8081, 8082, 608, 505, 508	DB-CLP2					ZB-CLP2
Rtx-1614 (p. 66)	Brominated flame retardants	Restek innovation					
Rtx-PCB (p. 70)	Polychlorinated biphenyl -	Restek innovation					
Rxi-XLB (p. 71)	EPA Methods 8082, 608, PCB congeners	DB-XLB, VF-XMS					MR1, ZB-XLB
Rtx-OPPesticides (p. 74)	Organophosphorus pesticides -	Restek innovation					
Rtx-OPPesticides2 (p. 74)	EPA Method 8141	Restek innovation					
Rtx-Dioxin2 (p. 68)	Dioxin & Furans - EPA Methods	Restek innovation					
Rtx-Mineral Oil (p. 69)	DIN EN ISO 9377-2	Select Mineral Oil					
Rxi-17Sil MS (p. 76)	Polycyclic aromatic hydrocarbons	DB-17ms, VF-17ms		OPTIMA 17 MS	BPX50		
Foods, Flavors, & Fragrances							
Rt-2560 (p. 83)	cis/trans FAMES	HP-88	SPB-2560				
FAMEWAX (p. 83)	Marine oils	Select FAME	Omegawax			AT-AquaWax, AT-FAME	
Rxi-PAH (p. 84)	PAHs	Restek innovation					
Rtx-65 TG (p. 89)	Triglycerides	Restek innovation					
Petroleum & Petrochemical							
Rt-Alumina BOND/CFC (p. 123)	Chlorinated fluorocarbons (CFCs)						
Rtx-DHA (p. 92)	Detailed hydrocarbon analysis	HP-PONA, DB-Petro, CP-Sil PONA CB	Petrocol DH		BP1PONA		
Rtx-2887 (p. 93)	Hydrocarbons - ASTM 2887	DB-2887	Petrocol 2887, Petrocol EX2887			AT-2887	
MXT-2887 (p. 113)							
D3606 (p. 138)	Ethanol - ASTM 3606	Restek innovation					
Rt-TCEP (p. 90)		CP-TCEP	TCEP				
MXT-1HT SimDist (p. 114)	Simulated distillation	DB-HT-SimDis, CP-SimDist, CP-SimDst Ultimet			BPX1	AT-3710	ZB-1XT SimDist
MXT-500 SimDist (p. 115)		Restek innovation					
Rtx-Biodiesel TG (p. 91)	Triglycerides in biodiesel	Biodiesel, Select Biodiesel		OPTIMA Biodiesel			ZB-Bioethanol
MXT-Biodiesel TG (p. 113)							
Pharmaceutical							
Rtx-G27 w/IntegraGuard (p. 97)	Organic volatile impurities (OVI) - USP 467						
Rtx-G43 w/IntegraGuard (p. 97)			OVI-G43				
Rxi-624Sil MS (p. 94)		DB-624, VF-624ms, CP-Select 624 CB		OPTIMA 624 LB	BP624		
Rtx-5 (G27) (p. 96)		HP-5, DB-5, CP-Sil 8 CB	SPB-5	OPTIMA 5	BP5	EC-5, AT-5	ZB-5
Stabilwax (G16) (p. 95)		HP-INNOWax, CP-Wax 52 CB, VF-WAX MS	Supelcowax-10	OPTIMA WAXplus		AT-WAX	ZB-WAXplus
Specially deactivated phases							
Rtx-Volatile Amine (p. 99)	Volatile amines	CP-Volamine					
Rtx-5Amine (p. 100)	Amines			OPTIMA 5 Amine			
Rtx-35Amine (p. 101)		Restek innovation					
Stabilwax-DB (p. 102)		CAM, CP-WAX 51 for Amines	Carbowax Amine	FS-CW 20 M-AM		AT-CAM	
Stabilwax-DA (p. 98)	Free fatty acids	HP-FFAP, DB-FFAP, CP-WAX 58 FFAP CB	NUKOL	PERMABOND FFAP, OPTIMA FFAP, OPTIMA FFAP Plus	BP-21	AT-AquaWax DA, AT-1000	ZB-FFAP