


The GC Racer is a Fast Temperature Programmer. It overcomes the power limitations of conventional GC ovens by adding a second convection heater to the oven. The GC Racer is controlled by the "host" GC and can boost ramp rates up to 120 °C/min.

## Advantages of the GC Racer

<b>Saves Time</b>	Faster runs mean higher efficiency, quick turnaround and greater profitability.
<b>Chromatographic Integrity</b>	Does not require specialized columns or interfaces to connect to the injector and detector.
<b>Easy to Implement</b>	Installs in minutes and uses standard columns from any vendor.
<b>Fully Integrated</b>	It is fully controlled by the GC keypad or operating software. No operator training required.
<b>Rugged</b>	Employs the same type of convective heater used in the GC oven.
<b>Utility</b>	One system can be used for dual column GCs.
<b>Versatile</b>	Can run on 110 - 240V, 50/60Hz and coupled with 110 - 240V GCs.
<b>Cost Effective</b>	One-fifth the purchase price compared to resistively heated accessories and no expensive specialized columns.

## Learn more about the GC Racer's


- Performance Data: [Heating Power](#) > [Heating Performance](#) > [RT Precision](#) > [Fast Analysis](#)



**Heating Power:** The GC Racer can more than double the heating power of the GC. Power can be measured in Volt-Amps (VA). An Agilent GC oven runs on 120V, 13 amp = 1560VA. The GC Racer supplies 1800 VA (120V X 15 amp) additional power for a total of 3360VA heating power. By contrast a "fast-ramp" oven (240 V X 10 amp) = 2400VA.

**Q: Why is the GC Racer better than a "fast oven"?**

**A: More heating power and you don't need 240V circuits.**




## GC Power Ratings with and without GC Racer

GC Oven			GC Racer			Total Power (VA)
Volts	Amps	Power (VA)	Volts	Amps	Power (VA)	
120	13	1560	120	15	1800	3360
240	10	2400	Without GC Racer			2400
240	10	2400	240	13	3120	5520

Learn more about the GC Racer's


# Zip Scientific Fast GC

Good People, Top-notch Products and Excellent Service



## ROI Calculator

**Retention Time Precision for 900 runs over 7 days:** The GC Racer performance can be evaluated by measuring the relative standard deviation (RSD) of the RT for each peak eluted at 70°C/min ramp rate. This graph plots the retention time for each peak from a sample containing 16 hydrocarbons from C10 – C25.



Temp Program:  
50°C (0.33 min hold)  
70°C/min to 300°C  
(hold 0.1 min)

Splitless injection with  
purge on at 0.66 min

Column: Rtx-5  
15m 0.25mm X 0.25um

Inj = 300°C  
FID = 325°C

## Learn more about the GC Racer's

- Performance Data: [Heating Power](#) > [Heating Performance](#) > [Fast Analysis](#)




www.chromtech.net.au E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel : +61 3 9762 2034 Fax : +61 3 9761 1169




# Zip Scientific Fast GC

Good People, Top-notch Products and Excellent Service



**FL PRO:** Here's a chromatogram and report from a Fast GC separation of hydrocarbons ranging from C8 to C40. Run time was less than 6 min. Use the zoom feature on your browser to inspect the peaks closely.



```

=====
Injection Date : 2/1/06 12:18:29 PM          Seq. Line : 1
Sample Name   : FL TPN 8000MS              Vial : 4
Acq. Operator : Steve MacDonald              Inj : 1
                                           Inj Volume : 1 µl

Acq. Method   : C:\NPGHEM\1\METHODS\FID_MC.M
Last changed  : 2/1/06 12:00:13 PM by Steve MacDonald
Analysis Method : C:\NPGHEM\1\METHODS\FID_MC.M
Last changed  : 2/1/06 1:40:29 PM by Steve MacDonald
              (modified after loading)
=====
    
```

## Height Percent Report

```

=====
Sorted By      : Signal
Multiplier     : 3.0500
Dilution       : 1.0000
Sample Amount  : 20.00000 [µl] (not used in calc.)
    
```

Signal 1: FID1 8,

Peak #	RetTime [min]	Type	Width [min]	Area counts*s	Height [counts]	Height %
1	0.407	EV	5.73e-3	8410.89160	2.18155e4	4.38534
2	0.584	EV	5.81e-3	7.51584e4	2.16044e4	43.42828
3	1.012	EV	6.48e-3	7.45015e4	1.90339e4	38.24298
4	1.540	EV	7.46e-3	7.72130e4	1.64946e4	33.15743
5	1.968	EV	7.29e-3	7.05807e4	1.64587e4	33.00529
6	2.324	EV	7.91e-3	7.91226e4	1.56335e4	31.42507
7	2.674	EV	8.73e-3	8.18097e4	1.47186e4	29.58733
8	2.996	EV	7.92e-3	8.26824e4	1.62395e4	32.78298
9	3.272	EV	8.05e-3	8.43944e4	1.34886e4	27.11526
10	3.535	EV	8.94e-3	8.65440e4	1.35675e4	27.27335
11	3.777	EV	8.92e-3	8.85245e4	1.53595e4	30.98467
12	4.008	EV	8.87e-3	8.82699e4	1.46037e4	29.22500
13	4.219	EV	8.01e-3	8.82755e4	1.52317e4	30.61882
14	4.424	EV	8.09e-3	8.62029e4	1.46940e4	29.53784
15	4.623	EV	8.88e-3	9.04315e4	1.41490e4	28.44237
16	4.812	EV	0.5103	4.91539e4	1.24980e4	21.10207
17	5.005	EV	0.5107	5.98795e4	8.54587e4	17.17893
18	5.231	EV	0.5154	4.00752e4	4.70954e4	9.48714

Totals : 1.33723e6 2.48731e6

Results obtained with enhanced integrator?

\*\*\* End of Report \*\*\*

Learn more about the GC Racer's

- Performance Data: [Heating Power](#) > [Heating Performance](#) > [RT Precision](#)




www.chromtech.net.au E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel : +61 3 9762 2034 Fax : +61 3 9761 1169






# Zip Scientific Fast GC

Good People, Top-notch Products and Excellent Service



**Interface:** Our Fast GC Accessories interface with the host GC by sensing the current applied to the oven heater. A current transformer snaps around one of the heater wires of the GC. That's it! Simple, Easy, non-Intrusive.

The GC Racer applies a proportionate amount of current to the accessory heater whenever oven heater current is sensed. Both heaters work together to achieve fast temperature programming up to 120°C/min.



The GC Chaser determines the end of the temperature program when the oven current is stopped. The blower turns on to start the fast oven cool-down cycle and stops automatically when the oven reaches the start temperature.

**HROMalytic Chromatography**  
**Products '08**  
**Australian Distributors ECH nology**

www.chromtech.net.au E-mail : info@chromtech.net.au Tel : +61 3 9762 2034 Fax : +61 3 9761 1169

**RESTEK**