

13.0 Product case specifications



14.0 Calibration and service

The Restek ProFLOW 6000 comes factory calibrated and carries a one year warranty from time of purchase. All units are calibrated to NIST traceable standards.

Recommended schedule for recalibration is once every year from time of purchase. Customers will need to return the unit to Restek for recalibration. At that time, preventative maintenance services can also be performed. A fee will be charged for recalibration and servicing of the unit. Prolonged failure to recalibrate the instrument may result in increased error.

**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4 (or your Restek representative)
if you have any questions about this product or any other Restek product.**

Please have the serial number available when calling Restek with any concerns you may have.
Additional charges may apply if the warranty is expired or the damage is due to misuse.

This manual is also available in electronic format at **www.restek.com**.



www.restek.com



Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au