



Operating Instructions

1.0 Introduction

Restek's portable Leak Detector is specifically designed for use with gas chromatography (GC) systems. It detects minute leaks of any gas with a thermal conductivity different from air. The reference gas inlet (Figure 1) draws in ambient air for comparison to air drawn into the sample probe. A leak is indicated by both an LED light display and an audible alarm.

Restek's Leak Detector is manufactured by Restek so you are assured of the same Restek quality and service you have come to recognize with the Restek name.

Should you require assistance at anytime regarding our Leak Detector, please contact Restek Technical Service at 1-800-356-1688 or 814-353-1300, ext. 4.