



## DYNAMIC DATA EXCHANGE

Link PeakSimple to your DDE compatible spreadsheet or word processor (Excel, Word, 123, etc.). Analytical results are automatically transferred after every run, or can be accumulated within PeakSimple and copied as a block of data. Use the Copy Picture option to paste the chromatogram itself into Excel, etc. along with the results.




## MERGE RESULTS FROM MULTIPLE CHANNELS

PeakSimple lets you merge the results from any channel or all channels into one report. This feature is handy when you're combining results from different detectors into a single report for export to Excel or other data analysis programs.



## PRINT CHROMATOGRAMS IN COLOR



Use any Windows™ supported color printer to create impressive hardcopies of your chromatographic data.

Print multiple chromatograms per page for easy detector-to-detector comparisons and paperwork consolidation. Print overlaid data in contrasting colors with adjustable line weight.