

For these columns:

MXT®-35, MXT®-50, MXT®-65, MXT®-65TG, MXT®-200, MXT®-624, MXT®-1301, MXT®-1701, MXT®-502.2, MXT®-Volatiles

Rtx®-17, Rtx®-20, Rtx®-35, Rtx®-35MS, Rtx®-35 Amine, Rtx®-50, Rtx®-65, Rtx®-65TG, Rtx®-200, Rtx®-200MS, Rtx®-624, Rtx®-1301, Rtx®-1701, Rtx®-BAC1, Rtx®-BAC2, Rtx®-CLPesticides, Rtx®-CLPesticides2, Rtx®-Dioxin, Rtx®-Dioxin2, Rtx®-OPPesticides, Rtx®-OPPesticides2, Rtx®-G43, Rtx®-502.2, Rtx®-VMS, Rtx®-VGC, Rtx®-VRX, Rtx®-Volatiles

Stx -CLPesticides, Stx -CLPesticides2

Column Length (meters)	Film Thickness (µm)	Approx. Time	
		(min.)	(hr.)
≤15	0.1 - 0.25	20	0.3
	0.5 - 1.0	40	0.7
	1.0 - 1.5	60	1
	1.5 - 3.0	80	1.3
30	0.1 - 0.25	40	0.7
	0.5 - 1.0	60	1
	1.0 - 1.5	80	1.3
	1.5 - 3.0	100	1.7
≥60	0.1 - 0.25	80	1.3
	0.5 - 1.0	120	2
	1.0 - 1.5	160	2.7
	1.5 - 3.0	200	3.3

For these columns:

DEX chiral phases

FAMEWAX

MXT®-WAX

Rt-CW20M F&F, Rt-2560, Rt-TCEP

Rtx®-225, Rtx®-2330, Rtx®-WAX

Stabilwax®, Stabilwax®-DA, Stabilwax®-DB

Column Length (meters)	Film Thickness (µm)	Approx. Time	
		(min.)	(hr.)
≤15	0.1 - 0.25	30	0.5
	0.5 - 1.0	45	0.75
	1.0 - 1.5	60	1
	1.5 - 3.0	90	1.5
30	0.1 - 0.25	60	1
	0.5 - 1.0	90	1.5
	1.0 - 1.5	120	2
	1.5 - 3.0	150	2.5
≥60	0.1 - 0.25	80	1.3
	0.5 - 1.0	120	2
	1.0 - 1.5	160	2.7
	1.5 - 3.0	200	3.3