

ISO 80369 Compliant Connectors

Nordson MEDICAL now offers compliant connectors for the new ISO 80369-3 standards for enteral and gastric small bore connectors as well as the ISO 80369-5 limb cuff inflation standards for blood pressure cuffs. **PG-12**



Nordson MEDICAL's Quick Connect Line is Incredible

Incredibly broad - We have over 170 quick connect configurations, arguably the most extensive line of quick connect couplers developed exclusively for the medical market.

Incredibly versatile - Our quick connector series now offer a full range of fluid path diameters from 1/8" to 1/2" in push-to-connect, quarter-turn, thumb latch, and side-latch connections. Also, the available pressure range is staggering, with a variety of sealing surfaces to choose from. The connectors also offer resin choices that afford multiple sterilization options, MRI compatibility, animal-free options, as well as bioburden, cytotoxicity, and extractable testing data available upon request. **PG-16**



PQC Series Small Bore Quick Connects

PQC Series products are produced in acetal and polypropylene in four configurations – threaded, inline straight, elbows and panel mounts. These quick connectors bring performance and convenience to your tubing applications. Offered in small tubing sizes to fit 1/8" (3.2 mm) and 1/4" (6.4 mm) ID tubing, the 3200 Series triple barb style works well with semi-rigid tubing. **PG-17**



High Flow Check Valve

Featuring a flowrate $\geq 3,000$ ml/min with water at a pressure of 1 bar, our high flow check valve fits 8.8 O.D. tubing sizes and has a back pressure of 1 bar/14.5 psig. The check valve is designed for flushing and irrigation applications and for use with certain insufflation tubing sets used in minimally invasive surgeries. **PG-31**



FC Series Stopcocks

Nordson MEDICAL's stopcocks are designed to easily control the flow of a liquid or gas through the use of an on/off diverter, and are available in three basic configurations – 1-way, 3-way, and 4-way. These chemically resistant stopcocks are available as standard or high flow and are offered in polycarbonate, PEI and nylon with HDPE diverters. **PG-36**



RENOLIT SOLMED Granuflex 4301 Resin

Nordson MEDICAL is now offering our popular bag ports made with RENOLIT SOLMED Granuflex 4301 resin. These products demonstrate lower leachable levels and high purity. **PG-39**



PharmaLok™ Clamps

Nordson MEDICAL's line of PharmaLok™ clamps feature a state-of-the-art design that simplifies the assembly process by making it effortless to securely attach sanitary fittings to bag ports and filters. All PharmaLok™ clamps are ASME BPE compliant and operate at a maximum pressure of 80 psi. Single-use and easy-open clamp sizes include .5" to .75" and 1" to 1.5", and the tube clamp supports tubing sizes from .5" to 1.25". The tube clamp is also designed for use with silicone and TPE tubing. **PG-45**



RQX Series Quick Connects

RQX Series quick connects feature unsurpassed ergonomics with a large locking button for easier handling and connection as well as a single barb open flow design. Available in widely used polypropylene and polycarbonate materials, these connectors are compatible with industry standard quick connects to simplify qualification and fit 3/8" and 1/2" tubing sizes. **PG-49**



TBA Series Tuohy Borst Adapters

Tuohy Borst adapters are designed to facilitate catheter introduction and are intended for interventional and diagnostic procedures where fluid backflow is a concern. **PG-56**



Need Samples for Evaluation?

You can call to request most of our standard product in quantities of 5, or less, for free application evaluation.

+1 (970) 267-5200 Toll Free: (888) 404-5837