

**ENTME050 (-8003)**

Enteral Male Connector with 5.0 FR Tubing Bond-In Pocket in Natural ABS



**ENTME064 (-8003)**

Enteral Male Connector with 6.4 FR Tubing Bond-In Pocket in Natural ABS



**ENTME078 (-8003)**

Enteral Male Connector with 7.8 FR Tubing Bond-In Pocket in Natural ABS



**ENTME098 (-8003)**

Enteral Male Connector with 9.8 FR Tubing Bond-In Pocket in Natural ABS



**ENTM230 (-6005)**

Enteral Male Connector Barb for 1/8" (3.2 mm) ID Tubing in Animal-Free Polypropylene



**ENTFE098 (-8003)**

Enteral Female Connector with 10 FR Tubing Bond-In Pocket in Natural ABS



**ISO 80369-3 Accessories:**

Tethered caps and plugs are not covered under the ISO 80369-3 standard, but are compatible with our enteral male and female connectors.

**ENTFTC (-VP1)**

Tethered Cap for Enteral Female Connector in Animal-Free Polyethylene



**ENTMTP (-VP1)**

Tethered Plug for Enteral Male Connector in Animal-Free Polyethylene



# ISO 80369-5

## Blood Pressure

Nordson MEDICAL continues to offer a wide range of connectors for non-invasive blood pressure (NIBP) monitoring equipment. With the publication of the new ISO 80369-5 limb cuff inflation standard for blood pressure cuffs, many organizations will move towards this new worldwide standard – bayonet-style connectors. The tight dimensional specifications for the compliant connectors will drastically reduce the probability of a misconnection. Nordson MEDICAL has already modified designs for parts to be 100% compliant with the blood pressure male connectors that fit singly or in pairs on blood pressure cuffs.

These industry-standard male connectors work well with the Nordson MEDICAL female blood pressure connectors, Parker/Rectus connectors, and other bayonet-style females, both compliant and non-compliant. And the Nordson MEDICAL single-barb designs offer an excellent secure connection with better tensile pull results and hydrostatic burst testing results than traditional multi-barb designs. The toughness of the DuPont Zytel nylon is well proven, and Nordson MEDICAL has been making these connectors for more than a decade.

### ISO 80369-5 COMPLIANT PARTS:

#### Bayonet Male Cuff Connectors

**BPM013 (-1)**

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 80369-5 Compliant, 1/8" (3.2 mm) ID Tubing



**BPM240 (-1)**

Blood Pressure Monitor Bayonet Connector to 200 Series Barb, 80369-5 Compliant, 5/32" (4.0 mm) ID Tubing



**BPM025 (-1)**

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 80369-5 Compliant, 5/32" (4.0 mm) ID Tubing



**BPM035 (-1)**

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 80369-5 Compliant, 3/16" (4.8 mm) ID Tubing

