

QUICK CONNECT

Fittings



MQC Series Female Fittings
PG-16



MQC Series Male Fittings
PG-17



PQC Series Female Fittings
PG-17



PQC Series Male Fittings
PG-20



RQC Series Female Fittings
PG-24



RQC Series Male Fittings
PG-24



RQX Series Female Fittings
PG-25



RQX Series Male Fittings
PG-25



SBL Series Female Fittings
PG-26



SBL Series Male Fittings
PG-26



XQ Series Female Fittings
PG-27



XQ Series Male Fittings
PG-28



Our medical quick connect line is incredible -

Nordson MEDICAL's quick connect line offers six distinct families of quick connect fittings (quick disconnect couplings), the MQC, PQC, RQC, RQX, SBL, and XQ Series products, for applications covering 1/8" (3.2 mm) through 1/2" (12.7 mm) ID tubing.

MQC Series

Quick Connects

The MQC Series is an easy to use Push-to-Connect design and arguably one of the most versatile quick connect solutions on the market. With its large, ergonomic buttons and wide selection of colors, both designers and end users will embrace its unique features. The MQC Series of connectors are available in nylon and fit 1/8" (3.2 mm), 5/32" (4.0 mm) and 3/16" (4.8 mm) ID tubing. Between its versatility, one-of-a-kind design, and loud "click" indicating a secure fit, misconnects will be a problem of the past!

MQC SERIES: Female Open Flow

MQCF230 (-1 -2 -3 -4 -5 -8)

Dual Latch Quick Connect Female 200 Series Barb 1/8" (3.2 mm) Tubing ID



MQCF240 (-1 -2 -3 -4 -5 -8)

Dual Latch Quick Connect Female 200 Series Barb 5/32" (4.0 mm) Tubing ID



MQCF250 (-1 -2 -3 -4 -5 -8)

Dual Latch Quick Connect Female 200 Series Barb 3/16" (4.8 mm) Tubing ID

