

MQC SERIES: Male Open Flow

MQCRM230 (-001 -002 -003 -004 -005 -08B)

Nylon, Buna-N O-Ring

Dual Latch Quick Connect Male 200 Series Barb 1/8" (3.2 mm) Tubing ID



-001 White Nylon
-005 Blue Nylon

MQCRM240 (-001 -002 -003 -004 -005 -08B)

Nylon, Buna-N O-Ring

Dual Latch Quick Connect Male 200 Series Barb 5/32" (4.0 mm) Tubing ID



-001 White Nylon
-005 Blue Nylon

MQCRM250 (-001 -002 -003 -004 -005 -08B)

Nylon, Buna-N O-Ring

Dual Latch Quick Connect Male 200 Series Barb 3/16" (4.8 mm) Tubing ID



-001 White Nylon
-005 Blue Nylon

MQC SERIES: Male and Female Sets

MQCMF230-901

Bag of 50 Dual Latch Quick Connect Male and Female Sets, 200 Series Barb 1/8" (3.2 mm) Tubing ID, Blue Female and White Male



MQCMF240-901

Bag of 50 Dual Latch Quick Connect Male and Female Sets, 200 Series Barb 5/32" (4.0 mm) Tubing ID, Blue Female and White Male



MQCMF250-901

Bag of 50 Dual Latch Quick Connect Male and Female Sets, 200 Series Barb 3/16" (4.8 mm) Tubing ID, Blue Female and White Male



PQC Series

Quick Connects

PQC Series products are produced in acetal and polypropylene in four configurations – threaded, inline straight, elbows and panel mounts. These quick connectors bring performance and convenience to your tubing applications. Offered in small tubing sizes to fit 1/8" (3.2 mm) and 1/4" (6.4 mm) ID tubing, the 3200 Series triple barb style works well with semi-rigid tubing.

PQC SERIES: Female Open Flow Threaded

PQCF18-AC-001

Female Open Flow Connector, 1/8-27 NPT Thread, Natural Acetal



PQCF18-PP-001

Female Open Flow Connector, 1/8-27 NPT Thread, Almond Polypropylene



PQCF18-PP-002

Female Open Flow Connector, 1/8-27 NPT Thread, Grey Polypropylene



PQCF14-AC-001

Female Open Flow Connector, 1/4-18 NPT Thread, Natural Acetal



PQCF14-PP-001

Female Open Flow Connector, 1/4-18 NPT Thread, Almond Polypropylene



PQCF14-PP-002

Female Open Flow Connector, 1/4-18 NPT Thread, Grey Polypropylene



QUICK
CONNECTS