

PQC SERIES:
Male Valved Threaded

PQCV18-AC-001

Male Valved Connector, 1/8-27 NPT Thread, Natural Acetal, Buna-N O-Ring



PQCV18-PP-001

Male Valved Connector, 1/8-27 NPT Thread, Almond Polypropylene, EPDM O-Ring



PQCV18-PP-002

Male Valved Connector, 1/8-27 NPT Thread, Grey Polypropylene, EPDM O-Ring



PQCV14-AC-001

Male Valved Connector, 1/4-18 NPT Thread, Natural Acetal, Buna-N O-Ring



PQC SERIES:
Male Open Flow Inline Straight (Barbed)

PQCM3230-AC-001

Male Open Flow Connector, 3200 Series Barb, 1/8" (3.2 mm) Tubing ID, Natural Acetal, Buna-N O-Ring



PQCM3230-PP-001

Male Open Flow Connector, 3200 Series Barb, 1/8" (3.2 mm) Tubing ID, Almond Polypropylene, EPDM O-Ring



PQCM3230-PP-002

Male Open Flow Connector, 3200 Series Barb, 1/8" (3.2 mm) Tubing ID, Grey Polypropylene, EPDM O-Ring



PQCM3260-AC-001

Male Open Flow Connector, 3200 Series Barb, 1/4" (6.4 mm) Tubing ID, Natural Acetal, Buna-N O-Ring



PQCM3260-PP-001

Male Open Flow Connector, 3200 Series Barb, 1/4" (6.4 mm) Tubing ID, Almond Polypropylene, EPDM O-Ring



PQCM3260-PP-002

Male Open Flow Connector, 3200 Series Barb, 1/4" (6.4 mm) Tubing ID, Grey Polypropylene, EPDM O-Ring



QUICK
CONNECTS

Need Samples for Evaluation?

You can call to request most of our standard product in quantities of 5, or less, for free application evaluation.

+1 (970) 267-5200 Toll Free: (888) 404-5837

Need More Details?

Check out the MQC, PQC, RQ, SBL, and XQ Datasheets at www.nordsonmedical.com or contact Sales at sales@nordsonmedical.com.