

GASKETS

Material options include Class VI compliant PTFE/EPDM and platinum-cured silicone in gasket sizes that fit both the 1-1/2" (Maxi) and the 3/4" (Mini) Flange Series. Keeping customer requirements in mind, these components were manufactured for applications requiring the highest sanitary standards and chemical purity.

GMN450 (-PX)

Molded Gasket for Mini Flange SFMN Series Fittings



G42MPGR-TF075 (-E)

Molded Tuf-Flex® Gasket for Mini Flange SFMN Series Fittings



GMX455 (-PX)

Molded Gasket for Maxi Flange SFMX Series Fittings



G40MPGR-TF150 (-E)

Molded Tuf-Flex® Gasket for Maxi Flange SFMX Series Fittings



G424TRX-0200 Platinum-Cured Silicone
Molded Gasket for 2" Flange Fittings



G424TRX-0300 Platinum-Cured Silicone
Molded Gasket for 3" Flange Fittings



Need More Details?

Check out all of our bioprocessing components at www.bioprocesscomponents.com or contact Sales at sales@nordsonmedical.com.

Custom Assembly

We produce standard components to fit your current assemblies and devices, but we also provide custom solutions by applying a variety of assembly options, such as joining fittings, connectors, filters, and other components to tubing. From systems using a single fitting to complex configurations utilizing various components, we can design complete assemblies to your requirements.