

# TUOHY BORST

## Adapters



Adult Without Side Port  
**PG-57**



Pediatric Without Side Port  
**PG-57**



Adult With Side Port  
**PG-57**



## TBASeries

### Tuohy Borst Adapters

Tuohy Borst adapters are designed to facilitate catheter introduction and are intended for interventional and diagnostic procedures where fluid backflow is a concern.

Tuohy Borst adapters are available in convenient bag sizes of 100 and 25 with larger sizes attainable on request and lower minimum order quantities to suit your needs.

#### Features & Benefits:

- Male luer terminations that comply with EN/ISO standards
- With and without side ports for delivering medication
- Cap does not remove at fully opened position
- Leak proof seal during interventional and diagnostic procedures
- Materials: clear PVC, silicone adult gaskets, and Kraton pediatric gaskets
- Gasket through hole ID from 0.025" to 0.100"
- Different cap color options

#### Part Drawings Online!

Nordson MEDICAL provides dimensional drawings available for download for all of its fittings.  
[www.nordsonmedical.com](http://www.nordsonmedical.com)

#### Need it Fast?

You can order any of our 4,000+ parts online now with no order \$ minimums, 24 hours a day, 7 days a week!  
[www.nordsonmedical.com](http://www.nordsonmedical.com)