

TBA SERIES: Adult Without Side Port


TBA2025-001

Tuohy Borst Adapter without Side Port, 9 Fr Clear PVC Body, Blue Triax Cap, 0.025" ID Silicone Gasket



TBA2025-002

Tuohy Borst Adapter without Side Port, 9 Fr Clear PVC Body, Clear PVC Cap, 0.025" ID Silicone Gasket



TBA2040-001

Tuohy Borst Adapter without Side Port, Clear PVC Body, Blue Triax Cap, 0.040" ID Silicone Gasket



TBA2052-001

Tuohy Borst Adapter without Side Port, Clear PVC Body, Blue Triax Cap, 0.052" ID Silicone Gasket



TBA2062-001

Tuohy Borst Adapter without Side Port, Clear PVC Body, Blue Triax Cap, 0.062" ID Silicone Gasket



TBA2081-001

Tuohy Borst Adapter without Side Port, Clear PVC Body, Blue Triax Cap, 0.081" ID Silicone Gasket



TBA2100-001

Tuohy Borst Adapter without Side Port, Clear PVC Body, Clear PVC Cap, 0.100" ID Silicone Gasket



TBA SERIES: Pediatric Without Side Port

TBA3084-001

Pediatric Tuohy Borst Adapter without Side Port, Clear PVC Body, Red Cap, 0.084" ID Natural Kraton Gasket



TBA SERIES: Adult With Side Port

TBA4030P-001

Tuohy Borst Adapter with Side Port, Clear PVC Body, Blue Triax Cap, Pellethane Tubing, 0.030" ID Silicone Gasket, Acetal Mini-Clamp, Clear PVC Female Luer, and SS Washer



TBA4030P-002

Tuohy Borst Adapter with Side Port, Clear PVC Body, Blue Triax Cap, 0.030" ID Silicone Gasket, Pellethane Tubing, Clear PVC Female Luer, Acetal Mini Clamp, HDPE and Teflon Washers



TUOHY BORST
ADAPTERS

Looking for Something Different?

Most of our best product ideas have come from you, our customers. While we pride ourselves on our extensive product offering, we recognize that your needs are ever-changing and that our products, too, must adapt to your needs.

If there's a new product you think Nordson MEDICAL should make, please send us an Email at designcenter@nordsonmedical.com describing your product concept in as much detail as you desire. We will contact you with any questions, and who knows - the next catalog edition could include products YOU told us to build.

See page 14 for more information on the Design Center!

