

# LUER

## Fittings

Male Luers  
**PG-59**



Female Luers  
**PG-61**



Adapters  
**PG-63**



Couplers  
**PG-63**



Tees  
**PG-63**



Elbows  
**PG-63**



Plugs and Caps  
**PG-64**



Luer to Thread  
**PG-65**



Panel Mounts & Accessories  
**PG-66**



### Our Luer Tapers

Nordson MEDICAL's luer tapers are made to the ISO 594-1 and ISO 594-2 standards, and may be connected to other luers meeting these ISO standards. Please test your application to verify the performance of the interconnection.

Nordson MEDICAL's luer fittings come in over 100 configurations and several materials for the device designer, joining tubing from 1/16" (1.6 mm) ID to 1/4" (6.4 mm) ID.

### Features & Benefits:

#### Published ISO Design Standard

- Well understood worldwide
- Often assures interconnections without need to retrofit

#### Rotating Lock Permits Tube Positioning Prior To Tightening

- Minimizes twisting and kinking of tubing

#### Tapered Cone Seal

- Reliable
- Inexpensive
- Uncomplicated

#### Available In Many Materials

- Facilitates optimal material selection for your application

#### Inexpensive

- Facilitates single-use

#### Easy To Understand

- Ease of use in field

#### Many Configuration Options

- Adaptable to wide range of applications

#### Security Options (Fixed Lock, Rotating Lock, No Lock)

- Tailorable to your specific application

### Custom Assembly

We produce standard components to fit your current assemblies and devices, but we also provide custom solutions by applying a variety of assembly options, such as joining fittings, connectors, filters, and other components to tubing. From systems using a single fitting to complex configurations utilizing various components, we can design complete assemblies to your requirements.