

FEMALE LUER:  
Thread Style

**FTLL240 (-1 -2 -3 -5 -6005 -6012 -9 -9002 -J1A)**  
Female Luer Thread Style to 200 Series Barb, 5/32" (4.0 mm) ID Tubing



**FTLL025 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002)**  
Female Luer Thread Style to 500 Series Barb, 5/32" (4.0 mm) ID Tubing



**FTLL445 (-1 -6005 -9 -9002 -J1A)**  
Female Luer Thread Style to 400 Series Barb, .170" (4.3 mm) ID Tubing



**FTLL250 (-1 -3 -5 -6005 -9 -9002 -J1A)**  
Female Luer Thread Style to 200 Series Barb, 3/16" (4.8 mm) ID Tubing



**FTLL450 (-1 -6005 -9 -9002 -J1A)**  
Female Luer Thread Style to 400 Series Barb, 3/16" (4.8 mm) ID Tubing



**FTLL035 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002 -J1A)**  
Female Luer Thread Style to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



**FTLL360 (-1 -3 -5 -6005 -6012 -9 -9002 -J1A)**  
Female Luer Thread Style to 300 Series Barb, 1/4" (6.4 mm) ID Tubing



**FTLL055 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002)**  
Female Luer Thread Style to 500 Series Barb, 1/4" (6.4 mm) ID Tubing



FEMALE LUER:  
Smooth Bore Thread Style

**FTLLSB004 (-1 -6005 -9)**  
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/16" (1.6 mm) ID Tubing



**FTLLSB007 (-1 -4 -6005 -9)**  
Female Luer Thread Style Smooth Bore to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



**FTLLSB013 (-1 -6005 -9)**  
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



Why Choose a Smooth Bore?

Our Smooth Bore Luers have a smooth inner transition providing for easy insertion of guide wires, probes and other devices.

FEMALE LUER:  
Lug Style

**FTL10 (-1 -6005 -9 -J1A)**  
Female Luer Lug Style to Classic Series Barb, 1/16" (1.6 mm) ID Tubing



**FTL210 (-1 -6005 -9 -J1A)**  
Female Luer Lug Style to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



**FTL20 (-1 -6005 -07 -81 -9 -J1A)**  
Female Luer Lug Style to Classic Series Barb, 3/32" (2.4 mm) ID Tubing



**FTL220 (-1 -3 -5 -6005 -6012 -9 -J1A)**  
Female Luer Lug Style to 200 Series Barb, 3/32" (2.4 mm) ID Tubing



**FTL30 (-1 -6005 -6012 -9 -J1A)**  
Female Luer Lug Style to Classic Series Barb, 1/8" (3.2 mm) ID Tubing

