

## LUER TO THREAD:

Male

### **ABSML-6005 UNDER (-1 -2 -3 -4 -5 -06 -07 -6012)**

Male Luer with Rotating Lock Ring and 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread



### **ABSML-J1A UNDER (-1 -2 -3 -4 -5 -06 -07)**

Male Luer with Rotating Lock Ring and 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread



### **XMTLL (-1 -6005 -9 -J1A)**

Male Luer Integral Lock Ring to 10-32 Special Tapered Thread



### **SMTLL (-1 -6005 -6012 -9 -J1A)**

Male Luer Integral Lock Ring to 1/4-28 UNF Thread



### **18MTLL (-1 -2 -6005 -6012 -J1A)**

Male Luer Integral Lock Ring with 7/16" Hex to 1/8-27 NPT Thread



### **14MTLL (-1 -2 -6005 -6012)**

Male Luer Integral Lock Ring with 9/16" Hex to 1/4-18 NPT Thread



### **XMLRL (-1 -6005 -9 -9002)**

Male Luer with 1/4" Hex to 10-32 Special Tapered Thread (May be used with separate rotating lock ring; FSLLR)



### **SMLRL (-1 -6005 -9 -J1A)**

Male Luer with 5/16" Hex to 1/4-28 UNF Thread (May be used with separate rotating lock ring; FSLLR)



### **BSMLRL (-6005 -J1A)**

Male Luer with 5/16" Hex to 1/4-28 Bottom-Sealing Thread (May be used with separate rotating lock ring; FSLLR)



LUER  
FITTINGS

## Looking for Something Different?

Most of our best product ideas have come from you, our customers. While we pride ourselves on our extensive product offering, we recognize that your needs are ever-changing and that our products, too, must adapt to your needs.

If there's a new product you think Nordson MEDICAL should make, please send us an Email at **designcenter@nordsonmedical.com** describing your product concept in as much detail as you desire. We will contact you with any questions, and who knows - the next catalog edition could include products YOU told us to build.

See page 14 for more information on the Design Center!



## Need More Details?

Part drawings, minimum order quantities, and CAD requests can be found online. Visit **www.nordsonmedical.com** or contact Sales at **sales@nordsonmedical.com**.