

# BLOOD PRESSURE

## Fittings



### BPF Series Female Bayonet Style Quick Connects PG-68



### Bayonet Male Cuff Connectors PG-69



### Bayonet Panel Mount Connectors PG-69



### Dinamap® Monitor Adapters PG-69



### Dinamap® Cuff Connectors PG-70



### BPL Series 1/4 Turn Cuff Connectors PG-71



## BPF Series

### Blood Pressure Fittings

The BPF Series features an ergonomic design for simple one-handed operation. Easy to connect and disconnect to existing tubing, many locations are replacing their metal bayonet connectors with these light weight connectors. The BPF Series also works well with the new ISO 80369-5 limb cuff inflation standard compliant parts. The standard defines the bayonet as the only approved design.

#### BPF SERIES:

#### Female Bayonet Style Quick Connect

##### BPF220 (-001) Black ABS Housing

Bayonet-Style Female Quick Connect Fitting, 200 Series Barb, 3/32" (2.4 mm) ID Tubing



works with...

<b>BPM013</b>	<b>BPMB230</b>
<b>BPM240</b>	<b>BPMB240</b>
<b>BPM025</b>	<b>BPMB250</b>
<b>BPM035</b>	<b>BPAM</b>
<b>BPMMTLL</b>	

##### BPF230 (-001) Black ABS Housing

Bayonet-Style Female Quick Connect Fitting, 200 Series Barb, 1/8" (3.2 mm) ID Tubing



works with...

<b>BPM013</b>	<b>BPMB230</b>
<b>BPM240</b>	<b>BPMB240</b>
<b>BPM025</b>	<b>BPMB250</b>
<b>BPM035</b>	<b>BPAM</b>
<b>BPMMTLL</b>	

##### BPF240 (-001) Black ABS Housing

Bayonet-Style Female Quick Connect Fitting, 200 Series Barb, 5/32" (4.0 mm) ID Tubing



works with...

<b>BPM013</b>	<b>BPMB230</b>
<b>BPM240</b>	<b>BPMB240</b>
<b>BPM025</b>	<b>BPMB250</b>
<b>BPM035</b>	<b>BPAM</b>
<b>BPMMTLL</b>	

##### BPF250 (-001) Black ABS Housing

Bayonet-Style Female Quick Connect Fitting, 200 Series Barb, 3/16" (4.8 mm) ID Tubing



works with...

<b>BPM013</b>	<b>BPMB230</b>
<b>BPM240</b>	<b>BPMB240</b>
<b>BPM025</b>	<b>BPMB250</b>
<b>BPM035</b>	<b>BPAM</b>
<b>BPMMTLL</b>	