

**LARGE TUBE-TO-TUBE:**  
Elbow Reducer



Reference part shown: **L070/055**

Due to size we are not able to display all large Tube-To-Tube connectors. Please see [www.nordsonmedical.com](http://www.nordsonmedical.com) for part reference. Minimum package quantities vary by tube IDs, please verify minimum package sizes when placing orders.

**L070/055 (-0 -6005)**

Elbow Reduction Connector with 500 Series Barbs, 3/8" (9.5 mm) and 1/4" (6.4 mm) ID Tubing

**L080/070 (-0 -6005)**

Elbow Reduction Connector with 500 Series Barbs, 1/2" (12.7 mm) and 3/8" (9.5 mm) ID Tubing

**L100/080 (-0 -6005)**

Elbow Reduction Connector with 500 Series Barbs, 3/4" (19.0 mm) and 1/2" (12.7 mm) ID Tubing

**L110/100 (-0 -6005)**

Elbow Reduction Connector with 500 Series Barbs, 1" (25.4 mm) and 3/4" (19.0 mm) ID Tubing

**TUBE-TO-TUBE:**  
Tee

**T10 (-1 -2 -6005 -J1A)**

Tee Connector with Classic Series Barbs, 1/16" (1.6 mm) ID Tubing



**T210 (-1 -2 -6005 -9 -J1A)**

Tee Connector with 200 Series Barbs, 1/16" (1.6 mm) ID Tubing



**T410 (-1 -2 -6005 -J1A)**

Tee Connector with 400 Series Barbs, 1/16" (1.6 mm) ID Tubing



**T20 (-1 -2 -6005 -J1A)**

Tee Connector with Classic Series Barbs, 3/32" (2.4 mm) ID Tubing



**T220 (-1 -2 -6005 -6012 -9 -J1A)**

Tee Connector with 200 Series Barbs, 3/32" (2.4 mm) ID Tubing



**TUBE-TO-TUBE:**  
Tee

**T420 (-1 -2 -6005 -J1A)**

Tee Connector with 400 Series Barbs, 3/32" (2.4 mm) ID Tubing



**T30 (-1 -2 -6005 -9 -J1A)**

Tee Connector with Classic Series Barbs, 1/8" (3.2 mm) ID Tubing



**T230 (-1 -2 -5 -6005 -6012 -9 -J1A)**

Tee Connector with 200 Series Barbs, 1/8" (3.2 mm) ID Tubing



**T430 (-1 -2 -6005 -9002 -J1A)**

Tee Connector with 400 Series Barbs, 1/8" (3.2 mm) ID Tubing



**T013 (-6005)**

Tee Connector with 500 Series Barbs, 1/8" (3.2 mm) ID Tubing



**T40 (-1 -2 -6005 -J1A)**

Tee Connector with Classic Series Barbs, 5/32" (4 mm) ID Tubing



**T250 (-1 -2 -6005 -9 -J1A)**

Tee Connector with 200 Series Barbs, 3/16" (4.8 mm) ID Tubing



**T60 (-1 -2 -6005 -9 -J1A)**

Tee Connector with Classic Series Barbs, 1/4" (6.4 mm) ID Tubing



TUBE-TO-TUBE  
FITTINGS