

BUILD-A-PART®:
Accessories, Barbs and Elbows

RB01 (-X2) SB01 (-X0)

Check Valve Ball 5/32" Diameter



BDMRP (-40 -81)

Bondable Plug



BDMD (-40 -81)

Bondable Mounting Device



BDMFTLL (-40 -81 -9)

Bondable Female Luer Thread Style




BDMMTL (-40 -81)

Bondable Male Luer (May be used with separate rotating lock ring; RMLLR)



RMLLR (-1 -2 -3 -4 -5 -06 -07 -6005 -9 -9002 -J1A)

Rotating Male Luer Lock Ring (For use with BDMMTL)



BDMR210 (-40 -81 -9 -9002)

Bondable 200 Series Barb, 1/16" (1.6 mm) ID Tubing



BDMR410 (-40 -81)

Bondable 400 Series Barb, 1/16" (1.6 mm) ID Tubing



BDMR220 (-40 -81)

Bondable 200 Series Barb, 3/32" (2.4 mm) ID Tubing



BDMR420 (-40 -81)

Bondable 400 Series Barb, 3/32" (2.4 mm) ID Tubing



BDMR230 (-40 -81)

Bondable 200 Series Barb, 1/8" (3.2 mm) ID Tubing



BDMR430 (-40 -81)

Bondable 400 Series Barb, 1/8" (3.2 mm) ID Tubing



BDMR240 (-40 -6005 -81)

Bondable 200 Series Barb, 5/32" (4.0 mm) ID Tubing



BDMR445 (-40 -81)

Bondable 400 Series Barb, .170" (4.3 mm) ID Tubing



BDML210 (-40 -81 -9)

Bondable Elbow 200 Series Barb, 1/16" (1.6 mm) ID Tubing



BDML410 (-40 -81)

Bondable Elbow 400 Series Barb, 1/16" (1.6 mm) ID Tubing



BDML220 (-40 -81 -9)

Bondable Elbow 200 Series Barb, 3/32" (2.4 mm) ID Tubing



BDML420 (-40 -81)

Bondable Elbow 400 Series Barb, 3/32" (2.4 mm) ID Tubing



BDML230 (-40 -81)

Bondable Elbow 200 Series Barb, 1/8" (3.2 mm) ID Tubing



BDML430 (-40 -81)

Bondable Elbow 400 Series Barb, 1/8" (3.2 mm) ID Tubing



BUILD-A-PART®