

## SMALL BORE CONNECTORS:

Female

### TSC2F210 (-6005)

Female Taper Seal Non-Luer Connector to 200 Series Barb,  
1/16" (1.6 mm) ID Tubing



### TSC2F007 (-6005)

Female Taper Seal Non-Luer Connector to 500 Series Barb,  
3/32" (2.4 mm) ID Tubing



### TSC2F013 (-6005)

Female Taper Seal Non-Luer Connector to 500 Series Barb,  
1/8" (3.2 mm) ID Tubing



### TSC2F035 (-6005)

Female Taper Seal Non-Luer Connector to 500 Series Barb,  
3/16" (4.8 mm) ID Tubing



### TSC2FB230 (-6005)

Female Taper Seal Non-Luer with 1/4-28 UNF Panel Mount to  
200 Series Barb, 1/8" (3.2 mm) ID Tubing



## Small BORE CONNECTORS:

Plugs and Caps

### TSC2MP (-6005) 1,000 piece minimum

Male Taper Seal Non-Luer Plug with Tether Recess



## LARGE BORE CONNECTORS:

Male

### MLB035 (-1 -8012 -9)

Male Large Bore Non-Luer Connector to 500 Series Barb,  
3/16" (4.8 mm) ID Tubing




### MLB055 (-1 -8012 -9)

Male Large Bore Non-Luer Connector to 500 Series Barb,  
1/4" (6.4 mm) ID Tubing



### MLB065 (-1 -8012 -9)

Male Large Bore Non-Luer Connector to 500 Series Barb,  
5/16" (8.0 mm) ID Tubing




### MLBC7 (-9002)

Male Large Bore Non-Luer Connector with Integral Lock Ring to Male Slip Luer  
(May be used with separate rotating lock ring; FSLLR)



### MLBC8 (-9002)

Male Large Bore Non-Luer Connector with Integral Lock Ring to Female  
Thread-Style Locking Luer



## Need it Fast?

You can order any of our 4,000+ parts online now with no  
order \$ minimums, 24 hours a day, 7 days a week!

[www.nordsonmedical.com](http://www.nordsonmedical.com)

OTHER  
COMPONENTS