

Material - Manufacturer & Grade		VP Code	Regulatory Compliance of Base Resin*		RoHS Compliant	BPA Free	DEHP Free	Animal Derivative Free*
ABS (Lace White) - INEOS Lustran 348-000000		-81		USP Class VI Certified; ISO 10993 Part I Compliant	Yes	Yes	Yes	No
MABS (Methyl Methacrylate) - BASF Terlux 2802 TR		-8003, SCV Series		USP Class VI Certified	Yes	Yes	Yes	Yes
(Sno White) - INEOS Lustran 348-012002		-8012		USP Class VI Certified; ISO 10993 Part I Compliant	Yes	Yes	Yes	No
tal - DuPont Delrin 500P		-10		21 CFR 177.2480 FDA Approved	Yes	Yes	Yes	No
tal Copolymer - BASF Ultraform N 23020 003		-1006, -1007		21 CFR 177.2470 FDA Approved	Yes	Yes	Yes	No
ilic - Cyro Cyrolite Med 2		-50		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1010 FDA Approved	Yes	Yes	Yes	Yes
n (PA 66) - DuPont Zytel 101F		-0, -1, -2, -3, -4, -5, -06, -07, -08B, -8, -8B, -CM030		21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No
n (PA 66) - DuPont Zytel SC310 NC010		-U1, -U5, -U01, -U05		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No
n (Glass Reinforced) - DuPont Zytel FGF5171		-U71, -U75		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No
n (Glass Reinforced) - DuPont Zytel 70G33L		-7, -72		21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No
carbonate (Clear) - Bayer Makrolon 2558 550115		-9		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1580 FDA Approved	Yes	No	Yes	No
carbonate (Clear) - Bayer Makrolon Rx 1805 451118		-9010, SCV Series		USP Class VI Certified; ISO 10993 Compliant	Yes	No	Yes	Yes
carbonate (Tinted) - Sabic Lexan HPS1		-9024		USP Class VI Certified; ISO 10993 Compliant	Yes	No	Yes	No
carbonate (Tinted) - Styron Calibre Megarad 2081-15		-9002, DCV Series		ISO 10993 Compliant	Yes	No	Yes	No
ethylene - LyondellBasell Purell PE 1840 H		SCV Series		Please call for details	Yes	Yes	Yes	Yes
ethylene - Trilliant HC5420-0002 LL Natural Resin		-VP1		USP Class VI Certified	Yes	Yes	Yes	Yes
propylene - Flint Hills P5M6K-080		-6005, -6006		USP Class VI Certified; 21 CFR 177.1520 FDA Approved	Yes	Yes	Yes	Yes
Polysulfone - Solvay Udel P1700 QL2611		-40		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1655 FDA Approved	Yes	Yes	Yes	Yes
PVDF - Arkema Kynar 1000 HD		-J1A		USP Class VI Certified; 21 CFR 177.2510 FDA Approved	Yes	Yes	Yes	Yes
SAN - BASF 388S TR 77770		SCV Series		ISO 10993 Compliant	Yes	Yes	Yes	No
Silicone - Wacker Elastosil LR 3003		SCV Series		USP Class VI Certified	Yes	Yes	Yes	Yes

*As determined by resin manufacturer