



| Cat. No. | Filter Medium | Pore Size(μm) | Diameter (mm) | Gamma Sterile | Qty/ pack |
|---------------------------------------|---------------------|---------------|---------------|---------------|-----------|
| Sterile Syringe Filter With Prefilter | | | | | |
| SFPES013022SG | PES / GF Prefilter | 0.22 | 13 | Yes | 100 |
| SFPES013045SG | PES / GF Prefilter | 0.45 | 13 | Yes | 100 |
| SFPES025022SG | PES / GF Prefilter | 0.22 | 25 | Yes | 100 |
| SFPES025045SG | PES / GF Prefilter | 0.45 | 25 | Yes | 100 |
| SFPVDF013022SG | PVDF/ GF Prefilter | 0.22 | 13 | Yes | 100 |
| SFPVDF013045SG | PVDF/ GF Prefilter | 0.45 | 13 | Yes | 100 |
| SFPVDF025022SG | PVDF/ GF Prefilter | 0.22 | 25 | Yes | 100 |
| SFPVDF025045SG | PVDF/ GF Prefilter | 0.45 | 25 | Yes | 100 |
| SFNY013022SG | Nylon/ GF Prefilter | 0.22 | 13 | Yes | 100 |
| SFNY013045SG | Nylon/ GF Prefilter | 0.45 | 13 | Yes | 100 |
| SFNY025022SG | Nylon/ GF Prefilter | 0.22 | 25 | Yes | 100 |
| SFNY025045SG | Nylon/ GF Prefilter | 0.45 | 25 | Yes | 100 |

MS® Nylon Syringe Filter (Polyamide)

- Hydrophilic property
- No need to moist beforehand
- Uniform aperture
- Strong tenacity and adsorbability
- Compatible with aqueous and alcoholic solutions and solvents
- Suitable for HPLC



| Cat. No. | Filter Medium | Pore Size(μm) | Diameter (mm) | Gamma Sterile | Qty/ pack |
|-------------|---------------|---------------|---------------|---------------|-----------|
| SFNY013022N | Nylon | 0.22 | 13 | No | 100 |
| SFNY013045N | Nylon | 0.45 | 13 | No | 100 |
| SFNY025022N | Nylon | 0.22 | 25 | No | 100 |
| SFNY025045N | Nylon | 0.45 | 25 | No | 100 |
| SFNY033022N | Nylon | 0.22 | 33 | No | 100 |
| SFNY033045N | Nylon | 0.45 | 33 | No | 100 |

Note: Membrane Solutions offer Customized Pore Sizes: 0.1μm, 0.8μm, 1.0μm, 3.0μm, 5.0μm

MS® MCE Syringe Filter (Mixed Cellulose Ester)

- Uniform aperture
- No fiber migration
- Naturally hydrophilic



| Cat. No. | Filter Medium | Pore Size(μm) | Diameter (mm) | Gamma Sterile | Qty/ pack |
|--------------|---------------|---------------|---------------|---------------|-----------|
| SFMCE013022N | MCE | 0.22 | 13 | No | 100 |
| SFMCE013045N | MCE | 0.45 | 13 | No | 100 |
| SFMCE025022N | MCE | 0.22 | 25 | No | 100 |
| SFMCE025045N | MCE | 0.45 | 25 | No | 100 |
| SFMCE033022N | MCE | 0.22 | 33 | No | 100 |
| SFMCE033045N | MCE | 0.45 | 33 | No | 100 |

Note: Membrane Solutions offer Customized Pore Sizes: 0.8μm, 1.0μm, 3.0μm, 5.0μm