

Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFNY037300	Nylon	3.0	37	100
MFNY047300	Nylon	3.0	47	50
MFNY090300	Nylon	3.0	90	25
MFNY142300	Nylon	3.0	142	25
MFNY293300	Nylon	3.0	293	25
MFNY013500	Nylon	5.0	13	200
MFNY025500	Nylon	5.0	25	100
MFNY037500	Nylon	5.0	37	100
MFNY047500	Nylon	5.0	47	50
MFNY090500	Nylon	5.0	90	25
MFNY142500	Nylon	5.0	142	25
MFNY293500	Nylon	5.0	293	25

MS® PTFE Membrane Filter

- PTFE membrane with supporting layer polyester or polypropylene
- Suitable for applications involving aggressive organic solvents, strong acids, and alkalis
- Hydrophobic nature of the membrane has applications for air and gas sterilization
- High temperature resistance



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFPTFE013010B	Hydrophobic PTFE	0.1	13	200
MFPTFE025010B	Hydrophobic PTFE	0.1	25	100
MFPTFE037010B	Hydrophobic PTFE	0.1	37	100
MFPTFE047010B	Hydrophobic PTFE	0.1	47	50
MFPTFE090010B	Hydrophobic PTFE	0.1	90	25
MFPTFE142010B	Hydrophobic PTFE	0.1	142	25
MFPTFE293010B	Hydrophobic PTFE	0.1	293	25
MFPTFE013022B	Hydrophobic PTFE	0.22	13	200
MFPTFE025022B	Hydrophobic PTFE	0.22	25	100
MFPTFE037022B	Hydrophobic PTFE	0.22	37	100
MFPTFE047022B	Hydrophobic PTFE	0.22	47	50
MFPTFE090022B	Hydrophobic PTFE	0.22	90	25
MFPTFE142022B	Hydrophobic PTFE	0.22	142	25
MFPTFE293022B	Hydrophobic PTFE	0.22	293	25
MFPTFE013045B	Hydrophobic PTFE	0.45	13	200
MFPTFE025045B	Hydrophobic PTFE	0.45	25	100
MFPTFE037045B	Hydrophobic PTFE	0.45	37	100
MFPTFE047045B	Hydrophobic PTFE	0.45	47	50
MFPTFE090045B	Hydrophobic PTFE	0.45	90	25
MFPTFE142045B	Hydrophobic PTFE	0.45	142	25
MFPTFE293045B	Hydrophobic PTFE	0.45	293	25
MFPTFE013100B	Hydrophobic PTFE	1.0	13	200
MFPTFE025100B	Hydrophobic PTFE	1.0	25	100