



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFPTFE037100B	Hydrophobic PTFE	1.0	37	100
MFPTFE047100B	Hydrophobic PTFE	1.0	47	50
MFPTFE090100B	Hydrophobic PTFE	1.0	90	25
MFPTFE142100B	Hydrophobic PTFE	1.0	142	25
MFPTFE293100B	Hydrophobic PTFE	1.0	293	25
MFPTFE013300B	Hydrophobic PTFE	3.0	13	200
MFPTFE025300B	Hydrophobic PTFE	3.0	25	100
MFPTFE037300B	Hydrophobic PTFE	3.0	37	100
MFPTFE047300B	Hydrophobic PTFE	3.0	47	50
MFPTFE090300B	Hydrophobic PTFE	3.0	90	25
MFPTFE142300B	Hydrophobic PTFE	3.0	142	25
MFPTFE293300B	Hydrophobic PTFE	3.0	293	25
MFPTFE013500B	Hydrophobic PTFE	5.0	13	200
MFPTFE025500B	Hydrophobic PTFE	5.0	25	100
MFPTFE037500B	Hydrophobic PTFE	5.0	37	100
MFPTFE047500B	Hydrophobic PTFE	5.0	47	50
MFPTFE090500B	Hydrophobic PTFE	5.0	90	25
MFPTFE142500B	Hydrophobic PTFE	5.0	142	25
MFPTFE293500B	Hydrophobic PTFE	5.0	293	25

MS® PES Membrane Filter

- Inherently hydrophilic
- Low protein binding
- Good chemical compatibility
- Superior thermo stability



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFPE013005	PES	0.05	13	200
MFPE025005	PES	0.05	25	100
MFPE037005	PES	0.05	37	100
MFPE047005	PES	0.05	47	50
MFPE090005	PES	0.05	90	25
MFPE142005	PES	0.05	142	25
MFPE013010	PES	0.1	13	200
MFPE025010	PES	0.1	25	100
MFPE037010	PES	0.1	37	100
MFPE047010	PES	0.1	47	50
MFPE090010	PES	0.1	90	25
MFPE142010	PES	0.1	142	25