



# Lab Filtration

## MS® Microbiology Test Membrane Filter(MCE)

- Available in white and black
- Gridded or non-gridded
- Individual pack, pre-sterilized
- Standard for microbiological analysis of water, waste water, and beverages.



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
<b>MCE Membrane Filter, White, Gridded, Sterile, Individually packed</b>				
MFMC025022GWS	MCE	0.22	25	100
MFMC037022GWS	MCE	0.22	37	100
MFMC047022GWS	MCE	0.22	47	100
MFMC025045GWS	MCE	0.45	25	100
MFMC037045GWS	MCE	0.45	37	100
MFMC047045GWS	MCE	0.45	47	100
MFMC025080GWS	MCE	0.8	25	100
MFMC037080GWS	MCE	0.8	37	100
MFMC047080GWS	MCE	0.8	47	100
MFMC025100GWS	MCE	1.0	25	100
MFMC037100GWS	MCE	1.0	37	100
MFMC047100GWS	MCE	1.0	47	100
MFMC025300GWS	MCE	3.0	25	100
MFMC037300GWS	MCE	3.0	37	100
MFMC047300GWS	MCE	3.0	47	100
<b>MCE Membrane Filter, Black, Gridded, Sterile, Individually packed</b>				
MFMC025045GBS	MCE	0.45	25	100
MFMC037045GBS	MCE	0.45	37	100
MFMC047045GBS	MCE	0.45	47	100
MFMC025080GBS	MCE	0.8	25	100
MFMC037080GBS	MCE	0.8	37	100
MFMC047080GBS	MCE	0.8	47	100