

Lab Filtration

MS® Absorbent Pad

- Pure cellulose pad will not inhibit any bacterial grow.
- Available non-sterile or sterilized by gamma irradiation.
- Individually packed for pre-sterilized
- Customized diameter is available



Cat. No.	Description	Diameter (mm)	Qty/ pack
MFAP047N	Non-sterile Absorbent Pad	47	100
MFAP047S	Sterile Absorbent Pad, individually packed	47	100
MFAP050N	Non-sterile Absorbent Pad	50	100
MFAP050S	Sterile Absorbent Pad, individually packed	50	100
MFAP015050D	Petri-Pad dish, 55×15mm, with absorbent pad, sterile	50	200

MS® PVDF Membrane Filter

- Wide chemical compatibility
- Excellent mechanical properties
- High temperature capabilities
- Low extractable levels



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFPVDF013022	PVDF	0.22	13	200
MFPVDF025022	PVDF	0.22	25	100
MFPVDF037022	PVDF	0.22	37	100
MFPVDF047022	PVDF	0.22	47	50
MFPVDF090022	PVDF	0.22	90	25
MFPVDF142022	PVDF	0.22	142	25
MFPVDF293022	PVDF	0.22	293	25
MFPVDF013045	PVDF	0.45	13	200
MFPVDF025045	PVDF	0.45	25	100
MFPVDF037045	PVDF	0.45	37	100
MFPVDF047045	PVDF	0.45	47	50
MFPVDF090045	PVDF	0.45	90	25
MFPVDF142045	PVDF	0.45	142	25
MFPVDF293045	PVDF	0.45	293	25