


MS® Polycarbonate (PC) Membrane Filter

- Precise pore sizes and pore distribution for absolute filtration and separation
- Excellent chemical resistance and thermal stability
- Smooth glass-like surface with cylindrical pores for maximum particulate capture
- Optically transparent in most pore sizes



Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFPC013005	PC	0.05	13	100
MFPC025005	PC	0.05	25	100
MFPC037005	PC	0.05	37	100
MFPC047005	PC	0.05	47	100
MFPC090005	PC	0.05	90	30
MFPC142005	PC	0.05	142	25
MFPC293005	PC	0.05	293	25
MFPC013010	PC	0.1	13	100
MFPC025010	PC	0.1	25	100
MFPC037010	PC	0.1	37	100
MFPC047010	PC	0.1	47	100
MFPC090010	PC	0.1	90	30
MFPC142010	PC	0.1	142	25
MFPC293010	PC	0.1	293	25
MFPC013020	PC	0.2	13	100
MFPC025020	PC	0.2	25	100
MFPC037020	PC	0.2	37	100
MFPC047020	PC	0.2	47	100
MFPC090020	PC	0.2	90	30
MFPC142020	PC	0.2	142	25
MFPC293020	PC	0.2	293	25
MFPC013040	PC	0.4	13	100
MFPC025040	PC	0.4	25	100
MFPC037040	PC	0.4	37	100
MFPC047040	PC	0.4	47	100
MFPC090040	PC	0.4	90	30
MFPC142040	PC	0.4	142	25
MFPC293040	PC	0.4	293	25
MFPC013080	PC	0.8	13	100
MFPC025080	PC	0.8	25	100
MFPC037080	PC	0.8	37	100
MFPC047080	PC	0.8	47	100