


## MS® Regenerated Cellulose(RC) Membrane Filter

- Hydrophilic
- Easily wettable
- Resistant to almost all solvents and aqueous solutions in pH range 3-12
- Low non-specific of adsorption



### Note:

- Other sizes available upon request

Cat. No.	Filter Medium	Pore Size (µm)	Diameter (mm)	Qty/pack
MFRC013022	RC	0.22	13	100
MFRC025022	RC	0.22	25	100
MFRC037022	RC	0.22	37	100
MFRC047022	RC	0.22	47	100
MFRC090022	RC	0.22	90	25
MFRC142022	RC	0.22	142	25
MFRC013045	RC	0.45	13	100
MFRC025045	RC	0.45	25	100
MFRC037045	RC	0.45	37	100
MFRC047045	RC	0.45	47	100
MFRC090045	RC	0.45	90	25
MFRC142045	RC	0.45	142	25