



Cat. No.	Description	Qty/pack
VFG030010	1) 300ml Glass Funnel with cover	1set
	2) 1000ml Glass Solvent Collection	
	3) Stainless Steel Clamp	
	4) Solid Glass Frit	
	5) Packaging: 30cmX30cmX25cm	
VFG025020	1) 250ml Glass Funnel with cover	1set
	2) 2000ml Glass Solvent Collection	
	3) Stainless Steel Clamp	
	4) Solid Glass Frit	
	5) Packaging: 30cmX30cmX25cm	
VFG030020	1) 300ml Glass Funnel with cover	1set
	2) 2000ml Glass Solvent Collection	
	3) Stainless Steel Clamp	
	4) Solid Glass Frit	
	5) Packaging: 30cmX30cmX25cm	
VFG147SS	Pyrex glass support screen	1set
VFG047SS	Pyrex glass support screen	1set
VFG010RB	Solvent Collection Bottle, 1000ml	1set
VFG020RB	Solvent Collection Bottle, 2000ml	1set

MS® Disposable Plastic Vacuum Filters

- Funnel is manufactured from optically clear polystyrene and graduated.
- Receiver bottle and filter adapter are manufactured from 100% virgin polypropylene.
- Hose connector is designed to accept multiple hose diameter and features an easy gripping collar to simplify tightening/loosening and adjustment.
- Individually wrapped sterile, certified RNase-free, DNase-free, non-pyrogenic, and DNA free.



Funnel Capacity	Filter Diameter	Process Volume	Hold-up Volume after purge	Maximum Operating Temperature	Fitting Outlet (Thread)	Full Unit Overall Height	Housing Material
150mL	50mm	150mL	3mL	45	45mm	156mm	ABS
250mL	50mm	250mL	3mL	45	45mm	156mm	ABS
500mL	50mm	500mL	3mL	45	45mm	156mm	ABS