

2ml 12*32mm Vials, flat base, clear and amber - Nominal sizes Only - designed to fit common chromatography Autosamplers	
Standard Opening Crimp Top 8-425	<ol style="list-style-type: none"> 1. Narrow neck crimp top vials fit many older GC autosamplers. 2. Compatible with all 11mm crimp seals and 5mm nominal diameter inserts
Wide Opening Crimp Top 9-425 eg Agilent Autosamplers	<ol style="list-style-type: none"> 1. Easier to fill 2. Used with either 11mm aluminum seals
Standard Opening Screw Thread 8-425	Standard neck vial—8-425 thread finish, 12x32mm outer profile
Wide opening Screw Thread 9-425	Wide neck vial is easy to fill
11mm Snap	<ol style="list-style-type: none"> 1. Wide opening vial is easier to fill and provides a larger target area for the autosampler needle 2. Used with snap caps

4ml 15*45mm Vials, flat base

Screw Thread, Crimp Top, Snap Top Vials

Headspace Vials

Beveled Edge, Square Rim

Also Available with 10ml, 20ml, and 40ml vials.

More information please visit our website : www.chromtech.net.au .

Or mail : Chromalytic Technology Pty Ltd in AUSTRALIA at : info@chromtech.net.au

