

MS® Syringe Filters

MS® syringe filters are simply quality filters, well packaged, and offered at a fair and competitive price. The Classic range is available in all of the major membranes including Nylon, PTFE, PES, MCE and PVDF which are supplied in 13mm, 25mm and 33mm formats in virgin polypropylene housings.

The emphasis is very much on quality. Membrane materials are supplied by the best names in the industry and the ISO9000 certified manufacturing is carried out to the highest standards, in certified clean room conditions, using the latest manufacturing technology to ensure a high quality, consistent product.

All items are quality tests for filter efficacy and housing integrity. The housing is pressure tested for use with up to 75 psig (5.0 bar) of pressure. Designed with a Female Luer-Lok inlet and Male Luer slip outlets. Some Filters are individually wrapped sterile, certified RNase-free, DNase-free, Non-pyrogenic, and DNA-free.



Colour Coded:

MS® syringe filters are colour coded, providing **easy identification of the membrane type and porosity.**

Click on any of the filters below to view the range we stock with that membrane:


White -- PES
0.22µm, 0.45µm
13mm, 25mm, 33mm


Green -- MCE
0.22µm, 0.45µm
13mm, 25mm, 33mm


Purple -- PTFE
0.22µm, 0.45µm
13mm, 25mm, 33mm


Yellow -- Nylon
0.22µm, 0.45µm
13mm, 25mm, 33mm


Black -- PVDF
0.22µm, 0.45µm
13mm, 25mm, 33mm


Blue -- MCE
0.22µm, 0.45µm
13mm, 25mm, 33mm


Orange -- PTFE
0.22µm, 0.45µm
13mm, 25mm, 33mm


Pink -- Nylon
0.22µm, 0.45µm
13mm, 25mm, 33mm


Red -- PVDF
0.22µm, 0.45µm
13mm, 25mm, 33mm