


MS® PTFE Syringe Filter



Product Description:

MS® Syringe filters are purpose-built with features designed to bring the highest levels of performance and purity to your research. We incorporate a variety of membranes to offer separation and purification solutions for the majority of your laboratory needs.

Features and Benefits:

- Broad chemical compatibility
- Strong chemical stability and inertia
- Strong hydrophobicity
- Syringe Filters for Cell Culture provide effective filtration for a wide variety of sample types. They are available in two pore sizes (0.22µm and 0.45µm) and four different membrane types.
- All items are quality tested for filter efficacy and housing integrity. The housing is pressure tested for use with up to 75 psig (5.0 bar) of pressure
- Designed with a Female Luer-Lok inlet and Male Luer slip outlets.
- Some Filters are individually wrapped sterile, certified RNase-free, DNase-free,
- Non-pyrogenic, and DNA-free.

Application:

- Organic solvent with strong chemical causticity filtration
- Strong acid solvent filtration
- Alkali solvent filtration

Technical Specification:

Parameters	13mm		25mm		33mm	
Membrane material	PTFE		PTFE		PTFE	
Housing material	PP		PP		PP	
Filter diameter (mm)	13mm		25mm		33mm	
Filtration area (cm ²)	0.65		3.90		4.60	
Pore Size (µm)	0.22	0.45	0.22	0.45	0.22	0.45
Holdup volume (µl)	<10		<30		<55	
Sample volume (ml)	<12		<100		<140	
Maximum Operating Temperature	130°C		130°C		130°C	
Maximum Operating Pressure (psi)	130		130		130	
Applicable pH value	1-14		1-14		1-14	