

MS® PES Syringe Filter



Product Description:

MS® Syringe filters are purpose-built with features designed to bring the highest levels of performance and purity to your research. We incorporate a variety of membranes to offer separation and purification solutions for the majority of your laboratory needs. PES(Polyethersulfone) – low affinity for proteins and extractable with substantially faster flow rates than PVDF; suitable for pre-filtration and filtration of buffers and culture media.

Features and Benefits:

- High filtration speed
- Lowest protein binding
- Syringe Filters for Cell Culture provide effective filtration for a wide variety of sample types. They are available in two pore sizes(0.22µm and 0.45µm)
- All items are quality tests for filter efficacy and housing integrity. The housing is pressure tested for use with up to 75 psig (5.0 bar) of pressure
- Designed with a Female Luer-Lok inlet and Male Luer slip outlets.
- some Filters are individually wrapped sterile, certified RNase-free, DNase- free,
- Non-pyrogenic, and DNA –free.
- Low exeractables

Application:

- Sterile filtering protein solution
- Tissue culture media filtration
- Tissue culture additive filtration

Technical Specification:

Parameters	13mm		25mm		33mm	
Membrane material	PES		PES		PES	
Housing material	PP		PP		PP	
Filter diameter (mm)	13mm		25mm		33mm	
Filtration area (cm ²)	0.65		3.90		4.60	
Pore Size(µm)	0.22	0.45	0.22	0.45	0.22	0.45
Holdup volume (µl)	<10		<30		<55	
Sample volume (ml)	<12		<100		<140	
Maximum Operating Temperature	90°C		90°C		90°C	
Maximum Operating Pressure (psi)	50		95		120	
Applicable pH value	1-14		1-14		1-14	