

## *MS<sup>®</sup> SteriBio Syringe Filters*

*----- Ideal for Proteinaceous Samples and Tissue Culture Work*



MS<sup>®</sup> SteriBio syringe filters are available with Polyethersulphone (PES) and Cellulose Acetate (CA) membranes. Each filter is individually packed and sterilized by Gama Radiation. Every box is printed with a Batch Number and Expiry Date for quick and easy QC tracking.

SteriBio Cellulose Acetate and PES syringe filters have particularly low adsorption which ensures minimal loss of proteins and preservatives.

Membrane Solutions offer Syringes (sterile) which are suitable for sterilization filtration together with MS<sup>®</sup> SteriBio syringe filters.



- All the Syringes are sterilized by Ethylene Oxide,
- Individually packaged.
- No-toxic
- Pyrogen free