

[!\[\]\(a3ea015cc5581cad732d1eb81613fe7b_img.jpg\) PTFE Syringe Filter](#)

[!\[\]\(c8d96c8885d3000a912c2582004aed63_img.jpg\) PES Syringe Filter](#)

[!\[\]\(919a2cb85b99741a73c0c31a427236a8_img.jpg\) MCE Syringe Filter](#)

[!\[\]\(666e09182d4cd268646ea700ea60dcdf_img.jpg\) PVDF Syringe Filter](#)

[!\[\]\(c3d993ca47bfe2a953c700506ce31fa0_img.jpg\) CA Syringe Filter](#)

[!\[\]\(d66ff64371a51729ac8c1cdaa685ba6f_img.jpg\) PP Syringe Filter](#)

[!\[\]\(e3f8612927870f2e0f9f5989e6dd3064_img.jpg\) Glass Fiber Syringe Filter](#)

[!\[\]\(003082e50e3009141f59bd5df831749f_img.jpg\) Reusable Syringe Filter Holders](#)

[!\[\]\(17413706fd4997a1a4bdf85c6864eee1_img.jpg\) Syringe Filter with Prefilter](#)

[!\[\]\(faf942dc3e59ce8eb64b4ac481eca7e0_img.jpg\) 50mm Syringe Filter](#)

[!\[\]\(cf531ed27e91483460120fcc057b3901_img.jpg\) Cartridge Filter](#)

[!\[\]\(d3102649f02e825ddb76dc3de0190154_img.jpg\) Filter Bag](#)

- Compatibility Chemical resistance as for polycarbonate and silicone
- Dead volume Less than 0.3 ml after bubble point
- Filter Requires 13, 25, 50 mm diameter membrane filter
- 25 mm diameter Filtration area 3 cm²
- Flow rates• Typical values for water at 1 bar (100 kPa), 70 ml/min with 0.2µm, 110 ml/min with 0.45 µm pore size filters
- Materials• polypropylene top and bottom part. Silicone gasket (20.5 x 26.5mm, replacement for a pack of 10)
- Pressure limit• Max. operating pressure, 7 bar (700 kPa)
- Sterilisation By autoclaving (121°C)

Application

- Gas particulate and bacteria filtration and then inspect them
- Oil particulate and bacteria filtration and then inspect them
- Alcohol particulate and bacteria filtration and then inspect them
- Other solvent particulate and bacteria filtration and then inspect them

Order Information

Item#	Description	Pcs per box		
SFRPP13	Reusable Syringe Filter, Polypropylene Housing,Diameter:13 (mm)	10		
SFRPP25	Reusable Syringe Filter, Polypropylene Housing,Diameter:25 (mm)	10		
SFRPP35	Reusable Syringe Filter, Polypropylene Housing,Diameter:35 (mm)	10		
SFRPP50	Reusable Syringe Filter, Polypropylene Housing,Diameter:50 (mm)	10		
SFRSS13	Reusable Syringe Filter, Stainless Steel Housing,Diameter:13 (mm)	10		