

SimplePure Syringe Filter Product List 2014

Membrane Solutions

- 2) All the syringe filters are printed with **SimplePure** brand.
 4) 0.1µm, 0.8µm, 3.0µm, 5.0µm is available, but at \$15/pk - Quote required
 5) **Min. order qty is 5 packs per item.**
 6) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

Item#	Description	Packaging
Nylon66 Syringe Filter (Polyamide-66, 13mm/25mm)		
SFNY013022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFNY013045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFNY025022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFNY025045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
PTFE Syringe Filter (Polytetrafluoroethylene, Hydrophobic, 13mm/25mm)		
SFPTFE013022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFPTFE013045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFPTFE025022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFPTFE025045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
PTFE Syringe Filter (Polytetrafluoroethylene, Hydrophilic, 13mm/25mm)		
SFPTFE013022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFPTFE013045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFPTFE025022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFPTFE025045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
PES Syringe Filter (13mm/25mm)		
SFPES013022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFPES013045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFPES025022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFPES025045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
PVDF Syringe Filter (Polyvinylidene Fluoride, 13mm/25mm)		
SFPVDF013022N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFPVDF013045N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFPVDF025022N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFPVDF025045N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
GF Syringe Filter (Glass Fiber, 13mm/25mm)		
SFGF013070N	Nonsterile GF Syringe Filters, Pore:0.7(µm), Housing Diameter:13(mm)	100
SFGF013100N	Nonsterile GF Syringe Filters, Pore:1.0(µm), Housing Diameter:13(mm)	100
SFGF025070N	Nonsterile GF Syringe Filters, Pore:0.7(µm), Housing Diameter:25(mm)	100
SFGF025100N	Nonsterile GF Syringe Filters, Pore:1.0(µm), Housing Diameter:25(mm)	100
CA Syringe Filter (Cellulose Acetate, 13mm/25mm)		
SFCA013022N	Nonsterile CA Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFCA013045N	Nonsterile CA Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFCA025022N	Nonsterile CA Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFCA025045N	Nonsterile CA Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
MCE Syringe Filter (Mixed Cellulose Ester, 13mm/25mm)		
SFMCE013022N	Nonsterile MCE Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFMCE013045N	Nonsterile MCE Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFMCE025022N	Nonsterile MCE Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFMCE025045N	Nonsterile MCE Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100
PP Syringe Filter (Polypropylene, 13mm/25mm)		
SFPP013022N	Nonsterile PP Syringe Filters, Pore:0.22(µm), Housing Diameter:13(mm)	100
SFPP013045N	Nonsterile PP Syringe Filters, Pore:0.45(µm), Housing Diameter:13(mm)	100
SFPP025022N	Nonsterile PP Syringe Filters, Pore:0.22(µm), Housing Diameter:25(mm)	100
SFPP025045N	Nonsterile PP Syringe Filters, Pore:0.45(µm), Housing Diameter:25(mm)	100

Membrane Filters Product List 2014

Item#	Type	Description	Packaging
Nylon Membrane Filter			
MFNY013022	Pore size: 0.22um	Nylon Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFNY025022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFNY047022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFNY090022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFNY142022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFNY293022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:293(mm)	25
MFNY013045	Pore size: 0.45um	Nylon Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFNY025045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFNY047045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFNY090045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFNY142045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFNY293045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:293(mm)	25
MFNY013080	Pore size: 0.8um	Nylon Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFNY025080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFNY047080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFNY090080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	200
MFNY142080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	200
MFNY293080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:293(mm)	200
MFNY013100	Pore size: 1.2um	Nylon Membrane Filter, Pore:1.2(μm), Diameter:13(mm)	400
MFNY025100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFNY047100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFNY090100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFNY142100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFNY293100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:293(mm)	25
MFNY013300	Pore size: 3.0um	Nylon Membrane Filter, Pore:3.0(μm), Diameter:13(mm)	400
MFNY025300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFNY047300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFNY090300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFNY142300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MFNY293300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:293(mm)	25
MFNY013500	Pore size: 5.0um	Nylon Membrane Filter, Pore:5.0(μm), Diameter:13(mm)	400
MFNY025500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:25(mm)	200
MFNY047500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:47(mm)	200
MFNY090500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:90(mm)	100
MFNY142500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:142(mm)	50
MFNY293500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:293(mm)	25
Hydrophobic PTFE Membrane Filter			
MFPTFE013022B	Pore size: 0.22um	Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPTFE025022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPTFE047022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPTFE090022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPTFE142022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPTFE013045B	Pore size: 0.45um	Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPTFE025045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200

MFPTFE047045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPTFE090045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPTFE142045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPTFE025010B	Pore size: 0.1μm	Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:25(mm)	200
MFPTFE047010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:47(mm)	200
MFPTFE090010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:90(mm)	100
MFPTFE142010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:142(mm)	50
MFPTFE025100B	Pore size: 1.2μm	Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFPTFE047100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFPTFE090100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFPTFE142100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFPTFE025300B	Pore size: 3.0μm	Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFPTFE047300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFPTFE090300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFPTFE142300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MFPTFE025300B	Pore size: 5.0μm	Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFPTFE047500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:47(mm)	200
MFPTFE090500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:90(mm)	100
MFPTFE142500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:142(mm)	50
Hydrophilic PTFE Membrane Filter			
MFPTFE013022L	Pore size: 0.22μm	Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPTFE025022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPTFE047022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPTFE090022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPTFE013045L	Pore size: 0.45μm	Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPTFE025045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPTFE047045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPTFE090045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
CA Membrane Filter			
MFCA013022	Pore size: 0.22μm	CA Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFCA025022		CA Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFCA047022		CA Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFCA090022		CA Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFCA142022		CA Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFCA013045	Pore size: 0.45μm	CA Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFCA025045		CA Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFCA047045		CA Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFCA090045		CA Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFCA142045		CA Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFCA013080	Pore size: 0.8μm	CA Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFCA025080		CA Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFCA047080		CA Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFCA090080		CA Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	100
MFCA142080		CA Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	50
CN Membrane Filter			

Membrane Filters Product List 2014

MFCN047022	Pore size: 0.22um	CN Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFCN047045		CN Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFCN047080		CN Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MCE Membrane Filter			
MFMC013022	Pore size: 0.22um	MCE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFMC025022		MCE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFMC037022		MCE Membrane Filter, Pore:0.22(μm), Diameter:37(mm)	200
MFMC047022		MCE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFMC090022		MCE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFMC142022		MCE Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFMC013045	Pore size: 0.45um	MCE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFMC025045		MCE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFMC037045		MCE Membrane Filter, Pore:0.45(μm), Diameter:37(mm)	200
MFMC047045		MCE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFMC090045		MCE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFMC142045		MCE Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFMC013080	Pore size: 0.8um	MCE Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFMC025080		MCE Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFMC037080		MCE Membrane Filter, Pore:0.8(μm), Diameter:37(mm)	200
MFMC047080		MCE Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFMC090080		MCE Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	100
MFMC142080		MCE Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	50
MFMC013100	Pore size: 1.2um	MCE Membrane Filter, Pore:1.2(μm), Diameter:13(mm)	400
MFMC025100		MCE Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFMC037100		MCE Membrane Filter, Pore:1.2(μm), Diameter:37(mm)	200
MFMC047100		MCE Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFMC090100		MCE Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFMC142100		MCE Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFMC013300	Pore size: 3.0um	MCE Membrane Filter, Pore:3.0(μm), Diameter:13(mm)	400
MFMC025300		MCE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFMC037300		MCE Membrane Filter, Pore:3.0(μm), Diameter:37(mm)	200
MFMC047300		MCE Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFMC090300		MCE Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFMC142300		MCE Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MCE Membrane Filter, White, Gridded, Non-sterile			
MFMC047022GWN	Pore size: 0.22um	MCE Gridded Membrane Filter, White, 0.22(μm), 47(mm)	200
MFMC047045GWN	Pore size: 0.45um	MCE Gridded Membrane Filter, White, 0.45(μm), 47(mm)	200
MFMC047080GWN	Pore size: 0.80um	MCE Gridded Membrane Filter, White, 0.80(μm), 47(mm)	200
MCE Membrane Filter, Black, Gridded, Non-sterile			
MFMC047022GBN	Pore size: 0.22um	MCE Gridded Membrane Filter, Black, 0.22(μm), 47(mm)	200
MFMC047045GBN	Pore size: 0.45um	MCE Gridded Membrane Filter, Black, 0.45(μm), 47(mm)	200
MFMC047080GBN	Pore size: 0.80um	MCE Gridded Membrane Filter, Black, 0.80(μm), 47(mm)	200
PES Membrane Filter			
MFPE013022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPE025022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPE047022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200

Membrane Filters Product List 2014

MFPE090022	Pore size: 0.22m	PES Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPE0142022	Pore size: 0.22m	PES Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPE013045	Pore size: 0.45m	PES Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPE025045		PES Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPE047045		PES Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPE090045		PES Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPE0142045		PES Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPE013010	Pore size: 0.1m	PES Membrane Filter, Pore:0.1(μm), Diameter:13(mm)	400
MFPE025010		PES Membrane Filter, Pore:0.1(μm), Diameter:25(mm)	200
MFPE047010		PES Membrane Filter, Pore: 0.1(μm), Diameter: 47(mm)	200
MFPE090010		PES Membrane Filter, Pore:0.1(μm), Diameter:90(mm)	100
PVDF Membrane Filter			
MFPVDF013022	Pore size: 0.22m	PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 13(mm)	400
MFPVDF025022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 25(mm)	200
MFPVDF047022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 47(mm)	200
MFPVDF090022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 90(mm)	100
MFPVDF142022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 142(mm)	50
MFPVDF013045	Pore size: 0.45m	PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 13(mm)	400
MFPVDF025045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 25(mm)	200
MFPVDF047045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 47(mm)	200
MFPVDF090045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 90(mm)	100
MFPVDF142045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 142(mm)	50
PP Membrane Filter			
MFPP025022	Pore size: 0.22m	PP Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPP047022		PP Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPP090022		PP Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPP142022		PP Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPP250022		PP Membrane Filter, Pore:0.22(μm), Diameter:250(mm)	50
MFPP293022		PP Membrane Filter, Pore:0.22(μm), Diameter:293(mm)	25
MFPP025045	Pore size: 0.45m	PP Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPP047045		PP Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPP090045		PP Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPP142045		PP Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPP250045		PP Membrane Filter, Pore:0.45(μm), Diameter:250(mm)	50
MFPP293045		PP Membrane Filter, Pore:0.45(μm), Diameter:293(mm)	25

MS Plastic Vacuum Filter System Product List 2013

Note: 1) All of the price listed below are **Ex.work price**.

4) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty(pcs/cs)	Case Size(cm)
Filter System (Membrane Diameter: 50mm, Funnel Capacity : 250ml,				
Receiver Bottle Capacity : 250ml)				
VFPPE122250	PES	0.22	12	33*33*30
VFPPE145250		0.45	12	
VFPMCE122250	MCE	0.22	12	
VFPMCE145250		0.45	12	
VFPCA122250	CA	0.22	12	
VFPCA145250		0.45	12	
VFPNY122250	Nylon	0.22	12	
VFPNY145250		0.45	12	
VFPVDF122250	PVDF	0.22	12	
VFPVDF145250		0.45	12	
Filter System (Membrane Diameter : 50mm, Funnel Capacity : 250ml,				
Receiver Bottle Capacity : 500ml)				
VFPPE122500	PES	0.22	12	33*33*30
VFPPE145500		0.45	12	
VFPMCE122500	MCE	0.22	12	
VFPMCE145500		0.45	12	
VFPCA122500	CA	0.22	12	
VFPCA145500		0.45	12	
VFPNY122500	Nylon	0.22	12	
VFPNY145500		0.45	12	
VFPVDF122500	PVDF	0.22	12	
VFPVDF145500		0.45	12	
Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty (pcs/cs)	Case Size(cm)
Filter System (Membrane Diameter : 90mm, Funnel Capacity : 500ml,				
Receiver Bottle Capacity : 500ml)				
VFPPE222500	PES	0.22	12	52.5*38.5*29.5
VFPPE245500		0.45	12	
VFPMCE222500	MCE	0.22	12	
VFPMCE245500		0.45	12	
VFPCA222500	CA	0.22	12	
VFPCA245500		0.45	12	
VFPNY222500	Nylon	0.22	12	
VFPNY245500		0.45	12	
VFPVDF222500	PVDF	0.22	12	
VFPVDF245500		0.45	12	

MS Plastic Vacuum Filter System Product List 2013

Filter System (Membrane Diameter : 90mm, Funnel Capacity : 500ml, Receiver Bottle Capacity:1000ml)				
VFPES2221000	PES	0.22	12	52.5*38.5*29.5
VFPES2451000		0.45	12	
VFPMCE2221000	MCE	0.22	12	
VFPMCE2451000		0.45	12	
VFPCA2221000	CA	0.22	12	
VFPCA2451000		0.45	12	
VFPNY2221000	Nylon	0.22	12	
VFPNY2451000		0.45	12	
VFPVDF2221000	PVDF	0.22	12	
VFPVDF2451000		0.45	12	
Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty (pcs/cs)	Case Size(cm)
Funnel Capacity : 250ml, Membrane Diameter : 50mm				
VFPES122250F	PES	0.22	24	33*33*30
VFPES145250F		0.45	24	
VFPMCE122250F	MCE	0.22	24	
VFPMCE145250F		0.45	24	
VFPCA122250F	CA	0.22	24	
VFPCA145250F		0.45	24	
VFPNY122250F	Nylon	0.22	24	
VFPNY145250F		0.45	24	
VFPVDF122250F	PVDF	0.22	24	
VFPVDF145250F		0.45	24	
Funnel Capacity : 500ml, Membrane Diameter : 90mm				
VFPES222500F	PES	0.22	24	52.5*38.5*29.5
VFPES245500F		0.45	24	
VFPMCE222500F	MCE	0.22	24	
VFPMCE245500F		0.45	24	
VFPCA222500F	CA	0.22	24	
VFPCA245500F		0.45	24	
VFPNY222500F	Nylon	0.22	24	
VFPNY245500F		0.45	24	
VFPVDF222500F	PVDF	0.22	24	
VFPVDF245500F		0.45	24	
Cat.No.	Description	Qty(pcs/cs)	Case Size(cm)	
Receiver Bottle				
VFP250B	250ml polystyrene receiver bottle, sterile	24	33*33*30	
VFP500B	500ml polystyrene receiver bottle, sterile	24	52.5*38.5*29.5	
VFP1000B	1000ml polystyrene receiver bottle, sterile	24	52.5*38.5*29.5	

MS Glass Filtration System Product List 2013

Note: 1) All of the price listed below are **Ex work price** below are Ex work price.

2) All the price listed below are for Membrane Solutions' Brand label.

3) **Order Qty** means **total items** in one order.

4) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

	Item#	Description	Packaging
Glass Vacuum Filter	VFG030010	1) 300ml Glass Funnel with cover	1set
		2) 1000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Glass Vacuum Filter	VFG030020	1) 300ml Glass Funnel with cover	1set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Glass Vacuum Filter	VFG050020	1) 500ml Glass Funnel with cover	1set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Support Screen	VFG147SS	Pyrex glass support screen (1000ML)	1pcs
	VFG047SS	Pyrex glass support screen (2000ML)	1pcs
	VFG010RB	Solvent Collection Bottle, 1000ml	1pcs
	VFG020RB	Solvent Collection Bottle, 2000ml	1pcs

Vacuum Pump list 2013

	Item#	Description	Packaging
Vacuun Pump	VPJ0201	Vacuum Pump, 300mbar, 12L/min GM-0.20((Oil free))	1set
Vacuun Pump	VPJ0332	Vacuum Pump, 250mbar, 20L/min GM-0.33A(Oil free),	1set
Vacuun Pump	VPJ0333	Vacuum Pump, 50mbar, 20L/min	1set
Vacuun Pump	VPJ0501	Vacuum Pump, 250mbar, 30L/min	1set
Vacuun Pump	VPJ0502	Vacuum Pump ,50mbar 30L/min	1set
Vacuun Pump	VPJ1001	Vacuum Pump, 250mbar, 60L/min GM-1.0A(Oil free),	1set