



# VICI Valco-Cheminert PRODUCT LIST 2012

Part No.	Description 1	Specs
.2FR.5-5	Frit 1/32in 0.2 micron	.25mm thick 316SS 5/pkg
.5FR.5-5	Frit 1/32in 0.5 micron	.25mm thick 316SS 5/pkg
.5FR1-10	Frit 1/16in 0.5 micron	.75mm thick 316SS 10/pkg
.5FR1HC-10	Frit 1/16in 0.5 micron	.75mm thick HC 10/pkg
.5FR2-10	Frit 1/8in 0.5 micron	1mm thick 316SS 10/pkg
.5FR4-10	Frit 1/4in 0.5 micron	1mm thick 316SS 10/pkg
.5SR.5-10	Screen 1/32in OD 0.5 micron 316SS 10/pkg	
.5SR1-10	Screen 1/16in OD 0.5 micron 316SS 10/pkg	
.5SR2-10	Screen 1/8in OD 0.5 micron 316SS 10/pkg	
.5SR4-10	Screen 1/4in 0.5 micron 316SS 10/pkg	
10FR1-10	Frit 1/16in 10 micron	.75mm thick 316SS 10/pkg
10FR2-10	Frit 1/8in 10 micron	1mm thick 316SS 10/pkg
10FR4-10	Frit 1/4in 10 micron	1mm thick 316SS 10/pkg
10FR4HC-10	Frit 1/4in 10 micron	1mm thick HC 10/pkg
10SR.5-10	Screen 1/32in 10 micron 316SS 10/pkg	
10SR1-10	Screen 1/16in 10 micron 316SS 10/pkg	
10SR2-10	Screen 1/8in 10 micron 316SS 10/pkg	
10SR4-10	Screen 1/4in 10 micron 316SS 10/pkg	
1SR.5-10	Screen 1/32in OD 1 micron 316SS 10/pkg	
1SR1-10	Screen 1/16in OD 1 micron 316SS 10/pkg	
1SR2-10	Screen 1/8in OD 1 micron 316SS 10/pkg	
1SR4-10	Screen 1/4in 1 micron 316SS 10/pkg	
210UW	10 port 2-pos valve manual/2in so	1/8 x .75mm 75C/5000psi liq N60/H
24UW	4 port 2-pos valve manual/2in so	1/8 x .75mm 75C/5000psi liq N60/H
26UW	6 port 2-pos valve manual/2in so	1/8 x .75mm 75C/5000psi liq N60/H
28UW	8 port 2-pos valve manual/2in so	1/8 x .75mm 75C/5000psi liq N60/H
2C10UW	10 port 2-pos valve manual/2in so	1/16 x .75mm 75C/5000psi liq N60/H
2C4UW	4 port 2-pos valve manual/2in so	1/16 x .75mm 75C/5000psi liq N60/H
2C6UW	6 port 2-pos valve manual/2in so	1/16 x .75mm 75C/5000psi liq N60/H
2C8UW	8 port 2-pos valve manual/2in so	1/16 x .75mm 75C/5000psi liq N60/H
2Ci4UW.2	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 75C/5000psi liq .2ul N60/H
2Ci4UW.5	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 75C/5000psi liq .5ul N60/H
2Ci4UW1	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 75C/5000psi liq 1ul N60/H
2Ci4UW2	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 75C/5000psi liq 2ul N60/H
2Ci4UWE.2	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 175C/1000psi liq .2ul N60/E
2Ci4UWE.5	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 175C/1000psi liq .5ul N60/E
2Ci4UWE1	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 175C/1000psi liq 1ul N60/E
2Ci4UWE2	4 port 2-pos int vol valve manual/2in so	1/16 x .75mm 175C/1000psi liq 2ul N60/E
2Ci4WE.06	4 port 2-pos int vol valve manual/2in so	1/16 x .4mm 175C/1000psi liq .06ul N60/E
2Ci4WE.1	4 port 2-pos int vol valve manual/2in so	1/16 x .4mm 175C/1000psi liq .1ul N60/E
2Ci4WE.2	4 port 2-pos int vol valve manual/2in so	1/16 x .4mm 175C/1000psi liq .2ul N60/E
2Ci4WE.5	4 port 2-pos int vol valve manual/2in so	1/16 x .4mm 175C/1000psi liq .5ul N60/E
2CSC10MWE	10 mpos common outlet path MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSC12MWE	12 mpos common outlet path MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSC16MWE	16 mpos common outlet path MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSC6MWE	6 mpos common outlet path MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSD10MWE	10 mpos dead-end path MW manual/2in so	1/16 x .75mm 200C/400psi gas N60/E
2CSD12MWE	12 mpos dead-end path MW manual/2in so	1/16 x .75mm 200C/400psi gas N60/E
2CSD16MWE	16 mpos dead-end path MW manual/2in so	1/16 x .75mm 200C/400psi gas N60/E
2CSD6MWE	6 mpos dead-end path MW manual/2in so	1/16 x .75mm 200C/400psi gas N60/E
2CSF10MWE	10 mpos flow-through MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSF12MWE	12 mpos flow-through MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CSF16MWE	16 mpos flow-through MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E




# VICI Valco-Cheminert PRODUCT LIST 2012




2CSF6MWE	6 mpos flow-through MW manual/2in so	1/16 x 1mm 200C/200psi gas N60/E
2CST10MWE	10 mpos trapping path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CST12MWE	12 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CST16MWE	16 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CST6MWE	6 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CSTF10MWE	10 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CSTF12MWE	12 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CSTF16MWE	16 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2CSTF6MWE	6 mpos trap/flow path MW manual/2in so	1/16 x .75mm 200C/200psi gas N60/E
2FR1-10	Frit 1/16in 2 micron	.75mm thick 316SS 10/pkg
2FR1HC-10	Frit 1/16in 2 micron	.75mm thick HC 10/pkg
2FR1K	Frit 1in 2 micron	1.5mm thick 316SS
2FR1KHC	Frit 1in 2 micron	1.5mm thick HC
2FR1KTI	Frit 1in 2 micron	1.5mm thick Titanium
2FR1TI-10	Frit 1/16in 2 micron	.75mm thick Titanium 10/pkg
2FR2-10	Frit 1/8in 2 micron	1mm thick 316SS 10/pkg
2FR2HC-10	Frit 1/8in 2 micron	1mm thick HC 10/pkg
2FR2TI-10	Frit 1/8in 2 micron	1mm thick Titanium 10/pkg
2FR4-10	Frit 1/4in 2 micron	1mm thick 316SS 10/pkg
2FR4HC-10	Frit 1/4in 2 micron	1mm thick HC 10/pkg
2FR4TI-10	Frit 1/4in 2 micron	1mm thick Titanium 10/pkg
2FR6	Frit 3/8in 2 micron	1mm thick 316SS
2FR6HC	Frit 3/8in 2 micron	1mm thick HC
2FR6TI	Frit 3/8in 2 micron	1mm thick Titanium
2FR8	Frit 1/2in 2 micron	1mm thick 316SS
2FR8HC	Frit 1/2in 2 micron	1mm thick HC
2FR8TI	Frit 1/2in 2 micron	1mm thick Titanium
2I4UWE.2	4 port 2-pos int vol valve manual/2in so	1/8 x .75mm 175C/1000psi liq .2ul N60/E
2I4UWE.5	4 port 2-pos int vol valve manual/2in so	1/8 x .75mm 175C/1000psi liq .5ul N60/E
2I4UWE1	4 port 2-pos int vol valve manual/2in so	1/8 x .75mm 175C/1000psi liq 1ul N60/E
2I4UWE2	4 port 2-pos int vol valve manual/2in so	1/8 x .75mm 175C/1000psi liq 2ul N60/E
2L10UW	10 port 2-pos valve manual/2in so	1/8 x 1.0mm 75C/5000psi liq N60/H
2L4UW	4 port 2-pos valve manual/2in so	1/8 x 1.7mm 75C/5000psi liq N60/H
2L6UW	6 port 2-pos valve manual/2in so	1/8 x 1.7mm 75C/5000psi liq N60/H
2L8UW	8 port 2-pos valve manual/2in so	1/8 x 1.3mm 75C/5000psi liq N60/H
2NI4WE.06	4 port 2-pos int vol valve manual/2in so	1/32 x .25mm 175C/1000psi liq .06ulN60/E
2NI4WE.1	4 port 2-pos int vol valve manual/2in so	1/32 x .25mm 175C/1000psi liq .1ulN60/E
2NI4WE.2	4 port 2-pos int vol valve manual/2in so	1/32 x .25mm 175C/1000psi liq .2ulN60/E
2NI4WE.5	4 port 2-pos int vol valve manual/2in so	1/32 x .25mm 175C/1000psi liq .5ulN60/E
2RAD	Right angle drive 2in standoff	
2SC10MWE	10 mpos common outlet pathMWmanual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SC12MWE	12 mpos common outlet pathMWmanual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SC16MWE	16 mpos common outlet pathMWmanual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SC6MWE	6 mpos common outlet pathMWmanual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SD10MWE	10 mpos dead-end path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SD12MWE	12 mpos dead-end path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SD16MWE	16 mpos dead-end path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SD6MWE	6 mpos dead-end path MW manual/2in so	1/8 x 1mm 200C/400psi gas N60/E
2SF10MWE	10 mpos flow-through MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SF12MWE	12 mpos flow-through MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SF16MWE	16 mpos flow-through MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SF6MWE	6 mpos flow-through MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2SOA	Standoff assembly 2in	for 2 pos actuator
2SOAMMP	Standoff assembly 2in keyed	for m-electric actuator
2SOAMP	Standoff assembly 2in	for mpos actuator
2SR.5-10	Screen 1/32in 2 micron 316SS 10/pkg	
2SR1-10	Screen 1/16in 2 micron 316SS 10/pkg	
2SR2-10	Screen 1/8in 2 micron 316SS 10/pkg	
2SR4-10	Screen 1/4in 2 micron 316SS 10/pkg	
2ST10MWE	10 mpos trapping path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2ST12MWE	12 mpos trapping path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




2ST16MWE	16 mpos trapping path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2ST6MWE	6 mpos trapping path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2STF10MWE	10 mpos trap/flow path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2STF12MWE	12 mpos trap/flow path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2STF16MWE	16 mpos trap/flow path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
2STF6MWE	6 mpos trap/flow path MW manual/2in so	1/8 x 1mm 200C/200psi gas N60/E
410UWE	10 port 2-pos valve manual/4in so	1/8 x .75mm 225C/400psi gas N60/E
410UWT	10 port 2-pos valve manual/4in so	1/8 x .75mm 330C/300psi gas N60/T
41E1-120VAC	4-way solenoid air valve 120Vac	Non-CE !!!
41E1-12VDC	4-way solenoid air valve 12Vdc	Non-CE !!!
41E1-220VAC	4-way solenoid air valve 220Vac	Non-CE !!!
41E1-24VAC	4-way solenoid air valve 24Vac	Non-CE !!!
41E1-24VDC	4-way solenoid air valve 24Vdc	Non-CE !!!
44UWE	4 port 2-pos valve manual/4in so	1/8 x .75mm 225C/400psi gas N60/E
44UWT	4 port 2-pos valve manual/4in so	1/8 x .75mm 330C/300psi gas N60/T
46UWE	6 port 2-pos valve manual/4in so	1/8 x .75mm 225C/400psi gas N60/E
46UWT	6 port 2-pos valve manual/4in so	1/8 x .75mm 330C/300psi gas N60/T
48UWE	8 port 2-pos valve manual/4in so	1/8 x .75mm 225C/400psi gas N60/E
48UWT	8 port 2-pos valve manual/4in so	1/8 x .75mm 330C/300psi gas N60/T
4C10UWE	10 port 2-pos valve manual/4in so	1/16 x .75mm 225C/400psi gas N60/E
4C10UWT	10 port 2-pos valve manual/4in so	1/16 x .75mm 330C/300psi gas N60/T
4C10WE	10 port 2-pos valve manual/4in so	1/16 x .4mm 225C/400psi gas N60/E
4C10WT	10 port 2-pos valve manual/4in so	1/16 x .4mm 350C/300psi gas N60/T
4C4UWE	4 port 2-pos valve manual/4in so	1/16 x .75mm 225C/400psi gas N60/E
4C4UWT	4 port 2-pos valve manual/4in so	1/16 x .75mm 330C/300psi gas N60/T
4C4WE	4 port 2-pos valve manual/4in so	1/16 x .4mm 225C/400psi gas N60/E
4C4WT	4 port 2-pos valve manual/4in so	1/16 x .4mm 350C/300psi gas N60/T
4C6UWE	6 port 2-pos valve manual/4in so	1/16 x .75mm 225C/400psi gas N60/E
4C6UWT	6 port 2-pos valve manual/4in so	1/16 x .75mm 330C/300psi gas N60/T
4C6WE	6 port 2-pos valve manual/4in so	1/16 x .4mm 225C/400psi gas N60/E
4C6WT	6 port 2-pos valve manual/4in so	1/16 x .4mm 350C/300psi gas N60/T
4C8UWE	8 port 2-pos valve manual/4in so	1/16 x .75mm 225C/400psi gas N60/E
4C8UWT	8 port 2-pos valve manual/4in so	1/16 x .75mm 330C/300psi gas N60/T
4C8WE	8 port 2-pos valve manual/4in so	1/16 x .4mm 225C/400psi gas N60/E
4C8WT	8 port 2-pos valve manual/4in so	1/16 x .4mm 350C/300psi gas N60/T
4N10WE	10 port 2-pos valve manual/4in so	1/32 x .25mm 225C/400psi gas N60/E
4N10WT	10 port 2-pos valve manual/4in so	1/32 x .25mm 350C/300psi gas N60/T
4N4WE	4 port 2-pos valve manual/4in so	1/32 x .25mm 225C/400psi gas N60/E
4N4WT	4 port 2-pos valve manual/4in so	1/32 x .25mm 350C/300psi gas N60/T
4N6WE	6 port 2-pos valve manual/4in so	1/32 x .25mm 225C/400psi gas N60/E
4N6WT	6 port 2-pos valve manual/4in so	1/32 x .25mm 350C/300psi gas N60/T
4N8WE	8 port 2-pos valve manual/4in so	1/32 x .25mm 225C/400psi gas N60/E
4N8WT	8 port 2-pos valve manual/4in so	1/32 x .25mm 350C/300psi gas N60/T
4SOA	Standoff assembly 4in	for 2 pos actuator
4SOAMP	Standoff assembly 4in	for mpos actuator
4SOUTH	Handle with 4in standoff Valco UW-type	2-pos valve with 1 or 2 mounting holes
4SOUTHMP	Handle with 4in standoff Valco UW-type	2-pos valve with no mounting holes
4SOWK	4in standoff assy for 2 pos vlv W/WK	with 1 or 2 mounting holes
4SOWKMP	4in Standoff assy for 2 pos vlv W/WK	with no mounting holesand Cheminert vlv
4UW	!!! RECOMMEND 24UW	SUPPLIED WITH LARGER HANDLE
4UWE	4 port 2-pos valve manual	1/8 x .75mm 225C/400psi gas N60/E
4VL4MWE2	4 port 2-pos valve MW manual/4in so	1/4 x 4mm 75C/100psi gas N60/E2
4VL6MWE2	6 port 2-pos valve MW manual/4in so	1/4 x 4mm 75C/100psi gas N60/E2
4VL8MWE2	8 port 2-pos valve MW manual/4in so	1/4 x 4mm 75C/100psi gas N60/E2
6UW	!!! RECOMMEND 26UW	SUPPLIED WITH LARGER HANDLE
6UWE	6 port 2-pos valve manual	1/8 x .75mm 225C/400psi gas N60/E
8UW	!!! RECOMMEND 28UW	SUPPLIED WITH LARGER HANDLE
8UWE	8 port 2-pos valve manual	1/8 x .75mm 225C/400psi gas N60/E
A10	Air actuator 10 mpos close mount	36 deg rotation
A102	Air actuator 10 mpos 2in standoff	36 deg rotation
A10S	Air actuator 10 mpos for standoff	36 deg rotation

# VICI Valco-Cheminert PRODUCT LIST 2012




A10UW	10 port 2-pos valve air actuated	1/8 x .75mm 75C/5000psi liq N60/H
A12	Air actuator 12 mpos close mount	30 deg rotation
A122	Air actuator 12 mpos 2in standoff	30 deg rotation
A12S	Air actuator 12 mpos for standoff	30 deg rotation
A16	Air actuator 16 mpos close mount	22.5 deg rotation
A162	Air actuator 16 mpos 2in standoff	22.5 deg rotation
A16S	Air actuator 16 mpos for standoff	22.5 deg rotation
A2CI4UWE.2	4 port 2-pos int vol valve air act/2in so	1/16 x .75mm 175C/1000psi liq .2ul N60/E
A2CI4UWE.5	4 port 2-pos int vol valve air act/2in so	1/16 x .75mm 175C/1000psi liq .5ul N60/E
A2CI4UWE1	4 port 2-pos int vol valve air act/2in so	1/16 x .75mm 175C/1000psi liq 1ul N60/E
A2CI4UWE2	4 port 2-pos int vol valve air act/2in so	1/16 x .75mm 175C/1000psi liq 2ul N60/E
A2CI4WE.06	4 port 2-pos int vol valve air act/2in so	1/16 x .4mm 175C/1000psi liq .06ul N60/E
A2CI4WE.1	4 port 2-pos int vol valve air act/2in so	1/16 x .4mm 175C/1000psi liq .1ul N60/E
A2CI4WE.2	4 port 2-pos int vol valve air act/2in so	1/16 x .4mm 175C/1000psi liq .2ul N60/E
A2CI4WE.5	4 port 2-pos int vol valve air act/2in so	1/16 x .4mm 175C/1000psi liq .5ul N60/E
A2CSC10MWE	10 mpos common outlet path MWair act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSC12MWE	12 mpos common outlet path MWair act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSC16MWE	16 mpos common outlet path MWair act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSC6MWE	OBSOLETE USE AH2CSC6MWE	
A2CSD10MWE	10 mpos dead-end path MW air act/2in so	1/16 x .75mm 200C/400psi gas N60/E
A2CSD12MWE	12 mpos dead-end path MW air act/2in so	1/16 x .75mm 200C/400psi gas N60/E
A2CSD16MWE	16 mpos dead-end path MW air act/2in so	1/16 x .75mm 200C/400psi gas N60/E
A2CSD6MWE	6 mpos dead-end path MW air act/2in so	1/16 x .75mm 200C/400psi gas N60/E
A2CSF10MWE	10 mpos flow-through MW air act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSF12MWE	12 mpos flow-through MW air act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSF16MWE	16 mpos flow-through MW air act/2in so	1/16 x 1mm 200C/200psi gas N60/E
A2CSF6MWE	OBSOLETE USE AH2CSF6MWE	
A2CST10MWE	10 mpos trapping path MW air act high torq./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CST12MWE	12 mpos trapping path MW air act high torq./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CST16MWE	16 mpos trapping path MW air act high torq./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CST6MWE	OBSOLETE USE AH2CST6MWE	
A2CSTF10MWE	10 mpos trap/flow path MW air act high t./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CSTF12MWE	12 mpos trap/flow path MW air act high t./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CSTF16MWE	16 mpos trap/flow path MW air act high t./2in so	1/16 x .75mm 200C/200psi gas N60/E
A2CSTF6MWE	OBSOLETE USE AH2CSTF6MWE	
A2I4UWE.2	4 port 2-pos int vol valve air act/2in so	1/8 x .75mm 175C/1000psi liq .2ul N60/E
A2I4UWE.5	4 port 2-pos int vol valve air act/2in so	1/8 x .75mm 175C/1000psi liq .5ul N60/E
A2I4UWE1	4 port 2-pos int vol valve air act/2in so	1/8 x .75mm 175C/1000psi liq 1ul N60/E
A2I4UWE2	4 port 2-pos int vol valve air act/2in so	1/8 x .75mm 175C/1000psi liq 2ul N60/E
A2NI4WE.06	4 port 2-pos int vol valve air act/2in so	1/32 x .25mm 175C/1000psi liq .06ulN60/E
A2NI4WE.1	4 port 2-pos int vol valve air act/2in so	1/32 x .25mm 175C/1000psi liq .1ulN60/E
A2NI4WE.2	4 port 2-pos int vol valve air act/2in so	1/32 x .25mm 175C/1000psi liq .2ulN60/E
A2NI4WE.5	4 port 2-pos int vol valve air act/2in so	1/32 x .25mm 175C/1000psi liq .5ulN60/E
A2SC10MWE	10 mpos common outlet path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SC12MWE	12 mpos common outlet path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SC16MWE	16 mpos common outlet path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SC6MWE	OBSOLETE USE AH2SC6MWE	
A2SD10MWE	10 mpos dead-end path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SD12MWE	12 mpos dead-end path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SD16MWE	16 mpos dead-end path MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SD6MWE	6 mpos dead-end path MW air act/2in so	1/8 x 1mm 200C/400psi gas N60/E
A2SF10MWE	10 mpos flow-through MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SF12MWE	12 mpos flow-through MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SF16MWE	16 mpos flow-through MW air act/2in so	1/8 x 1mm 200C/200psi gas N60/E
A2SF6MWE	OBSOLETE USE AH2SF6MWE	
A2ST10MWE	10 mpos trapping path MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
A2ST12MWE	12 mpos trapping path MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
A2ST16MWE	16 mpos trapping path MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
A2ST6MWE	OBSOLETE USE AH2ST6MWE	
A2STF10MWE	10 mpos trap/flow path MW air act high t./2in so	1/8 x 1mm 200C/200psi gas N60/E
A2STF12MWE	12 mpos trap/flow path MW air act high t./2in so	1/8 x 1mm 200C/200psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




A2STF16MWE	16 mpos trap/flow path MW air act high t./2in so	1/8 x 1mm 200C/200psi gas N60/E
A2STF6MWE	OBSOLETE USE AH2STF6MWE	
A2VLSC8MWE2	OBSOLETE USE AH2VLSC8MWE2	
A30	Air actuator 2-pos close mount	30 deg rotation
A304	Air actuator 2-pos 4in standoff	30 deg rotation
A30S	Air actuator 2-pos for standoff	30 deg rotation
A36	Air actuator 2-pos close mount	36 deg rotation
A364	Air actuator 2-pos 4in standoff	36 deg rotation
A36S	Air actuator 2-pos for standoff	36 deg rotation
A4	Air actuator 4 mpos close mount	90 deg rotation
A410UWE	10 port 2-pos valve air actuated/4in so	1/8 x .75mm 225C/400psi gas N60/E
A410UWT	10 port 2-pos valve air actuated/4in so	1/8 x .75mm 330C/300psi gas N60/T
A42	Air actuator 4 mpos 2in standoff	90 deg rotation
A44UWE	4 port 2-pos valve air actuated/4in so	1/8 x .75mm 225C/400psi gas N60/E
A44UWT	4 port 2-pos valve air actuated/4in so	1/8 x .75mm 330C/300psi gas N60/T
A45	Air actuator 2-pos close mount	45 deg rotation
A454	Air actuator 2-pos 4in standoff	45 deg rotation
A45S	Air actuator 2-pos for standoff	45 deg rotation
A46UWE	6 port 2-pos valve air actuated/4in so	1/8 x .75mm 225C/400psi gas N60/E
A46UWT	6 port 2-pos valve air actuated/4in so	1/8 x .75mm 330C/300psi gas N60/T
A48UWE	8 port 2-pos valve air actuated/4in so	1/8 x .75mm 225C/400psi gas N60/E
A48UWT	8 port 2-pos valve air actuated/4in so	1/8 x .75mm 330C/300psi gas N60/T
A4C10UWE	10 port 2-pos valve air actuated/4in so	1/16 x .75mm 225C/400psi gas N60/E
A4C10UWT	10 port 2-pos valve air actuated/4in so	1/16 x .75mm 330C/300psi gas N60/T
A4C10WE	10 port 2-pos valve air actuated/4in so	1/16 x .4mm 225C/400psi gas N60/E
A4C10WT	10 port 2-pos valve air actuated/4in so	1/16 x .4mm 350C/300psi gas N60/T
A4C4UWE	4 port 2-pos valve air actuated/4in so	1/16 x .75mm 225C/400psi gas N60/E
A4C4UWT	4 port 2-pos valve air actuated/4in so	1/16 x .75mm 330C/300psi gas N60/T
A4C4WE	4 port 2-pos valve air actuated/4in so	1/16 x .4mm 225C/400psi gas N60/E
A4C4WT	4 port 2-pos valve air actuated/4in so	1/16 x .4mm 350C/300psi gas N60/T
A4C6UWE	6 port 2-pos valve air actuated/4in so	1/16 x .75mm 225C/400psi gas N60/E
A4C6UWT	6 port 2-pos valve air actuated/4in so	1/16 x .75mm 330C/300psi gas N60/T
A4C6WE	6 port 2-pos valve air actuated/4in so	1/16 x .4mm 225C/400psi gas N60/E
A4C6WT	6 port 2-pos valve air actuated/4in so	1/16 x .4mm 350C/300psi gas N60/T
A4C8UWE	8 port 2-pos valve air actuated/4in so	1/16 x .75mm 225C/400psi gas N60/E
A4C8UWT	8 port 2-pos valve air actuated/4in so	1/16 x .75mm 330C/300psi gas N60/T
A4C8WE	8 port 2-pos valve air actuated/4in so	1/16 x .4mm 225C/400psi gas N60/E
A4C8WT	8 port 2-pos valve air actuated/4in so	1/16 x .4mm 350C/300psi gas N60/T
A4N10WE	10 port 2-pos valve air actuated/4in so	1/32 x .25mm 225C/400psi gas N60/E
A4N10WT	10 port 2-pos valve air actuated/4in so	1/32 x .25mm 350C/300psi gas N60/T
A4N4WE	4 port 2-pos valve air actuated/4in so	1/32 x .25mm 225C/400psi gas N60/E
A4N4WT	4 port 2-pos valve air actuated/4in so	1/32 x .25mm 350C/300psi gas N60/T
A4N6WE	6 port 2-pos valve air actuated/4in so	1/32 x .25mm 225C/400psi gas N60/E
A4N6WT	6 port 2-pos valve air actuated/4in so	1/32 x .25mm 350C/300psi gas N60/T
A4N8WE	8 port 2-pos valve air actuated/4in so	1/32 x .25mm 225C/400psi gas N60/E
A4N8WT	8 port 2-pos valve air actuated/4in so	1/32 x .25mm 350C/300psi gas N60/T
A4S	Air actuator 4 mpos for standoff	90 deg rotation
A4UW	4 port 2-pos valve air actuated	1/8 x .75mm 75C/5000psi liq N60/H
A4VL4MWE2	4 port 2-pos valve MW air act/4in so	1/4 x 4mm 75C/100psi gas N60/E2
A4VL6MWE2	6 port 2-pos valve MW air act/4in so	1/4 x 4mm 75C/100psi gas N60/E2
A4VL8MWE2	8 port 2-pos valve MW air act/4in so	1/4 x 4mm 75C/100psi gas N60/E2
A6	Air actuator 6 mpos close mount	60 deg rotation
A60	Air actuator 2-pos close mount	60 deg rotation
A604	Air actuator 2-pos 4in standoff	60 deg rotation
A60S	Air actuator 2-pos for standoff	60 deg rotation
A62	Air actuator 6 mpos 2in standoff	60 deg rotation
A6S	Air actuator 6 mpos for standoff	60 deg rotation
A6UW	6 port 2-pos valve air actuated	1/8 x .75mm 75C/5000psi liq N60/H
A8	Air actuator 8 mpos close mount	45 deg rotation
A82	Air actuator 8 mpos 2in standoff	45 deg rotation
A8S	Air actuator 8 mpos for standoff	45 deg rotation

# VICI Valco-Cheminert PRODUCT LIST 2012




A8UW	8 port 2-pos valve air actuated	1/8 x .75mm 75C/5000psi liq N60/H
A90	Air actuator 2-pos close mount	90 deg rotation
A904	Air actuator 2-pos 4in standoff	90 deg rotation
A90S	Air actuator 2-pos for standoff	90 deg rotation
AC10UW	10 port 2-pos valve air actuated	1/16 x .75mm 75C/5000psi liq N60/H
AC10W	10 port 2-pos valve air actuated	1/16 x .4mm 75C/5000psi liq N60/H
AC4UW	4 port 2-pos valve air actuated	1/16 x .75mm 75C/5000psi liq N60/H
AC4W	4 port 2-pos valve air actuated	1/16 x .4mm 75C/5000psi liq N60/H
AC6UW	6 port 2-pos valve air actuated	1/16 x .75mm 75C/5000psi liq N60/H
AC6W	6 port 2-pos valve air actuated	1/16 x .4mm 75C/5000psi liq N60/H
AC8UW	8 port 2-pos valve air actuated	1/16 x .75mm 75C/5000psi liq N60/H
AC8W	8 port 2-pos valve air actuated	1/16 x .4mm 75C/5000psi liq N60/H
ACI4UW.2	4 port 2-pos int vol valve air actuated	1/16 x .75mm 75C/5000psi liq .2ul N60/H
ACI4UW.5	4 port 2-pos int vol valve air actuated	1/16 x .75mm 75C/5000psi liq .5ul N60/H
ACI4UW1	4 port 2-pos int vol valve air actuated	1/16 x .75mm 75C/5000psi liq 1ul N60/H
ACI4UW2	4 port 2-pos int vol valve air actuated	1/16 x .75mm 75C/5000psi liq 2ul N60/H
ACI4W.06	4 port 2-pos int vol valve air actuated	1/16 x .4mm 75C/5000psi liq .06ul N60/H
ACI4W.1	4 port 2-pos int vol valve air actuated	1/16 x .4mm 75C/5000psi liq .1ul N60/H
ACI4W.2	4 port 2-pos int vol valve air actuated	1/16 x .4mm 75C/5000psi liq .2ul N60/H
ACI4W.5	4 port 2-pos int vol valve air actuated	1/16 x .4mm 75C/5000psi liq .5ul N60/H
ACSD10UW	10 mpos dead-end path air actuated	1/16 x .4mm 75C/5000psi liq N60/E
ACSD4UW	OBSOLETE USE AHCS4UW	1/16 x .4mm 75C/5000psi liq N60/E
ACSD6UW	OBSOLETE USE AHCS6UW	1/16 x 1mm 200C/200psi gas N60/E
ACST4UW	4 mpos trapping path air actuated	1/16 x 1mm 200C/200psi gas N60/E
ACST6UW	6 mpos trapping path air actuated	1/16 x .75mm 200C/200psi gas N60/E
AH2CSC6MWE	6 mpos common outlet pathMW/air act.high tor/2in so	1/16 x .75mm 200C/200psi gas N60/E
AH2CSF6MWE	6 mpos flow-through MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2CST6MWE	6 mpos trapping path MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2CSTF6MWE	6 mpos trap/flow path MW air act high t./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2SC6MWE	6 mpos common outlet pathMWair act high tor/2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2SF6MWE	6 mpos flow-through MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2ST6MWE	6 mpos trapping path MW air act high torq./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2STF6MWE	6 mpos trap/flow path MW air act high t./2in so	1/8 x 1mm 200C/200psi gas N60/E
AH2VLS4MWE2	4 mpos common outlet path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS6MWE2	6 mpos common outlet path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS8MWE2	8 mpos common outlet path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS10MWE2	10 mpos dead-end path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS4MWE2	4 mpos dead-end path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS6MWE2	6 mpos dead-end path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLS8MWE2	8 mpos dead-end path MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLSF4MWE2	4 mpos flow-through MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLSF6MWE2	6 mpos flow-through MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AH2VLSF8MWE2	8 mpos flow-through MW air act/2in so	1/4 x 4mm 75C/100psi gas N60/E2
AHCS4UW	4 mpos dead-end path air actuated high torque	1/16 x .4mm 75C/5000psi liq N60/E
AHCS6UW	6 mpos dead-end path air actuated hi torque	1/16 x .4mm 75C/5000psi liq N60/E
AHSD4UW	4 mpos dead-end path air actuated high torque	1/8 x .75mm 75C/5000psi liq N60/E
AHSD6UW	6 mpos dead-end path air actuated high torque	1/8 x .75mm 75C/5000psi liq N60/E
AL10UW	10 port 2-pos valve air actuated	1/8 x 1.0mm 75C/5000psi liq N60/H
AL4UW	4 port 2-pos valve air actuated	1/8 x 1.7mm 75C/5000psi liq N60/H
AL6UW	6 port 2-pos valve air actuated	1/8 x 1.7mm 75C/5000psi liq N60/H
AL8UW	8 port 2-pos valve air actuated	1/8 x 1.3mm 75C/5000psi liq N60/H
ASD4UW	OBSOLETE USE AHSD4UW	
ASD6UW	OBSOLETE USE AHSD6UW	
ASD8UW	8 mpos dead-end path air actuated	1/8 x .75mm 75C/5000psi liq N60/E
ASFV	Prime/Purge valve air actuated/1in so	3 x 1/16in .50 mm 100C/10'000psi SS
ASFV2HT	Prime/Purge valve air actuated/1in so	3 x 1/8in 1.50 mm 300C/2'000psi SS
ASFV2HT4	Prime/Purge valve air actuated/4in so	3 x 1/8in 1.50 mm 300C/2'000psi SS
ASFV40K.5	Prime/purge valve air actuated/1.625in so	3 x 1/32in .15 mm bore 50C/40'000 psi SS
ASFV40K1	Prime/purge valve air actuated/1.625in so	3 x 1/16" .15 mm bore 50C/40'000 psi SS
ASFV40K360	Prime/purge valve air actuated/1.625in so	3 x 360um .15 mm bore 50C/40'000 psi SS
ASFVHT	Prime/Purge valve air actuated/1in so	3 x 1/16in .50 mm 300C/6'000psi SS

# VICI Valco-Cheminert PRODUCT LIST 2012




ASFVHT4	Prime/Purge valve air actuated/4in so	3 x 1/16in .50 mm 300C/6'000psi SS
ASFVL	Prime/Purge valve air actuated/1in so	3 x 1/16in .75 mm 100C/10'000psi SS
ASFVLHT	Prime/Purge valve air actuated/1in so	3 x 1/16in .75 mm 300C/6'000psi SS
ASFVLHT4	Prime/Purge valve air actuated/4in so	3 x 1/16in .75 mm 300C/6'000psi SS
ASFVO	On/Off valve air actuated/1in so	2 x 1/16in .50 mm bore 100C/10'000 psi SS
ASFVO2HT	On/Off valve air actuated/1in so	2 x 1/8in 1.50 mm bore 300C/2'000 psi SS
ASFVO2HT4	On/Off valve air actuated 4in so	2 x 1/8in 1.50 mm bore 300C/2'000 psi SS
ASFVO40K.5	On/Off valve air actuated/1.625in so	2 x 1/32in .15 mm bore 50C/40'000 psi SS
ASFVO40K1	On/Off valve air actuated/1.625in so	2 x 1/16" .15 mm bore 50C/40'000 psi SS
ASFVO40K360	On/Off valve air actuated/1.625in so	2 x 360um .15 mm bore 50C/40'000 psi SS
ASFVOD40K.5	Dual On/Off valve air actuated/1.625in so	3 x 1/32in .15 mm bore 50C/40'000 psi SS
ASFVOD40K1	Dual On/Off valve air actuated/1.625in so	3 x 1/16in .15 mm bore 50C/40'000 psi SS
ASFVOD40K360	Dual On/Off valve air actuated/1.625in so	3 x 360um .15 mm bore 50C/40'000 psi SS
ASFVOHT	On/Off valve air actuated/1in so	2 x 1/16in .50 mm bore 300C/6'000 psi SS
ASFVOHT4	On/Off valve air actuated 4in so	2 x 1/16in .50 mm bore 300C/6'000psi SS
ASFVOL	On/Off valve air actuated/1in so	2 x 1/16in .75 mm bore 100C/10'000 psi SS
ASFVOLHT	On/Off valve air actuated/1in so	2 x 1/16in .75 mm bore 300C/6'000 psi SS
ASFVOLHT4	On/Off valve air actuated 4in so	2 x 1/16in .75 mm bore 300C/6'000 psi SS
AT10	Air actuator 10 mpos high T close mount	36 deg rotation
AT104	Air actuator 10 mpos high T 4in standoff	36 deg rotation
AT10S	Air actuator 10 mpos high T for standoff	36 deg rotation
AT12	Air actuator 12 mpos high T close mount	30 deg rotation
AT124	Air actuator 12 mpos high T 4in standoff	30 deg rotation
AT12S	Air actuator 12 mpos high T for standoff	30 deg rotation
AT16	Air actuator 16 mpos high T close mount	22.5 deg rotation
AT164	Air actuator 16 mpos high T 4in standoff	22.5 deg rotation
AT16S	Air actuator 16 mpos high T for standoff	22.5 deg rotation
AT30	Air actuator 2-pos high T close mount	30 deg rotation
AT302	Air actuator 2-pos high T 2in standoff	30 deg rotation
AT30S	Air actuator 2-pos high T for standoff	30 deg rotation
AT36	Air actuator 2-pos high T close mount	36 deg rotation
AT362	Air actuator 2-pos high T 2in standoff	36 deg rotation
AT36S	Air actuator 2-pos high T for standoff	36 deg rotation
AT4	Air actuator 4 mpos high T close mount	90 deg rotation
AT44	Air actuator 4 mpos high T 4in standoff	90 deg rotation
AT45	Air actuator 2-pos high T close mount	45 deg rotation
AT452	Air actuator 2-pos high T 2in standoff	45 deg rotation
AT45S	Air actuator 2-pos high T for standoff	45 deg rotation
AT4S	Air actuator 4 mpos high T for standoff	90 deg rotation
AT6	Air actuator 6 mpos high T close mount	60 deg rotation
AT60	Air actuator 2-pos high T close mount	60 deg rotation
AT602	Air actuator 2-pos high T 2in standoff	60 deg rotation
AT60S	Air actuator 2-pos high T for standoff	60 deg rotation
AT64	Air actuator 6 mpos high T 4in standoff	60 deg rotation
AT6S	Air actuator 6 mpos high T for standoff	60 deg rotation
AT8	Air actuator 8 mpos high T close mount	45 deg rotation
AT84	Air actuator 8 mpos high T 4in standoff	45 deg rotation
AT8S	Air actuator 8 mpos high T for standoff	45 deg rotation
AT90	Air actuator 2-pos high T close mount	90 deg rotation
AT902	Air actuator 2-pos high T 2in standoff	90 deg rotation
AT90S	Air actuator 2-pos high T for standoff	90 deg rotation
BNV1	Micrometering needle valve bulkhead 2-225 ml/min	1/16in tube 18in /45cm2um frit 225Clubr. SS/Viton
BNV1-D	Micrometering needle valve bulkhead 2-225 ml/min	1/16in tube 18in /45cmnon-lubricated225C SS/Viton
BNV1-KZ	Micrometering needle valve bulkhead 2-225 ml/min	1/16in tube 18in /45cmnon-lubricated315CSS/Kalrez
BNV1LF	Micrometering needle valve bulkhead 2-90 ml/min	1/16in tube 18in /45cm2um frit225Clubr. SS/Viton
BNV1LF-D	Micrometering needle valve bulkhead 2-90 ml/min	1/16in tube 18in /45cm2um frit225Cnon-lubr.SS/Viton
BNV1LF-KZ	Micrometering needle valve bulkhead 2-90 ml/min	1/16in tube 18in /45cm2um frit315Cnon-lub.SS/Kalrez
C10UW	!!! RECOMMEND 2C10UW	SUPPLIED WITH LARGER HANDLE
C10UWE	10 port 2-pos valve manual	1/16 x .75mm 225C/400psi gas N60/E
C10W	10 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C10WE	10 port 2-pos valve manual	1/16 x .4mm 225C/400psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




C1-1006	6 port 2-pos vlv man/front loading+PFC	1/16 x .25mm 75C/5000psi liq N60/H
C1-10R6	Rotor seal 6 port 2-pos valve	1/16 x .25mm 75C/5000psi liq H
C1-1346	6 port 2-pos vlv man/front loading+PFC	1/16 x .25mm 50C/5000psi liq PAEK/E
C1-13R6	Rotor seal 6 port 2-pos valve	1/16 x .25mm 75C/5000psi liq E
C1-2006	6 port 2-pos vlv man/front loading+PFC	1/16 x .4mm 75C/5000psi liq N60/H
C12-310	Rotor seal 10 port 2-pos C22/C22Z-valve	1/16 x .75mm tube 75C/250psi liq E2
C12-314	Rotor seal 4 port 2-pos C22/C22Z-valve	1/16 x .75mm tube 75C/250psi liq E2
C12-316	Rotor seal 6 port 2-pos C22/C22Z-valve	1/16 x .75mm tube 75C/250psi liq E2
C12-318	Rotor seal 8 port 2-pos C22/C22Z-valve	1/16 x .75mm tube 75C/250psi liq E2
C1-2346	6 port 2-pos vlv man/front loading+PFC	1/16 x .4mm 50C/5000psi liq PAEK/E
C1-23R6	Rotor seal 6 port 2-pos valve PAEK	1/16 x .4mm 50C/5000psi liq E
C15-310	Rotor seal mpos dead-end path C15/C25Z-v	1/16 x .75mm 75C/250psi liq E2
C-1C00H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .25mm N60
C-1C04	Stator 4 port 2-pos C2-valve	1/16 x .25mm N60
C-1C06	Stator 6 port 2-pos C2-valve	1/16 x .25mm N60
C-1C08H	Stator 8 port 2-pos (.180BC) C2-valve	1/16 x .25mm N60
C-1C30H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .25mm Ti
C-1C34	Stator 4 port 2-pos C2-valve	1/16 x .25mm Ti
C-1C36	Stator 6 port 2-pos C1/C2-valve	1/16 x .25mm Ti
C-1C38H	Stator 8 port 2-pos (.180BC) C2-valve	1/16 x .25mm Ti
C-1C40H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .25mm PAEK
C-1C44	Stator 4 port 2-pos C2-valve	1/16 x .25mm PAEK
C-1C46	Stator 6 port 2-pos C1/C2-valve	1/16 x .25mm PAEK
C-1C48H	Stator 8 port 2-pos (.180) C2-valve	1/16 x .25mm PAEK
C1CF-1006	6 port 2-pos vlv man/front loading cont flow+PFC	1/16 x .25mm 75C/5000psi liq N60/H
C1CF-1346	6 port 2-pos vlv man/front loading cont flow+PFC	1/16 x .25mm 50C/5000psi liq PAEK/E
C1CF-2006	6 port 2-pos vlv man/front loading cont flow+PFC	1/16 x .4mm 75C/5000psi liq N60/H
C1CF-2346	6 port 2-pos vlv man/front loading cont flow+PFC	1/16 x .4mm 75C/5000psi liq PAEK/E
C2-1004	4 port 2-pos valve manual	1/16 x .25mm 75C/5000psi liq N60/H
C2-1004A	4 port 2-pos valve air actuated	1/16 x .25mm 75C/5000psi liq N60/H
C2-1004D	4 port 2-pos valve replacement	1/16 x .25mm 75C/5000psi liq N60/H
C2-1004E	4 port 2-pos valve std-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2-1004EH	4 port 2-pos valve m-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2-1006	6 port 2-pos valve manual	1/16 x .25mm 75C/5000psi liq N60/H
C2-1006A	6 port 2-pos valve air actuated	1/16 x .25mm 75C/5000psi liq N60/H
C2-1006D	6 port 2-pos valve replacement	1/16 x .25mm 75C/5000psi liq N60/H
C2-1006E	6 port 2-pos valve std-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2-1006EH	6 port 2-pos valve m-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2-1034	4 port 2-pos valve manual	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1034A	4 port 2-pos valve air actuated	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1034D	4 port 2-pos valve replacement	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1034E	4 port 2-pos valve std-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1034EH	4 port 2-pos valve m-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1036	6 port 2-pos valve manual	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1036A	6 port 2-pos valve air actuated	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1036D	6 port 2-pos valve replacement	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1036E	6 port 2-pos valve std-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2-1036EH	6 port 2-pos valve m-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2-10R0H	Rotor seal 10 port 2-pos (.180BC) C2-valve	1/16 x .25mm 75C/5000psi liq H
C2-10R4	Rotor seal 4 port 2-pos C2-valve	1/16 x .25mm 75C/5000psi liq H
C2-10R6	Rotor seal 6 port 2-pos C2/C3-valve	1/16 x .25mm 75C/5000psi liq H
C2-10R8H	Rotor seal 8 port 2-pos (.180BC) C2-valve	1/16 x .25mm 75C/5000psi liq H
C2-1344	4 port 2-pos valve manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1344A	4 port 2-pos valve air actuated	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1344D	4 port 2-pos valve replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1344E	4 port 2-pos valve std-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1344EH	4 port 2-pos valve m-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1346	6 port 2-pos valve manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1346A	6 port 2-pos valve air actuated	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1346D	6 port 2-pos valve replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-1346E	6 port 2-pos valve std-electric	1/16 x .25mm 50C/5000psi liq PAEK/E

# VICI Valco-Cheminert PRODUCT LIST 2012




C2-1346EH	6 port 2-pos valve m-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2-13R0H	Rotor seal 10 port 2-pos C2-valve (B.C.180)	1/16 x .25mm 50C/5000psi liq E
C2-13R4	Rotor seal 4 port 2-pos C2-valve	1/16 x .25mm 50C/5000psi liq E
C2-13R6	Rotor seal 6 port 2-pos C2/C3-valve	1/16 x .25mm 50C/5000psi liq E
C2-13R8H	Rotor seal 8 port 2-pos (.180BC) C2-valve	1/16 x .25mm 50C/5000psi liq E
C2-2004	4 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C2-2004A	4 port 2-pos valve air act	1/16 x .4mm 75C/5000psi liq N60/H
C2-2004D	4 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
C2-2004E	4 port 2-pos valve std-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2-2004EH	4 port 2-pos valve m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2-2006	6 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C2-2006A	6 port 2-pos valve air act	1/16 x .4mm 75C/5000psi liq N60/H
C2-2006D	6 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
C2-2006E	6 port 2-pos valve std-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2-2006EH	6 port 2-pos valve m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2-2034	4 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2034A	4 port 2-pos valve air act	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2034D	4 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2034E	4 port 2-pos valve std-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2034EH	4 port 2-pos valve m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2036	6 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2036A	6 port 2-pos valve air act	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2036D	6 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2036E	6 port 2-pos valve std-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C2-2036EH	6 port 2-pos valve m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C2-20R0H	Rotor seal 10 port 2-pos (.180BC) C2-valve	1/16 x .4mm 75C/5000psi liq H
C2-20R4	Rotor seal 4 port 2-pos C2-valve	1/16 x .4mm 75C/5000psi liq H
C2-20R6	Rotor seal 6 port 2-pos C2/C3-valve	1/16 x .4mm 75C/5000psi liq H
C2-20R8H	Rotor seal 8 port 2-pos (.180BC) C2-valve	1/16 x .4mm (.180 BC) 75C/5000psi liq H
C22-310	Rotor seal 10 port 2-pos C12/C22Z-valve	1/4in -28 x .75mm 75C/250psi liq E2
C22-314	Rotor seal 4 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq E2
C22-316	Rotor seal 6 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq E2
C22-318	Rotor seal 8 port 2-pos C12/C22Z-valve	1/4in -28 x .75mm 75C/250psi liq E2
C22-3180	10 port 2-pos valve manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3180A	10 port 2-pos valve air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3180D	10 port 2-pos valve replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3180E	10 port 2-pos valve std-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3180EH	10 port 2-pos valve m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3184	4 port 2-pos valve manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3184A	4 port 2-pos valve air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3184D	4 port 2-pos valve replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3184E	4 port 2-pos valve std-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3184EH	4 port 2-pos valve m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3186	6 port 2-pos valve manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3186A	6 port 2-pos valve air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3186D	6 port 2-pos valve replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3186E	6 port 2-pos valve std-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3186EH	6 port 2-pos valve m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3188	8 port 2-pos valve manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3188A	8 port 2-pos valve air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3188D	8 port 2-pos valve replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3188E	8 port 2-pos valve std-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C22-3188EH	8 port 2-pos valve m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C2-2344	4 port 2-pos valve manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2344A	4 port 2-pos valve air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2344D	4 port 2-pos valve replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2344E	4 port 2-pos valve std-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2344EH	4 port 2-pos valve m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2346	6 port 2-pos valve manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2346A	6 port 2-pos valve air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2346D	6 port 2-pos valve replacement	1/16 x .4mm 50C/5000psi liq PAEK/E

# VICI Valco-Cheminert PRODUCT LIST 2012




C2-2346E	6 port 2-pos valve std-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C2-2346EH	6 port 2-pos valve m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C22-380	Stator 10 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C22-384	Stator 4 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C22-386	Stator 6 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C22-388	Stator 8 port 2-pos C22-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C2-23R0H	Rotor seal 10 port 2-pos (.180BC) C2-valve	1/16 x .4mm 50C/5000psi liq E
C2-23R4	Rotor seal 4 port 2-pos C2-valve	1/16 x .4mm 50C/5000psi liq E
C2-23R6	Rotor seal 6 port 2-pos C2/C3-valve	1/16 x .4mm 50C/5000psi liq E
C2-23R8H	Rotor seal 8 port 2-pos (.180BC) C2-valve	1/16 x .4mm(.180 BC) 50C/5000psi liq E
C22-610	Rotor seal 10 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C22-614	Rotor seal 4 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C22-616	Rotor seal 6 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C22-618	Rotor seal 8 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C22-6180	10 port 2-pos valve manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6180A	10 port 2-pos valve air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6180D	10 port 2-pos valve replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6180E	10 port 2-pos valve std-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6180EH	10 port 2-pos valve m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6184	4 port 2-pos valve manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6184A	4 port 2-pos valve air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6184D	4 port 2-pos valve replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6184E	4 port 2-pos valve std-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6184EH	4 port 2-pos valve m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6186	6 port 2-pos valve manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6186A	6 port 2-pos valve air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6186D	6 port 2-pos valve replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6186E	6 port 2-pos valve std-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6186EH	6 port 2-pos valve m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6188	8 port 2-pos valve manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6188A	8 port 2-pos valve air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6188D	8 port 2-pos valve replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6188E	8 port 2-pos valve std-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-6188EH	8 port 2-pos valve m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C22-680	Stator 10 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C22-684	Stator 4 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C22-686	Stator 6 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C22-688	Stator 8 port 2-pos C22-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C22Z-3180	10 port 2-pos valve manual	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3180A	10 port 2-pos valve air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3180D	10 port 2-pos valve replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3180E	10 port 2-pos valve std-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3180EH	10 port 2-pos valve m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3184	4 port 2-pos valve manual	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3184A	4 port 2-pos valve air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3184D	4 port 2-pos valve replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3184E	4 port 2-pos valve std-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3184EH	4 port 2-pos valve m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3186	6 port 2-pos valve manual	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3186A	6 port 2-pos valve air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3186D	6 port 2-pos valve replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3186E	6 port 2-pos valve std-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3186EH	6 port 2-pos valve m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3188	8 port 2-pos valve manual	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3188A	8 port 2-pos valve air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3188D	8 port 2-pos valve replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3188E	8 port 2-pos valve std-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-3188EH	8 port 2-pos valve m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C22Z-380	Stator 10 port 2-pos C22Z-valve	1/16 x .75mm 75C/250psi liq PPS
C22Z-384	Stator 4 port 2-pos C22Z-valve	1/16 x .75mm 75C/250psi liq PPS
C22Z-386	Stator 6 port 2-pos C22Z-valve	1/16 x .75mm 75C/250psi liq PPS

# VICI Valco-Cheminert PRODUCT LIST 2012




C22Z-388	Stator 8 port 2-pos C22Z-valve	1/16 x .75mm 75C/250psi liq PPS
C24-10R-.2	Rotor seal 4port 2-pos int vol C24-valve	1/4in -28 x .5mm 75C/250psi liq .2ul E2
C24-10R-.5	Rotor seal 4port 2-pos int vol C24-valve	1/4in -28 x .5mm 75C/250psi liq .5ul E2
C24-10R-1	Rotor seal 4port 2-pos int vol C24-valve	1/4in -28 x .5mm 75C/250psi liq 1ul E2
C24-10R-2	Rotor seal 4port 2-pos int vol C24-valve	1/4in -28 x .5mm 75C/250psi liq 2ul E2
C24-1C8	Stator 4 port 2-pos int vol C24-valve	1/4in -28 x .5mm 75C/250psi liq PPS
C24-2184-.5	4 port 2-pos int vol valve manual	1/4in -28 x .5mm 75C/250psi liq.5ulPPS/E2
C24-2184-.5A	4 port 2-pos int vol valve air act	1/4in -28 x .5mm 75C/250psi liq.5ulPPS/E2
C24-2184-.5D	4 port 2-pos int vol valve replacement	1/4in -28 x .5mm 75C/250psi liq.5ulPPS/E2
C24-2184-.5E	4 port 2-pos int vol valve std-elect	1/4in -28 x .5mm 75C/250psi liq.5ulPPS/E2
C24-2184-.5EH	4 port 2-pos int vol valve m-elect	1/4in -28 x .5mm 75C/250psi liq.5ulPPS/E2
C24-2184-1	4 port 2-pos int vol valve manual	1/4in -28 x .5mm 75C/250psi liq 1ulPPS/E2
C24-2184-1A	4 port 2-pos int vol valve air act	1/4in -28 x .5mm 75C/250psi liq 1ulPPS/E2
C24-2184-1D	4 port 2-pos int vol valve replacement	1/4in -28 x .5mm 75C/250psi liq 1ulPPS/E2
C24-2184-1E	4 port 2-pos int vol valve std-elect	1/4in -28 x .5mm 75C/250psi liq 1ulPPS/E2
C24-2184-1EH	4 port 2-pos int vol valve m-elect	1/4in -28 x .5mm 75C/250psi liq 1ulPPS/E2
C24-2184-2	4 port 2-pos int vol valve manual	1/4in -28 x .5mm 75C/250psi liq 2ulPPS/E2
C24-2184-2A	4 port 2-pos int vol valve air act	1/4in -28 x .5mm 75C/250psi liq 2ulPPS/E2
C24-2184-2D	4 port 2-pos int vol valve replacement	1/4in -28 x .5mm 75C/250psi liq 2ulPPS/E2
C24-2184-2E	4 port 2-pos int vol valve std-elect	1/4in -28 x .5mm 75C/250psi liq 2ulPPS/E2
C24-2184-2EH	4 port 2-pos int vol valve m-elect	1/4in -28 x .5mm 75C/250psi liq 2ulPPS/E2
C24Z-1C8	Stator 4 port 2-pos int vol C24Z-valve	1/16 x .4mm 75C/250psi liq PPS
C24Z-2184-.2	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/250psi liq .2ul PPS/E2
C24Z-2184-.2A	4 port 2-pos int vol valve air act	1/16 x .4mm 75C/250psi liq .2ul PPS/E2
C24Z-2184-.2D	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/250psi liq .2ul PPS/E2
C24Z-2184-.2E	4 port 2-pos int vol valve std-elect	1/16 x .4mm 75C/250psi liq .2ul PPS/E2
C24Z-2184-.2EH	4 port 2-pos int vol valve m-elect	1/16 x .4mm 75C/250psi liq .2ul PPS/E2
C24Z-2184-.5	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/250psi liq .5ul PPS/E2
C24Z-2184-.5A	4 port 2-pos int vol valve air act	1/16 x .4mm 75C/250psi liq .5ul PPS/E2
C24Z-2184-.5D	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/250psi liq .5ul PPS/E2
C24Z-2184-.5E	4 port 2-pos int vol valve std-elect	1/16 x .4mm 75C/250psi liq .5ul PPS/E2
C24Z-2184-.5EH	4 port 2-pos int vol valve m-elect	1/16 x .4mm 75C/250psi liq .5ul PPS/E2
C24Z-2184-1	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/250psi liq 1ul PPS/E2
C24Z-2184-1A	4 port 2-pos int vol valve air act	1/16 x .4mm 75C/250psi liq 1ul PPS/E2
C24Z-2184-1D	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/250psi liq 1ul PPS/E2
C24Z-2184-1E	4 port 2-pos int vol valve std-elect	1/16 x .4mm 75C/250psi liq 1ul PPS/E2
C24Z-2184-1EH	4 port 2-pos int vol valve m-elect	1/16 x .4mm 75C/250psi liq 1ul PPS/E2
C25-310	Rotor seal 8/10 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq E2
C25-314	Rotor seal 4 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq E2
C25-316	Rotor seal 6 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq E2
C25-318	OBSOLETE USE PART-NO. C25-310	
C25-3180	10 mpos dead-end path manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3180A	10 mpos dead-end path air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3180D	10 mpos dead-end path replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3180EMH	10 mpos dead-end path m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3184	4 mpos dead-end path manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3184A	4 mpos dead-end path air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3184D	4 mpos dead-end path replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3184E	OBSOLETE USE C25-3184EMH	
C25-3184EMH	4 mpos dead-end path m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3186	6 mpos dead-end path manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3186A	6 mpos dead-end path air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3186D	6 mpos dead-end path replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3186E	OBSOLETE USE C25-3186EMH	
C25-3186EMH	6 mpos dead-end path m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3188	8 mpos dead-end path manual	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3188A	8 mpos dead-end path air actuated	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3188D	8 mpos dead-end path replacement	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-3188E	OBSOLETE USE C25-3188EMH	
C25-3188EMH	8 mpos dead-end path m-electric	1/4in -28 x .75mm 75C/250psi liq PPS/E2
C25-380	Stator 10 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq PPS

# VICI Valco-Cheminert PRODUCT LIST 2012




C25-384	Stator 4 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C25-386	Stator 6 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C25-388	Stator 8 mpos dead-end path C25-valve	1/4in -28 x .75mm 75C/250psi liq PPS
C25-610	Rotor seal 10 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C25-614	Rotor seal 4 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C25-616	Rotor seal 6 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C25-618	Rotor seal 8 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq E2
C25-6180	10 mpos dead-end path manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6180A	10 mpos dead-end path air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6180D	10 mpos dead-end path replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6180E	OBSOLETE USE C25-6180EMH	
C25-6180EMH	10 mpos dead-end path m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6184	4 mpos dead-end path manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6184A	4 mpos dead-end path air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6184D	4 mpos dead-end path replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6184E	OBSOLETE USE C25-6184EMH	
C25-6184EMH	4 mpos dead-end path m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6186	6 mpos dead-end path manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6186A	6 mpos dead-end path air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6186D	6 mpos dead-end path replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6186E	OBSOLETE USE C25-6186EMH	
C25-6186EMH	6 mpos dead-end path m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6188	8 mpos dead-end path manual	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6188A	8 mpos dead-end path air actuated	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6188D	8 mpos dead-end path replacement	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-6188E	OBSOLETE USE C25-6188EMH	
C25-6188EMH	8 mpos dead-end path m-electric	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C25-680	Stator 10 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C25-684	Stator 4 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C25-686	Stator 6 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C25-688	Stator 8 mpos dead-end path C25-valve	1/4in -28 x 1.5mm 75C/250psi liq PPS
C25Z-3180	10 mpos dead-end path manual	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3180A	10 mpos dead-end path air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3180D	10 mpos dead-end path replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3180E	OBSOLETE USE C25Z-3180EMH	
C25Z-3180EMH	10 mpos dead-end path m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-31814	14 mpos dead-end path manual	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-31814A	OBSOLETE	
C25Z-31814D	14 mpos dead-end path replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-31814E	OBSOLETE USE C25Z-31814EMH	
C25Z-31814EMH	14 mpos dead-end path m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3186	6 mpos dead-end path manual	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3186A	6 mpos dead-end path air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3186D	6 mpos dead-end path replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3186E	OBSOLETE USE C25Z-3186EMH	
C25Z-3186EMH	6 mpos dead-end path m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3188	8 mpos dead-end path manual	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3188A	8 mpos dead-end path air actuated	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3188D	8 mpos dead-end path replacement	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-3188E	OBSOLETE USE C25Z-3188EMH	
C25Z-3188EMH	8 mpos dead-end path m-electric	1/16 x .75mm 75C/250psi liq PPS/E2
C25Z-380	Stator 10 mpos dead-end path valve	1/16 x .75mm 75C/250psi liq PPS
C25Z-3814	Stator 14 mpos dead-end path valve	1/16 x .75mm 75C/250psi liq PPS
C25Z-38-14	Stator 14 mpos dead-end path valve	1/16 x .75mm 75C/250psi liq PPS
C25Z-386	Stator 6 mpos dead-end path valve	1/16 x .75mm 75C/250psi liq PPS
C25Z-388	Stator 8 mpos dead-end path valve	1/16 x .75mm 75C/250psi liq PPS
C-261	Syringe fitting assy. head knurled	C1 PEEK/TFE
C-2C00H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .4mm N60
C-2C04	Stator 4 port 2-pos C2-valve	1/16 x .4mm N60
C-2C06	Stator 6 port 2-pos C2-valve	1/16 x .4mm N60
C-2C08H	Stator 8 port 2-pos (.180BC) C2-valve	1/16 x .4mm N60

# VICI Valco-Cheminert PRODUCT LIST 2012




C-2C30H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .4mm Titanium
C-2C34	Stator 4 port 2-pos C2-valve	1/16 x .4mm Titanium
C-2C36	Stator 6 port 2-pos C1/C2-valve	1/16 x .4mm Titanium
C-2C38H	Stator 8 port 2-pos (.180BC) C2-valve	1/16 x .4mm (.180 BC) Titanium
C-2C40H	Stator 10 port 2-pos (.180BC) C2-valve	1/16 x .4mm PAEK
C-2C44	Stator 4 port 2-pos C2-valve	1/16 x .4mm PAEK
C-2C46	Stator 6 port 2-pos C1/C2-valve	1/16 x .4mm PAEK
C-2C48H	Stator 8 port 2-pos (.180BC) C2-valve	1/16 x .4mm (.180 BC) PAEK
C2H-1000	10 port 2-pos valve (.180BC) manual	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1000A	10 port 2-pos valve (.180BC) air actuated	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1000D	10 port 2-pos valve (.180BC) replacement	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1000E	10 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1000EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1008	8 port 2-pos valve (.180BC) manual	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1008A	8 port 2-pos valve (.180BC) air actuated	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1008D	8 port 2-pos valve (.180BC) replacement	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1008E	8 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1008EH	8 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 75C/5000psi liq N60/H
C2H-1030	10 port 2-pos valve (.180BC) manual	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1030A	10 port 2-pos valve (.180BC) air actuated	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1030D	10 port 2-pos valve (.180BC) replacement	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1030E	10 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1030EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1038	8 port 2-pos valve (.180BC) manual	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1038A	8 port 2-pos valve (.180BC) air actuated	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1038D	8 port 2-pos valve (.180BC) replacement	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1038E	8 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1038EH	8 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 75C/5000psi liq Ti/H
C2H-1340	10 port 2-pos valve (.180BC) manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1340A	10 port 2-pos valve (.180BC) air actuator	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1340D	10 port 2-pos valve (.180BC) replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1340E	10 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1340EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1348	8 port 2-pos valve (.180BC) manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1348A	8 port 2-pos valve (.180BC) air actuated	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1348D	8 port 2-pos valve (.180BC) replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1348E	8 port 2-pos valve (.180BC) std-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-1348EH	8 port 2-pos valve (.180BC) m-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C2H-2000	10 port 2-pos valve (.180BC) manual	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2000A	10 port 2-pos valve (.180BC) air act	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2000D	10 port 2-pos valve (.180BC) replacement	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2000E	10 port 2-pos valve (.180BC) std-electric	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2000EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2008	8 port 2-pos valve (.180BC) manual	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2008A	8 port 2-pos valve (.180BC) air act	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2008D	8 port 2-pos valve (.180BC) replacement	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2008E	8 port 2-pos valve (.180BC) std-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2008EH	8 port 2-pos valve (.180BC) m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C2H-2030	10 port 2-pos valve (.180BC) manual	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2030A	10 port 2-pos valve (.180BC) air act	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2030D	10 port 2-pos valve (.180BC) replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2030E	10 port 2-pos valve (.180BC) std-electric	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2030EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2038	8 port 2-pos valve (.180BC) manual	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2038A	8 port 2-pos valve (.180BC) air act	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2038D	8 port 2-pos valve (.180BC) replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2038E	8 port 2-pos valve (.180BC) std-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2038EH	8 port 2-pos valve (.180BC) m-electric	1/16 x .4mm 75C/5000psi liq Ti/H
C2H-2340	10 port 2-pos valve (.180BC) manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2340A	10 port 2-pos valve (.180BC) air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2340D	10 port 2-pos valve (.180BC) replacement	1/16 x .4mm 50C/5000psi liq PAEK/E

# VICI Valco-Cheminert PRODUCT LIST 2012




C2H-2340E	10 port 2-pos valve (.180BC) std-electric	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2340EH	10 port 2-pos valve (.180BC) m-electric	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2348	8 port 2-pos valve (.180BC) manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2348A	8 port 2-pos valve (.180BC) air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2348D	8 port 2-pos valve (.180BC) replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2348E	8 port 2-pos valve (.180BC) std-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C2H-2348EH	8 port 2-pos valve (.180BC) m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-1006	6 port 2-pos vlv manual vert. port injection	1/16 x .25mm 75C/5000psi liq N60/H
C2V-1006A	6 port 2-pos vlv air act. vert. port injection	1/16 x .25mm 75C/5000psi liq N60/H
C2V-1006D	6 port 2-pos vlv replacement vert. port injection	1/16 x .25mm 75C/5000psi liq N60/H
C2V-1006E	6 port 2-pos vlv std-act. vert. port injection	1/16 x .25mm 75C/5000psi liq N60/H
C2V-1006EH	6 port 2-pos vlv m-elect vert. port injection	1/16 x .25mm 75C/5000psi liq N60/H
C2V-1346	6 port 2-pos vlv manual vert. port injector	1/16 x .25mm 50C/5000psi liq PAEK/E
C2V-1346A	6 port 2-pos vlv air act vert. port injector	1/16 x .25mm 50C/5000psi liq PAEK/E
C2V-1346D	6 port 2-pos vlv replacement vert. port injector	1/16 x .25mm 50C/5000psi liq PAEK/E
C2V-1346E	6 port 2-pos vlv std-elect vert. port injector	1/16 x .25mm 50C/5000psi liq PAEK/E
C2V-1346EH	6 port 2-pos vlv m-elect vert. port injector	1/16 x .25mm 50C/5000psi liq PAEK/E
C2V-1C06	Stator 6 port 2-pos C2-valve vertical	1/16 x .25mm N60
C2V-1C46	Stator 6 port 2-pos C2-valve vertical	1/16 x .25mm PAEK
C2V-2006	6 port 2-pos vlv manual vert. port injector	1/16 x .4mm 75C/5000psi liq N60/H
C2V-2006A	6 port 2-pos vlv air. act. vert. port injector	1/16 x .4mm 75C/5000psi liq N60/H
C2V-2006D	6 port 2-pos vlv replacement vert. port injector	1/16 x .4mm 75C/5000psi liq N60/H
C2V-2006E	6 port 2-pos vlv std-elect vert. port injector	1/16 x .4mm 75C/5000psi liq N60/H
C2V-2006EH	6 port 2-pos vlv m-elect vert. port injector	1/16 x .4mm 75C/5000psi liq N60/H
C2V-2346	6 port 2-pos vlv manual vert. port injector	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-2346A	6 port 2-pos vlv air act. vert. port injector	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-2346D	6 port 2-pos vlv replacement vert. port injector	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-2346E	6 port 2-pos vlv std-elect vert. port injector	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-2346EH	6 port 2-pos vlv m-elect vert. port injector	1/16 x .4mm 50C/5000psi liq PAEK/E
C2V-2C06	Stator 6 port 2-pos C2-valve vertical	1/16 x .4mm N60
C2V-2C46	Stator 6 port 2-pos C2-valve vertical	1/16 x .4mm PAEK
C3-1006	6 port 2-pos valve manual	1/16 x .25mm 75C/5000psi liq N60/H
C3-1006A	6 port 2-pos valve air act	1/16 x .25mm 75C/5000psi liq N60/H
C3-1006D	6 port 2-pos valve replacement	1/16 x .25mm 75C/5000psi liq N60/H
C3-1006E	6 port 2-pos valve std-elect	1/16 x .25mm 75C/5000psi liq N60/H
C3-1006EH	6 port 2-pos valve m-elect	1/16 x .25mm 75C/5000psi liq N60/H
C3-1346	6 port 2-pos valve manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C3-1346A	6 port 2-pos valve air act	1/16 x .25mm 50C/5000psi liq PAEK/E
C3-1346D	6 port 2-pos valve replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C3-1346E	6 port 2-pos valve std-elect	1/16 x .25mm 50C/5000psi liq PAEK/E
C3-1346EH	6 port 2-pos valve m-elect	1/16 x .25mm 50C/5000psi liq PAEK/E
C3-1C06	Stator 6 port 2-pos C3-valve	1/16 x .25mm N60
C3-1C46	Stator 6 port 2-pos C3-valve	1/16 x .25mm PAEK
C3-2006	6 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C3-2006A	6 port 2-pos valve air actuated	1/16 x .4mm 75C/5000psi liq N60/H
C3-2006D	6 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
C3-2006E	6 port 2-pos valve std-elect	1/16 x .4mm 75C/5000psi liq N60/H
C3-2006EH	6 port 2-pos valve m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C3-2346	6 port 2-pos valve manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C3-2346A	6 port 2-pos valve air actuated	1/16 x .4mm 50C/5000psi liq PAEK/E
C3-2346D	6 port 2-pos valve replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C3-2346E	6 port 2-pos valve std-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C3-2346EH	6 port 2-pos valve m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C3-2C06	Stator 6 port 2-pos C3-valve	1/16 x .4mm N60
C3-2C46	Stator 6 port 2-pos C3-valve	1/16 x .4mm PAEK
C35Z-31820D	20 mpos dead-end path replacement	1/16 x .75mm 50C/100psi liq PPS/E2
C35Z-31820EMH	!!! RECOMMEND C35Z-31820EMT	
C35Z-31820EMT	20 mpos dead-end path m-electric	1/16 x .75mm 50C/100psi liq PPS/E2
C35Z-31824D	24 mpos dead-end path replacement	1/16 x .75mm 50C/100psi liq PPS/E2
C35Z-31824EMH	!!! RECOMMEND C35Z-31824EMT	
C35Z-31824EMT	24 mpos dead-end path m-electric	1/16 x .75mm 50C/100psi liq PPS/E2

# VICI Valco-Cheminert PRODUCT LIST 2012




C35Z-31826D	26 mpos dead-end path replacement	1/16 x .75mm 50C/100psi liq PPS/E2
C35Z-31826EMH	!!! RECOMMEND C35Z-31826EMT	
C35Z-31826EMT	26 mpos dead-end path m-electric	1/16 x .75mm 50C/100psi liq PPS/E2
C35Z-31R20	Rotor seal 20mpos dead-end path C35-vlv	1/16 x .75mm 50C/100psi liq E2
C360C	Cap 360um for metal tubing	SS nut included
C360CFS	Cap 360um for fused silica tubing	SS nut included
C360CPK	Cap 360um	PEEK
C360CPKG	OBSOLETE replaced by C360CPK	
C360ET	Make up extender tool	for 360um nuts
C360IZR.5S6	Internal reducer 1/32in x 360um	.38mm bore SS w/o frit
C360IZR.5S6PKG	Internal reducer 1/32in x 360um	.38mm bore SS with PEEK nut w/o frit
C360NFFS	Nut 360um No Twist	for FS tubing (for > 10 kpsi)
C360NFPKG	Nut 360um No Twist PEEK	with glass-filled PEEK ferrule
C360NFS6	Nut 360um No Twist SS	with 316 ferrule
C360PPK	Plug 360um	
C360QTPK2	Tee 360um quick mount 50um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QTPK4	Tee 360um quick mount 100um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QTPK6	Tee 360um quick mount 150um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QTPKG2	OBSOLETE replaced by C360QTPK2	
C360QTPKG4	OBSOLETE replaced by C360QTPK4	
C360QTPKG6	OBSOLETE replaced by C360QTPK6	
C360QUPK2	Union 360um quick mount 50um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QUPK4	Union 360um quick mount 100um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QUPK6	Union 360um quick mount 150um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QUPKG2	OBSOLETE replaced by C360QUPK2	
C360QUPKG4	OBSOLETE replaced by C360QUPK4	
C360QUPKG6	OBSOLETE replaced by C360QUPK6	
C360QXPK2	Cross 360um quick mount 50um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QXPK4	Cross 360um quick mount 100um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QXPK6	Cross 360um quick mount 150um bore	for PEEK and FS Tubing (for < 10 kpsi)
C360QXPKG2	OBSOLETE replaced by C360QXPK2	
C360QXPKG4	OBSOLETE replaced by C360QXPK4	
C360QXPKG6	OBSOLETE replaced by C360QXPK6	
C360RU.5FS2	Reducing union 1/32in to 360um	50um bore for FS tubing (for > 10 kpsi)
C360RU.5FS4	Reducing union 1/32in to 360um	100um bore for FS tubing (for > 10 kpsi)
C360RU.5FS6	Reducing union 1/32in to 360um	150um bore for FS tubing (for > 10 kpsi)
C360RU.5S62	Reducing union 1/32in to 360um	50um bore for metal tubing (up to 20 kpsi)
C360RU.5S64	Reducing union 1/32in to 360um	100um bore for metal tubing (up to 20 kpsi)
C360RU.5S66	Reducing union 1/32in to 360um	150um bore for metal tubing (up to 20 kpsi)
C360RU1FS6	Reducing union 1/16in to 360um	150um bore for FS tubing (for > 10 kpsi)
C360RU1S66	Reducing union 1/16in to 360um	150um bore for metal tubing (up to 20 kpsi)
C360UFS2	Union 360um 50um bore	for FS tubing (for > 10 kpsi)
C360UFS4	Union 360um 100um bore	for FS tubing (for > 10 kpsi)
C360UFS6	Union 360um 150um bore	for FS tubing (for > 10 kpsi)
C360UPK2	Union 360um 50um bore	for PEEK and FS tubing (for < 10 kpsi)
C360UPK4	Union 360um 100um bore	for PEEK and FS tubing (for < 10 kpsi)
C360UPK6	Union 360um 150um bore	for PEEK and FS tubing (for < 10 kpsi)
C360UPKG2	OBSOLETE replaced by C360UPK2	
C360UPKG4	OBSOLETE replaced by C360UPK4	
C360UPKG6	OBSOLETE replaced by C360UPK6	
C360US62	Union 360um 50um bore	for metal tubing (up to 20 kpsi)
C360US64	Union 360um 100um bore	for metal tubing (up to 20 kpsi)
C360US66	Union 360um 150um bore	for metal tubing (up to 20 kpsi)
C4-0004-.01	4 port 2-pos int vol valve manual	1/16 x .15mm 75C/5000psi liq .01ul N60/H
C4-0004-.01A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 75C/5000psi liq .01ul N60/H
C4-0004-.01D	4 port 2-pos int vol valve replacement	1/16 x .15mm 75C/5000psi liq .01ul N60/H
C4-0004-.01E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 75C/5000psi liq .01ul N60/H
C4-0004-.01EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 75C/5000psi liq .01ul N60/H
C4-0004-.02	4 port 2-pos int vol valve manual	1/16 x .15mm 75C/5000psi liq .02ul N60/H
C4-0004-.02A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 75C/5000psi liq .02ul N60/H
C4-0004-.02D	4 port 2-pos int vol valve replacement	1/16 x .15mm 75C/5000psi liq .02ul N60/H

# VICI Valco-Cheminert PRODUCT LIST 2012




C4-0004-.02E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 75C/5000psi liq .02ul N60/H
C4-0004-.02EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 75C/5000psi liq .02ul N60/H
C4-0004-.05	4 port 2-pos int vol valve manual	1/16 x .15mm 75C/5000psi liq .05ul N60/H
C4-0004-.05A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 75C/5000psi liq .05ul N60/H
C4-0004-.05D	4 port 2-pos int vol valve replacement	1/16 x .15mm 75C/5000psi liq .05ul N60/H
C4-0004-.05E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 75C/5000psi liq .05ul N60/H
C4-0004-.05EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 75C/5000psi liq .05ul N60/H
C4-00R-.01	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 75C/5000psi liq .01ul H
C4-00R-.02	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 75C/5000psi liq .02ul H
C4-00R-.05	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 75C/5000psi liq .05ul H
C4-0344-.01	4 port 2-pos int vol valve manual	1/16 x .15mm 50C/5000psi liq .01ul PAEK/E
C4-0344-.01A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 50C/5000psi liq .01ul PAEK/E
C4-0344-.01D	4 port 2-pos int vol valve replacement	1/16 x .15mm 50C/5000psi liq .01ul PAEK/E
C4-0344-.01E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 50C/5000psi liq .01ul PAEK/E
C4-0344-.01EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 50C/5000psi liq .01ul PAEK/E
C4-0344-.02	4 port 2-pos int vol valve manual	1/16 x .15mm 50C/5000psi liq .02ul PAEK/E
C4-0344-.02A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 50C/5000psi liq .02ul PAEK/E
C4-0344-.02D	4 port 2-pos int vol valve replacement	1/16 x .15mm 50C/5000psi liq .02ul PAEK/E
C4-0344-.02E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 50C/5000psi liq .02ul PAEK/E
C4-0344-.02EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 50C/5000psi liq .02ul PAEK/E
C4-0344-.05	4 port 2-pos int vol valve manual	1/16 x .15mm 50C/5000psi liq .05ul PAEK/E
C4-0344-.05A	4 port 2-pos int vol valve air actuated	1/16 x .15mm 50C/5000psi liq .05ul PAEK/E
C4-0344-.05D	4 port 2-pos int vol valve replacement	1/16 x .15mm 50C/5000psi liq .05ul PAEK/E
C4-0344-.05E	4 port 2-pos int vol valve std-electric	1/16 x .15mm 50C/5000psi liq .05ul PAEK/E
C4-0344-.05EH	4 port 2-pos int vol valve m-electric	1/16 x .15mm 50C/5000psi liq .05ul PAEK/E
C4-03R-.01	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 50C/5000psi liq .01ul E
C4-03R-.02	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 50C/5000psi liq .02ul E
C4-03R-.05	Rotor seal 4 port 2-pos int vol valve	1/16 x .15mm 50C/5000psi liq .05ul E
C4-0C0	Stator 4 port 2-pos int vol C4-valve	1/16 x .15mm N60
C4-0C4	Stator 4 port 2-pos int vol C4-valve	1/16 x .15mm PAEK
C4-1004-.1	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .1ul N60/H
C4-1004-.1A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .1ul N60/H
C4-1004-.1D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .1ul N60/H
C4-1004-.1E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .1ul N60/H
C4-1004-.1EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .1ul N60/H
C4-1004-.2	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .2ul N60/H
C4-1004-.2A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .2ul N60/H
C4-1004-.2D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .2ul N60/H
C4-1004-.2E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .2ul N60/H
C4-1004-.2EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .2ul N60/H
C4-1004-.5	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .5ul N60/H
C4-1004-.5A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .5ul N60/H
C4-1004-.5D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .5ul N60/H
C4-1004-.5E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .5ul N60/H
C4-1004-.5EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .5ul N60/H
C4-1034-.1	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .1ul Ti/H
C4-1034-.1A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .1ul Ti/H
C4-1034-.1D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .1ul Ti/H
C4-1034-.1E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .1ul Ti/H
C4-1034-.1EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .1ul Ti/H
C4-1034-.2	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .2ul Ti/H
C4-1034-.2A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .2ul Ti/H
C4-1034-.2D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .2ul Ti/H
C4-1034-.2E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .2ul Ti/H
C4-1034-.2EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .2ul Ti/H
C4-1034-.5	4 port 2-pos int vol valve manual	1/16 x .25mm 75C/5000psi liq .5ul Ti/H
C4-1034-.5A	4 port 2-pos int vol valve air act	1/16 x .25mm 75C/5000psi liq .5ul Ti/H
C4-1034-.5D	4 port 2-pos int vol valve replacement	1/16 x .25mm 75C/5000psi liq .5ul Ti/H
C4-1034-.5E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 75C/5000psi liq .5ul Ti/H
C4-1034-.5EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 75C/5000psi liq .5ul Ti/H
C4-10R-.1	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 75C/5000psi liq .1ul H

# VICI Valco-Cheminert PRODUCT LIST 2012




C4-10R-.2	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 75C/5000psi liq .2ul H
C4-10R-.5	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 75C/5000psi liq .5ul H
C4-1344-.1	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/5000psi liq.1ul PAEK/E
C4-1344-.1A	4 port 2-pos int vol valve air act	1/16 x .25mm 50C/5000psi liq.1ul PAEK/E
C4-1344-.1D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/5000psi liq.1ul PAEK/E
C4-1344-.1E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 50C/5000psi liq.1ul PAEK/E
C4-1344-.1EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 50C/5000psi liq.1ul PAEK/E
C4-1344-.2	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/5000psi liq.2ul PAEK/E
C4-1344-.2A	4 port 2-pos int vol valve air act	1/16 x .25mm 50C/5000psi liq.2ul PAEK/E
C4-1344-.2D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/5000psi liq.2ul PAEK/E
C4-1344-.2E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 50C/5000psi liq.2ul PAEK/E
C4-1344-.2EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 50C/5000psi liq.2ul PAEK/E
C4-1344-.5	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/5000psi liq.5ul PAEK/E
C4-1344-.5A	4 port 2-pos int vol valve air act	1/16 x .25mm 50C/5000psi liq.5ul PAEK/E
C4-1344-.5D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/5000psi liq.5ul PAEK/E
C4-1344-.5E	4 port 2-pos int vol valve std-elect	1/16 x .25mm 50C/5000psi liq.5ul PAEK/E
C4-1344-.5EH	4 port 2-pos int vol valve m-elect	1/16 x .25mm 50C/5000psi liq.5ul PAEK/E
C4-13R-.1	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 50C/5000psi liq .1ul E
C4-13R-.2	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 50C/5000psi liq .2ul E
C4-13R-.5	Rotor seal 4 port 2-pos int vol valve	1/16 x .25mm 50C/5000psi liq .5ul E
C4-1C0	Stator 4 port 2-pos int vol C4-valve	1/16 x .25mm N60
C4-1C3	Stator 4 port 2-pos int vol C4-valve	1/16 x .25mm Titanium
C4-1C4	Stator 4 port 2-pos int vol C4-valve	1/16 x .25mm PAEK
C45-9784A	4 mpos dead-end path air act	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9784D	4 mpos dead-end path replacement	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9784E	4 mpos dead-end path m-elect	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9784EMT	4 mpos dead-end path m-elect	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9786A	OBSOLETE USE C45-9786AH	
C45-9786D	6 mpos dead-end path replacement	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9786E	6 mpos dead-end path m-elect	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-9786EMT	6 mpos dead-end path m-elect	1/2in -20 x 4.6mm 50C/100psi liq PPS/PTFE
C45-97R4	Rotor 4 mpos C45-valve	1/2in -20 x 4.6mm 50C/100psi liq E2
C45-97R6	OBSOLETE USE C-316	
C4UW	!!! RECOMMEND 2C4UW	
C4UWE	4 port 2-pos valve manual	SUPPLIED WITH LARGER HANDLE
C4W	4 port 2-pos valve manual	1/16 x .75mm 225C/400psi gas N60/E
C4WE	4 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C5-2004	4 mpos dead-end path manual	1/16 x .4mm 225C/400psi gas N60/E
C5-2004A	4 mpos dead-end path air act	1/16 x .4mm 75C/5000psi liq N60/H
C5-2004D	4 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq N60/H
C5-2004E	OBSOLETE USE C5-2004EMH	
C5-2004EMH	4 mpos dead-end path m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C5-2006	6 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq N60/H
C5-2006A	6 mpos dead-end path air act	1/16 x .4mm 75C/5000psi liq N60/H
C5-2006D	6 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq N60/H
C5-2006E	OBSOLETE USE C5-2006EMH	
C5-2006EMH	6 mpos dead-end path m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C5-2006EMTD	6 mpos dead-end path m-elect dual drive	1/16 x .4mm 75C/5000psi liq N60/H
C5-2034	4 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2034A	4 mpos dead-end path air act	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2034D	4 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2034E	OBSOLETE USE C5-2034EMH	
C5-2034EMH	4 mpos dead-end path m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2036	6 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2036A	6 mpos dead-end path air act	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2036D	6 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C5-2036E	OBSOLETE USE C5-2036EMH	
C5-2036EMH	6 mpos dead-end path m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C5-20R0H	Rotor seal 10 mpos (.180BC) dead-end path valve	1/16 x .4mm 75C/5000psi liq H
C5-20R4	Rotor seal 4 mpos dead-end path valve	1/16 x .4mm 75C/5000psi liq H
C5-20R6	Rotor seal 6 mpos dead-end path valve	1/16 x .4mm 75C/5000psi liq H

# VICI Valco-Cheminert PRODUCT LIST 2012




C5-20R8H	Rotor seal 8 mpos (.180BC) dead-end path valve	1/16 x .4mm 75C/5000psi liq H
C52-1000I	10 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq N60/H
C52-1000I-S	10 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq N60/H
C52-1000IX	10 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq N60/H
C52-1004I	4 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq N60/H
C52-1004I-S	4 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq N60/H
C52-1004IX	4 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq N60/H
C52-1006I	6 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq N60/H
C52-1006I-S	6 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq N60/H
C52-1006IX	6 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq N60/H
C52-1008I	8 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq N60/H
C52-1008I-S	8 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq N60/H
C52-1008IX	8 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq N60/H
C52-1340I	10 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1340I-S	10 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1340IX	10 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1344I	4 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1344I-S	4 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1344IX	4 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1346I	6 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1346I-S	6 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1346IX	6 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1348I	8 port 2-pos valve motor sensor controller	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1348I-S	8 port 2-pos valve motor sensor	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1348IX	8 port 2-pos valve motor only	1/16 x .25mm 40C/5000psi liq PAEK/E
C52-1C00	Stator 10 port 2-pos C52-valve	1/16 x .25mm (.180) N60
C52-1C04	Stator 4 port 2-pos C52-valve	1/16 x .25mm (.180) N60
C52-1C06	Stator 6 port 2-pos C52-valve	1/16 x .25mm N60
C52-1C08	Stator 8 port 2-pos C52-valve	1/16 x .25mm (.180) N60
C52-1C40	Stator 10 port 2-pos C52-valve	1/16 x .25mm (.100) E3
C52-1C44	Stator 4 port 2-pos C52-valve	1/16 x .25mm (.100) E3
C52-1C46	Stator 6 port 2-pos C52-valve	1/16 x .25mm (.100) E3
C52-1C48	Stator 8 port 2-pos C52-valve	1/16 x .25mm (.100) E3
C52-2000I	10 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq N60/H
C52-2000I-S	10 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq N60/H
C52-2000IX	10 port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq N60/H
C52-2004I	4 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq N60/H
C52-2004I-S	4port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq N60/H
C52-2004IX	4port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq N60/H
C52-2006I	6 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq N60/H
C52-2006I-S	6 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq N60/H
C52-2006IX	6port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq N60/H
C52-2008I	8 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq N60/H
C52-2008I-S	8 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq N60/H
C52-2008IX	8port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq N60/H
C52-2340I	10 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2340I-S	10 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2340IX	10 port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2344I	4 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2344I-S	4 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2344IX	4 port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2346I	6 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2346I-S	6 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2346IX	6 port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2348I	8 port 2-pos valve motor sensor controller	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2348I-S	6 port 2-pos valve motor sensor	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2348IX	8 port 2-pos valve motor only	1/16 x .4mm 40C/5000psi liq PAEK/E
C52-2C00	Stator 10 port 2-pos C52-valve	1/16 x .4mm N60
C52-2C04	Stator 4 port 2-pos C52-valve	1/16 x .4mm N60
C52-2C06	Stator 6 port 2-pos C52-valve	1/16 x .4mm N60
C52-2C08	Stator 8 port 2-pos C52-valve	1/16 x .4mm N60

# VICI Valco-Cheminert PRODUCT LIST 2012




C52-2C40	Stator 10 port 2-pos C52-valve	1/16 x .4mm PAEK
C52-2C44	Stator 4 port 2-pos C52-valve	1/16 x .4mm PAEK
C52-2C46	Stator 6 port 2-pos C52-valve	1/16 x .4mm PAEK
C52-2C48	Stator 8 port 2-pos C52-valve	1/16 x .4mm PAEK
C5-2344	4 mpos dead-end path manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2344A	4 mpos dead-end path air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2344D	4 mpos dead-end path replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2344E	OBSOLETE USE C5-2344EMH	
C5-2344EMH	4 mpos dead-end path m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2346	6 mpos dead-end path manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2346A	6 mpos dead-end path air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2346D	6 mpos dead-end path replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2346E	OBSOLETE USE C5-2346EMH	
C5-2346EMH	6 mpos dead-end path m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-2346EMHD	!!! RECOMMEND C5-2346EMTD	
C5-2346EMTD	6 mpos dead-end path m-elect dual drive	1/16 x .4mm 50C/5000psi liq PAEK/E
C5-23R0H	Rotor seal 10 mpos (.180BC) dead-end path valve	1/16 x .4mm 50C/5000psi liq E
C5-23R4	Rotor seal 4 mpos dead-end path valve	1/16 x .4mm 50C/5000psi liq E
C5-23R6	Rotor seal 6 mpos dead-end path valve	1/16 x .4mm 50C/5000psi liq E
C5-23R8H	Rotor seal 8 mpos (.180BC) dead-end path valve	1/16 x .4mm 50C/5000psi liq E
C5-2C00H	Stator 10 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm N60
C5-2C04	Stator 4 mpos dead-end path C5-valve	1/16 x .4mm N60
C5-2C06	Stator 6 mpos dead-end path C5-valve	1/16 x .4mm N60
C5-2C08H	Stator 8 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm N60
C5-2C30H	Stator 10 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm Titanium
C5-2C34	Stator 4 mpos dead-end path C5-valve	1/16 x .4mm Titanium
C5-2C36	Stator 6 mpos dead-end path C5-valve	1/16 x .4mm Titanium
C5-2C38H	Stator 8 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm Titanium
C5-2C40H	Stator 10 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm PEAK
C5-2C44	Stator 4 mpos dead-end path C5-valve	1/16 x .4mm PAEK
C5-2C46	Stator 6 mpos dead-end path C5-valve	1/16 x .4mm PAEK
C5-2C48H	Stator 8 mpos (.180BC) dead-end path C5-valve	1/16 x .4mm PAEK
C55-2000I	10 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq N60/H
C55-2000I-S	10 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq N60/H
C55-2000IX	10 mpos valve motor only	1/16 x .4mm 50C/5000psi liq N60/H
C55-2004I	4 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq N60/H
C55-2004I-S	4 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq N60/H
C55-2004IX	4 mpos valve motor only	1/16 x .4mm 50C/5000psi liq N60/H
C55-2006I	6 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq N60/H
C55-2006I-S	6 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq N60/H
C55-2006IX	6 mpos valve motor only	1/16 x .4mm 50C/5000psi liq N60/H
C55-2008I	8 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq N60/H
C55-2008I-S	8 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq N60/H
C55-2008IX	8 mpos valve motor only	1/16 x .4mm 50C/5000psi liq N60/H
C55-2340I	10 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2340I-S	10 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2340IX	10 mpos valve motor only	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2344I	4 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2344I-S	4 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2344IX	4 mpos valve motor only	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2346I	6 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2346I-S	6 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2346IX	6 mpos valve motor only	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2348I	8 mpos valve motor sensor controller	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2348I-S	8 mpos valve motor sensor	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2348IX	8 mpos valve motor only	1/16 x .4mm 50C/5000psi liq PAEK/E
C55-2C00	Stator 10 mpos dead-end path C55-valve	1/16 x .4mm N60
C55-2C04	Stator 4 mpos dead-end path C55-valve	1/16 x .4mm N60
C55-2C06	Stator 6 mpos dead-end path C55-valve	1/16 x .4mm N60
C55-2C08	Stator 8 mpos dead-end path C55-valve	1/16 x .4mm N60
C55-2C40	Stator 10 mpos dead-end path C55-valve	1/16 x .4mm PAEK

# VICI Valco-Cheminert PRODUCT LIST 2012




C55-2C44	Stator 4 mpos dead-end path C55-valve	1/16 x .4mm PAEK
C55-2C46	Stator 6 mpos dead-end path C55-valve	1/16 x .4mm PAEK
C55-2C48	Stator 8 mpos dead-end path C55-valve	1/16 x .4mm PAEK
C5H-2000	10 mpos dead-end path (.180BC) manual	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2000A	10 mpos dead-end path (.180BC) air act	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2000D	10 mpos dead-end path (.180BC) replacement	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2000E	OBSOLETE USE C5H-2000EMT	
C5H-2000EMT	10 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2000EMTD	10 mpos dead-end path(.180BC) m-elect dual drive	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2008	8 mpos dead-end path (.180BC) manual	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2008A	8 mpos dead-end path (.180BC) air act	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2008D	8 mpos dead-end path (.180BC) replacement	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2008E	OBSOLETE USE C5H-2008EMT	
C5H-2008EMT	8 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2008EMTD	8 mpos dead-end path (.180BC) m-elect dual drive	1/16 x .4mm 75C/5000psi liq N60/H
C5H-2030	10 mpos dead-end path (.180BC) manual	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2030A	10 mpos dead-end path (.180BC) air act	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2030D	10 mpos dead-end path (.180BC) replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2030E	OBSOLETE USE C5H-2030EMT	
C5H-2030EMT	10 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2038	8 mpos dead-end path (.180BC) manual	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2038A	8 mpos dead-end path (.180BC) air act	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2038D	8 mpos dead-end path (.180BC) replacement	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2038E	OBSOLETE USE C5H-2038EMT	
C5H-2038EMT	8 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 75C/5000psi liq Ti/H
C5H-2340	10 mpos dead-end path (.180BC) manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2340A	10 mpos dead-end path (.180BC) air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2340D	10 mpos dead-end path (.180BC) replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2340E	OBSOLETE USE C5H-2340EMT	
C5H-2340EMT	10 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2340EMTD	10 mpos dead-end path(.180BC) m-elect dual drive	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2348	8 mpos dead-end path (.180BC) manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2348A	8 mpos dead-end path (.180BC) air act	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2348D	8 mpos dead-end path (.180BC) replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2348E	OBSOLETE USE C5H-2348EMT	
C5H-2348EMT	8 mpos dead-end path (.180BC) m-elect	1/16 x .4mm 50C/5000psi liq PAEK/E
C5H-2348EMTD	8 mpos dead-end path (.180BC) m-elect dual drive	1/16 x .4mm 50C/5000psi liq PAEK/E
C5M1KF	Cheminert manifold 1/4in -28 1/16in OD tubing	4 inlets 1 outlet .75mm bore CTFE
C5M1PK	Cheminert manifold 1/4in -28 1/16in OD tubing	4 inlets 1 outlet .75mm bore PEEK
C5M2KF	Cheminert manifold 1/4in -28 1/8in OD tubing	4 inlets 1 outlet 1.5mm bore CTFE
C5M2PK	Cheminert manifold 1/4in -28 1/8in OD tubing	4 inlets 1 outlet 1.5mm bore PEEK
C6-1006	6 port 2-pos valve cont. flow manual	1/16 x .25mm 75C/5000psi liq N60/H
C6-1006A	6 port 2-pos valve cont. flow air actuator	1/16 x .25mm 75C/5000psi liq N60/H
C6-1006D	6 port 2-pos valve cont. flow replacement	1/16 x .25mm 75C/5000psi liq N60/H
C6-1006E	6 port 2-pos valve cont. flow std. electric	1/16 x .25mm 75C/5000psi liq N60/H
C6-1006EH	6 port 2-pos valve cont. flow m-electric	1/16 x .25mm 75C/5000psi liq N60/H
C6-1346	6 port 2-pos valve cont. flow manual	1/16 x .25mm 50C/5000psi liq PAEK/E
C6-1346A	6 port 2-pos valve cont. flow air actuator	1/16 x .25mm 50C/5000psi liq PAEK/E
C6-1346D	6 port 2-pos valve cont. flow replacement	1/16 x .25mm 50C/5000psi liq PAEK/E
C6-1346E	6 port 2-pos valve cont. flow std. electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C6-1346EH	6 port 2-pos valve cont. flow m-electric	1/16 x .25mm 50C/5000psi liq PAEK/E
C6-1C06	Stator 6 port 2-pos C6-valve	1/16 x .25mm N60
C6-1C46	Stator 6 port 2-pos cont. flow valve	1/16 x .25mm PAEK
C6-2006	6 port 2-pos valve cont. flow manual	1/16 x .4mm 75C/5000psi liq N60/H
C6-2006A	6 port 2-pos valve cont. flow air actuator	1/16 x .4mm 75C/5000psi liq N60/H
C6-2006D	6 port 2-pos valve cont. flow replacement	1/16 x .4mm 75C/5000psi liq N60/H
C6-2006E	6 port 2-pos valve cont. flow std. electric	1/16 x .4mm 75C/5000psi liq N60/H
C6-2006EH	6 port 2-pos valve cont. flow m-electric	1/16 x .4mm 75C/5000psi liq N60/H
C62-3180I	10 port 2-pos valve motor sensor controller	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C62-3180I-S	10 port 2-pos valve motor sensor	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C62-3184I	4 port 2-pos valve motor sensor controller	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2

# VICI Valco-Cheminert PRODUCT LIST 2012




C62-3184I-S	4 port 2-pos valve motor sensor	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C62-3186I	6 port 2-pos valve motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C62-3186I-S	6 port 2-pos valve motor sensor	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C62-3188I	8 port 2-pos valve motor sensor controller	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C62-3188I-S	8 port 2-pos valve motor sensor	1/4in -28 x 0.75mm 75C/250psi liq PPS/E2
C6-2346	6 port 2-pos valve cont. flow manual	1/16 x .4mm 50C/5000psi liq PAEK/E
C6-2346A	6 port 2-pos valve cont. flow air actuator	1/16 x .4mm 50C/5000psi liq PAEK/E
C6-2346D	6 port 2-pos valve cont. flow replacement	1/16 x .4mm 50C/5000psi liq PAEK/E
C6-2346E	6 port 2-pos valve cont. flow std. electric	1/16 x .4mm 50C/5000psi liq PAEK/E
C6-2346EH	6 port 2-pos valve cont. flow m-electric	1/16 x .4mm 50C/5000psi liq PAEK/E
C62-6180I	10 port 2-pos valve motor sensor controller	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6180I-S	10 port 2-pos valve motor sensor	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6184I	4 port 2-pos valve motor sensor controller	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6184I-S	4 port 2-pos valve motor sensor	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6186I	6 port 2-pos valve motor sensor controller	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6186I-S	6 port 2-pos valve motor sensor	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C62-6188I	8 port 2-pos valve motor sensor controller	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C62-6188I-S	6 port 2-pos valve motor sensor	1/4in -28 x 1.5mm 75C/250psi liq PPS/E2
C6-2C06	Stator 6 port 2-pos cont. flow valve	1/16 x .4mm N60
C6-2C46	Stator 6 port 2-pos cont. flow valve	1/16 x .4mm PAEK
C62Z-3180I	10 port 2-pos valve motor sensor controller	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3180I-S	10 port 2-pos valve motor sensor	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3184I	4 port 2-pos valve motor sensor controller	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3184I-S	4 port 2-pos valve motor sensor	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3186I	6 port 2-pos valve motor sensor controller	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3186I-S	6 port 2-pos valve motor sensor	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3188I	8 port 2-pos valve motor sensor controller	1/16in x .75mm 50C/250psi liq PPS/E2
C62Z-3188I-S	8 port 2-pos valve motor sensor	1/16in x .75mm 50C/250psi liq PPS/E2
C65-3180I	10 mpos dead-end path motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3180I-S	10 mpos dead-end path motor sensor	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3184I	4 mpos dead-end path motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3184I-S	4 mpos dead-end path motor sensor	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3186I	6 mpos dead-end path motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3186I-S	6 mpos dead-end path motor sensor	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3188I	8 mpos dead-end path motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-3188I-S	8 mpos dead-end path motor sensor	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-6180I	10 mpos dead-end path motor sensor controller	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-6180I-S	10 mpos dead-end path motor sensor	1/4in -28 x 0.75mm 50C/250psi liq PPS/E2
C65-6184I	4 mpos dead-end path motor sensor controller	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65-6184I-S	4 mpos dead-end path motor sensor	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65-6186I	6 mpos dead-end path motor sensor controller	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65-6186I-S	6 mpos dead-end path motor sensor	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65-6188I	8 mpos dead-end path motor sensor controller	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65-6188I-S	8 mpos dead-end path motor sensor	1/4in -28 x 1.5mm 50C/250psi liq PPS/E2
C65Z-3180I	10 mpos dead-end path motorsensorw/o controller	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3180I-S	10 mpos dead-end pathmotor sensorw/o controller	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3184I	4 mpos dead-end path integrated actuator	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3184I-S	4 mpos dead-end path motor sensorw/o controller	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3186I	6 mpos dead-end path integrated actuator	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3186I-S	6 mpos dead-end path motor sensorw/o controller	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3188I	8 mpos dead-end path integrated actuator	1/16in x .75mm 50C/250psi liq PPS/E2
C65Z-3188I-S	8 mpos dead-end path motor sensorw/o controller	1/16in x .75mm 50C/250psi liq PPS/E2
C6UW	!!! RECOMMEND 2C6UW	SUPPLIED WITH LARGER HANDLE
C6UWE	6 port 2-pos valve manual	1/16 x .75mm 225C/400psi gas N60/E
C6W	6 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C6WE	6 port 2-pos valve manual	1/16 x .4mm 225C/400psi gas N60/E
C72-16R0	Rotor seal 10 port 2-pos C72-valve	1/16 x .25mm 50C/all pressures liqE3
C72-16R4	Rotor seal 4 port 2-pos C72-valve	1/16 x .25mm 50C/all pressures liqE3
C72-16R6	Rotor seal 6 port 2-pos C72-valve	1/16 x .25mm 50C/all pressures liqE3
C72-16R8	Rotor seal 8 port 2-pos C72-valve	1/16 x .25mm 50C/all pressures liqE3
C72-1C90	Stator 10 port 2-pos C72-valve	1/16 x .25mm N60C

# VICI Valco-Cheminert PRODUCT LIST 2012




C72-1C94	Stator 4 port 2-pos C72-valve	1/16 x .25mm N60C
C72-1C96	Stator 6 port 2-pos C72-valve	1/16 x .25mm N60C
C72-1C98	Stator 8 port 2-pos C72-valve	1/16 x .25mm N60C
C72H-1690	10 port 2-pos valve	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1690A	10 port 2-pos valve air actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1690D	10 port 2-pos valve replacement	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1690E	10 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1690ED	10 port 2-pos valve m-electric	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1694	4 port 2-pos valve	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1694A	4 port 2-pos valve air actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1694D	4 port 2-pos valve replacement	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1694E	4 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1694EH	4 port 2-pos valve m-electric	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1696	6 port 2-pos valve manual	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1696A	6 port 2-pos valve air actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1696D	6 port 2-pos valve replacement	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1696E	6 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1696EH	6 port 2-pos valve m-electric	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1698	8 port 2-pos valve manual	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1698A	8 port 2-pos valve air actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1698D	8 port 2-pos valve replacement	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1698E	8 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/10000psi liq N60C/E3
C72H-1698ED	8 port 2-pos valve m-electric	1/16 x .25mm 50C/10000psi liq N60C/E3
C72M-46R6	Rotor seal 6 port 2-pos C72M-valve	360um .10mm 50C/all pressures liqE3
C72M-4C94	Stator 4 port 2-pos C72M-valve	360um .10mm N60C
C72M-4C96	Stator 6 port 2-pos C72M-valve	360um .10mm N60C
C72M-66R0	Rotor seal 10 port 2-pos C72M-valve	360um .15mm 50C/all pressures liqE3
C72M-66R6	Rotor seal 6 port 2-pos C72M-valve	360um .15mm 50C/all pressures liqE3
C72M-6C90	Stator 10 port 2-pos C72M-valve	360um .15mm N60C
C72M-6C96	Stator 6 port 2-pos C72M-valve	360um .15mm N60C
C72MH-4696D	6 port 2-pos valve replacement	360um x .10mm 50C/10000psi liq N60C/E3
C72MH-4696EH	6 port 2-pos valve m-electric	360um x .10mm 50C/10000psi liq N60C/E3
C72MU-4694	4 port 2-pos valve manual	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4694A	4 port 2-pos valve air act	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4694D	4 port 2-pos valve replacement	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4694E	4 port 2-pos valve std-elect	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4694EH	4 port 2-pos valve m-elect	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4696	6 port 2-pos valve manual	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4696A	6 port 2-pos valve air act	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4696D	6 port 2-pos valve replacement	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4696E	6 port 2-pos valve std-elect	360um x .10mm 50C/20000psi liq N60C/E3
C72MU-4696EH	6 port 2-pos valve m-elect	360um x .10mm 50C/20000psi liq N60C/E3
C72MX-6690	10 port 2-pos valve manual	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6690A	10 port 2-pos valve air act.	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6690D	10 port 2-pos valve replacement	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6690E	10 port 2-pos valve std-elect.	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6690ED	10 port 2-pos valve m-elect.	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6696	6 port 2-pos valve manual	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6696A	6 port 2-pos valve air act.	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6696D	6 port 2-pos valve replacement	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6696E	6 port 2-pos valve std-elect.	360um x .15mm 50C/15000psi liq N60C/E3
C72MX-6696EH	6 port 2-pos valve m-elect.	360um x .15mm 50C/15000psi liq N60C/E3
C72N-66R0	Rotor seal 10 port 2-pos C72N-valve	1/32 x .15mm 50C/all pressures liqE3
C72N-66R6	Rotor seal 6 port 2-pos C72N-valve	1/32 x .15mm 50C/all pressures liqE3
C72N-6C90	Stator 10 port 2-pos C2N-valve	1/32 x .15mm N60C
C72N-6C96	Stator 6 port 2-pos C2N-valve	1/32 x .15mm N60C
C72NX-6690	10 port 2-pos valve manual	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6690A	10 port 2-pos valve air actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6690D	10 port 2-pos valve replacement	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6690E	10 port 2-pos valve standard electric actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6690ED	10 port 2-pos valve m-electric	1/32 x .15mm 50C/15000psi liq N60C/E3

# VICI Valco-Cheminert PRODUCT LIST 2012




C72NX-6696	6 port 2-pos valve manual	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6696A	6 port 2-pos valve air actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6696D	6 port 2-pos valve replacement	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6696E	6 port 2-pos valve standard electric actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C72NX-6696EH	6 port 2-pos valve m-electric	1/32 x .15mm 50C/15000psi liq N60C/E3
C72X-1690	10 port 2-pos valvemanual	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1690A	10 port 2-pos valve air actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1690D	10 port 2-pos valve replacement	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1690E	10 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1690ED	10 port 2-pos valve m-electric	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1694	4 port 2-pos valvemanual	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1694A	4 port 2-pos valve air actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1694D	4 port 2-pos valve replacement	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1694E	4 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1694EH	4 port 2-pos valve m-electric	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1696	6 port 2-pos valve manual	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1696A	6 port 2-pos valve air actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1696D	6 port 2-pos valve replacement	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1696E	6 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1696EH	6 port 2-pos valve m-electric	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1698	8 port 2-pos valvemanual	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1698A	8 port 2-pos valve air actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1698D	8 port 2-pos valve replacement	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1698E	8 port 2-pos valve standard electric actuator	1/16 x .25mm 50C/15000psi liq N60C/E3
C72X-1698ED	8 port 2-pos valve m-electric	1/16 x .25mm 50C/15000psi liq N60C/E3
C74-16R-.01	Rotor seal 4 port 2-pos int vol valve C74-valve	1/16 x .25mm 50C/all pressures liq .01ul E3
C74-16R-.02	Rotor seal 4 port 2-pos int vol valve C74-valve	1/16 x .25mm 50C/all pressures liq .02ul E3
C74-16R-.05	Rotor seal 4 port 2-pos int vol valve C74-valve	1/16 x .25mm 50C/all pressures liq .05ul E3
C74-1C9	Stator 4 port 2-pos int vol C74-valve	1/16 x .25mm N60C
C74H-1694-.01	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/10000psi liq .01ul N60C/E3
C74H-1694-.01A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/10000psi liq .01ul N60C/E3
C74H-1694-.01D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/10000psi liq .01ul N60C/E3
C74H-1694-.01E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/10000psi liq .01ul N60C/E3
C74H-1694-.01EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/10000psi liq .01ul N60C/E3
C74H-1694-.02	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/10000psi liq .02ul N60C/E3
C74H-1694-.02A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/10000psi liq .02ul N60C/E3
C74H-1694-.02D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/10000psi liq .02ul N60C/E3
C74H-1694-.02E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/10000psi liq .02ul N60C/E3
C74H-1694-.02EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/10000psi liq .02ul N60C/E3
C74H-1694-.05	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/10000psi liq .05ul N60C/E3
C74H-1694-.05A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/10000psi liq .05ul N60C/E3
C74H-1694-.05D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/10000psi liq .05ul N60C/E3
C74H-1694-.05E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/10000psi liq .05ul N60C/E3
C74H-1694-.05EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/10000psi liq .05ul N60C/E3
C74N-66R-.004	Rotor seal 4 port 2-pos int vol valve C74N-valve	1/32 x .15mm 50C/all pressures liq .004ul E3
C74N-66R-.01	Rotor seal 4 port 2-pos int vol valve C74N-valve	1/32 x .15mm 50C/all pressures liq .01ul E3
C74N-66R-.02	Rotor seal 4 port 2-pos int vol valve C74N-valve	1/32 x .15mm 50C/all pressures liq .02ul E3
C74N-6C9	Stator 4 port 2-pos int vol C74N-valve	1/32 x .15mm N60C
C74NX-6694-.01	4 port 2-pos int vol valve manual	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.004	4 port 2-pos int vol valve manual	1/32 x .15mm 50C/15000psi liq .004ul N60C/E3
C74NX-6694-.004A	4 port 2-pos int vol valve air actuator	1/32 x .15mm 50C/15000psi liq .004ul N60C/E3
C74NX-6694-.004D	4 port 2-pos int vol valve replacement valve	1/32 x .15mm 50C/15000psi liq .004ul N60C/E3
C74NX-6694-.004E	4 port 2-pos int vol valve std. electric actuator	1/32 x .15mm 50C/15000psi liq .004ul N60C/E3
C74NX-6694-.004EH	4 port 2-pos int vol valve microelectric actuator	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.01	4 port 2-pos int vol valve manual	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.01A	4 port 2-pos int vol valve air actuator	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.01D	4 port 2-pos int vol valve replacement valve	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.01E	4 port 2-pos int vol valve std. electric actuator	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.01EH	4 port 2-pos int vol valve microelectric actuator	1/32 x .15mm 50C/15000psi liq .01ul N60C/E3
C74NX-6694-.02	4 port 2-pos int vol valve manual	1/32 x .15mm 50C/15000psi liq .02ul N60C/E3
C74NX-6694-.02A	4 port 2-pos int vol valve air actuator	1/32 x .15mm 50C/15000psi liq .02ul N60C/E3

# VICI Valco-Cheminert PRODUCT LIST 2012




C74NX-6694-.02D	4 port 2-pos int vol valve replacement valve	1/32 x .15mm 50C/15000psi liq .02ul N60C
C74NX-6694-.02E	4 port 2-pos int vol valve std. electric actuator	1/32 x .15mm 50C/15000psi liq .02ul N60C
C74NX-6694-.02EH	4 port 2-pos int vol valve microelectric actuator	1/32 x .15mm 50C/15000psi liq .02ul N60C
C74X-1694-.01	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/15000psi liq .01ul N60C
C74X-1694-.01A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/15000psi liq .01ul N60C
C74X-1694-.01D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/15000psi liq .01ul N60C
C74X-1694-.01E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/15000psi liq .01ul N60C
C74X-1694-.01EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/15000psi liq .01ul N60C/E3
C74X-1694-.02	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/15000psi liq .02ul N60C/E3
C74X-1694-.02A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/15000psi liq .02ul N60C/E3
C74X-1694-.02D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/15000psi liq .02ul N60C/E3
C74X-1694-.02E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/15000psi liq .02ul N60C/E3
C74X-1694-.02EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/15000psi liq .02ul N60C/E3
C74X-1694-.05	4 port 2-pos int vol valve manual	1/16 x .25mm 50C/15000psi liq .05ul N60C/E3
C74X-1694-.05A	4 port 2-pos int vol valve air actuator	1/16 x .25mm 50C/15000psi liq .05ul N60C/E3
C74X-1694-.05D	4 port 2-pos int vol valve replacement	1/16 x .25mm 50C/15000psi liq .05ul N60C/E3
C74X-1694-.05E	4 port 2-pos int vol valve std-electric	1/16 x .25mm 50C/15000psi liq .05ul N60C/E3
C74X-1694-.05EH	4 port 2-pos int vol valve m-electric	1/16 x .25mm 50C/15000psi liq .05ul N60C/E3
C75-16R0	Rotor seal 11 port 10-pos C75-valve	1/16 x .25mm 50C/all pressures liqE3
C75-16R6	Rotor seal 7 port 6-pos C75-valve	1/16 x .25mm 50C/all pressures liqE3
C75-16R8	Rotor seal 9 port 8-pos C75-valve	1/16 x .25mm 50C/all pressures liqE3
C75-1C90	Stator 11 port 10-pos C75-valve	1/16 x .25mm N60C
C75-1C96	Stator 7 port 6-pos C75-valve	1/16 x .25mm N60C
C75-1C98	Stator 7 port 6-pos C75-valve	1/16 x .25mm N60C
C75H-1690A	10 mpos dead-end path air actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1690D	10 mpos dead-end path replacement	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1690E	10 mpos dead-end path standard electric actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1690EMT	10 mpos dead-end path m-electric	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1696A	6 mpos dead-end path air actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1696D	6 mpos dead-end path replacement	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1696E	6 mpos dead-end path standard electric actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1696EMH	6 mpos dead-end path m-electric	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1698A	8 mpos dead-end path air actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1698D	8 mpos dead-end path replacement	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1698E	8 mpos dead-end path standard electric actuator	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75H-1698EMT	8 mpos dead-end path m-electric	1/16 x .25mm 50C/10'000psi liq N60C/E3
C75N-66R0	Rotor seal 10 mpos C75N-valve dead-end path	1/32 x .15mm 50C/all pressures liqE3
C75N-66R6	Rotor seal 6 mpos C75N-valve dead-end path	1/32 x .15mm 50C/all pressures liqE3
C75N-66R8	Rotor seal 8 mpos C75N-valve dead-end path	1/32 x .15mm 50C/all pressures liqE3
C75N-6C90	Stator 10 mpos dead-end path C75N-valve	1/32 x .15mm N60C
C75N-6C96	Stator 6 mpos dead-end path C75N-valve	1/32 x .15mm N60C
C75N-6C98	Stator 8 mpos dead-end path C75N-valve	1/32 x .15mm N60C
C75NX-6690A	10 mpos dead-end path air actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6690D	10 mpos dead-end path replacement	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6690E	10 mpos dead-end path standard electric actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6690EMT	10 mpos dead-end path m-electric	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6696A	6 mpos dead-end path air actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6696D	6 mpos dead-end path replacement	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6696E	6 mpos dead-end path standard electric actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6696EMH	6 mpos dead-end path m-electric	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6698A	8 mpos dead-end path air actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6698D	8 mpos dead-end path replacement	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6698E	8 mpos dead-end path standard electric actuator	1/32 x .15mm 50C/15000psi liq N60C/E3
C75NX-6698EMT	8 mpos dead-end path m-electric	1/32 x .15mm 50C/15000psi liq N60C/E3
C8UW	!!! RECOMMEND 2C8UW	SUPPLIED WITH LARGER HANDLE
C8UWE	8 port 2-pos valve manual	1/16 x .75mm 225C/400psi gas N60/E
C8W	8 port 2-pos valve manual	1/16 x .4mm 75C/5000psi liq N60/H
C8WE	8 port 2-pos valve manual	1/16 x .4mm 225C/400psi gas N60/E
C9M1KF	Cheminert manifold 1/4in -28 1/16in OD tubing	8 inlets 1 outlet .75mm bore CTFE
C9M1PK	Cheminert manifold 1/4in -28 1/16in OD tubing	8 inlets 1 outlet .75mm bore PEEK
C9M2KF	Cheminert manifold 1/4in -28 1/8in OD tubing	8 inlets 1 outlet 1.5mm bore CTFE

# VICI Valco-Cheminert PRODUCT LIST 2012




C9M2PK	Cheminert manifold 1/4in -28 1/8in OD tubing	8 inlets 1 outlet 1.5mm bore PEEK
C-BOT16	Cheminert solvent bottle 500ml	with plain cap PE
C-BOT32	Cheminert solvent bottle 1000ml	with plain cap PE
CBUCKF	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.25mm bore CTFE
CBUCPK	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.25mm bore PEEK
CBUCS6	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.25mm bore 316SS
CBUFLKF	Cheminert luer adapter bulkhead union	female luer to 1/4-28 1.5mm bore CTFE
CBUFLPK	Cheminert luer adapter bulkhead union	female luer to 1/4-28 1.5mm bore PEEK
CBUKF	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.5mm bore CTFE
CBULKF	Cheminert bulkhead union 1/4in -28 1/8in OD tubing	1.5mm bore CTFE
CBULPK	Cheminert bulkhead union 1/4in -28 1/8in OD tubing	1.5mm bore PEEK
CBULS6	Cheminert bulkhead union 1/4in -28 1/8in OD tubing	1.5mm bore 316SS
CBUMKF	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.75mm bore CTFE
CBUMPK	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.75mm bore PEEK
CBUMS6	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.75mm bore 316SS
CBUPK	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.5mm bore PEEK
CBUS6	Cheminert bulkhead union 1/4in -28 1/16in OD tubing	.5mm bore 316SS
CCKF-5	Cheminert cap 1/4in -28	CTFE 5/pkg
CCPK-5	Cheminert cap 1/4in -28	PEEK 5/pkg
CD10-4841-M1A	M10 diluter/dispenser	10 nl - 10 ml
CD50-8572-M2A	M50 diluter/dispenser	50 ul - 50 ml
CEF1	Column end fitting 1/16in	.25mm bore 1mm cone without frit 316SS
CEF1.5	Column end fitting 1/16in x 1/32in	.25mm bore 1mm cone without frit 316SS
CEF1.5F	Column end fitting 1/16in x 1/32in	.25mm bore 1mm cone removable 2um frit 316SS
CEF1F	Column end fitting 1/16in	.25mm bore 1mm cone removable 2um frit 316SS
CEF2.51.0	Column end fitting 1/8in x 1/32in	.25mm bore 1mm cone without frit 316SS
CEF2.51.0F	Column end fitting 1/8in x 1/32in	.25mm bore 1mm cone removable 2um frit 316SS
CEF211.0	Column end fitting 1/8in x 1/16in	.25mm bore 1mm cone without frit 316SS
CEF211.0F	Column end fitting 1/8in x 1/16in	.25mm bore 1mm cone removable 2um frit 316SS
CEF212.0	Column end fitting 1/8in x 1/16in	.25mm bore 2mm cone without frit 316SS
CEF212.0F	Column end fitting 1/8in x 1/16in	.25mm bore 2mm cone removable 2um frit 316SS
CEF412.1	Column end fitting 1/4in x 1/16in	.4mm bore 2.1mm cone without frit 316SS
CEF412.1F	Column end fitting 1/4in x 1/16in	.4mm bore 2.1mm cone removable 2um frit 316SS
CEF413.0	Column end fitting 1/4in x 1/16in	.4mm bore 3mm cone without frit 316SS
CEF413.0F	Column end fitting 1/4in x 1/16in	.4mm bore 3mm cone removable 2um frit 316SS
CEF414.0	Column end fitting 1/4in x 1/16in	.4mm bore 4mm cone without frit 316SS
CEF414.0F	Column end fitting 1/4in x 1/16in	.4mm bore 4mm cone removable 2um frit 316SS
CEF414.6	Column end fitting 1/4in x 1/16in	.4mm bore 4.6mm cone without frit 316SS
CEF414.6F	Column end fitting 1/4in x 1/16in	.4mm bore 4.6mm cone removable 2um frit 316SS
C-EN.5FPKB	Nut External 1/32in Cheminert	PEEK black
CEN1KF	Cheminert external nut 1/4in -28 tube OD 1/16in	CTFE 5/pkg
CEN1PK	Cheminert external nut 1/4in -28 tube OD 1/16in	PEEK 5/pkg
CEN2KF	Cheminert external nut 1/4in -28 tube OD 1/8in	CTFE 5/pkg
CEN2PK	Cheminert external nut 1/4in -28 tube OD 1/8in	PEEK 5/pkg
CEPA2KF	Cheminert pipe connector 1/8in NPT	1/4in -28male x M-NPT 1.5mm bore CTFE
CEPA2PK	Cheminert pipe connector 1/8in NPT	1/4in -28male x M-NPT 1.5mm bore PEEK
CEPA4KF	Cheminert pipe connector 1/4in NPT	1/4in -28male x M-NPT 1.5mm bore CTFE
CEPA4PK	Cheminert pipe connector 1/4in NPT	1/4in -28male x M-NPT 1.5mm bore PEEK
C-F1.5TI	Frit Titanium PEEK-encapsulated	0.5 micron
CF-1A	Cheminert nut 1/4in -28 for 1/16in OD tubing	Assorted incl. washer for flange 12/pkg
CF-1BE	Cheminert nut 1/4in -28 for 1/16in OD tubing	Blue incl. washer for flange 10/pkg
CF-1BK	Cheminert nut 1/4in -28 for 1/16in OD tubing	Black incl. washer for flange 10/pkg
CF-1BR	Cheminert nut 1/4in -28 for 1/16in OD tubing	Brown incl. washer for flange 10/pkg
CF-1DG	Cheminert nut 1/4in -28 for 1/16in OD tubing	Dark gray incl. washer for flange 10/pkg
CF-1E	Cheminert nut 1/4in -28 for 1/16in OD tubing	Orange incl. washer for flange 10/pkg
CF-1G	Cheminert nut 1/4in -28 for 1/16in OD tubing	Green incl. washer for flange 10/pkg
CF-1L	Cheminert nut 1/4in -28 for 1/16in OD tubing	Lavender incl. washer for flange 10/pkg
CF-1N	Cheminert nut 1/4in -28 for 1/16in OD tubing	Natural incl. washer for flange 10/pkg
CF-1P	Cheminert nut 1/4in -28 for 1/16in OD tubing	Purple incl. washer for flange 10/pkg
CF-1R	Cheminert nut 1/4in -28 for 1/16in OD tubing	Red incl. washer for flange 10/pkg
CF-1W	Cheminert nut 1/4in -28 for 1/16in OD tubing	White incl. washer for flange 10/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




CF-1Y	Cheminert nut 1/4in -28 for 1/16in OD tubing	Yellow incl. washer for flange 10/pkg
CF-2A	Cheminert nut 1/4in -28 for 1/8in OD tubing	Assorted incl. washer for flange 12/pkg
CF-2BE	Cheminert nut 1/4in -28 for 1/8in OD tubing	Blue incl. washer for flange 10/pkg
CF-2BK	Cheminert nut 1/4in -28 for 1/8in OD tubing	Black incl. washer for flange 10/pkg
CF-2BR	Cheminert nut 1/4in -28 for 1/8in OD tubing	Brown incl. washer for flange 10/pkg
CF-2DG	Cheminert nut 1/4in -28 for 1/8in OD tubing	Dark gray incl. washer for flange 10/pkg
CF-2E	Cheminert nut 1/4in -28 for 1/8in OD tubing	Orange incl. washer for flange 10/pkg
CF-2G	Cheminert nut 1/4in -28 for 1/8in OD tubing	Green incl. washer for flange 10/pkg
CF-2L	Cheminert nut 1/4in -28 for 1/8in OD tubing	Lavender incl. washer for flange 10/pkg
CF-2N	Cheminert nut 1/4in -28 for 1/8in OD tubing	Natural incl. washer for flange 10/pkg
CF-2P	Cheminert nut 1/4in -28 for 1/8in OD tubing	Purple incl. washer for flange 10/pkg
CF-2R	Cheminert nut 1/4in -28 for 1/8in OD tubing	Red incl. washer for flange 10/pkg
CF-2W	Cheminert nut 1/4in -28 for 1/8in OD tubing	White incl. washer for flange 10/pkg
CF-2Y	Cheminert nut 1/4in -28 for 1/8in OD tubing	Yellow incl. washer for flange 10/pkg
CFE-S10	Cheminert in-line filter	1/4in -28 10um screen SS316 PTFE/CTFE
CFE-S2	Cheminert in-line filter	1/4in -28 2um screen SS316 PTFE/CTFE
CFE-S75	Cheminert in-line filter	1/4in -28 75um screen SS316 PTFE/CTFE
CFL-1A	Cheminert nut 1/4in -28 for 1/16in OD tubing	Assort collapsible ferrule CTFE 12/pkg
CFL-1A-PK	Cheminert nut 1/4in -28 for 1/16in OD tubing	Assort collapsible ferrule PEEK 12/pkg
CFL-1BE	Cheminert nut 1/4in -28 for 1/16in OD tubing	Blue collapsible ferrule 5/pkg
CFL-1BK	Cheminert nut 1/4in -28 for 1/16in OD tubing	Black collapsible ferrule 5/pkg
CFL-1BR	Cheminert nut 1/4in -28 for 1/16in OD tubing	Brown collapsible ferrule 5/pkg
CFL-1DG	Cheminert nut 1/4in -28 for 1/16in OD tubing	Dark gray collapsible ferrule 5/pkg
CFL-1E	Cheminert nut 1/4in -28 for 1/16in OD tubing	Orange collapsible ferrule 5/pkg
CFL-1G	Cheminert nut 1/4in -28 for 1/16in OD tubing	Green collapsible ferrule 5/pkg
CFL-1L	Cheminert nut 1/4in -28 for 1/16in OD tubing	Lavender collapsible ferrule 5/pkg
CFL-1N	Cheminert nut 1/4in -28 for 1/16in OD tubing	Natural collapsible ferrule 5/pkg
CFL-1P	Cheminert nut 1/4in -28 for 1/16in OD tubing	Purple collapsible ferrule 5/pkg
CFL-1PK	Cheminert nut 1/4in -28 for 1/16in OD tubing	collapsible ferrule PEEK 5/pkg
CFL-1R	Cheminert nut 1/4in -28 for 1/16in OD tubing	Red collapsible ferrule 5/pkg
CFL-1W	Cheminert nut 1/4in -28 for 1/16in OD tubing	White collapsible ferrule 5/pkg
CFL-1Y	Cheminert nut 1/4in -28 for 1/16in OD tubing	Yellow collapsible ferrule 5/pkg
CFL-2A	Cheminert nut 1/4in -28 for 1/8in OD tubing	Assort collapsible ferrule CTFE 12/pkg
CFL-2A-PK	Cheminert nut 1/4in -28 for 1/8in OD tubing	Assort collapsible ferrule PEEK 12/pkg
CFL-2BE	Cheminert nut 1/4in -28 for 1/8in OD tubing	Blue collapsible ferrule 5/pkg
CFL-2BK	Cheminert nut 1/4in -28 for 1/8in OD tubing	Black collapsible ferrule 5/pkg
CFL-2BR	Cheminert nut 1/4in -28 for 1/8in OD tubing	Brown collapsible ferrule 5/pkg
CFL-2DG	Cheminert nut 1/4in -28 for 1/8in OD tubing	Dark gray collapsible ferrule 5/pkg
CFL-2E	Cheminert nut 1/4in -28 for 1/8in OD tubing	Orange collapsible ferrule 5/pkg
CFL-2G	Cheminert nut 1/4in -28 for 1/8in OD tubing	Green collapsible ferrule 5/pkg
CFL-2L	Cheminert nut 1/4in -28 for 1/8in OD tubing	Lavender collapsible ferrule 5/pkg
CFL-2N	Cheminert nut 1/4in -28 for 1/8in OD tubing	Natural collapsible ferrule 5/pkg
CFL-2P	Cheminert nut 1/4in -28 for 1/8in OD tubing	Purple collapsible ferrule 5/pkg
CFL-2PK	Cheminert nut 1/4in -28 for 1/8in OD tubing	collapsible ferrule PEEK 5/pkg
CFL-2R	Cheminert nut 1/4in -28 for 1/8in OD tubing	Red collapsible ferrule 5/pkg
CFL-2W	Cheminert nut 1/4in -28 for 1/8in OD tubing	White collapsible ferrule 5/pkg
CFL-2Y	Cheminert nut 1/4in -28 for 1/8in OD tubing	Yellow collapsible ferrule 5/pkg
CFL-4D	Cheminert nut 1/2in -20 for 1/4in OD tubing	Delrin
CFL-4KF	Cheminert nut 1/2in -20 for 1/4in OD tubing	CTFE
CFL-4PPS	Cheminert nut 1/2in -20 for 1/4in OD tubing	PPS
CFLAKF	Cheminert female luer adapter to 1/4-28	1.5mm bore CTFE
CFLAPFA	Cheminert female luer adapter to 1/4-28	1.5mm bore PFA
CFLAPK	Cheminert female luer adapter to 1/4-28	1.5mm bore PEEK
CFL-CB1KF	Cheminert ferrule collapsible for 1/16in OD tubing	CTFE 10/pkg
CFL-CB1PK	Cheminert ferrule collapsible for 1/16in OD tubing	PEEK 10/pkg
CFL-CB2KF	Cheminert ferrule collapsible for 1/8in OD tubing	CTFE 10/pkg
CFL-CB2PK	Cheminert ferrule collapsible for 1/8in OD tubing	PEEK 10/pkg
CFL-CB4KF-S	Cheminert bushing collapsible 1/4in OD tubing	CTFE
CFPA2KF	Cheminert pipe connector 1/8in NPT	1/4in -28female x F-NPT 1.5mm bore CTFE
CFPA2PK	Cheminert pipe connector 1/8in NPT	1/4in -28female x F-NPT 1.5mm bore PEEK
CFPA4KF	Cheminert pipe connector 1/4in NPT	1/4in -28female x F-NPT 1.5mm bore CTFE

# VICI Valco-Cheminert PRODUCT LIST 2012




CFPA4PK	Cheminert pipe connector 1/4in NPT	1/4in -28female x F-NPT 1.5mm bore PEE
CFS-A01010-010B	GC-Column FS ValcoBond VB-1	10m length ID 0.10mm df 0.10um
CFS-A01010-020B	GC-Column FS ValcoBond VB-1	10m length ID 0.10mm df 0.20um
CFS-A01018-010B	GC-Column FS ValcoBond VB-1	10m length ID 0.18mm df 0.10um
CFS-A01018-018B	GC-Column FS ValcoBond VB-1	10m length ID 0.18mm df 0.18um
CFS-A01018-040B	GC-Column FS ValcoBond VB-1	10m length ID 0.18mm df 0.40um
CFS-A01525-010B	GC-Column FS ValcoBond VB-1	15m length ID 0.25mm df 0.10um
CFS-A01525-025B	GC-Column FS ValcoBond VB-1	15m length ID 0.25mm df 0.25um
CFS-A01525-050B	GC-Column FS ValcoBond VB-1	15m length ID 0.25mm df 0.50um
CFS-A01525-100B	GC-Column FS ValcoBond VB-1	15m length ID 0.25mm df 1.00um
CFS-A01525-150B	GC-Column FS ValcoBond VB-1	15m length ID 0.25mm df 1.50um
CFS-A01532-010B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 0.10um
CFS-A01532-025B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 0.25um
CFS-A01532-032B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 0.32um
CFS-A01532-050B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 0.50um
CFS-A01532-100B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 1.00um
CFS-A01532-200B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 2.00um
CFS-A01532-300B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 3.00um
CFS-A01532-400B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 4.00um
CFS-A01532-500B	GC-Column FS ValcoBond VB-1	15m length ID 0.32mm df 5.00um
CFS-A01553-015B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 0.15um
CFS-A01553-050B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 0.50um
CFS-A01553-100B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 1.00um
CFS-A01553-150B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 1.50um
CFS-A01553-300B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 3.00um
CFS-A01553-500B	GC-Column FS ValcoBond VB-1	15m length ID 0.53mm df 5.00um
CFS-A02010-010B	GC-Column FS ValcoBond VB-1	20m length ID 0.10mm df 0.10um
CFS-A02010-020B	GC-Column FS ValcoBond VB-1	20m length ID 0.10mm df 0.20um
CFS-A02010-040B	GC-Column FS ValcoBond VB-1	20m length ID 0.10mm df 0.40um
CFS-A03018-010B	GC-Column FS ValcoBond VB-1	30m length ID 0.18mm df 0.10um
CFS-A03018-018B	GC-Column FS ValcoBond VB-1	30m length ID 0.18mm df 0.18um
CFS-A03018-040B	GC-Column FS ValcoBond VB-1	30m length ID 0.18mm df 0.40um
CFS-A03018-100B	GC-Column FS ValcoBond VB-1	30m length ID 0.18mm df 1.00um
CFS-A03025-010B	GC-Column FS ValcoBond VB-1	30m length ID 0.25mm df 0.10um
CFS-A03025-025B	GC-Column FS ValcoBond VB-1	30m length ID 0.25mm df 0.25um
CFS-A03025-050B	GC-Column FS ValcoBond VB-1	30m length ID 0.25mm df 0.50um
CFS-A03025-100B	GC-Column FS ValcoBond VB-1	30m length ID 0.25mm df 1.00um
CFS-A03025-150B	GC-Column FS ValcoBond VB-1	30m length ID 0.25mm df 1.50um
CFS-A03032-010B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 0.10um
CFS-A03032-025B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 0.25um
CFS-A03032-032B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 0.32um
CFS-A03032-050B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 0.50um
CFS-A03032-100B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 1.00um
CFS-A03032-200B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 2.00um
CFS-A03032-300B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 3.00um
CFS-A03032-400B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 4.00um
CFS-A03032-500B	GC-Column FS ValcoBond VB-1	30m length ID 0.32mm df 5.00um
CFS-A03053-015B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 0.15um
CFS-A03053-050B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 0.50um
CFS-A03053-100B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 1.00um
CFS-A03053-150B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 1.50um
CFS-A03053-300B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 3.00um
CFS-A03053-500B	GC-Column FS ValcoBond VB-1	30m length ID 0.53mm df 5.00um
CFS-A04018-018B	GC-Column FS ValcoBond VB-1	40m length ID 0.18mm df 0.18um
CFS-A04018-040B	GC-Column FS ValcoBond VB-1	40m length ID 0.18mm df 0.40um
CFS-A06025-010B	GC-Column FS ValcoBond VB-1	60m length ID 0.25mm df 0.10um
CFS-A06025-025B	GC-Column FS ValcoBond VB-1	60m length ID 0.25mm df 0.25um
CFS-A06025-050B	GC-Column FS ValcoBond VB-1	60m length ID 0.25mm df 0.50um
CFS-A06025-100B	GC-Column FS ValcoBond VB-1	60m length ID 0.25mm df 1.00um
CFS-A06025-150B	GC-Column FS ValcoBond VB-1	60m length ID 0.25mm df 1.50um
CFS-A06032-010B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 0.10um

# VICI Valco-Cheminert PRODUCT LIST 2012




CFS-A06032-032B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 0.32um
CFS-A06032-050B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 0.50um
CFS-A06032-100B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 1.00um
CFS-A06032-200B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 2.00um
CFS-A06032-300B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 3.00um
CFS-A06032-400B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 4.00um
CFS-A06032-500B	GC-Column FS ValcoBond VB-1	60m length ID 0.32mm df 5.00um
CFS-A06053-100B	GC-Column FS ValcoBond VB-1	60m length ID 0.53mm df 1.00um
CFS-A06053-150B	GC-Column FS ValcoBond VB-1	60m length ID 0.53mm df 1.50um
CFS-A06053-300B	GC-Column FS ValcoBond VB-1	60m length ID 0.53mm df 3.00um
CFS-A06053-500B	GC-Column FS ValcoBond VB-1	60m length ID 0.53mm df 5.00um
CFS-B01010-010B	GC-Column FS ValcoBond VB-5	10m length ID 0.10mm df 0.10um
CFS-B01010-020B	GC-Column FS ValcoBond VB-5	10m length ID 0.10mm df 0.20um
CFS-B01018-018B	GC-Column FS ValcoBond VB-5	10m length ID 0.18mm df 0.18um
CFS-B01018-040B	GC-Column FS ValcoBond VB-5	10m length ID 0.18mm df 0.40um
CFS-B01525-010B	GC-Column FS ValcoBond VB-5	15m length ID 0.25mm df 0.10um
CFS-B01525-025B	GC-Column FS ValcoBond VB-5	15m length ID 0.25mm df 0.25um
CFS-B01525-050B	GC-Column FS ValcoBond VB-5	15m length ID 0.25mm df 0.50um
CFS-B01525-100B	GC-Column FS ValcoBond VB-5	15m length ID 0.25mm df 1.00um
CFS-B01532-010B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 0.10um
CFS-B01532-025B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 0.25um
CFS-B01532-050B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 0.50um
CFS-B01532-100B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 1.00um
CFS-B01532-200B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 2.00um
CFS-B01532-300B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 3.00um
CFS-B01532-500B	GC-Column FS ValcoBond VB-5	15m length ID 0.32mm df 5.00um
CFS-B01553-050B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 0.50um
CFS-B01553-100B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 1.00um
CFS-B01553-150B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 1.50um
CFS-B01553-200B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 2.00um
CFS-B01553-265B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 2.65um
CFS-B01553-300B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 3.00um
CFS-B01553-500B	GC-Column FS ValcoBond VB-5	15m length ID 0.53mm df 5.00um
CFS-B02010-010B	GC-Column FS ValcoBond VB-5	20m length ID 0.10mm df 0.10um
CFS-B02010-020B	GC-Column FS ValcoBond VB-5	20m length ID 0.10mm df 0.20um
CFS-B02018-018B	GC-Column FS ValcoBond VB-5	20m length ID 0.18mm df 0.18um
CFS-B02018-040B	GC-Column FS ValcoBond VB-5	20m length ID 0.18mm df 0.40um
CFS-B03018-018B	GC-Column FS ValcoBond VB-5	30m length ID 0.18mm df 0.18um
CFS-B03018-040B	GC-Column FS ValcoBond VB-5	30m length ID 0.18mm df 0.40um
CFS-B03025-010B	GC-Column FS ValcoBond VB-5	30m length ID 0.25mm df 0.25um
CFS-B03025-025B	GC-Column FS ValcoBond VB-5	30m length ID 0.25mm df 0.25um
CFS-B03025-050B	GC-Column FS ValcoBond VB-5	30m length ID 0.25mm df 0.50um
CFS-B03025-100B	GC-Column FS ValcoBond VB-5	30m length ID 0.25mm df 1.00um
CFS-B03032-010B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 0.10um
CFS-B03032-025B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 0.25um
CFS-B03032-050B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 0.50um
CFS-B03032-100B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 1um
CFS-B03032-200B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 2.00um
CFS-B03032-300B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 3.00um
CFS-B03032-500B	GC-Column FS ValcoBond VB-5	30m length ID 0.32mm df 5.00um
CFS-B03053-050B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df .50um
CFS-B03053-100B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df 1.00um
CFS-B03053-150B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df 1.50um
CFS-B03053-265B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df 2.65um
CFS-B03053-300B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df 3.00um
CFS-B03053-500B	GC-Column FS ValcoBond VB-5	30m length ID 0.53mm df 5.00um
CFS-B04018-018B	GC-Column FS ValcoBond VB-5	40m length ID 0.18mm df 0.18um
CFS-B04018-040B	GC-Column FS ValcoBond VB-5	40m length ID 0.18mm df 0.40um
CFS-B06025-010B	GC-Column FS ValcoBond VB-5	60m length ID 0.25mm df 0.10um
CFS-B06025-025B	GC-Column FS ValcoBond VB-5	60m length ID 0.25mm df 0.25um
CFS-B06025-050B	GC-Column FS ValcoBond VB-5	60m length ID 0.25mm df 0.50um

# VICI Valco-Cheminert PRODUCT LIST 2012




CFS-B06025-100B	GC-Column FS ValcoBond VB-5	60m length ID 0.25mm df 1.00um
CFS-B06032-010B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 0.10um
CFS-B06032-025B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 0.25um
CFS-B06032-050B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 0.50um
CFS-B06032-100B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 1.00um
CFS-B06032-200B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 2.00um
CFS-B06032-300B	GC-Column FS ValcoBond VB-5	60m length ID 0.32mm df 3.00um
CFS-B06032-500B	GC-Column FS ValcoBond VB-7	60m x 0.32mm ID x 5.00df
CFS-B06053-100B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 1.00um
CFS-B06053-150B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 1.00um
CFS-B06053-200B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 2.00um
CFS-B06053-265B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 2.65um
CFS-B06053-300B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 3.00um
CFS-B06053-500B	GC-Column FS ValcoBond VB-5	60m length ID 0.53mm df 5.00um
CFS-C01525-025B	GC-Column FS ValcoBond VB-35	15m length ID 0.25mm df 0.25um
CFS-C01525-050B	GC-Column FS ValcoBond VB-35	15m length ID 0.25mm df 0.50um
CFS-C01532-025B	GC-Column FS ValcoBond VB-35	15m length ID 0.32mm df 0.25um
CFS-C01532-050B	GC-Column FS ValcoBond VB-35	15m length ID 0.32mm df 0.50um
CFS-C01553-050B	GC-Column FS ValcoBond VB-35	15m length ID 0.53mm df 0.50um
CFS-C01553-100B	GC-Column FS ValcoBond VB-35	15m length ID 0.53mm df 1.00um
CFS-C03025-025B	GC-Column FS ValcoBond VB-35	30m length ID 0.25mm df 0.25um
CFS-C03025-050B	GC-Column FS ValcoBond VB-35	30m length ID 0.25mm df 0.50um
CFS-C03032-025B	GC-Column FS ValcoBond VB-35	30m length ID 0.32mm df 0.25um
CFS-C03032-050B	GC-Column FS ValcoBond VB-35	30m length ID 0.32mm df 0.50um
CFS-C03053-050B	GC-Column FS ValcoBond VB-35	30m length ID 0.53mm df 0.50um
CFS-C03053-100B	GC-Column FS ValcoBond VB-35	30m length ID 0.53mm df 1.00um
CFS-C06025-025B	GC-Column FS ValcoBond VB-35	60m length ID 0.25mm df 0.25um
CFS-C06025-050B	GC-Column FS ValcoBond VB-35	60m length ID 0.25mm df 0.50um
CFS-C06032-025B	GC-Column FS ValcoBond VB-35	60m length ID 0.32mm df 0.25um
CFS-C06032-050B	GC-Column FS ValcoBond VB-35	60m length ID 0.32mm df 0.50um
CFS-C06053-050B	GC-Column FS ValcoBond VB-35	60m length ID 0.53mm df 0.50um
CFS-C06053-100B	GC-Column FS ValcoBond VB-35	60m length ID 0.53mm df 1.00um
CFS-D01525-025B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.25mm df 0.25um
CFS-D01525-050B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.25mm df 0.50um
CFS-D01532-025B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.32mm df 0.25um
CFS-D01532-050B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.32mm df 0.50um
CFS-D01532-100B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.32mm df 1.0um
CFS-D01553-050B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.53mm df 0.50um
CFS-D01553-083B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.53mm df 0.83um
CFS-D01553-100B	GC-Column FS ValcoBond VB-50/608	15m length ID 0.53mm df 1.00um
CFS-D03025-025B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.25mm df 0.25um
CFS-D03025-050B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.25mm df 0.50um
CFS-D03032-025B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.32mm df 0.25um
CFS-D03032-050B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.32mm df 0.50um
CFS-D03032-100B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.32mm df 1.00um
CFS-D03053-050B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.53mm df 0.50um
CFS-D03053-083B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.53mm df 0.83um
CFS-D03053-100B	GC-Column FS ValcoBond VB-50/608	30m length ID 0.53mm df 1.00um
CFS-D06025-025B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.25mm df 0.25um
CFS-D06025-050B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.25mm df 0.50um
CFS-D06032-025B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.32mm df 0.25um
CFS-D06032-050B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.32mm df 0.50um
CFS-D06032-100B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.32mm df 1.00um
CFS-D06053-050B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.53mm df 0.50um
CFS-D06053-100B	GC-Column FS ValcoBond VB-50/608	60m length ID 0.53mm df 1.00um
CFS-E01018-100A	GC-Column FS ValcoBond VB-624	10m length ID 0.18mm df 1.00um
CFS-E01525-140A	GC-Column FS ValcoBond VB-624	15m length ID 0.25mm df 1.40um
CFS-E01532-180A	GC-Column FS ValcoBond VB-624	15m length ID 0.32mm df 1.80um
CFS-E01553-300A	GC-Column FS ValcoBond VB-624	15m length ID 0.53mm df 3.00um
CFS-E02018-100A	GC-Column FS ValcoBond VB-624	20m length ID 0.18mm df 1.00um
CFS-E02520-112A	GC-Column FS ValcoBond VB-624	25m length ID 0.20mm df 1.12um

# VICI Valco-Cheminert PRODUCT LIST 2012




CFS-E03018-100A	GC-Column FS ValcoBond VB-624	30m length ID 0.18mm df 1.00um
CFS-E03025-140A	GC-Column FS ValcoBond VB-624	30m length ID 0.25mm df 1.40um
CFS-E03032-180A	GC-Column FS ValcoBond VB-624	30m length ID 0.32mm df 1.80um
CFS-E03053-300A	GC-Column FS ValcoBond VB-624	30m length ID 0.53mm df 3.00um
CFS-E04018-100A	GC-Column FS ValcoBond VB-624	40m length ID 0.18mm df 1.00um
CFS-E06025-140A	GC-Column FS ValcoBond VB-624	60m length ID 0.25mm df 1.40um
CFS-E06032-180A	GC-Column FS ValcoBond VB-624	60m length ID 0.32mm df 1.80um
CFS-E06053-300A	GC-Column FS ValcoBond VB-624	60m length ID 0.53mm df 3.00um
CFS-E07553-300A	GC-Column FS ValcoBond VB-624	75m length ID 0.53mm df 3.00um
CFS-F01525-025A	GC-Column FS ValcoBond VB-1701	15m length ID 0.25mm df 0.25um
CFS-F01525-050A	GC-Column FS ValcoBond VB-1701	15m length ID 0.25mm df 0.50um
CFS-F01532-025A	GC-Column FS ValcoBond VB-1701	15m length ID 0.32mm df 0.25um
CFS-F01532-050A	GC-Column FS ValcoBond VB-1701	15m length ID 0.32mm df 0.50um
CFS-F01532-100A	GC-Column FS ValcoBond VB-1701	15m length ID 0.32mm df 1.00um
CFS-F01553-050A	GC-Column FS ValcoBond VB-1701	15m length ID 0.53mm df 0.50um
CFS-F01553-100A	GC-Column FS ValcoBond VB-1701	15m length ID 0.53mm df 1.00um
CFS-F03025-025A	GC-Column FS ValcoBond VB-1701	30m length ID 0.25mm df 0.25um
CFS-F03025-050A	GC-Column FS ValcoBond VB-1701	30m length ID 0.25mm df 0.50um
CFS-F03032-025A	GC-Column FS ValcoBond VB-1701	30m length ID 0.32mm df 0.25um
CFS-F03032-050A	GC-Column FS ValcoBond VB-1701	30m length ID 0.32mm df 0.50um
CFS-F03032-100A	GC-Column FS ValcoBond VB-1701	30m length ID 0.32mm df 1.00um
CFS-F03053-050A	GC-Column FS ValcoBond VB-1701	30m length ID 0.53mm df 0.50um
CFS-F03053-100A	GC-Column FS ValcoBond VB-1701	30m length ID 0.53mm df 1.00um
CFS-F06025-025A	GC-Column FS ValcoBond VB-1701	60m length ID 0.25mm df 0.25um
CFS-F06025-050A	GC-Column FS ValcoBond VB-1701	60m length ID 0.25mm df 0.50um
CFS-F06032-025A	GC-Column FS ValcoBond VB-1701	60m length ID 0.32mm df 0.25um
CFS-F06032-050A	GC-Column FS ValcoBond VB-1701	60m length ID 0.32mm df 0.50um
CFS-F06032-100A	GC-Column FS ValcoBond VB-1701	60 length ID 0.32mm df 1.00um
CFS-F06053-050A	GC-Column FS ValcoBond VB-1701	60m length ID 0.53mm df 0.50um
CFS-F06053-100A	GC-Column FS ValcoBond VB-1701	60 length ID 0.53mm df 1.00um
CFS-G01010-010A	GC-Column FS ValcoBond VB-Wax	10m length ID 0.10mm df 0.10um
CFS-G01018-018A	GC-Column FS ValcoBond VB-WAX	10m length ID 0.18mm df 0.18um
CFS-G01525-025A	GC-Column FS ValcoBond VB-Wax	VB-Wax column 15m x 0.25mm ID x 0.25df
CFS-G01532-025A	GC-Column FS ValcoBond VB-WAX	15m length ID 0.32mm df 0.25um
CFS-G01532-050A	GC-Column FS ValcoBond VB-WAX	15m length ID 0.32mm df 0.50um
CFS-G01532-100A	GC-Column FS ValcoBond VB-WAX	15m length ID 0.32mm df 1.00um
CFS-G01553-050A	GC-Column FS ValcoBond VB-WAX	15m length ID 0.53mm df 0.50um
CFS-G01553-100A	GC-Column FS ValcoBond VB-Wax	VB-Wax column 15m x 0.53mm ID x 1.00df
CFS-G02010-010A	GC-Column FS ValcoBond VB-WAX	20m length ID 0.10mm df 0.10um
CFS-G02018-018A	GC-Column FS ValcoBond VB-WAX	20m length ID 0.18mm df 0.18um
CFS-G03025-025A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.25mm df 0.25um
CFS-G03032-025A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.32mm df 0.25um
CFS-G03032-050A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.32mm df 0.50um
CFS-G03032-100A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.32mm df 1.00um
CFS-G03053-050A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.53mm df 0.50um
CFS-G03053-100A	GC-Column FS ValcoBond VB-Wax	30m length ID 0.53mm df 1.00um
CFS-G06025-025A	GC-Column FS ValcoBond VB-Wax	60m length ID 0.25mm df 0.25um
CFS-G06032-025A	GC-Column FS ValcoBond VB-Wax	60m length ID 0.32mm df 0.25um
CFS-G06032-050A	GC-Column FS ValcoBond VB-Wax	60m length ID 0.32mm df 0.50um
CFS-G06053-100A	GC-Column FS ValcoBond VB-Wax	60m length ID 0.53mm df 1um
CFSL100FEP	Cheminert sample loop 100ul	flangeless connections 1/4in -28 FEP
CFSL100PK	Cheminert sample loop 100ul	flangeless connections 1/4in -28 PEEK
CFSL100TF	Cheminert sample loop 100ul	flangeless connections 1/4in -28 PTFE
CFSL1KFEP	Cheminert sample loop 1ml	flangeless connections 1/4in -28 FEP
CFSL1KPK	Cheminert sample loop 1ml	flangeless connections 1/4in -28 PEEK
CFSL1KTF	Cheminert sample loop 1ml	flangeless connections 1/4in -28 PTFE
CFSL20FEP	Cheminert sample loop 20ul	flangeless connections 1/4in -28 FEP
CFSL20PK	Cheminert sample loop 20ul	flangeless connections 1/4in -28 PEEK
CFSL20TF	Cheminert sample loop 20ul	flangeless connections 1/4in -28 PTFE
CFSL250FEP	Cheminert sample loop 250ul	flangeless connections 1/4in -28 FEP
CFSL250PK	Cheminert sample loop 250ul	flangeless connections 1/4in -28 PEEK

# VICI Valco-Cheminert PRODUCT LIST 2012




CFSL250TF	Cheminert sample loop 250ul	flangeless connections 1/4in -28 PTFE
CFSL2KFEP	Cheminert sample loop 2ml	flangeless connections 1/4in -28 FEP
CFSL2KPK	Cheminert sample loop 2ml	flangeless connections 1/4in -28 PEEK
CFSL2KTF	Cheminert sample loop 2ml	flangeless connections 1/4in -28 PTFE
CFSL500FEP	Cheminert sample loop 500ul	flangeless connections 1/4in -28 FEP
CFSL500PK	Cheminert sample loop 500ul	flangeless connections 1/4in -28 PEEK
CFSL500TF	Cheminert sample loop 500ul	flangeless connections 1/4in -28 PTFE
CFSL50FEP	Cheminert sample loop 50ul	flangeless connections 1/4in -28 FEP
CFSL50PK	Cheminert sample loop 50ul	flangeless connections 1/4in -28 PEEK
CFSL50TF	Cheminert sample loop 50ul	flangeless connections 1/4in -28 PTFE
CFS-N03025-025	GC-Column FS ValcoBond VB-Fluoro	30m length ID 0.25mm df 0.25um
CFS-N03053-100	GC-Column FS ValcoBond VB-Fluoro	30m length ID 0.53mm ID df 1.00um
CFS-PA1532-100	GC-Column FS ValcoPLOT HayeSep A	15m length ID 0.32mm df 10um
CFS-PA1553-200	GC-Column FS ValcoPLOT HayeSep A	15m length ID 0.53mm df 20um
CFS-PA3032-100	GC-Column FS ValcoPLOT HayeSep A	30m length ID 0.32mm df 10um
CFS-PA3053-200	GC-Column FS ValcoPLOT HayeSep A	30m length ID 0.53mm df 20um
CFS-PB1532-100	VPHS-B Column FS ValcoBond HayeSep B	15m length ID 0.32mm df 0.10um
CFS-PB1553-200	VPHS-B Column FS ValcoBond HayeSep B	15m length ID 0.53mm df 0.20um
CFS-PB3032-100	VPHS-B Column FS ValcoBond HayeSep B	30m length ID 0.32mm df 0.10um
CFS-PB3053-200	GC-Column FS ValcoPLOT HayeSep B	30m length ID 0.53mm df 20um
CFS-PC1532-100	VPHS-C Column FS ValcoBond HayeSep C	15m length ID 0.32mm df 0.10um
CFS-PC1553-200	GC-Column FS ValcoPLOT HayeSep C	15m length ID 0.53mm df 20um
CFS-PC3032-100	VPHS-C Column FS ValcoBond HayeSep C	30m length ID 0.32mm df 0.10um
CFS-PC3053-200	GC-Column FS ValcoPLOT HayeSep C	30m length ID 0.53mm df 20um
CFS-PD1532-100	GC-Column FS ValcoPLOT HayeSep D	15m length ID 0.32mm df 10um
CFS-PD1553-200	GC-Column FS ValcoPLOT HayeSep D	15m length ID 0.53mm df 20um
CFS-PD3032-100	GC-Column FS ValcoPLOT HayeSep D	30m length ID 0.32mm df 10um
CFS-PD3053-200	GC-Column FS ValcoPLOT HayeSep D	30m length ID 0.53mm df 20um
CFS-PN1532-100	VPHS-N Column FS ValcoBond HayeSep N	15m length ID 0.32mm df 0.10um
CFS-PN1553-200	VPHS-N Column FS ValcoBond HayeSep N	15m length ID 0.53mm df 0.20um
CFS-PN3032-100	VPHS-N Column FS ValcoBond HayeSep N	30m length ID 0.32mm df 0.10um
CFS-PN3053-200	GC-Column FS ValcoPLOT HayeSep N	30m length ID 0.53mm df 20um
CFS-PP1532-100	VPHS-P Column FS ValcoBond HayeSep P	15m length ID 0.32mm df 0.10um
CFS-PP1553-200	VPHS-P Column FS ValcoBond HayeSep P	15m length ID 0.53mm df 0.20um
CFS-PP3032-100	VPHS-P Column FS ValcoBond HayeSep P	30m length ID 0.32mm df 0.10um
CFS-PP3053-200	GC-Column FS ValcoPLOT HayeSep P	30m length ID 0.53mm df 20um
CFS-PQ1532-100	GC-Column FS ValcoPLOT HayeSep Q	15m length ID 0.32mm df 10um
CFS-PQ1553-200	GC-Column FS ValcoPLOT HayeSep Q	15m length ID 0.53mm df 20um
CFS-PQ3032-100	VPHS-Q Column FS ValcoBond HayeSep Q	30m length ID 0.32mm df 0.10um
CFS-PQ3053-200	GC-Column FS ValcoPLOT HayeSep Q	30m length ID 0.53mm df 20um
CFS-PR1532-100	VPHS-R Column FS ValcoBond HayeSep R	15m length ID 0.32mm df 0.10um
CFS-PR1553-200	VPHS-R Column FS ValcoBond HayeSep R	15m length ID 0.53mm df 0.20um
CFS-PR3032-100	VPHS-R Column FS ValcoBond HayeSep R	30m length ID 0.32mm df 0.10um
CFS-PR3053-200	GC-Column FS ValcoPLOT HayeSep R	30m length ID 0.53mm df 20um
CFS-PS1532-100	VPHS-S Column FS ValcoBond HayeSep S	15m length ID 0.32mm df 0.10um
CFS-PS1553-200	VPHS-S Column FS ValcoBond HayeSep S	15m length ID 0.53mm df 0.20um
CFS-PS3032-100	GC-Column FS ValcoPLOT HayeSep S	30m length ID 0.32mm df 10um
CFS-PS3053-200	VPHS-S Column FS ValcoBond HayeSep S	30m length ID 0.53mm df 0.20um
CFS-X1553-200	VP-Molesieve ColumnFS ValcoBond 5A Fused Silica	15m length ID 0.53mm df 0.20um
CFS-X1553-500	VP-Molesieve ColumnFS ValcoBond 5A Fused Silica	15m length ID 0.53mm df 0.50um
CFS-X3053-200	GC-Column FS ValcoPLOT Molesieve 5A	30m length ID 0.53mm df 20um
CFS-X3053-500	GC-Column FS ValcoPLOT Molesieve 5A	30m length ID 0.53mm df 50um
CFS-Y3053-100	GC-Column FS ValcoPlot Alumina	30m length ID 0.53mm df 10um
CFS-Y5053-100	GC-Column FS ValcoPlot Alumina	50m length ID 0.53mm df 10um
CFS-Z1553-100	VP-Alumina ColumnFS ValcoBond VP-Alumina Na2SO4	15m length ID 0.53mm df 0.10um
CFS-Z3053-100	GC-Column FS ValcoPlot Alumina	30m length ID 0.53mm df 10um
CFS-Z5053-100	GC-Column FS ValcoPlot Alumina	50m length ID 0.53mm df 10um
CFT-110	Cheminert flanging tool	110Vac non CE !!!
CFT1K-110	Cheminert starter kit for 1/16in tubing	110Vac non CE !!!
CFT1K-220	Cheminert starter kit for 1/16in tubing	220Vac non CE !!!
CFT-220	Cheminert flanging tool	220Vac non CE !!!

# VICI Valco-Cheminert PRODUCT LIST 2012




CFT2K-110	Cheminert starter kit for 1/8in tubing	110Vac non CE !!!
CFT2K-220	Cheminert starter kit for 1/8in tubing	220Vac non CE !!!
CFT-H	Cheminert holder flanging tool	
CFT-R	Cheminert razor blade 10/pkg	
CFT-TC	Cheminert tip flanging tool	for tubing ID <075 mm
CFT-TL	Cheminert tip flanging tool	for tubing ID <150 mm
CFT-TM	Cheminert tip flanging tool	for tubing ID <100 mm
CFT-TXC	Cheminert tip flanging tool	for tubing ID <0.25 mm
CFT-TXL	Cheminert tip flanging tool	for tubing ID <200 mm
CF-W1	Cheminert washers for 1/16in OD tubing	10/pkg
CF-W2	Cheminert washers for 1/8in OD tubing	10/pkg
CGC41	Cheminert glass connector	1mm bore Acetal
CGC41KF	Cheminert glass connector	1mm bore CTFE
CGC42	Cheminert glass connector	2mm bore Acetal
CGC42KF	Cheminert glass connector	2mm bore CTFE
CGCB	Cheminert glass bushing	Acetal
CGCBKF	Cheminert glass bushing	CTFE
CGCG41	Cheminert glass rod 3-1/4in long	1mm bore Acetal
CGCG42	Cheminert glass rod 3-1/4in long	2mm bore Acetal
CGCN	Cheminert nut glass connector	Acetal
CGCNKF	Cheminert nut glass connector	CTFE
CGCR	Cheminert split ring	Acetal 5/pkg
CGCRKF	Cheminert split ring	CTFE 5/pkg
CI4UW.2	!!! RECOMMEND 2CI4UW.2	SUPPLIED WITH LARGER HANDLE
CI4UW.5	!!! RECOMMEND 2CI4UW.5	SUPPLIED WITH LARGER HANDLE
CI4UW1	!!! RECOMMEND 2CI4UW1	SUPPLIED WITH LARGER HANDLE
CI4UW2	!!! RECOMMEND 2CI4UW2	SUPPLIED WITH LARGER HANDLE
CI4UWE.2	4 port 2-pos int vol valve manual	1/16 x .75mm 175C/1000psi liq .2ul N60/E
CI4UWE.5	4 port 2-pos int vol valve manual	1/16 x .75mm 175C/1000psi liq .5ul N60/E
CI4UWE1	4 port 2-pos int vol valve manual	1/16 x .75mm 175C/1000psi liq 1ul N60/E
CI4UWE2	4 port 2-pos int vol valve manual	1/16 x .75mm 175C/1000psi liq 2ul N60/E
CI4W.06	4 port 2-pos int/vol valve manual	1/16 x .4mm 75C/5000psi liq .06ul N60/H
CI4W.1	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/5000psi liq .1ul N60/H
CI4W.2	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/5000psi liq .2ul N60/H
CI4W.5	4 port 2-pos int vol valve manual	1/16 x .4mm 75C/5000psi liq .5ul N60/H
CI4WE.06	4 port 2-pos int vol valve manual	1/16 x .4mm 175C/1000psi liq .06ulN60/E
CI4WE.1	4 port 2-pos int vol valve manual	1/16 x .4mm 175C/1000psi liq .1ul N60/E
CI4WE.2	4 port 2-pos int vol valve manual	1/16 x .4mm 175C/1000psi liq .2ulN60/E
CI4WE.5	4 port 2-pos int vol valve manual	1/16 x .4mm 175C/1000psi liq .5ul N60/E
C-LFP	Loop fill port assy.	
C-LZP1PK	Plug 1/16in long hex high pressure	Cheminert Polymeric vlv series C1-C5 PEEK
CM	Column mandrel	
CM1XKF	Cheminert mixing tee 1/4in -28	for 1/16in tubing 0.75mm bore CTFE
CM1XPK	Cheminert mixing tee 1/4in -28	for 1/16in tubing 0.75mm bore PEEK
CM2XKF	Cheminert mixing tee 1/4in -28	for 1/8in tubing 1.5mm bore CTFE
CM2XPK	Cheminert mixing tee 1/4in -28	for 1/8in tubing 1.5mm bore PEEK
CMH	Close mount hardware 2-pos valve	1 or 2 holes air or elec. act.
CMH11H	Close mount hardware	Cheminert h/p air or elec. act.
CMH11L	Close mount hardware	Cheminert l/p air or elec. act.
CMH12H	Close mount hardware 2-pos. valve	m-electric actuator
CMH12L	Close mount hardware 2-pos. valve	Cheminert l/p m-electric actuator
CMH13	OBSOLETE USE CMH13H	
CMH13H	Close mount hardware m-pos. valve	Cheminert h/p m-electric actuator
CMH13L	Close mount hardware m-pos. valve	Cheminert l/p m-electric actuator
CMHMP	Close mount hardware 2-pos/m-pos valve	no mounting holes air or elec. actuator
CMLAKF	Cheminert male luer adapter to 1/4-28	1.5mm bore CTFE
CMLAPFA	Cheminert male luer adapter to 1/4-28	1.5mm bore PFA
CMLAPK	Cheminert male luer adapter to 1/4-28	1.5mm bore PEEK
C-MPFTI10	Cheminert mobile phase filter	for 1/8in OD tubing 10um Titanium frit PTFE
C-MPFTI2	Cheminert mobile phase filter	for 1/8in OD tubing 2um Titanium frit PTFE
C-MPR2	Cheminert solvent reservoir 500ml	with 2-hole cap 1/4in -28 fittings PE

# VICI Valco-Cheminert PRODUCT LIST 2012




C-MPR3	Cheminert solvent reservoir 500ml	with 3-hole cap 1/4in -28 fittings PE
C-MZP1PK	Plug 1/16in hex high pressure	Cheminert Polymeric vlv series C1-C5 PEEK
CN2-4340	10 port 2-pos valve manual	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-4340D	10 port 2-pos valve replacement	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-4340EH	10 port 2-pos valve m-electric	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-4346	6 port 2-pos valve manual	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-4346D	6 port 2-pos valve replacement	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-4346EH	6 port 2-pos valve m-elect	1/32 x .10mm 50C/5000psi liq PAEK/E
CN2-43R0	Rotor seal 10 port CN2-valve	1/32 x .10mm 50C/5000psi liq E
CN2-43R6	Rotor seal 6 port CN2-valve	1/32 x .10mm 50C/5000psi liq E
CN2-4C40I	Stator/Insert 10 port CN2-valve	0.10mm bore (.004in) .080in bolt circle E3
CN2-4C46I	Stator/Insert 6 port CN2-valve	0.10mm bore (.004in) .080in bolt circle E3
CN4-4344-.004	4 port 2-pos int vol valve manual	1/32 x .10mm 50C/5000psi liq .004ul PAEK/E
CN4-4344-.004D	4 port 2-pos int vol valve replacement	1/32 x .10mm 50C/5000psi liq .004ul PAEK/E
CN4-4344-.004EH	4 port 2-pos int vol valve m-elect	1/32 x .10mm 50C/5000psi liq .004ul PAEK/E
CN4-4344-.01	4 port 2-pos int vol valve manual	1/32 x .10mm 50C/5000psi liq .01ul PAEK/E
CN4-4344-.01D	4 port 2-pos int vol valve replacement	1/32 x .10mm 50C/5000psi liq .01ul PAEK/E
CN4-4344-.01EH	4 port 2-pos int vol valve m-elect	1/32 x .10mm 50C/5000psi liq .01ul PAEK/E
CN4-4344-.02	4 port 2-pos int vol valve manual	1/32 x .10mm 50C/5000psi liq .02ul PAEK/E
CN4-4344-.02D	4 port 2-pos int vol valve replacement	1/32 x .10mm 50C/5000psi liq .02ul PAEK/E
CN4-4344-.02EH	4 port 2-pos int vol valve m-elect	1/32 x .10mm 50C/5000psi liq .02ul PAEK/E
CN4-43R-.004	Rotor seal 4 port CN4-valve int .004ul	1/32 x .10mm 50C/5000psi liq E
CN4-43R-.01	Rotor seal 4 port CN4-valve int .01ul	1/32 x .10mm 50C/5000psi liq E
CN4-43R-.02	Rotor seal 4 port CN4-valve int .02ul	1/32 x .10mm 50C/5000psi liq E
CN4-4C4I	Stator/Insert 4 port CN4-valve	0.10mm bore (.004in) .080in int. E3
C-NEF.5FPK.15S1	Column end fitting Nano 1/32in Cheminert 1um scr	150um PEEK with liners for FS 125-175um
C-NEF.5FPK.20S1	Column end fitting Nano 1/32in Cheminert 1um scr	150um PEEK with liners for FS 175-225um
C-NEF.5FPK.25S1	Column end fitting Nano 1/32in Cheminert 1um scr	150um PEEK with liners for FS 225-275um
C-NEF.5FPK.30S1	Column end fitting Nano 1/32in Cheminert 1um scr	150um PEEK with liners for FS 275-325um
C-NEF.5FPK.35S1	Column end fitting Nano 1/32in Cheminert 1um scr	150um PEEK with liners for FS 325-375um
C-NEF.5XFPK.15S1	Column end fitting Nano 1/32in Cheminert 1um scr	100um PEEK with liners for FS 125-175um
C-NEF.5XFPK.20S1	Column end fitting Nano 1/32in Cheminert 1um scr	100um PEEK with liners for FS 175-225um
C-NEF.5XFPK.25S1	Column end fitting Nano 1/32in Cheminert 1um scr	100um PEEK with liners for FS 225-275um
C-NEF.5XFPK.30S1	Column end fitting Nano 1/32in Cheminert 1um scr	100um PEEK with liners for FS 275-325um
C-NEF.5XFPK.35S1	Column end fitting Nano 1/32in Cheminert 1um scr	100um PEEK with liners for FS 325-375um
C-NEU.5FPK	Union Int/Ext 1/32in Cheminert	150um PEEK without liners
C-NEU.5FPK.15	Union Int/Ext Nano 1/32in Cheminert	150um PEEK with liners for FS 125-175um
C-NEU.5FPK.20	Union Int/Ext Nano 1/32in Cheminert	150um PEEK with liners for FS 175-225um
C-NEU.5FPK.25	Union Int/Ext Nano 1/32in Cheminert	150um PEEK with liners for FS 225-275um
C-NEU.5FPK.30	Union Int/Ext Nano 1/32in Cheminert	150um PEEK with liners for FS 275-325um
C-NEU.5FPK.35	Union Int/Ext Nano 1/32in Cheminert	150um PEEK with liners for FS 325-375um
C-NEU.5XFPK	Union Int/Ext Nano 1/32in Cheminert	100um PEEK without liners
C-NEU.5XFPK.15	Union Int/Ext Nano 1/32in Cheminert	100um PEEK with liners for FS 125-175um
C-NEU.5XFPK.20	Union Int/Ext Nano 1/32in Cheminert	100um PEEK with liners for FS 175-225um
C-NEU.5XFPK.25	Union Int/Ext Nano 1/32in Cheminert	100um PEEK with liners for FS 225-275um
C-NEU.5XFPK.30	Union Int/Ext Nano 1/32in Cheminert	100um PEEK with liners for FS 275-325um
C-NEU.5XFPK.35	Union Int/Ext Nano 1/32in Cheminert	100um PEEK with liners for FS 325-375um
CNFT	Tightening Tool for PEEK hex-head	1/32 nuts green
C-NL.15L-5	Liner long Nano 1/32in Cheminert 27 mm	PEEK for FS OD 125 - 175um 5/pkg
C-NL.15S-5	Liner short Nano 1/32in Cheminert 10 mm	PEEK for FS OD 125 - 175um 5/pkg
C-NL.20L-5	Liner long Nano 1/32in Cheminert 27 mm	PEEK for FS OD 175 - 225um 5/pkg
C-NL.20S-5	Liner long Nano 1/32in Cheminert 10 mm	PEEK for FS OD 175 - 225um 5/pkg
C-NL.25L-5	Liner long Nano 1/32in Cheminert 27 mm	PEEK for FS OD 225 - 275um 5/pkg
C-NL.25S-5	Liner short Nano 1/32in Cheminert 10 mm	PEEK for FS OD 225 - 275um 5/pkg
C-NL.30L-5	Liner long Nano 1/32in Cheminert 27 mm	PEEK for FS OD 275 - 325um 5/pkg
C-NL.30S-5	Liner short Nano 1/32in Cheminert 10 mm	PEEK for FS OD 275 - 325um 5/pkg
C-NL.35L-5	Liner long Nano 1/32in Cheminert 27 mm	PEEK for FS OD 325 - 375um 5/pkg
C-NL.35S-5	Liner short Nano 1/32in Cheminert 10 mm	PEEK for FS OD 325 - 375um 5/pkg
C-NLS1.15	Liner Nano 1/32in Cheminert screen 1um 27 mm	PEEK for FS OD 125-175um
C-NLS1.20	Liner Nano 1/32in Cheminert screen 1um 27 mm	PEEK for FS OD 175 - 225um
C-NLS1.25	Liner Nano 1/32in Cheminert screen 1um 27 mm	PEEK for FS OD 225-275um

# VICI Valco-Cheminert PRODUCT LIST 2012




C-NLS1.30	Liner Nano 1/32in Cheminert screen 1um 27 mm	PEEK for FS OD 275-325um
C-NLS1.35	Liner Nano 1/32in Cheminert screen 1um 27 mm	PEEK for FS OD 325-375um
C-NNFFPK	Nut/Ferrule One-piece 1/32in Cheminert	
C-NNFLFPK	Nut/Ferrule One-piece 1/32in Cheminert	long
C-NPFPK	Plug 1/32 Cheminert	
CNSL1KPK	Cheminert sample loop 1ul (1000nl)	1/32in ends with nuts PEEK
CNSL250PK	Cheminert sample loop 250nl	1/32in ends with nuts PEEK
CNSL2KPK	Cheminert sample loop 2ul (2000nl)	1/32in ends with nuts PEEK
CNSL500PK	Cheminert sample loop 500nl	1/32in ends with nuts PEEK
CNSL5KPK	Cheminert sample loop 5ul (5000nl)	1/32in ends with nuts PEEK
C-NTFPK	Tee Nano 1/32in Cheminert	150um PEEK
C-NTXFPK	Tee Nano 1/32in Cheminert	100um PEEK
CNV1A10S1	Combo needle valve silver colored knob	0-10ml N2 at 40psi 1/16in int Al/Viton
CNV1A150S1	Combo needle valve silver colored knob	150ml N2 at 40psi 1/16in int Al/Viton
CNV1A250S1	Combo needle valve silver colored knob	250ml N2 at 40psi 1/16in int Al/Viton
CNV1A500S1	Combo needle valve silver colored knob	500ml N2 at 40psi 1/16in int Al/Viton
CNV1A50S1	Combo needle valve silver colored knob	50ml N2 at 40psi 1/16in int Al/Viton
CNV1S10S1	Combo needle valve silver colored knob	0-10ml N2 at 40psi 1/16in int SS/Viton
CNV1S150S1	Combo needle valve silver colored knob	150ml N2 at 40psi 1/16in int SS/Viton
CNV1S250S1	Combo needle valve silver colored knob	250ml N2 at 40psi 1/16in int SS/Viton
CNV1S500S1	Combo needle valve silver colored knob	500ml N2 at 40psi 1/16in int SS/Viton
CNV1S50S1	Combo needle valve silver colored knob	50ml N2 at 40psi 1/16in int SS/Viton
CNVEKB	Combo needle valve standard knob black	length .62in /16mm long
CNVEKG	Combo needle valve standard knob green	length .62in /16mm long
CNVEKLB	Combo needle valve long knob black	length 1.25in /32mm long
CNVEKLG	Combo needle valve long knob green	length 1.25in /32mm long
CNVEKLR	Combo needle valve long long knob red	length 1.25in /32mm long
CNVEKLS	Combo needle valve long knob silver	length 1.25in /32mm long
CNVEKLU	Combo needle valve long knob blue	length 1.25in /32mm long
CNVEKR	Combo needle valve standard knob red	length .62in /16mm long
CNVEKS	Combo needle valve standard knob silver	length .62in /16mm long
CNVEKU	Combo needle valve standard knob blue	length .62in /16mm long
C-NVISF	Syringe Fill Port 1/32in Cheminert	for 0.85in /22mm long needle 26 Gauge OD
C-NXFPK	Cross Nano 1/32in Cheminert	150um PEEK
C-NXXFPK	Cross Nano 1/32in Cheminert	100um PEEK
C-NYFPK	Y-assy Nano 1/32in Cheminert	150um PEEK
C-NYXFPK	Y-assy Nano 1/32in Cheminert	100um PEEK
CP2-4841-100D	M6 Pump only replacement	vol. 100ul Sapphire/TF E3/P VICI Logo
CP2-4841-100M1	M6 Pump with controller and stepper motor	vol. 100ul Sapphire/TF E3/P Viton VICI Logo
CP2-4841-100SM	M6 Pump with stepper motor (no controller)	vol. 100ul Sapphire/TF E3/P VICI Logo
CP3-8182-625D	M50 Pump only	vol. 625ul Sapphire/TF ceramic rotor VICI Logo
CP3-8182-625M2	M50 Pump with controller and stepper motor	vol. 625ul Sapphire/TF ceramic rotor VICI Logo
CP3-8182-625SM2	M50 Pump and stepper motor	vol. 625ul Sapphire/TF ceramic rotor VICI Logo
CPA2KF	Cheminert pipe connector 1/8in NPT	1/4in -28female x M-NPT 1.5mm bore CTFE
CPA2PK	Cheminert pipe connector 1/8in NPT	1/4in -28female x M-NPT 1.5mm bore PEEK
CPA4KF	Cheminert pipe connector 1/4in NPT	1/4in -28female x M-NPT 1.5mm bore CTFE
CPA4PK	Cheminert pipe connector 1/4in NPT	1/4in -28female x M-NPT 1.5mm bore PEEK
CP-CM1-P	Controller MicroLynx-4	
CP-DSM	Driver Pump M6	
CP-DSM2	Driver Pump M50	
C-PF	Cheminert perfit 1/4in -28 kit of 6 pcs	for peristaltic tubing .5 to 2.29mm ID
C-PFL	Cheminert perfit 1/4in -28	for peristaltic tubing 1.85 to 2.29mm ID
C-PFM	Cheminert perfit 1/4in -28	for peristaltic tubing 1.12 to 1.65mm ID
C-PFS	Cheminert perfit 1/4in -28	for peristaltic tubing .5 to 1.02mm ID
CPKF	Cheminert plug 1/4in -28	CTFE 5/pkg
CPPK	Cheminert plug 1/4in -28	PEEK 5/pkg
CR10	Clamp ring for standoff assembly	x SOAMP m-pos m-electric actuator
CR11	Clamp ring closomount 2-pos/m-pos vlv.	2-pos/m-pos std. electric/air actuator
CR3	Clamp ring for standoff assembly	x SOA/ x SOAMP
CR4	Clamp ring for close mount assembly	2-pos air/electric actuator and DV valve
CR8	Clamp ring for standoff assembly	x SOA/ x SOAMP 2-pos m-electric actuator

# VICI Valco-Cheminert PRODUCT LIST 2012




CSD10UW	10 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq N60/E
CSD4UW	4 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq N60/E
CSD6UW	6 mpos dead-end path manual	1/16 x .4mm 75C/5000psi liq N60/E
CSL10	Cheminert sample loop 10ul	1/16in ends with nuts & ferrules 316SS
CSL100	Cheminert sample loop 100ul	1/16in ends with nuts & ferrules 316SS
CSL100TI	Cheminert sample loop 100ul	1/16in ends with nuts & ferrules Ti
CSL10K	Cheminert sample loop 10ml	1/16in ends with nuts & ferrules 316SS
CSL10TI	Cheminert sample loop 10ul	1/16in ends with nuts & ferrules Ti
CSL1K	Cheminert sample loop 1ml	1/16in ends with nuts & ferrules 316SS
CSL1KTi	Cheminert sample loop 1ml	1/16in ends with nuts & ferrules Ti
CSL2	Cheminert sample loop 2ul	1/16in ends with nuts & ferrules 316SS
CSL20	Cheminert sample loop 20ul	1/16in ends with nuts & ferrules 316SS
CSL20TI	Cheminert sample loop 20ul	1/16in ends with nuts & ferrules Ti
CSL250	Cheminert sample loop 250ul	1/16in ends with nuts & ferrules 316SS
CSL250TI	Cheminert sample loop 250ul	1/16in ends with nuts & ferrules Ti
CSL2K	Cheminert sample loop 2ml	1/16in ends with nuts & ferrules 316SS
CSL5	Cheminert sample loop 5ul	1/16in ends with nuts & ferrules 316SS
CSL50	Cheminert sample loop 50ul	1/16in ends with nuts & ferrules 316SS
CSL500	Cheminert sample loop 500ul	1/16in ends with nuts & ferrules 316SS
CSL500TI	Cheminert sample loop 500ul	1/16in ends with nuts & ferrules Ti
CSL50TI	Cheminert sample loop 50ul	1/16in ends with nuts & ferrules Ti
CSL5K	Cheminert sample loop 5ml	1/16in ends with nuts & ferrules 316SS
CSLN1K	Cheminert sample loop 1ul (1000nl)	1/32in ends with nuts SS
CSLN2K	Cheminert sample loop 2ul (2000nl)	1/32in ends with nuts SS
CSLN5K	Cheminert sample loop 5ul (5000nl)	1/32in ends with nuts SS
CSS-X1553-200	GC-Column SS ValcoPLOT Molesieve 5A	15m length ID 0.53mm df 20um
CSS-X3053-200	GC-Column SS ValcoPLOT Molesieve 5A	30m length ID 0.53mm df 20um
CST	Cheminert setting tool	for 1/4in -28 flangeless fitting
CTA2KF	Cheminert adapter for tubing 1/8in OD	1.5mm bore CTFE
CTA2PK	Cheminert adapter for tubing 1/8in OD	1.5mm bore PEEK
CTA2S6	Cheminert adapter for tubing 1/8in OD	1.5mm bore 316SS
CTA4KF	Cheminert adapter for tubing 1/4in OD	1.5mm bore CTFE
CTA4PK	Cheminert adapter for tubing 1/4in OD	1.5mm bore PEEK
CTA4S6	Cheminert adapter for tubing 1/4in OD	1.5mm bore 316SS
CTCKF	Cheminert tee 1/4in -28 1/16in OD tubing	.25mm bore CTFE
CTCPK	Cheminert tee 1/4in -28 1/16in OD tubing	.25mm bore PEEK
CTKF	Cheminert tee 1/4in -28 1/16in OD tubing	.5mm bore CTFE
CTLKF	Cheminert tee 1/4in -28 1/8in OD tubing	1.5mm bore CTFE
CTLPK	Cheminert tee 1/4in -28 1/8in OD tubing	1.5mm bore PEEK
CTMKF	Cheminert tee 1/4in -28 1/16in OD tubing	.75mm bore CTFE
CTMPK	Cheminert tee 1/4in -28 1/16in OD tubing	.75mm bore PEEK
CTPK	Cheminert tee 1/4in -28 1/16in OD tubing	.5mm bore PEEK
CU4LKF	Cheminert union 1/4in -28 to 1/2-20	1.5mm bore 1/8in to 1/4in tube CTFE
CU4LPK	Cheminert union 1/4in -28 to 1/2-20	1.6mm bore 1/8in to 1/4in tube PEEK
CUCKF	Cheminert union 1/4in -28 1/16in OD tubing	.25mm bore CTFE
CUCPK	Cheminert union 1/4in -28 1/16in OD tubing	.25mm bore PEEK
CUKF	Cheminert union 1/4in -28 1/16in OD tubing	.5mm bore CTFE
CULKF	Cheminert union 1/4in -28 1/8in OD tubing	1.5mm bore CTFE
CULPK	Cheminert union 1/4in -28 1/8in OD tubing	1.5mm bore PEEK
CUMKF	Cheminert union 1/4in -28 1/16 OD tubing	.75mm bore CTFE
CUMPK	Cheminert union 1/4in -28 1/16in OD tubing	.75mm bore PEEK
CUPK	Cheminert union 1/4in -28 1/16in OD tubing	.5mm bore PEEK
CUTKF	Cheminert union 1/4in -28 1/8in OD tubing	butt connection CTFE
CUTPK	Cheminert union 1/4in -28 OD 1/8in tubing	butt connection PEEK
CUTPP	Cheminert union 1/4in -28 w/o fitting	butt connection Polypropylene 5/pkg
CVA10GS1	Combo needle valve short green knob	10ml/min 40psi He/N2 1/16in int Al/Viton
CVA10GS2	Combo needle valve short green knob	10ml/min 40psi He/N2 1/8in int Al/Viton
CVA150US1	Combo needle valve short blue knob	150ml/min 40psi He/N2 1/16in int Al/Viton
CVA150US2	Combo needle valve short blue knob	150ml/min 40psi He/N2 1/8in int Al/Viton
CVA1KYS1	Combo needle valve short yellow knob	1 liter/min 40psi He/N2 1/16in int Al/Viton
CVA1KYS2	Combo needle valve short yellow knob	1 liter/min 40psi He/N2 1/8in int Al/Viton

# VICI Valco-Cheminert PRODUCT LIST 2012




CVA500BS1	Combo needle valve short black knob	500ml/min 40psi He/N2 1/16in int Al/Viton
CVA500BS2	Combo needle valve short black knob	500ml/min 40psi He/N2 1/8in int Al/Viton
CVA50RS1	Combo needle valve short red knob	50ml/min 40psi He/N2 1/16in int Al/Viton
CVA50RS2	Combo needle valve short red knob	50ml/min 40psi He/N2 1/8in int Al/Viton
C-VISF-1	Syringe fillport Cheminert C22Z polymeric valve	for .75in /19mm long needle 22Gauge OD
C-VISF-1H	Syringe fillport Cheminert C2/C3/C4 polymeric vlv	for .75in /19mm long needle 22Gauge OD
CXCKF	Cheminert cross 1/4in -28 1/16in OD tubing	.25mm bore CTFE
CXCPK	Cheminert cross 1/4in -28 1/16in OD tubing	.25mm bore PEEK
CXKF	Cheminert cross 1/4in -28 1/16in OD tubing	.5mm bore CTFE
CXLKF	Cheminert cross 1/4in -28 1/8in OD tubing	1.5mm bore CTFE
CXLPK	Cheminert cross 1/4in -28 1/8in OD tubing	1.5mm bore PEEK
CXMKF	Cheminert cross 1/4in -28 1/16in OD tubing	.75 mm bore CTFE
CXMPK	Cheminert cross 1/4in -28 1/16in OD tubing	.75mm bore PEEK
CXPK	Cheminert cross 1/4in -28 1/16in OD tubing	.5mm bore PEEK
CZBUCKF	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.25mm bore CTFE
CZBUCPK	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.25mm bore PEEK
CZBUCS6	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.25mm bore 316SS
CZBUKF	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.5mm bore CTFE
CZBUMKF	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.75mm bore CTFE
CZBUMPK	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.75mm bore PEEK
CZBUMS6	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.75mm bore 316SS
CZBUPK	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.5mm bore PEEK
CZBUS6	Cheminert bulkhead union 1/4in -28 x 1/16in OD tubing	.5mm bore 316SS
C-ZP1FPK	Plug 1/16in fingertight high pressure	Cheminert Polymeric vlv series C1-C5 PEEK
CZSL100FEP	Cheminert sample loop 100ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL100PK	Cheminert sample loop 100ul	1/16in PEEK nuts & ferrules PEEK
CZSL100TF	Cheminert sample loop 100ul	1/16in PEEK nuts & ferrules PTFE
CZSL10FEP	Cheminert sample loop 10ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL10PK	Cheminert sample loop 10ul	1/16in PEEK nuts & ferrules PEEK
CZSL10TF	Cheminert sample loop 10ul	1/16in PEEK nuts & ferrules PTFE
CZSL1KFEP	Cheminert sample loop 1ml	1/16in PEEK nuts & CTFE ferrules FEP
CZSL1KPK	Cheminert sample loop 1ml	1/16in PEEK nuts & ferrules PEEK
CZSL1KTF	Cheminert sample loop 1ml	1/16in PEEK nuts & ferrules PTFE
CZSL20FEP	Cheminert sample loop 20ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL20PK	Cheminert sample loop 20ul	1/16in PEEK nuts & ferrules PEEK
CZSL20TF	Cheminert sample loop 20ul	1/16in PEEK nuts & ferrules PTFE
CZSL250FEP	Cheminert sample loop 250ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL250PK	Cheminert sample loop 250ul	1/16in PEEK nuts & ferrules PEEK
CZSL250TF	Cheminert sample loop 250ul	1/16in PEEK nuts & ferrules PTFE
CZSL2KFEP	Cheminert sample loop 2ml	1/16in PEEK nuts & CTFE ferrules FEP
CZSL2KPK	Cheminert sample loop 2ml	1/16in PEEK nuts & ferrules PEEK
CZSL2KTF	Cheminert sample loop 2ml	1/16in PEEK nuts & ferrules PTFE
CZSL2PK	Cheminert sample loop 2ul	1/16in PEEK nuts & ferrules PEEK
CZSL500FEP	Cheminert sample loop 500ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL500PK	Cheminert sample loop 500ul	1/16in PEEK nuts & ferrules PEEK
CZSL500TF	Cheminert sample loop 500ul	1/16in PEEK nuts & ferrules PTFE
CZSL50FEP	Cheminert sample loop 50ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL50PK	Cheminert sample loop 50ul	1/16in PEEK nuts & ferrules PEEK
CZSL50TF	Cheminert sample loop 50ul	1/16in PEEK nuts & ferrules PTFE
CZSL5FEP	Cheminert sample loop 5ul	1/16in PEEK nuts & CTFE ferrules FEP
CZSL5KPK	Cheminert sample loop 5ml	1/16in PEEK nuts & ferrules PEEK
CZSL5PK	Cheminert sample loop 5ul	1/16in PEEK nuts & ferrules PEEK
CZSL5TF	Cheminert sample loop 5ul	1/16in PEEK nuts & ferrules PTFE
CZUCKF	Cheminert union 1/4in -28 x 1/16in OD tubing	.25mm bore CTFE
CZUCPK	Cheminert union 1/4in -28 x 1/16in OD tubing	.25mm bore PEEK
CZUCS6	Cheminert union 1/4in -28 x 1/16in OD tubing	.25mm bore 316SS
CZUKF	Cheminert union 1/4in -28 x 1/16in OD tubing	.5mm bore CTFE
CZUMKF	Cheminert union 1/4in -28 x 1/16in OD tubing	.75mm bore CTFE
CZUMPK	Cheminert union 1/4in -28 x 1/16in OD tubing	.75mm bore PEEK
CZUMS6	Cheminert union 1/4in -28 x 1/16in OD tubing	.75mm bore 316SS
CZUPK	Cheminert union 1/4in -28 x 1/16in OD tubing	.5mm bore PEEK

# VICI Valco-Cheminert PRODUCT LIST 2012




CZUS6	Cheminert union 1/4in -28 x 1/16in OD tubing	.5mm bore 316SS
D10UW	10 port 2-pos valve replacement	1/8 x .75mm 75C/5000psi liq N60/H
D10UWE	10 port 2-pos valve replacement	1/8 x .75mm 225C/400psi gas N60/E
D10UWT	10 port 2-pos valve replacement	1/8 x .75mm 330C/300psi gas N60/T
D-2-110	Pulsed discharge detector PDECD/HID/PID	48V detector controller HP2 110Vac CE
D-2-220	Pulsed discharge detector PDECD/HID/PID	48V detector controller HP2 220Vac CE
D-2-I-110	Pulsed discharge detector optimized HID	48V detector controller HP2 110Vac CE
D-2-I-220	Pulsed discharge detector optimized HID	48V detector controller HP2 220Vac CE
D-3-I-7890-110	Pulsed discharge detector HID/PID HP7890	pulser HP electrometer HP2-110 HPM-110
D-3-I-7890-220	Pulsed discharge detector HID/PID HP7890	pulser HP electrometer HP2-220 HPM-220
D-3-I-HP-110	Pulsed discharge detector HID/PID HP6890	pulser HP electrometer HP2-110 HPM-110
D-3-I-HP-220	Pulsed discharge detector HID/PID HP6890	pulser HP electrometer HP2-220 HPM-220
D-4-I-110	Pulsed discharge detector HeID/PID	48V detector controller HP2 110Vac CE
D-4-I-220	Pulsed discharge detector HeID/PID	48V detector controller HP2 220Vac CE
D-4-I-HP58-110	Pulsed discharge detector HeID/PID HP5890	48V detector controller HP2 110Vac CE
D-4-I-HP58-220	Pulsed discharge detector HeID/PID HP5890	48V detector controller HP2 220Vac CE
D-4-I-SH14-R-110	Pulsed discharge detector HeID/PID Shimadzu 14	w/o controller HP2-110
D-4-I-SH14-R-220	Pulsed discharge detector HeID/PID Shimadzu 14	w/o controller HP2-220
D-4-I-SH17-R-110	Pulsed discharge detector HeID/PID Shimadzu 17	w/o controller HP2-110 suitable for 2010 & 2014
D-4-I-SH17-R-220	Pulsed discharge detector HeID/PID Shimadzu 17	w/o controller HP2-220 suitable for 2010 & 2014
D-4-I-TQ-R-110	Pulsed discharge detector HeID/PID for Trace GC	pulser HP2-110 ThermoQuest
D-4-I-TQ-R-220	Pulsed discharge detector HeID/PID for Trace GC	pulser HP2-220 ThermoQuest
D-4-I-VA38-R-110	Pulsed discharge detector HeID/PID Varian 3800	w/o controller with HP2-110
D-4-I-VA38-R-220	Pulsed discharge detector HeID/PID Varian 3800	w/o controller with HP2-220
D4UW	4 port 2-pos valve replacement	1/8 x .75mm 75C/5000psi liq N60/H
D4UWE	4 port 2-pos valve replacement	1/8 x .75mm 225C/400psi gas N60/E
D4UWT	4 port 2-pos valve replacement	1/8 x .75mm 330C/300psi gas N60/T
D-5-6890-110	Pulsed discharge detector PDECD mode HP 6890	pulser HP electrometer HP2-110 HPM-110
D-5-6890-220	Pulsed discharge detector PDECD mode HP 6890	pulser HP electrometer HP2-220 HPM-220
D6UW	6 port 2-pos valve replacement	1/8 x .75mm 75C/5000psi liq N60/H
D6UWE	6 port 2-pos valve replacement	1/8 x .75mm 225C/400psi gas N60/E
D6UWT	6 port 2-pos valve replacement	1/8 x .75mm 330C/300psi gas N60/T
D8UW	8 port 2-pos valve replacement	1/8 x .75mm 75C/5000psi liq N60/H
D8UWE	8 port 2-pos valve replacement	1/8 x .75mm 225C/400psi gas N60/E
D8UWT	8 port 2-pos valve replacement	1/8 x .75mm 330C/300psi gas N60/T
DC10UW	10 port 2-pos valve replacement	1/16 x .75mm 75C/5000psi liq N60/H
DC10UWE	10 port 2-pos valve replacement	1/16 x .75mm 225C/400psi gas N60/E
DC10UWT	10 port 2-pos valve replacement	1/16 x .75mm 330C/300psi gas N60/T
DC10W	10 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
DC10WE	10 port 2-pos valve replacement	1/16 x .4mm 225C/400psi gas N60/E
DC10WT	10 port 2-pos valve replacement	1/16 x .4mm 350C/300psi gas N60/T
DC4UW	4 port 2-pos valve replacement	1/16 x .75mm 75C/5000psi liq N60/H
DC4UWE	4 port 2-pos valve replacement	1/16 x .75mm 225C/400psi gas N60/E
DC4UWT	4 port 2-pos valve replacement	1/16 x .75mm 330C/300psi gas N60/T
DC4W	4 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
DC4WE	4 port 2-pos valve replacement	1/16 x .4mm 225C/400psi gas N60/E
DC4WT	4 port 2-pos valve replacement	1/16 x .4mm 350C/300psi gas N60/T
DC6UW	6 port 2-pos valve replacement	1/16 x .75mm 75C/5000psi liq N60/H
DC6UWE	6 port 2-pos valve replacement	1/16 x .75mm 225C/400psi gas N60/E
DC6UWT	6 port 2-pos valve replacement	1/16 x .75mm 330C/300psi gas N60/T
DC6W	6 port 2-pos valve K-hole replacement	1/16 x .4mm 75C/5000psi liq N60/H
DC6WE	6 port 2-pos valve replacement	1/16 x .4mm 225C/400psi gas N60/E
DC6WT	6 port 2-pos valve replacement	1/16 x .4mm 350C/300psi gas N60/T
DC8UW	8 port 2-pos valve replacement	1/16 x .75mm 75C/5000psi liq N60/H
DC8UWE	8 port 2-pos valve replacement	1/16 x .75mm 225C/400psi gas N60/E
DC8UWT	8 port 2-pos valve replacement	1/16 x .75mm 330C/300psi gas N60/T
DC8W	8 port 2-pos valve replacement	1/16 x .4mm 75C/5000psi liq N60/H
DC8WE	8 port 2-pos valve replacement	1/16 x .4mm 225C/400psi gas N60/E
DC8WT	8 port 2-pos valve replacement	1/16 x .4mm 350C/300psi gas N60/T
DCI4UW.2	4 port 2-pos int vol valve replacement	1/16 x .75mm 75C/5000psi liq .2ul N60/H
DCI4UW.5	4 port 2-pos int vol valve replacement	1/16 x .75mm 75C/5000psi liq .5ul N60/H

# VICI Valco-Cheminert PRODUCT LIST 2012




DCI4UW1	4 port 2-pos int vol valve replacement	1/16 x .75mm 75C/5000psi liq 1ul N60/H
DCI4UW2	4 port 2-pos int vol valve replacement	1/16 x .75mm 75C/5000psi liq 2ul N60/H
DCI4UWE.2	4 port 2-pos int vol valve replacement	1/16 x .75mm 175C/1000psi liq .2ul N60/E
DCI4UWE.5	4 port 2-pos int vol valve replacement	1/16 x .75mm 175C/1000psi liq .5ul N60/E
DCI4UWE1	4 port 2-pos int vol valve replacement	1/16 x .75mm 175C/1000psi liq 1ul N60/E
DCI4UWE2	4 port 2-pos int vol valve replacement	1/16 x .75mm 175C/1000psi liq 2ul N60/E
DCI4W.06	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/5000psi liq .06ul N60/H
DCI4W.1	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/5000psi liq .1ul N60/H
DCI4W.2	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/5000psi liq .2ul N60/H
DCI4W.5	4 port 2-pos int vol valve replacement	1/16 x .4mm 75C/5000psi liq .5ul N60/H
DCI4WE.06	4 port 2-pos int vol valve replacement	1/16 x .4mm 175C/1000psi liq .06ulN60/E
DCI4WE.1	4 port 2-pos int vol valve replacement	1/16 x .4mm 175C/1000psi liq .1ul N60/E
DCI4WE.2	4 port 2-pos int vol valve replacement	1/16 x .4mm 175C/1000psi liq .2ul N60/E
DCI4WE.5	4 port 2-pos int vol valve replacement	1/16 x .4mm 175C/1000psi liq .5ul N60/E
DCSC10MWE	10 mpos common outlet path MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSC12MWE	12 mpos common outlet path MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSC16MWE	16 mpos common outlet path MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSC6MWE	6 mpos common outlet path MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSD10MWE	10 mpos dead-end path MW replacement	1/16 x .75mm 200C/400psi gas N60/E
DCSD10UW	10 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq N60/E
DCSD12MWE	12 mpos dead-end path MW replacement	1/16 x .75mm 200C/400psi gas N60/E
DCSD16MWE	16 mpos dead-end path MW replacement	1/16 x .75mm 200C/400psi gas N60/E
DCSD4UW	4 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq N60/E
DCSD6MWE	6 mpos dead-end path MW replacement	1/16 x .75mm 200C/400psi gas N60/E
DCSD6UW	6 mpos dead-end path replacement	1/16 x .4mm 75C/5000psi liq N60/E
DCSF10MWE	10 mpos flow-through MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSF12MWE	12 mpos flow-through MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSF16MWE	16 mpos flow-through MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCSF6MWE	6 mpos flow-through MW replacement	1/16 x 1mm 200C/200psi gas N60/E
DCST10MWE	10 mpos trapping path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCST12MWE	12 mpos trapping path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCST16MWE	16 mpos trapping path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCST4UW	4 mpos trapping path replacement	1/16 x .4mm 75C/5000psi liq N60/E
DCST6MWE	6 mpos trapping path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCST6UW	6 mpos trapping path replacement	1/16 x .4mm 75C/5000psi liq N60/E
DCSTF10MWE	10 mpos trap/flow path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCSTF12MWE	12 mpos trap/flow path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCSTF16MWE	16 mpos trap/flow path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DCSTF6MWE	6 mpos trap/flow path MW replacement	1/16 x .75mm 200C/200psi gas N60/E
DI4UWE.2	4 port 2-pos int vol valve replacement	1/8 x .75mm 175C/1000psi liq .2ul N60/E
DI4UWE.5	4 port 2-pos int vol valve replacement	1/8 x .75mm 175C/1000psi liq .5ul N60/E
DI4UWE1	4 port 2-pos int vol valve replacement	1/8 x .75mm 175C/1000psi liq 1ul N60/E
DI4UWE2	4 port 2-pos int vol valve replacement	1/8 x .75mm 175C/1000psi liq 2ul N60/E
DL10UW	10 port 2-pos valve replacement	1/8 x 1.0mm 75C/5000psi liq N60/H
DL4UW	4 port 2-pos valve replacement	1/8 x 1.7mm 75C/5000psi liq N60/H
DL6UW	6 port 2-pos valve replacement	1/8 x 1.7mm 75C/5000psi liq N60/H
DL8UW	8 port 2-pos valve replacement	1/8 x 1.3mm 75C/5000psi liq N60/H
DN10WE	10 port 2-pos valve replacement	1/32 x .25mm 225C/400psi gas N60/E
DN10WT	10 port 2-pos valve replacement	1/32 x .25mm 350C/300psi gas N60/T
DN4WE	4 port 2-pos valve replacement	1/32 x .25mm 225C/400psi gas N60/E
DN4WT	4 port 2-pos valve replacement	1/32 x .25mm 350C/300psi gas N60/T
DN6WE	6 port 2-pos valve replacement	1/32 x .25mm 225C/400psi gas N60/E
DN6WT	6 port 2-pos valve replacement	1/32 x .25mm 350C/300psi gas N60/T
DN8WE	8 port 2-pos valve replacement	1/32 x .25mm 225C/400psi gas N60/E
DN8WT	8 port 2-pos valve replacement	1/32 x .25mm 350C/300psi gas N60/T
DNI4WE.06	4 port 2-pos int vol valve replacement	1/32 x .25mm 175C/1000psi liq .06ulN60/E
DNI4WE.1	4 port 2-pos int vol valve replacement	1/32 x .25mm 175C/1000psi liq .1ulN60/E
DNI4WE.2	4 port 2-pos int vol valve replacement	1/32 x .25mm 175C/1000psi liq .2ulN60/E
DNI4WE.5	4 port 2-pos int vol valve replacement	1/32 x .25mm 175C/1000psi liq .5ulN60/E
DSC10MWE	10 mpos common outlet path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSC12MWE	12 mpos common outlet path MW replacement	1/8 x 1mm 200C/200psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




DSC16MWE	16 mpos common outlet path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSC6MWE	6 mpos common outlet path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSD10MWE	10 mpos dead-end path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSD12MWE	12 mpos dead-end path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSD16MWE	16 mpos dead-end path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSD4UW	4 mpos dead-end path replacement	1/8 x .75mm 75C/5000psi liq N60/E
DSD6MWE	6 mpos dead-end path MW replacement	1/8 x 1mm 200C/400psi gas N60/E
DSD6UW	6 mpos dead-end path replacement	1/8 x .75mm 75C/5000psi liq N60/E
DSD8UW	8 mpos dead-end path replacement	1/8 x .75mm 75C/5000psi liq N60/E
DSF10MWE	10 mpos flow-through MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSF12MWE	12 mpos flow-through MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSF16MWE	16 mpos flow-through MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSF6MWE	6 mpos flow-through MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DST10MWE	10 mpos trapping path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DST12MWE	12 mpos trapping path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DST16MWE	16 mpos trapping path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DST6MWE	6 mpos trapping path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSTF10MWE	10 mpos trap/flow path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSTF12MWE	12 mpos trap/flow path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSTF16MWE	16 mpos trap/flow path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DSTF6MWE	6 mpos trap/flow path MW replacement	1/8 x 1mm 200C/200psi gas N60/E
DV12-1110	Mini diaphragm 10 port 2-pos valve	1/32 x .25mm 200C/300psi liq N60/Polyimide
DV12-1114-.5	Mini diaphragm 4 port 2-pos int/vol valve	1/32 x .25mm 50C/750psi liq. N60/Polyimide .5ul
DV12-1114-1	Mini diaphragm 4 port 2-pos int/vol valve	1/32 x .25mm 50C/750psi liq. N60/Polyimide 1ul
DV12-1116	Mini diaphragm 6 port 2-pos valve	1/32 x .25mm 200C/300psi liq N60/Polyimide
DV22-2110	Mini diaphragm 10 port 2-pos valve	1/16 x .4mm 200C/300psi liq N60/Polyimide
DV22-2114-.5	Mini diaphragm 4 port 2-pos int/vol valve	1/16 x .4mm 50C/750psi liq N60/Polyimide .5ul
DV22-2114-1	Mini diaphragm 4 port 2-pos int/vol valve	1/16 x .4mm 50C/750psi liq N60/Polyimide 1ul
DV22-2116	Mini diaphragm 6 port 2-pos valve	1/16 x .4mm 200C/300psi liq N60/Polyimide
DV22-21D	Diaphragm replacement 6/10 port Polyimide	for valve bore 0.25 and 0.40 mm
DV22-22D	Diaphragm replacement 6/10 port Teflon	
DV22-3110	Mini diaphragm 10 port 2-pos valve	1/16 x .75mm 200C/300psi liq N60/Polyimide
DV22-3114-.5	Mini diaphragm 4 port 2-pos int/vol valve	1/16 x .75mm 50C/750psi liq N60/Polyimide .5ul
DV22-3114-1	Mini diaphragm 4 port 2-pos int/vol valve	1/16 x .75mm 50C/750psi liq N60/Polyimide 1ul
DV22-3116	Mini diaphragm 6 port 2-pos valve	1/16 x .75mm 200C/300psi liq N60/Polyimide
DV22-31D	Diaphragm replacement 6/10 port Polyimide	for valve bore 0.75 mm
DV22-PURGE	Purge ring for DV22	
DVI-110	Digital valve interface 110Vac	Non-CE !!!
DVI-220	Digital valve interface 220Vac	Non-CE !!!
DVL4MWE2	4 port 2-pos valve MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVL6MWE2	6 port 2-pos valve MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVL8MWE2	8 port 2-pos valve MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSC4MWE2	4 mpos common outlet path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSC6MWE2	6 mpos common outlet path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSC8MWE2	8 mpos common outlet path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSD10MWE2	10 mpos dead-end path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSD4MWE2	4 mpos dead-end path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSD6MWE2	6 mpos dead-end path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSD8MWE2	8 mpos dead-end path MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSF4MWE2	4 mpos flow-through MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSF6MWE2	6 mpos flow-through MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVLSF8MWE2	8 mpos flow-through MW replacement	1/4 x 4mm 75C/100psi gas N60/E2
DVSP4-220	Digital valve sequence programmer	4 intervals 220Vac Non-CE !!!
E10-110	Electric actuator mpos close mount	10 position 110Vac
E102-110	Electric actuator mpos 2in standoff	10 position 110Vac
E10S-110	Electric actuator mpos for standoff	10 position 110Vac
E10UW	10 port 2-pos valve std-electric	1/8 x .75mm 75C/5000psi liq N60/H
E12-110	Electric actuator mpos close mount	12 position 110Vac
E122-110	Electric actuator mpos 2in standoff	12 position 110 Vac
E12S-110	Electric actuator mpos for standoff	12 position 110Vac
E16-110	Electric actuator mpos close mount	16 position 110Vac

# VICI Valco-Cheminert PRODUCT LIST 2012




E162-110	Electric actuator mpos 2in standoff	16 position 110Vac
E16S-110	Electric actuator mpos for standoff	16 position 110Vac
E2CI4UWE.2	4 port 2-pos int vol valve std-elect/2in so	1/16 x .75mm 175C/1000psi liq .2ul N60/E
E2CI4UWE.5	4 port 2-pos int vol valve std-elect/2in so	1/16 x .75mm 175C/1000psi liq .5ul N60/E
E2CI4UWE1	4 port 2-pos int vol valve std-elect/2in so	1/16 x .75mm 175C/1000psi liq 1ul N60/E
E2CI4UWE2	4 port 2-pos int vol valve std-elect/2in so	1/16 x .75mm 175C/1000psi liq 2ul N60/E
E2CI4WE.06	4 port 2-pos int vol valve std-elect/2in so	1/16 x .4mm 175C/1000psi liq .06ul N60/E
E2CI4WE.1	4 port 2-pos int vol valve std-elect/2in so	1/16 x .4mm 175C/1000psi liq .1ul N60/E
E2CI4WE.2	4 port 2-pos int vol valve std-elect/2in so	1/16 x .4mm 175C/1000psi liq .2ul N60/E
E2CI4WE.5	4 port 2-pos int vol valve std-elect/2in so	1/16 x .4mm 175C/1000psi liq .5ul N60/E
E2CSC10MWE	10 mpos common outlet path MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSC12MWE	12 mpos common outlet path MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSC16MWE	16 mpos common outlet path MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSC6MWE	6 mpos common outlet path MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSD10MWE	10 mpos dead-end path MW std-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
E2CSD12MWE	12 mpos dead-end path MW std-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
E2CSD16MWE	16 mpos dead-end path MW std-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
E2CSD6MWE	6 mpos dead-end path MW std-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
E2CSF10MWE	10 mpos flow-through MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSF12MWE	12 mpos flow-through MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSF16MWE	16 mpos flow-through MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CSF6MWE	6 mpos flow-through MW std-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
E2CST10MWE	10 mpos trapping path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CST12MWE	12 mpos trapping path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CST16MWE	16 mpos trapping path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CST6MWE	6 mpos trapping path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CSTF10MWE	10 mpos trap/flow path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CSTF12MWE	12 mpos trap/flow path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CSTF16MWE	16 mpos trap/flow path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2CSTF6MWE	6 mpos trap/flow path MW std-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
E2I4UWE.2	4 port 2-pos int vol valve std-elect/2in so	1/8 x .75mm 175C/1000psi liq .2ul N60/E
E2I4UWE.5	4 port 2-pos int vol valve std-elect/2in so	1/8 x .75mm 175C/1000psi liq .5ul N60/E
E2I4UWE1	4 port 2-pos int vol valve std-elect/2in so	1/8 x .75mm 175C/1000psi liq 1ul N60/E
E2I4UWE2	4 port 2-pos int vol valve std-elect/2in so	1/8 x .75mm 175C/1000psi liq 2ul N60/E
E2NI4WE.06	4 port 2-pos int vol valve std-elect/2in so	1/32 x .25mm 175C/1000psi liq .06ulN60/E
E2NI4WE.1	4 port 2-pos int vol valve std-elect/2in so	1/32 x .25mm 175C/1000psi liq .1ulN60/E
E2NI4WE.2	4 port 2-pos int vol valve std-elect/2in so	1/32 x .25mm 175C/1000psi liq .2ulN60/E
E2NI4WE.5	4 port 2-pos int vol valve std-elect/2in so	1/32 x .25mm 175C/1000psi liq .5ulN60/E
E2SC10MWE	10 mpos common outlet path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SC12MWE	12 mpos common outlet path MW std-elect/4in so	1/8 x 1mm 200C/200psi gas N60/E
E2SC16MWE	16 mpos common outlet path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SC6MWE	6 mpos common outlet path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SD10MWE	10 mpos dead-end path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SD12MWE	12 mpos dead-end path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SD16MWE	16 mpos dead-end path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SD6MWE	6 mpos dead-end path MW std-elect/2in so	1/8 x 1mm 200C/400psi gas N60/E
E2SF10MWE	10 mpos flow-through MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SF12MWE	12 mpos flow-through MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SF16MWE	16 mpos flow-through MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2SF6MWE	6 mpos flow-through MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2ST10MWE	10 mpos trapping path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2ST12MWE	12 mpos trapping path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2ST16MWE	16 mpos trapping path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2ST6MWE	6 mpos trapping path MW std-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2STF10MWE	10 mpos trap/flow path MW std-elec/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2STF12MWE	12 mpos trap/flow path MW std-elec/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2STF16MWE	16 mpos trap/flow path MW std-elec/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2STF6MWE	6 mpos trap/flow path MW std-elec/2in so	1/8 x 1mm 200C/200psi gas N60/E
E2VLSC4MWE2	4 mpos common outlet path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLSC6MWE2	6 mpos common outlet path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLSC8MWE2	8 mpos common outlet path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2

# VICI Valco-Cheminert PRODUCT LIST 2012




E2VLS10MWE2	10 mpos dead-end path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLS4MWE2	4 mpos dead-end path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLS6MWE2	6 mpos dead-end path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLS8MWE2	8 mpos dead-end path MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLSF4MWE2	4 mpos flow-through MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLSF6MWE2	6 mpos flow-through MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E2VLSF8MWE2	8 mpos flow-through MW std-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
E30-110	Electric actuator 2-pos close mount	30 deg rotation 110Vac
E302-110	Electric actuator 2-pos 2in standoff	30 deg rotation 110Vac
E302-CE	Electric actuator 2-pos 2in standoff	30 deg rot. 24Vdcw.110/220Vac p.supply
E30-CE	Electric actuator 2-pos close mount	30 deg rot. 24Vdcw.110/220Vac p.supply
E30S-110	Electric actuator 2-pos for standoff	30 deg rotation 110Vac
E30S-CE	Electric actuator 2-pos for standoff	30 deg rot. 24Vdcw.110/220Vac p.supply
E36-110	Electric actuator 2-pos close mount	36 deg rotation 110Vac
E362-110	Electric actuator 2-pos 2in standoff	36 deg rotation 110Vac
E362-CE	Electric actuator 2-pos 2in standoff	36 deg rot. 24Vdcw.110/220Vac p.supply
E36-CE	Electric actuator 2-pos close mount	36 deg rot. 24Vdcw.110/220Vac p.supply
E36S-110	Electric actuator 2-pos for standoff	36 deg rotation 110Vac
E36S-CE	Electric actuator 2-pos for standoff	36 deg rot. 24Vdcw.110/220Vac p.supply
E410UWE	10 port 2-pos valve std-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
E410UWT	10 port 2-pos valve std-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
E4-110	Electric actuator mpos close mount	4 position 110Vac
E42-110	Electric actuator mpos 2in standoff	4 position 110Vac
E44UWE	4 port 2-pos valve std-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
E44UWT	4 port 2-pos valve std-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
E45-110	Electric actuator 2-pos close mount	45 deg rotation 110Vac
E452-110	Electric actuator 2-pos 2in standoff	45 deg rotation 110Vac
E452-CE	Electric actuator 2-pos 2in standoff	45 deg rot. 24Vdcw.110/220Vac p.supply
E45-CE	Electric actuator 2-pos close mount	45 deg rot.24Vdcw. 110/220Vac p.supply
E45S-110	Electric actuator 2-pos for standoff	45 deg rotation 110Vac
E45S-CE	Electric actuator 2-pos for standoff	45 deg rot. 24Vdcw.110/220Vac p.supply
E46UWE	6 port 2-pos valve std-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
E46UWT	6 port 2-pos valve std-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
E48UWE	8 port 2-pos valve std-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
E48UWT	8 port 2-pos valve std-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
E4C10UWE	10 port 2-pos valve std-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
E4C10UWT	10 port 2-pos valve std-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
E4C10WE	10 port 2-pos valve std-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
E4C10WT	10 port 2-pos valve std-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
E4C4UWE	4 port 2-pos valve std-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
E4C4UWT	4 port 2-pos valve std-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
E4C4WE	4 port 2-pos valve std-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
E4C4WT	4 port 2-pos valve std-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
E4C6UWE	6 port 2-pos valve std-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
E4C6UWT	6 port 2-pos valve std-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
E4C6WE	6 port 2-pos valve std-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
E4C6WT	6 port 2-pos valve std-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
E4C8UWE	8 port 2-pos valve std-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
E4C8UWT	8 port 2-pos valve std-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
E4C8WE	8 port 2-pos valve std-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
E4C8WT	8 port 2-pos valve std-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
E4N10WE	10 port 2-pos valve std-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
E4N10WT	10 port 2-pos valve std-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
E4N4WE	4 port 2-pos valve std-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
E4N4WT	4 port 2-pos valve std-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
E4N6WE	6 port 2-pos valve std-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
E4N6WT	6 port 2-pos valve std-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
E4N8WE	8 port 2-pos valve std-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
E4N8WT	8 port 2-pos valve std-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
E4S-110	Electric actuator mpos for standoff	4 position 110Vac
E4UW	4 port 2-pos valve std-electric	1/8 x .75mm 75C/5000psi liq N60/H

# VICI Valco-Cheminert PRODUCT LIST 2012




E4VL4MWE2	4 port 2-pos valve MW std-elect/4in so	1/4 x 4mm 75C/100psi gas N60/E2
E4VL6MWE2	6 port 2-pos valve MW std-elect/4in so	1/4 x 4mm 75C/100psi gas N60/E2
E4VL8MWE2	8 port 2-pos valve MW std-elect/4in so	1/4 x 4mm 75C/100psi gas N60/E
E4X2-110	Electric actuator mpos close mount	4x2 position 110Vac
E4X22-110	Electric actuator mpos 2in standoff	4x2 position 110Vac
E4X2S-110	Electric actuator mpos for standoff	4x2 position 110Vac
E60-110	Electric actuator 2-pos close mount	60 deg rotation 110Vac
E602-110	Electric actuator 2-pos 2in standoff	60 deg rotation 110Vac
E602-CE	Electric actuator 2-pos 2in standoff	60 deg rot. 24Vdcw.110/220Vac p.supply
E60-CE	Electric actuator 2-pos close mount	60 deg rot. 24Vdcw.110/220Vac p.supply
E60S-110	Electric actuator 2-pos for standoff	60 deg rotation 110Vac
E60S-CE	Electric actuator 2-pos for standoff	60 deg rot. 24Vdcw.110/220Vac p.supply
E6-110	Electric actuator mpos close mount	6 position 110Vac
E62-110	Electric actuator mpos 2in standoff	6 position 110Vac
E6S-110	Electric actuator mpos for standoff	6 position 110Vac
E6UW	6 port 2-pos valve std-electric	1/8 x .75mm 75C/5000psi liq N60/H
E6X2-110	Electric actuator mpos close mount	6x2 position 110Vac
E6X22-110	Electric actuator mpos 2in standoff	6x2 position 110Vac
E6X2S-110	Electric actuator mpos for standoff	6x2 position 110Vac
E8-110	Electric actuator mpos close mount	8 position 110Vac
E82-110	Electric actuator mpos 2in standoff	8 position 110Vac
E8S-110	Electric actuator mpos for standoff	8 position 110Vac
E8UW	8 port 2-pos valve std-electric	1/8 x .75mm 75C/5000psi liq N60/H
E90-110	Electric actuator 2-pos close mount	90 deg rotation 110Vac
E902-110	Electric actuator 2-pos 2in standoff	90 deg rotation 110Vac
E902-CE	Electric actuator 2-pos 2in standoff	90 deg rot. 24Vdcw.110/220Vac p.supply
E90-CE	Electric actuator 2-pos close mount	90 deg rot. 24Vdcw.110/220Vac p.supply
E90S-110	Electric actuator 2-pos for standoff	90 deg rotation 110Vac
E90S-CE	Electric actuator 2-pos for standoff	90 deg rot. 24Vdcw.110/220Vac p.supply
EAOR21	External adaptor 1/8in to 10-32 thread	O-ring seal SS
EAOR22	External adaptor 1/8in to 5/16-24 thread	O-ring seal SS
EBRU12L	Bulkhead external red. union 1/8in x 1/16in	1mm bore SS
EBRU12T	Bulkhead external red. union 1/8in x 1/16in	1/16in trough bore SS
EBRU14L	Bulkhead external red. union 1/4in x 1/16in	1mm bore SS
EBRU14T	Bulkhead external red. union 1/4in x 1/16in	1/16in trough bore SS
EBRU24L	Bulkhead external red. union 1/4in x 1/8in	2mm bore SS
EBU2L	Bulkhead external union 1/8in	2mm bore SS
EBU2T	Bulkhead external union 1/8in	1/8in through bore SS
EBU4	Bulkhead external union 1/4in	2mm bore SS
EBU4L	Bulkhead external union 1/4in	4.6mm bore SS
EBU4T	Bulkhead external union 1/4in	1/4in through bore SS
EC10UW	10 port 2-pos valve std-electric	1/16 x .75mm 75C/5000psi liq N60/H
EC10W	10 port 2-pos valve std-electric	1/16 x .4mm 75C/5000psi liq N60/H
EC4UW	4 port 2-pos valve std-electric	1/16 x .75mm 75C/5000psi liq N60/H
EC4W	4 port 2-pos valve std-electric	1/16 x .4mm 75C/5000psi liq N60/H
EC6UW	6 port 2-pos valve std-electric	1/16 x .75mm 75C/5000psi liq N60/H
EC6W	6 port 2-pos valve std-electric	1/16 x .4mm 75C/5000psi liq N60/H
EC8UW	8 port 2-pos valve std-electric	1/16 x .75mm 75C/5000psi liq N60/H
EC8W	8 port 2-pos valve std-electric	1/16 x .4mm 75C/5000psi liq N60/H
ECEF111.0	External column end fitting 1/16in	.25mm bore 1mm cone without frit 316SS
ECEF111.0F	External column end fitting 1/16in	.25mm bore 1mm cone removable 2um frit 316SS
ECEF1K1	External column end fitting 1in x 1/16in	.75mm bore 20mm cone without frit 316SS
ECEF1K1F	External column end fitting 1in x 1/16in	.75mm bore 20mm cone removable 2um frit 316SS
ECEF211.0	External column end fitting 1/8in x 1/16in	.25mm bore 1mm cone without frit 316SS
ECEF211.0F	External column end fitting 1/8in x 1/16in	.25mm bore 1mm cone removable 2um frit 316SS
ECEF412.1	External column end fitting 1/4in x 1/16in	.4mm bore 2.1mm conewithout frit 316SS
ECEF412.1F	External column end fitting 1/4in x 1/16in	.4mm bore 2.1mm coneremovable 2um frit 316SS
ECEF413.0	External column end fitting 1/4in x 1/16in	.4mm bore 3mm cone without frit 316SS
ECEF413.0F	External column end fitting 1/4in x 1/16in	.4mm bore 3mm cone removable 2um frit 316SS
ECEF414.0	External column end fitting 1/4in x 1/16in	.4mm bore 4mm cone without frit 316SS
ECEF414.0F	External column end fitting 1/4in x 1/16in	.4mm bore 4mm cone removable 2um frit 316SS

# VICI Valco-Cheminert PRODUCT LIST 2012




ECEF414.6	External column end fitting 1/4in x 1/16in	.4mm bore 4.6mm conewithout frit 316SS
ECEF414.6F	External column end fitting 1/4in x 1/16in	.4mm bore 4.6mm coneremovable 2um frit 316SS
ECEF616.0	External column end fitting 3/8in x 1/16in	.35mm bore 6mm cone without frit 316SS
ECEF616.0F	External column end fitting 3/8in x 1/16in	.35mm bore 6mm cone removable 2um frit 316SS
ECEF617.0	External column end fitting 3/8in x 1/16in	.4mm bore 7mm cone without frit 316SS
ECEF617.0F	External column end fitting 3/8in x 1/16in	.4mm bore 7mm cone removable 2um frit 316SS
ECEF8110.0	External column end fitting 1/2in x 1/16in	.75mm bore 10mm cone without frit 316SS
ECEF8110.0F	External column end fitting 1/2in x 1/16in	.75mm bore 10mm cone removable 2um frit 316SS
ECEF819.0	External column end fitting 1/2in x 1/16in	.75mm bore 9mm cone without frit 316SS
ECEF819.0F	External column end fitting 1/2in x 1/16in	.75mm bore 9mm cone removable 2um frit 316SS
ECI4UW.2	4 port 2-pos int vol valve std-electric	1/16 x .75mm 75C/5000psi liq .2ul N60/H
ECI4UW.5	4 port 2-pos int vol valve std-electric	1/16 x .75mm 75C/5000psi liq .5ul N60/H
ECI4UW1	4 port 2-pos int vol valve std-electric	1/16 x .75mm 75C/5000psi liq 1ul N60/H
ECI4UW2	4 port 2-pos int vol valve std-electric	1/16 x .75mm 75C/5000psi liq 2ul N60/H
ECI4W.06	4 port 2-pos int vol valve std-electric	1/16 x .4mm 75C/5000psi liq .06ul N60/H
ECI4W.1	4 port 2-pos int vol valve std-electric	1/16 x .4mm 75C/5000psi liq .1ul N60/H
ECI4W.2	4 port 2-pos int vol valve std-electric	1/16 x .4mm 75C/5000psi liq .2ul N60/H
ECI4W.5	4 port 2-pos int vol valve std-electric	1/16 x .4mm 75C/5000psi liq .5ul N60/H
ECMH	OBSOLETE! Please ask for alternatives!	
ECMH2	OBSOLETE! Please ask for alternatives!	
ECMH2-CE	OBSOLETE! Please ask for alternatives!	
ECMH-CE	OBSOLETE! Please ask for alternatives!	
ECMHS	OBSOLETE! Please ask for alternatives!	
ECMHS-CE	OBSOLETE! Please ask for alternatives!	
ECMT	OBSOLETE! Please ask for alternatives!	
ECMT2	OBSOLETE! Please ask for alternatives!	
ECMT2-CE	OBSOLETE! Please ask for alternatives!	
ECMT-CE	OBSOLETE! Please ask for alternatives!	
ECMTS	OBSOLETE! Please ask for alternatives!	
ECMTS-CE	OBSOLETE! Please ask for alternatives!	
ECSD10UW	10 mpos dead-end path std-electric	1/16 x .4mm 75C/5000psi liq N60/E
ECSD4UW	4 mpos dead-end path std-electric	1/16 x .4mm 75C/5000psi liq N60/E
ECSD6UW	6 mpos dead-end path std-electric	1/16 x .4mm 75C/5000psi liq N60/E
ECST4UW	4 mpos trapping path std-electric	1/16 x .4mm 75C/5000psi liq N60/E
ECST6UW	6 mpos trapping path std-electric	1/16 x .4mm 75C/5000psi liq N60/E
ED	Micro electric actuator 2-pos closemount	medium torque 24Vdc w.110/220Vac p-supply
ED10UW	10 port 2-pos valve m-electric	1/8 x .75mm 75C/5000psi liq N60/H
ED2	Micro electric actuator 2-pos/2in so	high torque 24Vdc w.110/220Vac p-supply
ED2CI4UWE.2	4 port 2-pos int vol valve m-elect/2in so	1/16 x .75mm 175C/1000psi liq .2ul N60/E
ED2CI4UWE.5	4 port 2-pos int vol valve m-elect/2in so	1/16 x .75mm 175C/1000psi liq .5ul N60/E
ED2CI4UWE1	4 port 2-pos int vol valve m-elect/2in so	1/16 x .75mm 175C/1000psi liq 1ul N60/E
ED2CI4UWE2	4 port 2-pos int vol valve m-elect/2in so	1/16 x .75mm 175C/1000psi liq 2ul N60/E
ED2I4UWE.2	4 port 2-pos int vol valve m-elect/2in so	1/8 x .75mm 175C/1000psi liq .2ul N60/E
ED2I4UWE.5	4 port 2-pos int vol valve m-elect/2in so	1/8 x .75mm 175C/1000psi liq .5ul N60/E
ED2I4UWE1	4 port 2-pos int vol valve m-elect/2in so	1/8 x .75mm 175C/1000psi liq 1ul N60/E
ED2I4UWE2	4 port 2-pos int vol valve m-elect/2in so	1/8 x .75mm 175C/1000psi liq 2ul N60/E
ED410UWE	10 port 2-pos valve m-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
ED410UWT	10 port 2-pos valve m-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
ED44UWE	4 port 2-pos valve m-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
ED44UWT	4 port 2-pos valve m-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
ED46UWE	6 port 2-pos valve m-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
ED46UWT	6 port 2-pos valve m-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
ED48UWE	8 port 2-pos valve m-electric/4in so	1/8 x .75mm 225C/400psi gas N60/E
ED48UWT	8 port 2-pos valve m-electric/4in so	1/8 x .75mm 330C/300psi gas N60/T
ED4C10UWE	10 port 2-pos valve m-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
ED4C10UWT	10 port 2-pos valve m-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
ED4C4UWE	4 port 2-pos valve m-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
ED4C4UWT	4 port 2-pos valve m-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
ED4C6UWE	6 port 2-pos valve m-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E
ED4C6UWT	6 port 2-pos valve m-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
ED4C8UWE	8 port 2-pos valve m-electric/4in so	1/16 x .75mm 225C/400psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




ED4C8UWT	8 port 2-pos valve m-electric/4in so	1/16 x .75mm 330C/300psi gas N60/T
ED4UW	4 port 2-pos valve m-electric	1/8 x .75mm 75C/5000psi liq N60/H
ED6UW	6 port 2-pos valve m-electric	1/8 x .75mm 75C/5000psi liq N60/H
ED8UW	8 port 2-pos valve m-electric	1/8 x .75mm 75C/5000psi liq N60/H
EDC10UW	10 port 2-pos valve m-electric	1/16 x .75mm 75C/5000psi liq N60/H
EDC4UW	4 port 2-pos valve m-electric	1/16 x .75mm 75C/5000psi liq N60/H
EDC6UW	6 port 2-pos valve m-electric	1/16 x .75mm 75C/5000psi liq N60/H
EDC8UW	8 port 2-pos valve m-electric	1/16 x .75mm 75C/5000psi liq N60/H
EDCI4UW.2	4 port 2-pos int vol valve m-electric	1/16 x .75mm 75C/5000psi liq .2ul N60/H
EDCI4UW.5	4 port 2-pos int vol valve m-electric	1/16 x .75mm 75C/5000psi liq .5ul N60/H
EDCI4UW1	4 port 2-pos int vol valve m-electric	1/16 x .75mm 75C/5000psi liq 1ul N60/H
EDCI4UW2	4 port 2-pos int vol valve m-electric	1/16 x .75mm 75C/5000psi liq 2ul N60/H
EDL10UW	10 port 2-pos valve m-electric	1/8 x 1.0mm 75C/5000psi liq N60/H
EDL4UW	4 port 2-pos valve m-electric	1/8 x 1.7mm 75C/5000psi liq N60/H
EDL6UW	6 port 2-pos valve m-electric	1/8 x 1.7mm 75C/5000psi liq N60/H
EDL8UW	8 port 2-pos valve m-electric	1/8 x 1.3mm 75C/5000psi liq N60/H
EDS	Micro electric actuator 2-pos for so	high torque 24Vdc w.110/220Vac p-supply
EH	Micro electric actuator 2-pos closelount	high speed 24Vdc w.110/220Vac p-supply
EH2	Micro electric actuator 2-pos/2in so	high speed 24Vdc w.110/220Vac p-supply
EH4C10WE	10 port 2-pos valve m-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
EH4C10WT	10 port 2-pos valve m-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
EH4C4WE	4 port 2-pos valve m-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
EH4C4WT	4 port 2-pos valve m-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
EH4C6WE	6 port 2-pos valve m-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
EH4C6WT	6 port 2-pos valve m-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
EH4C8WE	8 port 2-pos valve m-electric/4in so	1/16 x .4mm 225C/400psi gas N60/E
EH4C8WT	8 port 2-pos valve m-electric/4in so	1/16 x .4mm 350C/300psi gas N60/T
EH4N10WE	10 port 2-pos valve m-electric act/4in so	1/32 x .25mm 225C/400psi gas N60/E
EH4N10WT	10 port 2-pos valve m-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
EH4N4WE	4 port 2-pos valve m-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
EH4N4WT	4 port 2-pos valve m-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
EH4N6WE	6 port 2-pos valve m-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
EH4N6WT	6 port 2-pos valve m-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
EH4N8WE	8 port 2-pos valve m-electric/4in so	1/32 x .25mm 225C/400psi gas N60/E
EH4N8WT	8 port 2-pos valve m-electric/4in so	1/32 x .25mm 350C/300psi gas N60/T
EHS	Micro electric actuator 2-pos for so	high speed 24Vdc w.110/220Vac p-supply
EL10UW	10 port 2-pos valve std-electric	1/8 x 1.0mm 75C/5000psi liq N60/H
EL4UW	4 port 2-pos valve std-electric	1/8 x 1.7mm 75C/5000psi liq N60/H
EL6UW	6 port 2-pos valve std-electric	1/8 x 1.7mm 75C/5000psi liq N60/H
EL8UW	8 port 2-pos valve std-electric	1/8 x 1.3mm 75C/5000psi liq N60/H
EMH	!!!RECOMMENDED EMH-CE	
EMH2	!!!RECOMMENDED EMH2-CE	
EMH2-CE	Micro electric actuator mpos 2in so	high speed 24Vdc w. 110/220Vac p-supply
EMH-CE	Micro electric actuator mpos closelount	high speed 24Vdc w. 110/220Vac p-supply
EMHS	!!!RECOMMENDED EMHS-CE	
EMHS-CE	Micro electric actuator mpos for so	high speed 24Vdc w. 110/220Vac p-supply
EMT	!!!RECOMMENDED EMT-CE	
EMT2	!!!RECOMMENDED EMT2-CE	
EMT2-CE	Micro electric actuator mpos 2in so	high torque 24Vdc w.110/220Vac p-supply
EMT2CSC10MWE	10 mpos common outlet path MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSC12MWE	12 mpos common outlet path MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSC16MWE	16 mpos common outlet path MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSC6MWE	6 mpos common outlet path MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSD10MWE	10 mpos dead-end path MW m-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
EMT2CSD12MWE	12 mpos dead-end path MW m-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
EMT2CSD16MWE	16 mpos dead-end path MW m-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
EMT2CSD6MWE	6 mpos dead-end path MW m-elect/2in so	1/16 x .75mm 200C/400psi gas N60/E
EMT2CSF10MWE	10 mpos flow-through MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSF12MWE	12 mpos flow-through MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSF16MWE	16 mpos flow-through MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E
EMT2CSF6MWE	6 mpos flow-through MW m-elect/2in so	1/16 x 1mm 200C/200psi gas N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




EMT2CST10MWE	10 mpos trapping path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CST12MWE	12 mpos trapping path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CST16MWE	16 mpos trapping path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CST6MWE	6 mpos trapping path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CSTF10MWE	10 mpos trap/flow path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CSTF12MWE	12 mpos trap/flow path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CSTF16MWE	16 mpos trap/flow path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2CSTF6MWE	6 mpos trap/flow path MW m-elect/2in so	1/16 x .75mm 200C/200psi gas N60/E
EMT2SC10MWE	10 mpos common outlet path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SC12MWE	12 mpos common outlet path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SC16MWE	16 mpos common outlet path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SC6MWE	6 mpos common outlet path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SD10MWE	10 mpos dead-end path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SD12MWE	12 mpos dead-end path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SD16MWE	16 mpos dead-end path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SD6MWE	6 mpos dead-end path MW m-elect/2in so	1/8 x 1mm 200C/400psi gas N60/E
EMT2SF10MWE	10 mpos flow-through MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SF12MWE	12 mpos flow-through MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SF16MWE	16 mpos flow-through MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2SF6MWE	6 mpos flow-through MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2ST10MWE	10 mpos trapping path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2ST12MWE	12 mpos trapping path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2ST16MWE	16 mpos trapping path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2ST6MWE	6 mpos trapping path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2STF10MWE	10 mpos trap/flow path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2STF12MWE	12 mpos trap/flow path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2STF16MWE	16 mpos trap/flow path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2STF6MWE	6 mpos trap/flow path MW m-elect/2in so	1/8 x 1mm 200C/200psi gas N60/E
EMT2VLS4MWE2	4 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS6MWE2	6 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS8MWE2	8 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS10MWE2	10 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS12MWE2	12 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS16MWE2	16 mpos common outlet path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS6MWE2	6 mpos dead-end path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS8MWE2	8 mpos dead-end path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS10MWE2	10 mpos dead-end path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS12MWE2	12 mpos dead-end path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLS16MWE2	16 mpos dead-end path MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLSF4MWE2	4 mpos flow-through MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLSF6MWE2	6 mpos flow-through MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT2VLSF8MWE2	8 mpos flow-through MW m-elect/2in so	1/4 x 4mm 75C/100psi gas N60/E2
EMT-CE	Micro electric actuator mpos closemount	high torque 24Vdc w.110/220Vac p-supply
EMTCS10UW	10 mpos dead-end path m-electric	1/16 x .4mm 75C/5000psi liq N60/E
EMTCS4UW	4 mpos dead-end path m-electric	1/16 x .4mm 75C/5000psi liq N60/E
EMTCS6UW	6 mpos dead-end path m-electric	1/16 x .4mm 75C/5000psi liq N60/E
EMTCS4UW	4 mpos trapping path m-electric	1/16 x .4mm 75C/5000psi liq N60/E
EMTCS6UW	6 mpos trapping path m-electric	1/16 x .4mm 75C/5000psi liq N60/E
EMTS	RECOMMENDED EMTS-CE	
EMTS-CE	Micro electric actuator mpos for so	high torque 24Vdc w.110/220Vac p-supply
EMTSD4UW	4 mpos dead-end path m-electric	1/8 x .75mm 75C/5000psi liq N60/E
EMTSD6UW	6 mpos dead-end path m-electric	1/8 x .75mm 75C/5000psi liq N60/E
EMTSD8UW	8 mpos dead-end path m-electric	1/8 x .75mm 75C/5000psi liq N60/E
EN.5	External nut 1/32in hex	SS
EN.5KN	External nut 1/32in knurled	SS
EN1	External nut 1/16in hex	SS
EN1K	External nut 1in hex	SS
EN2	External nut 1/8in hex	SS
EN2SI	Septum injector nut 1/8in connection	with polyimide support and 1/4in septum
EN4	External nut 1/4in hex	SS
EN6	External nut 3/8in hex	SS
EN8	External nut 1/2in hex	SS
EP	Micro electric actuator 2-pos closemount	medium torque 24Vdc w.110/220Vac p-supply
EP2	Micro electric actuator 2-pos/2in so	medium torque 24Vdc w.110/220Vac p-supply
EP2CI4WE.06	4 port 2-pos int vol valve m-elect/2in so	1/16 x .4mm 175C/1000psi liq .06ul N60/E

# VICI Valco-Cheminert PRODUCT LIST 2012




EP2CI4WE.1	4 port 2-pos int vol valve m-elect/2in so	1/16 x .4mm 175C/1000psi liq .1ul N60/E
EP2CI4WE.2	4 port 2-pos int vol valve m-elect/2in so	1/16 x .4mm 175C/1000psi liq .2ul N60/E
EP2CI4WE.5	4 port 2-pos int vol valve m-elect/2in so	1/16 x .4mm 175C/1000psi liq .5ul N60/E
EP2NI4WE.06	4 port 2-pos int vol valve m-elect/2in so	1/32 x .25mm 175C/1000psi liq .06ulN60/E
EP2NI4WE.1	4 port 2-pos int vol valve m-elect/2in so	1/32 x .25mm 175C/1000psi liq .1ulN60/E
EP2NI4WE.2	4 port 2-pos int vol valve m-elect/2in so	1/32 x .25mm 175C/1000psi liq .2ulN60/E
EP2NI4WE.5	4 port 2-pos int vol valve m-elect/2in so	1/32 x .25mm 175C/1000psi liq .5ulN60/E
EPC10W	10 port 2-pos valve m-electric hi-psi	1/16 x .4mm 75C/5000psi liq N60/H
EPC4W	4 port 2-pos valve m-electric hi-psi	1/16 x .4mm 75C/5000psi liq N60/H
EPC6W	6 port 2-pos valve m-electric hi-psi	1/16 x .4mm 75C/5000psi liq N60/H
EPC8W	8 port 2-pos valve m-electric hi-psi	1/16 x .4mm 75C/5000psi liq N60/H
EPCI4W.06	4 port 2-pos int vol valve m-elect hi-psi	1/16 x .4mm 75C/5000psi liq .06ul N60/H
EPCI4W.1	4 port 2-pos int vol valve m-elect hi-psi	1/16 x .4mm 75C/5000psi liq .1ul N60/H
EPCI4W.2	4 port 2-pos int vol valve m-elect hi-psi	1/16 x .4mm 75C/5000psi liq .2ul N60/H
EPCI4W.5	4 port 2-pos int vol valve m-elect hi-psi	1/16 x .4mm 75C/5000psi liq .5ul N60/H
EPS	Micro electric actuator 2-pos for so	medium torque 24Vdc w.110/220Vac p-supply
EQ	Micro electric actuator 2-pos closemount	highest speed 24Vdc w.110/220Vac p-supply
EQ2	Micro electric actuator 2-pos/2in so	highest speed 24Vdc w.110/220Vac p-supply
EQS	Micro electric actuator 2-pos for so	highest speed 24Vdc w.110/220Vac p-supply
ERU21	External reducing union 1/8in x 1/16in	.75mm bore SS
ERU21L	External reducing union 1/8in x 1/16in	1mm bore SS
ERU21T	External reducing union 1/8in x 1/16in	1/16in through bore SS
ERU41	External reducing union 1/4in x 1/16in	.75mm bore SS
ERU41T	External reducing union 1/4in x 1/16in	1/16in through bore SS
ERU42	External reducing union 1/4in x 1/8in	0.75mm bore SS
ERU42L	External reducing union 1/4in x 1/8in	2mm bore SS
ERU42T	External reducing union 1/4in x 1/8in	1/8in through bore SS
ESD4UW	4 mpos dead-end path std-electric	1/8 x .75mm 75C/5000psi liq N60/E
ESD6UW	6 mpos dead-end path std-electric	1/8 x .75mm 75C/5000psi liq N60/E
ESD8UW	8 mpos dead-end path std-electric	1/8 x .75mm 75C/5000psi liq N60/E
ET	Micro electric actuator 2-pos closemount	highest torque 24Vdc w.110/220Vac p-supply
ET2	Micro electric actuator 2-pos/2in so	highest torque 24Vdc w.110/220Vac p-supply
ET4VL4MWE2	4 port 2-pos valve MW m-electric/4in so	1/4 x 4mm 75C/100psi gas N60/E2
ET4VL6MWE2	6 port 2-pos valve MW m-electric/4in so	1/4 x 4mm 75C/100psi gas N60/E2
ET4VL8MWE2	8 port 2-pos valve MW m-electric/4in so	1/4 x 4mm 75C/100psi gas N60/E2
ETS	Micro electric actuator 2-pos for so	highest torque 24Vdc w.110/220Vac p-supply
EU.5	External union 1/32in	.25mm bore SS
EU.5L	External union 1/32in	.50mm bore SS
EU.5T	External union 1/32in	1/32in through bore SS
EU2	External union 1/8in	1mm bore SS
EU2L	External union 1/8in	2mm bore SS
EU2T	External union 1/8in	1/8in through bore SS
EU4	External union 1/4in	2mm bore SS
EU4L	External union 1/4in	4.6mm bore SS
EU4T	External union 1/4in	1/4in through bore SS
EUD	Universal actuator no interface no clamp ring	medium torque 24Vdc w.110/220Vac p-supply
EUDA	Universal actuator RS-232/485 no clamp ring	medium torque 24Vdc w.110/220Vac p-supply
EUDB	Universal actuator USB no clamp ring	medium torque 24Vdc w.110/220Vac p-supply
EUDC	Universal actuator BCD no clamp ring	medium torque 24Vdc w.110/220Vac p-supply
EUH	Universal actuator no interface no clamp ring	high speed 24Vdc w.110/220Vac p-supply
EUHA	Universal actuator RS-232/485 no clamp ring	high speed 24Vdc w.110/220Vac p-supply
EUHB	Universal actuator USB no clamp ring	high speed 24Vdc w.110/220Vac p-supply
EUHC	Universal actuator BCD no clamp ring	high speed 24Vdc w.110/220Vac p-supply
EUT	Universal actuator no interface no clamp ring	high torque 24Vdc w.110/220Vac p-supply
EUTA	Universal actuator RS-232/485 no clamp ring	high torque 24Vdc w.110/220Vac p-supply
EUTB	Universal actuator USB no clamp ring	high torque 24Vdc w.110/220Vac p-supply
EUTC	Universal actuator BCD no clamp ring	high torque 24Vdc w.110/220Vac p-supply
EZA1	Ext./int. adapter 1/16in x 1/16in SS	
EZA2	Ext./int. adapter 1/8in x 1/8in SS	
EZBRU.51	Bulkhead ext/int red. union 1/16in x 1/32in	.25mm bore SS
EZBRU.51L	Bulkhead ext/int red. union 1/16in x 1/32in	.5mm bore SS

# VICI Valco-Cheminert PRODUCT LIST 2012




EZBRU.51T	Bulkhead ext/int red. union 1/16in x 1/32in	1/32in through bore SS
EZBRU1.5L	Bulkhead ext/int red. union 1/16in x 1/32in	.5mm bore SS
EZBRU1.5T	Bulkhead ext/int red. union 1/16in x 1/32	1/32in through bore SS
EZBRU2.5L	Bulkhead ext/int red. union 1/8in x 1/32in	.5mm bore SS
EZBRU2.5T	Bulkhead ext/int red. union 1/8in x 1/32in	1/32in through bore SS
EZBRU21	Bulkhead ext/int red. union 1/8in x 1/16in	.75mm bore SS
EZBRU21T	Bulkhead ext/int red. union 1/8in x 1/16in	1/16in through bore SS
EZBRU41	Bulkhead ext/int red. union 1/4in x 1/16in	.75mm bore SS
EZBRU41T	Bulkhead ext/int red. union 1/4in x 1/16in	1/16in through bore SS
EZBRU42	Bulkhead ext/int red. union 1/4in x 1/8in	1mm bore SS
EZBRU42L	Bulkhead ext/int red. union 1/4in x 1/8in	2mm bore SS
EZBRU42T	Bulkhead ext/int red. union 1/4in x 1/8in	1/8in through bore SS
EZBU1	Bulkhead ext./int union 1/16in	.75mm bore SS
EZBU1C	Bulkhead ext./int union 1/16in	.25mm bore SS
EZBU1M	Bulkhead ext./int union 1/16in	.5mm bore SS
EZBU1T	Bulkhead ext./int union 1/16in	1/16in through bore SS
EZBU2	Bulkhead ext./int union 1/8in	1mm bore SS
EZBU2L	Bulkhead ext./int union 1/8in	2mm bore SS
EZBU2T	Bulkhead ext./int union 1/8in	1/8in through bore SS
EZR1.5	Ext./int. reducing union 1/16in x 1/32in	.25mm bore SS
EZR1.5T	Ext./int. reducing union 1/16in x 1/32in	1/32in through bore SS
EZR2.5	Ext./int. reducing union 1/8in x 1/32in	.25mm bore SS
EZR21	Ext./int. reducing union 1/8in x 1/16in	.50mm bore SS
EZR21T	Ext./int. reducing union 1/8in x 1/16in	1/16in through bore SS
EZR41	Ext./int. reducing union 1/4in x 1/16in	1mm bore SS
EZR41T	Ext./int. reducing union 1/4in x 1/16in	1/16in through bore SS
EZR42	Ext./int. reducing union 1/4in x 1/8in	1mm bore SS
EZR42T	Ext./int. reducing union 1/4in x 1/8in	1/8in through bore SS
EZRF1.5KF-5	External reducing ferrule 1/16in x 1/32in	CTFE 5/pkg
EZRF1.5PK-5	External reducing ferrule 1/16in x 1/32in	PEEK 5/pkg
EZRF1.5TFG-5	External reducing ferrule 1/16in x 1/32in	PTFE Glass-filled 25% 5/pkg
EZRF1.5V-5	External reducing ferrule 1/16in x 1/32in	Polyimide Valcon 5/pkg
EZRF2.5KF-5	External reducing ferrule 1/8in x 1/32in	CTFE 5/pkg
EZRF2.5PK-5	External reducing ferrule 1/8in x 1/32in	PEEK 5/pkg
EZRF2.5TFG-5	External reducing ferrule 1/8in x 1/32in	PTFE Glass-filled 25% 5/pkg
EZRF2.5V-5	External reducing ferrule 1/8in x 1/32in	Polyimide Valcon 5/pkg
EZRF21KF-5	External reducing ferrule 1/8in x 1/16in	CTFE 5/pkg
EZRF21PK-5	External reducing ferrule 1/8in x 1/16in	PEEK 5/pkg
EZRF21TFG-5	External reducing ferrule 1/8in x 1/16in	PTFE Glass-filled 25% 5/pkg
EZRF21V-5	External reducing ferrule 1/8in x 1/16in	Polyimide Valcon 5/pkg
EZRF41KF-5	External reducing ferrule 1/4in x 1/16in	CTFE 5/pkg
EZRF41PK-5	External reducing ferrule 1/4in x 1/16in	PEEK 5/pkg
EZRF41TFG-5	External reducing ferrule 1/4in x 1/16in	PTFE Glass-filled 25% 5/pkg
EZRF41V-5	External reducing ferrule 1/4in x 1/16in	Poyimide Valcon 5/pkg
EZRF42KF-5	External reducing ferrule 1/4in x 1/8in	CTFE 5/pkg
EZRF42PK-5	External reducing ferrule 1/4in x 1/8in	PEEK 5/pkg
EZRF42TFG-5	External reducing ferrule 1/4in x 1/8in	PTFE Glass-filled 25% 5/pkg
EZRF42V-5	External reducing ferrule 1/4in x 1/8in	Polyimide Valcon 5/pkg
EZRU.51	Int./ext. reducing union 1/16in x 1/32in	.25mm bore SS
EZRU.51L	Int./ext. reducing union 1/16in x 1/32in	.5mm bore SS
EZRU.51T	Int./ext. reducing union 1/16in x 1/32in	1/32in through bore SS
EZRU1.5	Ext./int. reducing union 1/16in x 1/32in	.25mm bore SS
EZRU1.5L	Ext./int. reducing union 1/16in x 1/32in	.5mm bore SS
EZRU1.5T	Ext./int. reducing union 1/16in x 1/32in	1/32in through bore SS
EZRU2.5	Ext./int. reducing union 1/8in x 1/32in	.25mm bore SS
EZRU2.5L	Ext./int. reducing union 1/8in x 1/32in	.5mm bore SS
EZRU2.5T	Ext./int. reducing union 1/8in x 1/32in	1/32in through bore SS
EZRU21	Ext./int. reducing union 1/8in x 1/16in	.75mm bore SS
EZRU21C	Ext./int. reducing union 1/8in x 1/16in	.25mm bore SS
EZRU21T	Ext./int. reducing union 1/8in x 1/16in	1/16in through bore SS
EZRU41	Ext./Int. reducing unions 1/4in x 1/16in	.75mm bore SS

# VICI Valco-Cheminert PRODUCT LIST 2012




EZRU41C	Ext./Int. reducing unions 1/4in x 1/16in	.25mm bore SS
EZRU41T	Ext./Int. reducing unions 1/4in x 1/16in	1/16in through bore SS
EZRU42	Ext./int. reducing union 1/4in x 1/8in	1mm bore SS
EZRU42L	Ext./int. reducing union 1/4in x 1/8in	2mm bore SS
EZRU42T	Ext./int. reducing union 1/4in x 1/8in	1/8in through bore SS
EZU.5	Ext./int union 1/32in	.25mm bore SS
EZU.5L	Ext./int union 1/32in	.5mm bore SS
EZU1	Ext./int union 1/16in	.75mm bore SS
EZU1C	Ext./int union 1/16in	.25mm bore SS
EZU1M	Ext./int union 1/16in	.5mm bore SS
EZU1T	Ext./int union 1/16in	1/16in through bore SS
EZU2	Ext./int. union 1/8in	1mm bore SS
EZU2L	Ext./int. union 1/8in	2mm bore SS
EZU2T	Ext./int. union 1/8in	1/8in through bore SS
FC10AS1K	Flow controller M100 standard knob	1/8in ext 150ml N2 40psi Aluminum/SS
FC10AS1S3	Flow controller M100 Spectrol 3-dial	1/8in ext 150ml N2 40psi Aluminum/SS
FC10AS1S4	Flow controller M100 Spectrol 4-dial	1/8in ext 150ml N2 40psi Aluminum/SS
FC10AS2K	Flow controller M100 standard knob	1/8in ext 250ml N2 40psi Aluminum/SS
FC10AS2S3	Flow controller M100 Spectrol 3-dial	1/8in ext 250ml N2 40psi Aluminum/SS
FC10AS2S4	Flow controller M100 Spectrol 4-dial	1/8in ext 250ml N2 40psi Aluminum/SS
FC10AS3K	Flow controller M100 standard knob	1/8in ext 850ml N2 40psi Aluminum/SS
FC10AS3S3	Flow controller M100 Spectrol 3-dial	1/8in ext 850ml N2 40psi Aluminum/SS
FC10AS3S4	Flow controller M100 Spectrol 4-dial	1/8in ext 850ml N2 40psi Aluminum/SS
FC10AS4K	Flow controller M100 standard knob	1/8in ext 1.2l N2 40psi Aluminum/SS
FC10AS4S3	Flow controller M100 Spectrol 3-dial	1/8in ext 1.2l N2 40psi Aluminum/SS
FC10AS4S4	Flow controller M100 Spectrol 4-dial	1/8in ext 1.2l N2 40psi Aluminum/SS
FC10AS5K	Flow controller M100 standard knob	1/8in ext 4.5l N2 40psi Aluminum/SS
FC10AS5S3	Flow controller M100 Spectrol 3-dial	1/8in ext 4.5l N2 40psi Aluminum/SS
FC10AS5S4	Flow controller M100 Spectrol 4-dial	1/8in ext 4.5l N2 40psi Aluminum/SS
FC10AS6K	Flow controller M100 standard knob	1/8in ext 10.0l N2 40psi Aluminum/SS
FC10AS6S3	Flow controller M100 Spectrol 3-dial	1/8in ext 10.0l N2 40psi Aluminum/SS
FC10AS6S4	Flow controller M100 Spectrol 4-dial	1/8in ext 10.0l N2 40psi Aluminum/SS
FC10AV1K	Flow controller M100 standard knob	1/8in ext 150ml N2 40psi Aluminum/Viton
FC10AV1S3	Flow controller M100 Spectrol 3-dial	1/8in ext 150ml N2 40psi Aluminum/Viton
FC10AV1S4	Flow controller M100 Spectrol 4-dial	1/8in ext 150ml N2 40psi Aluminum/Viton
FC10AV2K	Flow controller M100 standard knob	1/8in ext 250ml N2 40psi Aluminum/Viton
FC10AV2S3	Flow controller M100 Spectrol 3-dial	1/8in ext 250ml N2 40psi Aluminum/Viton
FC10AV2S4	Flow controller M100 Spectrol 4-dial	1/8in ext 250ml N2 40psi Aluminum/Viton
FC10AV3K	Flow controller M100 standard knob	1/8in ext 850ml N2 40psi Aluminum/Viton
FC10AV3S3	Flow controller M100 Spectrol 3-dial	1/8in ext 850ml N2 40psi Aluminum/Viton
FC10AV3S4	Flow controller M100 Spectrol 4-dial	1/8in ext 850ml N2 40psi Aluminum/Viton
FC10AV4K	Flow controller M100 standard knob	1/8in ext 1.2l N2 40psi Aluminum/Viton
FC10AV4S3	Flow controller M100 Spectrol 3-dial	1/8in ext 1.2l N2 40psi Aluminum/Viton
FC10AV4S4	Flow controller M100 Spectrol 4-dial	1/8in ext 1.2l N2 40psi Aluminum/Viton
FC10AV5K	Flow controller M100 standard knob	1/8in ext 4.5l N2 40psi Aluminum/Viton
FC10AV5S3	Flow controller M100 Spectrol 3-dial	1/8in ext 4.5l N2 40psi Aluminum/Viton
FC10AV5S4	Flow controller M100 Spectrol 4-dial	1/8in ext 4.5l N2 40psi Aluminum/Viton
FC10AV6K	Flow controller M100 standard knob	1/8in ext 10.0l N2 40psi Aluminum/Viton
FC10AV6S3	Flow controller M100 Spectrol 3-dial	1/8in ext 10.0l N2 40psi Aluminum/Viton
FC10AV6S4	Flow controller M100 Spectrol 4-dial	1/8in ext 10.0l N2 40psi Aluminum/Viton
FC10SS1K	Flow controller M100 standard knob	1/8in ext 150ml N2 40psi SS/SS
FC10SS1S3	Flow controller M100 Spectrol 3-dial	1/8in ext 150ml N2 40psi SS/SS
FC10SS1S4	Flow controller M100 Spectrol 4-dial	1/8in ext 150ml N2 40psi SS/SS
FC10SS2K	Flow controller M100 standard knob	1/8in ext 250ml N2 40psi SS/SS
FC10SS2S3	Flow controller M100 Spectrol 3-dial	1/8in ext 250ml N2 40psi SS/SS
FC10SS2S4	Flow controller M100 Spectrol 4-dial	1/8in ext 250ml N2 40psi SS/SS
FC10SS3K	Flow controller M100 standard knob	1/8in ext 850ml N2 40psi SS/SS
FC10SS3S3	Flow controller M100 Spectrol 3-dial	1/8in ext 850ml N2 40psi SS/SS
FC10SS3S4	Flow controller M100 Spectrol 4-dial	1/8in ext 850ml N2 40psi SS/SS
FC10SS4K	Flow controller M100 standard knob	1/8in ext 1.2l N2 40psi SS/SS
FC10SS4S3	Flow controller M100 Spectrol 3-dial	1/8in ext 1.2l N2 40psi SS/SS

# VICI Valco-Cheminert PRODUCT LIST 2012




FC10SS4S4	Flow controller M100 Spectrol 4-dial	1/8in ext 1.2l N2 40psi SS/SS
FC10SS5K	Flow controller M100 standard knob	1/8in ext 4.5l N2 40psi SS/SS
FC10SS5S3	Flow controller M100 Spectrol 3-dial	1/8in ext 4.5l N2 40psi SS/SS
FC10SS5S4	Flow controller M100 Spectrol 4-dial	1/8in ext 4.5l N2 40psi SS/SS
FC10SS6K	Flow controller M100 standard knob	1/8in ext 10.0l N2 40psi SS/SS
FC10SS6S3	Flow controller M100 Spectrol 3-dial	1/8in ext 10.0l N2 40psi SS/SS
FC10SS6S4	Flow controller M100 Spectrol 4-dial	1/8in ext 10.0l N2 40psi SS/SS
FC10SV1K	Flow controller M100 standard knob	1/8in ext 150ml N2 40psi SS/Viton
FC10SV1S3	Flow controller M100 Spectrol 3-dial	1/8in ext 150ml N2 40psi SS/Viton
FC10SV1S4	Flow controller M100 Spectrol 4-dial	1/8in ext 150ml N2 40psi SS/Viton
FC10SV2K	Flow controller M100 standard knob	1/8in ext 250ml N2 40psi SS/Viton
FC10SV2S3	Flow controller M100 Spectrol 3-dial	1/8in ext 250ml N2 40psi SS/Viton
FC10SV2S4	Flow controller M100 Spectrol 4-dial	1/8in ext 250ml N2 40psi SS/Viton
FC10SV3K	Flow controller M100 standard knob	1/8in ext 850ml N2 40psi SS/Viton
FC10SV3S3	Flow controller M100 Spectrol 3-dial	1/8in ext 850ml N2 40psi SS/Viton
FC10SV3S4	Flow controller M100 Spectrol 4-dial	1/8in ext 850ml N2 40psi SS/Viton
FC10SV4K	Flow controller M100 standard knob	1/8in ext 1.2l N2 40psi SS/Viton
FC10SV4S3	Flow controller M100 Spectrol 3-dial	1/8in ext 1.2l N2 40psi SS/Viton
FC10SV4S4	Flow controller M100 Spectrol 4-dial	1/8in ext 1.2l N2 40psi SS/Viton
FC10SV5K	Flow controller M100 standard knob	1/8in ext 4.5l N2 40psi SS/Viton
FC10SV5S3	Flow controller M100 Spectrol 3-dial	1/8in ext 4.5l N2 40psi SS/Viton
FC10SV5S4	Flow controller M100 Spectrol 4-dial	1/8in ext 4.5l N2 40psi SS/Viton
FC10SV6K	Flow controller M100 standard knob	1/8in ext 10.0l N2 40psi SS/Viton
FC10SV6S3	Flow controller M100 Spectrol 3-dial	1/8in ext 10.0l N2 40psi SS/Viton
FC10SV6S4	Flow controller M100 Spectrol 4-dial	1/8in ext 10.0l N2 40psi SS/Viton
FC22AS1K	Flow controller M202 standard knob	1/8in NPT adjustable span Aluminum/SS
FC22AS1S3	Flow controller M202 Spectrol 3-dial	1/8in NPT adjustable span Aluminum/SS
FC22AS1S4	Flow controller M202 Spectrol 4-dial	1/8in NPT adjustable span Aluminum/SS
FC22AV1K	Flow controller M202 standard knob	1/8in NPT adjustable span Aluminum/Viton
FC22AV1S3	Flow controller M202 Spectrol 3-dial	1/8in NPT adjustable span Aluminum/Viton
FC22AV1S4	Flow controller M202 Spectrol 4-dial	1/8in NPT adjustable span Aluminum/Viton
FC22SS1K	Flow controller M202 standard knob	1/8in NPT adjustable span SS/SS
FC22SS1S3	Flow controller M202 Spectrol 3-dial	1/8in NPT adjustable span SS/SS
FC22SS1S4	Flow controller M202 Spectrol 4-dial	1/8in NPT adjustable span SS/SS
FC22SV1K	Flow controller M202 standard knob	1/8in NPT adjustable span SS/Viton
FC22SV1S3	Flow controller M202 Spectrol 3-dial	1/8in NPT adjustable span SS/Viton
FC22SV1S4	Flow controller M202 Spectrol 4-dial	1/8in NPT adjustable span SS/Viton
FC30AS1K	Flow controller M300 standard knob	1/8in ext 200ml N2 40psi Aluminum/SS
FC30AS1S3	Flow controller M300 Spectrol 3-dial	1/8in ext 200ml N2 40psi Aluminum/SS
FC30AS1S4	Flow controller M300 Spectrol 4-dial	1/8in ext 200ml N2 40psi Aluminum/SS
FC30AS2K	Flow controller M300 standard knob	1/8in ext 300ml N2 40psi Aluminum/SS
FC30AS2S3	Flow controller M300 Spectrol 3-dial	1/8in ext 300ml N2 40psi Aluminum/SS
FC30AS2S4	Flow controller M300 Spectrol 4-dial	1/8in ext 300ml N2 40psi Aluminum/SS
FC30AS3K	Flow controller M300 standard knob	1/8in ext 800ml N2 40psi Aluminum/SS
FC30AS3S3	Flow controller M300 Spectrol 3-dial	1/8in ext 800ml N2 40psi Aluminum/SS
FC30AS3S4	Flow controller M300 Spectrol 4-dial	1/8in ext 800ml N2 40psi Aluminum/SS
FC30AS4K	Flow controller M300 standard knob	1/8in ext 1.6l N2 40psi Aluminum/SS
FC30AS4S3	Flow controller M300 Spectrol 3-dial	1/8in ext 1.6l N2 40psi Aluminum/SS
FC30AS4S4	Flow controller M300 Spectrol 4-dial	1/8in ext 1.6l N2 40psi Aluminum/SS
FC30AV1K	Flow controller M300 standard knob	1/8in ext 200ml N2 40psi Aluminum/Viton
FC30AV1S3	Flow controller M300 Spectrol 3-dial	1/8in ext 200ml N2 40psi Aluminum/Viton
FC30AV1S4	Flow controller M300 Spectrol 4-dial	1/8in ext 200ml N2 40psi Aluminum/Viton
FC30AV2K	Flow controller M300 standard knob	1/8in ext 300ml N2 40psi Aluminum/Viton
FC30AV2S3	Flow controller M300 Spectrol 3-dial	1/8in ext 300ml N2 40psi Aluminum/Viton
FC30AV2S4	Flow controller M300 Spectrol 4-dial	1/8in ext 300ml N2 40psi Aluminum/Viton
FC30AV3K	Flow controller M300 standard knob	1/8in ext 800ml N2 40psi Aluminum/Viton
FC30AV3S3	Flow controller M300 Spectrol 3-dial	1/8in ext 800ml N2 40psi Aluminum/Viton
FC30AV3S4	Flow controller M300 Spectrol 4-dial	1/8in ext 800ml N2 40psi Aluminum/Viton
FC30AV4K	Flow controller M300 standard knob	1/8in ext 1.6l N2 40psi Aluminum/Viton
FC30AV4S3	Flow controller M300 Spectrol 3-dial	1/8in ext 1.6l N2 40psi Aluminum/Viton
FC30AV4S4	Flow controller M300 Spectrol 4-dial	1/8in ext 1.6l N2 40psi Aluminum/Viton

# VICI Valco-Cheminert PRODUCT LIST 2012




FC30SS1K	Flow controller M300 standard knob	1/8in ext 200ml N2 40psi SS/SS
FC30SS1S3	Flow controller M300 Spectrol 3-dial	1/8in ext 200ml N2 40psi SS/SS
FC30SS1S4	Flow controller M300 Spectrol 4-dial	1/8in ext 200ml N2 40psi SS/SS
FC30SS2K	Flow controller M300 standard knob	1/8in ext 300ml N2 40psi SS/SS
FC30SS2S3	Flow controller M300 Spectrol 3-dial	1/8in ext 300ml N2 40psi SS/SS
FC30SS2S4	Flow controller M300 Spectrol 4-dial	1/8in ext 300ml N2 40psi SS/SS
FC30SS3K	Flow controller M300 standard knob	1/8in ext 800ml N2 40psi SS/SS
FC30SS3S3	Flow controller M300 Spectrol 3-dial	1/8in ext 800ml N2 40psi SS/SS
FC30SS3S4	Flow controller M300 Spectrol 4-dial	1/8in ext 800ml N2 40psi SS/SS
FC30SS4K	Flow controller M300 standard knob	1/8in ext 1.6l N2 40psi SS/SS
FC30SS4S3	Flow controller M300 Spectrol 3-dial	1/8in ext 1.6l N2 40psi SS/SS
FC30SS4S4	Flow controller M300 Spectrol 4-dial	1/8in ext 1.6l N2 40psi SS/SS
FC30SV1K	Flow controller M300 standard knob	1/8in ext 200ml N2 40psi SS/Viton
FC30SV1S3	Flow controller M300 Spectrol 3-dial	1/8in ext 200ml N2 40psi SS/Viton
FC30SV1S4	Flow controller M300 Spectrol 4-dial	1/8in ext 200ml N2 40psi SS/Viton
FC30SV2K	Flow controller M300 standard knob	1/8in ext 300ml N2 40psi SS/Viton
FC30SV2S3	Flow controller M300 Spectrol 3-dial	1/8in ext 300ml N2 40psi SS/Viton
FC30SV2S4	Flow controller M300 Spectrol 4-dial	1/8in ext 300ml N2 40psi SS/Viton
FC30SV3K	Flow controller M300 standard knob	1/8in ext 800ml N2 40psi SS/Viton
FC30SV3S3	Flow controller M300 Spectrol 3-dial	1/8in ext 800ml N2 40psi SS/Viton
FC30SV3S4	Flow controller M300 Spectrol 4-dial	1/8in ext 800ml N2 40psi SS/Viton
FC30SV4K	Flow controller M300 standard knob	1/8in ext 1.6l N2 40psi SS/Viton
FC30SV4S3	Flow controller M300 Spectrol 3-dial	1/8in ext 1.6l N2 40psi SS/Viton
FC30SV4S4	Flow controller M300 Spectrol 4-dial	1/8in ext 1.6l N2 40psi SS/Viton
FC31AV1	Flow controller M300 screwdriver adjustment	1/8in ext 750ml N2 40psi Aluminium/Viton
FP1Z3M21	Manifold pipe adaptor	one 1/8in F-NPT/three 1/16in 1mm bore SS
FP1Z3M22	Manifold pipe adaptor	one 1/8in F-NPT/three 1/8in 2mm bore SS
FP1Z3M24	Manifold pipe adaptor	one 1/8in F-NPT/three 1/4in 4.6mm bore SS
FP1Z3M41	Manifold pipe adaptor	one 1/4in F-NPT/three 1/16in 1mm bore SS
FP1Z3M42	Manifold pipe adaptor	one 1/4in F-NPT/three 1/8in 2mm bore SS
FP1Z3M44	Manifold pipe adaptor	one 1/4in F-NPT/three 1/4in 4.6mm bore SS
FP2Z3M21	Manifold pipe adaptor	two 1/8in F-NPT/three 1/16in 1mm bore SS
FP2Z3M22	Manifold pipe adaptor	two 1/8in F-NPT/three 1/8in 2mm bore SS
FP2Z3M24	Manifold pipe adaptor	two 1/8in F-NPT/three 1/4in 4.6mm bore SS
FP2Z3M41	Manifold pipe adaptor	two 1/4in F-NPT/three 1/16in 1mm bore SS
FP2Z3M42	Manifold pipe adaptor	two 1/4in F-NPT/three 1/8in 2mm bore SS
FP2Z3M44	Manifold pipe adaptor	two 1/4in F-NPT/three 1/4in 4.6mm bore SS
FPEA22	Female pipe adaptor 1/8in F-NPT x 1/8in ext	2mm bore SS
FPEA24	Female pipe adaptor 1/8in F-NPT x 1/4in ext	4.6mm bore SS
FPEA42	Female pipe adaptor 1/4in F-NPT x 1/8in ext	2mm bore SS
FPEA44	Female pipe adaptor 1/4in F-NPT x 1/4in ext	4.6mm bore SS
FPEA82	Female pipe adaptor 1/2in F-NPT x 1/8in ext	2mm bore SS
FPEA84	Female pipe adaptor 1/2in F-NPT x 1/4in ext	4.6mm bore SS
FPZA21	Female pipe adaptor 1/8in F-NPT x 1/16in int	1mm bore SS
FPZA22	Female pipe adaptor 1/8in F-NPT x 1/8in int	1mm bore SS
FPZA22L	Female pipe adaptor 1/8in F-NPT x 1/8in int	2mm bore SS
FPZA41	Female pipe adaptor 1/4in F-NPT x 1/16in int	1mm bore SS
FPZA42	Female pipe adaptor 1/4in F-NPT x 1/8in int	1mm bore SS
FPZA42L	Female pipe adaptor 1/4in F-NPT x 1/8in int	2mm bore SS
FPZA44L	Female pipe adaptor 1/4in F-NPT x 1/4in int	4.6mm bore SS
FPZA81	Female pipe adaptor 1/2in F-NPT x 1/16in int	1mm bore SS
FPZA82	Female pipe adaptor 1/2in F-NPT x 1/8in int	1mm bore SS
FPZA82L	Female pipe adaptor 1/2in F-NPT x 1/8in int	2mm bore SS
FPZA84L	Female pipe adaptor 1/2in F-NPT x 1/4in int	4.6mm bore SS
FRK1	Ferrule removal kit	
FS.25-5	Fused silica adaptor 1/32in	for tubing .25mm OD Valcon T 5/pkg
FS.25PK-5	Fused silica adaptor 1/32in	for tubing .25mm OD PEEK 5/pkg
FS.25V1-5	Fused silica adaptor 1/32in	for tubing .25mm OD Polyimide Virgin 5/pkg
FS.4-5	Fused silica adaptor 1/32in	for tubing .4mm OD Valcon T 5/pkg
FS.4PK-5	Fused silica adaptor 1/32in	for tubing .4mm OD PEEK 5/pkg
FS.4V1-5	Fused silica adaptor 1/32in	for tubing .4mm OD Polyimide Virgin 5/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




FS.5-5	Fused silica adaptor 1/32in	for tubing .5mm OD Valcon T 5/pkg
FS.5PK-5	Fused silica adaptor 1/32in	for tubing .5mm OD PEEK 5/pkg
FS.5V1-5	Fused silica adaptor 1/32in	for tubing .5mm OD Polyimide Virgin 5/pkg
FS1.2-5	Fused silica adaptor 1/16in	for tubing .2mm OD Valcon T 5/pkg
FS1.25-5	Fused silica adaptor 1/16in	for tubing .25mm OD Valcon T 5/pkg
FS1.25PK-5	Fused silica adaptor 1/16in	for tubing .25mm OD PEEK 5/pkg
FS1.25V1-5	Fused silica adaptor 1/16in	for tubing .25mm OD Polyimide Virgin 5/pkg
FS1.2PK-5	Fused silica adaptor 1/16in	for tubing .2mm OD PEEK 5/pkg
FS1.2V1-5	Fused silica adaptor 1/16in	for tubing .2mm OD Polyimide Virgin 5/pkg
FS1.3-5	Fused silica adaptor 1/16in	for tubing .3mm OD Valcon T 5/pkg
FS1.35-R10	Ferrule for fused silica union	for .25mm ID tubing Valcon T 10/pkg
FS1.3PK-5	Fused silica adaptor 1/16in	for tubing .3mm OD PEEK 5/pkg
FS1.3V1-5	Fused silica adaptor 1/16in	for tubing .3mm OD Polyimide Virgin 5/pkg
FS1.4-5	Fused silica adaptor 1/16in	for tubing .4mm OD Valcon T 5/pkg
FS1.45-R10	Ferrule for fused silica union	for .32mm ID tubing Valcon T 10/pkg
FS1.4PK-5	Fused silica adaptor 1/16in	for tubing .4mm OD PEEK 5/pkg
FS1.4V1-5	Fused silica adaptor 1/16in	for tubing .4mm OD Polyimide Virgin 5/pkg
FS1.5-5	Fused silica adaptor 1/16in	for tubing .5mm OD Valcon T 5/pkg
FS1.5PK-5	Fused silica adaptor 1/16in	for tubing .5mm OD PEEK 5/pkg
FS1.5V1-5	Fused silica adaptor 1/16in	for tubing .5mm OD Polyimide Virgin 5/pkg
FS1.75-R10	Ferrule for fused silica union	for .53mm ID tubing Valcon T 10/pkg
FS1.8-5	Fused silica adaptor 1/16in	for tubing .8mm OD Valcon T 5/pkg
FS1.8PK-5	Fused silica adaptor 1/16in	for tubing .8mm OD PEEK 5/pkg
FS1.8V1-5	Fused silica adaptor 1/16in	for tubing .8mm OD Polyimide Virgin 5/pkg
FS1.9-5	Fused silica adaptor 1/16in	for tubing .9mm OD Valcon T 5/pkg
FS1.9PK-5	Fused silica adaptor 1/16in	for tubing .9mm OD PEEK 5/pkg
FS1.9V1-5	Fused silica adaptor 1/16in	for tubing .9mm OD Polyimide Virgin 5/pkg
FS11.0-5	Fused silica adaptor 1/16in	for tubing 1mm OD Valcon T 5/pkg
FS11.0PK-5	Fused silica adaptor 1/16in	for tubing 1mm OD PEEK 5/pkg
FS11.0V1-5	Fused silica adaptor 1/16in	for tubing 1mm OD Polyimide Virgin 5/pkg
FS1L.15PK-5	Liner fused silica adaptor 1/16in	for tubing .15mm OD PEEK 5/pkg
FS1L.2-5	Liner fused silica adaptor 1/16in	for tubing .2mm OD Polyimide Valcon 5/pkg
FS1L.2PK-5	Liner fused silica adaptor 1/16in	for tubing .2mm OD PEEK 5/pkg
FS1L.4-5	Liner fused silica adaptor 1/16in	for tubing .4mm OD Polyimide Valcon 5/pkg
FS1L.4PK-5	Liner fused silica adaptor 1/16in	for tubing .4mm OD PEEK 5/pkg
FS1L.5-5	Liner fused silica adaptor 1/16in	for tubing .5mm OD Polyimide Valcon 5/pkg
FS1L.5PK-5	Liner fused silica adaptor 1/16in	for tubing .5mm OD PEEK 5/pkg
FS1L.8-5	Liner fused silica adaptor 1/16in	for tubing .8mm OD Polyimide Valcon 5/pkg
FS1L.8PK-5	Liner fused silica adaptor 1/16in	for tubing .8mm OD PEEK 5/pkg
FS1L1.0-5	Liner fused silica adaptor 1/16in	for tubing 1mm OD Polyimide Valcon 5/pkg
FS1L1.0PK-5	Liner fused silica adaptor 1/16in	for tubing 1mm OD PEEK 5/pkg
FS1R.15PK-5	Fused silica adaptor 1/16in removable	for tubing .15mm OD PEEK 5/pkg
FS1R.2-5	Fused silica adaptor 1/16in removable	for tubing .2mm OD Polyimide Valcon 5/pkg
FS1R.2PK-5	Fused silica adaptor 1/16in removable	for tubing .2mm OD PEEK 5/pkg
FS1R.4-5	Fused silica adaptor 1/16in removable	for tubing .4mm OD Polyimide Valcon 5/pkg
FS1R.4PK-5	Fused silica adaptor 1/16in removable	for tubing .4mm OD PEEK 5/pkg
FS1R.5-5	Fused silica adaptor 1/16in removable	for tubing .5mm OD Polyimide Valcon 5/pkg
FS1R.5PK-5	Fused silica adaptor 1/16in removable	for tubing .5mm OD PEEK 5/pkg
FS1R.8-5	Fused silica adaptor 1/16in removable	for tubing .8mm OD Polyimide Valcon 5/pkg
FS1R.8PK-5	Fused silica adaptor 1/16in removable	for tubing .8mm OD PEEK 5/pkg
FS1R1.0-5	Fused silica adaptor 1/16in removable	for tubing 1mm OD Polyimide Valcon 5/pkg
FS1R1.0PK-5	Fused silica adaptor 1/16in removable	for tubing 1mm OD PEEK 5/pkg
FSKZU1	Fused silica union 1/16in fingertight	.5mm ID quartz liner
FSL.25-5	Liner fused silica adaptor 1/32in	for tubing .25mm OD Polyimide Valcon 5/pkg
FSL.4-5	Liner fused silica adaptor 1/32in	for tubing .4mm OD Polyimide Valcon 5/pkg
FSL.5-5	Liner fused silica adaptor 1/32in	for tubing .5mm OD Polyimide Valcon 5/pkg
FSMUA1.5	Fused silica make-up adapter 1/16in x 1/32in	89mm long .75mm bore SS
FSMUAS1.5	Fused silica make-up adapter 1/16in x 1/32in	38mm long .75mm bore SS
FSMUAS1.5L	Fused silica make-up adapter 1/16in x 1/32in	38mm long 1.0mm bore SS
FSMUAS1.5M	Fused silica make-up adapter 1/16in x 1/32in	38mm long .5mm bore SS
FSQ1	Replacement liner for fused silica union	.5mm ID quartz

# VICI Valco-Cheminert PRODUCT LIST 2012




FSR.25-5	Fused silica adaptor 1/32in	for tubing .25mm OD Polyimide Valcon 5/pkg
FSR.3-5	Fused silica adaptor 1/32in removable	for tubing .3mm OD PolyimideValcon 5/pkg
FSR.4-5	Fused silica adaptor 1/32in	for tubing .4mm OD Polyimide Valcon 5/pkg
FSR.5-5	Fused silica adaptor 1/32in	for tubing .5mm OD Polyimide Valcon 5/pkg
FSZA-HP	Nut self compensating for HP5890/6890	with split/splitless seal disk
FSZN1	Replacement nut for fused silica union	self compensating SS
FSZNA-HP	Nut self compensating for HP5890/6890	SS
F-TCF	Tube connector for air actuator	compressed fitting 1/8 ftg to 10-32 thread
GH-LNR-CIS4-B-5	Reduced breakdown liner Gerstel CIS-4/PTV	Baffled 5/pkg
GH-LNR-DGS4-5	Reduced breakdown liner Agilent/Finnigan	4mm double gooseneck 5/pkg
GH-LNR-GS2-5	Reduced breakdown liner Agilent/Finnigan/Varian CP-1177	2mm gooseneck 5/pkg
GH-LNR-GS4-5	Reduced breakdown liner Agilent/Finnigan/Varian CP-1177	4mm gooseneck 5/pkg
GH-LNR-HP2-5	Reduced breakdown liner Agilent/Finnigan	2mm straight splitless 5/pkg
GH-LNR-HP4-5	Reduced breakdown liner Agilent/Finnigan	4mm straight splitless 5/pkg
GH-LNR-VAR2-5	replaced by GH-LNR-GS2-5	
GH-LNR-VAR3.4-5	Reduced breakdown liner Varian 1078/9	3.4mm gooseneck 5/pkg
GH-LNR-VAR4-5	replaced by GH-LNR-GS4-5	
GH-LNR-VARG2-5	Reduced breakdown liner Varian 1078/9	2mm gooseneck 5/pkg
HA1-110	Heater assembly for single valve	65W/110Vac
HA1-220	Heater assembly for single valve	65W/220Vac
HA1T-110	Heater assembly for single valve	100W/110Vac
HA1T-220	Heater assembly for single valve	100W/220Vac
HA2-110	Heater assembly for two valves	125W/110Vac
HA2-220	Heater assembly for two valves	125W/220Vac
HA3-110	Heater assembly for three valves	150W/110Vac
HA3-220	Heater assembly for three valves	150W/220Vac
HA6-110	Heater assembly for six valves	300W/110Vac
HA6-220	Heater assembly for six valves	300W/220Vac
HB	Heater block standard for single valve	w/o heater cartridge
HB1N	Heater block for single NW type valve	w/o heater cartridge
HBS	Heater block low mass for single valve	w/o heater cartridge
HCE1-110	Heated column enclosure 1 col. mandrel	4in x4.25in x3.63in d 65W/110Vac
HCE1-220	Heated column enclosure 1 col. mandrel	4in x4.25in x3.63in d 65W/220Vac
HCE2-110	Heated column enclosure 1 col. mandrel	8in x5.25in x4in d 65W/110Vac
HCE2-220	Heated column enclosure 1 col. mandrel	8in x5.25in x4in d 65W/220Vac
HCEB-110	Heated column enclosure 1 col. mandrel	4.25in x5.13in x3.63in d 65W/110Vac
HCEB-220	Heated column enclosure 1 col. mandrel	4.25in x5.13in x3.63in d 65W/220Vac
HCEC-110	Heated column enclosure 1 col. mandrel	8in x8in x6in d 65W/110Vac
HCEC-220	Heated column enclosure 1 col. mandrel	8in x8in x6in d 65W/220Vac
HKS	Hex key set	
HP2-110	Helium purifier Model 2	24VDC input 110/230Vac US power cable
HP2-220	Helium purifier Model 2	24VDC input 110/230Vac EUR power cable
HPM-110	Helium purifier Mini version	24VDC input 110/230Vac US power cable
HPM-220	Helium purifier Mini version	24VDC input 110/230Vac EUR power cable
HSSA	High speed switching assembly for DVI	
HVE2-110	Heated valve enclosure for 2 valves	8in x5.25in x4in 125W/110Vac; Non-CE !!!
HVE2-220	Heated valve enclosure for 2 valves	8in x5.25in x4in 125W/220Vac; Non-CE !!!
HVE3-110	Heated valve enclosure for 3 valves	13.5in x5.75in x4in 150W/110Vac; Non-CE !!!
HVE3-220	Heated valve enclosure for 3 valves	13.5in x5.75in x4in 150W/220Vac; Non-CE !!!
HVE6-110	Heated valve enclosure for 6 valves	13.75in x8in x6in 300W/110Vac; Non-CE !!!
HVE6-220	Heated valve enclosure for 6 valves	13.75in x8in x6in 300W/220Vac; Non-CE !!!
HVEA-110	Heated valve enclosure for 1 valve	4in x4.25in x3.63in 65W/110Vac; Non-CE !!!
HVEA-220	Heated valve enclosure for 1 valve	4in x4.25in x3.63in 65W/220Vac; Non-CE !!!
HVEAN-110	Heated valve enclosure for NW valve	4in x4.25in x3.63in 65W/110Vac Non-CE !!!
HVEAN-220	Heated valve enclosure for NW valve	4in x4.25in x3.63in 65W/220Vac; Non-CE !!!
HVEB-110	Heated valve enclosure for 1 valve	4.25in x5.13in x3.63in 65W/110Vac; Non-CE !!!
HVEB-220	Heated valve enclosure for 1 valve	4.25in x5.13in x3.63in 65W/220Vac; Non-CE !!!
HVEBN-110	Heated valve enclosure for NW valve	4.25in x5.13in x3.63in 65W/110Vac Non-CE !!!
HVEBN-220	Heated valve enclosure for NW valve	4.25in x5.13in x3.63in 65W/220Vac; Non-CE !!!
HVEC-110	Heated valve enclosure for 1 valve	8in x8in x6in 100W/110Vac; Non-CE !!!
HVEC-220	Heated valve enclosure for 1 valve	8in x8in x6in 100W/220Vac; Non-CE !!!

# VICI Valco-Cheminert PRODUCT LIST 2012




HWSC-SC8-6	Screw 8-32 x 3/8in long	counter sunk SS
HWSC-SC8-8B	Screw 8-32 x 1/2in long	counter sunk black oxide SS.
HWSC-SC8-8TDH	Screw 8-32 x 1/2in long	counter sunk turned down head SS
I-21014-01	Thermocouple wire assy	ITC typ in Kin
I-21208-05	Heater cartridge repl. for HVEC	100W/110Vac
I-21208-06	Heater cartridge repl. for HVEC	100W/220Vac
I-21208-32	Heater cartridge for HVEA/HVEB	65W/120Vac
I-21208-33	Heater cartridge repl. for HVEA/HVEB	65W/220Vac
I-22239	Interface cable DVI/SVI	
I-22697	Interface cable RS-232 6 feet	for 2-pos/mpos micro electric actuator
I-22897-02	Interface cable RS-232 2 output	for 2-pos/mpos micro electric actuator
I-22897-03	Interface cable RS-232 3 output	for 2-pos/mpos micro electric actuator
I-22897-04	Interface cable RS-232 4 output	for 2-pos/mpos micro electric actuator
I-22897-05	Interface cable RS-232 5 output	for 2-pos/mpos micro electric actuator
I-22897-06	Interface cable RS-232 6 output	for 2-pos/mpos micro electric actuator
I-22897-08	Interface cable RS-232 8 output	for 2-pos/mpos micro electric actuator
I-23572HP2	Getter assembly replacement for HP2	helium purifier
I-23572NP2	Getter assembly replacement for NP2	nitrogen purifier
I4UWE.2	4 port 2-pos int vol valve manual	1/8 x .75mm 175C/1000psi liq .2ul N60/E
I4UWE.5	4 port 2-pos int vol valve manual	1/8 x .75mm 175C/1000psi liq .5ul N60/E
I4UWE1	4 port 2-pos int vol valve manual	1/8 x .75mm 175C/1000psi liq 1ul N60/E
I4UWE2	4 port 2-pos int vol valve manual	1/8 x .75mm 175C/1000psi liq 2ul N60/E
ITC10399-110	Instrumentation temperature controller	0-399 deg. C 10A/110Vac
ITC10399-220	Instrumentation temperature controller	0-399 deg. C 10A/220Vac
IZERA1.5	Int./Ext. reducer/adaptor 1/16in x 1/32in	1mm bore SS
IZERA1.5C	Int./Ext. reducer/adaptor 1/16in x 1/32in	.25mm bore SS
IZERA1.5M	Int./Ext. reducer/adaptor 1/16in x 1/32in	.5mm bore SS
IZR1.5	Internal reducer 1/16in x 1/32in	.25mm bore SS
IZR1.5F	Internal reducer 1/16in x 1/32in	.25mm bore 2 um frit SS
IZR1.5L	Internal reducer 1/16in x 1/32in	.50mm bore SS
IZR1.5LF	Internal reducer 1/16in x 1/32in	.50mm bore 2 um frit SS
IZR1.5T	Internal reducer 1/16in x 1/32in	1/32in through bore SS
IZR21	Internal reducer 1/8in x 1/16in	.50mm bore SS
IZR21C	Internal reducer 1/8in x 1/16in	.25mm bore SS
IZR21CF	Internal reducer 1/8in x 1/16in	.25mm bore 2 um frit SS
IZR21F	Internal reducer 1/8in x 1/16in	.50mm bore 2 um frit SS
IZR21L	Internal reducer 1/8in x 1/16in	1mm bore SS
IZR21LF	Internal reducer 1/8in x 1/16in	1mm bore 2 um frit SS
IZR21T	Internal reducer 1/8in x 1/16in	1/16 trough bore SS
IZR41	Internal reducer 1/4in x 1/16in	1mm bore SS
IZR41F	Internal reducer 1/4in x 1/16in	1mm bore 2 um frit SS
IZR42	Internal reducer 1/4in x 1/8in	1mm bore SS
IZR42F	Internal reducer 1/4in x 1/8in	1mm bore 2 um frit SS
IZR42L	Internal reducer 1/4in x 1/8in	2mm bore SS
IZR42LF	Internal reducer 1/4in x 1/8in	2mm bore 2 um frit SS
JR-201540	Easy-Flange kit in plastic case	
JR-26501	Coupler PEEK universal 0.13 mm ID red	
JR-26502	Coupler PEEK universal 0.17 mm ID yellow	
JR-26503	Coupler PEEK universal 0.25 mm ID blue	
JR-26504	Coupler PEEK universal 0.50 mm ID orange	
JR-32178	Filter PE mobile phase No-Met 20 um 1/8"	
JR-32179	Filter PE mobile phase replacement No-Met 20 um	
JR-357090-5	Tubing Elbow 90° 5/pkg	
JR-357180-5	Tubing Elbow 180° 5/pkg	
JR-35P	PEEK Starter Kit	
JR-367008-10	Filter SS mobile phase 10 um	fittings for 1/8" tubing
JR-367008-2	Filter SS mobile phase 2 um	fittings for 1/8" tubing
JR-367008-20	Filter SS mobile phase 20 um	fittings for 1/8" tubing
JR-367016-10	Filter SS mobile phase 10 um	fittings for 1/16" tubing
JR-367016-2	Filter SS mobile phase 2 um	fittings for 1/16" tubing
JR-367016-20	Filter SS mobile phase 20 um	fittings for 1/16" tubing

# VICI Valco-Cheminert PRODUCT LIST 2012




JR-55010-5	Color-It fingertight adapter blue	5/pkg
JR-55011-5	Color-It fingertight adapter yellow	5/pkg
JR-55012-5	Color-It fingertight adapter green	5/pkg
JR-55013-5	Color-It fingertight adapter black	5/pkg
JR-55014-5	Color-It fingertight adapter white	5/pkg
JR-55015-5	Color-It fingertight adapter red	5/pkg
JR-55016-12	Color-It fingertight adapter multi-color	assorted 2 pcs. each color 12/pkg
JR-55020-5	Fitting PEEK one-piece natural 10-32	5/pkg
JR-55021-5	Fitting PEEK one-piece black 10-32	5/pkg
JR-55022-5	Fitting PEEK one-piece red 10-32	5/pkg
JR-55023-5	Fitting PEEK one-piece yellow 10-32	5/pkg
JR-55024-5	Fitting PEEK one-piece blue 10-32	5/pkg
JR-55025-5	Fitting PEEK one-piece green 10-32	5/pkg
JR-5508-5	Fitting PEEK 1/16in OD tubing	one-piece fingertight hex-head 10-32 5/pkg
JR-797	Tool Clean-Cut Tubing Cutter	
JR-798	Tool replacement blade for JR-797	
JR-9000-0001	Cap PTFE VICI-Cap GL-45	
JR-9000-0006	Cap PTFE VICI-Cap S40	
JR-9000-0520	Filter PTFE Last Drop 2.5 um	tripod stepped tubing connector
JR-9000-0530	Filter SS Last Drop 2 um	tripod stepped tubing connector
JR-9000-0602	Filter/Sparger PTFE Last Drop 2.5 um filter	10 um sparger
JR-9000-0640	Filter/Sparger SS Last Drop 2 um filter	10 um sparger
JR-9001-5	Tubing Clip 5/pkg	
JR-T-5804-M10	Tubing PEEK 1/16" x 0.100 mm ID black striped	10 m/pkg
JR-T-5804-M25	Tubing PEEK 1/16" x 0.100 mm ID black striped	25 m/pkg
JR-T-5804-M3	Tubing PEEK 1/16" x 0.100 mm ID black striped	3 m/pkg
JR-T-5804-M5	Tubing PEEK 1/16" x 0.100 mm ID black striped	5 m/pkg
JR-T-5999-M10	Tubing PEEK 1/16" x 0.13 mm ID red striped	10m/pkg
JR-T-5999-M25	Tubing PEEK 1/16" x 0.13 mm ID red stripe	25m/pkg
JR-T-5999-M3	Tubing PEEK 1/16" x 0.13 mm ID red striped	3m/pkg
JR-T-5999-M5	Tubing PEEK 1/16" x 0.13 mm ID red stripe	5m/pkg
JR-T-6000-M10	Tubing PEEK 1/16" x 0.18 mm ID yellow striped	10m/pkg
JR-T-6000-M25	Tubing PEEK 1/16" x 0.18 mm ID yellow striped	25 m/pkg
JR-T-6000-M3	Tubing PEEK 1/16" x 0.18 mm ID yellow striped	3m/pkg
JR-T-6000-M5	Tubing PEEK 1/16" x 0.18 mm ID yellow striped	5 m/pkg
JR-T-6001-M10	Tubing PEEK 1/16" x 0.25 mm ID blue striped	10m/pkg
JR-T-6001-M25	Tubing PEEK 1/16" x 0.25 mm ID blue striped	25m/pkg
JR-T-6001-M3	Tubing PEEK 1/16" x 0.25 mm ID blue striped	3m/pkg
JR-T-6001-M5	Tubing PEEK 1/16" x 0.25 mm ID blue stripe	5m/pkg
JR-T-6002-M10	Tubing PEEK 1/16" x 0.50 mm ID orange striped	10m/pkg
JR-T-6002-M25	Tubing PEEK 1/16" x 0.50 mm ID orange stripe	25m/pkg
JR-T-6002-M3	Tubing PEEK 1/16" x 0.50 mm ID orange striped	3m/pkg
JR-T-6002-M5	Tubing PEEK 1/16" x 0.50 mm ID orange stripe	5m/pkg
JR-T-60031-M10	Tubing PEEK 1/16" x 1.00 mm ID grey striped	10m/pkg
JR-T-60031-M25	Tubing PEEK 1/16" x 1.00 mm ID grey striped	25 m/pkg
JR-T-60031-M3	Tubing PEEK 1/16" x 1.00 mm ID grey striped	3m/pkg
JR-T-60031-M5	Tubing PEEK 1/16" x 1.00 mm ID grey striped	5 m/pkg
JR-T-6003-M10	Tubing PEEK 1/16" x 0.75 mm ID green striped	10m/pkg
JR-T-6003-M25	Tubing PEEK 1/16" x 0.75 mm ID green stripe	25m/pkg
JR-T-6003-M3	Tubing PEEK 1/16" x 0.75 mm ID green striped	3m/pkg
JR-T-6003-M5	Tubing PEEK 1/16" x 0.75 mm ID green stripe	5m/pkg
JR-T-98013	Tubing PEEK precut kit 1/16" x 0.13 mm	
JR-T-98017	Tubing PEEK precut kit 1/16" x 0.18 mm	
JR-T-98025	Tubing PEEK precut kit 1/16" x 0.25 mm	
JR-T-98050	Tubing PEEK precut kit 1/16" x 0.50 mm	
JR-T-98075	Tubing PEEK precut kit 1/16" x 0.75 mm	
JR-T-98200	Tubing PEEK precut super kit 1/16"	
JR-TD-5804-M10	Tubing Dual Layer PEEK 1/16" x 100 um ID	solid black outside natural inside 10m/pkg
JR-TD-5804-M25	Tubing Dual Layer PEEK 1/16" x 100 um ID	solid black outside natural inside 25m/pkg
JR-TD-5804-M3	Tubing Dual Layer PEEK 1/16" x 100 um ID	solid black outside natural inside 3m/pkg
JR-TD-6007-M10	Tubing Dual Layer PEEK 1/16" x 0.13 mm ID	solid red outside natural inside 10m/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




JR-TD-6007-M25	Tubing Dual Layer PEEK 1/16" x 0.13 mm ID	solid red outside natural inside 25m/pkg
JR-TD-6007-M3	Tubing Dual Layer PEEK 1/16" x 0.13 mm ID	solid red outside natural inside 3m/pkg
JR-TD-6008-M10	Tubing Dual Layer PEEK 1/16" x 0.18 mm ID	solid yellow outside natural inside 10m/pkg
JR-TD-6008-M25	Tubing Dual Layer PEEK 1/16" x 0.18 mm ID	solid yellow outside natural inside 25m/pkg
JR-TD-6008-M3	Tubing Dual Layer PEEK 1/16" x 0.18 mm ID	solid yellow outside natural inside 3m/pkg
JR-TD-6009-M10	Tubing Dual Layer PEEK 1/16" x 0.25 mm ID	solid blue outside natural inside 10m/pkg
JR-TD-6009-M25	Tubing Dual Layer PEEK 1/16" x 0.25 mm ID	solid blue outside natural inside 25m/pkg
JR-TD-6009-M3	Tubing Dual Layer PEEK 1/16" x 0.25 mm ID	solid blue outside natural inside 3m/pkg
JR-TD-6010-M10	Tubing Dual Layer PEEK 1/16" x 0.50 mm ID	solid orange outside natural inside 10m/pkg
JR-TD-6010-M25	Tubing Dual Layer PEEK 1/16" x 0.50 mm ID	solid orange outside natural inside 25m/pkg
JR-TD-6010-M3	Tubing Dual Layer PEEK 1/16" x 0.50 mm ID	solid orange outside natural inside 3m/pkg
JR-TD-6011-M10	Tubing Dual Layer PEEK 1/16" x 0.75 mm ID	solid green outside natural inside 10m/pkg
JR-TD-6011-M25	Tubing Dual Layer PEEK 1/16" x 0.75 mm ID	solid green outside natural inside 25m/pkg
JR-TD-6011-M3	Tubing Dual Layer PEEK 1/16" x 0.75 mm ID	solid green outside natural inside 3m/pkg
L4UW	!!! RECOMMEND 2L4UW	SUPPLIED WITH LARGER HANDLE
L6UW	!!! RECOMMEND 2L6UW	SUPPLIED WITH LARGER HANDLE
L8UW	!!! RECOMMEND 2L8UW	SUPPLIED WITH LARGER HANDLE
LZCN1-10	Counter bored nut 1/16in long	SS 10/pkg
LZN.5-10	Nut 1/32in 11mm long SS 10/pkg	
LZN.5FPK-10	Nut 1/32in fingertight 13.7mm long	PEEK 10/pkg
LZN1-10	Nut 1/16in 19mm long SS 10/pkg	
LZN1PK-10	Nut 1/16in hex 22mm long	PEEK 10/pkg
LZN2-10	Nut 1/8in 21mm long SS 10/pkg	
LZN4-10	Nut 1/4in 28mm long SS 10/pkg	
LZNF1KF-5	Nut 1/16in No Twist hex long PEEK	with CTFE ferrule 5/pkg
LZNF1PK-5	Nut 1/16in No Twist hex long PEEK	with PEEK ferrule 5/pkg
LZNF1PKG-5	Nut 1/16in No Twist hex long PEEK	with glass-filled PEEK ferrule 5/pkg
LZP.5FPK	Plug 1/32in fingertight PEEK	
LZP1	Plug 1/16in 19mm long	SS
LZP1H	Plug 1/16in 19mm long high pressure	SS
LZP1PK	Plug 1/16in long hex PEEK	
LZP2	Plug 1/8in 21mm long	SS
LZP2H	Plug 1/8in 21mm long high pressure	SS
MRX.5PK	Micro volume replacement ring 1/32in	for union tee cross PEEK
MRX.5S6	Micro volume replacement ring 1/32in	for union tee cross 316SS
MRX1PK	Micro volume replacement ring 1/16in	for union tee cross PEEK
MRX1S6	Micro volume replacement ring 1/16in	for union tee cross 316SS
MRY.5PK	Micro volume replacement ring 1/32in	for Y PEEK
MRY.5S6	Micro volume replacement ring 1/32in	for Y 316SS
MRY1PK	Micro volume replacement ring 1/16in	for Y PEEK
MRY1S6	Micro volume replacement ring 1/16in	for Y 316SS
MS-P100-1	Helium Purifier	1/8in brass caps
MS-P100-2	Helium Purifier	1/4in brass caps
MS-P200-1	Hydrogen Purifier	1/8in brass caps
MS-P200-2	Hydrogen Purifier	1/4in brass caps
MS-P300-1	Nitrogen Purifier	1/8in brass caps
MS-P300-2	Nitrogen Purifier	1/4in brass caps
MS-P310-1	Nitrogen Purifier for LC/MS applications	1/8in brass caps
MS-P310-2	Nitrogen Purifier for LC/MS applications	1/4in brass caps
MS-P350-1	Nitrogen Purifier for nitrogen generators	1/8in brass caps
MS-P350-2	Nitrogen Purifier for nitrogen generators	1/4in brass caps
MS-P400-1	Air Purifier	1/8in brass caps
MS-P400-2	Air Purifier	1/4in brass caps
MS-P500-1	Methane Purifier/CI	1/8in brass caps
MS-P500-2	Methane Purifier/CI	1/4in brass caps
MS-P600-1	Carbon dioxide	1/8in brass caps
MS-P600-2	Carbon dioxide	1/4in brass caps
MS-P700-1	CO2 Purifier	1/8in brass caps
MS-P700-2	CO2 Purifier	1/4in brass caps
MS-T100-1	Moisture Trap	1/8in brass caps
MS-T100-2	Moisture Trap	1/4in brass caps

# VICI Valco-Cheminert PRODUCT LIST 2012



MS-T200-1	Hydrocarbon Trap	1/8in brass caps
MS-T200-2	Hydrocarbon Trap	1/4in brass caps
MS-T300-1	Oxygen Trap	1/8in brass caps
MS-T300-2	Oxygen Trap	1/4in brass caps
MS-T400-1	Sulfur/Halocarbon Trap	1/8in brass caps
MS-T400-2	Sulfur/Halocarbon Trap	1/4in brass caps
MS-T401-1	Sulfur/Halocarbon Trap	1/8in brass caps
MS-T401-2	Sulfur/Halocarbon Trap	1/4in brass caps
MS-T700-1	Mercury Trap	1/8in brass caps
MS-T700-2	Mercury Trap	1/4in brass caps
MSVA-110VAC	3-way manifold solenoid valve assy. 110Vac	Non-CE !!!
MSVA-12VDC	3-way manifold solenoid valve assy. 12Vdc	Non-CE !!!
MSVA-220VAC	3-way manifold solenoid valve assy. 220Vac	Non-CE !!!
MSVA-24VAC	3-way manifold solenoid valve assy. 24Vac	Non-CE !!!
MSVA-24VDC	3-way manifold solenoid valve assy. 24Vdc	Non-CE !!!
MT.5CHC	Micro volume tee-assy 1/32in	.25mm bore HC
MT.5CKF	Micro volume tee-assy 1/32in	.25mm bore CTFE
MT.5CPK	Micro volume tee-assy 1/32in	.25mm bore PEEK
MT.5CS6	Micro volume tee-assy 1/32in	.25mm bore 316SS
MT.5CTI	Micro volume tee-assy 1/32in	.25mm bore Titanium
MT.5XCHC	Micro volume tee-assy 1/32in	.15mm bore HC
MT.5XCKF	Micro volume tee-assy 1/32in	.15mm bore CTFE
MT.5XCPK	Micro volume tee-assy 1/32in	.15mm bore PEEK
MT.5XCS6	Micro volume tee-assy 1/32in	.15mm bore 316SS
MT.5XCTI	Micro volume tee-assy 1/32in	.15mm bore Titanium
MT1CHC	Micro volume tee-assy 1/16in	.25mm bore HC
MT1CKF	Micro volume tee-assy 1/16in	.25mm bore CTFE
MT1CPK	Micro volume tee-assy 1/16in	.25mm bore PEEK
MT1CS6	Micro volume tee-assy 1/16in	.25mm bore 316SS
MT1CTI	Micro volume tee-assy 1/16in	.25mm bore Titanium
MT1XCHC	Micro volume tee-assy 1/16in	.15mm bore HC
MT1XCKF	Micro volume tee-assy 1/16in	.15mm bore CTFE
MT1XCPK	Micro volume tee-assy 1/16in	.15mm bore PEEK
MT1XCS6	Micro volume tee-assy 1/16in	.15mm bore 316SS
MT1XCTI	Micro volume tee-assy 1/16in	.15mm bore Titanium
MU.5CHC	Micro volume union-assy 1/32in	.25mm bore HC
MU.5CKF	Micro volume union-assy 1/32in	.25mm bore CTFE
MU.5CPK	Micro volume union-assy 1/32in	.25mm bore PEEK
MU.5CS6	Micro volume union-assy 1/32in	.25mm bore 316SS
MU.5CTI	Micro volume union-assy 1/32in	.25mm bore Titanium
MU.5XCHC	Micro volume union-assy 1/32in	.15mm bore HC
MU.5XCKF	Micro volume union-assy 1/32in	.15mm bore CTFE
MU.5XCPK	Micro volume union-assy 1/32in	.15mm bore PEEK
MU.5XCS6	Micro volume union-assy 1/32in	.15mm bore 316SS
MU.5XCTI	Micro volume union-assy 1/32in	.15mm bore Titanium
MU1CHC	Micro volume union-assy 1/16in	.25mm bore HC
MU1CKF	Micro volume union-assy 1/16in	.25mm bore CTFE
MU1CPK	Micro volume union-assy 1/16in	.25mm bore PEEK
MU1CS6	Micro volume union-assy 1/16in	.25mm bore 316SS
MU1CTI	Micro volume union-assy 1/16in	.25mm bore Titanium
MU1XCHC	Micro volume union-assy 1/16in	.15mm bore HC
MU1XCKF	Micro volume union-assy 1/16in	.15mm bore CTFE
MU1XCPK	Micro volume union-assy 1/16in	.15mm bore PEEK
MU1XCS6	Micro volume union-assy 1/16in	.15mm bore 316SS
MU1XCTI	Micro volume union-assy 1/16in	.15mm bore Titanium
MX.5CHC	Micro volume cross-assy 1/32in	.25mm bore HC
MX.5CKF	Micro volume cross-assy 1/32in	.25mm bore CTFE
MX.5CPK	Micro volume cross-assy 1/32in	.25mm bore PEEK
MX.5CS6	Micro volume cross-assy 1/32in	.25mm bore 316SS
MX.5CTI	Micro volume cross-assy 1/32in	.25mm bore Titanium
MX.5XCHC	Micro volume cross-assy 1/32in	.15mm bore HC


# VICI Valco-Cheminert PRODUCT LIST 2012




MX.5XCKF	Micro volume cross-assy 1/32in	.15mm bore CTFE
MX.5XCPK	Micro volume cross-assy 1/32in	.15mm bore PEEK
MX.5XCS6	Micro volume cross-assy 1/32in	.15mm bore 316SS
MX.5XCTI	Micro volume cross-assy 1/32in	.15mm bore Titanium
MX1CHC	Micro volume cross-assy 1/16in	.25mm bore HC
MX1CKF	Micro volume cross-assy 1/16in	.25mm bore CTFE
MX1CPK	Micro volume cross-assy 1/16in	.25mm bore PEEK
MX1CS6	Micro volume cross-assy 1/16in	.25mm bore 316SS
MX1CTI	Micro volume cross-assy 1/16in	.25mm bore Titanium
MX1XCHC	Micro volume cross-assy 1/16in	.15mm bore HC
MX1XCKF	Micro volume cross-assy 1/16in	.15mm bore CTFE
MX1XCPK	Micro volume cross-assy 1/16in	.15mm bore PEEK
MX1XCS6	Micro volume cross-assy 1/16in	.15mm bore 316SS
MX1XCTI	Micro volume cross-assy 1/16in	.15mm bore Titanium
MY.5CHC	Micro volume Y-assy 1/32in	.25mm bore HC
MY.5CKF	Micro volume Y-assy 1/32in	.25mm bore CTFE
MY.5CPK	Micro volume Y-assy 1/32in	.25mm bore PEEK
MY.5CS6	Micro volume Y-assy 1/32in	.25mm bore 316SS
MY.5CTI	Micro volume Y-assy 1/32in	.25mm bore Titanium
MY.5XCHC	Micro volume Y-assy 1/32in	.15mm bore HC
MY.5XCKF	Micro volume Y-assy 1/32in	.15mm bore CTFE
MY.5XCPK	Micro volume Y-assy 1/32in	.15mm bore PEEK
MY.5XCS6	Micro volume Y-assy 1/32in	.15mm bore 316SS
MY.5XCTI	Micro volume Y-assy 1/32in	.15mm bore Titanium
MY1CHC	Micro volume Y-assy 1/16in	.25mm bore HC
MY1CKF	Micro volume Y-assy 1/16in	.25mm bore CTFE
MY1CPK	Micro volume Y-assy 1/16in	.25mm bore PEEK
MY1CS6	Micro volume Y-assy 1/16in	.25mm bore 316SS
MY1CTI	Micro volume Y-assy 1/16in	.25mm bore Titanium
MY1XCHC	Micro volume Y-assy 1/16in	.15mm bore HC
MY1XCKF	Micro volume Y-assy 1/16in	.15mm bore CTFE
MY1XCPK	Micro volume Y-assy 1/16in	.15mm bore PEEK
MY1XCS6	Micro volume Y-assy 1/16in	.15mm bore 316SS
MY1XCTI	Micro volume Y-assy 1/16in	.15mm bore Titanium
MZN1-10	Nut 1/16in 13mm long SS 10/pkg	
MZN1PK-10	Nut 1/16in hex 15.7mm long	PEEK 10/pkg
MZNF1KF-5	Nut 1/16in No Twist hex medium PEEK	with CTFE ferrule 5/pkg
MZNF1PK-5	Nut 1/16in No Twist hex medium PEEK	with PEEK ferrule 5/pkg
MZNF1PKG-5	Nut 1/16in No Twist hex medium PEEK	with glass-filled PEEK ferrule 5/pkg
MZP1PK	Plug 1/16in hex PEEK	
NP2-110	Nitrogen purifier Model 2	24VDC input 110/230Vac US power cable
NP2-220	Nitrogen purifier Model 2	24VDC input 110/230Vac EUR power cable
NPM-110	Nitrogen purifier Mini version	24VDC input 110/230Vac US power cable
NPM-220	Nitrogen purifier Mini version	24VDC input 110/230Vac EUR power cable
OEW-1	Open end wrench 3/16in x 1/4in	for 1/32in and 1/16in fittings
OEW-2	Open end wrench 3/8in x 7/16in	for 1/8in fittings
OEW-3	Open end wrench 1/2in x 9/16in	for 1/4in fittings
OR	O-ring set incl. lubricant	standard 2-pos air actuator
ORMP	O-ring set incl. lubricant	standard mpos air actuator
ORT	O-ring set incl. lubricant	high temperature 2-pos air actuator
ORTMP	O-ring set incl. lubricant	high temperature mpos air actuator
PCL412	Precolumn 1/4in x 2cm 4.6mm ID	incl. insert nut and ferrule 2um frit SS
PCL415	Precolumn 1/4in x 5cm 4.6mm ID	incl. insert nut and ferrule 2um frit SS
PCS412F	Precolumn system 1/4in x 2cm	incl. insert body and 2um frit SS
PCS415F	Precolumn system 1/4in x 5cm	incl. insert body and 2um frit SS
PEA22	Male pipe adaptor 1/8in M-NPT x 1/8in ext	2mm bore SS
PEA24	Male pipe adaptor 1/8in M-NPT x 1/4in ext	4.6mm bore SS
PEA42	Male pipe adaptor 1/4in M-NPT x 1/8in ext	2mm bore SS
PEA44	Male pipe adaptor 1/4in M-NPT x 1/4in ext	4.6mm bore SS
PEA82	Male pipe adaptor 1/2in M-NPT x 1/8in ext	2mm bore SS
PEA84	Male pipe adaptor 1/2in M-NPT x 1/4in ext	4.6mm bore SS

# VICI Valco-Cheminert PRODUCT LIST 2012

PFAF	Position feedback adapter	for air actuator
PFC2	Position feedback	Cheminert C2 and C3 series 6/10 port
PFC4	Position feedback	to flat plate valve C4
PFW36	Position feedback for 2-pos manual valve	W type 8 and 10 port 36 deg
PFW60	Position feedback for 2-pos manual valve	W type 6 port and int. 4 port 60 deg
PFW90	Position feedback for 2-pos manual valve	W type 4 port 90 deg
PM	Pencil magnet	for W and UW rotor seals
PR50A15E2	Combo pressure regulator 0-15psi	1/8in ext Aluminum/SS
PR50A15Z1	Combo pressure regulator 0-15psi	1/16in int Aluminum/SS
PR50A15Z2	Combo pressure regulator 0-15psi	1/8in int Aluminum/SS
PR50A30E2	Combo pressure regulator 0-30psi	1/8in ext Aluminum/SS
PR50A30Z1	Combo pressure regulator 0-30psi	1/16in int Aluminum/SS
PR50A30Z2	Combo pressure regulator 0-30psi	1/8in int Aluminium/SS
PR50A60E2	Combo pressure regulator 0-60psi	1/8in ext Aluminum/SS
PR50A60Z1	Combo pressure regulator 0-60psi	1/16in int Aluminum/SS
PR50A60Z2	Combo pressure regulator 0-60psi	1/8in int Aluminium/SS
PS-050023	Gas syringe A-2 25 ul	RN .028in x.005in x2in bevel open end
PS-050024	Gas syringe A-2 50 ul	RN .029in x.005in x2in bevel open end
PS-050025	Gas syringe A-2 100 ul	RN .028in x.005in x2in bevel open end
PS-050031	Gas syringe A-2 250 ul	RN .029in x.012in x2in bevel open end
PS-050032	Gas syringe A-2 500 ul	RN .029in x.012in x2in bevel open end
PS-050033	Gas syringe A-2 1 ml	RN .029in x.012in x2in bevel open end
PS-050034	Gas syringe A-2 2 ml	RN .029in x.012in x2in bevel open end
PS-050035	Gas syringe A-2 5 ml	RN .029in x.012in x2in bevel open end
PS-050036	Gas syringe A-2 10 ml	RN .029in x.012in x2in bevel open end
PS-050043	Gas syringe A-2 25 ul	RN .028in x.005in x2in bevel open end Luer
PS-050044	Gas syringe A-2 50 ul	RN .028in x.005in x2in bevel open end Luer
PS-050045	Gas syringe A-2 100 ul	RN .028in x.005in x2in bevel open end Luer
PS-050051	Gas syringe A-2 250 ul	RN .029in x.012in x2in bevel open end Luer
PS-050052	Gas syringe A-2 500 ul	RN .029in x.012in x2in bevel open end Luer
PS-050053	Gas syringe A-2 1 ml	RN .029in x.012in x2in bevel open end Luer
PS-050054	Gas syringe A-2 2 ml	RN .029in x.012in x2in bevel open end Luer
PS-050055	Gas syringe A-2 5 ml	RN .029in x.012in x2in bevel open end Luer
PS-050056	Gas syringe A-2 10 ml	RN .029in x.012in x2in bevel open end Luer
PS-123050	Needles removable C-160 CG-130	.019in x.005in x2.25in bevel open end 3/pkg
PS-123050R	Needles removable C-160	22GAX2in blunt tip open end V/C/Rh 3/pkg
PS-160021	Liquid syringe C-160 5 ul	FN .019in x.005in x2in bevel open end
PS-160021R	Liquid syringe C-160 5 ul	FN 22GAX2in blunt tip open end V/C/Rh
PS-160022	Liquid syringe C-160 10 ul	FN .019in x.005in x2in bevel open end
PS-160022R	Liquid syringe C-160 10 ul	FN 22GAX2in blunt tip open end V/C/Rh
PS-160023	Liquid syringe C-160 25 ul	FN .019in x.005in x2in bevel open end
PS-160023R	Liquid syringe C-160 25 ul	FN 22GAX2in blunt tip open end V/C/Rh
PS-160024	Liquid syringe C-160 50 ul	FN .019in x.005in x2in bevel open end
PS-160024R	Liquid syringe C-160 50 ul	FN 22GAX2in blunt tip open end V/C/Rh
PS-160025	Liquid syringe C-160 100 ul	FN 019in x.005in x2in bevel open end
PS-160025R	Liquid syringe C-160 100 ul	FN 22GAX2in blunt tip open end V/C/Rh
PS-160221	Liquid syringe C-160 5 ul	RN .019in x.005in x2.25in bevel open end
PS-160221R	Liquid syringe C-160 5 ul	RN 22GAX2in blunt tip open end V/C/Rh
PS-160222	Liquid syringe C-160 10 ul	RN .019in x.005in x2.25in bevel open end
PS-160222R	Liquid syringe C-160 10 ul	RN 22GAX2in blunt tip open end V/C/Rh
PS-160223	Liquid syringe C-160 25 ul	RN .019in x.005in x2.25in bevel open end
PS-160223R	Liquid syringe C-160 25 ul	RN 22GAX2in blunt tip open end V/C/Rh
PS-160224	Liquid syringe C-160 50 ul	RN .019in x.005in x2.25in bevel open end
PS-160224R	Liquid syringe C-160 50 ul	RN 22GAX2in blunt tip open end V/C/Rh
PS-160225	Liquid syringe C-160 100 ul	RN .019in x.005in x2.25in bevel open end
PS-160225R	Liquid syringe C-160 100 ul	RN 22GAX2in blunt tip open end V/C/Rh
PS-614158	Mininert valve 13mm-425 12/pkg	
PS-614160	Mininert valve 15mm-425 12/pkg	
PS-614161	Mininert valve 18mm-400 12/pkg	
PS-614163	Mininert valve 24mm-400 12/pkg	
PS-614170	Mininert valve 20mm-400 12/pkg	




# VICI Valco-Cheminert PRODUCT LIST 2012



PS-614250	Mininert valve crimp-top 12/pkg	for 13 mm ID glassware
PS-631201	Mininert valve termination	1/4in -28 male to Mininert
PS-631202	Mininert valve termination	1/8in male NPT to Mininert
PS-631203	Mininert valve termination	1/4in -28 female to Mininert
PS-631204	Mininert valve termination	1/8in female NPT to Mininert
PS-631205	Mininert valve in-line	1/4in -28 male to 1/4in -28 male
PS-644350	Mininert needle seal septa 50/pkg	
PS-644850	Septum installation tool	
PS-654050	Mininert syringe valve	for C or D syringe
PS-654051	Mininert syringe valve	for luer-tip
PS-660100	Valve 3-port manual T-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660110	Valve 4-port manual T-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660200	Valve 2-port manual 180-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660210	Valve 4-port manual 180-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660300	Valve 2-port manual 90-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660310	Valve 3-port manual 90-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-660320	Valve 4-port manual 90-flow w/o nutsferrules	200psi .06" bore 1/4-28 fitting detail
PS-943050	Needles RN p-lok A-2 C	.028in x.005in x2in bevel open end 3/pkg
PS-943051	Needles RN p-lok A-2 C D	.029in x.012in x2in bevel open end 3/pkg
PS-943052	Needles RN p-lok A-2 C D	.029in x.012in x2in side-port taper 3/pkg
PS-943060	Needles RN luer A-2 C Magnum	.028in x.006in x2in bevel open end 3/pkg
PS-943061	Needles RN luer A-2 D	.028in x.012in x2in bevel open end 3/pkg
PS-943062	Needles RN luer A-2 C D Magnum	.028in x.012in x2in side port taper 3/pkg
PV	Pin vise and drill index .34-1mm	
PZA21	Male pipe adaptor 1/8in M-NPT x 1/16in int	1mm bore SS
PZA21T	Male pipe adaptor 1/8in M-NPT x 1/16in int	1/16in trough bore SS
PZA22	Male pipe adaptor 1/8in M-NPT x 1/8in int	1mm bore SS
PZA41	Male pipe adaptor 1/4in M-NPT x 1/16in int	1mm bore SS
PZA42	Male pipe adaptor 1/4in M-NPT x 1/8in int	1mm bore SS
PZA42L	Male pipe adaptor 1/4in M-NPT x 1/8in int	2mm bore SS
PZA44L	Male pipe adaptor 1/4in M-NPT x 1/4in int	4.6mm bore SS
PZA81	Male pipe adaptor 1/2in M-NPT x 1/16in int	1mm bore SS
PZA82	Male pipe adaptor 1/2in M-NPT x 1/8in int	1mm bore SS
PZA82L	Male pipe adaptor 1/2in M-NPT x 1/8in int	2mm bore SS
PZA84L	Male pipe adaptor 1/2in M-NPT x 1/4in int	4.6mm bore SS
RAD	Right angle drive close mount	
RF1.5KF-5	Standard reducing ferrule 1/16in x 1/32in	CTFE 5/pkg
RF1.5PK-5	Standard reducing ferrule 1/16in x 1/32in	PEEK 5/pkg
RF1.5TFG-5	Standard reducing ferrule 1/16in x 1/32in	PTFE Glass-filled 25% 5/pkg
RF1.5V-5	Standard reducing ferrule 1/16in x 1/32in	Polyimide Valcon 5/pkg
RF2.5KF-5	Standard reducing ferrule 1/8in x 1/32in	CTFE 5/pkg
RF2.5PK-5	Standard reducing ferrule 1/8in x 1/32in	PEEK 5/pkg
RF2.5TFG-5	Standard reducing ferrule 1/8in x 1/32in	PTFE Glass-filled 25% 5/pkg
RF2.5V-5	Standard reducing ferrule 1/8in x 1/32in	Polyimide Valcon 5/pkg
RF21KF-5	Standard reducing ferrule 1/8in x 1/16in	CTFE 5/pkg
RF21PK-5	Standard reducing ferrule 1/8in x 1/16in	PEEK 5/pkg
RF21TFG-5	Standard reducing ferrule 1/8in x 1/16in	PTFE Glass-filled 25% 5/pkg
RF21V-5	Standard reducing ferrule 1/8in x 1/16in	Polyimide Valcon 5/pkg
RF41KF-5	Standard reducing ferrule 1/4in x 1/16in	CTFE 5/pkg
RF41PK-5	Standard reducing ferrule 1/4in x 1/16in	PEEK 5/pkg
RF41TFG-5	Standard reducing ferrule 1/4in x 1/16in	PTFE Glass-filled 25% 5/pkg
RF41V-5	Standard reducing ferrule 1/4in x 1/16in	Polyimide Valcon 5/pkg
RF42KF-5	Standard reducing ferrule 1/4in x 1/8in	CTFE 5/pkg
RF42PK-5	Standard reducing ferrule 1/4in x 1/8in	PEEK 5/pkg
RF42TFG-5	Standard reducing ferrule 1/4in x 1/8in	PTFE Glass-filled 25% 5/pkg
RF42V-5	Standard reducing ferrule 1/4in x 1/8in	Polyimide Valcon 5/pkg
SD4UW	4 mpos dead-end path manual	1/8 x .75mm 75C/5000psi liq N60/E
SD6UW	6 mpos dead-end path manual	1/8 x .75mm 75C/5000psi liq N60/E
SD8UW	8 mpos dead-end path manual	1/8 x .75mm 75C/5000psi liq N60/E
SEAL1-HP	Seal disk split/splitless for HP5890	SS
SFV	Prime/Purge valve manual	3 x 1/16in .50 mm bore 100C/10'000psi SS


# VICI Valco-Cheminert PRODUCT LIST 2012




SFVL	Prime/Purge valve manual	3 x 1/16in .75 mm bore 100C/10'000psi SS
SFVO	On/Off valve manual	2 x 1/16in .50 mm bore 100C/10'000psi SS
SFVOHT	On/Off valve manual high temperature/pressure	2 x 1/16in .50 mm bore 300C/6000psi SS
SFVOHT4	On/Off valve manual high temp/pressure/4in so	2 x 1/16in .50 mm bore 100C/10'000psi SS
SFVOL	On/Off valve manual	2 x 1/16in .75 mm bore 100C/10'000psi SS
SI4G	Replacement septa for septum injector 350C	1/4in OD low bleed 10/pkg
SL100CSTP	Sample loop 100ul mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL100CSTUW	Sample loop 100ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL100CUW	Sample loop 100ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL100CW	Sample loop 100ul W type valve	1/16in ends with nuts & ferrules SS
SL100NW	Sample loop 100ul W type valve	1/32in ends with nuts & ferrules 316SS
SL100STP	Sample loop 100ul mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL100UW	Sample loop 100ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL10CSTUW	Sample loop 10ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL10CUW	Sample loop 10ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL10CW	Sample loop 10ul W type valve	1/16in ends with nuts & ferrules 316SS
SL10KCSTP	Sample loop 10ml mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL10KCSTUW	Sample loop 10ml mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL10KCUW	Sample loop 10ml UW type valve	1/16in ends with nuts & ferrules 316SS
SL10KCW	Sample loop 10ml W type valve	1/16in ends with nuts & ferrules 316SS
SL10KSTP	Sample loop 10ml mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL10KUW	Sample loop 10ml UW type valve	1/8in ends with nuts & ferrules 316SS
SL10NW	Sample loop 10ul W type valve	1/32in ends with nuts & ferrules 316SS
SL10UW	Sample loop 10ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL15CSTUW	Sample loop 15ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL15CUW	Sample loop 15ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL15CW	Sample loop 15ul W type valve	1/16in ends with nuts & ferrules 316SS
SL15NW	Sample loop 15ul W type valve	1/32in ends with nuts & ferrules 316SS
SL15UW	Sample loop 15ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL1KCSTP	Sample loop 1ml mpos ST/MT P-type valve	1/16in ends with nuts & ferrules 316SS
SL1KCSTUW	Sample loop 1ml mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL1KCUW	Sample loop 1ml UW type valve	1/16in ends with nuts & ferrules 316SS
SL1KCW	Sample loop 1ml W type valve	1/16in ends with nuts & ferrules 316SS
SL1KSTP	Sample loop 1ml mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL1KUW	Sample loop 1ml UW type valve	1/8in ends with nuts & ferrules 316SS
SL20CSTUW	Sample loop 20ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL20CUW	Sample loop 20ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL20CW	Sample loop 20ul W type valve	1/16in ends with nuts & ferrules 316SS
SL20KUW	Sample loop 20 ml UW type valve	1/8in ends with nuts & ferrules 316SS
SL20NW	Sample loop 20ul W type valve	1/32in ends with nuts & ferrules 316SS
SL20UW	Sample loop 20ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL250CSTP	Sample loop 250ul mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL250CSTUW	Sample loop 250ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL250CUW	Sample loop 250ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL250CW	Sample loop 250ul W type valve	1/16in ends with nuts & ferrules 316SS
SL250NW	Sample loop 250ul W type valve	1/32in ends with nuts & ferrules 316SS
SL250STP	Sample loop 250ul mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL250UW	Sample loop 250ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL25CSTUW	Sample loop 25ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL25CUW	Sample loop 25ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL25CW	Sample loop 25ul W type valve	1/16in ends with nuts & ferrules 316SS
SL25NW	Sample loop 25ul W type valve	1/32in ends with nuts & ferrules 316SS
SL25UW	Sample loop 25ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL2CW	Sample loop 2ul W type valve	1/16in ends with nuts & ferrules 316SS
SL2KCSTP	Sample loop 2ml mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL2KCSTUW	Sample loop 2ml mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL2KCUW	Sample loop 2ml UW type valve	1/16in ends with nuts & ferrules 316SS
SL2KCW	Sample loop 2ml W type valve	1/16in ends with nuts & ferrules 316SS
SL2KSTP	Sample loop 2ml mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL2KUW	Sample loop 2ml UW type valve	1/8in ends with nuts & ferrules 316SS
SL2NW	Sample loop 2ul W type valve	1/32in ends with nuts & ferrules 316SS

# VICI Valco-Cheminert PRODUCT LIST 2012

SL3CUW	Sample loop 3ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL500CSTP	Sample loop 500ul mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL500CSTUW	Sample loop 500ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL500CUW	Sample loop 500ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL500CW	Sample loop 500ul W type valve	1/16in ends with nuts & ferrules 316SS
SL500NW	Sample loop 500ul W type valve	1/32in ends with nuts & ferrules 316SS
SL500STP	Sample loop 500ul mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL500UW	Sample loop 500ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL50CSTP	Sample loop 50ul mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL50CSTUW	Sample loop 50ul mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL50CUW	Sample loop 50ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL50CW	Sample loop 50ul W type valve	1/16in ends with nuts & ferrules 316SS
SL50NW	Sample loop 50ul W type valve	1/32in ends with nuts & ferrules 316SS
SL50UW	Sample loop 50ul UW type valve	1/8in ends with nuts & ferrules 316SS
SL5CUW	Sample loop 5ul UW type valve	1/16in ends with nuts & ferrules 316SS
SL5CW	Sample loop 5ul W type valve	1/16in ends with nuts & ferrules 316SS
SL5KCSTP	Sample loop 5ml mpos ST/MW P-type valve	1/16in ends with nuts & ferrules 316SS
SL5KCSTUW	Sample loop 5ml mpos ST UW-type valve	1/16in ends with nuts & ferrules 316SS
SL5KCUW	Sample loop 5ml UW type valve	1/16in ends with nuts & ferrules 316SS
SL5KCW	Sample loop 5ml W type valve	1/16in ends with nuts & ferrules 316SS
SL5KSTP	Sample loop 5ml mpos ST/MW P-type valve	1/8in ends with nuts & ferrules 316SS
SL5KUW	Sample loop 5ml UW type valve	1/8in ends with nuts & ferrules 316SS
SL5NW	Sample loop 5ul W type valve	1/32in ends with nuts & ferrules 316SS
SOMAP	Manual adaptor for standoff	incl. MHA and SVH
SOMAPPF	Manual adaptor with Position feedback	for standoff incl. MHA and SVH
SOMAW	Manual adaptor for standoff	incl. MHA and WK
SOUTH	Handle for use with standoff	2-pos. UW-type valve
SPSS40	Ultra High Pressure Injection System	40K psi liq
SSA10UW	Rotor seal 10 port 2-pos UW-valve	1/8 x .75mm 75C/5000psi liq H
SSA10UWE	Rotor seal 10 port 2-pos UW-valve	1/8 x .75mm 225C/400psi gas E
SSA10UWT	Rotor seal 10 port 2-pos UW-valve	1/8 x .75mm 330C/300psi gas T
SSA4UW	Rotor seal 4 port 2-pos UW-valve	1/8 x .75mm 75C/5000psi liq H
SSA4UWE	Rotor seal 4 port 2-pos UW-valve	1/8 x .75mm 225C/400psi gas E
SSA4UWT	Rotor seal 4 port 2-pos UW-valve	1/8 x .75mm 330C/300psi gas T
SSA6UW	Rotor seal 6 port 2-pos UW-valve	1/8 x .75mm 75C/5000psi liq H
SSA6UWE	Rotor seal 6 port 2-pos UW-valve	1/8 x .75mm 225C/400psi gas E
SSA6UWT	Rotor seal 6 port 2-pos UW-valve	1/8 x .75mm 330C/300psi gas T
SSA8UW	Rotor seal 8 port 2-pos UW-valve	1/8 x .75mm 75C/5000psi liq H
SSA8UWE	Rotor seal 8 port 2-pos UW-valve	1/8 x .75mm 225C/400psi gas E
SSA8UWT	Rotor seal 8 port 2-pos UW-valve	1/8 x .75mm 330C/300psi gas T
SSAC10UW	Rotor seal 10 port 2-pos UW-valve	1/16 x .75mm 75C/5000psi liq H
SSAC10UWE	Rotor seal 10 port 2-pos UW-valve	1/16 x .75mm 200C/400psi gas E
SSAC10UWT	Rotor seal 10 port 2-pos UW-valve	1/16 x .75mm 330C/300psi gas T
SSAC10W	Rotor seal 10 port 2-pos W-valve	1/16 x .4mm 75C/5000psi liq H
SSAC10WE	Rotor seal 10 port 2-pos W-valve	1/16 x .4mm 225C/400psi gas E
SSAC10WT	Rotor seal 10 port 2-pos W-valve	1/16 x .4mm 350C/300psi gas T
SSAC4UW	Rotor seal 4 port 2-pos UW-valve	1/16 x .75mm 75C/5000psi liq H
SSAC4UWE	Rotor seal 4 port 2-pos UW-valve	1/16 x .75mm 225C/400psi gas E
SSAC4UWT	Rotor seal 4 port 2-pos UW-valve	1/16 x .75mm 330C/300psi gas T
SSAC4W	Rotor seal 4 port 2-pos W-valve	1/16 x .4mm 75C/5000psi liq H
SSAC4WE	Rotor seal 4 port 2-pos W-valve	1/16 x .4mm 225C/400psi gas E
SSAC4WT	Rotor seal 4 port 2-pos W-valve	1/16 x .4mm 350C/300psi gas T
SSAC6UW	Rotor seal 6 port 2-pos UW-valve	1/16 x .75mm 75C/5000psi liq H
SSAC6UWE	Rotor seal 6 port 2-pos UW-valve	1/16 x .75mm 225C/400psi gas E
SSAC6UWT	Rotor seal 6 port 2-pos UW-valve	1/16 x .75mm 330C/300psi gas T
SSAC6W	Rotor seal 6 port 2-pos W-valve	1/16 x .4mm 75C/5000psi liq H
SSAC6WE	Rotor seal 6 port 2-pos W-valve	1/16 x .4mm 225C/400psi gas E
SSAC6WT	Rotor seal 6 port 2-pos W-valve	1/16 x .4mm 350C/300psi gas T
SSAC8UW	Rotor seal 8 port 2-pos UW-valve	1/16 x .75mm 75C/5000psi liq H
SSAC8UWE	Rotor seal 8 port 2-pos UW-valve	1/16 x .75mm 225C/400psi gas E
SSAC8UWT	Rotor seal 8 port 2-pos UW-valve	1/16 x .75mm 330C/300psi gas T




# VICI Valco-Cheminert PRODUCT LIST 2012




SSAC8W	Rotor seal 8 port 2-pos W-valve	1/16 x .4mm 75C/5000psi liq H
SSAC8WE	Rotor seal 8 port 2-pos W-valve	1/16 x .4mm 225C/400psi gas E
SSAC8WT	Rotor seal 8 port 2-pos valve	1/16 x .4mm 350C/300psi gas T
SSACI4UW.2	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 75C/5000psi liq .2ul H
SSACI4UW.5	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 75C/5000psi liq .5ul H
SSACI4UW1	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 75C/5000psi liq 1ul H
SSACI4UW2	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 75C/5000psi liq 2ul H
SSACI4UWE.2	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 175C/1000psi liq .2ul E
SSACI4UWE.5	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 175C/1000psi liq .5ul E
SSACI4UWE1	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 175C/1000psi liq 1ul E
SSACI4UWE2	Rotor seal 4 port 2-pos int vol UW-valve	1/16 x .75mm 175C/1000psi liq 2ul N60/E
SSACI4W.06	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 75C/5000psi liq .06ul H
SSACI4W.1	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 75C/5000psi liq .1ul H
SSACI4W.2	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 75C/5000psi liq .2ul H
SSACI4W.5	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 75C/5000psi liq .5ul H
SSACI4WE.06	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 175C/1000psi liq .06ul E
SSACI4WE.1	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 175C/1000psi liq .1ul E
SSACI4WE.2	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 175C/1000psi liq .2ul E
SSACI4WE.5	Rotor seal 4 port 2-pos int vol W-valve	1/16 x .4mm 175C/1000psi liq .5ul E
SSACSC10MWE	Rotor seal 10 mpos common outlet MW valve	1/16 x 1mm 200C/200psi gas N60/E
SSACSC12MWE	Rotor seal 12 mpos common outlet MW valve	1/16 x 1mm 200C/200psi gas N60/E
SSACSC16MWE	Rotor seal 16 mpos common outlet MW valve	1/16 x 1mm 200C/200psi gas N60/E
SSACSC6MWE	Rotor seal 6 mpos common outlet MW valve	1/16 x 1mm 200C/200psi gas N60/E
SSACSD10MWE	Rotor seal 10 mpos dead-end path MW-valve	1/16 x .75mm 200C/400psi gas E
SSACSD10UW	Rotor seal 10 mpos dead-end path UW-valve	1/16 x .4mm 75C/5000psi liq E
SSACSD12MWE	Rotor seal 12 mpos dead-end path MW-valve	1/16 x .75mm 200C/400psi gas E
SSACSD16MWE	Rotor seal 16 mpos dead-end path MW-valve	1/16 x .75mm 200C/400psi gas E
SSACSD4UW	Rotor seal 4 mpos dead-end path UW-valve	1/16 x .4mm 75C/5000psi liq E
SSACSD6MWE	Rotor seal 6 mpos dead-end path MW-valve	1/16 x .75mm 200C/400psi gas E
SSACSD6UW	Rotor seal 6 mpos dead-end path UW-valve	1/16 x .4mm 75C/5000psi liq E
SSACSF10MWE	Rotor seal 10 mpos flow-through MW-valve	1/16 x 1mm 200C/200psi gas E
SSACSF12MWE	Rotor seal 12 mpos flow-through MW-valve	1/16 x 1mm 200C/200psi gas E
SSACSF16MWE	Rotor seal 16 mpos flow-through MW-valve	1/16 x 1mm 200C/200psi gas E
SSACSF6MWE	Rotor seal 6 mpos flow-through MW-valve	1/16 x 1mm 200C/200psi gas E
SSACST10MWE	Rotor seal 10 mpos trapping path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACST12MWE	Rotor seal 12 mpos trapping path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACST16MWE	Rotor seal 16 mpos trapping path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACST4UW	Rotor seal 4 mpos trapping-path UW-valve	1/16 x .4mm 75C/5000psi liq E
SSACST6MWE	Rotor seal 6 mpos trapping path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACST6UW	Rotor seal 6 mpos trapping-path UW-valve	1/16 x .4mm 75C/5000psi liq E
SSACSTF10MWE	Rotor seal 10 mpos trap/flow path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACSTF12MWE	Rotor seal 12 mpos trap/flow path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACSTF16MWE	Rotor seal 16 mpos trap/flow path MW-valve	1/16 x .75mm 200C/200psi gas E
SSACSTF6MWE	Rotor seal 6 mpos trap/flow path MW-valve	1/16 x .75mm 200C/200psi gas E
SSAI4UWE.2	Rotor seal 4 port 2-pos int vol UW-valve	1/8 x .75mm 175C/1000psi liq .2ul E
SSAI4UWE.5	Rotor seal 4 port 2-pos int vol UW-valve	1/8 x .75mm 175C/1000psi liq .5ul E
SSAI4UWE1	Rotor seal 4 port 2-pos int vol UW-valve	1/8 x .75mm 175C/1000psi liq 1ul E
SSAI4UWE2	Rotor seal 4 port 2-pos int vol UW-valve	1/8 x .75mm 175C/1000psi liq 2ul E
SSAL10UW	Rotor seal 10 port 2-pos UW-valve	1/8 x 1.0mm 75C/5000psi liq H
SSAL4UW	Rotor seal 4 port 2-pos UW-valve	1/8 x 1.7mm 75C/5000psi liq H
SSAL6UW	Rotor seal 6 port 2-pos UW-valve	1/8 x 1.7mm 75C/5000psi liq H
SSAL8UW	Rotor seal 8 port 2-pos UW-valve	1/8 x 1.3mm 75C/5000psi liq H
SSAN10WE	Rotor seal 10 port 2-pos W-valve	1/32 x .25mm 225C/400psi gas E
SSAN10WT	Rotor seal 10 port 2-pos W-valve	1/32 x .25mm 300C/300psi gas T
SSAN4WE	Rotor seal 4 port 2-pos W-valve	1/32 x .25mm 225C/400psi gas E
SSAN4WT	Rotor seal 4 port 2-pos W-valve	1/32 x .25mm 350C/300psi gas T
SSAN6WE	Rotor seal 6 port 2-pos W-valve	1/32 x .25mm 250C/400psi gas E
SSAN6WT	Rotor seal 6 port 2-pos W-valve	1/32 x .25mm 350C/300psi gas T
SSAN8WE	Rotor seal 8 port 2-pos W-valve	1/32 x .25mm 225C/400psi gas E
SSAN8WT	Rotor seal 8 port 2-pos W-valve	1/32 x .25mm 350C/300psi gas T
SSANI4WE.06	Rotor seal 4 port 2-pos int vol W-valve	1/32 x .25mm 175C/1000psi liq .06ul E

# VICI Valco-Cheminert PRODUCT LIST 2012

SSANI4WE.1	Rotor seal 4 port 2-pos int vol W-valve	1/32 x .25mm 175C/1000psi liq .1ul E
SSANI4WE.2	Rotor seal 4 port 2-pos int vol W-valve	1/32 x .25mm 175C/1000psi liq .2ul E
SSANI4WE.5	Rotor seal 4 port 2-pos int vol W-valve	1/32 x .25mm 175C/1000psi liq .5 ul E
SSASC10MWE	Rotor seal 10 mpos common outlet MW valve	1/8 x 1mm 200C/200psi gas E
SSASC12MWE	Rotor seal 12 mpos common outlet MW valve	1/8 x 1mm 200C/200psi gas E
SSASC16MWE	Rotor seal 16 mpos common outlet MW valve	1/8 x 1mm 200C/200psi gas E
SSASC6MWE	Rotor seal 6 mpos common outlet MW valve	1/8 x 1mm 200C/200psi gas E
SSASD10MWE	Rotor seal 10 mpos dead-end path MW-valve	1/8 x 1mm 200C/200psi gas E
SSASD12MWE	Rotor seal 12 mpos dead-end path MW-valve	1/8 x 1mm 200C/200psi gas E
SSASD16MWE	Rotor seal 16 mpos dead-end path MW-valve	1/8 x 1mm 200C/200psi gas E
SSASD4UW	Rotor seal 4 mpos dead-end path UW-valve	1/8 x .75mm 75C/5000psi liq E
SSASD6MWE	Rotor seal 6 mpos dead-end path MW-valve	1/8 x 1mm 200C/400psi gas E
SSASD6UW	Rotor seal 6 mpos dead-end path UW-valve	1/8 x .75mm 75C/5000psi liq E
SSASD8UW	Rotor seal 8 mpos dead-end path UW-valve	1/8 x .75mm 75C/5000psi liq E
SSASF10MWE	Rotor seal 10 mpos flow-through MW-valve	1/8 x 1mm 200C/200psi gas E
SSASF12MWE	Rotor seal 12 mpos flow-through MW-valve	1/8 x 1mm 200C/200psi gas E
SSASF16MWE	Rotor seal 16 mpos flow-through MW-valve	1/8 x 1mm 200C/200psi gas E
SSASF6MWE	Rotor seal 6 mpos flow-through MW-valve	1/8 x 1mm 200C/200psi gas E
SSAST10MWE	Rotor seal 10 mpos trapping path MW-valve	1/8 x 1mm 200C/200psi gas E
SSAST12MWE	Rotor seal 12 mpos trapping path MW-valve	1/8 x 1mm 200C/200psi gas E
SSAST16MWE	Rotor seal 16 mpos trapping path MW-valve	1/8 x 1mm 200C/200psi gas E
SSAST6MWE	Rotor seal 6 mpos trapping path MW-valve	1/8 x 1mm 200C/200psi gas E
SSASTF10MWE	Rotor seal 10 mpos trap/flow path MWvalve	1/8 x 1mm 200C/200psi gas E
SSASTF12MWE	Rotor seal 12 mpos trap/flow path MWvalve	1/8 x 1mm 200C/200psi gas E
SSASTF16MWE	Rotor seal 16 mpos trap/flow path MWvalve	1/8 x 1mm 200C/200psi gas E
SSASTF6MWE	Rotor seal 6 mpos trap/flow path MW-valve	1/8 x 1mm 200C/200psi gas E
SSAVL4MWE2	Rotor seal 4 port 2 pos MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVL6MWE2	Rotor seal 6 port 2 pos MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVL8MWE2	Rotor seal 8 port 2 pos MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSC4MWE2	Rotor seal 4 mpos common outlet MW valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSC6MWE2	Rotor seal 6 mpos common outlet MW valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSC8MWE2	Rotor seal 8 mpos common outlet MW valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSD10MWE2	Rotor seal 10 mpos dead-end path MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSD4MWE2	Rotor seal 4 mpos dead-end path MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSD6MWE2	Rotor seal 6 mpos dead-end path MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSD8MWE2	Rotor seal 8 mpos dead-end path MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSF4MWE2	Rotor seal 4 mpos flow-through MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSF6MWE2	Rotor seal 6 mpos flow-through MW-valve	1/4 x 4mm 75C/100psi gas E2
SSAVLSF8MWE2	Rotor seal 8 mpos flow-through MW-valve	1/4 x 4mm 75C/100psi gas E2
SVI-110	Serial valve interface 110Vac	Non-CE !!!
SVI-220	Serial valve interface 220Vac	Non-CE !!!
SWH4	Socket wrench handle 1/4in	for 1/16in nuts
T100C10-10	Tube-kit 1/16in .010in /.25mm	100 cm long 316SS 10/pkg
T100C10-100	Tube-kit 1/16in .010in /.25mm	100 cm long 316SS 100/pkg
T100C10-50	Tube-kit 1/16in .010in /.25mm	100 cm long 316SS 50/pkg
T100C10D	Tube-kit 1/16in .010in /.25mm	100 cm long 316SS 2/pkg
T100C20-10	Tube-kit 1/16in .020in /.50mm	100 cm long 316SS 10/pkg
T100C20-100	Tube-kit 1/16in .020in /.50mm	100 cm long 316SS 100/pkg
T100C20-50	Tube-kit 1/16in .020in /.50mm	100 cm long 316SS 50/pkg
T100C20D	Tube-kit 1/16in .020in /.50mm	100 cm long 316SS 2/pkg
T100C30-10	Tube-kit 1/16in .030in /.75mm	100 cm long 316SS 10/pkg
T100C30-100	Tube-kit 1/16in .030in /.75mm	100 cm long 316SS 100/pkg
T100C30-50	Tube-kit 1/16in .030in /.75mm	100 cm long 316SS 50/pkg
T100C30D	Tube-kit 1/16in .030in /.75mm	100 cm long 316SS 2/pkg
T100C40-10	Tube-kit 1/16in .040in /1.0mm	100 cm long 316SS 10/pkg
T100C40-100	Tube-kit 1/16in .040in /1.0mm	100 cm long 316SS 100/pkg
T100C40-50	Tube-kit 1/16in .040in /1.0mm	100 cm long 316SS 50/pkg
T100C40D	Tube-kit 1/16in .040in /1.0mm	100 cm long 316SS 2/pkg
T100N10-10	Tube-kit 1/32in .010in /.25mm	100 cm long 316SS 10/pkg
T100N10-100	Tube-kit 1/32in .010in /.25mm	100 cm long 316SS 100/pkg
T100N10-50	Tube-kit 1/32in .010in /.25mm	100 cm long 316SS 50/pkg




# VICI Valco-Cheminert PRODUCT LIST 2012




T100N10D	Tube-kit 1/32in .010in /.25mm	100 cm long 316SS 2/pkg
T100N20-10	Tube-kit 1/32in .020in /.50mm	100 cm long 316SS 10/pkg
T100N20-100	Tube-kit 1/32in .020in /.50mm	100 cm long 316SS 100/pkg
T100N20-50	Tube-kit 1/32in .020in /.50mm	100 cm long 316SS 50/pkg
T100N20D	Tube-kit 1/32in .020in /.50mm	100 cm long 316SS 2/pkg
T10C10-10	Tube-kit 1/16in .010in /.25mm	10 cm long 316SS 10/pkg
T10C10-100	Tube-kit 1/16in .010in /.25mm	10 cm long 316SS 100/pkg
T10C10-50	Tube-kit 1/16in .010in /.25mm	10 cm long 316SS 50/pkg
T10C10D	Tube-kit 1/16in .010in /.25mm	10 cm long 316SS 2/pkg
T10C20-10	Tube-kit 1/16in .020in /.50mm	10 cm long 316SS 10/pkg
T10C20-100	Tube-kit 1/16in .020in /.50mm	10 cm long 316SS 100/pkg
T10C20-50	Tube-kit 1/16in .020in /.50mm	10 cm long 316SS 50/pkg
T10C20D	Tube-kit 1/16in .020in /.50mm	10 cm long 316SS 2/pkg
T10C30-10	Tube-kit 1/16in .030in /.75mm	10 cm long 316SS 10/pkg
T10C30-100	Tube-kit 1/16in .030in /.75mm	10 cm long 316SS 100/pkg
T10C30-50	Tube-kit 1/16in .030in /.75mm	10 cm long 316SS 50/pkg
T10C30D	Tube-kit 1/16in .030in /.75mm	10 cm long 316SS 2/pkg
T10C40-10	Tube-kit 1/16in .040in /1.0mm	10 cm long 316SS 10/pkg
T10C40-100	Tube-kit 1/16in .040in /1.0mm	10 cm long 316SS 100/pkg
T10C40-50	Tube-kit 1/16in .040in /1.0mm	10 cm long 316SS 50/pkg
T10C40D	Tube-kit 1/16in .040in /1.0mm	10 cm long 316SS 2/pkg
T10C5-10	Tube-kit 1/16in .005in /.12mm	10 cm long 316SS 10/pkg
T10C5-100	Tube-kit 1/16in .005in /.12mm	10 cm long 316SS 100/pkg
T10C5-50	Tube-kit 1/16in .005in /.12mm	10 cm long 316SS 50/pkg
T10C5D	Tube-kit 1/16in .005in /.12mm	10 cm long 316SS 2/pkg
T10N10-10	Tube-kit 1/32in .010in /.25mm	10 cm long 316SS 10/pkg
T10N10-100	Tube-kit 1/32in .010in /.25mm	10 cm long 316SS 100/pkg
T10N10-50	Tube-kit 1/32in .010in /.25mm	10 cm long 316SS 50/pkg
T10N10D	Tube-kit 1/32in .010in /.25mm	10 cm long 316SS 2/pkg
T10N20-10	Tube-kit 1/32in .020in /.50mm	10 cm long 316SS 10/pkg
T10N20-100	Tube-kit 1/32in .020in /.50mm	10 cm long 316SS 100/pkg
T10N20-50	Tube-kit 1/32in .020in /.50mm	10 cm long 316SS 50/pkg
T10N20D	Tube-kit 1/32in .020in /.50mm	10 cm long 316SS 2/pkg
T10N5-10	Tube-kit 1/32in .005in /.12mm	10 cm long 316SS 10/pkg
T10N5-100	Tube-kit 1/32in .005in /.12mm	10 cm long 316SS 100/pkg
T10N5-50	Tube-kit 1/32in .005in /.12mm	10 cm long 316SS 50/pkg
T10N5D	Tube-kit 1/32in .005in /.12mm	10 cm long 316SS 2/pkg
T20C10-10	Tube-kit 1/16in .010in /.25mm	20 cm long 316SS 10/pkg
T20C10-100	Tube-kit 1/16in .010in /.25mm	20 cm long 316SS 100/pkg
T20C10-50	Tube-kit 1/16in .010in /.25mm	20 cm long 316SS 50/pkg
T20C10D	Tube-kit 1/16in .010in /.25mm	20 cm long 316SS 2/pkg
T20C20-10	Tube-kit 1/16in .020in /.50mm	20 cm long 316SS 10/pkg
T20C20-100	Tube-kit 1/16in .020in /.50mm	20 cm long 316SS 100/pkg
T20C20-50	Tube-kit 1/16in .020in /.50mm	20 cm long 316SS 50/pkg
T20C20D	Tube-kit 1/16in .020in /.50mm	20 cm long 316SS 2/pkg
T20C30-10	Tube-kit 1/16in .030in /.75mm	20 cm long 316SS 10/pkg
T20C30-100	Tube-kit 1/16in .030in /.75mm	20 cm long 316SS 100/pkg
T20C30-50	Tube-kit 1/16in .030in /.75mm	20 cm long 316SS 50/pkg
T20C30D	Tube-kit 1/16in .030in /.75mm	20 cm long 316SS 2/pkg
T20C40-10	Tube-kit 1/16in .040in /1.0mm	20 cm long 316SS 10/pkg
T20C40-100	Tube-kit 1/16in .040in /1.0mm	20 cm long 316SS 100/pkg
T20C40-50	Tube-kit 1/16in .040in /1.0mm	20 cm long 316SS 50/pkg
T20C40D	Tube-kit 1/16in .040in /1.0mm	20 cm long 316SS 2/pkg
T20C5-10	Tube-kit 1/16in .005in /.12mm	20 cm long 316SS 10/pkg
T20C5-100	Tube-kit 1/16in .005in /.12mm	20 cm long 316SS 100/pkg
T20C5-50	Tube-kit 1/16in .005in /.12mm	20 cm long 316SS 50/pkg
T20C5D	Tube-kit 1/16in .005in /.12mm	20 cm long 316SS 2/pkg
T20N10-10	Tube-kit 1/32in .010in /.25mm	20 cm long 316SS 10/pkg
T20N10-100	Tube-kit 1/32in .010in /.25mm	20 cm long 316SS 100/pkg
T20N10-50	Tube-kit 1/32in .010in /.25mm	20 cm long 316SS 50/pkg
T20N10D	Tube-kit 1/32in .010in /.25mm	20 cm long 316SS 2/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




T20N20-10	Tube-kit 1/32in .020in /.50mm	20 cm long 316SS 10/pkg
T20N20-100	Tube-kit 1/32in .020in /.50mm	20 cm long 316SS 100/pkg
T20N20-50	Tube-kit 1/32in .020in /.50mm	20 cm long 316SS 50/pkg
T20N20D	Tube-kit 1/32in .020in /.50mm	20 cm long 316SS 2/pkg
T20N5-10	Tube-kit 1/32in .005in /.12mm	20 cm long 316SS 10/pkg
T20N5-100	Tube-kit 1/32in .005in /.12mm	20 cm long 316SS 100/pkg
T20N5-50	Tube-kit 1/32in .005in /.12mm	20 cm long 316SS 50/pkg
T20N5D	Tube-kit 1/32in .005in /.12mm	20 cm long 316SS 2/pkg
T30C10-10	Tube-kit 1/16in .010in /.25mm	30 cm long 316SS 10/pkg
T30C10-100	Tube-kit 1/16in .010in /.25mm	30 cm long 316SS 100/pkg
T30C10-50	Tube-kit 1/16in .010in /.25mm	30 cm long 316SS 50/pkg
T30C10D	Tube-kit 1/16in .010in /.25mm	30 cm long 316SS 2/pkg
T30C20-10	Tube-kit 1/16in .020in /.50mm	30 cm long 316SS 10/pkg
T30C20-100	Tube-kit 1/16in .020in /.50mm	30 cm long 316SS 100/pkg
T30C20-50	Tube-kit 1/16in .020in /.50mm	30 cm long 316SS 50/pkg
T30C20D	Tube-kit 1/16in .020in /.50mm	30 cm long 316SS 2/pkg
T30C30-10	Tube-kit 1/16in .030in /.75mm	30 cm long 316SS 10/pkg
T30C30-100	Tube-kit 1/16in .030in /.75mm	30 cm long 316SS 100/pkg
T30C30-50	Tube-kit 1/16in .030in /.75mm	30 cm long 316SS 50/pkg
T30C30D	Tube-kit 1/16in .030in /.75mm	30 cm long 316SS 2/pkg
T30C40-10	Tube-kit 1/16in .040in /1.0mm	30 cm long 316SS 10/pkg
T30C40-100	Tube-kit 1/16in .040in /1.0mm	30 cm long 316SS 100/pkg
T30C40-50	Tube-kit 1/16in .040in /1.0mm	30 cm long 316SS 50/pkg
T30C40D	Tube-kit 1/16in .040in /1.0mm	30 cm long 316SS 2/pkg
T30C5-10	Tube-kit 1/16in .005in /.12mm	30 cm long 316SS 10/pkg
T30C5-100	Tube-kit 1/16in .005in /.12mm	30 cm long 316SS 100/pkg
T30C5-50	Tube-kit 1/16in .005in /.12mm	30 cm long 316SS 50/pkg
T30C5D	Tube-kit 1/16in .005in /.12mm	30 cm long 316SS 2/pkg
T30N10-10	Tube-kit 1/32in .010in /.25mm	30 cm long 316SS 10/pkg
T30N10-100	Tube-kit 1/32in .010in /.25mm	30 cm long 316SS 100/pkg
T30N10-50	Tube-kit 1/32in .010in /.25mm	30 cm long 316SS 50/pkg
T30N10D	Tube-kit 1/32in .010in /.25mm	30 cm long 316SS 2/pkg
T30N20-10	Tube-kit 1/32in .020in /.50mm	30 cm long 316SS 10/pkg
T30N20-100	Tube-kit 1/32in .020in /.50mm	30 cm long 316SS 100/pkg
T30N20-50	Tube-kit 1/32in .020in /.50mm	30 cm long 316SS 50/pkg
T30N20D	Tube-kit 1/32in .020in /.50mm	30 cm long 316SS 2/pkg
T30N5-10	Tube-kit 1/32in .005in /.12mm	30 cm long 316SS 10/pkg
T30N5-100	Tube-kit 1/32in .005in /.12mm	30 cm long 316SS 100/pkg
T30N5-50	Tube-kit 1/32in .005in /.12mm	30 cm long 316SS 50/pkg
T30N5D	Tube-kit 1/32in .005in /.12mm	30 cm long 316SS 2/pkg
T50C10-10	Tube-kit 1/16in .010in /.25mm	50 cm long 316SS 10/pkg
T50C10-100	Tube-kit 1/16in .010in /.25mm	50 cm long 316SS 100/pkg
T50C10-50	Tube-kit 1/16in .010in /.25mm	50 cm long 316SS 50/pkg
T50C10D	Tube-kit 1/16in .010in /.25mm	50 cm long 316SS 2/pkg
T50C20-10	Tube-kit 1/16in .020in /.50mm	50 cm long 316SS 10/pkg
T50C20-100	Tube-kit 1/16in .020in /.50mm	50 cm long 316SS 100/pkg
T50C20-50	Tube-kit 1/16in .020in /.50mm	50 cm long 316SS 50/pkg
T50C20D	Tube-kit 1/16in .020in /.50mm	50 cm long 316SS 2/pkg
T50C30-10	Tube-kit 1/16in .030in /.75mm	50 cm long 316SS 10/pkg
T50C30-100	Tube-kit 1/16in .030in /.75mm	50 cm long 316SS 100/pkg
T50C30-50	Tube-kit 1/16in .030in /.75mm	50 cm long 316SS 50/pkg
T50C30D	Tube-kit 1/16in .030in /.75mm	50 cm long 316SS 2/pkg
T50C40-10	Tube-kit 1/16in .040in /1.0mm	50 cm long 316SS 10/pkg
T50C40-100	Tube-kit 1/16in .040in /1.0mm	50 cm long 316SS 100/pkg
T50C40-50	Tube-kit 1/16in .040in /1.0mm	50 cm long 316SS 50/pkg
T50C40D	Tube-kit 1/16in .040in /1.0mm	50 cm long 316SS 2/pkg
T50C5-10	Tube-kit 1/16in .005in /.12mm	50 cm long 316SS 10/pkg
T50C5-100	Tube-kit 1/16in .005in /.12mm	50 cm long 316SS 100/pkg
T50C5-50	Tube-kit 1/16in .005in /.12mm	50 cm long 316SS 50/pkg
T50C5D	Tube-kit 1/16in .005in /.12mm	50 cm long 316SS 2/pkg
T50N10-10	Tube-kit 1/32in .010in /.25mm	50 cm long 316SS 10/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




T50N10-100	Tube-kit 1/32in .010in /.25mm	50 cm long 316SS 100/pkg
T50N10-50	Tube-kit 1/32in .010in /.25mm	50 cm long 316SS 50/pkg
T50N10D	Tube-kit 1/32in .010in /.25mm	50 cm long 316SS 2/pkg
T50N20-10	Tube-kit 1/32in .020in /.50mm	50 cm long 316SS 10/pkg
T50N20-100	Tube-kit 1/32in .020in /.50mm	50 cm long 316SS 100/pkg
T50N20-50	Tube-kit 1/32in .020in /.50mm	50 cm long 316SS 50/pkg
T50N20D	Tube-kit 1/32in .020in /.50mm	50 cm long 316SS 2/pkg
T50N5-10	Tube-kit 1/32in .005in /.12mm	50 cm long 316SS 10/pkg
T50N5-100	Tube-kit 1/32in .005in /.12mm	50 cm long 316SS 100/pkg
T50N5-50	Tube-kit 1/32in .005in /.12mm	50 cm long 316SS 50/pkg
T50N5D	Tube-kit 1/32in .005in /.12mm	50 cm long 316SS 2/pkg
T5C10-10	Tube-kit 1/16in .010in /.25mm	5 cm long 316SS 10/pkg
T5C10-100	Tube-kit 1/16in .010in /.25mm	5 cm long 316SS 100/pkg
T5C10-50	Tube-kit 1/16in .010in /.25mm	5 cm long 316SS 50/pkg
T5C10D	Tube-kit 1/16in .010in /.25mm	5 cm long 316SS 2/pkg
T5C20-10	Tube-kit 1/16in .020in /.50mm	5 cm long 316SS 10/pkg
T5C20-100	Tube-kit 1/16in .020in /.50mm	5 cm long 316SS 100/pkg
T5C20-50	Tube-kit 1/16in .020in /.50mm	5 cm long 316SS 50/pkg
T5C20D	Tube-kit 1/16in .020in /.50mm	5 cm long 316SS 2/pkg
T5C30-10	Tube-kit 1/16in .030in /.75mm	5cm long 316SS 10/pkg
T5C30-100	Tube-kit 1/16in .030in /.75mm	5 cm long 316SS 100/pkg
T5C30-50	Tube-kit 1/16in .030in /.75mm	5 cm long 316SS 50/pkg
T5C30D	Tube-kit 1/16in .030in /.75mm	5 cm long 316SS 2/pkg
T5C40-10	Tube-kit 1/16in .040in /1.0mm	5 cm long 316SS 10/pkg
T5C40-100	Tube-kit 1/16in .040in /1.0mm	5 cm long 316SS 100/pkg
T5C40-50	Tube-kit 1/16in .040in /1.0mm	5 cm long 316SS 50/pkg
T5C40D	Tube-kit 1/16in .040in /1.0mm	5 cm long 316SS 2/pkg
T5C5-10	Tube-kit 1/16in .005in /.12mm	5 cm long 316SS 10/pkg
T5C5-100	Tube-kit 1/16in .005in /.12mm	5 cm long 316SS 100/pkg
T5C5-50	Tube-kit 1/16in .005in /.12mm	5 cm long 316SS 50/pkg
T5C5D	Tube-kit 1/16in .005in /.12mm	5 cm long 316SS 2/pkg
T5N10-10	Tube-kit 1/32in .010in /.25mm	5 cm long 316SS 10/pkg
T5N10-100	Tube-kit 1/32in .010in /.25mm	5 cm long 316SS 100/pkg
T5N10-50	Tube-kit 1/32in .010in /.25mm	5 cm long 316SS 50/pkg
T5N10D	Tube-kit 1/32in .010in /.25mm	5 cm long 316SS 2/pkg
T5N20-10	Tube-kit 1/32in .020in /.50mm	5 cm long 316SS 10/pkg
T5N20-100	Tube-kit 1/32in .020in /.50mm	5 cm long 316SS 100/pkg
T5N20-50	Tube-kit 1/32in .020in /.50mm	5 cm long 316SS 50/pkg
T5N20D	Tube-kit 1/32in .020in /.50mm	5 cm long 316SS 2/pkg
T5N5-10	Tube-kit 1/32in .005in /.12mm	5 cm long 316SS 10/pkg
T5N5-100	Tube-kit 1/32in .005in /.12mm	5 cm long 316SS 100/pkg
T5N5-50	Tube-kit 1/32in .005in /.12mm	5 cm long 316SS 50/pkg
T5N5D	Tube-kit 1/32in .005in /.12mm	5 cm long 316SS 2/pkg
TCD2-C-110	TCD2 controller	110Vac
TCD2-C-220	TCD2 controller	220Vac
TCD2-NIFE-110	TCD2 detector with controller	Nickel-Iron filaments 110Vac
TCD2-NIFE-220	TCD2 detector with controller	Nickel-Iron filaments 220Vac
TCD2-NIFED-110	TCD2 detector w/o controller	NIFE filaments 110Vac
TCD2-NIFED-220	TCD2 detector w/o controller	NIFE filaments 220Vac
TCD2-NIFE-FA	TCD2 filament assembly set of 2	Nickel/Iron NOT AVAILABLE see bemerkung
TCEF1.5.5C	Tee column end fitting 1/16in x 1/32in x 1/32in	.25mm in/out .25mm side 1mm cone SS
TCEF1.5.5T	Tee column end fitting 1/16in x 1/32in x 1/32in	.9mm in/out .25mm side 1mm cone SS
TCEF111	Tee column end fitting 1/16in x 1/16in x 1/16in	.25mm in/out .25mm side 1mm cone SS
TCEF211	Tee column end fitting 1/8in x 1/16in x 1/16in	.75mm in/out .75mm side 1mm cone SS
TCEF211T	Tee column end fitting 1/8in x 1/16in x 1/16in	1.65mm in/out .4mm side 1mm cone SS
TCEF411	Tee column end fitting 1/4in x 1/16in x 1/16in	.75mm in/out .75mm side 4.6mm cone SS
TCEF411C	Tee column end fitting 1/4in x 1/16in x 1/16in	.25mm in/out .25mm side 4.6mm cone SS
TCEF411T	Tee column end fitting 1/4in x 1/16in x 1/16in	1.65mm in/out .75mm side 4.6mm cone SS
TCEF421	Tee column end fitting 1/4in x 1/8in x 1/16in	.75mm in/out .75mm side 4.6mm cone SS
TCEF611	Tee column end fitting 3/8in x 1/16in x 1/16in	.75mm in/out .75mm side 6mm cone SS
TCEF611T	Tee column end fitting 3/8in x 1/16in x 1/16in	1.65mm in/out .75mm side 6mm cone SS

# VICI Valco-Cheminert PRODUCT LIST 2012




TCEF811	Tee column end fitting 1/2in x 1/16in x 1/16in	.75mm in/out .75mm side 9mm cone SS
TCEF811T	Tee column end fitting 1/2in x 1/16in x 1/16in	1.65mm in/out .75mm side 9mm cone SS
TEFNI.101	Tubing electroformed.360um OD .001in /25um ID	Nickel
TEFNI.102	Tubing electroformed.360um OD	.002in /50um ID Nickel
TEFNI.104	Tubing electroformed.360um OD	.004in /100um ID Nickel
TEFNI.105	Tubing electroformed.360um OD	.005in /125um ID Nickel
TEFNI.110	Tubing electroformed.360um OD	.007in /175um ID Nickel
TEFNI.502	Tubing electroformed.1/32in OD.002in /.05mm ID	Nickel
TEFNI.504	Tubing electroformed.1/32in OD.004in /.1mm ID	Nickel
TEFNI.505	Tubing electroformed.1/32in OD.005in /.12mm ID	Nickel
TEFNI.507	Tubing electroformed.1/32in OD.007in /.18mm ID	Nickel
TEFNI.510	Tubing electroformed.1/32in OD.010in /.25mm ID	Nickel
TEFNI.512	Tubing electroformed.1/32in OD	.0126in /.32mm ID Nickel
TEFNI.515	Tubing electroformed.1/32in OD.015in /.38mm ID	Nickel
TEFNI.520	Tubing electroformed.1/32in OD.020in /.5mm ID	Nickel
TEFNI120	Tubing electroformed.1/16in OD.020in /.5mm ID	Nickel
TEFNI130	Tubing electroformed.1/16in OD.030in /.75mm ID	Nickel
TEFNI140	Tubing electroformed. 1/16in OD.040in /1mm ID	Nickel
TEMPLATE1	Template for flow diagrams	
TFEP110-10M	Tubing 1/16in OD .010in /.25mm ID	FEP 10m/pkg
TFEP110-25M	Tubing 1/16in OD .010in /.25mm ID	FEP 25m/pkg
TFEP110-5M	Tubing 1/16in OD .010in /.25mm ID	FEP 5m/pkg
TFEP120-10M	Tubing 1/16in OD .020in /.5mm ID	FEP 10m/pkg
TFEP120-25M	Tubing 1/16in OD .020in /.5mm ID	FEP 25m/pkg
TFEP120-5M	Tubing 1/16in OD .020in /.5mm ID	FEP 5m/pkg
TFEP130-10M	Tubing 1/16in OD .030in /.75mm ID	FEP 10m/pkg
TFEP130-25M	Tubing 1/16in OD .030in /.75mm ID	FEP 25m/pkg
TFEP130-5M	Tubing 1/16in OD .030in /.75mm ID	FEP 5m/pkg
TFEP260-10M	Tubing 1/8in OD .060in /1.52mm ID	FEP 10m/pkg
TFEP260-25M	Tubing 1/8in OD .060in /1.52mm ID	FEP 25m/pkg
TFEP260-5M	Tubing 1/8in OD .060in /1.52mm ID	FEP 5m/pkg
THC130	Tubing 1/16in OD .030in /.75mm ID HC	
THC270	Tubing 1/8in OD .070in /1.77mm ID HC	
TINCO261	Tubing 1/8in OD .061in /1.55mm ID	Inconel 600
TINCO282	Tubing 1/8in OD .082in /2.1mm ID	Inconel 600
TNI.507	Tubing 1/32in OD .007in /.18 mm ID Nickel	
TNI.510	OBSOLETE NOT AVAILABLE ANYMORE	
TNI.520	Tubing 1/32in OD .020in /.5mm ID Nickel	
TNI105	Tubing 1/16in OD .005in /.12mm ID Nickel	
TNI110	Tubing 1/16in OD .010in /.25mm ID Nickel	
TNI120	Tubing 1/16in OD .020in /.5mm ID Nickel	
TNI130	Tubing 1/16in OD .030in /.75mm ID Nickel	
TNI140	Tubing 1/16in OD .040in /1mm ID Nickel	
TNI285	Tubing 1/8in OD .085in /2.15mm ID Nickel	
TPK.102	Tubing 360um OD .002in /50um ID	PEEK
TPK.104	Tubing 360um OD .004in /100um ID	PEEK
TPK.105	Tubing 360um OD .005in /125um ID	PEEK
TPK.106	Tubing 360um OD .006in /150um ID	PEEK
TPK.502-10M	Tubing 1/32in OD .0025in /63um ID	PEEK 10m/pkg
TPK.502-25M	Tubing 1/32in OD .0025in /63um ID	PEEK 25m/pkg
TPK.502-5M	Tubing 1/32in OD .0025in /63um ID	PEEK 5m/pkg
TPK.505-10M	Tubing 1/32in OD .005in /125um	PEEK 10m/pkg
TPK.505-25M	Tubing 1/32in OD .005in /125um	PEEK 25m/pkg
TPK.505-5M	Tubing 1/32in OD .005in /125um	PEEK 5m/pkg
TPK.510-10M	Tubing 1/32in OD .010in /.25mm	PEEK 10m/pkg
TPK.510-25M	Tubing 1/32in OD .010in /.25mm	PEEK 25m/pkg
TPK.510-5M	Tubing 1/32in OD .010in /.25mm	PEEK 5m/pkg
TPK.515-10M	Tubing 1/32in OD .015in /.38mm	PEEK 10m/pkg
TPK.515-25M	Tubing 1/32in OD .015in /.38mm	PEEK 25m/pkg
TPK.515-5M	Tubing 1/32in OD .015in /.38mm	PEEK 5m/pkg
TPK106-10M	Tubing 1/16in OD .006in /.15mm ID	PEEK 10m/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




TPK106-25M	Tubing 1/16in OD .006in /.15mm ID	PEEK 25m/pkg
TPK106-5M	Tubing 1/16in OD .006in /.15mm ID	PEEK 5m/pkg
TPK110-10M	Tubing 1/16in OD .010in /.25mm ID	PEEK 10m/pkg
TPK110-25M	Tubing 1/16in OD .010in /.25mm ID	PEEK 25m/pkg
TPK110-5M	Tubing 1/16in OD .010in /.25mm ID	PEEK 5m/pkg
TPK120-10M	Tubing 1/16in OD .020in /.5mm ID	PEEK 10m/pkg
TPK120-25M	Tubing 1/16in OD .020in /.5mm ID	PEEK 25m/pkg
TPK120-5M	Tubing 1/16in OD .020in /.5mm ID	PEEK 5m/pkg
TPK130-10M	Tubing 1/16in OD .030in /.75mm ID	PEEK 10m/pkg
TPK130-25M	Tubing 1/16in OD .030in /.75mm ID	PEEK 25m/pkg
TPK130-5M	Tubing 1/16in OD .030in /.75mm ID	PEEK 5m/pkg
TPK260-10M	Tubing 1/8in OD .060in /1.52mm ID	PEEK 10m/pkg
TPK260-25M	Tubing 1/8in OD .060in /1.52mm ID	PEEK 25m/pkg
TPK260-5M	Tubing 1/8in OD .060in /1.52mm ID	PEEK 5m/pkg
TSS.505	Tubing 1/32in OD .005in /.12mm 316SS	
TSS.507	Tubing 1/32in OD .007in /.18mm 316SS	
TSS.510	Tubing 1/32in OD .010in /.25mm 316SS	
TSS.520	Tubing 1/32in OD .020in /.5mm ID 316SS	
TSS105	Tubing 1/16in OD .005in /.12mm ID 316SS	
TSS110	Tubing 1/16in OD .010in /.25mm ID 316SS	
TSS115	Tubing 1/16in OD .015in /.38mm ID 316SS	
TSS120	Tubing 1/16in OD .020in /.5mm ID 316SS	
TSS126	Tubing 1/16in OD .026in /.66mm ID 316SS	
TSS130	Tubing 1/16in OD .030in /.75mm ID 316SS	
TSS140	Tubing 1/16in OD .040in /1mm ID 316SS	
TSS230	Tubing 1/8in OD .030in /.75mm ID 316SS	
TSS240	Tubing 1/8in OD .040in /1mm ID 316SS	
TSS260	Tubing 1/8in OD .060in /1.50mm ID 316SS	
TSS267	Tubing 1/8in OD .067in /1.7mm ID 316SS	
TSS285	Tubing 1/8in OD .085in /2.15mm ID 316SS	
TTF106-10M	Tubing 1/16in OD .006in /.15mm ID	PTFE 10m/pkg
TTF106-25M	Tubing 1/16in OD .006in /.15mm ID	PTFE 25m/pkg
TTF106-5M	Tubing 1/16in OD .006in /.15mm ID	PTFE 5m/pkg
TTF110-10M	Tubing 1/16in OD .010in /.25mm ID	PTFE 10m/pkg
TTF110-25M	Tubing 1/16in OD .010in /.25mm ID	PTFE 25m/pkg
TTF110-5M	Tubing 1/16in OD .010in /.25mm ID	PTFE 5m/pkg
TTF115-10M	Tubing 1/16in OD .015in /.38mm ID	PTFE 10m/pkg
TTF115-25M	Tubing 1/16in OD .015in /.38mm ID	PTFE 25m/pkg
TTF115-5M	Tubing 1/16in OD .015in /.38mm ID	PTFE 5m/pkg
TTF120-10M	Tubing 1/16in OD .020in /.5mm ID	PTFE 10m/pkg
TTF120-25M	Tubing 1/16in OD .020in /.5mm ID	PTFE 25m/pkg
TTF120-5M	Tubing 1/16in OD .020in /.5mm ID	PTFE 5m/pkg
TTF130-10M	Tubing 1/16in OD .030in /.75mm ID	PTFE 10m/pkg
TTF130-25M	Tubing 1/16in OD .030in /.75mm ID	PTFE 25m/pkg
TTF130-5M	Tubing 1/16in OD .030in /.75mm ID	PTFE 5m/pkg
TTF230-10M	Tubing 1/8in OD .030in /.75mm ID	PTFE 10m/pkg
TTF230-25M	Tubing 1/8in OD .030in /.75mm ID	PTFE 25m/pkg
TTF230-5M	Tubing 1/8in OD .030in /.75mm ID	PTFE 5m/pkg
TTF260-10M	Tubing 1/8in OD .060in /1.52mm ID	PTFE 10m/pkg
TTF260-25M	Tubing 1/8in OD .060in /1.52mm ID	PTFE 25m/pkg
TTF260-5M	Tubing 1/8in OD .060in /1.52mm ID	PTFE 5m/pkg
TTF285-10M	Tubing 1/8in OD .085in /2.15mm ID	PTFE 10m/pkg
TTF285-25M	Tubing 1/8in OD .085in /2.15mm ID	PTFE 25m/pkg
TTF285-5M	Tubing 1/8in OD .085in /2.15mm ID	PTFE 5m/pkg
TTZ110-10M	Tubing 1/16in OD .010in /.25mm ID	ETFE 10m/pkg
TTZ110-25M	Tubing 1/16in OD .010in /.25mm ID	ETFE 25m/pkg
TTZ110-5M	Tubing 1/16in OD .010in /.25mm ID	ETFE 5m/pkg
TTZ120-10M	Tubing 1/16in OD .020in /.5mm ID	ETFE 10m/pkg
TTZ120-25M	Tubing 1/16in OD .020in /.5mm ID	ETFE 25m/pkg
TTZ120-5M	Tubing 1/16in OD .020in /.5mm ID	ETFE 5m/pkg
TTZ130-10M	Tubing 1/16in OD .030in /.75mm ID	ETFE 10m/pkg

# VICI Valco-Cheminert PRODUCT LIST 2012




TTZ130-25M	Tubing 1/16in OD .030in / .75mm ID	ETFE 25m/pkg
TTZ130-5M	Tubing 1/16in OD .030in / .75mm ID	ETFE 5m/pkg
TTZ260-10M	Tubing 1/8in OD .060in / 1.52mm ID	ETFE 10m/pkg
TTZ260-25M	Tubing 1/8in OD .060in / 1.52mm ID	ETFE 25m/pkg
TTZ260-5M	Tubing 1/8in OD .060in / 1.52mm ID	ETFE 5m/pkg
V-EMPMCR-HP1100	Remote cable VICI micro E mpos. actuator	4/6 pos to HP1100 G1351A interface board
V-EMPMCR-HP1100-10	Remote cable VICI micro E mpos. actuator	8/10 pos to HP1100 G1351A interface board
V-EMPMCR-HP6890	BCD cable VICI micro E mpos. actuator	(26 pin connector) to HP 6890/7890
V-EMPMCR-HP6890N	BCD cable VICI micro E mpos. actuator	(26 pin connector) to HP 6890N/7890N
V-EMPMCR-WA2690	Remote cable VICI micro E mpos. actuator	4/6pos to Waters 2690 Separation Module
V-EMPMCR-WA2690-1C	Remote cable VICI micro E mpos. actuator	10pos to Waters 2690 Separation Module
VISF.5FPK	Syringe fillport for 1/32in connection fingertight	for .75in /19mm long needle 26Gauge OD PEEK
VISF-1	Syringe fillport for 1/16in connection	for .75in /19mm long needle 22Gauge OD
VISF-2	Syringe fillport for 1/16in connection	for 2.0in /51mm long needle 22Gauge OD
VISF-A	Syringe fillport for 1/8in fittings	for 2.0in /51mm long needle 22Gauge OD
VISL-1	Liner replacement for VISF-1 or C-VISF-1L	PTFE
VISL-1H	Liner replacement for C-VISF-1H	PTFE
VISL-2	Liner replacement for VISF-2 or VISF-A	PTFE
V-RA-24VDC-HP6890	Relay assy. for 2pos micro actuator	to HP6890/7890 24VDC
V-RA-24VDC-VA3800	Relay assy. for 2pos micro actuator	to VA3800 24VDC
V-RA-5VDC-HP1100	Relay assy. for 2pos micro actuator	to HP1100 external event 5VDC from micr act.
V-RA-5VDC-WA2690	Relay assy. for 2pos micro actuator	to WA2690 external event 5VDC from micro act.
VSH	Valve spanner handle for mpos valve	
V-SV-D32-110VAC	3-way manifold solenoid valve assy.	110Va CE
V-SV-D32-220VAC	3-way manifold solenoid valve assy.	220Vac CE
V-SV-D32-24VUC	3-way manifold solenoid valve assy.	24VacVdc CE
V-SV-S42-120VAC	4-way solenoid air valve	120Vac CE
V-SV-S42-220VAC	4-way solenoid air valve	220Vac CE
V-SV-S42-24VUC	4-way solenoid air valve	24VacVdc CE
WMMA	Close mount hdw. manual W/UW	2-pos 3468 port valves
WMMA10	Close mount hdw. manual	valves with unthreaded mounting holes
WMMASO	Manual W-knob for use with standoff	2-pos W-type valves
XLZN1-10	Nut 1/16in 25mm long SS 10/pkg	
XLZN2-10	Nut 1/8in 27mm long SS 10/pkg	
Z10M1	Manifold 1/16in 10 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z10M1PK	Manifold 1/16in 10 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z10M2	Manifold 1/8in 10 inlets 1 outlet	inlet 2mm/outlet 2mm bore SS
Z12M1	Manifold 1/16in 12 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z12M1PK	Manifold 1/16in 12 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z12M2	Manifold 1/8in 12 inlets 1 outlet	inlet 2mm/outlet 2mm bore SS
Z14M1	Manifold 1/16in 14 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z14M1PK	Manifold 1/16in 14 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z16M1PK	Manifold 1/16in 16 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z4M1	Manifold 1/16in 4 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z4M1PK	Manifold 1/16in 4 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z4M2	Manifold 1/8in 4 inlets 1 outlet	inlet 2mm/outlet 2mm bore SS
Z6M1	Manifold 1/16in 6 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z6M1PK	Manifold 1/16in 6 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z6M2	Manifold 1/8in 6 inlets 1 outlet	inlet 2mm/outlet 2mm bore SS
Z8M1	Manifold 1/16in 8 inlets 1 outlet	inlet .25mm/outlet .75mm bore SS
Z8M1PK	Manifold 1/16in 8 inlets 1 outlet	inlet .25mm/outlet .75mm bore PEEK
Z8M2	Manifold 1/8in 8 inlets 1 outlet	inlet 2mm/outlet 2mm bore SS
ZAOR11	Internal adaptor 1/16in to 10-32 thread	O-ring seal SS
ZAOR12	Internal adaptor 1/16in to 5/16-24 thread	O-ring seal SS
ZAOR22	Internal adaptor 1/8in to 5/16-24 thread	O-ring seal SS
ZBAA1	Aerosol adaptor bulkhead union 1/16in SS	Bore 0.75 mm
ZBNV1	Micrometering needle valve bulkhead 2-225 ml/min	1/16in int 2um screenlubricated 225CSS/Viton
ZBNV1-D	Micrometering needle valve bulkhead 2-225 ml/min	1/16in int 2um screen225Cnon-lubricatedSS/Viton
ZBNV1F	Micrometering needle valve bulkhead 2-175 ml/min	1/16in int 2um screenlubricated225CSS/Viton
ZBNV1F-D	Micrometering needle valve bulkhead 2-175 ml/min	1/16in int 2um screen225Cnon-lubricatedSS/Viton
ZBNV1F-KZ	Micrometering needle valve bulkhead 2-175 ml/min	1/16in int 2um screen315Cnon-lubricatSS/Kalrez

# VICI Valco-Cheminert PRODUCT LIST 2012




ZBNV1-KZ	Micrometering needle valve bulkhead 2-225 ml/min	1/16in int 2um screen315Cnon-lubricat.SS
ZBNV1LF	Micrometering needle valve bulkhead 2-90 ml/min	1/16in int2um screen225C lubricatedSS/Vi
ZBNV1LF-D	Micrometering needle valve bulkhead 2-90 ml/min	1/16in int2um screen225Cnon-lubricatedSS
ZBNV1LF-KZ	Micrometering needle valve bulkhead 2-90 ml/min	1/16in int2um screen315Cnon-lubricatedSS
ZBRU1.5	Bulkhead reducing union 1/16in x 1/32in	.25mm bore SS
ZBRU1.5FPK	Bulkhead reducing union 1/16in x 1/32in	.25mm bore PEEK
ZBRU1.5L	Bulkhead reducing union 1/16in x 1/32in	.5mm bore SS
ZBRU1.5LFPK	Bulkhead reducing union 1/16in x 1/32in	.5mm bore PEEK
ZBRU1.5T	Bulkhead reducing union 1/16in x 1/32in	1/32in through bore SS
ZBRU1.5TFPK	Bulkhead reducing union 1/16in x 1/32in	1/32in trough bore PEEK
ZBRU2.5	Bulkhead reducing union 1/8in x 1/32in	.25mm bore SS
ZBRU2.5FPK	Bulkhead reducing union 1/8in x 1/32in	.25mm bore PEEK
ZBRU2.5L	Bulkhead reducing union 1/8in x 1/32in	.5mm bore SS
ZBRU2.5LFPK	Bulkhead reducing union 1/8in x 1/32in	.5mm bore PEEK
ZBRU2.5T	Bulkhead reducing union 1/8in x 1/32in	1/32in through bore SS
ZBRU2.5TFPK	Bulkhead reducing union 1/8in x 1/32in	1/32in trough bore PEEK
ZBRU21	Bulkhead reducing union 1/8in x 1/16in	.75mm bore SS
ZBRU21C	Bulkhead reducing union 1/8in x 1/16in	.25mm bore SS
ZBRU21CFPK	Bulkhead reducing union 1/8in x 1/16in	.25mm bore PEEK
ZBRU21FPK	Bulkhead reducing union 1/8in x 1/16in	.75mm bore PEEK
ZBRU21LFPK	Bulkhead reducing union 1/8in x 1/16in	1mm bore PEEK
ZBRU21T	Bulkhead reducing union 1/8in x 1/16in	1/16in through bore SS
ZBRU21TFPK	Bulkhead reducing union 1/8in x 1/16in	1/16in trough bore PEEK
ZBRU41	Bulkhead reducing union 1/4in x 1/16in	.75mm bore SS
ZBRU41C	Bulkhead reducing union 1/4in x 1/16in	.25mm bore SS
ZBRU41T	Bulkhead reducing union 1/4in x 1/16in	1/16in through bore SS
ZBRU42	Bulkhead reducing union 1/4in x 1/8in	0.75mm bore SS
ZBRU42L	Bulkhead reducing union 1/4in x 1/8in	2mm bore 316SS
ZBRU42T	Bulkhead reducing union 1/4in x 1/8in	1/8in through bore 316SS
ZBRUF1.5	Bulkhead reducing union 1/16in x 1/32in	2 um frit pressed SS
ZBRUF21	Bulkhead reducing union 1/8in x 1/16in	.75mm bore pressed 2 um frit SS
ZBRUF42	Bulkhead reducing union 1/4in x 1/8in	2mm bore pressed 2 um frit SS
ZBRUFR1.5	Bulkhead reducing union 1/16in x 1/32in	.25mm bore removable 2 um screen SS
ZBRUFR1.5F	Bulkhead reducing union 1/16in x 1/32in	.25mm bore removable 2um frit 316SS
ZBRUFR21	Bulkhead reducing union 1/8in x 1/16in	.75mm bore removable 2 um screen SS
ZBRUFR21F	Bulkhead reducing union 1/8in x 1/16in	.75mm bore removable 2 um frit 316SS
ZBRUFR41	Bulkhead reducing union 1/4in x 1/16in	1mm bore removable 2 um screen SS
ZBRUFR41F	Bulkhead reducing union 1/4in x 1/16in	1mm bore removable 2 um frit 316SS
ZBRUFR42	Bulkhead reducing union 1/4in x 1/8in	2mm bore removable 2 um screen SS
ZBRUFR42F	Bulkhead reducing union 1/4in x 1/8in	2mm bore removable 2 um frit 316SS
ZBU.5	Bulkhead union 1/32in	.25mm bore SS
ZBU.5FPK	Bulkhead union 1/32in fingertight	.25mm bore PEEK
ZBU.5L	Bulkhead union 1/32in	.5mm bore SS
ZBU.5LFPK	Bulkhead union 1/32in fingertight	.5mm bore PEEK
ZBU.5T	Bulkhead union 1/32in	1/32in through bore SS
ZBU.5TFPK	Bulkhead union 1/32in fingertight	1/32in trough bore PEEK
ZBU.5XC	Bulkhead union 1/32in	.15mm bore SS
ZBU1	Bulkhead union 1/16in	.75mm bore SS
ZBU1C	Bulkhead union 1/16in	.25mm bore SS
ZBU1CFPK	Bulkhead union 1/16in fingertight	.25mm bore PEEK
ZBU1CPK	Bulkhead union 1/16in hex	.25mm bore PEEK
ZBU1FPK	Bulkhead union 1/16in fingertight	.75mm bore PEEK
ZBU1L	Bulkhead union 1/16in	1mm bore SS
ZBU1M	Bulkhead union 1/16in	.5mm bore SS
ZBU1MFPK	Bulkhead union 1/16in fingertight	.5mm bore PEEK
ZBU1MPK	Bulkhead union 1/16in hex	.5mm bore PEEK
ZBU1PK	Bulkhead union 1/16in hex	.75mm bore PEEK
ZBU1T	Bulkhead union 1/16in	1/16in through bore SS
ZBU1TFPK	Bulkhead union 1/16in fingertight	1/16in trough bore PEEK
ZBU1TPK	Bulkhead union 1/16in hex	1/16in trough bore PEEK
ZBU1XC	Bulkhead union 1/16in	.15mm bore SS

# VICI Valco-Cheminert PRODUCT LIST 2012




ZBU2	Bulkhead union 1/8in	.75mm bore SS
ZBU2L	Bulkhead union 1/8in	2mm bore SS
ZBU2LPK	Bulkhead union 1/8in hex	2mm bore PEEK
ZBU2PK	Bulkhead union 1/8in hex	.75mm bore PEEK
ZBU2T	Bulkhead union 1/8in	1/8in through bore SS
ZBU2TPK	Bulkhead union 1/8in hex	1/8in trough bore PEEK
ZBU4	Bulkhead union 1/4in	.75mm bore SS
ZBU4L	Bulkhead union 1/4in	4.6mm bore SS
ZBU4T	Bulkhead union 1/4in	1/4in through bore SS
ZBUF1	Bulkhead union 1/16in	.75mm bore pressed 2 um frit SS
ZBUF2	Bulkhead union 1/8in	.75mm bore pressed 2 um frit SS
ZBUF4	Bulkhead union 1/4in	4.6mm bore pressed 2 um frit SS
ZBUFLKF	Cheminert luer adapter bulkhead union	female luer to 10-32 1mm bore CTFE
ZBUFLPK	Cheminert luer adapter bulkhead union	female luer to 10-32 1 mm bore PEEK
ZBUFR.5	Bulkhead union 1/32in	.25mm bore removable 2 um screen SS
ZBUFR.5F	Bulkhead union 1/32in	.25mm bore removable 2 um frit 316SS
ZBUFR1	Bulkhead union 1/16in	.50mm bore removable 2 um screen SS
ZBUFR1C	Bulkhead union 1/16in	.25mm bore removable 2 um screen SS
ZBUFR1CF	Bulkhead union 1/16in	.25mm bore removable 2 um frit 316SS
ZBUFR1F	Bulkhead union 1/16in	.5mm bore removable 2 um frit 316SS
ZBUFR2	Bulkhead union 1/8in	2mm bore removable 2 um screen SS
ZBUFR2F	Bulkhead union 1/8in	2mm bore removable 2 um frit 316SS
ZBUMLKF	Cheminert luer adapter bulkhead union	luer male to 10-32 ZDV 1mm bore CTFE
ZBUMLPK	Cheminert luer adapter bulkhead union	luer male to 10-32 ZDV 1mm bore PEEK
ZC.5	Cap 1/32in	SS
ZC.5FPK	Cap 1/32in fingertight PEEK	
ZC1	Cap 1/16in	SS
ZC1FPK	Cap 1/16in fingertight PEEK	
ZC1PK	Cap 1/16in hex PEEK	
ZC2	Cap 1/8in	SS
ZC2PK	Cap 1/8in hex PEEK	
ZC4	Cap 1/4in	SS
ZCN1-10	Counter bored nut 1/16in	SS 10/pkg
ZERU1.5FPK	Union Reducing Int/Ext 1/16in 1/32in 20mm	200um PEEK
ZF.5-10	Ferrule 1/32in 303SS 10/pkg	
ZF.5B-10	Ferrule 1/32in Brass 10/pkg	
ZF.5FEP-10	Ferrule 1/32in FEP 10/pkg	
ZF.5GP-10	Ferrule 1/32in SS gold plated 10/pkg	
ZF.5GV-5	Ferrule 1/32in	Polyimide Graphite 5/pkg
ZF.5HC	Ferrule 1/32in HC	
ZF.5KF-10	Ferrule 1/32in CTFE 10/pkg	
ZF.5NI	Ferrule 1/32in Nickel	
ZF.5PFA-10	Ferrule 1/32in PFA 10/pkg	
ZF.5PK-10	Ferrule 1/32in PEEK 10/pkg	
ZF.5PK-5	Ferrule 1/32in PEEK 5/pkg	
ZF.5S6-10	Ferrule 1/32in 316SS 10/pkg	
ZF.5TF-10	Ferrule 1/32in	PTFE Virgin 10/pkg
ZF.5TFG-10	Ferrule 1/32in	PTFE Glass-filled 25% 10/pkg
ZF.5TI	Ferrule 1/32in Titanium	
ZF.5V1-5	Ferrule 1/32in	Polyimide Virgin 5/pkg
ZF.5V-5	Ferrule 1/32in	Polyimide Valcon 5/pkg
ZF1-10	Ferrule 1/16in 303SS 10/pkg	
ZF1B-10	Ferrule 1/16in Brass 10/pkg	
ZF1FEP-10	Ferrule 1/16in FEP 10/pkg	
ZF1GP-10	Ferrule 1/16in SS gold plated 10/pkg	
ZF1GV-5	Ferrule 1/16in	Polyimide Graphite 5/pkg
ZF1HC	Ferrule 1/16in HC	
ZF1KF-10	Ferrule 1/16in CTFE 10/pkg	
ZF1NI	Ferrule 1/16in Nickel	
ZF1PFA-10	Ferrule 1/16in PFA 10/pkg	
ZF1PK-10	Ferrule 1/16in PEEK 10/pkg	

# VICI Valco-Cheminert PRODUCT LIST 2012




ZF1S6-10	Ferrule 1/16in 316SS 10/pkg	
ZF1TF-10	Ferrule 1/16in	PTFE Virgin 10/pkg
ZF1TFG-10	Ferrule 1/16in	PTFE Glass-filled 25% 10/pkg
ZF1TI	Ferrule 1/16in Titanium	
ZF1V1-5	Ferrule 1/16in	Polyimide Virgin 5/pkg
ZF1V-5	Ferrule 1/16in	Polyimide Valcon 5/pkg
ZF1VISF	Ferrule replacement for VISF-1 or VISF-2	or C-VISF-1 (L or H) CTFE
ZF2-10	Ferrule 1/8in 303SS 10/pkg	
ZF2B-10	Ferrule 1/8in Brass 10/pkg	
ZF2FEP-10	Ferrule 1/8in FEP 10/pkg	
ZF2GP-10	Ferrule 1/8in SS gold plated 10/pkg	
ZF2GV-5	Ferrule 1/8in	Polyimide Graphite 5/pkg
ZF2HC	Ferrule 1/8in HC	
ZF2KF-10	Ferrule 1/8in CTFE 10/pkg	
ZF2NI	Ferrule 1/8in Nickel	
ZF2PFA-10	Ferrule 1/8in PFA 10/pkg	
ZF2PK-10	Ferrule 1/8in PEEK 10/pkg	
ZF2S6-10	Ferrule 1/8in 316SS 10/pkg	
ZF2SI	Replacement support for septum injector	Valcon T polyimide
ZF2TF-10	Ferrule 1/8in	PTFE Virgin 10/pkg
ZF2TFG-10	Ferrule 1/8in	PTFE Glass-filled 25% 10/pkg
ZF2TI	Ferrule 1/8in Titanium	
ZF2V1-5	Ferrule 1/8in	Polyimide Virgin 5/pkg
ZF2V-5	Ferrule 1/8in	Polyimide Valcon 5/pkg
ZF4B-10	Ferrule 1/4in Brass 10/pkg	
ZF4FEP-10	Ferrule 1/4in FEP 10/pkg	
ZF4GP-10	Ferrule 1/4in SS gold plated 10/pkg	
ZF4GV-5	Ferrule 1/4in	Polyimide Graphite 5/pkg
ZF4HC	Ferrule 1/4in HC	
ZF4KF-10	Ferrule 1/4in CTFE 10/pkg	
ZF4NI	Ferrule 1/4in Nickel	
ZF4PFA-10	Ferrule 1/4in PFA 10/pkg	
ZF4PK-10	Ferrule 1/4in PEEK 10/pkg	
ZF4S6-10	Ferrule 1/4in 316SS 10/pkg	
ZF4TF-10	Ferrule 1/4in	PTFE Virgin 10/pkg
ZF4TFG-10	Ferrule 1/4in	PTFE Glass-filled 25% 10/pkg
ZF4TI	Ferrule 1/4in Titanium	
ZF4V1-5	Ferrule 1/4in	Polyimide Virgin 5/pkg
ZF4V-5	Ferrule 1/4in	Polyimide Valcon 5/pkg
ZF6FEP-10	Ferrule 3/8in FEP 10/pkg	
ZF6GV-5	Ferrule 3/8in	Polyimide Graphite 5/pkg
ZF6KF-10	Ferrule 3/8in CTFE 10/pkg	
ZF6PFA-10	Ferrule 3/8in PFA 10/pkg	
ZF6PK-10	Ferrule 3/8in PEEK 10/pkg	
ZF6TF-10	Ferrule 3/8in	PTFE Virgin 10/pkg
ZF6TFG-10	Ferrule 3/8in	PTFE Glass-filled 25% 10/pkg
ZF6V1-5	Ferrule 3/8in	Polyimide Virgin 5/pkg
ZF6V-5	Ferrule 3/8in	Polyimide Valcon 5/pkg
ZF8FEP-10	Ferrule 1/2in FEP 10/pkg	
ZF8GV-5	Ferrule 1/2in	Polyimide Graphite 5/pkg
ZF8KF-10	Ferrule 1/2in CTFE 10/pkg	
ZF8PFA-10	Ferrule 1/2in PFA 10/pkg	
ZF8PK-10	Ferrule 1/2in PEEK 10/pkg	
ZF8TF-10	Ferrule 1/2in	PTFE Virgin 10/pkg
ZF8TFG-10	Ferrule 1/2in	PTFE Glass-filled 25% 10/pkg
ZF8V1-5	Ferrule 1/2in	Polyimide Virgin 5/pkg
ZF8V-5	Ferrule 1/2in	Polyimide Valcon 5/pkg
ZGF.5PK	Ferrule double grip 1/32in	PEEK
ZGF.5PK-10	Ferrule double grip 1/32in PEEK 10/pkg	
ZGF1PK-10	Ferrule double grip 1/16in PEEK 10/pkg	
ZLA-.5	Female luer adapter	to 1/32in fittings

# VICI Valco-Cheminert PRODUCT LIST 2012




ZLA-1	Female luer adapter	to 1/16in fittings
ZLA-2	Female luer adapter	to 1/8in fittings
ZLCA1PK	Cheminert adaptor 1/4-28 to 10-32 ZDV	0.5mm bore PEEK
ZLTA41	Tube adaptor 1/4in x 1/16in 2.075in (52.70 mm) long	1/16in bore SS
ZN.5	Nut 1/32in standard	SS
ZN.5-10	Nut 1/32in standard SS 10/pkg	
ZN.5FPK	Nut 1/32 fingertight 10.7mm long	PEEK
ZN.5FPK-10	Nut 1/32in fingertight 10.7 mm long	PEEK 10/pkg
ZN1	Nut 1/16in standard	SS SS
ZN1-10	Nut 1/16in standard SS 10/pkg	
ZN1FPK	Nut 1/16in fingertight 22.4mm long	PEEK
ZN1FPK-10	Nut 1/16in fingertight 22.4mm long	PEEK 10/pkg
ZN1PK-10	Nut 1/16in hex 11.4mm long	PEEK 10/pkg
ZN1Q	Nut 1/16in length .43in	SS
ZN2-10	Nut 1/8in standard SS 10/pkg	
ZN2PK-10	Nut 1/8in hex 15.7mm long	PEEK 10/pkg
ZN4-10	Nut 1/4in standard SS 10/pkg	
ZNF1FKF-5	Nut 1/16in No Twist fingertight PEEK	with CTFE ferrule 5/pkg
ZNF1FPK-5	Nut 1/16in No Twist fingertight PEEK	with PEEK ferrule 5/pkg
ZNF1FPKG-5	Nut 1/16in No Twist fingertight PEEK	with glass-filled PEEK ferrule 5/pkg
ZNF1KF-5	Nut 1/16in No Twist hex short PEEK	with CTFE ferrule 5/pkg
ZNF1PK-5	Nut 1/16in No Twist hex short PEEK	with PEEK ferrule 5/pkg
ZNF1PKG-5	Nut 1/16in No Twist hex short PEEK	with glass-filled PEEK ferrule 5/pkg
ZNFT	Tightening Tool for PEEK hex-head	1/16in nuts blue
ZP.5	Plug 1/32in	SS
ZP.5FPK	Plug 1/32in fingertight PEEK	
ZP.5H	Plug 1/32in high pressure	SS
ZP1	Plug 1/16in	SS
ZP1FPK	Plug 1/16in fingertight PEEK	
ZP1H	Plug 1/16in high pressure	SS
ZP2	Plug 1/8in	SS
ZP2H	Plug 1/8in high pressure	SS
ZP2PK	Plug 1/8in hex PEEK	
ZP4	Plug 1/4in	SS
ZRF1.5KF-5	Internal reducing ferrule 1/16in x 1/32in	CTFE 5/pkg
ZRF1.5PK-5	Internal reducing ferrule 1/16in x 1/32in	PEEK 5/pkg
ZRF1.5TFG-5	Internal reducing ferrule 1/16in x 1/32in	PTFE Glass-filled 25% 5/pkg
ZRF1.5V1-5	Internal reducing ferrule 1/16in x 1/32in	Polyimide Virgin 5/pkg
ZRF1.5V-5	Internal reducing ferrule 1/16in x 1/32in	Polyimide Valcon 5/pkg
ZRF2.5KF-5	Internal reducing ferrule 1/8in x 1/32in	CTFE 5/pkg
ZRF2.5PK-5	Internal reducing ferrule 1/8in x 1/32in	PEEK 5/pkg
ZRF2.5TFG-5	Internal reducing ferrule 1/8in x 1/32in	PTFE Glass-filled 25% 5/pkg
ZRF2.5V1-5	Internal reducing ferrule 1/8in x 1/32in	Polyimide Virgin 5/pkg
ZRF2.5V-5	Internal reducing ferrule 1/8in x 1/32in	Polyimide Valcon 5/pkg
ZRF21KF-5	Internal reducing ferrule 1/8in x 1/16in	CTFE 5/pkg
ZRF21PK-5	Internal reducing ferrule 1/8in x 1/16in	PEEK 5/pkg
ZRF21TFG-5	Internal reducing ferrule 1/8in x 1/16in	PTFE Glass-filled 25% 5/pkg
ZRF21V1-5	Internal reducing ferrule 1/8in x 1/16in	Polyimide Virgin 5/pkg
ZRF21V-5	Internal reducing ferrule 1/8in x 1/16in	Polyimide Valcon 5/pkg
ZRF41KF-5	Internal reducing ferrule 1/4in x 1/16in	CTFE 5/pkg
ZRF41PK-5	Internal reducing ferrule 1/4in x 1/16in	PEEK 5/pkg
ZRF41TFG-5	Internal reducing ferrule 1/4in x 1/16in	PTFE Glass-filled 25% 5/pkg
ZRF41V1-5	Internal reducing ferrule 1/4in x 1/16in	Polyimide Virgin 5/pkg
ZRF41V-5	Internal reducing ferrule 1/4in x 1/16in	Polyimide Valcon 5/pkg
ZRF42KF-5	Internal reducing ferrule 1/4in x 1/8in	CTFE 5/pkg
ZRF42PK-5	Internal reducing ferrule 1/4in x 1/8in	PEEK 5/pkg
ZRF42TFG-5	Internal reducing ferrule 1/4in x 1/8in	PTFE Glass-filled 25% 5/pkg
ZRF42V1-5	Internal reducing ferrule 1/4in x 1/8in	Polyimide Virgin 5/pkg
ZRF42V-5	Internal reducing ferrule 1/4in x 1/8in	Polyimide Valcon 5/pkg
ZRU1.5	Reducing union 1/16in x 1/32in	.25mm bore SS
ZRU1.5FPK	Reducing union 1/16in x 1/32in	.25mm bore PEEK

# VICI Valco-Cheminert PRODUCT LIST 2012



ZRU1.5L	Reducing union 1/16in x 1/32in	.5mm bore SS
ZRU1.5LFPK	Reducing union 1/16in x 1/32in	.5mm bore PEEK
ZRU1.5T	Reducing union 1/16in x 1/32in	1/32in through bore SS
ZRU1.5TFPK	Reducing union 1/16in x 1/32in	1/32in trough bore PEEK
ZRU1.5XC	Reducing union 1/16in x 1/32in	.15mm bore SS
ZRU2.5	Reducing union 1/8in x 1/32in	.25mm bore SS
ZRU2.5FPK	Reducing union 1/8in x 1/32in	.25mm bore PEEK
ZRU2.5L	Reducing union 1/8in x 1/32in	.5mm bore SS
ZRU2.5LFPK	Reducing union 1/8in x 1/32in	.5mm bore PEEK
ZRU2.5T	Reducing union 1/8in x 1/32in	1/32in through bore SS
ZRU2.5TFPK	Reducing union 1/8in x 1/32in	1/32in trough bore PEEK
ZRU21	Reducing union 1/8in x 1/16in	.75mm bore SS
ZRU21C	Reducing union 1/8in x 1/16in	.25mm bore SS
ZRU21CFPK	Reducing union 1/8in x 1/16in	.25mm bore PEEK
ZRU21FPK	Reducing union 1/8in x 1/16in	.75mm bore PEEK
ZRU21LFPK	Reducing union 1/8in x 1/16in	1mm bore PEEK
ZRU21T	Reducing union 1/8in x 1/16in	1/16in through bore SS
ZRU21TFPK	Reducing union 1/8in x 1/16in	1/16in trough bore PEEK
ZRU41	Reducing union 1/4in x 1/16in	.75mm bore SS
ZRU41C	Reducing union 1/4in x 1/16in	.25mm bore SS
ZRU41T	Reducing union 1/4in x 1/16in	1/16in through bore SS
ZRU42	Reducing union 1/4in x 1/8in	0.75 mm bore SS
ZRU42L	Reducing union 1/4in x 1/8in	2mm bore SS
ZRU42T	Reducing union 1/4in x 1/8in	1/8in through bore SS
ZRUF1.5	Reducing union 1/16in x 1/32in	.25mm bore pressed 2 um frit SS
ZRUF21	Reducing union 1/8in x 1/16in	.75mm bore pressed 2 um frit SS
ZRUF42	Reducing union 1/4in x 1/8in	2mm bore pressed 2 um frit SS
ZRUFR1.5	Reducing filter union 1/16in x 1/32in	.25mm bore removable 2 um screen SS
ZRUFR1.5F	Reducing union 1/16in x 1/32in	.25mm bore removable 2 um frit 316SS
ZRUFR21	Reducing filter union 1/8in x 1/16in	.75mm bore removable 2 um screen SS
ZRUFR21F	Reducing filter union 1/8in x 1/16in	.75mm bore removable 2 um frit 316SS
ZRUFR41	Reducing filter union 1/4in x 1/16in	1mm bore removable 2 um screen SS
ZRUFR41F	Reducing filter union 1/4in x 1/16in	1mm bore removable 2 um frit 316SS
ZRUFR42	Reducing filter union 1/4in x 1/8in	2mm bore removable 2 um screen SS
ZRUFR42F	Reducing filter union 1/4in x 1/8in	2mm bore removable 2 um frit 316SS
ZSN1	Nut 1/16in for purge housing	SS
ZT.5	Tee 1/32in .25mm bore SS	
ZT.5FPK	Tee 1/32in fingertight	.25mm bore PEEK
ZT.5L	Tee 1/32in .50mm bore SS	
ZT.5LFPK	Tee 1/32in fingertight	.5mm bore PEEK
ZT1	Tee 1/16in .75mm bore SS	
ZT1C	Tee 1/16in .25mm bore SS	
ZT1CFPK	Tee 1/16in fingertight	.25mm bore PEEK
ZT1FPK	Tee 1/16in fingertight	.75mm bore PEEK
ZT1L	Tee 1/16in 1.0mm bore SS	
ZT1LFPK	Tee 1/16in fingertight	1mm bore PEEK
ZT1M	Tee 1/16in .50mm bore SS	
ZT1MFPK	Tee 1/16in fingertight	.5mm bore PEEK
ZT2	Tee 1/8in .75mm bore SS	
ZT2L	Tee 1/8in 2.0mm bore SS	
ZT2LPK	Tee 1/8in hex	2mm bore PEEK
ZT2PK	Tee 1/8in hex	.75mm bore PEEK
ZT4	Tee 1/4in 1.0mm bore SS	
ZT4L	Tee 1/4in 4.6mm bore SS	
ZTA41	Tube adaptor 1/4in x 1/16in 0.975in (24.76 mm) long	1/16in bore SS
ZU.5	Union 1/32in	.25mm bore SS
ZU.5FPK	Union 1/32in fingertight	.25mm bore PEEK
ZU.5L	Union 1/32in	.5mm bore SS
ZU.5LFPK	Union 1/32in fingertight	.5mm bore PEEK
ZU.5T	Union 1/32in	1/32in through bore SS
ZU.5TFPK	Union 1/32in fingertight	1/32in trough bore PEEK

# VICI Valco-Cheminert PRODUCT LIST 2012



ZU.5XC	Union 1/32in	.15mm bore SS
ZU1	Union 1/16in	.75mm bore SS
ZU1C	Union 1/16in	.25mm bore SS
ZU1CFPK	Union 1/16in fingertight	.25mm bore PEEK
ZU1CPK	Union 1/16in hex	.25mm bore PEEK
ZU1FPK	Union 1/16in fingertight	.75mm bore PEEK
ZU1FPK.5	Filter 1/16in fingertight 0.5 micron filter Ti	.5mm bore PEEK
ZU1L	Union 1/16in	1mm bore SS
ZU1M	Union 1/16in	.5mm bore SS
ZU1MFPK	Union 1/16in fingertight	.5mm bore PEEK
ZU1MPK	Union 1/16in hex	.5mm bore PEEK
ZU1PK	Union 1/16in hex	.75mm bore PEEK
ZU1T	Union 1/16in	1/16in through bore SS
ZU1TFPK	Union 1/16in fingertight	1/16in through bore PEEK
ZU1TPK	Union 1/16in hex	1/16in trough bore PEEK
ZU1XC	Union 1/16in	.15mm bore SS
ZU2	Union 1/8in	.75mm bore 316SS
ZU2L	Union 1/8in	2mm bore 316SS
ZU2LPK	Union 1/8in hex	2mm bore PEEK
ZU2PK	Union 1/8in hex	.75mm bore PEEK
ZU2T	Union 1/8in	1/8in through bore SS
ZU2TPK	Union 1/8in hex	1/8in trough bore PEEK
ZU4	Union 1/4in	.75mm bore SS
ZU4L	Union 1/4in	4.6mm bore SS
ZU4T	Union 1/4in	1/4in through bore SS
ZUF1	Filter union 1/16in	.75mm bore pressed 2 um frit SS
ZUF2	Filter union 1/8in	.75mm bore pressed 2 um frit SS
ZUF4	Filter union 1/4in	4.6mm bore pressed 2 um frit SS
ZUFLKF	Cheminert female luer adapter to 10-32	0.75mm bore CTFE
ZUFLPK	Cheminert female luer adapter to 10-32	0.75mm bore PEEK
ZUFR.5	Filter union 1/32in	.25mm bore removable 2 um screen SS
ZUFR.5F	Filter union 1/32in	.25mm bore removable 2 um frit 316SS
ZUFR1	Filter union 1/16in	.50mm bore removable 2 um screen SS
ZUFR1C	Filter union 1/16in	.25mm bore removable 2 um screen SS
ZUFR1CF	Filter union 1/16in	.25mm bore removable 2 um frit 316SS
ZUFR1F	Filter union 1/16in	.5mm bore removable 2um frit 316SS
ZUFR2	Filter 1/8in	2mm bore removable 2 um screen SS
ZUFR2F	Filter union 1/8in	2mm bore removable 2 um frit 316SS
ZVISF-1PE	Syringe fillport for 1/16in high throughput HPE	for minimum 35mm long needle 22Gauge OD
ZVISF-1PEH	Syringe fillport for 1/16in high throughput HPE	for minimum 35mm long needle 22Gauge OD
ZVISF-1PFA	Fillport ZDV for 1/16in PEEK with PFA bushing	for minimum 35mm long needle 22 Gauge
ZVISF-1PFAH	Fillport ZDV for 1/16in PEEK with PFA bushing HP	for minimum 35mm long needle 22 Gauge
ZX.5	Cross 1/32in .25mm bore SS	
ZX.5FPK	Cross 1/32in fingertight	.25mm bore PEEK
ZX.5L	Cross 1/32in .50mm bore SS	
ZX.5LFPK	Cross 1/32in fingertight	.5mm bore PEEK
ZX1	Cross 1/16in .75mm bore SS	
ZX1C	Cross 1/16in .25mm bore SS	
ZX1CFPK	Cross 1/16in fingertight	.25mm bore PEEK
ZX1FPK	Cross 1/16in fingertight	.75mm bore PEEK
ZX1L	Cross 1/16in 1.0mm bore SS	
ZX1LFPK	Cross 1/16in fingertight	1mm bore PEEK
ZX1M	Cross 1/16in .50mm bore SS	
ZX1MFPK	Cross 1/16in fingertight	.5mm bore PEEK
ZX2	Cross 1/8in .75mm bore SS	
ZX2L	Cross 1/8in 2.0mm bore SS	
ZX2LPK	Cross 1/8in hex	2mm bore PEEK
ZX2PK	Cross 1/8in hex	.75mm bore PEEK
ZX4	Cross 1/4in 1.0mm bore SS	
ZX4L	Cross 1/4in 4.6mm bore SS	
ZXLT4A1	Tube adaptor 1/4in x 1/16in 2.800in (71.12 mm) long	1/16in bore SS