



25259



25270



25271

Preventive Maintenance (PM) Kits for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Autosampler PM Kit Includes: rotor seal, needle seat, needle assembly, seat cap	1050	01078-68721	kit	25259
Autosampler PM Kit Includes: rotor seal, piston seals (2), needle assembly, needle seat, finger caps (3)	1100, 1200	G1313-68709	kit	25271
Pump PM Kit Includes: PTFE frits (2), outlet cap, gold disk seal, active inlet cartridge, piston seals (4), glass solvent filters (2)	1050, 1100, 1200	G1311-68710	kit	25270
Pump Maintenance Kit Includes: sapphire plunger (2), piston seals (2), outlet ball check valve, active inlet cartridge, PTFE frits (5)	1100, 1200	5065-4499	kit	25915



25915



25258



25916, 25278



25272



25275




25349



26424

For even more options, visit
www.restek.com/LCacc



25265

25058

Injector & Autosampler Supplies for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Needle Seat	1050, 1090	79846-67101	ea.	25258
Needle Seat Assembly	1100, 1200, 1220	G1329-87017	ea.	25265
Needle Assembly	1100, 1200	G1313-87201	ea.	25278
Needle Assembly, 900 µL	1050, 1100, 1200	G1313-87202	ea.	25916
Rotor Seal, 2-Groove (not for use with 7125 injection valve)	1050	0101-0626	ea.	25272
Rotor Seal, 2-Groove	1100, 1200	0100-1853	ea.	25275
Rotor Seal (Rheodyne-Style), 3-Groove	1090	1535-4048	ea.	25349
Stator Face Assembly	1100, 1200	0100-1851	ea.	26424
Connecting Tube	1050, 1100	G1311-67304	ea.	25058
Injector Lubricant, PTFE, 10 g	1050, 1100	79841-65501	ea.	26520
Solvent Inlet Filter, SS	1050, 1090, 1100, 1200	01018-60025	ea.	26423
Isolation Seal for Injection Valve	1100, 1200	0100-1852	ea.	25917



26520



26423



25917

NEW!