

Instrument Replacement Parts for PerkinElmer HPLC Systems

Injector & Autosampler Supplies for PerkinElmer HPLC Systems

Description	Model #	Similar to PE Part #	qty.	cat.#
Needle	PE Series 200 Pumps	N2930023	ea.	25432



25432



25434



26425



25438



25435



25437

Pump Supplies for PerkinElmer HPLC Systems

Description	Model #	Similar to PE Part #	qty.	cat.#
Low-Pressure Seal, Gold*	Int. 4000, 1, 10, 2, 200 Series, 250, 3, 3B, 4, 400, 410, 620	09907339	ea.	25435
Standard High-Pressure Seal, Gold*	Int. 4000, 1, 10, 2, 200 Series, 250, 3, 3B, 4, 400, 410, 620	09907324	ea.	25434
Standard Pump Seal Kit	200 Series	N2910383	kit	26425
High-Pressure Seal Kit Includes: pump seals (4), backup rings (4), O-rings (4)	Int. 4000, 1, 10, 2, 200 Series, 250, 3, 3B, 4, 400, 410, 620	0254-0275	kit	25433
Inlet Check Valve	Int. 4000, 1, 10, 2, 200 Series, 250, 3, 3B, 4, 400, 410, 620, Series 200 Micropump	02540177	ea.	25438
Outlet Check Valve	PE 200, 4, 250, 400, 410, 620	02540197	ea.	25437

*Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion; best for buffers.



25433



25436



25431

Deuterium Lamps for PerkinElmer HPLC Systems

Description	Model #	Similar to PE Part #	qty.	cat.#
Deuterium Lamp	Lambda-2, 5, 7, 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 25, 40, 800, 900	B016-0917	ea.	25436
Deuterium Lamp	200 Series, 785A	N2920149	ea.	25431