

Instrument Replacement Parts for Waters HPLC & ACQUITY UPLC® Systems



25143



25145



26519



25144



25146



26430

Preventive Maintenance (PM) Kits for Waters HPLC Systems

Description	Model #	Similar to Waters Part #	qty.	cat.#
2690/2695 Pump & Autosampler PM Kit Includes: sapphire plungers (2), seal wash plunger seals (4), head plunger seals (4), wash tube seals (4), battery, sparge diffuser, filter insert, face seals (4), solvent reservoir 10 µm filters (4), 250 µL syringe, check valve cartridges (2), wash tube seal, seal wash tube, PTFE washer, filter retainer, lower wash seal frit, needle wash frit, TFE washer, needle assembly, gold injector seals (2), stainless steel ferrule, compression screw	Alliance 2690, 2695	WAT270944	kit	25143
717 Autosampler PM Kit Includes: seal pack assembly, tube assembly (0.020" ID), needle, needle compression screw, 0.062 stainless steel ferrule, pre-column filter assembly, filter insert, 250 µL WISP syringe	717 Autosampler	WAT052669	kit	25145
600 Pump PM Kit Includes: PerformancePLUS cartridges (4), sparge diffusers (4), Super Seals (2), solvent reservoir 10 µm filters (4), sapphire plungers (2), reference valve button, valve disk spacer, valve disk, TFE ball plug, TFE seat, ruby ball, inlet tube body assembly manifold insert, insert seal, Belleville washers (2), flat washer	600 Pump	WAT052675	kit	25144
616 Pump PM Kit Includes: sapphire plungers (2), sparge diffusers (4), solvent reservoir 10 µm filters (4), plunger seals (2), plunger wash seals (2), O-rings (2), backup rings (2), check valve cartridges (4)	616 Pump	WAT052672	kit	25146
515 Pump PM Kit Includes: PerformancePLUS check valves (4), sparge diffuser, solvent reservoir 10 µm filter, sapphire plungers (2), plunger seals (2), pivot inserts (2), pivot guides (2), washer (2), plunger springs (2), retaining rings (2)	515 Pump	WAT052587	kit	26519
1525 Pump PM Kit Includes: sapphire plungers (4), check valve cartridges (8), plunger seals (4), solvent reservoir 10 µm filters (2), reference valve button, valve disk spacer, valve disk	1525 Pump	201000114	kit	26430

Simplify Routine LC Maintenance
Preventive Maintenance Kits

- Significant savings over instrument manufacturer prices.
- High-quality components in every kit.
- Wide range of options for HPLC systems and pumps.

www.restek.com/LC-Maintenance


26430