

## Instrument Replacement Parts for Waters HPLC &amp; ACQUITY UPLC® Systems

**Pump Supplies** for Waters HPLC & ACQUITY UPLC® Systems (Also see pages 332–333.)

Description	Model #	Similar to Waters Part #	qty.	cat.#
Head Plunger Seal Kit Includes: head plunger seals (2), back-up washers (2)	ACQUITY, nanoACQUITY	700002599	kit	26428
Wash Seal	ACQUITY, nanoACQUITY	700002598	2-pk.	26429
Pump O-Ring	ACQUITY, nanoACQUITY	WAT076152	ea.	25951
In-Line Filter Assembly, Stainless Steel Frit	ACQUITY, nanoACQUITY	289003547	ea.	25952
Filter Frit Cartridge, Stainless Steel	ACQUITY	700002913	ea.	25953
Solvent Bottle Filters	ACQUITY	700003616	7-pk.	25954
Tube Assembly, SSV to In-line Filter	nanoACQUITY	430001470	ea.	25948
Tube Assembly, Transducer to Check Valve, BSM	ACQUITY, nanoACQUITY	430001121	ea.	25949
Tube Assembly, Transducer to Check Valve, QSM	ACQUITY	430002357	ea.	25950
Tube Assembly, SSV to I2V	ACQUITY	430001443	ea.	25941
Super Seal for Analytical Heads	M6KA, 1515, 1525, 510, 515, 590, 600, 610, LC Module 1	WAT022946	ea.	25374
Plunger Seal, Gold**	1515, 1525, 510, 515, 600, 610, LC Module 1	—	ea.	25375
Plunger Seal, Black*	1515, 1525, 510, 515, 600, 610, LC Module 1	WAT026613	ea.	25378
Plunger Seal, Gold for EF Heads**	510, 515, 590, 600, 610, LC Module 1	700002282	ea.	25380
Pump Seal Kit Includes: plungers seals (2), wash seals (2), O-rings (2), back-up rings (2)	616, 1525 Micro	WAT034515	kit	25497
Seal Wash Plunger Seals	2690, 2695, 2790, 2795	WAT270160, WAT271018	2-pk.	25386
Head Plunger Seals, Gold	2690, 2695, 2790, 2795	WAT270789, WAT270938	2-pk.	25387
Head Plunger Seals, Black	2690, 2695, 2790, 2795	WAT271066	2-pk.	25388
Seal Wash Face Seal	2690, 2695, 2790, 2795	WAT271017	ea.	25428
Wash Tube Seal Kit	2690, 2695, 2790, 2795	WAT270940	4-pk.	25429
Face Seals Replacement Kit	2690, 2695, 2790, 2795	WAT270163, WAT270939	4-pk.	26547

\*Graphite-filled PTFE; best for organic solvents.

\*\*Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion; best for buffers.

