



Pump Supplies for Waters HPLC & ACQUITY UPLC® Systems (cont.)

Description	Model #	Similar to Waters Part #	qty.	cat.#
Primary Check Valve	ACQUITY, nanoACQUITY	700002596	2-pk.	25955
Accumulator Check Valve (Double Ball & Seat)	ACQUITY, nanoACQUITY	700002968	2-pk.	25956
I2 Check Valve Cartridge	ACQUITY	700005165	ea.	25957
Inlet Check Valve Assembly	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT033679, WAT025214	ea.	25360
Inlet Check Valve Housing	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT025203	ea.	25361
Inlet Check Valve Rebuild Kit Includes: blue washer (2), retainer gasket (2), Kel-F insert (2), ruby ball (2), TFE seat (2), TFE washer (2)	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT060495	kit	25362
Outlet Check Valve Assembly (Actuator Style)	M6KA, 510, 590	WAT025028	ea.	25363
Outlet Check Valve Housing (Actuator Style)	M6KA, 510, 590	WAT025212	ea.	25364
Outlet Check Valve Rebuild Kit (Actuator Style) Includes: 2 µm cup filter (2), compression spring (2), actuator (2), TFE washer actuator seat (2), TFE washer (2)	M6KA, 510, 590	WAT026016	kit	25365
Outlet Check Valve Assembly (Ball & Seat Style)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT025216	ea.	25366
Outlet Check Valve Housing (Ball & Seat Style)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT025207	ea.	25367
Outlet Check Valve Rebuild Kit (Ball & Seat Style) Includes: 2 µm cup filter (2), ruby ball (2), ball seat (2), Kel-F insert (2), blue washer (2), washer (2)	1525EF, 510, 515, 600, 610, LC Module 1	WAT026014	kit	25368
Inlet Check Valve Assembly, 225 µL (Extended Flow)	1525EF, 510, 515, 600, 610, LC Module 1	WAT060307	ea.	25369
Check Valve Rebuild Kit (Extended Flow) Includes: 3/16" EF ruby ball (2), EF ball seat (2), EF guide seal (2), EF ball guide (2), EF housing seal (2)	510, 590, 600, LC Module 1	WAT088223	kit	25371
PerformancePLUS Check Valve Housing	1515, 1525, 510, 515, 525, 600, 610	—	ea.	25372
PerformancePLUS Check Valve Cartridge	1515, 1525, 2690, 2695, 2795, 510, 515, 600, 610	700000254	2-pk.	25370
Check Valve Cartridges	2690, 2695, 2790, 2795	WAT270941	2-pk.	25373
Check Valve Cartridge, Stainless Steel	616	WAT024960	ea.	25498